

December 21, 2021

Terumo Clinical Supply Co., Ltd. c/o Vaibhav Sivaramakrishan, Terumo Medical Corporation Regulatory Affairs Specialist II 265 Davidson Ave., Suite 320 Somerset, New Jersey 08873

Re: K211078

Trade/Device Name: Progreat Lambda Regulation Number: 21 CFR 870.1200

Regulation Name: Diagnostic Intravascular Catheter

Regulatory Class: Class II Product Code: DQO, KRA Dated: November 30, 2021 Received: December 1, 2021

Dear Vaibhay Siyaramakrishan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Gregory O'Connell
Assistant Director
DHT2C: Division of Coronary
and Peripheral Intervention Devices
OHT2: Office of Cardiovascular Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

510(k) Number (if known)	
K211078	
Device Name Progreat Lambda	
Indications for Use (Describe) Progreat Lambda is intended for the infusion of contrast media, or vasculature, excluding the blood vessels belonging to the central cidrug infusion in intra-arterial therapy in the peripheral vasculature. coronary vessels.	rculatory system. Progreat Lambda is also indicated for
Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

CONTINUE ON A SEPARATE PAGE IF NEEDED.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(K) SUMMARY K211078

A. SUBMITTER INFORMATION (807.92(a)(1))

Prepared by: Vaibhav Sivaramakrishan

Sr. Regulatory Affairs Specialist

Terumo Medical Corporation

Phone#: (609) 613-3958

Fax#: (410) 398-6079

Prepared for: Owner/Operator

TERUMO CLINICAL SUPPLY CO., LTD.

3, Kawashima-Takehayamachi, Kakamigahara,

Gifu, Japan, 501-6024

Owner/Operator Number: 100 396 57

Manufacturer and Sterilization Facility (Applicant)

TERUMO CLINICAL SUPPLY CO., LTD.

3, Kawashima-Takehayamachi, Kakamigahara,

Gifu, Japan, 501-6024

Registration Number: 3009500972

Contact Person: Vaibhav Sivaramakrishan

Sr. Regulatory Affairs Specialist

Terumo Medical Corporation

265 Davidson Avenue, Suite 320

Somerset, NJ 08873

Phone#: (609) 613-3958

Fax#: (410) 398-6079

Email: vai.sivaramakrishan@terumomedical.com

Date prepared: December 21, 2021



B. DEVICE NAME (807.92(a)(2))

Proprietary Name: Progreat Lambda

Common Name: Micro Catheter System

Classification Name: Diagnostic Intravascular Catheter (DQO),

Continuous flush catheter (KRA)

Classification Panel: Cardiovascular

Regulation: 21 CFR 870.1200

Product Code: DQO, KRA

Classification: Class II

C. PREDICATE DEVICE (807.92(a)(3))

The legally marketed device to which substantial equivalence is claimed is:

Primary Predicate Device:

1. K033583: PROGREAT, manufactured by Ashitaka Factory of Terumo Corporation

Reference Devices:

- 1. K173548: Merit Pursue Microcatheter, Merit Medical Systems, Inc.
- 2. K201792: TRUSELECT Microcatheter, Boston Scientific Corp.

D. REASON FOR 510(k) SUBMISSION

This Traditional 510(k) is being submitted for Progreat Lambda for the purposes of establishing substantial equivalence to a legally marketed predicate device.

E. DEVICE DESCRIPTION (807.92(a)(4))

Principle of Operation Technology

Progreat Lambda submitted in this 510(k) and its predicate Progreat (K033583) are operated by a manual process.



Design/Construction

Progreat Lambda is a single use, ethylene oxide sterilized device that is intended for the infusion of contrast media, or embolic materials for hemostasis, into the peripheral vasculature, excluding the blood vessels belonging to the central circulatory system. Progreat Lambda is also indicated for drug infusion in intra-arterial therapy in the peripheral vasculature. Progreat Lambda should not be used in cerebral or coronary vessels.

The catheter consists of metal wire mesh reinforced multi-layer polymer tubing. The mesh is embedded in the catheter wall the entire length of the catheter with the exception of the distal tip section. This increases the flexibility, kink resistance, and pressure resistance of the catheter. The inner layer of the catheter is made of PTFE (polytetrafluoroethylene) to ensure smooth movement of devices such as the guide wire.

The outer surface of the catheter is coated with a hydrophilic polymer with the exception of the proximal end that is 60cm from the catheter hub. The coating becomes lubricious when wet with saline solution or blood.

The device is offered in effective lengths of 110, 130, 150, 165 and 175 cm. French size and shaft outer and inner diameter are given in Table 5.1.

Table 5.1: Catheter size

French Size	Shaft Outer Diameter (mm)		Shaft Inner Di	ameter (mm)
French Size	Distal part Proximal part		Distal part	Proximal part
1.7 Fr.	0.57	0.94	0.43	0.58
1.9 Fr.	0.64	0.94	0.48	0.60



Materials

The materials for Progreat Lambda are provided in Table 5.2.

Table 5.2: List of Materials

No.	Name of Component		Component	Raw material		
1*			Outer layer	Polyester Elastomer		
1		Shaft	Outer layer	Pigment		
2		Shart	Reinforcement wire	Tungsten		
3*	Catheter		Inner layer	Polytetrafluoroethylene		
4*	Cameter	D:-4-14:-	Outon lavon	Polyester Elastomer		
4			•	Distal tip Section Outer layer	Pigment	
5*		Section	Inner layer	Polytetrafluoroethylene		
6		Radiopaque marker		Pt-Ir alloy		
7*	Hydrophilic polymer coating		Hydrophilia polymer coating		olymer coeting	Dimethyl acrylamide - glycidyl
1			orymer coating	methacrylate copolymer		
8	Quick-drying glue		Quick-drying glue Cyanoacrylate			
9*	Catheter hub		Cosh oton hub		atar huh	Polyamide
9	Cameter nuo		etel muu	Silicone		
10	Catheter strain relief tube		ain relief tube	Polyester elastomer		
10			ani ichci tube	Pigment		

^{*}Blood contacting material.

Specifications

The specifications for Progreat Lambda are provided in Table 5.3.

Table 5.3: Progreat Lambda Specifications

French Size		1.7 Fr.	1.9 Fr.
Catheter I.D. (mm)	Distal	0.43	0.48
Cameter 1.D. (mm)	Proximal	0.58	0.60
Catheter O.D. (mm)	Distal	0.57	0.64
	proximal	0.94	0.94
Effective length (cm)*		110, 130, 15	50, 165, 175
Coating length(cm)		50, 70, 90, 105, 115	
Maximum guidewire outer diameter		0.0	16"

^{*}The length from the proximal catheter strain relief tube to the catheter distal tip.



F. INDICATIONS FOR USE (807.92(a)(5))

Progreat Lambda is intended for the infusion of contrast media, or embolic materials for hemostasis, into the peripheral vasculature, excluding the blood vessels belonging to the central circulatory system. Progreat Lambda is also indicated for drug infusion in intra-arterial therapy in the peripheral vasculature. Progreat Lambda should not be used in cerebral or coronary vessels.

The indications for use are equivalent to the predicate (K033583).

G. SUBSTANTIAL EQUIVALENCE COMPARISON (807.92(a)(6))

Progreat Lambda, the subject of this Traditional 510(k), is substantially equivalent in its intended use, technology/principle of operation, materials, and performance to the predicate, K033583 –Progreat, manufactured by Ashitaka Factory of Terumo Corporation.

In addition to the above-listed primary predicate, TERUMO CLINICAL SUPPLY CO., LTD. has identified the following reference devices. These are market leading devices with the same intended use and basic design as the subject device. Because these devices are frequently used in clinical practice, TERUMO CLINICAL SUPPLY CO., LTD. felt it was appropriate to use them as references when setting the acceptance criteria for Progreat Lambda performance testing.

- 1. Merit Medical Systems, Inc: Merit Pursue Microcatheter (K173548)
- 2. Boston Scientific Corp.: TRUSELECT Microcatheter (K201792)

A comparison of the technological characteristics is summarized in Table 5.4.

 Table 5.4: Summary of Comparative Information

	Subject Device:	Predicate Device:	Reference Device #1:	Reference Device #2:
Device Characteristic	Progreat Lambda	Progreat	Merit Pursue Microcatheter	TRUSELECT™ Microcatheter
16 0	menta to or nargar	K033583	K173548	K201792
Manufacturer	TERUMO CLINICAL	Ashitaka Factory of	Merit Medical Systems,	Boston Scientific
	SUPPLY CO., LTD.	Terumo Corporation	Inc.	Corporation
Intended Use	Progreat Lambda is	Progreat is intended for	The Microcatheter is	The TRUSELECT TM
/Indications for Use	intended for the infusion	the infusion of contrast	intended for general	Microcatheters are
	of contrast media, or	media into all peripheral	intravascular use,	intended for peripheral
	embolic materials for	vessels up to and	including peripheral and	vascular use. The
	hemostasis, into the	including the cervical	coronary vasculature.	microcatheter can be used
	peripheral vasculature,	vessels, all vessels in the	Once the subselective	for selective infusion of
	excluding the blood	lower and upper	region has been accessed,	diagnostic, embolic, or
	vessels belonging to the	extremities, and all	the Microcatheter can be	therapeutic materials into
	central circulatory	coronary vessels. Progreat	used for the controlled	the vessel.
	system. Progreat Lambda	is also intended for drug	and selective infusion of	
	is also indicated for drug	infusion in intra-arterial	diagnostic, embolic, or	
	infusion in intra-arterial	therapy and the infusion	therapeutic materials into	
	therapy in the peripheral	of embolic materials for	vessels. The catheter	
	vasculature. Progreat	hemostasis. Progreat	should not be used in the	
	Lambda should not be	should not be used in	cerebral vessels.	
	used in cerebral or	cerebral vessels.		
	coronary vessels.			
Operation Principle	Manual	Same	Same	Same

	Subject Device:	Predicate Device:	Reference Device	Reference Device #2:
Device Characteristic	Progreat Lambda	Progreat	#1: Merit Pursue Microcatheter	TRUSELECT™ Microcatheter
		K033583	K173548	K201792
Design/Construction	The catheter consists of inner layer, reinforcement wire, outer layer, radiopaque marker, hydrophilic polymer coating, catheter hub and catheter strain relief tube.	Same	Same	Same
Materials	Outer tube*: Polyester Elastomer Inner tube*: PTFE Reinforcement wire: Tungsten	Outer layer*: Polyester elastomer (distal) and Polyurethane elastomer (proximal) Inner layer*: Polytetrafluoroethylene Reinforcing coil: Tungsten	Information not publicly available.	Information not publicly available.
	Radiopaque marker: Pt-Ir alloy Coating*: Dimethyl acrylamide glycidyl methacrylate copolymer, Silicone Catheter hub*: Polyamide, Silicone Catheter strain relief tube: Polyester elastomer *: blood contacting material	Radiopaque maker: Pt-Ir alloy Hydrophilic coating*: Dimethyl acrylamide — glycidyl methacrylate copolymer Hub*: Nylon Anti-kink protector: Nylon elastomer *: blood contacting material		

Device Characteristic	Subject Device:	Predicate Device:	Reference Device #1:	Reference Device #2:
	Progreat Lambda	Progreat	Merit Pursue Microcatheter	TRUSELECT™ Microcatheter
		K033583	K173548	K201792
Package	 Individual package on which the product label and the peel-off labels are attached 1 unit per package 	Same	Same	Same
Specifications	 Effective lengths: 110, 130, 150, 165, 175 cm French size: 1.7Fr., 1.9Fr. O.D.(Distal/Proximal) 1.7Fr./2.8Fr.:0.57/0.94mm 1.9Fr./2.8Fr.:0.64/0.94mm 	 Effective lengths: 100, 110, 130, 150 cm French size: 2.0, 2.4, 2.7, 2.8 Fr. O.D.(Distal/Proximal) 2.0Fr./2.7Fr (0.67/0.90mm) 2.4/2.9Fr. (0.80/0.97mm) 2.7/2.9Fr. (0.90/0.97mm) 2.8/3.0Fr. (0.93/1.00mm) 	 Effective Length: 110, 130, 150cm French size: 1.7Fr., 2.0Fr. O.D.(Distal/Proximal) 1.7Fr./2.8Fr. 2.0Fr./2.9Fr. 	 Effective Length: 105, 130, 150, 175cm French size: 2.0Fr. O.D.(Distal/Proximal) 2.0Fr./2.8Fr.

Device Characteristic	Subject Device:	Predicate Device:	Reference Device #1:	Reference Device #2:
	Progreat Lambda	PROGREAT K033583	Merit Pursue Microcatheter K173548	TRUSELECT™ Microcatheter K201792
Specifications	• I.D. 1.7Fr.: 0.43mm 1.9Fr.: 0.48mm • Maximum Guide Wire outer diameter: 0.016"	• I.D. 2.0Fr.:0.019"/0.49mm 2.4Fr.:0.022"/0.57mm 2.7Fr.:0.025"/0.65mm 2.8Fr.:0.027"/0.70mm • Maximum Guide Wire outer diameter: 2.0Fr.type: 0.016" 2.4Fr.type: 0.018" 2.7Fr.type: 0.021" 2.8Fr.type: 0.021"	• I.D. 1.7Fr.type: 0.016" (0.40 mm) 2.0Fr.type: 0.020" (0.50 mm) •Maximum guidewire outer diameter: 1.7Fr.type: 0.014" 2.0Fr.type: 0.018"	• I.D. 0.021"(0.53mm) •Maximum guidewire outer diameter: 0.016" or 0.014"
Specifications	 Distal tip shape: Straight/Angle/Triple angle Maximum injection pressure: 900 psi 	 Distal tip shape: straight/angled Maximum injection pressure: 2.0Fr. : 750psi 2.4, 2.7Fr.: 750psi 2.8Fr. : 900 psi 	 Distal tip shape: straight/angled Maximum injection pressure: 800 psi 	 Distal tip shape: straight/angled(Bernshape) Maximum injection pressure: 800 psi
Sterilization	Ethylene oxide	Same	Same	Same
Shelf life	2 years	2 years	Information not publicly available.	Information not publicly available.

H. NON CLINICAL TESTS (807.92(b)(1))

Performance Testing

Performance testing was conducted to ensure the safety and effectiveness of Progreat Lambda throughout the shelf life, verify conformity to the applicable external and internal standards, and demonstrate substantial equivalence to the predicate device. With the exception of the Radio-detectability¹ test and Embolic device compatibility for coil and microsphere², the following performance tests were performed on nonaged and accelerated aged samples. Table 5.5 provides a list of performance tests that were performed on Progreat Lambda.

Table 5.5: Summary of Performance Testing

· · · · · · · · · · · · · · · · · · ·
Test Item
Radio-detectability
Surface
Peak tensile force
Freedom from leakage
Fluid leakage (Hub)
Sub-atmospheric pressure air leakage (Hub)
Stress cracking (Hub)
Resistance to separation from axial load (Hub)
Resistance to separation from unscrewing (Hub)
Resistance to overriding (Hub)
Power injection
Distal tip
Particulate evaluation
Coating integrity
Torque strength
Distal marker strength
Product dimension
Embolic device compatibility
Flexibility and kink test
Wire compatibility
Simulated Use

¹ Only non-aged sample was tested since the amount of metallic material contained in the product would not change over time.

² Only non-aged sample was tested since the inner diameter would not changed over time.

Performance testing met the predetermined acceptance criteria and is acceptable for clinical use throughout its shelf life.

Biocompatibility

In accordance with ISO 10993-1, Progreat Lambda is classified as: Externally Communicating Device, Circulating Blood, Limited Contact (<24 hours). The finished device's patient contacting parts were tested in accordance with the tests recommended in the FDA *Guidance for Industry and Food and Drug Administration Staff - Use of International Standard ISO-10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process."*Screening tests were performed on accelerated aged devices to show that the biocompatibility is maintained throughout the shelf life of the product. Table 5.6 provides a list of biocompatibility tests conducted on Progreat Lambda.

Table 5.6: Summary of ISO 10993 Biocompatibility Testing

Non-aged, sterile, whole device
Cytotoxicity
Sensitization
Intracutaneous Reactivity
Acute Systemic Toxicity
Pyrogenicity
Hemolysis
Thrombogenicity
Complement Activation (Immunology)
Physicochemical Profile (Physicochemical and FT-IR)
Accelerated-aged (3 years), sterile, whole device
Cytotoxicity
Hemolysis
Physicochemical Profile (Physicochemical and FT-IR)

Results of the testing demonstrate that the device is biocompatible throughout the shelf life of the product.

Sterilization

The sterility of the device is assured using a sterilization method validated in accordance with ISO 11135:2014, *Sterilization of Health Care Products – Ethylene Oxide – Requirements for the development, validation and routine control of a sterilization process for medical devices*, to provide a Sterility Assurance Level (SAL) of 10⁻⁶.

I. CLINICAL TESTS (807.92(b)(2))

This 510(k) does not include data from clinical tests.

J. CONCLUSION (807.92(b)(3))

In summary, Progreat Lambda, subject of this 510(k), is substantially equivalent in its intended use, technology/principle of operation, materials, and performance to the predicate, K033583–Progreat, manufactured by Ashitaka Factory of Terumo Corporation.