

March 24, 2022

Graftworx, Inc. dba Alio Allison Komiyama, PhD, RAC Principal Consultant RQM+ 2251 San Diego Ave. Suite B-257 San Diego, California 92110

Re: K211365

Trade/Device Name: Alio Medical Remote Monitoring System Regulation Number: 21 CFR 870.2910 Regulation Name: Radiofrequency Physiological Signal Transmitter And Receiver Regulatory Class: Class II Product Code: DRG, DQD Dated: March 21, 2022 Received: March 23, 2022

Dear Allison Komiyama:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jennifer Shih Kozen Assistant Director External Heart Rhythm and Rate Team Division of Cardiac Electrophysiology, Diagnostics and Monitoring Devices Office of Cardiovascular Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K211365

Device Name Alio Medical Remote Monitoring System

Indications for Use (Describe)

The Alio Medical Remote Monitoring System is a wireless remote monitoring system intended for use by healthcare professionals to intermittently collect physiological data in home use settings. The data includes skin temperature, auscultation sound data and heart rate. Data is transmitted wirelessly from the SmartPatch wearable sensor to a web-based portal for the healthcare provider's (HCP) review.

The Alio Medical RMS is intended for use on general care patients who are 18 years of age or older. The SmartPatch sensor is indicated to measure skin temperature and pulse rate where clinically indicated. The SmartPatch sensor is indicated to record and transmit auscultation sound data where clinically indicated.

The device is not intended for use in critical care or other high-acuity environments. The Alio Medical RMS is a secondary, adjunct patient monitor and is not intended to replace existing standard-of-care patient monitoring practices.

Type of Use (Select one or both, as applicable)		
Rescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)	

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary

1. General Information

510(k) Sponsor	Alio, Inc.	
Address	544B Bryant St	
	San Francisco, CA 94107	
Correspondence Person	Allison C. Komiyama, Ph.D., R.A.C.	
	RQM+	
Contact Information	Email: akomiyama@rqmplus.com	
	Phone: +1 (412) 816-8253	
Date Prepared	March 24, 2022	

2. Proposed Device

Proprietary Name	Alio Medical Remote Monitoring System	
Common Name	Alio Medical RMS	
Classification Name	Radiofrequency physiological signal transmitter and receiver	
Regulation Number	21 CFR 870.2910	
Product Code	DRG, DQD	
Regulatory Class	П	

3. Predicate Device

Proprietary Name	Vital Connect Platform
Premarket Notification	K152139
Classification Name	Radiofrequency physiological signal transmitter and receiver
Regulation Number	21 CFR 870.2910
Product Code	DRG, DSI, MHX
Regulatory Class	П

4. Reference Device

Proprietary Name	Eko Core
Premarket Notification	K200776
Classification Name	Stethoscope
Regulation Number	21 CFR 870.1875
Product Code	DQD
Regulatory Class	II



5. Reference Device

Proprietary Name	BB-613 WP	
Premarket Notification	K190792	
Classification Name	Oximeter, non-invasive blood pressure measurement system	
Regulation Number	21 C.F.R. 870.2700 Oximeter, 870.1130 Noninvasive blood pressure	
	measurement system	
Product Code	DQA, DXN, DRG	
Regulatory Class	11	

6. Device Description

Alio Medical Remote Monitoring System, or "Alio Medical RMS", utilizes a wearable device (SmartPatch) on the skin to gather physiological data and then transmits it to a device (Bedside Hub) located in the subject's home. The Bedside Hub then relays this raw data to the Alio Medical Cloud where it is processed and analyzed using Alio's proprietary algorithms. Data is accessible to Healthcare Professionals and the Alio clinical team via a web-based Clinician Portal. The SmartPatch and Bedside Hub are intended to be used on general care patients who are 18 years of age or older in a non-clinical environment. The web-based Clinical Portal is to be used by healthcare professionals in an office environment.

The Alio Medical Remote Monitoring System includes the following components:

- SmartPatch
- Bedside Hub
- Alio Medical Cloud (backend only not user facing)
- Clinician Portal

SmartPatch

A flexible, silicone-encased patch that can be worn where clinically indicated for up to seven days at a time. It houses numerous sensor technologies, which include a microphone, accelerometer, temperature sensors, and a PPG sensor. The sensors collect physiological data including skin temperature, auscultation sound data, and heart rate. Data is transmitted to the Cloud, via the Hub, where it is analyzed and sent to a Healthcare Professional via the Web Portal.

Bedside Hub

The Bedside Hub has the form and finish of an at-home device. It automatically communicates with the activated SmartPatch and uploads physiological data to the Alio Medical Cloud.



Alio Medical Cloud

The Cloud features a database that supports storage, analytics, system monitoring and visualization capabilities. The Alio Medical Cloud is encrypted and HIPAA compliant. All patient data is fully traceable to device and patient ID via the database.

Clinician Portal

The Clinician Portal is the interface tool between a user (healthcare professional users only) and the system that enables the user to visualize and interact with data being generated by the system.

Note - none of the above components/accessories of the Alio Medical RMS have received prior 510(k) clearance.

7. Intended Use/ Indications for Use

The Alio Medical Remote Monitoring System is a wireless remote monitoring system intended for use by healthcare professionals to intermittently collect physiological data in home use settings. The data includes skin temperature, auscultation sound data and heart rate. Data is transmitted wirelessly from the SmartPatch wearable sensor to a web-based portal for the healthcare provider's (HCP) review.

The Alio Medical RMS is intended for use on general care patients who are 18 years of age or older. The SmartPatch sensor is indicated to measure skin temperature and pulse rate where clinically indicated. The SmartPatch sensor is indicated to record and transmit auscultation sound data where clinically indicated.

The device is not intended for use in critical care or other high-acuity environments. The Alio Medical RMS is a secondary, adjunct patient monitor and is not intended to replace existing standard-of-care patient monitoring practices.



8. Substantial Equivalence

Feature/ Function	Device: Alio Medical RMS	Predicate Device: Vital Connect Platform (K152129)
	The Alio Medical Remote	The Vital Connect Platform is a wireless
	Monitoring System is a wireless	remote monitoring system intended for
	remote monitoring system	use by healthcare professionals for
	intended for use by healthcare	continuous collection of physiological
	professionals to intermittently	data in home and healthcare settings.
	collect physiological data in home	This can include heart rate,
	use settings. The data includes skin	electrocardiography (ECG), heart rate
	temperature, auscultation sound	variability, R-R interval, respiratory rate,
	data and heart rate. Data is	skin temperature, activity (including step
	transmitted wirelessly from the	count), and posture (body position
	SmartPatch wearable sensor to a	relative to gravity including fall). Data
	web-based portal for the	are transmitted wirelessly from the Vital
	healthcare provider's (HCP) review.	Connect Sensor for storage and analysis.
	The Alio Medical RMS is intended	The Vital Connect Platform can include
Indications	for use on general care patients	the ability to notify healthcare
for Use	who are 18 years of age or older.	professionals when physiological data
	The SmartPatch sensor is indicated	fall outside selected parameters. The
	to measure skin temperature and	device is intended for use on general
	pulse rate where clinically	care patients who are 18 years of age or
	indicated. The SmartPatch sensor is	older as a general patient monitor, to
	indicated to record and transmit	provide physiological information. The
	auscultation sound data where	data from the Vital Connect Platform are
	clinically indicated.	intended for use by healthcare
	The device is not intended for use	professionals as an aid to diagnosis and
	in critical care or other high-acuity	treatment. The device is not intended
	environments.	for use on critical care patients.
	The Alio Medical RMS is a	
	secondary, adjunct patient monitor	
	and is not intended to replace	
	existing standard-of-care patient	
	monitoring practices.	
Rx or OTC	Rx	Rx
Intended	Adult - 18+	Adult - 18+
Intended	Ноте	
Environment	- Home	Home & healthcare settings
Application	Patch - SmartPatch only	Patch
Single-use	Yes - SmartPatch only Yes - patch only	



		- heart rate
	Skin temperature	- ECG
		- heart rate variability
		- B-R interval
Data	Electronic stethoscope for sound	- respiratory rate
transmitted	auscultation data	- skin temperature
		- activity (including step count)
	Heart rate	- nosture (hody position relative to gravity
		including fall)
	Thermister	ECC electrodes to detect heart rate
	Microphone	- 2-avis MEMS accelerometer to detect
Sensor Type	Accelerometer	- 3-axis MENIS acceleronieter to detect
	Acceleronieter	Thermisters to detect hody temperature
Tananakana	PPG Claim to react unce	- memistors to detect body temperature
Temperature		
Data	15 C = 50 C	15 C - 50 C
	- neart	
Cound Date	- lungs	21/2
Sound Data	- bower	N/A
	- arteries	
Carried	- veins	
Sound	Yes	N/A
Amplification		
Record and	Vac (Dischards vie Clinician Dantel)	21/2
Раураск	Yes (Playback via Clinician Portal)	N/A
Operating		
Operating	Intermittent	Continuous
WIDUE	Bluetooth - SmartPatch only	
Data	Badio Frequency: 2.4 - 2.5 GHz	Bluetooth
Transmission	(Cellular) - SmartHub only	Bidetooth
	SmartDatch only hattory (Li ion)	Sonsor only
Energy Source	Hub - 5V Mains DC	hatteny (Li-jon)
	IEC 60601-1 3rd ed	Dattery (Linony
		10002 1:2000
		ISO 10993-1.2009
Standards	IEC 60086-4:2019	IEC 60601-2-25
	IEC 60601-1-2:2007/2014	IEC 60601-2-47
	IEC 62304:2006/A1:2016	IEC 60601-1-2
	IEC 62366-1:2007/2015	IEC/TS 62657-2
	FCC CRF47	FCC CRF47
	Part 15 Subpart B	Part 15 Subpart C
	ISO 10993-5:2009	
	ISO 10993-10:2010	



Feature/ Function	Device: Alio Medical RMS	Reference Device: BioBeat BB-613P (K190792)
Application	Skin - SmartPatch only	Wrist area and skin
Monitoring	Intermittent	Spot Check/Intermittent
Principle of	Pulse reflectance technology.	Pulse reflectance technology, Four LED (red
Operation	Light source is comprised of 6	+ IR) and photo diode absorbs reflected
	LEDs (2x green, 2x red, 2x IR),	light. Tracking changes of blood pressure is
	reflected light is detected by 4	done by pulse wave transit time (PWTT)
	photodiodes. An accelerometer	which is obtained utilizing pulse
	provides additional heart rate	measurements from the integrated skin
	data.	attached SpO2 sensor
Emitted light peak	530 nm (Green), 650 nm (Red),	880nm (IR), 650nm (Red)
wavelength	940 nm (IR)	
Measurement	30-200 bpm	40-250bpm
Range, HR		
A _{rms} , HR	5 bpm	±3%
Application Method	Biocompatible adhesive patch	Biocompatible adhesive patch

Technological Characteristic Comparison with Reference Devices

Feature/ Function	Device: Alio Medical RMS	Reference Device: EKO Device (K200776)
Stethoscope Type	Standalone electronic stethoscope	Attachment to an analog stethoscope
Auscultation Data	Yes	Yes
Measurement Frequency Range	55 - 2000 Hz	unknown
Sound Amplification	Yes	Yes
Connectivity	Bluetooth	Bluetooth
Record and Save	Yes	Yes

9. Performance Data

Nonclinical verification and validation test results established that the device meets its design requirements and intended use, that it is as safe, as effective, and performs as well as the predicate devices, and that no new issues of safety and effectiveness were raised. The Alio



Medical RMS was designed, verified, and validated according to the company's Design Control process and has been subjected to extensive safety and performance testing as shown in the test results provided in this submission. Verification and Validation testing data demonstrate that the device meets all of its specifications including compliance with the following standards:

Safety

- IEC 60601-1 3rd Ed.
- IEC 60601-1-11:2010
- IEC 80601-2-56:2017
- IEC 80601-2-61:2017
- IEC 62471:2008
- IEC 60529:2013
- IEC 60086-4:2019

EMC

- IEC 60601-1-2:2007/2014
- FCC Part 15 Radio Frequency Devices, Subpart B Unintentional Radiators

Software

- IEC 62304:2006/A1:2016
- FDA Guidance document, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" and "Content of Premarket Submission for Management of Cybersecurity in Medical Devices."

Usability

- IEC 62366-1:2007/2015
- FDA Guidance document, "Applying Human Factors and Usability Engineering to Medical Devices"

Biocompatibility

- ISO 10993-5:2009
- ISO 10993-10:2010

10. Conclusion

Based on the information submitted in this premarket notification, and based on the indications for use, technological characteristics and performance testing, Alio Medical Remote Monitoring System raises no new questions of safety and effectiveness and is substantially equivalent to the predicate devices in terms of safety, efficacy, and performance.