

July 25, 2023

DEKA Research and Development Paul Smolenski Regulatory Affairs 340 Commercial Street Manchester, New Hampshire 03101

Re: K213536

Trade/Device Name: DEKA ACE Pump System

Regulation Number: 21 CFR 880.5730

Regulation Name: Alternate Controller Enabled Infusion Pump

Regulatory Class: Class II Product Code: QFG, NDC Dated: November 5, 2021 Received: November 5, 2021

Dear Paul Smolenski:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

K213536 - Paul Smolenski Page 2

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Joshua Balsam -S

Joshua Balsam, Ph.D.
Branch Chief
Division of Chemistry
and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

K213536

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

CONTINUE ON A SEPARA	TE PAGE IF NEEDED.
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)
Type of Use (Select one or both, as applicable)	
The bolus calculator is indicated for use for aiding the user in dediabetes mellitus based on consumed carbohydrates, operator-escarbohydrate ratio, target glucose values, and current insulin on	ntered blood glucose, insulin sensitivity, insulin to
Indications for Use (Describe) The DEKA ACE Pump System is intended for the subcutaneou management of diabetes mellitus in persons requiring insulin, a securely communicate with compatible, digitally connected devereceive, execute, and confirm commands from these devices. The requires a prescription.	ges 13 and above. The pump is able to reliably and rices, including automated insulin dosing software, to
DEKA ACE Pump System	

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary: K213536

Submitter Information

510(k) Sponsor: DEKA Research & Development

340 Commercial Street Manchester, NH 03101

Contact Person: Paul Smolenski

Regulatory Affairs

DEKA Research & Development

Phone: (603) 669-5139 Fax: (603) 624-0573

psmolenski@dekaresearch.com

Date Prepared: 07/24/2023

Proposed Device

Common/Usual Name: ACE Pump

Trade/Proprietary Name: DEKA ACE Pump System

Classification Name: Alternate Controller Enabled ACE Pump; Calculator, Drug Dose

Device Classification: 880.5730; 868.1890

Product Code: QFG; NDC

Class: II

Device Panel: Clinical Chemistry

Predicate Device

The predicate device for this submission is the Tandem t:Slim X2 insulin pump with interoperable technology granted under De Novo DEN180058.

Device Description

The DEKA ACE Pump System is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The Pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The Pump is intended for single patient, home use and requires a prescription.

The system as described in this submission is able to be integrated with a Dexcom G6 interoperable Continuous Glycemic Controller (iCGM). This submission also details the integration process that can be used to incorporate an iAGC.

The DEKA ACE Pump System consists of the following components:

1. **Pump:** A durable pump that incorporates fluid delivery algorithms and interfaces to an DEKA ACE Pump cassette, Remote Interface, iCGM, and iAGC. The pump is powered by a rechargeable lithium ion battery.

- 2. Cassette: A single-use pumping cassette that combines microfluidic valves, a pump chamber, insulin reservoir, and Acoustic Volume Sensing (AVS) measurement chamber. The cassette interfaces to an DEKA ACE Pump and off-the-shelf infusion set.
- **3. Remote Interface (Controller):** A wireless controller that serves as the user interface to the DEKA ACE Pump system. This includes a large color touch display for ease of use.

Information is being supplied in this 510(k) premarket submission to demonstrate that the device is substantially equivalent in safety and effectiveness through comparison of indications for use and technological characteristics to the predicate Tandem t:Slim X2 insulin pump with interoperable technology granted on 12/03/2019 under De Novo DEN180058. As described throughout this submission, the subject DEKA ACE Pump system meets the product definition and all of the Special Controls defined in DEN180058 and 21 CFR 880.5730 for Alternate Controller Enabled Insulin Infusion Pumps, Product Code QFG.

Indications for Use

The DEKA ACE Pump System is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription.

The bolus calculator is indicated for use for aiding the user in determining the bolus insulin dosage for management of diabetes mellitus based on consumed carbohydrates, operator-entered blood glucose, insulin sensitivity, insulin to carbohydrate ratio, target glucose values, and current insulin on board.

Substantial Equivalence Discussion

Intended Use Comparison

The table below includes a summation matrix of the intended use between the new device and those of the current device:

Indications for Use The t:slim X2 insulin pump with interoperable technology (the Pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The Pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The Pump is intended for use with NovoLog or Humalog U-100 insulin. The Pump is indicated for use in dicated for use in the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The Pump is intended for single patient, home use and requires a prescription. The DEKA ACE Pump System is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing sutomated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription. The DEKA ACE Pump Subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing sutomated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription.	Characteristic	Predicate	Device
Indicated for use in the user in determining the		The t:slim X2 insulin pump with interoperable technology (the Pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The Pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The Pump is intended for single patient, home use and requires a prescription. The Pump is indicated for use with NovoLog or Humalog U-100 insulin. The Pump is	The DEKA ACE Pump System is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription. The bolus calculator is indicated for use for aiding

	individuals 6 years of age and greater.	management of diabetes mellitus based on consumed carbohydrates, operator-entered blood glucose, insulin sensitivity, insulin to carbohydrate ratio, target glucose values, and current insulin on board.
Prescription Only or Over the Counter	Prescription Only	Same
Intended Population	Persons with Diabetes Mellitus Ages 6 and up	Persons with Diabetes Mellitus Ages 13 and up
Environment of Use	In professional healthcare facility and home healthcare environments	Same

Discussions of differences in Indications of Use statement

The DEKA ACE Pump System is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin, ages 13 and above. The pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The pump is intended for single patient, home use and requires a prescription.

The bolus calculator is indicated for use for aiding the user in determining the bolus insulin dosage for management of diabetes mellitus based on consumed carbohydrates, operator-entered blood glucose, insulin sensitivity, insulin to carbohydrate ratio, target glucose values, and current insulin on board.

The predicate's indications for use from DEN180058 are as follows:

The t:slim X2 insulin pump with interoperable technology (the Pump) is intended for the subcutaneous delivery of insulin, at set and variable rates, for the management of diabetes mellitus in persons requiring insulin. The Pump is able to reliably and securely communicate with compatible, digitally connected devices, including automated insulin dosing software, to receive, execute, and confirm commands from these devices. The Pump is intended for single patient, home use and requires a prescription. The Pump is indicated for use with NovoLog or Humalog U-100 insulin. The Pump is indicated for use in individuals 6 years of age and greater.

Discussions of differences in intended population

Both the predicate and subject device are intended for use in persons with Diabetes Mellitus. The predicate device is indicated for ages 6 and older and the subject device is indicated for ages 13 and older. The difference in the lower age limit does not impact the safety and effectiveness of the device for its indicated use. Testing demonstrates equivalent safety and effectiveness for the

indicated population.

Discussions of differences in environment of use

Both the predicate and subject devices are intended to be used in professional healthcare facility and home healthcare environments. There are no differences in the environments of use.

Comparison of Technological Characteristics with the Predicate Device

The below table compares the characteristics of the subject device to the predicate, t:slim X2 insulin infusion pump, and includes an assessment of differences between them and why the differences between the subject device and the predicate device do not introduce new or different questions of safety or effectiveness.

Substantial Equivalence Discussion

The table below compares the intended use and technological characteristics of the subject device with that of the predicate device:

Characteristic	Predicate Device	Subject Device	Equivalence
Device	Alternate Controller Enabled	Alternate Controller Enabled	Same
Classification	Infusion Pump cleared under 21	Infusion Pump submitted under	
Regulation and	CFR 880.5730, Procode QFG	21 CFR 880.5730, Procode QFG	
Product Code			
Indications for Use	The t:slim X2 insulin pump with	The DEKA ACE Pump System is	Indications for use of the DEKA
	interoperable technology (the	intended for the subcutaneous	ACE Pump System are equivalent
	Pump) is intended for the	delivery of insulin, at set and	to the predicate.
	subcutaneous delivery of insulin,	variable rates, for the management	
	at set and variable rates, for the	of diabetes mellitus in persons	Indications for use of the bolus
	management of diabetes mellitus	requiring insulin, ages 13 and	calculator are equivalent to other
	in persons requiring insulin. The	above.	devices with this product code.
	Pump is able to reliably and		
	securely communicate with	The pump is able to reliably and	
	compatible, digitally connected	securely communicate with	
	devices, including automated	compatible, digitally connected	
	insulin dosing software, to receive,	devices, including automated	
	execute, and confirm commands	insulin dosing software, to receive,	
	from these devices. The Pump is	execute, and confirm commands	
	intended for single patient, home	from these devices. The pump is	
	use and requires a prescription.	intended for single patient, home	
	The Pump is indicated for use with	use and requires a prescription.	
	NovoLog or Humalog U-100		
	insulin. The Pump is indicated for	The bolus calculator is indicated	
	use in individuals 6 years of age	for use for aiding the user in	
	and greater.	determining the bolus insulin	
		dosage for management of	
		diabetes mellitus based on	
		consumed carbohydrates,	
		operator-	

Characteristic	Predicate Device	Subject Device	Equivalence
		entered blood glucose, insulin	
		sensitivity, insulin to carbohydrate	
		ratio, target glucose values, and	
		current insulin on	
		board.	
Prescription Use	Yes	Yes	Same
Intended Population	6 years and older	13 years and older	The difference in the lower age limit does not impact the safety and effectiveness of the device for its indicated use. Testing demonstrates equivalent safety and effectiveness for the indicated population.
Patient	On-body wearable ambulatory	On-body wearable ambulatory	Same
Environment	pump	pump	
Environment of Use	In professional healthcare facilities and home environments	In professional healthcare facilities and home environments	Same
Delivery Method	Micro-dosing threaded cartridge pump	Microprocessor controlled Micro- dosing pump mechanism supplemented with acoustic volume sensor (AVS) feedback for monitoring delivery accuracy	No impact to safety or effectiveness. Subject device meets all Special Controls requirements. Delivery method has been found SE to ambulatory pumps through clearance of reference device.
Insulin Basal Rate Delivery Range	0 units /hour- 15units/hour	0 units/hour - 30 units/hour	No impact on safety or effectiveness. Subject device meets all Special Controls requirements across the entire delivery rate range.
Insulin Bolus Delivery Range	0.01 U at volumes greater than 0.05 U units, Max Bolus Volume 25 U	Programmable from 0.05 - 25.00 Units in 0.01 Unit increments.	Same

Characteristic	Predicate Device	Subject Device	Equivalence
Basal Accuracy	See Delivery Accuracy	See Delivery Accuracy	No impact on safety or
	Comparison Below	Comparison Below	effectiveness. Subject device
			meets all Special Controls
			requirements across the entire delivery rate range.
Bolus Accuracy	See Delivery Accuracy	See Delivery Accuracy	No impact to safety or
Bolds Accuracy	Comparison Below	Comparison Below	effectiveness. Subject device
	Comparison Below	Companison Below	meets all Special Controls
			requirements.
Bolus Volume after	Less than 3 Units	No more than 0.74 units	No impact to safety or
Occlusion Release			effectiveness. Subject device
			meets all Special Controls
Time to occlusion	2 min (Dalus): 2 haves (Dass) 2	10 min (Dalva), 2 havra (Dagal 1	requirements.
alarm	3 min (Bolus); 2 hours (Basal, 2 U/hr); 36 hours (Basal, 0.1 U/hr):	10 min (Bolus); 3 hours (Basal, 1 U/h); 6 hours (Basal, 0.1 U/hr)	No impact to safety or effectiveness. Subject device
aiaiiii	C/III), 30 Hours (Basar, 0.1 C/III).	C/II), G Hours (Basar, 0.1 C/III)	meets all Special Controls
			requirements.
Material	Compliant with ISO-10993	Compliant with ISO-10993	Same
Biocompatibility			
Cartridge/Cassette	2 years	1 year	Performance testing over the one
Shelf Life			year shelf life of the cassette
			indicates that the cassette remains
Ingress Protection	IPX7: Watertight to a depth of 3	IP28, indicating protection from	safe and effective. The higher level of ingress
lligless i folection	feet (0.91 meters) for up to 30	continuous immersion in water.	protection than predicate meets
	minutes	The Pump can tolerate immersion	use model requirements.
		to depths of up to 12 feet (3.7 m)	1
		for 1 hour.	
Applicable Safety	• IEC 60601-1	• IEC 60601-1	Same
Standards	• IEC 60601-1-2	• IEC 60601-1-2	
	• IEC 60601-1-8	• IEC 60601-1-8	

Characteristic	Predicate Device	Subject Device	Equivalence
	• IEC 60601-1-11	• IEC 60601-1-11	
	• IEC 60601-2-24	• IEC 60601-2-24	
	• ISO 11137-1 (Sterilized via	• ISO 11137-1 (Sterilized via	
	Gamma Radiation)	Gamma Radiation)	
	• ISO 10993-1	• ISO 10993-1	
	• ISO 14971	• ISO 14971	
Power Source	Rechargeable Lithium Polymer Battery	Rechargeable Lithium Ion Battery	Same
Storage Conditions	Temperature: -4°F (-20°C) to 140°F (60°C)	Temperatures of -25 °C (-13 °F) to 70 °C (158 °F)	No impact to safety or effectiveness. Subject device
	Humidity: 20% to 90% RH non- condensing	Non-condensing humidity 15% to 90%	meets all Special Controls requirements.
Operating Conditions	Temperature: 41°F (5°C) to 98.6°F (37°C) Humidity: 20% to 90% RH non-condensing	Temperatures of 5 °C (41 °F) to 40 °C (104 °F) Non-condensing humidity of 15% to 90%	No impact to safety or effectiveness. Subject device meets all Special Controls requirements.
System User Feedback	Visual, audible, and vibratory	Visual, audio, and vibratory	Same
Battery Operating Time	4 – 7 days	72 hours	No impact to safety or effectiveness. Subject device meets all Special Controls requirements.

Basal and Bolus Accuracy Comparison

Below is a comparison of the basal and bolus accuracies of the subject and predicate devices, as reported in their respective User Guides per the ACE Pump Special controls.

Tandem (DEN180058) Basal Accuracy:

0.1 U/hr Basal Accuracy

Interval	Average Interval	Minimum (U)	Maximum (U)
	(U)		
1 hour	0.12	0.09	0.16
6 hours	0.67	0.56	0.76
12 hours	1.24	1.04	1.48

2.0 U/hr Basal Accuracy

Interval	Average Interval (U)	Minimum (U)	Maximum (U)
1 hour	2.1	2.1	2.2
6 hours	12.4	12.0	12.8
12 hours	24.3	22.0	24.9

15.0 U/hr Basal Accuracy

Interval	Average Interval (U)	Minimum (U)	Maximum (U)
1 hour	15.4	14.7	15.7
6 hours	90.4	86.6	93.0
12 hours	181	175	187

DEKA ACE Pump Basal Accuracy:

0.1 U/hr Basal Accuracy

Interval	Average Interval	Minimum (U)	Maximum (U)
	(U)		
1 hour	0.12	0.09	0.17
6 hours	0.62	0.57	0.66
12 hours	1.22	1.16	1.31

1 U/hr Basal Accuracy

Interval	Average Interval	Minimum (U)	Maximum (U)
	(U)		
1 hour	1.02	0.98	1.09

Interval	Average Interval	Minimum (U)	Maximum (U)
6 hours	6.05	5.84	6.22
12 hours	12.07	11.73	12.33

30 U/hr Basal Accuracy

Interval	Average Interval (U)	Minimum (U)	Maximum (U)
1 hour	30.16	29.80	30.61
6 hours	181.05	178.94	184.46

Tandem (DEN180058) Bolus Accuracy:

				0	.05U bo	lus accu	racy			
	<25%	25- 75%	75- 90%	90- 95%	95- 105%	105- 110%	110- 125%	125- 175%	175- 250%	>250%
Number of boluses	21 / 800	79 / 800	63/ 800	34 / 800	272 / 800	180 / 800	105 / 800	29 / 800	17 / 800	0 / 800
% of boluses	2.6%	9.9%	7.9%	4.3%	34.0%	22.5%	13.1%	3.6%	2.1%	0.0%

				2	2.5U bol	us accur	acy			
	<25%	25-	75-	90-	95-	105-	110-	125-	175-	>250%
		75%	90%	95%	105%	110%	125%	175%	250%	
Number	9 /	14 /	11 /	8 /	753 /	5 /	0 /	0 /	0 /	0 /
of boluses	800	800	800	800	800	800	800	800	800	800
% of	1.1%	1.8%	1.4%	1.0%	94.1%	0.6%	0.0%	0.0%	0.0%	0.0%
boluses	1.1 /0	1.070	1.470	1.070	74.170	0.070	0.070	0.070	0.070	0.070

				,	25U bolı	us accur	acy			
	<25%	25-	75-	90-	95- 1050/	105-	110-	125-	175-	>250%
		75%	90%	95%	105%	110%	125%	175%	250%	
Number	0 /	0 /	1 /	3 /	252 /	0 /	0 /	0 /	0 /	0 /
of boluses	256	256	256	256	256	256	256	256	256	256

				,	25U bolı	ıs accur	acy			
% of	0.0%	0.0%	0.4%	1.2%	98.4%	0.0%	0.0%	0.0%	0.0%	0.0%
boluses										

DEKA ACE Pump Bolus Accuracy:

				0	.05U bo	lus accu	racy			
	<25%	25- 75%	75- 90%	90- 95%	95- 105%	105- 110%	110- 125%	125- 175%	175- 250%	>250%
Number of boluses	0 / 800	53 / 800	202 / 800	107 / 800	278 / 800	80 / 800	69 / 800	11 / 800	0 / 800	0 / 800
% of boluses	0.0%	6.6%	25.3%	13.4%	34.8%	10.0%	8.6%	1.4%	0.0%	0.0%

					5U bolus	s accura	cy			
	<25%	25-	75-	90-	95-	105-	110-	125-	175-	>250%
		75%	90%	95%	105%	110%	125%	175%	250%	
Number	0 /	0 /	0 /	0 /	800 /	0 /	0 /	0 /	0 /	0 /
of	800	800	800	800	800	800	800	800	800	800
boluses										
% of	0.0%	0.0%	0.0%	0.0%	100.0%	0.0%	0.0%	0.0%	0.0%	0.0%
boluses										

				,	25U bolı	us accur	acy			
	<25%	25- 75%	75- 90%	90- 95%	95- 105%	105- 110%	110- 125%	125- 175%	175- 250%	>250%
Number of boluses	0 / 224	0 / 224	0 / 224	0 / 224	222 / 224	2 / 224	0 / 224	0 / 224	0 / 224	0 / 224
% of boluses	0.0%	0.0%	0.0%	0.0%	99.1%	0.9%	0.0%	0.0%	0.0%	0.0%

Non-Clinical/Performance Testing:

Performance testing was performed in order to establish substantial equivalence in terms of both safety and effectiveness, and to ensure the subject device met all applicable special controls. Performance testing was organized into the categories described below.

Nominal Basal Accuracy
Nominal Bolus Accuracy
Worst Case Accuracy
Occlusions
Fault Insertion
Sound Testing
Incidental Delivery
Reliability
Drug Compatibility and Particulate Testing
System Level Functionality
Battery Performance
Environmental Conditions

Testing was performed utilizing the following standards: IEC 60601-1, IEC 60601-1-2, IEC 60601-1-6, IEC 60601-1-8, IEC 60601-1-11, IEC 60601-2-24, IEC 62304.

Clinical Study

No clinical data was obtained in support of this premarket submission.

Design Control

The DEKA ACE Pump was specified and developed by DEKA. DEKA complies with the FDA Quality System Regulation as specified in 21 CFR 820, as well as to ISO 13485.

Conclusion

The DEKA ACE Pump System is substantially equivalent to the t:slim X2 insulin pump with interoperable technology. The differences summarized in this submission do not raise different questions of safety and effectiveness. The performance of the device is supported by DEKA's design control process which included non-clinical testing and risk management activities. The DEKA ACE Pump System complies with the ACE Pump Special Controls as established in DEN180058.