

June 12, 2023

Inofab Saglik Teknolojileri A.S % Ray Kelly Consultant Arazy Group Consultants Inc. 3422 Leonardo Lane New Smyrna Beach, Florida 32168

Re: K213754

Trade/Device Name: SpiroHome Ultrasonic Spirometer, SpiroHome Clinic, SpiroHome Personal Regulation Number: 21 CFR 868.1840 Regulation Name: Diagnostic Spirometer Regulatory Class: Class II Product Code: BZG Dated: January 17, 2023 Received: January 23, 2023

Dear Ray Kelly:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Rachana Visaria -S

Rachana Visaria, Ph.D.
Assistant Director
DHT1C: Division of Sleep Disordered
Breathing, Respiratory and
Anesthesia Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number *(if known)* K213754

#### Device Name

SpiroHome Ultrasonic Spirometer, SpiroHome Clinic, SpiroHome Personal

Indications for Use (Describe)

SpiroHome is intended to be used by adults and children over 5 years old in physician's offices, clinics and home setting to conduct basic lung function and spirometry testing.

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

### 510(k) Summary 21 CFR 807.92

In accordance with 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

### I. Submitter Information

Submitter Name:

	Inofab Health Sağlik Teknolojileri A.S Üniversiteler Mahallesi, İhsan Doğramacı Bulvarı, ODTÜ Teknokent Silikon Blok No:17/115, 06800 Çankaya/Ankara, Türkiye
Phone:	+90 312 988 03 08
Contact Person: Date Prepared:	Ray Kelly June 09, 2023

#### II. Device Information:

Name of Device:	SpiroHome Ultrasonic Spirometer,
	SpiroHome Clinic, SpiroHome Personal
Common Name:	Spirometer
Classification Name:	Diagnostic Spirometer (21 CFR 868.1840)
Regulatory Class:	II
Product Code:	BZG

### III. Predicate Device:

GoSpiro, K163249

### IV. Reference Device:

EasyOne Air Spirometer, K161536

### V. Device Description

The SpiroHome Ultrasonic Spirometer (SUS) is a portable spirometer designed to perform pulmonary function tests in patients over the age of 5 in office (clinical) and home settings. The SpiroHome spirometer is used together with a SpiroWay mouthpiece that is inserted into and lines the entire airway of the device. SpiroHome derives pulmonary function data from airflow measurements taken by its ultrasonic sensors during a spirometry test. All of the information recorded by the device is displayed on the relevant SpiroHome app running on a Bluetooth-connected device. The pulmonary function test (PFT) data recorded by the SpiroHome device during a spirometry test is also compared against the patient's predicted values which are

obtained from internationally accepted PFT equations. The user interfaces with the SpiroHome app during the entire use of the SpiroHome spirometer.

The associated accessories include: SpiroWay mouthpiece

### VI. Indication for Use:

SpiroHome is intended to be used by adults and children over 5 years old in physician's offices, clinics and home settings to conduct basic lung function and spirometry testing.

#### VII. Comparison of technological characteristics with the predicate and reference device:

The subject device has the same spirometry functions with the same test types, test parameters and utilize the same interpretation algorithms as the predicate device, however, the ultrasonic sensor technology used in the subject device is different to that of the predicate device which utilizes a turbine-based sensor. The subject device has the same technological characteristics with regards to spirometry testing as the reference device, in particular the same sensor technology is used for spirometry testing. The differences in sensor measurement technology do not raise concerns of safety and effectiveness for the subject device.

The subject device and predicate device both use Bluetooth-connected devices as displays in comparison to the reference device which has a touch-enabled display on the device. Each model of the subject device operates with its own respective application version whereas the predicate device has only one application with which it operates.

The subject device and predicate device both use off-the-shelf standard batteries whereas the reference device uses a rechargeable battery pack.

Each model of the subject device is used with its own respective mouthpiece model. The personal model of the subject device is used with a reusable mouthpiece and the clinic model is used with a disposable mouthpiece. The predicate device mouthpiece is constructed from both reusable and disposable components, and the reference device is used only with disposable mouthpieces. Design verification and validation demonstrates that the SpiroHome spirometer used in combination with the SpiroWay mouthpiece provides the same spirometry test results as the predicate and reference devices used in combination with their respective mouthpieces and is therefore substantially equivalent to the predicate and reference devices.

Table 1. Comparison of subject device attributes to predicate device and reference device

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
Indications for Use	Intended to be used by	Intended to be used by	The EasyOne Air	Identical for
	adults and children	adults and children	spirometer is intended	subject and
	over 5 years old	over 5 years old in	for prescription use only	predicate.
	in physician's offices,	physician's offices,	to conduct diagnostic	
	clinics and home	clinics and home	spirometry testing of	
	setting to conduct	settings to conduct	adults and pediatric	
	basic lung function and	basic lung function and	patients over 4 years	
	spirometry testing.	spirometry testing. It is	old. The EasyOne Air	
		a single-patient	spirometer is used by	
		use device.	general practitioners,	
			bospitals and clinics in	
			pharmacies and in	
			clinical settings in	
			occupational medicine.	
Common Name	Spirometer, Diagnostic	Spirometer, Diagnostic	Spirometer, Diagnostic	Identical
Regulation	868.1840	868.1840	868.1840	Identical
Product Code	BZG	BZG	BZG	Identical
Classification	II	II	11	Identical
Review Panel	Anesthesiology	Anesthesiology	Anesthesiology	Identical
Intended User	HCP or patient	HCP or patient	HCP	Identical for
				subject and
				predicate.
Environment of Use	Clinic or home	Clinic or home	Clinic	Identical for
				subject and
				predicate.
Principle of	Spirometer is designed	Spirometer is designed	Spirometer is designed	Identical for
Operation	for pulmonary function	for pulmonary function	for pulmonary function	subject and
	tests. It is a portable	tests. It is a portable	tests. It is a portable	predicate and
	spironneter capable of	spirometer capable of	spirometer capable of	similar to
	Eunction Test (DET)	Function Tost (PET)	Eurotion Tost (DET)	
	data based on the air	data based on the air	data based on the air	
	flow and volume. The	flow and volume. The	flow and volume. The	

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
	device measures a	device measures a	device measures a	
	patient's lung function	patient's lung function	patient's lung function	
	values and compares it	values and compares it	values and compares it	
	to their predicted	to their predicted	to their predicted values	
	values (obtained	values (obtained from	(obtained	
	from internationally	internationally accepted	from internationally	
	accepted lung function	lung function test	accepted lung function	
	test equations). The	equations). The	test equations). The	
	portable spirometer	portable spirometer	portable spirometer has	
	pairs to smart devices	pairs to smart devices	on-device screen	
	running the associated	running the associated	display.	
	apps.	apps.		
Conditions of Use	Spirometer is available	Spirometer is available	Spirometer is	Identical for
	for home use or clinical	for home use or clinical	available for only	subject and
	use. Clinic use is	use. Clinic use is	clinical use. Clinic use	predicate.
	indicated for healthcare	indicated for healthcare	is indicated	
	professionals (HCPs)	professionals (HCPs)	for healthcare	
	to use with patients	to use with patients	professionals (HCPs)	
	(pediatric between the	(pediatric between the	to use with patients	
	ages of 5-21 and adults	ages of 5-21 and adults	(pediatric between the	
	over the age of 21) who	over the age of 21) who	ages of 4-21 and	
	may have	may have been	adults over the age of	
	been diagnosed with a	diagnosed with a	21) who may have	
	chronic	chronic pulmonary	been diagnosed with	
	pulmonary disease.	disease. Clinic use is	a chronic pulmonary	
	Clinic use is intended	intended to be used in	disease. Clinic use is	
	to be used in a clinical	a clinical setting	intended to be used in	
	setting multiple times a	multiple times a day by	a clinical setting	
	day by multiple	multiple patients. The	multiple times a day	
	patients. The home	home user, after it is	by multiple patients.	
	user, after it is	prescribed by a		
	prescribed by a	clinical/professional,	Users need to interact	
	clinical/professional,	will be used by patients	with the device to	
	will be used by patients	at home. Home users	assemble the device by	
	at home. Home users	are (pediatric between	inserting batteries,	
	are (pediatric between	the ages of 5-21 and	turning on the device,	
	the ages of 5-21 and	adults over the age of	and go through the	
	adults over the age of	21) who may have	device set up	

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
	21) who may have	been diagnosed with a	procedure. Users then	
	been diagnosed with	chronic	perform a spirometry	
	a chronic pulmonary	pulmonary disease.	test by breathing into	
	disease. Patients or	Patients or their care	the device and	
	their care	providers/legal	receive feedback on	
	providers/legal	guardians will use	their performance	
	guardians will use the	the device multiple	though the app. Users	
	device multiple times	times a day by a	will clean the device	
	a day by a single	single patient. Users	as instructed.	
	patient. Users need	need to interact with		
	to interact with the	the device to assemble		
	device to assemble the	the device by		
	device by inserting	inserting batteries,		
	batteries, turning on	turning on the device,		
	the device, and go	and go through the		
	through the device set	device set up		
	up procedure in the	procedure in the app.		
	app. Users then	Users then perform		
	perform spirometry test	spirometry test		
	by breathing into the	by breathing into the		
	device and	device and		
	receive feedback on	receive feedback on		
	their performance	their performance		
	though the app. Users	though the app. Users		
	will clean the device	will clean the device		
	as instructed.	as instructed.		
Software	Mobile Medical App (MMA)	Mobile Medical App (MMA)	Mobile Medical App (MMA)	Identical
Software LOC	Moderate	Moderate	Moderate	Identical
Bluetooth	Yes	Yes	Yes	Identical
Used with PFT filter and / or mouthpiece	Mouthpiece Disposable	Mouthpiece Disposable	Mouthpiece Disposable	Similar for subject
	or Reusable	and Mouthport		and predicate.
		Reusable		Predicate
				mouthpiece has a
				reusable
				(mouthport) and
				disposable
				component and is

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
				for single- patient
				use. Subject
				device has
				reusable
				mouthpiece for
				home use and a
				disposable
				mouthpiece for
				clinic use, both are
				for single-patient
				use. Reference
				device has only a
				single-patient
				disposable
				mouthpiece. No
				concerns raised
				regarding subject
				device safety or
				effectiveness due
				to mouthpiece
				differences.
Configuration	Hand-held, portable	Hand-held, portable	Hand-held, portable	Identical
Patient	Over 5 years old	Over 5 years old	Over 4 years old	Identical for
Population				subject and
				predicate.
Sensor	Transit-time Ultrasound	Turbine	Transit-time Ultrasound	Identical for
Technology				subject and
				reference device.
Recalibration	Not required	Required	Not required	Identical for
				subject and
				reference.
Measured	FVC, FEV0.75 ,	FVC, FEV0.75 , FEV1 ,	*BEV, EOTV, FEF10,	Similar parameters
Parameters	FEV1 , FEV3 ,	FEV3 , FEV6	FEF25, FEF25-75,	for subject,
	FEV6 ,	, FEV0.75/FVC,	FEF25-75/FVC,	predicate and
	FEV0.75/FVC,	FEV1/FVC, FEV3/FVC	FEF40, FEF50,	reference devices.
	FEV1/FVC ,	, FEV6/FVC, PEF,	FEF50/FVC, FEF60,	Differences do
	FEV3/FVC ,	MMEF, FEF25 , FEF50	FEF75, FEF80, FET,	not affect the safe
	FEV6/FVC , PEF	,FEF75,FEF25-75,		and effective use

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome           , MMEF , FEF25 , FEF50 , FEF75           , FEF25-75 , MET25-75 , FEV0.75/FEV6           , FEV1/FEV6 , FEF50/FVC , MMEF/FVC , FET , BEV , FIV1 , FIVC , PIF , FIF25-75 , FIV1/FIVC , R50 (FEF50/FIF50), VC , VCin , VCex           , ERV , IRV , IC , Rf , VT , MVV , MVV6 , MVVtime	GoSpiro (K163249) FIV1 , FIVC , PIF, FIF25-75, FIF25 , FIF50 , FIF75, MET25- 75, FEV0.75/FEV6 , FEV1/FEV6 , FEF50/FVC , FIV1/FIVC, R50 (FEF50/FIF50) , FET, MVV	(K161536) FET25-75, FEV.25, FEV.5, FEV.5/FVC, FEV.75, FEV.75/FVC, FEV1, FEV1/FEV6, FEV1/FIV1, FEV1/FVC, FEV3, FEV3/FVC, FEV6, FIF25, FIF50, FIF50/FEF50, FIF75, FIV.25, FIV.5, FIV1, FVC6, MEF20, MEF25, MEF40, MEF50, MEF50, MEF60, MEF75, MEF90, MIF25, MIF50, MIF75, MMEF, PEF, PEFT, t0, MVV, MVV6, MVVtime, Rf, ERV, IC, IRV, Rf, VC, VCex, VCin, VCmax, VT	Differences of the subject device.
Power	2 x AAA Alkaline and Rechargeable Batteries	Rechargeable lithium- ion batteries	Rechargeable battery pack	The difference in power sources between devices does not raise concerns of safety or effectiveness of the subject device as validated by the electrical safety, performance, and EMC testing .
Air Resistance	48.54 Pa*s/L (Highest Expiratory Impedance)	137 Pa*s/L	*0.3 cm H2O/L/s at 16 L/s	Differences in dynamic air resistance between devices based on design

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
				differences and
				does not raise
				concerns over
				device safety and
				effectiveness of
				the subject device.
Volume Range & Accuracy	0-10 L,	0-8 L,	*±12 L,	Difference
	± 2.5% or ± 0.050 L	± 3% or ± 0.050 L	*±2% or 0.050 L	between volume
				range and
				accuracy
				between devices
				are minor, the
				subject device
				meets permissible
				margins given in
				ATS guidelines
				and ISO 286782.
Flow Range & Accuracy	0 - 14 L/s,	0 - 14 L/s,	*±16 L/s	Difference
	±10% or ± 170 mL/s	±5% or 200 mL/s	*±2% or 0.020 L/s	between flow
			(except PEF)	range and
			*±5% or 0.200 L/s PEF	accuracy
			accuracy	between devices
				are minor, the
				subject device
				conforms to
				permissible
				margins given in
				ISO 286782.
Display	Mobile app	Mobile app	Touch-enabled display	Identical for
Type/Size			on device	subject and
				predicate devices.
Flow-Volume Loop	During test and test	During test and test	*During test and test	Identical
	review	review	review	
Volume-Time Curve	During test and test	During test and test	*During test and test	Identical
	review	review	review	
Dimensions W x D x H	110mm x 40.8m x	88.9mm x 114.3mm x	*87mm x 155mm x	Subject device is
	63.3mm	50mm	36mm	smaller in overall
				volume than both

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
				predicate and
				reference devices.
				Difference in size
				does not raise
				concerns for
				subject device
				safety and
				effectiveness as
				demonstrated in
				usability and
				performance tests.
Weight	90g (with batteries)	300g (with batteries)	*356 g (with batteries)	Subject device
				weighs less than
				both predicate and
				reference devices.
				Difference in
				weight does not
				raise concerns for
				subject device
				safety and
				effectiveness as
				demonstrated in
				usability and
				performance tests.
Connection to patient	Mouthpiece	Mouthpiece	Mouthpiece	Identical
Storage	-20°C to 60°C	*-20°C to 70°C	*-20°C to 50°C	Minor differences.
Temperature				Differences do not
				affect the safety
				and effectiveness
				of the subject
				device as
				demonstrated in
	50( ) 050(	450( ) 050(	+====	product testing.
Storage	5% to 85%	15% to 95%	*5% to 90%	Minor differences.
Humidity (RH)				Differences do not
				attect the safety
				and effectiveness

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spironome	Gospiro (K163249)	(K161536)	Differences
				of the subject
				device as
				demonstrated in
Or section of		47%0 1. 05%0	*0°0 to 40°0	product testing.
	15°C to 35°C	17 °C to 35 °C	"0°C to 40°C	Minor differences.
Temperature				Differences do not
				anect the salety
				and enectiveness
				device as
				product testing
Operational	20% to 95%	20% to 75%	*5% to 00%	Minor differences
Humidity (PH)	30 % 10 83 %	30% 1075%	5 % 10 90 %	Differences do not
				affect the safety
				and effectiveness
				of the subject
				device as
				demonstrated in
				product testing
Expected	5 vears	Not Provided	*7 vears	Similar for subject
Service Life	- 5		,	and reference
				device, difference
				does not affect the
				safety or
				effectiveness of
				the subject
				device.
Water Ingress	IP22	IP22	Not Provided	Identical for
Protection				subject and
				predicate.
Electrical safety	Meets IEC 60601-1	Meets IEC 60601-1	Meets IEC 60601-1	Identical
EMC safety	Meets IEC 60601-1-2	Meets IEC 60601-1-2	Meets IEC 60601-1-2	Identical
Home use	Meets IEC 60601-1-11	Meets IEC 60601-1-11	Not Provided	Identical for
				subject and
				predicate device.
Biocompatibility	Meets ISO 10993-1,	Meets ISO 10993-1,	ISO 10993-1, ISO	Subject, predicate
	ISO 10993-5, ISO	ISO 10993-5,	10993-5, ISO 10993-10,	and reference

Attribute	Subject Device	Predicate Device	Reference Device	Similarities /
	Spirohome	GoSpiro (K163249)	(K161536)	Differences
	10993-10, ISO 10993- 11, ISO 18562- 1, ISO 18562-2, ISO 18562-3, USP 43-NF 38 (2020) <85> and EP 2.6.14 (EP 10.3)	ISO 10993-10, Volatile Organic Compounds (VOC), CO, CO2, Ozone, and PM2.5 testing	ISO 10993-18, ISO/FDIS 18562-2, ISO/FDIS 18562-3	devices meet requirements of ISO 10993-1 and ISO 18562-1. Subject device also meets pyrogen and LAL testing
				requirements.

## VII. PERFORMANCE DATA

The SpiroHome Ultrasonic Spirometer was tested against criteria for:

- Functional Requirements
  - ATS 2019 / ERS waveform simulator testing
  - o ISO 26782:2009
  - o ISO 23747:2015
  - High Altitude Performance
  - Flow Resistance
- Electrical Requirements
  - o AAMI ANSI ES 60601-1
  - o IEC 60601-1-11
  - IEC 60601-1-2
- Biocompatibility

Contact Type and Duration: Surface Contact, Mucosa, Externally Communicating, Tissue; permanent (>30 days)

- o ISO 10993-1
- o ISO 10993-3
- o ISO 10993-5
- o ISO 10993-10
- o ISO 10993-11
- o ISO/FDIS 18562-1
- o ISO/FDIS 18562-2
- o ISO/FDIS 18562-3
- Shipping Requirements, Packaging and Distribution
  - o ASTM D4332
  - ASTM D7386
  - o ASTM F1886/F1886M
  - o ISO 17664
  - o ISO 11737-1
- Cleaning
  - o ISO 17664
  - AAMI TIR30
  - AAMI TIR12
  - o ASTM E2314
  - o ISO 15883-1
  - o ISO 11737-1
  - o ISO/TS 15883-5
- Software and System Verification and Validation
  - Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices
  - IEC 62304
- Human Factors Study
  - o Guidance for Applying Human Factors and Usability Engineering to Medical Devices
  - o AAMI/ANSI HE75:2009

o ANSI/AAMI/IEC 62366-1:2015

#### VIII. CONCLUSIONS

Based upon the foregoing performance testing and comparison to the legally marketed predicate device, and reference device, for indications for use, technology, and performance, we believe we have demonstrated that the SpiroHome Spirometer is substantially equivalent to the predicate device.