## **QDOT MICRO™ Uni-Directional Navigation Catheter**

### **INSTRUCTIONS FOR USE**

DEVICE DESCRIPTION .....	4
INDICATIONS FOR USE .....	4
CONTRAINDICATIONS .....	5
WARNINGS AND PRECAUTIONS.....	5
ADVERSE EVENTS.....	8
SUMMARY OF CLINICAL STUDIES CONDUCTED FOR THE QDOT MICRO™ CATHETER .....	8
A. Objective .....	8
B. Study Design.....	9
C. Results.....	13
D. Study Conclusion.....	29
RF ABLATION.....	30
RF GENERATOR OPERATION .....	30
HOW SUPPLIED.....	30
PACKAGING.....	30
STORAGE.....	30
STERILIZATION AND USE-BY DATE.....	30
DISPOSAL .....	30
COMPATIBLE EP NAVIGATION SYSTEM .....	30
DIRECTIONS FOR USE .....	31
DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY .....	33

---

## Symbol Definitions

Following definitions are for reference only. Please refer to the product label for applicable usage.

**STERILE EO**

Sterilized Using Ethylene Oxide

Do Not Re-Use

Caution

Consult Instructions for Use

Do not use if package is damaged.

Do not use if package is opened.

Keep Away from Sunlight

Keep Dry

Use-By Date

**LOT**

Batch Code

**REF**

Catalog Number

Contents: 1

Manufacturer

Date of Manufacture

Pin Connector

**ELECTRODES EA**

Electrodes

**SPACING mm**

Spacing

Temperature Limit

---

**R<sub>x</sub>** only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Curve Type. Refer to label for colored circle containing applicable curve type.

**NAV**

Navigational Catheter

**C3**

Compatible with CARTO™ 3 EP Navigation System

## QDOT MICRO™ Uni-Directional Navigation Catheter

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize.
- Do not use if the package is open or damaged.

### DEVICE DESCRIPTION

The Biosense Webster QDOT MICRO™ Uni-Directional Catheter is a steerable multi-electrode luminal catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) energy to the catheter tip electrode for ablation purposes. The catheter shaft measures 7.5 F with 8 F ring electrodes. For ablation, the catheter is used in conjunction with a compatible RF generator and a dispersive pad (indifferent electrode). The catheter has force-sensing technology that provides a real-time measurement of contact force between the catheter tip and the heart wall.

The catheter has a high-torque shaft with a uni-directional deflectable tip section containing an array of electrodes which includes a 3.5 mm tip dome. All of the electrodes may be used for recording and stimulation purposes. The tip electrode serves to deliver RF energy from the generator to the desired ablation site. The tip electrode and ring electrodes are made from noble metals. The catheter incorporates six thermocouple temperature sensors and three ECG electrodes that are embedded in the 3.5 mm tip electrode. A thumb knob is used to deflect the tip. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Three curve configurations designated “D,” “F,” and “J” are available.

At the proximal end of the catheter, a saline input port with a standard Luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the tip electrode, to irrigate and cool the ablation site as well as the electrode tip. A compatible irrigation pump is used to control the saline irrigation. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

This catheter features a location sensor embedded in the tip section that transmits location and contact force information to the CARTO™ 3 Navigation System. An appropriate reference device is required for location reference position purposes. For information on using the catheter in mapping procedures and for information on appropriate reference devices, refer to the user manual for the CARTO™ 3 Navigation System. For further description of the operation of the irrigation pump, the generator, and the CARTO™ 3 Navigation System, refer to the applicable instructions for use.

### INDICATIONS FOR USE

The Biosense Webster QDOT MICRO™ Uni-Directional Navigation Catheter and related accessory devices are indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a compatible radiofrequency generator, for the treatment of:

- Type I atrial flutter in patients age 18 or older.
- Drug refractory recurrent symptomatic paroxysmal atrial fibrillation, when used with compatible three-dimensional electroanatomic mapping systems.

The Biosense Webster QDOT MICRO™ Uni-Directional Navigation Catheter provides a real-time measurement of contact force between the catheter tip and heart wall, as well as location information when used with CARTO™ 3 Navigation System.

## CONTRAINDICATIONS

Do not use this catheter:

1. If the patient has had a ventriculotomy or atriotomy within the preceding twelve weeks because the recent surgery may increase the risk of perforation.
2. In patients with a myxoma or an intracardiac thrombus as the catheter could precipitate an embolus.
3. In patients with prosthetic valves as the catheter may damage the prosthesis.
4. In the coronary arterial vasculature due to risk of damage to the coronary arterial vasculature.
5. In patients with an active systemic infection because this may increase the risk of cardiac infection.
6. Via the transseptal approach in a patient with an interatrial baffle or patch because the opening could persist and produce an iatrogenic atrial shunt.
7. Via the retrograde trans-aortic approach in patients who have had aortic valve replacement.
8. With a long sheath or short introducer < 8.5 F in order to avoid damage to the catheter shaft.

## WARNINGS AND PRECAUTIONS

1. Do not use excessive force to advance or withdraw the catheter when resistance is encountered during catheter manipulation through the sheath.
2. Do not manually pre-shape the distal shaft of the catheter by applying external forces intended to bend or affect the intended shape or curve of the catheter.
3. Prior to use, the catheter must be warmed up as specified in the Directions for Use section of this document. If the catheter has not reached a steady state condition, there is potential for a zero-offset drift to occur which could result in an inaccurate contact force reading.
4. Always zero the contact force reading following insertion into the patient or when moving the catheter from one chamber of the heart to another. Ensure the catheter is not in contact with heart tissue prior to zeroing. Refer to the user manual for your CARTO™ 3 System for instructions on how to zero the contact force reading.
5. The contact force reading might become inaccurate if the contact force sensor (located between the first and second ring electrode) comes into close proximity with a ferrous material, such as the braided shaft of another catheter. If extreme fluctuations in force are observed, ensure the catheter's contact force sensor is not in close proximity with another catheter's shaft, check zero on the catheter and, if necessary, remove and inspect the catheter. The contact force reading is for information only and is not intended to replace standard handling precautions.
6. To ensure proper operation of the contact force sensor, the tip electrode and the two distal ring electrodes must protrude from the distal tip of the guiding sheath.
7. When applying high lateral force during mapping and RF application, the user should monitor the contact force dashboard and vector display on the CARTO™ 3 System screen to ensure that contact force measurements remain within the accurate range. Refer to the Error Messages and Alerts section of the CARTO™ 3 System user manual for system-related alert messages and indications related to inaccurate force readings.
8. Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the tip section of the catheter does not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode. The temperature displayed on the generator is the temperature of the hottest of the six thermocouples in the electrode, not tissue temperature. Monitoring the temperature from the electrode during the application of RF energy ensures that the irrigation flow rate is being maintained.
9. At higher contact force settings, rapid power increases during ablation may result in steam pops, which increase the risk of perforation.
10. This catheter may damage the mechanical or prosthetic valve of a patient if the catheter is accidentally advanced through the valve.
11. The safety of discontinuing anticoagulation therapy following catheter ablation of atrial fibrillation has not been established. Anticoagulation therapy in such patients should be administered in accordance with the AHA/ACC/HRS 2014 Guideline for the Management of Patients With Atrial Fibrillation.

12. The safety and effectiveness of radiofrequency ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, substantial left atrial enlargement, and structural heart disease have not been established.
13. In accordance with your hospital's protocol, monitor the patient's fluid balance throughout the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible. Prior to the procedure, always identify the patient's risk of volume overload.
14. Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF energy. It is important to have temporary external sources of pacing and defibrillation available during ablation and to temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.
15. Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker. Patients who experience inadvertent complete AV block as a result of RF ablation may also require permanent pacing.
16. During the trans-aortic approach, adequate visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter, RF energy application, or both have been associated with myocardial infarction.
17. To prevent stenosis of the pulmonary veins, do not place the catheter in the pulmonary veins during the application of RF energy.
18. When ablating near adjacent anatomical structures, take precautions to minimize collateral damage to the adjacent structures.
19. When ablating near the phrenic nerve, take precautions to avoid injuring the phrenic nerve, including appropriately reducing RF power, pacing to identify the proximity of ablation electrode(s) to the nerve and/or fluoroscopic evaluation of the diaphragm post ablation.
20. When ablating near the esophagus (along the posterior wall of the left atrium), take precautions to avoid injuring the esophagus. These may include beginning the ablation with reduced RF power, reducing contact force, reducing application time, increasing the time interval between ablations, moving to new location if temperature rise observed, esophageal visualization, and/or intraluminal esophageal temperature monitoring.
21. Minimize or eliminate X-ray exposure during the procedure. Catheter ablation procedures using fluoroscopic imaging may present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of fluoroscopic imaging. Catheter ablation, when using fluoroscopic imaging, should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.
22. Do not expose the catheter to organic solvents such as alcohol.
23. Do not autoclave the catheter.
24. Do not immerse the proximal handle or cable connector in fluids. Electrical performance could be affected.
25. Do not scrub or twist the distal tip electrode during cleaning.
26. Inspect the saline within the irrigation tubing for air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause emboli.
27. Purge the catheter and irrigation tubing with heparinized normal saline prior to insertion of the catheter inside the patient body.
28. Electrophysiology catheters and systems are intended for use only in X-ray shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.
29. Do not attempt to operate the QDOT MICRO™ Uni-Directional Catheter or the generator prior to completely reading and understanding the applicable instructions for use.

30. Cardiac ablation procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory. Appropriate clinical instruction in the use of the QDOT MICRO™ Uni-Directional Catheter should also be completed.
31. The long-term risks of protracted fluoroscopy and creation of RF induced lesions have not been established. Furthermore, the risk/benefit in asymptomatic patients has not been studied.
32. To prevent thromboembolism, intravenous heparin (target ACT of  $\geq 350$  s) should be administered prior to or immediately following transeptal puncture during AF ablation procedures. The 2017 HRS/EHRA/ECAS/AOHRS/SOLAECCCE expert consensus statement on catheter and surgical ablation of atrial fibrillation recommends systemic anticoagulation with warfarin or a direct thrombin or factor Xa inhibitor for at least 2 months following an AF ablation procedure (Calkins H, Hindricks G, Cappato R et al. 2017 HRS/EHRA/ECAS/AOHRS/SOLAECCCE Expert Consensus Statement on Catheter and Surgical Ablation of Atrial Fibrillation. J Interv Card Electrophysiol. 2017 50:1-55).
33. When using the QDOT MICRO™ Uni-Directional Catheter with conventional systems (such as fluoroscopy or ultrasound imaging), or with the CARTO™ 3 Navigation System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be done under direct imaging guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. The firmness of the braided tip dictates that care must be taken to prevent perforation of the heart. The contact force reading is for information only and is not intended to replace standard handling precautions.
34. Always use the thumb knob to straighten the catheter tip before insertion or withdrawal of the catheter.
35. Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter.
36. When RF energy is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip should be inspected for char/coagulum that may be present on the tip. If present, do not continue the procedure with the same catheter and replace the catheter. If no char/coagulum is present, flush the tip to ensure the irrigation holes are not plugged prior to reinsertion of the catheter inside the patient body.
37. Apparent low power output, high impedance reading, or failure of the equipment to function correctly at normal settings may indicate faulty application of the indifferent electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication of the indifferent electrode or other electrical leads.
38. Read and follow the indifferent electrode manufacturer's instructions for use. The use of indifferent electrodes that meet or exceed ANSI/AAMI requirements (AAMI IEC 60601-2-2), is recommended.
39. The QDOT MICRO™ Uni-Directional Catheter is intended for use with a compatible RF generator, compatible irrigation pump, CARTO™ 3 Navigation System, Biosense Webster cables, and other appropriate interface cables and connectors. Use of a compatible irrigation pump is recommended to assure proper irrigation flow rate.
40. Care should be taken when ablating near structures such as the sino-atrial and atrioventricular nodes.
41. The sterile packaging and catheter should be inspected prior to use. Do not use if the packaging or catheter appears damaged.
42. The catheter is sterilized with ethylene oxide gas and should be used by the "Use By" date on the device package. Do not use the catheter if it is past the "Use By" date.
43. The QDOT MICRO™ Uni-Directional Catheter is intended for single patient use only.
44. Do not resterilize and reuse.
45. Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.
46. Use both direct imaging guidance (such as fluoroscopy or ultrasound) and electrogram data to monitor catheter advancement and reduce risk of tissue injury.
47. The QDOT MICRO™ Uni-Directional Catheter used in conjunction with a compatible RF generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and indifferent electrode, particularly when operating the catheter. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces.

48. The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
49. Electromagnetic interference (EMI) produced by the QDOT MICRO™ Uni-Directional Catheter, when used in conjunction with a compatible RF generator during normal operation, may adversely affect the performance of other equipment.
50. Electrodes and probes used for monitoring and stimulating devices can provide paths for high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes and probes as far away as possible from the ablation site and the indifferent electrode. Protective impedance may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.
51. If the generator does not display temperature, verify that the cables from the catheter to the CARTO™ System Patient Interface Unit (PIU) and the cable from the CARTO™ System Patient Interface Unit (PIU) to the generator are properly connected. If temperature still is not displayed, there may be a malfunction in the temperature sensing system that must be corrected prior to applying RF power.
52. The temperature measurement accuracy of the QDOT MICRO™ Uni-Directional Catheter, as with any temperature measurement electrophysiology catheter, is determined by the temperature accuracy of all of the connected devices. Please consult the appropriate user manuals for the connected devices for the temperature accuracy specification.
53. Before use, verify irrigation ports are patent by infusing heparinized normal saline through the catheter and tubing.
54. Regularly inspect and test reusable cables and accessories.

## **ADVERSE EVENTS**

---

**NOTE: The adverse events in the following summary were observed in clinical studies involving only the use of the QDOT Micro™ Catheter.**

---

### **Clinical trial for the QDOT Micro™ Catheter – Pivotal Study (Q-FFICIENCY)**

Of the 177 subjects in the Safety Analysis Set, 7 primary adverse events (AEs) were reported in 6 subjects. See the “Summary of Clinical Studies Conducted for the QDOT Micro™ Catheter” section below for a complete description of the AEs encountered during the study.

## **SUMMARY OF CLINICAL STUDIES CONDUCTED FOR THE QDOT MICRO™ CATHETER**

The following is a summary of the pivotal study (Q-FFICIENCY) performed under IDE for the QDOT Micro™ Catheter, using model numbers D-1394-XX-SI and D-1395-XX-SI. Data from this clinical study formed the basis for the US FDA PMA approval decision.

All patients underwent informed consent per the study protocol, and in compliance with the Code of Federal Regulations, 21 §50. The protocol and informed consent materials were reviewed and approved by the appropriate IRB prior to subject enrollment, and in compliance with the Code of Federal Regulations, 21 §56. Per the requirements of 21 CFR §812.20(b)(5), all participating investigators signed the Clinical Study Agreement prior to enrollment at their center, which included financial disclosure forms, current curricula vitae for the investigator and co-investigators, and written IRB approval for the investigational study. The clinical study was monitored in a manner consistent with 21CFR, Part 812, Subpart C, *Responsibilities of Sponsor*. Monitoring visits included, but were not limited to, verification of all study logs, verification that informed consent was being obtained in accordance with requirements described in the study protocol for all subjects participating in the study, verification of completeness of the Regulatory Binder, data source verification with the eCRFs, and identification and action to resolve any issues or problems with the study.

### **A. Objective**

The primary objective of this trial was to demonstrate the safety and effectiveness of the QDOT MICRO™ catheter when used with the nMARQ™ RF generator in the treatment of drug refractory symptomatic paroxysmal atrial fibrillation (PAF) during standard electrophysiology mapping and RF ablation procedures.

**B. Study Design**

The study was a prospective, single arm, unblinded, multicenter, pivotal clinical investigation conducted at 25 investigational sites in the US.

The study planned to enroll up to 185 subjects who are candidates for catheter ablation. Bayesian adaptive design was used to assess early success at one interim analysis performed 6 months after enrollment was completed and when 30% of subjects in the study completed 12 months of follow-up.

A CT/MRA imaging sub-study was integrated in the pivotal study. The first 40 consecutively enrolled subjects were included in the sub-study and underwent a 3-month CT/MRA in addition to the baseline CT/MRA to assess incidence of PV stenosis.

**B.1 - Study Objectives and Endpoints:**

The objectives and endpoints for this study are presented in Table 1. Failure modes for the primary effectiveness endpoint are listed in Table 2.

**Table 1 – Objectives and Endpoints for Q-EFFICIENCY Study**

Objectives	Endpoints
<b>Primary</b>	
<p>To demonstrate the safety and 12-month effectiveness of the QDOT MICRO™ catheter when used with the nMARQ™ RF generator for pulmonary vein isolation (PVI) in the treatment of subjects with drug refractory symptomatic paroxysmal atrial fibrillation (PAF).</p>	<ul style="list-style-type: none"> <li>• Primary safety endpoint: Incidence of any primary adverse event occurring within 7 days of the AF ablation procedure (initial or repeat procedure) using the QDOT MICRO™ catheter, except atrio-esophageal fistula and PV stenosis, which may also be considered as primary adverse events if occurring greater than 7 days and up to 90 days post the ablation procedure. Primary adverse events included the following: Death, myocardial infarction, PV stenosis, phrenic nerve injury/diaphragmatic paralysis, atrio-esophageal fistula, TIA, stroke/CVA, thromboembolism, pericarditis, cardiac tamponade/perforation, vagal nerve injury, major vascular access complications/bleeding, pulmonary edema (respiratory insufficiency), and heart block.</li> <li>• Primary effectiveness endpoint: Freedom from documented atrial fibrillation (AF), atrial tachycardia (AT) or atrial flutter (AFL) (hereinafter collectively referred to as “atrial tachyarrhythmias”) recurrence during the evaluation period (Day 91 through Day 365) and freedom from the following failure modes (defined in Table 2):                         <ul style="list-style-type: none"> <li>○ Acute procedure failure</li> <li>○ Repeat ablation failure</li> <li>○ Antiarrhythmic drug failure</li> </ul> </li> </ul>

Objectives	Endpoints
<b>Secondary</b>	
<p>To evaluate the incidence of (serious) adverse events during and after procedure up to 3 months following procedure</p>	<ul style="list-style-type: none"> <li>• Incidence of Unanticipated Adverse Device Effects (UADEs)</li> <li>• Incidence of Serious Adverse Events (SAEs) within 7 days (early onset), &gt; 7 to 30 days (peri-procedural) and &gt; 30 days (late onset) of initial ablation</li> <li>• Incidence of bleeding complication (ISTH definitions): a) major, b) clinically relevant non-major and c) minor bleeding</li> </ul>
<p>To evaluate Acute Procedural Success defined as confirmation of entrance block in all PVs</p>	<ul style="list-style-type: none"> <li>• Acute procedural success: Percent of subjects with electrical isolation of PVs (entrance block) at the end of the procedure                             <ul style="list-style-type: none"> <li>○ The percent of subjects with electrical isolation of PVs (entrance block) using QMODE+ as only ablation strategy</li> </ul> </li> <li>• Percent of subjects with electrical isolation of PVs (entrance block) after first encirclement (evaluated prior to the 20-minute waiting period and adenosine/isoproterenol challenge)</li> <li>• Percent of subjects with electrical isolation of all PVs (entrance block) after first encirclement without acute reconnection, after waiting period and adenosine/isoproterenol challenge</li> <li>• Percent of subjects and % of PVs with touch-up (i.e. touch-up is used to remove ablation of acute reconnection) among all targeted veins and touch-up location</li> <li>• Anatomical location of acute PV reconnection after first encirclement</li> </ul>
<p>To evaluate repeat ablation procedures during the 12-month period post-procedure</p>	<ul style="list-style-type: none"> <li>• Incidence (%) of repeat ablation procedures</li> <li>• Percent PVs re-isolated among all the targeted PVs at repeat procedure</li> <li>• Percent repeat ablation procedures requiring new linear lesions and/or identifying new foci outside of initially isolated area among the repeat ablation procedures</li> </ul>
<p>To evaluate the 12-month Single Procedure Success</p>	<ul style="list-style-type: none"> <li>• 12-month single procedure success: Freedom from documented AF/AFL/AT recurrence (episodes <math>\geq</math> 30 secs) during the Effectiveness Evaluation Period after a single ablation procedure and off anti-arrhythmia drugs (AADs) during the effectiveness evaluation period</li> </ul>

**Table 2 – Failure Modes for Primary Effectiveness Endpoint**

Failure Mode	Criteria
Acute procedural	<ul style="list-style-type: none"> <li>Failure to confirm entrance block in all pulmonary veins post-procedure.</li> <li>Use of a non-study catheter to treat left atrial ablation targets and cavo-tricuspid isthmus.</li> </ul>
Repeat ablation	<ul style="list-style-type: none"> <li>&gt; 2 repeat ablation procedures with the study catheter during the 3-Month blanking period (Day 0-90) after the index ablation procedure.</li> <li>Use of a non-study catheter to treat study arrhythmia ablation targets during the blanking period.</li> <li>Any repeat ablation procedure during the Evaluation Period.</li> </ul>
Antiarrhythmic drug	<ul style="list-style-type: none"> <li>Taking a new Class I or III AAD for AF or a previously failed Class I or III AAD at a greater than the highest ineffective historical dose for AF during the evaluation period.</li> </ul>

**B.2 - Safety and Effectiveness Performance Goals**

**B.2.1 - Safety Performance Goal**

The study was designed to compare the primary safety endpoint to a pre-determined safety performance goal of 14%. Data from recent PAF ablation studies for devices similar to the device in the current study were reviewed as a first step to deriving the performance goal for the primary safety endpoint. A meta-analysis approach was taken to estimate the average composite endpoint rate. Based on the meta-analysis, the upper bound of the 95% confidence interval was estimated to be equal to 9%. The pre-determined performance goal of 14% would reflect an approximately 50% increase in risk from the upper bound of the 95% CI.

**B.2.2 - Effectiveness Performance Goal**

The study was designed to compare the primary effectiveness endpoint to a pre-determined performance goal of 50%, which is indicated as the minimum acceptable success rate at 12 months for a paroxysmal AF population in the 2017 HRS/EHRA/ECAS/APHS/ SOLAECE Expert Consensus Statement on Catheter and Surgical Ablation of AF (Calkins, et al., 2017).

**B.3 - Subject Accountability:**

**Table 3 – Subject Accountability and Disposition**

Subject Enrollment	n
Enrolled Subjects	191
Excluded Subjects	14
Study Catheter Inserted	177
Discontinued Subjects	2
RF Energy Delivered with Study Catheter	175
Death	1
Withdrawn	7
Early Termination Subjects	0
Lost to Follow Up Subjects	5
Completed Subjects	162

Among the 191 enrolled subjects, 14 had no insertion of the study device and were excluded from the analysis sets for the following reasons:

- System connectivity related device deficiencies (n = 4)

- Not meeting study entrance criteria (n = 9)
- Occurrence of an adverse event (tamponade) before insertion of the study device (n = 1)

Two (2) subjects had the study catheter inserted but no RF applications were delivered due to system connectivity issues.

A total of 175 subjects received RF ablation using the study device. Among those, seven (7) subjects (4.0%) withdrew their consent and discontinued their study participation during follow-up from 48 to 384 days post ablation; five (5) subjects (2.9%) were lost to follow up between 104 days and 266 days post procedure; one (1) subject (0.6%) expired 7 months post procedure; and 162 subjects (92.6%) completed the 12-month follow-up visit.

### B.3.1 - Subject Disposition

The following definitions were used to classify subjects:

**Table 4 – Subject Classifications**

<b>Disposition</b>	<b>Definition</b>
Enrolled Subjects	Subjects who signed the informed consent.
Excluded Subjects	Subjects who were enrolled but never underwent insertion of the study catheter.
Discontinued Subjects	Subjects who had the investigational catheter inserted but did not undergo ablation (i.e., no RF energy was delivered via the study device). These discontinued subjects were followed up for 3 months.
Lost to Follow-up Subjects	Subjects who were enrolled, but contact was lost after the most recent follow-up visits (despite 3 documented attempts to contact the subject).
Withdrawn / Early Termination Subjects	Subjects who withdrew consent for study participation or were withdrawn by the investigator or were terminated from the study prior to completion of all follow-up visits.
Completed Subjects	Enrolled subjects who completed the 12-month follow-up visit and were not excluded, discontinued, withdrawn, early terminated, expired, or lost to follow-up from the study before the final study visit.

### B.3.2 - Analysis Sets

**Safety Analysis Set:** The Safety analysis set included all enrolled subjects who had the investigational device inserted, regardless of RF energy delivery.

**Modified Intent-To-Treat (mITT) Analysis Set:** The mITT analysis set consisted of all enrolled subjects who met all eligibility criteria and had the investigational device inserted.

**Per-Protocol (PP) Analysis Set:** The PP analysis set included subjects who satisfied the following criteria:

- Enrolled and met eligibility criteria
- Had undergone RF ablation
- Were treated with the study catheters and had been treated for the study related arrhythmia without major protocol deviations deemed to affect the scientific integrity of the study, included in the following list. No additional major protocol deviations beyond those below were identified that required exclusion of subjects from the Per Protocol analysis set
  - Missing all protocol-specified electronic effectiveness monitoring
  - Esophageal monitoring was not done per protocol and risk was not mitigated
  - The required waiting period of 20 minutes was not followed before pacing procedures and/or infusion of cardiac medications to induce AF/reconnection was not performed before the entrance block confirmation.

**Treatment Received Analysis Set:** The Treatment Received analysis set included all enrolled subjects who received RF ablation using the study catheter. This analysis set was created post-hoc because it was deemed relevant by the FDA.

**CT/MRA Analysis Set:** The CT/MRA set consisted of the first 40 subjects who are consecutively enrolled in the main arm of the study with 3-month CT/MRA assessment. The subject must have readable outcomes at baseline and 3 months.

**C. Results**

**C.1 - Subject Demographics:**

**Table 5 – Subject Demographics (Enrolled Subjects, N=191)**

Variables	All Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
Age (years) [1]				
n	191	177	167	166
Mean	63.5	63.3	63.1	63.2
SD	10.69	10.88	11.00	11.03
Min/Max	24 / 83	24 / 83	24 / 83	24 / 83
Sex, n/N (%)				
Male	116 / 191 (60.7)	106 / 177 (59.9)	101 / 167 (60.5)	101 / 166 (60.8)
Female	75 / 191 (39.3)	71 / 177 (40.1)	66 / 167 (39.5)	65 / 166 (39.2)
Race n (%)				
White, n/N (%)	166 / 191 (86.9)	157 / 177 (88.7)	149 / 167 (89.2)	148 / 166 (89.2)
Black or African American	7 / 191 (3.7)	6 / 177 (3.4)	6 / 167 (3.6)	6 / 166 (3.6)
Asian	1 / 191 (0.5)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
Race not reported	17 / 191 (8.9)	13 / 177 (7.3)	11 / 167 (6.6)	11 / 166 (6.6)
Ethnicity, n/N (%)				
Not Hispanic or Latino	164 / 191 (85.9)	154 / 177 (87.0)	146 / 167 (87.4)	145 / 166 (87.3)
Hispanic or Latino	8 / 191 (4.2)	6 / 177 (3.4)	6 / 167 (3.6)	6 / 166 (3.6)
Not reported	19 / 191 (9.9)	17 / 177 (9.6)	15 / 167 (9.0)	15 / 166 (9.0)

**C.2 - Baseline Characteristics:**

**Table 6 – Medical History (Enrolled Subjects, N=191)**

Medical History	Enrolled Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
<b>Cardiovascular</b>	156 / 191 (81.7)	143 / 177 (80.8)	135 / 167 (80.8)	134 / 166 (80.7)
Congestive heart failure	16 / 191 (8.4)	10 / 177 (5.6)	9 / 167 (5.4)	9 / 166 (5.4)

Medical History	Enrolled Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
Coronary artery disease	39 / 191 (20.4)	36 / 177 (20.3)	33 / 167 (19.8)	33 / 166 (19.9)
Vascular disease	25 / 191 (13.1)	23 / 177 (13.0)	22 / 167 (13.2)	22 / 166 (13.3)
Myocardial infarction	13 / 191 (6.8)	12 / 177 (6.8)	10 / 167 (6.0)	10 / 166 (6.0)
Hypertension	131 / 191 (68.6)	122 / 177 (68.9)	116 / 167 (69.5)	115 / 166 (69.3)
Cardiomyopathy	15 / 191 (7.9)	11 / 177 (6.2)	10 / 167 (6.0)	10 / 166 (6.0)
Left ventricular hypertrophy	1 / 191 (0.5)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
Significant valvular disease	3 / 191 (1.6)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
<b>Thromboembolic Events</b>	20 / 191 (10.5)	18 / 177 (10.2)	16 / 167 (9.6)	16 / 166 (9.6)
Transient ischemic attack	5 / 191 (2.6)	4 / 177 (2.3)	3 / 167 (1.8)	3 / 166 (1.8)
Stroke	6 / 191 (3.1)	5 / 177 (2.8)	5 / 167 (3.0)	5 / 166 (3.0)
Pulmonary embolus	2 / 191 (1.0)	2 / 177 (1.1)	2 / 167 (1.2)	2 / 166 (1.2)
<b>Diabetes</b>	37 / 191 (19.4)	36 / 177 (20.3)	33 / 167 (19.8)	33 / 166 (19.9)
<b>Obstructive Sleep Apnea</b>	50 / 191 (26.2)	46 / 177 (26.0)	44 / 167 (26.3)	44 / 166 (26.5)
CPAP Use	36 / 191 (18.8)	34 / 177 (19.2)	33 / 167 (19.8)	33 / 166 (19.9)
<b>CHA<sub>2</sub>DS<sub>2</sub>-VASc score</b>	2.4 ± 1.51	2.4 ± 1.53	2.4 ± 1.49	2.4 ± 1.50
<b>AF History (months)</b>	48.6 ± 73.63	49.6 ± 75.29	51.7 ± 76.92	51.9 ± 77.11
<b>PAF Episodes Within 12 Months</b>	24.9 ± 89.91	25.2 ± 92.80	26.0 ± 95.22	25.7 ± 95.38
<b>AAD Failed- Class I&amp;III</b>	1.2 ± 0.50	1.2 ± 0.50	1.2 ± 0.52	1.2 ± 0.52
<b>Arrhythmia other than Atrial Fibrillation</b>	87 / 191 (45.5)	81 / 177 (45.8)	73 / 167 (43.7)	73 / 166 (44.0)
Left/right Atrial Tachycardia (AT)	8 / 191 (4.2)	8 / 177 (4.5)	7 / 167 (4.2)	7 / 166 (4.2)
Supraventricular Tachycardia (SVT)	22 / 191 (11.5)	20 / 177 (11.3)	14 / 167 (8.4)	14 / 166 (8.4)
Atrial flutter (AFL)	43 / 191 (22.5)	41 / 177 (23.2)	40 / 167 (24.0)	40 / 166 (24.1)
AVNRT	0 / 191 (0.0)	0 / 177 (0.0)	0 / 167 (0.0)	0 / 166 (0.0)
Accessory pathway	0 / 191 (0.0)	0 / 177 (0.0)	0 / 167 (0.0)	0 / 166 (0.0)
Ventricular tachycardia	2 / 191 (1.0)	2 / 177 (1.1)	2 / 167 (1.2)	2 / 166 (1.2)
<b>LA Diameter (mm)</b>	38.06 ± 6.042	37.95 ± 5.917	38.09 ± 5.974	38.16 ± 5.917
<b>LVEF (%)</b>	59.4 ± 7.31	59.6 ± 7.02	59.7 ± 6.98	59.7 ± 6.99

Values in table represent n/N (%) or mean ± SD, as appropriate.

### C.3 - Index Ablation Procedure

#### **Ablation protocol**

A circumferential anatomical approach was used to isolate all PVs. Confirmation of entrance block in all PVs was required by the protocol with a 20-minute waiting period after the last PV encircling lesion, during isoproterenol infusion and/or after adenosine bolus.

The study ablation procedure used both QMODE and QMODE+ temperature control modes. The primary mode of ablation for PVI was QMODE+. QMODE was used for PVI once the investigator deemed

QMODE+ unable to complete PVI. QMODE was used primarily for RF application outside the PV ostia and for touch-up of the PVI, if necessary.

Table 7 and Table 8 present the ablation parameters for the QMODE+ mode and QMODE mode, respectively.

**Table 7 – QMODE+ RF and Flow Rate Settings during RF applications**

Power (W)	Target Temp (°C)		Cut-off Temp (°C)		Nominal Irrigation Flow rate (mL)
	Range	Maximum allowed	Range	Maximum allowed	
90†	40-60	60	60-70	70	8

† RF applications at this setting were limited to 4 sec. It was recommended to use lower target temperature setting for the posterior wall RF applications.

**Table 8 – QMODE Ablation Parameters\***

Power (W)	Target Temp (°C)		Cut-off Temp (°C)		Nominal Irrigation Flow rate (mL)
	Range	Recommended	Range	Recommended	
25-35	45-50	50	50-55	55	4
36-50	45-50	50	50-55	55	15

\* The study protocol recommended a maximum duration of 60 sec for RF applications using QMODE in general. For RF applications on the left atrial posterior wall, a maximum power of 35 W and a maximum duration of 30 sec were recommended.

Linear ablation lines were only required to treat documented macro-re-entry atrial tachycardias (ATs) and limited to the following targets only:

- LA roof line;
- Mitral Valve (MV) isthmus line;
- LA floor line;
- Cavo tricuspid isthmus (CTI).

A right atrial CTI linear ablation was required in cases with documented typical atrial flutter either before or during the procedure.

Prophylactic ablation of empirical sites was not allowed.

Ablation of spontaneous non-PV triggers was required.

Complex fractionated atrial electrogram (CFAE) ablation was not recommended.

All linear lesions required confirmation of bidirectional conduction block by pacing and/or mapping maneuvers.

The recommended contact force working range was 5-30 grams for both QMODE and QMODE+ modes.

**Procedural Data**

Among the 166 subjects in the PP analysis set, 59 subjects underwent the entire index procedure using QMODE+ mode only and 107 subjects using both QMODE+ & QMODE modes. QMODE+ only was used for PVI in 91 subjects in the PP cohort. Of the 91 subjects, 32 subjects received ablations outside PVs as well using QMODE.

Table 9-11 present the procedural data of the index procedures.

**Table 9 – Summary of Power, Temperature and Impedance Data Per Procedure (Per-Protocol Analysis Set, N=166)**

Description	QMODE+ Only n=59	QMODE+ and QMODE n=107
Average Power (W)	84.333 ± 0.4099	75.893 ± 7.5338
Average Temperature (°C)	39.415 ± 2.5587	39.555 ± 2.6122
Average Impedance Drop (ohms)	9.59 ± 1.696	9.47 ± 2.434

Values in table represent mean ± SD.

**Table 10 – Summary of Ablation Procedure Parameters (Per-Protocol Analysis Set, N=166)**

Variable	Per- Protocol Analysis Set, N=166	QMODE+ Only, N=59	QMODE + and QMODE Combined N=107
<b>Total mapping time (min)</b>			
n	166	59	107
Mean ± SD	10.49 ± 10.593	10.29 ± 7.575	10.61 ± 11.968
<b>Total procedure time (min)</b>			
n	166	59	107
Mean ± SD	144.26 ± 51.101	122.27 ± 28.463	156.38 ± 56.604
<b>Total fluoroscopy duration (min)</b>			
n	165	59	106
Mean ± SD	10.81 ± 10.600	8.65 ± 7.935	12.01 ± 11.690
<b>Total time to achieve PVI (min)</b>			
n	165	58	107
Mean ± SD	25.32 ± 33.948	16.38 ± 20.891	30.17 ± 38.483
<b>Total RF application time for ablating PVs (min)</b>			
n	165	58	107
Mean ± SD	11.66 ± 11.322	7.33 ± 5.273	14.01 ± 12.944
<b>Fluid delivered via the study catheter (mL)</b>			
n	154	52	102
Mean ± SD	546.5 ± 298.55	434.4 ± 190.79	603.7 ± 326.90
<b>Fluid delivered via intravenous fluids (mL)</b>			
n	159	56	103
Mean ± SD	861.2 ± 557.28	744.5 ± 472.39	924.7 ± 590.93

**Table 11 – Number of RF Applications, Average Duration of RF Application and Average Contact Force per Ablation Procedure (Per-Protocol Analysis Set, N=166)**

Variable	QMODE+ Only n=59	QMODE+ and QMODE Combined n=107
<b>Number of RF Applications</b>		
Mean ± SD	100.5 ± 39.36	132.7 ± 41.80
<b>Average Duration of RF Application (sec)*</b>		
Mean ± SD	3.870 ± 0.0393	7.384 ± 4.5099
<b>Average Contact Force (g)</b>		
Mean ± SD	15.800 ± 5.3346	15.194 ± 5.1492

\* Applications with RF energy >0 W and duration of RF application ≥1 sec were included for analysis.

### **Ablation Targets**

Table 12 summarizes ablation targets. As required by the protocol, PVs were targeted in all index procedures. Among the PP analysis set, about 40% of the subjects received ablations beyond PVI.

**Table 12 – Ablation Sites Targeted (Per-Protocol Analysis Set, N=166)**


Ablation Sites Targeted	Number of Subjects
<b>Total</b>	166
<b>PVI only</b>	100 / 166 (60.2%)
<b>Non-PV</b>	66 / 166 (39.8%)
Roof line	16 / 66 (24.2%)
Cavo-tricuspid isthmus	46 / 66 (69.7%)
Other linear lesion	15 / 66 (22.7%)
Other AF foci	16 / 66 (24.2%)
Other	3 / 66 (4.5%)

### **C.4 - Rhythm Monitoring Compliance**

During the 9-month evaluation period, rhythm monitoring included symptomatic and scheduled asymptomatic TTM transmissions (weekly from Month 3 visit through Month 5 post ablation and monthly starting from Month 6 to Month 12), 12-lead ECG at 3 and 12 months and at 6 months if subject returned for an in-clinic visit, and 24-hour Holter at Month 12.

Figure 1 shows the TTM, 12-lead ECG, and Holter compliance rates at each evaluation time point in the Per-Protocol Analysis Set. The overall TTM and Holter compliance rates were 70.3% and 58.5%, respectively. The ECG compliance rate was 93.3%, 83.7%, and 56.0% at Month 3, 6, and 12, respectively. The composite compliance rate was 97.4% and 91.2% at Month 6 and Month 12, respectively, where at least two rhythm monitoring methods were required per protocol.

**Figure 1 – TTM/ECG/HM Compliance by Time Point (Per-Protocol Analysis Set, N=166)**

**C.5 - Safety Results**

**C.5.1 - Primary Safety Endpoint**

*Primary Analysis*

The primary safety endpoint for this study was defined as the incidence of early onset Primary AEs for subjects undergoing a study ablation procedure. Per study protocol, the primary safety analysis was based on primary safety event rate using the mITT analysis set (N=166) as the primary analysis population.

Table 13 displays a summary for the primary safety endpoint in subjects with non-missing outcomes in the mITT analysis set. The raw incidence was 3.6%. The posterior mean of PAE rate was 4.2%, with a 95% Bayesian credible interval (BCI) ranging from 1.7% to 7.7%. This is well below the performance goal of 14%. The posterior probability that the PAE rate less than the performance goal of 0.14 was close to 1.0000 (i.e., > 0.99995), which exceeded the early success threshold of 0.975. The study met the early success threshold for the primary safety endpoint. There were no changes to the primary safety endpoint data between the early success interim analysis and the final analysis.

**Table 13 – Summary of Primary Safety Endpoint Analysis (mITT Analysis Set, N=167)**

Subjects with Primary Safety Events	6
Subjects without Primary Safety Events	160
Missing	1
Raw Incidence	3.6%
Posterior Mean*	4.2%
95% BCI	(1.7%, 7.7%)
Posterior Probability* that rate < 0.14	1.0000
Threshold to Pass	0.9750
Pass?	Yes

\* Using Beta (1,1) prior distribution

Table 14 presents the descriptive summaries by type of primary adverse events, which were adjudicated by CEC. The PAEs reported in this study were as follows: Three (3) major vascular access complication/bleeding events (reported in 3 subjects, 1.8%), two (2) cardiac tamponade/perforation events (2 subjects, 1.2%), and two (2) phrenic nerve injury/diaphragmatic paralysis events (1 subject, 0.6%). Both cardiac tamponade/perforation events required pericardiocentesis and one of them required blood transfusion; both were adjudicated as procedure related and possibly device related. The three (3) major vascular access complications were adjudicated as procedure related and included one femoral artery pseudoaneurysm requiring thrombin injection, one retroperitoneal bleed and one arteriovenous fistula. The events of phrenic nerve injury/diaphragmatic paralysis were adjudicated as both device and procedure related.

The subject adjudicated by CEC as phrenic nerve injury/diaphragmatic paralysis was re-evaluated at the 12-month visit to determine if the phrenic nerve injury is permanent and was deemed non-assessable by CEC. The sponsor has conservatively counted the adverse event as a PAE.

**Table 14 – Primary Safety Endpoint Outcome – Primary AEs (mITT Analysis Set, N=167)**

Variable	Number of Subjects with Event	Number of Events	Event Rate n/N (%) <sup>[1]</sup>
<b>Primary Adverse Event</b>	6	7	6 / 166 (3.6%)
<b>Type of Primary Adverse Event</b>			
Death (Device or procedure related)	0	0	0 / 166 (0.0%)
Atrio-Esophageal Fistula	0	0	0 / 166 (0.0%)
Cardiac Tamponade/Perforation	2	2	2 / 166 (1.2%)
Myocardial Infarction	0	0	0 / 166 (0.0%)
Stroke/Cardiovascular Accident (CVA)	0	0	0 / 166 (0.0%)
Thromboembolism	0	0	0 / 166 (0.0%)
Transient ischemic Attack (TIA)	0	0	0 / 166 (0.0%)
Phrenic Nerve Injury / Diaphragmatic Paralysis	1	2	1 / 166 (0.6%)
Heart Block	0	0	0 / 166 (0.0%)
Pulmonary Vein Stenosis	0	0	0 / 166 (0.0%)
Pulmonary Edema (Respiratory Insufficiency)	0	0	0 / 166 (0.0%)
Vagal Nerve Injury	0	0	0 / 166 (0.0%)
Pericarditis	0	0	0 / 166 (0.0%)
Major Vascular Access Complication / Bleeding	3	3*	3 / 166 (1.8%)

<sup>[1]</sup> One subject did not experience a PAE and exited the study prior to 3-month follow-up visit; subject was excluded from this analysis.

\* Major vascular access complication/bleeding events included one retroperitoneal hemorrhage, one vascular access site hemorrhage and one pseudoaneurysm.

### *Sensitivity Analysis*

The primary safety endpoint was also analyzed in the Safety Analysis set because it is considered by the FDA as a more clinically relevant analysis population for device safety evaluation. Of the 176 subjects in the Safety analysis set, 10 underwent the study ablation procedure and had the study catheter inserted but were identified post-procedure to not meet eligibility criteria for the study. These subjects were not included in the mITT analysis set but were included in the Safety analysis set. No PAEs were reported in these 10 subjects. One (1) subject withdrew from the study 48 days post the index procedure with no PAE, was considered missing and excluded from this analysis.

As summarized in Table 15, the rate of PAE was 3.4% (6/176) in the safety analysis set. The posterior mean of the PAE rate was 3.9%, with an 95% BCI ranging from 1.6% to 7.2%. With 10 more subjects in the analysis, the posterior probability of the PAE rate being less than 0.14 remain greater than 0.9999 which is well above the threshold of 0.975.

**Table 15 – Summary of Primary Safety Endpoint Analysis (Safety Analysis Set, N=177)**

Subjects with Primary Safety Events	6
Subjects without Primary Safety Events	170
Missing	1
Raw Incidence	3.4%
Posterior Mean*	3.9%
95% BCI*	(1.6%, 7.2%)
Posterior Probability* that rate < 0.14	1.0000
Threshold to Pass	0.9750

\* Using Beta (1,1) prior distribution

C.5.2 - Incidence of UADEs

There were no UADEs reported in the study.

C.5.3 - Incidence of Early Onset, Peri-Procedural and Late Onset SAEs

Table 16 Summaries the incidence of early onset, peri-procedural and late onset SAEs in the Safety Analysis set.

**Table 16 – Serious Adverse Events by Timing of Onset (Safety Analysis Set, N=177)**

Serious Adverse Events - by Timing of Onset <sup>[1]</sup>	Number of Events	Number of Subjects with Event	Event Rate <sup>[2]</sup> n/N (%)
<b>Total</b>	55	34	
0 - 7 Days	14	12	12 / 177 (6.8%)
8 - 30 Days	5	4	4 / 177 (2.3%)
> = 31 Days	36	25	25 / 177 (14.1%)

Note: This table summarizes all SAEs. Adjudication for seriousness is based on CEC for PAEs or site for non-PAEs.

<sup>[1]</sup> Including events on and after the day of procedure. Timing of AE based on onset in days post procedure (index or repeat).

<sup>[2]</sup> N denotes number of subjects with follow-up duration corresponding to each time interval.

Table 17 summarizes the SAEs (by causality and body system) occurring within 30 days of a study ablation procedure that were not classified as Primary AEs by protocol definition.

**Table 17 – Serious Non-Primary AEs Occurring within 30 Days Post Index/Repeat Ablation Procedure) by Causality and Body System (Safety Analysis Set, N=177)**

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Overall</b>	11	13	11 / 177 (6.2%)
<b>Causal Relationship to Device</b>	0	0	0 / 177 (0.0%)

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Probable Device Related</b>	0	0	0 / 177 (0.0%)
<b>Possible Device Related</b>	0	0	0 / 177 (0.0%)
<b>Unlikely Device Related</b>	0	0	0 / 177 (0.0%)
<b>Not Device Related</b>	11	13	11 / 177 (6.2%)
Cardiac disorders	7	8	7 / 177 (4.0%)
Atrial fibrillation	3	3	3 / 177 (1.7%)
Atrial flutter	2	2	2 / 177 (1.1%)
Congestive Cardiac failure	1	1	1 / 177 (0.6%)
Cardiac tamponade	1	1	1 / 177 (0.6%)
Pericardial effusion	1	1	1 / 177 (0.6%)
General disorders and administration site conditions	3	3	3 / 177 (1.7%)
Chest pain	2	2	2 / 177 (1.1%)
Pyrexia	1	1	1 / 177 (0.6%)
Respiratory, thoracic and mediastinal disorders	1	1	1 / 177 (0.6%)
Pulmonary embolism	1	1	1 / 177 (0.6%)
Vascular disorders	1	1	1 / 177 (0.6%)
Arteriovenous fistula	1	1	1 / 177 (0.6%)
<b>Causal Relationship to Procedure</b>	1	1	1 / 177 (0.6%)
General disorders and administration site conditions	1	1	1 / 177 (0.6%)
Chest pain	1	1	1 / 177 (0.6%)
<b>Probable Procedure Related</b>	3	5	3 / 177 (1.7%)
Cardiac disorders	2	3	2 / 177 (1.1%)
Congestive Cardiac failure	1	1	1 / 177 (0.6%)
Cardiac tamponade	1	1	1 / 177 (0.6%)
Pericardial effusion	1	1	1 / 177 (0.6%)
Respiratory, thoracic and mediastinal disorders	1	1	1 / 177 (0.6%)
Pulmonary embolism	1	1	1 / 177 (0.6%)
Vascular disorders	1	1	1 / 177 (0.6%)
Arteriovenous fistula	1	1	1 / 177 (0.6%)
<b>Possible Procedure Related</b>	0	0	0 / 177 (0.0%)
<b>Unlikely Procedure Related</b>	0	0	0 / 177 (0.0%)

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Not Procedure Related</b>	7	7	7 / 177 (4.0%)
Cardiac disorders	5	5	5 / 177 (2.8%)
Atrial fibrillation	3	3	3 / 177 (1.7%)
Atrial flutter	2	2	2 / 177 (1.1%)
General disorders and administration site conditions	2	2	2 / 177 (1.1%)
Chest pain	1	1	1 / 177 (0.6%)
Pyrexia	1	1	1 / 177 (0.6%)

<sup>[1]</sup> Timing of AE based on onset in days post the procedure (index or repeat).

<sup>[2]</sup> Event rate is the percentage of subjects with the event.

#### C.5.4 - Subject Death

There was one (1) death during the study. The subject was a 66-year-old male with history of PAF and multiple comorbidities (coronary artery disease, hypertension, non-ischemic cardiomyopathy, congestive heart failure, status post-biventricular pacer for cardiac resynchronization therapy). He underwent the study procedure without immediate complications. Five months following the ablation procedure, the subject suffered a witnessed cardiac arrest due to ventricular fibrillation and was resuscitated by immediate cardiopulmonary resuscitation and external defibrillations 10 minutes after the onset of ventricular fibrillation. He subsequently suffered from anoxic brain injury and pneumonia without clinical improvement after therapy. The subject expired about 7 months post ablation. The death was classified as not related to the procedure or the device by the investigator.

#### C.5.5 - CT/MRA Sub-Set Analysis

Table 18 presents the results of the CT/MRA subset analysis. A total of 40 subjects in the CT/MRA subset underwent the CT/MRA imaging of the PVs pre-procedure and 3 months post-procedure. No subjects showed moderate or severe pulmonary vein stenosis in the subset.

**Table 18 – PV Narrowing Post Ablation (CT/MRA Analysis Set, N=40)**

Pulmonary Vein Stenosis	Mild		Moderate	Severe	
	Less than 0%	0-20%	20-50%	50-70%	> 70%
Subject-Level Summary, n/N (%)	1 / 40 (2.5%)	15 / 40 (37.5%)	24 / 40 (60.0%)	0 / 40 (0.0%)	0 / 40 (0.0%)

### C.6 - Effectiveness Results

#### C.6.1 - Primary Effectiveness Endpoint

##### *Primary Analysis*

Per study protocol, the primary effectiveness analysis was based on primary effectiveness success using the PP Analysis Set as the primary analysis population.

An early success interim analysis was performed as pre-specified in the protocol with a data extract obtained on June 22nd, 2020 and repeated after the endpoint was met with an additional data extract obtained on August 14th, 2020 to incorporate FDA feedback.

At the planned interim analysis, the posterior mean of event-free rate at Month 12 was 74.5%, with a 95% Bayesian credible interval (BCI) ranging from 66.6% to 81.5%. The posterior probability that the

event-free rate at Month 12 was greater than the performance goal of 50% was close to 1.0000 (i.e., 0.99995), which exceeded the pre-specified threshold of 0.9975. The primary effectiveness endpoint was therefore met at the early success interim analysis.

The analysis of the primary effectiveness endpoint on the August 14, 2020 data extract yielded similar results. These results have been updated with complete follow-up on the full cohort using a beta-binomial model with a non-informative prior. The final primary effectiveness rate is nearly identical to the interim estimate.

Table 19 presents the primary effectiveness analysis based on interim analysis and final analysis in the PP Analysis set.

**Table 19 – Summary of Primary Effectiveness Endpoint Analysis based on the Interim Analysis and Final Analysis (Per-Protocol Analysis Set, N=166)**

	Jun-22, 2020 Data	Aug-14, 2020 Data	Final Data
	Based on 3-piece exponential model	Based on 3-piece exponential model	Based on Beta-Binomial Model
<b>Kaplan-Meier Estimate (39 weeks)</b>	0.765	0.778	0.767
<b>Posterior Mean</b>	0.745	0.755	0.759
<b>95% BCI</b>	(0.666, 0.815)	(0.684, 0.820)	(0.690, 0.823)
<b>Posterior Probability that <math>P_E &gt; 0.50</math></b>	1.0000	1.0000	1.0000
<b>Threshold to Pass</b>	0.9975	0.9975	
<b>Pass</b>	Yes	Yes	

The first failure reason for subjects in the Per-Protocol analysis set is presented in Table 20. Ten (10) subjects with missing effectiveness outcomes are excluded from the summary. Approximately seventy six percent (76.3%, 119/156) of the Per-Protocol population were free from documented atrial tachyarrhythmias and additional failure modes during their effectiveness evaluation period. Documented recurrence of atrial tachyarrhythmias (17.3%, 27/156) and repeat ablation failures (2.6%, 4/156) were the most common reasons for first failure.

**Table 20 – Primary Effectiveness Endpoint by Reason of The First Failure (Per-Protocol Analysis Set, N=166)**

Variable	Number of Subjects with Event	Event Rate n/N (%) <sup>[1]</sup>
Success	119	119/156 (76.3%)
First Failures <sup>[2]</sup>	37	37/156 (23.7%)
Atrial Tachyarrhythmia Recurrence	27	27/156 (17.3%)
AF	25	25/156 (16.0%)
AFL	2	2/156 (1.3%)
AT	0	0/156 (0.0%)
Acute Procedural Failure	3	3/156 (1.9%)
Acute Procedural Failure - Non-study Catheter	3	3/156 (1.9%)
Acute Procedural Failure - Entrance Block	0	0/156 (0.0%)
Repeat Ablation	4	4/156 (2.6%)
AAD Failure	3	3/156 (1.9%)
AAD Failure - New Drug	3	3/156 (1.9%)

Variable	Number of Subjects with Event	Event Rate n/N (%) <sup>[1]</sup>
AAD Failure - High Dose	0	0/156 (0.0%)
Missing	10	

<sup>[1]</sup> N is the number of subjects with non-missing primary effectiveness endpoint.

<sup>[2]</sup> First failures: If a subject has at least 1 failure event, only the earliest failure event is considered and ed by its failure mode.

*Sensitivity Analysis*

Sensitivity to Analysis Sets

Table 21 presents the posterior means of event-free rate at Month 12 using the beta binomial model in the mITT, Safety and Treatment Received Analysis Sets based on the final full follow up data. The Treatment Received Analysis Set was included in the sensitivity analyses per FDA’s feedback because it is considered by the FDA as the most clinically relevant analysis population for device effectiveness evaluation. The results suggested that the success rates in all three Analysis Sets met the effectiveness performance goal. Of note, of the two (2) discontinued subjects who had the study catheter inserted but received no RF ablation from the device, one was not included in the mITT Analysis Set due to study eligibility criteria violation, both were included in the Safety Analysis Set and neither was included in the Treatment Received Analysis Set. Both discontinued subjects were classified as primary effectiveness failures due to use of a non-study catheter for PVI in the respective analysis population.


**Table 21 – Primary Effectiveness Endpoint in mITT, Safety and Treatment Received Analysis Sets**

	mITT Analysis Set (N=167)	Safety Analysis Set (n = 177)	Treatment Received Analysis Set (n = 175)
<b>Posterior Mean</b>	0.755	0.746	0.754
<b>95% BCI</b>	(0.685, 0.818)	(0.678, 0.808)	(0.687, 0.817)


Kaplan-Meier analysis

Figures Figure 2 and Figure 3 present Kaplan-Meier freedom from effectiveness failure for subjects in the Per-Protocol Analysis Set and Treatment Received Analysis Set, respectively. The 12-month Kaplan-Meier estimate of freedom from effectiveness failure was 76.7% (97.5% LCB: 70.1%) in the Per-Protocol Analysis Set and 76.2% (97.5% LCB: 69.7%) in the Treatment Received Analysis Set, both of which are comparable to the posterior mean of 75.9% from the final data.

**Figure 2 – Kaplan-Meier Analysis of Time to First Primary Effectiveness Failures Post Procedure (Per-Protocol Analysis Set, N=166)**

**Figure 3 – Kaplan-Meier Analysis of Time to First Primary Effectiveness Failures Post Procedure (Treatment Received Analysis Set, N = 175)**

Worst-case scenario and tipping point analyses


In the worst-case scenario analysis, the 10 subjects missing primary effectiveness endpoint data were considered as having primary effectiveness failure, which resulted in the primary effectiveness success rate of 71.7% in the Per-Protocol Analysis Set, 71.3% in the mITT analysis set and 71.4% in Treatment Received Analysis Set.

Treating all of the 10 subjects with missing outcomes as effectiveness failures one by one, tipping point analyses using a posterior probability threshold of 0.9775 were performed for the primary effectiveness endpoint in the PP, mITT and Treatment Received Analysis Sets to assess the impact of missing outcomes on the primary effectiveness analysis. The posterior distribution was updated each time and the results showed that no tipping point was identified. Therefore, the primary effectiveness performance goal was met in the tipping point analysis.

*Impact of TTM Compliance*

To better understand the impact of TTM compliance on the primary effectiveness results, additional post-hoc Kaplan-Meier analysis on the time to first primary effectiveness failure was performed in subgroups dichotomized at TTM compliance of 84%, a compliance rate reported in a previous PAF ablation trial. The probability of freedom from primary effectiveness failures at 12 months was similar between subjects with TTM compliance rate of < 84% and those ≥ 84% (Figure 4).

**Figure 4 – Kaplan-Meier Analysis of Time to Primary Effectiveness Failures Post Procedure by TTM Compliance (≥ 84% vs. <84%) (Per-Protocol Analysis Set, N=166)**

C.6.2 - Acute Procedural Success

Acute procedural success was defined as confirmation of entrance block in all PVs. Off the 166 subjects included in the PP analysis set, two (2) subjects were treated with a non-study catheter for PVI due to CoolFlow Pump malfunction which resulted in replacement of the QDOT catheter with a non-study catheter to complete PVI. These two (2) subjects were excluded from the analysis of acute

procedural success. Acute procedural success was achieved in the remaining 164 subjects (100%, 97.5% LCB: 97.8%) in the PP Analysis Set.

Table 22 presents data on PVI after first encirclement at subject level and at ipsilateral PV level, respectively. Among a total of 164 subjects in the PP Analysis Set who received PVI ablation using the study catheter only (two subjects also treated with a non-study catheter for PVI and discussed above were excluded from the analysis), 94 (57.3%) achieved PVI with first encirclement without acute reconnections. For this analysis, both right-sided PVs and left-sided PVs for the subject had to achieve PVI without acute reconnections. At ipsilateral PV level, PVI was achieved in 72.6% (238/328) of the pulmonary veins (ipsilateral veins) after first encirclement without acute reconnections.

**Table 22 – PVI after First Encirclement at Subject Level and at Ipsilateral PV level (Per-Protocol Subjects Who Received PVI Ablation Using the Study Catheter Only, N=164)**

Variable	Number of Subjects Underwent Ablation Procedure	Number of Subjects	Rate
PVI after first encirclement	164	94	94 / 164 (57.3%)
Variable	Number of Targeted Ipsilateral Veins	Number of Targeted Ipsilateral Veins with first encirclement	Rate
PVI after first encirclement	328	238	238 / 328 (72.6%)

C.6.3. - Incidence of repeat ablation procedures


A total of 18 subjects (10.8%) in the Per-Protocol analysis set had one (1) repeat procedure and most (66.7%, 12/18) had their repeat procedure during the evaluation period (Day 91-365). Most subjects (88.9%, 16/18) were treated for AF during the repeat procedure.

C.6.4. - 12-month single procedure success

The 12-Month single procedure success was defined as freedom from documented AF/AFL/AT recurrence (episodes ≥ 30 secs) during the evaluation period after a single ablation and off AADs. Any repeat ablation procedure or Class I or III AAD use for AF/AFL/AT during the evaluation period was deemed effectiveness failure for this endpoint analysis.

Figure 5 presents the Kaplan- Meier analysis of 12-Month single procedure success off Class I/III AAD in the Per Protocol analysis set. The probability of 12-Month single procedure success off Class I/III AADs post blanking was 64.3% in the Per-Protocol analysis set.

**Figure 5 – 12-Month Single Procedure Success off Class I/III AADs Post Blanking (Per Protocol Analysis Set, N=166)**

**C.6.5 - Quality of Life**

The quality of life was measured using the Atrial Fibrillation Effect on Quality-of-Life (AFEQT) at baseline and the 3-, 6-, and 12-month visits.

Figure 6 presents the mean AFEQT composite scores at baseline, 3-, 6-, and 12-month follow up visits in the effectiveness evaluation period, indicating a sustained improvement on the AFEQT scores from 3 months to 12 months post ablation.

**Figure 6 – Mean AFEQT Composite Scores Throughout Study Period (Per-protocol Analysis Set, N=166)**



Figure 7 presents change in total AFEQT scores from baseline at 3, 6, and 12 months. The mean change of AFEQT score from baseline to 3 months (28.77) was sustained at 12 months (31.48) and was larger than a clinical important difference of 19 (Dorian et al., 2013) in the majority of subjects.

**Figure 7 – Atrial Fibrillation Effect on Quality of Life (AFEQT) Questionnaire Overall AFEQT Score (Per-Protocol Analysis Set, N=166)**

**C.7 - Subgroup Analyses**

Subgroup analyses for the primary endpoints were performed, including but not limited to: sex, age, CHA2DS2-VASc score, and ablation mode used for PVI. Per protocol, a test of independence by sex was performed for both primary endpoints, while for other subgroup analyses only descriptive statistics were presented. The PP Analysis Set and Safety Analysis Set were used for the subgroup analyses of the primary effectiveness endpoint and primary safety endpoint, respectively. The primary effectiveness success rates were similar by age (< 65 vs. ≥ 65 years) and CHA2DS2-VASc score (<=2 vs. >2).

i Gender analysis

The primary endpoints were assessed according to gender. Tables Table 23 and Table 24 summarize the primary effectiveness and safety endpoints by sex in the PP and Safety Analysis Sets, respectively. Effectiveness and safety outcomes were similar between male and female subjects.

**Table 23 – Primary Effectiveness Endpoint by Sex (Per-Protocol Analysis Set, N=166)**

Primary Effectiveness Endpoint	Male	Female	p-value <sup>[1]</sup>
Primary Effectiveness Endpoint Success			0.329
n/N	75 / 95	44 / 61	
%	78.9%	72.1%	

<sup>[1]</sup> Chi-square test

**Table 24 – Primary Safety Endpoint by Sex (Safety Analysis Set, N=177)**

Primary Safety Endpoint	Male	Female	p-value <sup>[1]</sup>
Primary AE			0.683
n/N	3 / 106	3 / 70	
%	2.8%	4.3%	

<sup>[1]</sup> Fisher's exact test

ii Primary Safety by Age (<65 vs. ≥65)

Table 25 summarizes the primary safety endpoint by age in the Safety Analysis Set. No PAEs were observed in subjects aged <65 years. All PAEs occurred in subjects who were aged ≥65 years. The big difference in the primary safety event rate between the two age groups was likely caused by chance in the study in which a very small number of primary safety events occurred. Previous AF ablation studies (Hao, et al., 2012 and Zado, et al., 2008) reported comparable major complication rates between patients 65 years of age or older and those under age 65.

**Table 25 – Primary Safety Endpoint by Age (Safety Analysis Set, N=177)**

Primary Safety Endpoint	Age < 65 years	Age ≥ 65 years
n/N	0 / 82	6 / 94
%	0.0%	6.4%

iii Subgroup analysis of primary endpoints by ablation mode used for PV isolation

Primary effectiveness success was evaluated by ablation mode used for isolating the PVs. Figure 8 presents the Kaplan- Meier analysis of time to primary effectiveness failures by ablation mode used for PVI in the PP Analysis Set. The probabilities of freedom from primary effectiveness failure at 12 months post-procedure were comparable in the 2 groups, QMODE+ for PVI (78.1%) vs. QMODE+ and QMODE for PVI (75.1%).

**Figure 8 – Kaplan-Meier Analysis of Time to First AF/AT/AFL Recurrence Through 12 Months Post Procedure by PV Ablation Mode (Per-Protocol Analysis Set, N=166)**


Table 26 summarizes the primary safety outcome by ablation mode used for PVI in the Safety Analysis Set. The PAE rate was similar between the QMODE+ only group and QMODE+ and QMODE group.

**Table 26 – Primary Safety Endpoint by PV Ablation Mode (Safety Analysis Set, N=177)**

Primary Safety Endpoint	QMODE+ Only	QMODE+ and QMODE
Primary AE		
n/N	3 / 95	3 / 79
% [1]	3.2%	3.8%

[1] The percentage of subjects with PAEs within each group. Subjects without at least 3 months of follow-up are excluded from the primary safety analysis unless they experienced a PAE.

**C.8 - Results in Subjects Who Received Cavo-tricuspid Isthmus Ablation Using Study Catheter**

Among a total of 175 subjects who received RF ablation using the study device in the pivotal trial, 47 underwent CTI ablation using the study device only (n = 45), a non-study catheter only (n = 1), or both (n = 1) during the index procedure per study protocol for typical atrial flutter documented either before or during the procedure. Protocol recommended QMODE mode was used to ablate the CTI in all 33 subjects with ablation mode data available and one of the 33 subjects also received ablations in the CTI using the QMODE+ mode.

Among the 46 subjects who received CTI ablation using the study catheter, 44 (95.7%) had bidirectional CTI block confirmed, and the remaining two subjects either did not achieve bidirectional CTI block (n=1) or had missing information about CTI block (n = 1).

Among the 46 subjects who received CTI ablation using the study catheter, none had a primary AE, three (3) had a total of three (3) serious adverse events and 11 had a total of 12 non-serious adverse events within 30 days after a study procedure. None of these serious adverse events or non-serious adverse events was adjudicated as study device-related; while 9 non-serious adverse events were adjudicated as ablation procedure related.

Among the 46 subjects who received CTI ablation using the study catheter, two (2) had documented atrial flutter recurrence during follow-up. The type of recurrent atrial flutter was atypical atrial flutter in one subject and undetermined in another because it was only documented on TTM.

**D. Study Conclusion**

The results of the Q-EFFICIENCY study demonstrated that there is a reasonable assurance of safety and effectiveness of the QDOT MICRO™ Catheter when used for the treatment of symptomatic drug refractory paroxysmal AF and typical atrial flutter.

## RF ABLATION

For RF ablation the catheter is connected to the CARTO™ 3 System Patient Interface Unit (PIU) via accessory cables, which connect to the RF generator. For setup procedures refer to the User Manual for your CARTO™ 3 Navigation System. For proper RF generator interface, use only a Biosense Webster or compatible interface cable. To complete the electrical circuit, an indifferent electrode must be connected to the indifferent electrode input on the RF generator. Verify that circuit impedance prior to RF ablation is within expected parameters. Verify that the RF generator displays a temperature not above 37°C after the catheter is inserted into the patient and before applying RF power.

## RF GENERATOR OPERATION

Refer to the applicable RF generator manual for proper connection of the catheter to the generator and for detailed instructions as to generator operation for RF ablation.

RF ablation application parameters will vary depending on the ablation site, the specific conditions present in each procedure and the RF generator control circuitry. Based on data obtained from prior animal and clinical studies, recommended RF application parameters are provided below in the “Directions for Use” and in Table 1. Always monitor temperature and impedance rise when using the QDOT MICRO™ Uni-Directional Navigation Catheter.

## HOW SUPPLIED

- The QDOT MICRO™ Uni-Directional Navigation Catheter is supplied STERILE (EtO).
- The catheter is supplied with a choice of three curve configurations: “D”, “F”, and “J”.
- Additional catheter accessory devices are provided separately.

## PACKAGING

The catheter is secured in a one-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed and placed in a unit carton. The sealed pouch forms the sterile barrier. Both the interior of the pouch and the thermoform tray are sterile unless the package is damaged or opened. The pouch and unit carton are labeled to indicate that the device contained within is sterile.

## STORAGE

Store in a cool, dry, dark place. Storage temperature should be between 5 and 25°C (41 and 77°F).

## STERILIZATION AND USE-BY DATE

This catheter has been sterilized with ethylene oxide gas. Product and package testing have been conducted to support the “Use By” date printed on the product labels. **DO NOT USE** after the “Use By” date.

This device is packaged and sterilized for single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

## DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

## COMPATIBLE EP NAVIGATION SYSTEM

The Biosense Webster QDOT MICRO™ Uni-Directional Navigation Catheter provides location and contact force information only when used with CARTO™ 3 Navigation System. Compatibility with the CARTO™ 3 System has been demonstrated via bench, animal, and clinical testing to confirm that the device is capable of providing accurate location and contact force information when used in accordance with the Instructions for Use.

## DIRECTIONS FOR USE

Please refer to the user manuals for the CARTO™ 3 Navigation System, the irrigation pump, the irrigation tubing, and the generator for instructions on connecting and operating these systems in conjunction with the QDOT MICRO™ Uni-Directional Catheter. Use appropriate Biosense Webster accessory cables to connect the QDOT MICRO™ Uni-Directional Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place it in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.
2. The catheter tip can be deflected to facilitate positioning by using the thumb knob to vary tip curvature. Pushing the thumb knob forward causes the catheter tip to deflect. When the thumb knob is pulled back, the tip straightens.
3. Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use and aspiration per its instructions for use.
4. To verify compatibility between the sheath and catheter, advance the catheter through the sheath prior to insertion. Any sheath < 8.5 F is contraindicated. For compatible sheaths, contact Customer Support or your Biosense Webster representative.
5. Connect the irrigation pump tubing to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. Open the stopcock on the end of the tubing set and fill the tubing set as slowly as possible. Remove any trapped air and then close the stopcock.
6. Load the irrigation tubing into the pump. Open the stopcock and flush per the irrigation pump instructions until the air is expelled through the open end of the tubing.
7. Connect the irrigation tubing to the Luer fitting of the catheter.
8. Flush the catheter and tubing to ensure purging of trapped air bubbles and to verify irrigation through the tip electrode.
9. Start continuous irrigation with the irrigation pump at the low flow rate.
10. Connect the catheter to the CARTO™ System Patient Interface Unit (PIU) with the appropriate Biosense Webster cables. Connect the PIU to the generator and to the appropriate recording and mapping systems, including the CARTO™ 3 Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.
11. Insert the QDOT MICRO™ Uni-Directional Catheter via the entrance site.
12. Advance the catheter to the desired area. Use both direct imaging guidance (such as fluoroscopy or ultrasound) and electrograms (EGM) to aid in proper positioning.
13. In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up for 2 minutes after connection to the CARTO™ 3 System, prior to use of the force feedback feature.
14. Zero the contact force reading following insertion into the patient. The tip electrode and the two distal ring electrodes on the catheter tip must be outside of the sheath so that the force sensor is inside the body. Ensure the catheter tip is not in contact with tissue by evaluating the location on fluoroscopy, the CARTO™ System or other direct imaging guidance along with the EGM amplitude and catheter movement. Variations in the force reading at the same rate as the cardiac or respiration cycle may indicate contact with cardiac structures. Once these markers indicate the tip is not in contact, the reading can be zeroed. Refer to the user manual for the CARTO™ 3 System for instructions on how to zero the contact force reading.
15. Zero the contact force reading when moving the catheter from one chamber of the heart to another or upon reinsertion.
16. Verify that the “QDOT MICRO” option is selected on the RF generator. When this option is chosen, the RF generator defaults to the safety parameters established for the QDOT MICRO™ Uni-Directional Navigation Catheter.

17. Table 1 below outlines the recommended ablation parameter settings when using the QDOT MICRO™ Uni-Directional Navigation Catheter for ablation in the atria:

**Table 1: Recommended Ablation Parameter Settings**

Power Mode	QMODE*	QMODE+	
Power (W)	25-35	36-50	90
Target Electrode Temperature (°C)	50	50	60
RF application time (s)	Up to 60	Up to 60	4
Nominal Irrigation flow rate (ml/min)**	4***	15***	8

\* RF applications on the left atrial posterior wall using QMODE should not exceed 35 W in power and 30 sec in duration.

\*\* A minimum flow rate of 2 mL/min during mapping is recommended.

\*\*\* The irrigation flow rate is set by the generator and automatically adjusted between 4ml/min or 15ml/min to reach and maintain the set maximum power within the target electrode temperature.

- For the treatment typical atrial flutter, the QMODE is recommended.
- For the treatment of paroxysmal atrial fibrillation, the QMODE+ is recommended as the primary ablation mode for creating pulmonary vein isolation lesions and the QMODE is recommended for applications outside the PV ostia and for touch-up ablations.
- The QMODE power setting allows a maximum power up to 50 W, and the clinical operator may select any combination of power and RF application time, guided by the clinical judgment and expertise of the clinical operator as well as the anatomical area of interest.
- The QMODE+ power setting allows a maximum power up to 90 W, with a maximum RF application time up to 4 seconds. Use caution when overlapping RF ablation points, particularly on the posterior wall. Maintain appropriate contact force within the recommended range. This combination of power and RF application time is optimized for atrial ablation and may be used independently or in conjunction with the QMODE power setting.
- For ablations with any power setting, the clinical operator should monitor commonly used ablation effectiveness parameters like electrogram amplitude reduction and/or impedance drop.

Contact Force (CF) ranges during ablation:

The QFFICIENCY IDE study recommended CF settings to operators for use with the study catheter in both QMODE and QMODE+ as follows:

- CF Working Range: 5-30g

During the study, the overall average CF recorded during a study ablation procedure for all 166 subjects in the Per-Protocol Analyss set was 15.4 ± 5.2 grams.

18. Do not use this catheter without irrigation flow. Do not alter the irrigation flow rate.
19. Monitor the catheter tip temperature response throughout the procedure to ensure adequate irrigation. If temperature increases very rapidly during RF application, power delivery should be interrupted. The irrigation system must be rechecked prior to restarting RF application. Note: The displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.
20. In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode inspected for coagulum before RF energy is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean. Do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode, or damage may also occur to the contact force sensor and affect measurement accuracy. Prior to reinsertion, flush the tip to ensure that the irrigation holes are not plugged.

If irrigation hole occlusion occurs:

- a. Fill a 1-2 ml syringe with sterile saline and attach it to the stopcock on the end of the tubing set.
- b. Inject the saline from the syringe into the catheter. A uniform flow of fluid should be visible from the tip of the catheter.
- c. Repeat steps a and b, if necessary, until the holes are cleared.
- d. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
- e. The catheter can now be introduced into the patient.
- f. Zero the catheter following reinsertion into the patient.

**WARNING: Do not continue use of the catheter if it is still occluded or if it is not functioning properly.**

## **DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY**

**THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.**

**WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.**

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

**Biosense Webster, Inc.**

31 Technology Drive, Suite 200  
Irvine, California 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010  
Fax: +1-909-468-2905  
[www.biosensewebster.com](http://www.biosensewebster.com)

BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, QDOT MICRO, THERMOCOOL SMARTTOUCH, THERMOCOOL, CARTO SMARTTOUCH, and NAVISTAR are trademarks of Biosense Webster, Inc.

The third-party trademarks used herein are trademarks of their respective owners.

Patent information is available at [www.biosensewebster.com/virtualpatentmarking](http://www.biosensewebster.com/virtualpatentmarking).

Revised: 2023-01  
© Biosense Webster, Inc. 2010-2023

## **QDOT MICRO™ Bi-Directional Navigation Catheter**

### **INSTRUCTIONS FOR USE**

DEVICE DESCRIPTION .....	4
INDICATIONS FOR USE .....	4
CONTRAINDICATIONS .....	5
WARNINGS AND PRECAUTIONS.....	5
ADVERSE EVENTS.....	8
SUMMARY OF CLINICAL STUDIES CONDUCTED FOR THE QDOT MICRO™ CATHETER .....	8
A. Objective .....	8
B. Study Design.....	9
C. Results.....	13
D. Study Conclusion.....	29
RF ABLATION.....	30
RF GENERATOR OPERATION .....	30
HOW SUPPLIED.....	30
PACKAGING.....	30
STORAGE.....	30
STERILIZATION AND USE-BY DATE.....	30
DISPOSAL .....	31
COMPATIBLE EP NAVIGATION SYSTEM .....	31
DIRECTIONS FOR USE .....	31
DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY .....	34

---

## Symbol Definitions

Following definitions are for reference only. Please refer to the product label for applicable usage.

**STERILE EO**

Sterilized Using Ethylene Oxide

Do Not Re-Use

Caution

Consult Instructions for Use

Do not use if package is damaged.

Do not use if package is opened.

Keep Away from Sunlight

Keep Dry

Use-By Date

**LOT**

Batch Code

**REF**

Catalog Number

Contents: 1

Manufacturer

Date of Manufacture

Pin Connector

**ELECTRODES EA**

Electrodes

**SPACING mm**

Spacing

Temperature Limit

---

**R<sub>x</sub>** only

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Curve Type. Refer to label for colored circle containing applicable curve type.

**NAV**

Navigational Catheter

**C3**

Compatible with CARTO™ 3 EP Navigation System

## QDOT MICRO™ Bi-Directional Navigation Catheter

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

- STERILE. Sterilized using ethylene oxide.
- For single use only.
- Do not resterilize.
- Do not use if the package is open or damaged.

### DEVICE DESCRIPTION

The Biosense Webster QDOT MICRO™ Bi-Directional Catheter is a steerable multi-electrode luminal catheter with a deflectable tip designed to facilitate electrophysiological mapping of the heart and to transmit radiofrequency (RF) energy to the catheter tip electrode for ablation purposes. The catheter shaft measures 7.5 F with 8 F ring electrodes. For ablation, the catheter is used in conjunction with a compatible RF generator and a dispersive pad (indifferent electrode). The catheter has force-sensing technology that provides a real-time measurement of contact force between the catheter tip and the heart wall.

The catheter has a high-torque shaft with a bi-directional deflectable tip section containing an array of electrodes which includes a 3.5 mm tip dome. All of the electrodes may be used for recording and stimulation purposes. The tip electrode serves to deliver RF energy from the generator to the desired ablation site. The tip electrode and ring electrodes are made from noble metals. The catheter incorporates six thermocouple temperature sensors and three ECG electrodes that are embedded in the 3.5 mm tip electrode. A Rocker Lever is used to deflect the tip. The high-torque shaft also allows the plane of the curved tip to be rotated to facilitate accurate positioning of the catheter tip at the desired site. Additionally, a variety of curve types are available in symmetric or asymmetric combinations, providing two 180° opposed, single-planed curves. Five curve configurations designated “DD”, “FF”, “JJ”, “DF”, and “FJ” are available.

At the proximal end of the catheter, a saline input port with a standard Luer fitting terminates from the open lumen. This saline port serves to permit the injection of normal saline to irrigate the tip electrode. During ablation, heparinized normal saline is passed through the internal lumen of the catheter and through the tip electrode, to irrigate and cool the ablation site as well as the electrode tip. A compatible irrigation pump is used to control the saline irrigation. The catheter interfaces with standard recording equipment and a compatible RF generator via accessory extension cables with the appropriate connectors.

This catheter features a location sensor embedded in the tip section that transmits location and contact force information to the CARTO™ 3 Navigation System. An appropriate reference device is required for location reference position purposes. For information on using the catheter in mapping procedures and for information on appropriate reference devices, refer to the user manual for the CARTO™ 3 Navigation System. For further description of the operation of the irrigation pump, the generator, and the CARTO™ 3 Navigation System, refer to the applicable instructions for use.

### INDICATIONS FOR USE

The Biosense Webster QDOT MICRO™ Bi-Directional Navigation Catheter and related accessory devices are indicated for catheter-based cardiac electrophysiological mapping (stimulating and recording) and, when used in conjunction with a compatible radiofrequency generator, for the treatment of:

- Type I atrial flutter in patients age 18 or older.
- Drug refractory recurrent symptomatic paroxysmal atrial fibrillation, when used with compatible three-dimensional electroanatomic mapping systems.

The Biosense Webster QDOT MICRO™ Bi-Directional Navigation Catheter provides a real-time measurement of contact force between the catheter tip and heart wall, as well as location information when used with CARTO™ 3 Navigation System.

## CONTRAINDICATIONS

Do not use this catheter:

1. If the patient has had a ventriculotomy or atriotomy within the preceding twelve weeks because the recent surgery may increase the risk of perforation.
2. In patients with a myxoma or an intracardiac thrombus as the catheter could precipitate an embolus.
3. In patients with prosthetic valves as the catheter may damage the prosthesis.
4. In the coronary arterial vasculature due to risk of damage to the coronary arterial vasculature.
5. In patients with an active systemic infection because this may increase the risk of cardiac infection.
6. Via the transseptal approach in a patient with an interatrial baffle or patch because the opening could persist and produce an iatrogenic atrial shunt.
7. Via the retrograde trans-aortic approach in patients who have had aortic valve replacement.
8. With a long sheath or short introducer < 8.5 F in order to avoid damage to the catheter shaft.

## WARNINGS AND PRECAUTIONS

1. Do not use excessive force to advance or withdraw the catheter when resistance is encountered during catheter manipulation through the sheath.
2. Do not manually pre-shape the distal shaft of the catheter by applying external forces intended to bend or affect the intended shape or curve of the catheter.
3. Prior to use, the catheter must be warmed up as specified in the Directions for Use section of this document. If the catheter has not reached a steady state condition, there is potential for a zero-offset drift to occur which could result in an inaccurate contact force reading.
4. Always zero the contact force reading following insertion into the patient or when moving the catheter from one chamber of the heart to another. Ensure the catheter is not in contact with heart tissue prior to zeroing. Refer to the user manual for your CARTO™ 3 System for instructions on how to zero the contact force reading.
5. The contact force reading might become inaccurate if the contact force sensor (located between the first and second ring electrode) comes into close proximity with a ferrous material, such as the braided shaft of another catheter. If extreme fluctuations in force are observed, ensure the catheter's contact force sensor is not in close proximity with another catheter's shaft, check zero on the catheter and, if necessary, remove and inspect the catheter. The contact force reading is for information only and is not intended to replace standard handling precautions.
6. To ensure proper operation of the contact force sensor, the tip electrode and the two distal ring electrodes must protrude from the distal tip of the guiding sheath.
7. When applying high lateral force during mapping and RF application, the user should monitor the contact force dashboard and vector display on the CARTO™ 3 System screen to ensure that contact force measurements remain within the accurate range. Refer to the Error Messages and Alerts section of the CARTO™ 3 System user manual for system-related alert messages and indications related to inaccurate force readings.
8. Do not use the temperature sensor to monitor tissue temperature. The temperature sensor located within the tip section of the catheter does not reflect either electrode-tissue interface or tissue temperature due to the cooling effects of the saline irrigation of the electrode. The temperature displayed on the generator is the temperature of the hottest of the six thermocouples in the electrode, not tissue temperature. Monitoring the temperature from the electrode during the application of RF energy ensures that the irrigation flow rate is being maintained.
9. At higher contact force settings, rapid power increases during ablation may result in steam pops, which increase the risk of perforation.
10. This catheter may damage the mechanical or prosthetic valve of a patient if the catheter is accidentally advanced through the valve.
11. The safety of discontinuing anticoagulation therapy following catheter ablation of atrial fibrillation has not been established. Anticoagulation therapy in such patients should be administered in accordance with the AHA/ACC/HRS 2014 Guideline for the Management of Patients With Atrial Fibrillation.

12. The safety and effectiveness of radiofrequency ablation for the treatment of atrial fibrillation in patients with significant left ventricular dysfunction, advanced heart failure, substantial left atrial enlargement, and structural heart disease have not been established.
13. In accordance with your hospital's protocol, monitor the patient's fluid balance throughout the procedure to avoid fluid volume overload. Some patients may have factors that reduce their ability to handle the volume overload, making them susceptible to developing pulmonary edema or heart failure during or after the procedure. Patients with congestive heart failure or renal insufficiency, and the elderly are particularly susceptible. Prior to the procedure, always identify the patient's risk of volume overload.
14. Implantable pacemakers and implantable cardioverter/defibrillator (ICDs) may be adversely affected by RF energy. It is important to have temporary external sources of pacing and defibrillation available during ablation and to temporarily reprogram the pacing system to minimum output or OFF mode to minimize the risk of inappropriate pacing. Exercise extreme caution during ablation when in close proximity to atrial or ventricular permanent leads; program the ICD to the OFF mode during the ablation procedure; and, perform complete implantable device analysis on all patients after ablation.
15. Patients undergoing septal accessory pathway ablation are at risk for complete AV block which requires the implantation of a permanent pacemaker. Patients who experience inadvertent complete AV block as a result of RF ablation may also require permanent pacing.
16. During the trans-aortic approach, adequate visualization is necessary to avoid placement of the catheter in the coronary vasculature. Intracoronary placement of the ablation catheter, RF energy application, or both have been associated with myocardial infarction.
17. To prevent stenosis of the pulmonary veins, do not place the catheter in the pulmonary veins during the application of RF energy.
18. When ablating near adjacent anatomical structures, take precautions to minimize collateral damage to the adjacent structures.
19. When ablating near the phrenic nerve, take precautions to avoid injuring the phrenic nerve, including appropriately reducing RF power, pacing to identify the proximity of ablation electrode(s) to the nerve and/or fluoroscopic evaluation of the diaphragm post ablation.
20. When ablating near the esophagus (along the posterior wall of the left atrium), take precautions to avoid injuring the esophagus. These may include beginning the ablation with reduced RF power, reducing contact force, reducing application time, increasing the time interval between ablations, moving to new location if temperature rise observed, esophageal visualization, and/or intraluminal esophageal temperature monitoring.
21. Minimize or eliminate X-ray exposure during the procedure. Catheter ablation procedures using fluoroscopic imaging may present the potential for significant X-ray exposure, which can result in acute radiation injury as well as increased risk for somatic and genetic effects, to both patients and laboratory staff due to the X-ray beam intensity and duration of fluoroscopic imaging. Catheter ablation, when using fluoroscopic imaging, should only be performed after adequate attention has been given to the potential radiation exposure associated with the procedure, and steps have been taken to minimize this exposure. Careful consideration must therefore be given for the use of the device in pregnant women.
22. Do not expose the catheter to organic solvents such as alcohol.
23. Do not autoclave the catheter.
24. Do not immerse the proximal handle or cable connector in fluids. Electrical performance could be affected.
25. Do not scrub or twist the distal tip electrode during cleaning.
26. Inspect the saline within the irrigation tubing for air bubbles prior to its use in the procedure. Air bubbles in the irrigation saline may cause emboli.
27. Purge the catheter and irrigation tubing with heparinized normal saline prior to insertion of the catheter inside the patient body.
28. Electrophysiology catheters and systems are intended for use only in X-ray shielded rooms due to electromagnetic compatibility requirements and other hospital safety guidelines.
29. Do not attempt to operate the QDOT MICRO™ Bi-Directional Catheter or the generator prior to completely reading and understanding the applicable instructions for use.

30. Cardiac ablation procedures should be performed by appropriately trained personnel in a fully equipped electrophysiology laboratory. Appropriate clinical instruction in the use of the QDOT MICRO™ Bi-Directional Catheter should also be completed.
31. The long-term risks of protracted fluoroscopy and creation of RF induced lesions have not been established. Furthermore, the risk/benefit in asymptomatic patients has not been studied.
32. To prevent thromboembolism, intravenous heparin (target ACT of  $\geq 350$  s) should be administered prior to or immediately following transeptal puncture during AF ablation procedures. The 2017 HRS/EHRA/ECAS/AOHRS/SOLAECCCE expert consensus statement on catheter and surgical ablation of atrial fibrillation recommends systemic anticoagulation with warfarin or a direct thrombin or factor Xa inhibitor for at least 2 months following an AF ablation procedure (Calkins H, Hindricks G, Cappato R et al. 2017 HRS/EHRA/ECAS/AOHRS/SOLAECCCE Expert Consensus Statement on Catheter and Surgical Ablation of Atrial Fibrillation. J Interv Card Electrophysiol. 2017 50:1-55).
33. When using the QDOT MICRO™ Bi-Directional Catheter with conventional systems (such as fluoroscopy or ultrasound imaging), or with the CARTO™ 3 Navigation System, careful catheter manipulation must be performed in order to avoid cardiac damage, perforation, or tamponade. Catheter advancement should be done under direct imaging guidance. Do not use excessive force to advance or withdraw the catheter when resistance is encountered. The firmness of the braided tip dictates that care must be taken to prevent perforation of the heart. The contact force reading is for information only and is not intended to replace standard handling precautions.
34. Always place the Rocker Lever in the neutral position to straighten the catheter tip before insertion or withdrawal of the catheter.
35. Always maintain a constant heparinized normal saline infusion to prevent coagulation within the lumen of the catheter.
36. When RF energy is interrupted for either a temperature or an impedance rise (the set limit is exceeded), the catheter should be removed, and the tip should be inspected for char/coagulum that may be present on the tip. If present, do not continue the procedure with the same catheter and replace the catheter. If no char/coagulum is present, flush the tip to ensure the irrigation holes are not plugged prior to reinsertion of the catheter inside the patient body.
37. Apparent low power output, high impedance reading, or failure of the equipment to function correctly at normal settings may indicate faulty application of the indifferent electrode(s) or failure of an electrical lead. Do not increase power before checking for obvious defects or misapplication of the indifferent electrode or other electrical leads.
38. Read and follow the indifferent electrode manufacturer's instructions for use. The use of indifferent electrodes that meet or exceed ANSI/AAMI requirements (AAMI IEC 60601-2-2), is recommended.
39. The QDOT MICRO™ Bi-Directional Catheter is intended for use with a compatible RF generator, compatible irrigation pump, CARTO™ 3 Navigation System, Biosense Webster cables, and other appropriate interface cables and connectors. Use of a compatible irrigation pump is recommended to assure proper irrigation flow rate.
40. Care should be taken when ablating near structures such as the sino-atrial and atrioventricular nodes.
41. The sterile packaging and catheter should be inspected prior to use. Do not use if the packaging or catheter appears damaged.
42. The catheter is sterilized with ethylene oxide gas and should be used by the "Use By" date on the device package. Do not use the catheter if it is past the "Use By" date.
43. The QDOT MICRO™ Bi-Directional Catheter is intended for single patient use only.
44. Do not resterilize and reuse.
45. Do not use near MRI equipment since movement or heating of the catheter may occur and the image on the display may become distorted.
46. Use both direct imaging guidance (such as fluoroscopy or ultrasound) and electrogram data to monitor catheter advancement and reduce risk of tissue injury.
47. The QDOT MICRO™ Bi-Directional Catheter used in conjunction with a compatible RF generator is capable of delivering significant electrical power. Patient or operator injury can result from improper handling of the catheter and indifferent electrode, particularly when operating the catheter. During energy delivery, the patient should not be allowed to come in contact with grounded metal surfaces.

48. The risk of igniting flammable gases or other materials is inherent in electrosurgery. Precautions must be taken to restrict flammable materials from the electrosurgical suite.
49. Electromagnetic interference (EMI) produced by the QDOT MICRO™ Bi-Directional Catheter, when used in conjunction with a compatible RF generator during normal operation, may adversely affect the performance of other equipment.
50. Electrodes and probes used for monitoring and stimulating devices can provide paths for high frequency current. The risk of burns can be reduced but not eliminated by placing the electrodes and probes as far away as possible from the ablation site and the indifferent electrode. Protective impedance may reduce the risk of burns, and permit continuous monitoring of the electrocardiogram during energy delivery.
51. If the generator does not display temperature, verify that the cables from the catheter to the CARTO™ System Patient Interface Unit (PIU) and the cable from the CARTO™ System Patient Interface Unit (PIU) to the generator are properly connected. If temperature still is not displayed, there may be a malfunction in the temperature sensing system that must be corrected prior to applying RF power.
52. The temperature measurement accuracy of the QDOT MICRO™ Bi-Directional Catheter, as with any temperature measurement electrophysiology catheter, is determined by the temperature accuracy of all of the connected devices. Please consult the appropriate user manuals for the connected devices for the temperature accuracy specification.
53. Before use, verify irrigation ports are patent by infusing heparinized normal saline through the catheter and tubing.
54. Regularly inspect and test reusable cables and accessories.

## ADVERSE EVENTS

---

**NOTE: The adverse events in the following summary were observed in clinical studies involving only the use of the QDOT Micro™ Catheter.**

---

### **Clinical trial for the QDOT Micro™ Catheter – Pivotal Study (Q-FFICIENCY)**

Of the 177 subjects in the Safety Analysis Set, 7 primary adverse events (AEs) were reported in 6 subjects. See the “Summary of Clinical Studies Conducted for the QDOT Micro™ Catheter” section below for a complete description of the AEs encountered during the study.

## SUMMARY OF CLINICAL STUDIES CONDUCTED FOR THE QDOT MICRO™ CATHETER

The following is a summary of the pivotal study (Q-FFICIENCY) performed under IDE for the QDOT Micro™ Catheter, using model numbers D-1394-XX-SI and D-1395-XX-SI. Data from this clinical study formed the basis for the US FDA PMA approval decision.

All patients underwent informed consent per the study protocol, and in compliance with the Code of Federal Regulations, 21 §50. The protocol and informed consent materials were reviewed and approved by the appropriate IRB prior to subject enrollment, and in compliance with the Code of Federal Regulations, 21 §56. Per the requirements of 21 CFR §812.20(b)(5), all participating investigators signed the Clinical Study Agreement prior to enrollment at their center, which included financial disclosure forms, current curricula vitae for the investigator and co-investigators, and written IRB approval for the investigational study. The clinical study was monitored in a manner consistent with 21CFR, Part 812, Subpart C, *Responsibilities of Sponsor*. Monitoring visits included, but were not limited to, verification of all study logs, verification that informed consent was being obtained in accordance with requirements described in the study protocol for all subjects participating in the study, verification of completeness of the Regulatory Binder, data source verification with the eCRFs, and identification and action to resolve any issues or problems with the study.

### **A. Objective**

The primary objective of this trial was to demonstrate the safety and effectiveness of the QDOT MICRO™ catheter when used with the nMARQ™ RF generator in the treatment of drug refractory symptomatic paroxysmal atrial fibrillation (PAF) during standard electrophysiology mapping and RF ablation procedures.

**B. Study Design**

The study was a prospective, single arm, unblinded, multicenter, pivotal clinical investigation conducted at 25 investigational sites in the US.

The study planned to enroll up to 185 subjects who are candidates for catheter ablation. Bayesian adaptive design was used to assess early success at one interim analysis performed 6 months after enrollment was completed and when 30% of subjects in the study completed 12 months of follow-up.

A CT/MRA imaging sub-study was integrated in the pivotal study. The first 40 consecutively enrolled subjects were included in the sub-study and underwent a 3-month CT/MRA in addition to the baseline CT/MRA to assess incidence of PV stenosis.

**B.1 - Study Objectives and Endpoints:**

The objectives and endpoints for this study are presented in Table 1. Failure modes for the primary effectiveness endpoint are listed in Table 2.

**Table 1 – Objectives and Endpoints for Q-EFFICIENCY Study**

Objectives	Endpoints
<b>Primary</b>	
<p>To demonstrate the safety and 12-month effectiveness of the QDOT MICRO™ catheter when used with the nMARQ™ RF generator for pulmonary vein isolation (PVI) in the treatment of subjects with drug refractory symptomatic paroxysmal atrial fibrillation (PAF).</p>	<ul style="list-style-type: none"> <li>• Primary safety endpoint: Incidence of any primary adverse event occurring within 7 days of the AF ablation procedure (initial or repeat procedure) using the QDOT MICRO™ catheter, except atrio-esophageal fistula and PV stenosis, which may also be considered as primary adverse events if occurring greater than 7 days and up to 90 days post the ablation procedure. Primary adverse events included the following: Death, myocardial infarction, PV stenosis, phrenic nerve injury/diaphragmatic paralysis, atrio-esophageal fistula, TIA, stroke/CVA, thromboembolism, pericarditis, cardiac tamponade/perforation, vagal nerve injury, major vascular access complications/bleeding, pulmonary edema (respiratory insufficiency), and heart block.</li> <li>• Primary effectiveness endpoint: Freedom from documented atrial fibrillation (AF), atrial tachycardia (AT) or atrial flutter (AFL) (hereinafter collectively referred to as “atrial tachyarrhythmias”) recurrence during the evaluation period (Day 91 through Day 365) and freedom from the following failure modes (defined in Table 2):                         <ul style="list-style-type: none"> <li>○ Acute procedure failure</li> <li>○ Repeat ablation failure</li> <li>○ Antiarrhythmic drug failure</li> </ul> </li> </ul>

Objectives	Endpoints
<b>Secondary</b>	
<p>To evaluate the incidence of (serious) adverse events during and after procedure up to 3 months following procedure</p>	<ul style="list-style-type: none"> <li>• Incidence of Unanticipated Adverse Device Effects (UADEs)</li> <li>• Incidence of Serious Adverse Events (SAEs) within 7 days (early onset), &gt; 7 to 30 days (peri-procedural) and &gt; 30 days (late onset) of initial ablation</li> <li>• Incidence of bleeding complication (ISTH definitions): a) major, b) clinically relevant non-major and c) minor bleeding</li> </ul>
<p>To evaluate Acute Procedural Success defined as confirmation of entrance block in all PVs</p>	<ul style="list-style-type: none"> <li>• Acute procedural success: Percent of subjects with electrical isolation of PVs (entrance block) at the end of the procedure                             <ul style="list-style-type: none"> <li>○ The percent of subjects with electrical isolation of PVs (entrance block) using QMODE+ as only ablation strategy</li> </ul> </li> <li>• Percent of subjects with electrical isolation of PVs (entrance block) after first encirclement (evaluated prior to the 20-minute waiting period and adenosine/isoproterenol challenge)</li> <li>• Percent of subjects with electrical isolation of all PVs (entrance block) after first encirclement without acute reconnection, after waiting period and adenosine/isoproterenol challenge</li> <li>• Percent of subjects and % of PVs with touch-up (i.e. touch-up is used to remove ablation of acute reconnection) among all targeted veins and touch-up location</li> <li>• Anatomical location of acute PV reconnection after first encirclement</li> </ul>
<p>To evaluate repeat ablation procedures during the 12-month period post-procedure</p>	<ul style="list-style-type: none"> <li>• Incidence (%) of repeat ablation procedures</li> <li>• Percent PVs re-isolated among all the targeted PVs at repeat procedure</li> <li>• Percent repeat ablation procedures requiring new linear lesions and/or identifying new foci outside of initially isolated area among the repeat ablation procedures</li> </ul>
<p>To evaluate the 12-month Single Procedure Success</p>	<ul style="list-style-type: none"> <li>• 12-month single procedure success: Freedom from documented AF/AFL/AT recurrence (episodes ≥ 30 secs) during the Effectiveness Evaluation Period after a single ablation procedure and off anti-arrhythmia drugs (AADs) during the effectiveness evaluation period</li> </ul>

**Table 2 – Failure Modes for Primary Effectiveness Endpoint**

Failure Mode	Criteria
Acute procedural	<ul style="list-style-type: none"> <li>Failure to confirm entrance block in all pulmonary veins post-procedure.</li> <li>Use of a non-study catheter to treat left atrial ablation targets and cavo-tricuspid isthmus.</li> </ul>
Repeat ablation	<ul style="list-style-type: none"> <li>&gt; 2 repeat ablation procedures with the study catheter during the 3-Month blanking period (Day 0-90) after the index ablation procedure.</li> <li>Use of a non-study catheter to treat study arrhythmia ablation targets during the blanking period.</li> <li>Any repeat ablation procedure during the Evaluation Period.</li> </ul>
Antiarrhythmic drug	<ul style="list-style-type: none"> <li>Taking a new Class I or III AAD for AF or a previously failed Class I or III AAD at a greater than the highest ineffective historical dose for AF during the evaluation period.</li> </ul>

**B.2 - Safety and Effectiveness Performance Goals**

**B.2.1 - Safety Performance Goal**

The study was designed to compare the primary safety endpoint to a pre-determined safety performance goal of 14%. Data from recent PAF ablation studies for devices similar to the device in the current study were reviewed as a first step to deriving the performance goal for the primary safety endpoint. A meta-analysis approach was taken to estimate the average composite endpoint rate. Based on the meta-analysis, the upper bound of the 95% confidence interval was estimated to be equal to 9%. The pre-determined performance goal of 14% would reflect an approximately 50% increase in risk from the upper bound of the 95% CI.

**B.2.2 - Effectiveness Performance Goal**

The study was designed to compare the primary effectiveness endpoint to a pre-determined performance goal of 50%, which is indicated as the minimum acceptable success rate at 12 months for a paroxysmal AF population in the 2017 HRS/EHRA/ECAS/APHS/ SOLAECE Expert Consensus Statement on Catheter and Surgical Ablation of AF (Calkins, et al., 2017).

**B.3 - Subject Accountability:**

**Table 3 – Subject Accountability and Disposition**

Subject Enrollment	n
Enrolled Subjects	191
Excluded Subjects	14
Study Catheter Inserted	177
Discontinued Subjects	2
RF Energy Delivered with Study Catheter	175
Death	1
Withdrawn	7
Early Termination Subjects	0
Lost to Follow Up Subjects	5
Completed Subjects	162

Among the 191 enrolled subjects, 14 had no insertion of the study device and were excluded from the analysis sets for the following reasons:

- System connectivity related device deficiencies (n = 4)

- Not meeting study entrance criteria (n = 9)
- Occurrence of an adverse event (tamponade) before insertion of the study device (n = 1)

Two (2) subjects had the study catheter inserted but no RF applications were delivered due to system connectivity issues.

A total of 175 subjects received RF ablation using the study device. Among those, seven (7) subjects (4.0%) withdrew their consent and discontinued their study participation during follow-up from 48 to 384 days post ablation; five (5) subjects (2.9%) were lost to follow up between 104 days and 266 days post procedure; one (1) subject (0.6%) expired 7 months post procedure; and 162 subjects (92.6%) completed the 12-month follow-up visit.

**B.3.1 - Subject Disposition**

The following definitions were used to classify subjects:

**Table 4 – Subject Classifications**

<b>Disposition</b>	<b>Definition</b>
Enrolled Subjects	Subjects who signed the informed consent.
Excluded Subjects	Subjects who were enrolled but never underwent insertion of the study catheter.
Discontinued Subjects	Subjects who had the investigational catheter inserted but did not undergo ablation (i.e., no RF energy was delivered via the study device). These discontinued subjects were followed up for 3 months.
Lost to Follow-up Subjects	Subjects who were enrolled, but contact was lost after the most recent follow-up visits (despite 3 documented attempts to contact the subject).
Withdrawn / Early Termination Subjects	Subjects who withdrew consent for study participation or were withdrawn by the investigator or were terminated from the study prior to completion of all follow-up visits.
Completed Subjects	Enrolled subjects who completed the 12-month follow-up visit and were not excluded, discontinued, withdrawn, early terminated, expired, or lost to follow-up from the study before the final study visit.

**B.3.2 - Analysis Sets**

**Safety Analysis Set:** The Safety analysis set included all enrolled subjects who had the investigational device inserted, regardless of RF energy delivery.

**Modified Intent-To-Treat (mITT) Analysis Set:** The mITT analysis set consisted of all enrolled subjects who met all eligibility criteria and had the investigational device inserted.

**Per-Protocol (PP) Analysis Set:** The PP analysis set included subjects who satisfied the following criteria:

- Enrolled and met eligibility criteria
- Had undergone RF ablation
- Were treated with the study catheters and had been treated for the study related arrhythmia without major protocol deviations deemed to affect the scientific integrity of the study, included in the following list. No additional major protocol deviations beyond those below were identified that required exclusion of subjects from the Per Protocol analysis set
  - Missing all protocol-specified electronic effectiveness monitoring
  - Esophageal monitoring was not done per protocol and risk was not mitigated
  - The required waiting period of 20 minutes was not followed before pacing procedures and/or infusion of cardiac medications to induce AF/reconnection was not performed before the entrance block confirmation.

**Treatment Received Analysis Set:** The Treatment Received analysis set included all enrolled subjects who received RF ablation using the study catheter. This analysis set was created post-hoc because it was deemed relevant by the FDA.

**CT/MRA Analysis Set:** The CT/MRA set consisted of the first 40 subjects who are consecutively enrolled in the main arm of the study with 3-month CT/MRA assessment. The subject must have readable outcomes at baseline and 3 months.

**C. Results**

**C.1 - Subject Demographics:**

**Table 5 – Subject Demographics (Enrolled Subjects, N=191)**

Variables	All Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
Age (years) [1]				
n	191	177	167	166
Mean	63.5	63.3	63.1	63.2
SD	10.69	10.88	11.00	11.03
Min/Max	24 / 83	24 / 83	24 / 83	24 / 83
Sex, n/N (%)				
Male	116 / 191 (60.7)	106 / 177 (59.9)	101 / 167 (60.5)	101 / 166 (60.8)
Female	75 / 191 (39.3)	71 / 177 (40.1)	66 / 167 (39.5)	65 / 166 (39.2)
Race n (%)				
White, n/N (%)	166 / 191 (86.9)	157 / 177 (88.7)	149 / 167 (89.2)	148 / 166 (89.2)
Black or African American	7 / 191 (3.7)	6 / 177 (3.4)	6 / 167 (3.6)	6 / 166 (3.6)
Asian	1 / 191 (0.5)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
Race not reported	17 / 191 (8.9)	13 / 177 (7.3)	11 / 167 (6.6)	11 / 166 (6.6)
Ethnicity, n/N (%)				
Not Hispanic or Latino	164 / 191 (85.9)	154 / 177 (87.0)	146 / 167 (87.4)	145 / 166 (87.3)
Hispanic or Latino	8 / 191 (4.2)	6 / 177 (3.4)	6 / 167 (3.6)	6 / 166 (3.6)
Not reported	19 / 191 (9.9)	17 / 177 (9.6)	15 / 167 (9.0)	15 / 166 (9.0)

**C.2 - Baseline Characteristics:**

**Table 6 – Medical History (Enrolled Subjects, N=191)**

Medical History	Enrolled Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
<b>Cardiovascular</b>	156 / 191 (81.7)	143 / 177 (80.8)	135 / 167 (80.8)	134 / 166 (80.7)
Congestive heart failure	16 / 191 (8.4)	10 / 177 (5.6)	9 / 167 (5.4)	9 / 166 (5.4)

Medical History	Enrolled Subjects (N=191)	Safety Analysis Set (N=177)	mITT Analysis Set (N=167)	Per-Protocol Analysis Set (N=166)
Coronary artery disease	39 / 191 (20.4)	36 / 177 (20.3)	33 / 167 (19.8)	33 / 166 (19.9)
Vascular disease	25 / 191 (13.1)	23 / 177 (13.0)	22 / 167 (13.2)	22 / 166 (13.3)
Myocardial infarction	13 / 191 (6.8)	12 / 177 (6.8)	10 / 167 (6.0)	10 / 166 (6.0)
Hypertension	131 / 191 (68.6)	122 / 177 (68.9)	116 / 167 (69.5)	115 / 166 (69.3)
Cardiomyopathy	15 / 191 (7.9)	11 / 177 (6.2)	10 / 167 (6.0)	10 / 166 (6.0)
Left ventricular hypertrophy	1 / 191 (0.5)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
Significant valvular disease	3 / 191 (1.6)	1 / 177 (0.6)	1 / 167 (0.6)	1 / 166 (0.6)
<b>Thromboembolic Events</b>	20 / 191 (10.5)	18 / 177 (10.2)	16 / 167 (9.6)	16 / 166 (9.6)
Transient ischemic attack	5 / 191 (2.6)	4 / 177 (2.3)	3 / 167 (1.8)	3 / 166 (1.8)
Stroke	6 / 191 (3.1)	5 / 177 (2.8)	5 / 167 (3.0)	5 / 166 (3.0)
Pulmonary embolus	2 / 191 (1.0)	2 / 177 (1.1)	2 / 167 (1.2)	2 / 166 (1.2)
<b>Diabetes</b>	37 / 191 (19.4)	36 / 177 (20.3)	33 / 167 (19.8)	33 / 166 (19.9)
<b>Obstructive Sleep Apnea</b>	50 / 191 (26.2)	46 / 177 (26.0)	44 / 167 (26.3)	44 / 166 (26.5)
CPAP Use	36 / 191 (18.8)	34 / 177 (19.2)	33 / 167 (19.8)	33 / 166 (19.9)
<b>CHA<sub>2</sub>DS<sub>2</sub>-VASc score</b>	2.4 ± 1.51	2.4 ± 1.53	2.4 ± 1.49	2.4 ± 1.50
<b>AF History (months)</b>	48.6 ± 73.63	49.6 ± 75.29	51.7 ± 76.92	51.9 ± 77.11
<b>PAF Episodes Within 12 Months</b>	24.9 ± 89.91	25.2 ± 92.80	26.0 ± 95.22	25.7 ± 95.38
<b>AAD Failed- Class I&amp;III</b>	1.2 ± 0.50	1.2 ± 0.50	1.2 ± 0.52	1.2 ± 0.52
<b>Arrhythmia other than Atrial Fibrillation</b>	87 / 191 (45.5)	81 / 177 (45.8)	73 / 167 (43.7)	73 / 166 (44.0)
Left/right Atrial Tachycardia (AT)	8 / 191 (4.2)	8 / 177 (4.5)	7 / 167 (4.2)	7 / 166 (4.2)
Supraventricular Tachycardia (SVT)	22 / 191 (11.5)	20 / 177 (11.3)	14 / 167 (8.4)	14 / 166 (8.4)
Atrial flutter (AFL)	43 / 191 (22.5)	41 / 177 (23.2)	40 / 167 (24.0)	40 / 166 (24.1)
AVNRT	0 / 191 (0.0)	0 / 177 (0.0)	0 / 167 (0.0)	0 / 166 (0.0)
Accessory pathway	0 / 191 (0.0)	0 / 177 (0.0)	0 / 167 (0.0)	0 / 166 (0.0)
Ventricular tachycardia	2 / 191 (1.0)	2 / 177 (1.1)	2 / 167 (1.2)	2 / 166 (1.2)
<b>LA Diameter (mm)</b>	38.06 ± 6.042	37.95 ± 5.917	38.09 ± 5.974	38.16 ± 5.917
<b>LVEF (%)</b>	59.4 ± 7.31	59.6 ± 7.02	59.7 ± 6.98	59.7 ± 6.99

Values in table represent n/N (%) or mean ± SD, as appropriate.

### C.3 - Index Ablation Procedure

#### **Ablation protocol**

A circumferential anatomical approach was used to isolate all PVs. Confirmation of entrance block in all PVs was required by the protocol with a 20-minute waiting period after the last PV encircling lesion, during isoproterenol infusion and/or after adenosine bolus.

The study ablation procedure used both QMODE and QMODE+ temperature control modes. The primary mode of ablation for PVI was QMODE+. QMODE was used for PVI once the investigator deemed

QMODE+ unable to complete PVI. QMODE was used primarily for RF application outside the PV ostia and for touch-up of the PVI, if necessary.

Table 7 and Table 8 present the ablation parameters for the QMODE+ mode and QMODE mode, respectively.

**Table 7 – QMODE+ RF and Flow Rate Settings during RF applications**

Power (W)	Target Temp (°C)		Cut-off Temp (°C)		Nominal Irrigation Flow rate (mL)
	Range	Maximum allowed	Range	Maximum allowed	
90†	40-60	60	60-70	70	8

† RF applications at this setting were limited to 4 sec. It was recommended to use lower target temperature setting for the posterior wall RF applications.

**Table 8 – QMODE Ablation Parameters\***

Power (W)	Target Temp (°C)		Cut-off Temp (°C)		Nominal Irrigation Flow rate (mL)
	Range	Recommended	Range	Recommended	
25-35	45-50	50	50-55	55	4
36-50	45-50	50	50-55	55	15

\* The study protocol recommended a maximum duration of 60 sec for RF applications using QMODE in general. For RF applications on the left atrial posterior wall, a maximum power of 35 W and a maximum duration of 30 sec were recommended.

Linear ablation lines were only required to treat documented macro-re-entry atrial tachycardias (ATs) and limited to the following targets only:

- LA roof line;
- Mitral Valve (MV) isthmus line;
- LA floor line;
- Cavo tricuspid isthmus (CTI).

A right atrial CTI linear ablation was required in cases with documented typical atrial flutter either before or during the procedure.

Prophylactic ablation of empirical sites was not allowed.

Ablation of spontaneous non-PV triggers was required.

Complex fractionated atrial electrogram (CFAE) ablation was not recommended.

All linear lesions required confirmation of bidirectional conduction block by pacing and/or mapping maneuvers.

The recommended contact force working range was 5-30 grams for both QMODE and QMODE+ modes.

**Procedural Data**

Among the 166 subjects in the PP analysis set, 59 subjects underwent the entire index procedure using QMODE+ mode only and 107 subjects using both QMODE+ & QMODE modes. QMODE+ only was used for PVI in 91 subjects in the PP cohort. Of the 91 subjects, 32 subjects received ablations outside PVs as well using QMODE.

Table 9-11 present the procedural data of the index procedures.

**Table 9 – Summary of Power, Temperature and Impedance Data Per Procedure (Per-Protocol Analysis Set, N=166)**

Description	QMODE+ Only n=59	QMODE+ and QMODE n=107
Average Power (W)	84.333 ± 0.4099	75.893 ± 7.5338
Average Temperature (°C)	39.415 ± 2.5587	39.555 ± 2.6122
Average Impedance Drop (ohms)	9.59 ± 1.696	9.47 ± 2.434

Values in table represent mean ± SD.

**Table 10 – Summary of Ablation Procedure Parameters (Per-Protocol Analysis Set, N=166)**

Variable	Per- Protocol Analysis Set, N=166	QMODE+ Only, N=59	QMODE + and QMODE Combined N=107
<b>Total mapping time (min)</b>			
n	166	59	107
Mean ± SD	10.49 ± 10.593	10.29 ± 7.575	10.61 ± 11.968
<b>Total procedure time (min)</b>			
n	166	59	107
Mean ± SD	144.26 ± 51.101	122.27 ± 28.463	156.38 ± 56.604
<b>Total fluoroscopy duration (min)</b>			
n	165	59	106
Mean ± SD	10.81 ± 10.600	8.65 ± 7.935	12.01 ± 11.690
<b>Total time to achieve PVI (min)</b>			
n	165	58	107
Mean ± SD	25.32 ± 33.948	16.38 ± 20.891	30.17 ± 38.483
<b>Total RF application time for ablating PVs (min)</b>			
n	165	58	107
Mean ± SD	11.66 ± 11.322	7.33 ± 5.273	14.01 ± 12.944
<b>Fluid delivered via the study catheter (mL)</b>			
n	154	52	102
Mean ± SD	546.5 ± 298.55	434.4 ± 190.79	603.7 ± 326.90
<b>Fluid delivered via intravenous fluids (mL)</b>			
n	159	56	103
Mean ± SD	861.2 ± 557.28	744.5 ± 472.39	924.7 ± 590.93

**Table 11 – Number of RF Applications, Average Duration of RF Application and Average Contact Force per Ablation Procedure (Per-Protocol Analysis Set, N=166)**

Variable	QMODE+ Only n=59	QMODE+ and QMODE Combined n=107
<b>Number of RF Applications</b>		
Mean ± SD	100.5 ± 39.36	132.7 ± 41.80
<b>Average Duration of RF Application (sec)*</b>		
Mean ± SD	3.870 ± 0.0393	7.384 ± 4.5099
<b>Average Contact Force (g)</b>		
Mean ± SD	15.800 ± 5.3346	15.194 ± 5.1492

\* Applications with RF energy >0 W and duration of RF application ≥1 sec were included for analysis.

**Ablation Targets**

Table 12 summarizes ablation targets. As required by the protocol, PVs were targeted in all index procedures. Among the PP analysis set, about 40% of the subjects received ablations beyond PVI.

**Table 12 – Ablation Sites Targeted (Per-Protocol Analysis Set, N=166)**


Ablation Sites Targeted	Number of Subjects
<b>Total</b>	166
<b>PVI only</b>	100 / 166 (60.2%)
<b>Non-PV</b>	66 / 166 (39.8%)
Roof line	16 / 66 (24.2%)
Cavo-tricuspid isthmus	46 / 66 (69.7%)
Other linear lesion	15 / 66 (22.7%)
Other AF foci	16 / 66 (24.2%)
Other	3 / 66 (4.5%)

**C.4 - Rhythm Monitoring Compliance**

During the 9-month evaluation period, rhythm monitoring included symptomatic and scheduled asymptomatic TTM transmissions (weekly from Month 3 visit through Month 5 post ablation and monthly starting from Month 6 to Month 12), 12-lead ECG at 3 and 12 months and at 6 months if subject returned for an in-clinic visit, and 24-hour Holter at Month 12.

Figure 1 shows the TTM, 12-lead ECG, and Holter compliance rates at each evaluation time point in the Per-Protocol Analysis Set. The overall TTM and Holter compliance rates were 70.3% and 58.5%, respectively. The ECG compliance rate was 93.3%, 83.7%, and 56.0% at Month 3, 6, and 12, respectively. The composite compliance rate was 97.4% and 91.2% at Month 6 and Month 12, respectively, where at least two rhythm monitoring methods were required per protocol.

**Figure 1 – TTM/ECG/HM Compliance by Time Point (Per-Protocol Analysis Set, N=166)**

**C.5 - Safety Results**

**C.5.1 - Primary Safety Endpoint**

*Primary Analysis*

The primary safety endpoint for this study was defined as the incidence of early onset Primary AEs for subjects undergoing a study ablation procedure. Per study protocol, the primary safety analysis was based on primary safety event rate using the mITT analysis set (N=166) as the primary analysis population.

Table 13 displays a summary for the primary safety endpoint in subjects with non-missing outcomes in the mITT analysis set. The raw incidence was 3.6%. The posterior mean of PAE rate was 4.2%, with a 95% Bayesian credible interval (BCI) ranging from 1.7% to 7.7%. This is well below the performance goal of 14%. The posterior probability that the PAE rate less than the performance goal of 0.14 was close to 1.0000 (i.e., > 0.99995), which exceeded the early success threshold of 0.975. The study met the early success threshold for the primary safety endpoint. There were no changes to the primary safety endpoint data between the early success interim analysis and the final analysis.

**Table 13 – Summary of Primary Safety Endpoint Analysis (mITT Analysis Set, N=167)**

Subjects with Primary Safety Events	6
Subjects without Primary Safety Events	160
Missing	1
Raw Incidence	3.6%
Posterior Mean*	4.2%
95% BCI	(1.7%, 7.7%)
Posterior Probability* that rate < 0.14	1.0000
Threshold to Pass	0.9750
Pass?	Yes

\* Using Beta (1,1) prior distribution

Table 14 presents the descriptive summaries by type of primary adverse events, which were adjudicated by CEC. The PAEs reported in this study were as follows: Three (3) major vascular access complication/bleeding events (reported in 3 subjects, 1.8%), two (2) cardiac tamponade/perforation events (2 subjects, 1.2%), and two (2) phrenic nerve injury/diaphragmatic paralysis events (1 subject, 0.6%). Both cardiac tamponade/perforation events required pericardiocentesis and one of them required blood transfusion; both were adjudicated as procedure related and possibly device related. The three (3) major vascular access complications were adjudicated as procedure related and included one femoral artery pseudoaneurysm requiring thrombin injection, one retroperitoneal bleed and one arteriovenous fistula. The events of phrenic nerve injury/diaphragmatic paralysis were adjudicated as both device and procedure related.

The subject adjudicated by CEC as phrenic nerve injury/diaphragmatic paralysis was re-evaluated at the 12-month visit to determine if the phrenic nerve injury is permanent and was deemed non-assessable by CEC. The sponsor has conservatively counted the adverse event as a PAE.

**Table 14 – Primary Safety Endpoint Outcome – Primary AEs (mITT Analysis Set, N=167)**

Variable	Number of Subjects with Event	Number of Events	Event Rate n/N (%) <sup>[1]</sup>
<b>Primary Adverse Event</b>	6	7	6 / 166 (3.6%)
<b>Type of Primary Adverse Event</b>			
Death (Device or procedure related)	0	0	0 / 166 (0.0%)
Atrio-Esophageal Fistula	0	0	0 / 166 (0.0%)
Cardiac Tamponade/Perforation	2	2	2 / 166 (1.2%)
Myocardial Infarction	0	0	0 / 166 (0.0%)
Stroke/Cardiovascular Accident (CVA)	0	0	0 / 166 (0.0%)
Thromboembolism	0	0	0 / 166 (0.0%)
Transient ischemic Attack (TIA)	0	0	0 / 166 (0.0%)
Phrenic Nerve Injury / Diaphragmatic Paralysis	1	2	1 / 166 (0.6%)
Heart Block	0	0	0 / 166 (0.0%)
Pulmonary Vein Stenosis	0	0	0 / 166 (0.0%)
Pulmonary Edema (Respiratory Insufficiency)	0	0	0 / 166 (0.0%)
Vagal Nerve Injury	0	0	0 / 166 (0.0%)
Pericarditis	0	0	0 / 166 (0.0%)
Major Vascular Access Complication / Bleeding	3	3*	3 / 166 (1.8%)

<sup>[1]</sup> One subject did not experience a PAE and exited the study prior to 3-month follow-up visit; subject was excluded from this analysis.

\* Major vascular access complication/bleeding events included one retroperitoneal hemorrhage, one vascular access site hemorrhage and one pseudoaneurysm.

*Sensitivity Analysis*

The primary safety endpoint was also analyzed in the Safety Analysis set because it is considered by the FDA as a more clinically relevant analysis population for device safety evaluation. Of the 176 subjects in the Safety analysis set, 10 underwent the study ablation procedure and had the study catheter inserted but were identified post-procedure to not meet eligibility criteria for the study. These subjects were not included in the mITT analysis set but were included in the Safety analysis set. No PAEs were reported in these 10 subjects. One (1) subject withdrew from the study 48 days post the index procedure with no PAE, was considered missing and excluded from this analysis.

As summarized in Table 15, the rate of PAE was 3.4% (6/176) in the safety analysis set. The posterior mean of the PAE rate was 3.9%, with an 95% BCI ranging from 1.6% to 7.2%. With 10 more subjects in the analysis, the posterior probability of the PAE rate being less than 0.14 remain greater than 0.9999 which is well above the threshold of 0.975.

**Table 15 – Summary of Primary Safety Endpoint Analysis (Safety Analysis Set, N=177)**

Subjects with Primary Safety Events	6
Subjects without Primary Safety Events	170
Missing	1
Raw Incidence	3.4%
Posterior Mean*	3.9%
95% BCI*	(1.6%, 7.2%)
Posterior Probability* that rate < 0.14	1.0000
Threshold to Pass	0.9750

\* Using Beta (1,1) prior distribution

C.5.2 - Incidence of UADEs

There were no UADEs reported in the study.

C.5.3 - Incidence of Early Onset, Peri-Procedural and Late Onset SAEs

Table 16 Summaries the incidence of early onset, peri-procedural and late onset SAEs in the Safety Analysis set.

**Table 16 – Serious Adverse Events by Timing of Onset (Safety Analysis Set, N=177)**

Serious Adverse Events - by Timing of Onset <sup>[1]</sup>	Number of Events	Number of Subjects with Event	Event Rate <sup>[2]</sup> n/N (%)
<b>Total</b>	55	34	
0 - 7 Days	14	12	12 / 177 (6.8%)
8 - 30 Days	5	4	4 / 177 (2.3%)
> = 31 Days	36	25	25 / 177 (14.1%)

Note: This table summarizes all SAEs. Adjudication for seriousness is based on CEC for PAEs or site for non-PAEs.

<sup>[1]</sup> Including events on and after the day of procedure. Timing of AE based on onset in days post procedure (index or repeat).

<sup>[2]</sup> N denotes number of subjects with follow-up duration corresponding to each time interval.

Table 17 summarizes the SAEs (by causality and body system) occurring within 30 days of a study ablation procedure that were not classified as Primary AEs by protocol definition.

**Table 17 – Serious Non-Primary AEs Occurring within 30 Days Post Index/Repeat Ablation Procedure) by Causality and Body System (Safety Analysis Set, N=177)**

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Overall</b>	11	13	11 / 177 (6.2%)
<b>Causal Relationship to Device</b>	0	0	0 / 177 (0.0%)

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Probable Device Related</b>	0	0	0 / 177 (0.0%)
<b>Possible Device Related</b>	0	0	0 / 177 (0.0%)
<b>Unlikely Device Related</b>	0	0	0 / 177 (0.0%)
<b>Not Device Related</b>	11	13	11 / 177 (6.2%)
Cardiac disorders	7	8	7 / 177 (4.0%)
Atrial fibrillation	3	3	3 / 177 (1.7%)
Atrial flutter	2	2	2 / 177 (1.1%)
Congestive Cardiac failure	1	1	1 / 177 (0.6%)
Cardiac tamponade	1	1	1 / 177 (0.6%)
Pericardial effusion	1	1	1 / 177 (0.6%)
General disorders and administration site conditions	3	3	3 / 177 (1.7%)
Chest pain	2	2	2 / 177 (1.1%)
Pyrexia	1	1	1 / 177 (0.6%)
Respiratory, thoracic and mediastinal disorders	1	1	1 / 177 (0.6%)
Pulmonary embolism	1	1	1 / 177 (0.6%)
Vascular disorders	1	1	1 / 177 (0.6%)
Arteriovenous fistula	1	1	1 / 177 (0.6%)
<b>Causal Relationship to Procedure</b>	1	1	1 / 177 (0.6%)
General disorders and administration site conditions	1	1	1 / 177 (0.6%)
Chest pain	1	1	1 / 177 (0.6%)
<b>Probable Procedure Related</b>	3	5	3 / 177 (1.7%)
Cardiac disorders	2	3	2 / 177 (1.1%)
Congestive Cardiac failure	1	1	1 / 177 (0.6%)
Cardiac tamponade	1	1	1 / 177 (0.6%)
Pericardial effusion	1	1	1 / 177 (0.6%)
Respiratory, thoracic and mediastinal disorders	1	1	1 / 177 (0.6%)
Pulmonary embolism	1	1	1 / 177 (0.6%)
Vascular disorders	1	1	1 / 177 (0.6%)
Arteriovenous fistula	1	1	1 / 177 (0.6%)
<b>Possible Procedure Related</b>	0	0	0 / 177 (0.0%)
<b>Unlikely Procedure Related</b>	0	0	0 / 177 (0.0%)

Relationship to the Device/Procedure by Body System	Number of Subjects with Event	Number of Events <sup>[1]</sup>	Event Rate <sup>[2]</sup>
<b>Not Procedure Related</b>	7	7	7 / 177 (4.0%)
Cardiac disorders	5	5	5 / 177 (2.8%)
Atrial fibrillation	3	3	3 / 177 (1.7%)
Atrial flutter	2	2	2 / 177 (1.1%)
General disorders and administration site conditions	2	2	2 / 177 (1.1%)
Chest pain	1	1	1 / 177 (0.6%)
Pyrexia	1	1	1 / 177 (0.6%)

<sup>[1]</sup> Timing of AE based on onset in days post the procedure (index or repeat).

<sup>[2]</sup> Event rate is the percentage of subjects with the event.

**C.5.4 - Subject Death**

There was one (1) death during the study. The subject was a 66-year-old male with history of PAF and multiple comorbidities (coronary artery disease, hypertension, non-ischemic cardiomyopathy, congestive heart failure, status post-biventricular pacer for cardiac resynchronization therapy). He underwent the study procedure without immediate complications. Five months following the ablation procedure, the subject suffered a witnessed cardiac arrest due to ventricular fibrillation and was resuscitated by immediate cardiopulmonary resuscitation and external defibrillations 10 minutes after the onset of ventricular fibrillation. He subsequently suffered from anoxic brain injury and pneumonia without clinical improvement after therapy. The subject expired about 7 months post ablation. The death was classified as not related to the procedure or the device by the investigator.

**C.5.5 - CT/MRA Sub-Set Analysis**

Table 18 presents the results of the CT/MRA subset analysis. A total of 40 subjects in the CT/MRA subset underwent the CT/MRA imaging of the PVs pre-procedure and 3 months post-procedure. No subjects showed moderate or severe pulmonary vein stenosis in the subset.

**Table 18 – PV Narrowing Post Ablation (CT/MRA Analysis Set, N=40)**

Pulmonary Vein Stenosis	Mild		Moderate	Severe	
	Less than 0%	0-20%	20-50%	50-70%	> 70%
Subject-Level Summary, n/N (%)	1 / 40 (2.5%)	15 / 40 (37.5%)	24 / 40 (60.0%)	0 / 40 (0.0%)	0 / 40 (0.0%)

**C.6 - Effectiveness Results**

**C.6.1 - Primary Effectiveness Endpoint**

*Primary Analysis*

Per study protocol, the primary effectiveness analysis was based on primary effectiveness success using the PP Analysis Set as the primary analysis population.

An early success interim analysis was performed as pre-specified in the protocol with a data extract obtained on June 22nd, 2020 and repeated after the endpoint was met with an additional data extract obtained on August 14th, 2020 to incorporate FDA feedback.

At the planned interim analysis, the posterior mean of event-free rate at Month 12 was 74.5%, with a 95% Bayesian credible interval (BCI) ranging from 66.6% to 81.5%. The posterior probability that the

event-free rate at Month 12 was greater than the performance goal of 50% was close to 1.0000 (i.e., 0.99995), which exceeded the pre-specified threshold of 0.9975. The primary effectiveness endpoint was therefore met at the early success interim analysis.

The analysis of the primary effectiveness endpoint on the August 14, 2020 data extract yielded similar results. These results have been updated with complete follow-up on the full cohort using a beta-binomial model with a non-informative prior. The final primary effectiveness rate is nearly identical to the interim estimate.

Table 19 presents the primary effectiveness analysis based on interim analysis and final analysis in the PP Analysis set.

**Table 19 – Summary of Primary Effectiveness Endpoint Analysis based on the Interim Analysis and Final Analysis (Per-Protocol Analysis Set, N=166)**

	Jun-22, 2020 Data	Aug-14, 2020 Data	Final Data
	Based on 3-piece exponential model	Based on 3-piece exponential model	Based on Beta-Binomial Model
<b>Kaplan-Meier Estimate (39 weeks)</b>	0.765	0.778	0.767
<b>Posterior Mean</b>	0.745	0.755	0.759
<b>95% BCI</b>	(0.666, 0.815)	(0.684, 0.820)	(0.690, 0.823)
<b>Posterior Probability that <math>P_E &gt; 0.50</math></b>	1.0000	1.0000	1.0000
<b>Threshold to Pass</b>	0.9975	0.9975	
<b>Pass</b>	Yes	Yes	

The first failure reason for subjects in the Per-Protocol analysis set is presented in Table 20. Ten (10) subjects with missing effectiveness outcomes are excluded from the summary. Approximately seventy six percent (76.3%, 119/156) of the Per-Protocol population were free from documented atrial tachyarrhythmias and additional failure modes during their effectiveness evaluation period. Documented recurrence of atrial tachyarrhythmias (17.3%, 27/156) and repeat ablation failures (2.6%, 4/156) were the most common reasons for first failure.

**Table 20 – Primary Effectiveness Endpoint by Reason of The First Failure (Per-Protocol Analysis Set, N=166)**

Variable	Number of Subjects with Event	Event Rate n/N (%) <sup>[1]</sup>
Success	119	119/156 (76.3%)
First Failures <sup>[2]</sup>	37	37/156 (23.7%)
Atrial Tachyarrhythmia Recurrence	27	27/156 (17.3%)
AF	25	25/156 (16.0%)
AFL	2	2/156 (1.3%)
AT	0	0/156 (0.0%)
Acute Procedural Failure	3	3/156 (1.9%)
Acute Procedural Failure - Non-study Catheter	3	3/156 (1.9%)
Acute Procedural Failure - Entrance Block	0	0/156 (0.0%)
Repeat Ablation	4	4/156 (2.6%)
AAD Failure	3	3/156 (1.9%)
AAD Failure - New Drug	3	3/156 (1.9%)

Variable	Number of Subjects with Event	Event Rate n/N (%) [1]
AAD Failure - High Dose	0	0/156 (0.0%)
Missing	10	

[1] N is the number of subjects with non-missing primary effectiveness endpoint.

[2] First failures: If a subject has at least 1 failure event, only the earliest failure event is considered and ed by its failure mode.

*Sensitivity Analysis*

Sensitivity to Analysis Sets

Table 21 presents the posterior means of event-free rate at Month 12 using the beta binomial model in the mITT, Safety and Treatment Received Analysis Sets based on the final full follow up data. The Treatment Received Analysis Set was included in the sensitivity analyses per FDA’s feedback because it is considered by the FDA as the most clinically relevant analysis population for device effectiveness evaluation. The results suggested that the success rates in all three Analysis Sets met the effectiveness performance goal. Of note, of the two (2) discontinued subjects who had the study catheter inserted but received no RF ablation from the device, one was not included in the mITT Analysis Set due to study eligibility criteria violation, both were included in the Safety Analysis Set and neither was included in the Treatment Received Analysis Set. Both discontinued subjects were classified as primary effectiveness failures due to use of a non-study catheter for PVI in the respective analysis population.


**Table 21 – Primary Effectiveness Endpoint in mITT, Safety and Treatment Received Analysis Sets**

	mITT Analysis Set (N=167)	Safety Analysis Set (n = 177)	Treatment Received Analysis Set (n = 175)
<b>Posterior Mean</b>	0.755	0.746	0.754
<b>95% BCI</b>	(0.685, 0.818)	(0.678, 0.808)	(0.687, 0.817)


Kaplan-Meier analysis

Figures Figure 2 and Figure 3 present Kaplan-Meier freedom from effectiveness failure for subjects in the Per-Protocol Analysis Set and Treatment Received Analysis Set, respectively. The 12-month Kaplan-Meier estimate of freedom from effectiveness failure was 76.7% (97.5% LCB: 70.1%) in the Per-Protocol Analysis Set and 76.2% (97.5% LCB: 69.7%) in the Treatment Received Analysis Set, both of which are comparable to the posterior mean of 75.9% from the final data.

**Figure 2 – Kaplan-Meier Analysis of Time to First Primary Effectiveness Failures Post Procedure (Per-Protocol Analysis Set, N=166)**

**Figure 3 – Kaplan-Meier Analysis of Time to First Primary Effectiveness Failures Post Procedure (Treatment Received Analysis Set, N = 175)**

Worst-case scenario and tipping point analyses


In the worst-case scenario analysis, the 10 subjects missing primary effectiveness endpoint data were considered as having primary effectiveness failure, which resulted in the primary effectiveness success rate of 71.7% in the Per-Protocol Analysis Set, 71.3% in the mITT analysis set and 71.4% in Treatment Received Analysis Set.

Treating all of the 10 subjects with missing outcomes as effectiveness failures one by one, tipping point analyses using a posterior probability threshold of 0.9775 were performed for the primary effectiveness endpoint in the PP, mITT and Treatment Received Analysis Sets to assess the impact of missing outcomes on the primary effectiveness analysis. The posterior distribution was updated each time and the results showed that no tipping point was identified. Therefore, the primary effectiveness performance goal was met in the tipping point analysis.

*Impact of TTM Compliance*

To better understand the impact of TTM compliance on the primary effectiveness results, additional post-hoc Kaplan-Meier analysis on the time to first primary effectiveness failure was performed in subgroups dichotomized at TTM compliance of 84%, a compliance rate reported in a previous PAF ablation trial. The probability of freedom from primary effectiveness failures at 12 months was similar between subjects with TTM compliance rate of < 84% and those ≥ 84% (Figure 4).

**Figure 4 – Kaplan-Meier Analysis of Time to Primary Effectiveness Failures Post Procedure by TTM Compliance (≥ 84% vs. <84%) (Per-Protocol Analysis Set, N=166)**

C.6.2 - Acute Procedural Success

Acute procedural success was defined as confirmation of entrance block in all PVs. Off the 166 subjects included in the PP analysis set, two (2) subjects were treated with a non-study catheter for PVI due to CoolFlow Pump malfunction which resulted in replacement of the QDOT catheter with a non-study catheter to complete PVI. These two (2) subjects were excluded from the analysis of acute

procedural success. Acute procedural success was achieved in the remaining 164 subjects (100%, 97.5% LCB: 97.8%) in the PP Analysis Set.

Table 22 presents data on PVI after first encirclement at subject level and at ipsilateral PV level, respectively. Among a total of 164 subjects in the PP Analysis Set who received PVI ablation using the study catheter only (two subjects also treated with a non-study catheter for PVI and discussed above were excluded from the analysis), 94 (57.3%) achieved PVI with first encirclement without acute reconnections. For this analysis, both right-sided PVs and left-sided PVs for the subject had to achieve PVI without acute reconnections. At ipsilateral PV level, PVI was achieved in 72.6% (238/328) of the pulmonary veins (ipsilateral veins) after first encirclement without acute reconnections.

**Table 22 – PVI after First Encirclement at Subject Level and at Ipsilateral PV level (Per-Protocol Subjects Who Received PVI Ablation Using the Study Catheter Only, N=164)**

Variable	Number of Subjects Underwent Ablation Procedure	Number of Subjects	Rate
PVI after first encirclement	164	94	94 / 164 (57.3%)
Variable	Number of Targeted Ipsilateral Veins	Number of Targeted Ipsilateral Veins with first encirclement	Rate
PVI after first encirclement	328	238	238 / 328 (72.6%)

C.6.3. - Incidence of repeat ablation procedures


A total of 18 subjects (10.8%) in the Per-Protocol analysis set had one (1) repeat procedure and most (66.7%, 12/18) had their repeat procedure during the evaluation period (Day 91-365). Most subjects (88.9%, 16/18) were treated for AF during the repeat procedure.

C.6.4. - 12-month single procedure success

The 12-Month single procedure success was defined as freedom from documented AF/AFL/AT recurrence (episodes ≥ 30 secs) during the evaluation period after a single ablation and off AADs. Any repeat ablation procedure or Class I or III AAD use for AF/AFL/AT during the evaluation period was deemed effectiveness failure for this endpoint analysis.

Figure 5 presents the Kaplan- Meier analysis of 12-Month single procedure success off Class I/III AAD in the Per Protocol analysis set. The probability of 12-Month single procedure success off Class I/III AADs post blanking was 64.3% in the Per-Protocol analysis set.

**Figure 5 – 12-Month Single Procedure Success off Class I/III AADs Post Blanking (Per Protocol Analysis Set, N=166)**

**C.6.5 - Quality of Life**

The quality of life was measured using the Atrial Fibrillation Effect on Quality-of-Life (AFEQT) at baseline and the 3-, 6-, and 12-month visits.

Figure 6 presents the mean AFEQT composite scores at baseline, 3-, 6-, and 12-month follow up visits in the effectiveness evaluation period, indicating a sustained improvement on the AFEQT scores from 3 months to 12 months post ablation.

**Figure 6 – Mean AFEQT Composite Scores Throughout Study Period (Per-protocol Analysis Set, N=166)**



Figure 7 presents change in total AFEQT scores from baseline at 3, 6, and 12 months. The mean change of AFEQT score from baseline to 3 months (28.77) was sustained at 12 months (31.48) and was larger than a clinical important difference of 19 (Dorian et al., 2013) in the majority of subjects.

**Figure 7 – Atrial Fibrillation Effect on Quality of Life (AFEQT) Questionnaire Overall AFEQT Score (Per-Protocol Analysis Set, N=166)**

**C.7 - Subgroup Analyses**

Subgroup analyses for the primary endpoints were performed, including but not limited to: sex, age, CHA2DS2-VASc score, and ablation mode used for PVI. Per protocol, a test of independence by sex was performed for both primary endpoints, while for other subgroup analyses only descriptive statistics were presented. The PP Analysis Set and Safety Analysis Set were used for the subgroup analyses of the primary effectiveness endpoint and primary safety endpoint, respectively. The primary effectiveness success rates were similar by age (< 65 vs. ≥ 65 years) and CHA2DS2-VASc score (<=2 vs. >2).

i Gender analysis

The primary endpoints were assessed according to gender. Tables Table 23 and Table 24 summarize the primary effectiveness and safety endpoints by sex in the PP and Safety Analysis Sets, respectively. Effectiveness and safety outcomes were similar between male and female subjects.

**Table 23 – Primary Effectiveness Endpoint by Sex (Per-Protocol Analysis Set, N=166)**

Primary Effectiveness Endpoint	Male	Female	p-value <sup>[1]</sup>
Primary Effectiveness Endpoint Success			0.329
n/N	75 / 95	44 / 61	
%	78.9%	72.1%	

<sup>[1]</sup> Chi-square test

**Table 24 – Primary Safety Endpoint by Sex (Safety Analysis Set, N=177)**

Primary Safety Endpoint	Male	Female	p-value <sup>[1]</sup>
Primary AE			0.683
n/N	3 / 106	3 / 70	
%	2.8%	4.3%	

<sup>[1]</sup> Fisher’s exact test

ii Primary Safety by Age (<65 vs. ≥65)

Table 25 summarizes the primary safety endpoint by age in the Safety Analysis Set. No PAEs were observed in subjects aged <65 years. All PAEs occurred in subjects who were aged ≥65 years. The big difference in the primary safety event rate between the two age groups was likely caused by chance in the study in which a very small number of primary safety events occurred. Previous AF ablation studies (Hao, et al., 2012 and Zado, et al., 2008) reported comparable major complication rates between patients 65 years of age or older and those under age 65.

**Table 25 – Primary Safety Endpoint by Age (Safety Analysis Set, N=177)**

Primary Safety Endpoint	Age < 65 years	Age ≥ 65 years
n/N	0 / 82	6 / 94
%	0.0%	6.4%

iii Subgroup analysis of primary endpoints by ablation mode used for PV isolation

Primary effectiveness success was evaluated by ablation mode used for isolating the PVs. Figure 8 presents the Kaplan- Meier analysis of time to primary effectiveness failures by ablation mode used for PVI in the PP Analysis Set. The probabilities of freedom from primary effectiveness failure at 12 months post-procedure were comparable in the 2 groups, QMODE+ for PVI (78.1%) vs. QMODE+ and QMODE for PVI (75.1%).

**Figure 8 – Kaplan-Meier Analysis of Time to First AF/AT/AFL Recurrence Through 12 Months Post Procedure by PV Ablation Mode (Per-Protocol Analysis Set, N=166)**


Table 26 summarizes the primary safety outcome by ablation mode used for PVI in the Safety Analysis Set. The PAE rate was similar between the QMODE+ only group and QMODE+ and QMODE group.

**Table 26 – Primary Safety Endpoint by PV Ablation Mode (Safety Analysis Set, N=177)**

Primary Safety Endpoint	QMODE+ Only	QMODE+ and QMODE
Primary AE		
n/N	3 / 95	3 / 79
% [1]	3.2%	3.8%

[1] The percentage of subjects with PAEs within each group. Subjects without at least 3 months of follow-up are excluded from the primary safety analysis unless they experienced a PAE.

**C.8 - Results in Subjects Who Received Cavo-tricuspid Isthmus Ablation Using Study Catheter**

Among a total of 175 subjects who received RF ablation using the study device in the pivotal trial, 47 underwent CTI ablation using the study device only (n = 45), a non-study catheter only (n = 1), or both (n = 1) during the index procedure per study protocol for typical atrial flutter documented either before or during the procedure. Protocol recommended QMODE mode was used to ablate the CTI in all 33 subjects with ablation mode data available and one of the 33 subjects also received ablations in the CTI using the QMODE+ mode.

Among the 46 subjects who received CTI ablation using the study catheter, 44 (95.7%) had bidirectional CTI block confirmed, and the remaining two subjects either did not achieve bidirectional CTI block (n=1) or had missing information about CTI block (n = 1).

Among the 46 subjects who received CTI ablation using the study catheter, none had a primary AE, three (3) had a total of three (3) serious adverse events and 11 had a total of 12 non-serious adverse events within 30 days after a study procedure. None of these serious adverse events or non-serious adverse events was adjudicated as study device-related; while 9 non-serious adverse events were adjudicated as ablation procedure related.

Among the 46 subjects who received CTI ablation using the study catheter, two (2) had documented atrial flutter recurrence during follow-up. The type of recurrent atrial flutter was atypical atrial flutter in one subject and undetermined in another because it was only documented on TTM.

**D. Study Conclusion**

The results of the Q-EFFICIENCY study demonstrated that there is a reasonable assurance of safety and effectiveness of the QDOT MICRO™ Catheter when used for the treatment of symptomatic drug refractory paroxysmal AF and typical atrial flutter.

## RF ABLATION

For RF ablation the catheter is connected to the CARTO™ 3 System Patient Interface Unit (PIU) via accessory cables, which connect to the RF generator. For setup procedures refer to the User Manual for your CARTO™ 3 Navigation System. For proper RF generator interface, use only a Biosense Webster or compatible interface cable. To complete the electrical circuit, an indifferent electrode must be connected to the indifferent electrode input on the RF generator. Verify that circuit impedance prior to RF ablation is within expected parameters. Verify that the RF generator displays a temperature not above 37°C after the catheter is inserted into the patient and before applying RF power.

## RF GENERATOR OPERATION


Refer to the applicable RF generator manual for proper connection of the catheter to the generator and for detailed instructions as to generator operation for RF ablation.

RF ablation application parameters will vary depending on the ablation site, the specific conditions present in each procedure and the RF generator control circuitry. Based on data obtained from prior animal and clinical studies, recommended RF application parameters are provided below in the “Directions for Use” and in Table 1. Always monitor temperature and impedance rise when using the QDOT MICRO™ Bi-Directional Navigation Catheter.

## HOW SUPPLIED

- The QDOT MICRO™ Bi-Directional Navigation Catheter is supplied STERILE (EtO).
- The catheter is supplied with a choice of five curve configurations: “DD”, “FF”, “JJ”, “DF”, and “FJ” (Figure 9).
- Additional catheter accessory devices are provided separately.

Figure 9 – Curve Types

## PACKAGING

The catheter is secured in a two-piece thermoform tray and placed into a Tyvek/Nylon film pouch, sealed, and placed in a box. Both the pouch and thermoform tray are sterile unless the package is damaged or opened.

## STORAGE

Store in a cool, dry, dark place. Storage temperature should be between 5 and 25°C (41 and 77°F).

## STERILIZATION AND USE-BY DATE

This catheter has been sterilized with ethylene oxide gas. Product and package testing have been conducted to support the “Use By” date printed on the product labels. **DO NOT USE** after the “Use By” date.

This device is packaged and sterilized for single use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/or lead to device failure that in turn may result in patient injury, illness, or death. Also, reprocessing or resterilization of single use devices may create a risk of contamination and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.

## DISPOSAL

Recycle components, or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

## COMPATIBLE EP NAVIGATION SYSTEM


The Biosense Webster QDOT MICRO™ Bi-Directional Navigation Catheter provides location and contact force information only when used with CARTO™ 3 Navigation System. Compatibility with the CARTO™ 3 System has been demonstrated via bench, animal, and clinical testing to confirm that the device is capable of providing accurate location and contact force information when used in accordance with the Instructions for Use.

## DIRECTIONS FOR USE

Please refer to the user manuals for the CARTO™ 3 Navigation System, the irrigation pump, the irrigation tubing, and the generator for instructions on connecting and operating these systems in conjunction with the QDOT MICRO™ Bi-Directional Catheter. Use appropriate Biosense Webster accessory cables to connect the QDOT MICRO™ Bi-Directional Catheter to the appropriate accessory equipment.

1. Using aseptic technique, remove the catheter from the package and place it in a sterile work area. Inspect the catheter carefully for electrode integrity and overall condition.
2. Use the Rocker Lever to deflect the catheter tip (Figure 10). When the lever is turned out of the neutral position, the tip deflects in the same direction as the lever. The amount of deflection is relative to the amount of lever rotation. To straighten the tip, return the Rocker Lever to neutral position.


**Figure 10 – Deflecting the Catheter Tip with the Rocker Lever**

The handle has an adjustable friction control that allows the operator to use the Rocker Lever and deflecting tip in a “free” state or adjust the friction to where the Rocker Lever and tip curve are “locked” in place (Figure 11). This knob is located on the opposite side of the Rocker Lever. Out of the package, the knob will be in the “off” position, which allows the freest movement for the lever and deflecting tip. The amount of friction increases as the Friction Control Knob is rotated clockwise until it reaches the fully “on” position.

**Figure 11 – Friction Control Knob**

Clockwise rotation from the “off” position increases the friction within the Deflection Mechanism.

Counter-Clockwise rotation from the “on” position decreases the friction within the Deflection Mechanism.

3. Follow standard practice for vessel puncture, guidewire insertion, and guiding sheath use and aspiration per its instructions for use.
4. To verify compatibility between the sheath and catheter, advance the catheter through the sheath prior to insertion. Any sheath < 8.5 F is contraindicated. For compatible sheaths, contact Customer Support or your Biosense Webster representative.
5. Connect the irrigation pump tubing to a room temperature, heparinized (1 IU heparin/ml) normal saline bag using standard safe hospital practices. Open the stopcock on the end of the tubing set and fill the tubing set as slowly as possible. Remove any trapped air and then close the stopcock.
6. Load the irrigation tubing into the pump. Open the stopcock and flush per the irrigation pump instructions until the air is expelled through the open end of the tubing.
7. Connect the irrigation tubing to the Luer fitting of the catheter.
8. Flush the catheter and tubing to ensure purging of trapped air bubbles and to verify irrigation through the tip electrode.
9. Start continuous irrigation with the irrigation pump at the low flow rate.
10. Connect the catheter to the CARTO™ System Patient Interface Unit (PIU) with the appropriate Biosense Webster cables. Connect the PIU to the generator and to the appropriate recording and mapping systems, including the CARTO™ 3 Navigation System, with appropriate interface cables. Use only Biosense Webster interface cables. To complete the electrical circuit, connect an indifferent electrode to the indifferent electrode input on the generator.
11. Insert the QDOT MICRO™ Bi-Directional Catheter via the entrance site.
12. Advance the catheter to the desired area. Use both direct imaging guidance (such as fluoroscopy or ultrasound) and electrograms (EGM) to aid in proper positioning.
13. In order to achieve optimal force reading accuracy and stability, allow the catheter to warm up for 2 minutes after connection to the CARTO™ 3 System, prior to use of the force feedback feature.
14. Zero the contact force reading following insertion into the patient. The tip electrode and the two distal ring electrodes on the catheter tip must be outside of the sheath so that the force sensor is inside the body. Ensure the catheter tip is not in contact with tissue by evaluating the location on fluoroscopy, the CARTO™ System or other direct imaging guidance along with the EGM amplitude and catheter movement. Variations in the force reading at the same rate as the cardiac or respiration cycle may indicate contact with cardiac structures. Once these markers indicate the tip is not in contact, the reading can be zeroed. Refer to the user manual for the CARTO™ 3 System for instructions on how to zero the contact force reading.
15. Zero the contact force reading when moving the catheter from one chamber of the heart to another or upon reinsertion.
16. Verify that the “QDOT MICRO” option is selected on the RF generator. When this option is chosen, the RF generator defaults to the safety parameters established for the QDOT MICRO™ Bi-Directional Navigation Catheter.
17. Table 1 below outlines the recommended ablation parameter settings when using the QDOT MICRO™ Bi-Directional Navigation Catheter for ablation in the atria:

**Table 1: Recommended Ablation Parameter Settings**

Power Mode	QMODE*		QMODE+
Power (W)	25-35	36-50	90
Target Electrode Temperature (°C)	50	50	60
RF application time (s)	Up to 60	Up to 60	4
Nominal Irrigation flow rate (ml/min)**	4***	15***	8

\* RF applications on the left atrial posterior wall using QMODE should not exceed 35 W in power and 30 sec in duration.

\*\* A minimum flow rate of 2 mL/min during mapping is recommended.

\*\*\* The irrigation flow rate is set by the generator and automatically adjusted between 4ml/min or 15ml/min to reach and maintain the set maximum power within the target electrode temperature.

- For the treatment typical atrial flutter, the QMODE is recommended.
- For the treatment of paroxysmal atrial fibrillation, the QMODE+ is recommended as the primary ablation mode for creating pulmonary vein isolation lesions and the QMODE is recommended for applications outside the PV ostia and for touch-up ablations.
- The QMODE power setting allows a maximum power up to 50 W, and the clinical operator may select any combination of power and RF application time, guided by the clinical judgment and expertise of the clinical operator as well as the anatomical area of interest.
- The QMODE+ power setting allows a maximum power up to 90 W, with a maximum RF application time up to 4 seconds. Use caution when overlapping RF ablation points, particularly on the posterior wall. Maintain appropriate contact force within the recommended range. This combination of power and RF application time is optimized for atrial ablation and may be used independently or in conjunction with the QMODE power setting.
- For ablations with any power setting, the clinical operator should monitor commonly used ablation effectiveness parameters like electrogram amplitude reduction and/or impedance drop.

Contact Force (CF) ranges during ablation:

The QFFICIENCY IDE study recommended CF settings to operators for use with the study catheter in both QMODE and QMODE+ as follows:

- CF Working Range: 5-30g

During the study, the overall average CF recorded during a study ablation procedure for all 166 subjects in the Per-Protocol Analysis set was  $15.4 \pm 5.2$  grams.

18. Do not use this catheter without irrigation flow. Do not alter the irrigation flow rate.
19. Monitor the catheter tip temperature response throughout the procedure to ensure adequate irrigation. If temperature increases very rapidly during RF application, power delivery should be interrupted. The irrigation system must be rechecked prior to restarting RF application. Note: The displayed temperature represents the temperature of the electrode only, not the temperature of the tissue.
20. In the event of a generator cutoff (impedance or temperature), the catheter must be withdrawn and the tip electrode inspected for coagulum before RF energy is reapplied. To remove any coagulum, if present, a sterile gauze pad dampened with sterile saline may be used to gently wipe the tip section clean. Do not scrub or twist the tip electrode as damage to the tip electrode bond may occur and loosen the tip electrode, or damage may also occur to the contact force sensor and affect measurement accuracy. Prior to reinsertion, flush the tip to ensure that the irrigation holes are not plugged.

If irrigation hole occlusion occurs:

- a. Fill a 1-2 ml syringe with sterile saline and attach it to the stopcock on the end of the tubing set.
- b. Inject the saline from the syringe into the catheter. A uniform flow of fluid should be visible from the tip of the catheter.
- c. Repeat steps a and b, if necessary, until the holes are cleared.
- d. Flush the catheter and tubing per standard technique to ensure purging of trapped air bubbles and to verify that the irrigation holes are patent.
- e. The catheter can now be introduced into the patient.
- f. Zero the catheter following reinsertion into the patient.

**WARNING: Do not continue use of the catheter if it is still occluded or if it is not functioning properly.**

**DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY**

**THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.**

**WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.**

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

---

**NOTES**

**Biosense Webster, Inc.**

31 Technology Drive, Suite 200  
Irvine, California 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010  
Fax: +1-909-468-2905  
[www.biosensewebster.com](http://www.biosensewebster.com)

BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, QDOT MICRO, THERMOCOOL SMARTTOUCH, THERMOCOOL, CARTO SMARTTOUCH, and NAVISTAR are trademarks of Biosense Webster, Inc.

The third-party trademarks used herein are trademarks of their respective owners.

Patent information is available at [www.biosensewebster.com/virtualpatentmarking](http://www.biosensewebster.com/virtualpatentmarking).

Revised: 2023-01  
© Biosense Webster, Inc. 2010-2023

# nGEN™ Generator User Manual

## Catheter Operation

### QDOT MICRO™ Catheter

**Caution: Federal (USA) law restricts  
this device to sale by or on the  
order of a licensed healthcare  
practitioner.**

EML-1384-02-IFU-02A

## Table of Contents

<b>1</b>	<b>DESCRIPTION OF THIS OPERATION MANUAL</b>	<b>1</b>
<b>2</b>	<b>ABLATION MODES FOR THE QDOT MICRO™ CATHETER</b>	<b>1</b>
<b>3</b>	<b>WARNINGS</b>	<b>1</b>
<b>4</b>	<b>GENERATOR COMPONENTS</b>	<b>1</b>
<b>5</b>	<b>CATHETER SETUP AND CONNECTIONS</b>	<b>2</b>
<b>6</b>	<b>TYPICAL WORKFLOW</b>	<b>3</b>
<b>7</b>	<b>GENERATOR SCREENS</b>	<b>4</b>
7.1	Connections Screen: Catheter Dialog Box .....	4
7.2	Presets Screen: QDOT MICRO™ Catheter.....	5
7.2.1	Preset Overview	6
7.2.2	Preset Settings	6
7.2.3	Selecting a Preset	9
7.2.4	Changing a Preset: Basic Instructions	10
7.2.5	Creating a New Preset	10
7.2.6	Deleting a Preset	11
7.3	Ablation Screen: QDOT MICRO™ Catheter.....	11
7.3.1	Features on the Ablation Screen	12
7.4	Irrigation Control Panel .....	13
7.5	Post Screen: QDOT MICRO™ Catheter.....	16
<b>8</b>	<b>TECHNICAL DATA</b>	<b>20</b>
8.1	Output Specifications: QDOT MICRO™ Catheter .....	20
<b>9</b>	<b>TROUBLESHOOTING</b>	<b>21</b>

## Table of Figures

Figure 1 - QDOT MICRO™ Catheter to CARTO™ System.....	2
Figure 2 - Connections Screen: QDOT MICRO™ Catheter .....	4
Figure 3 - Presets Screen .....	5
Figure 4 – QDOT MICRO™ Catheter, QMODE™ .....	6
Figure 5 – QDOT MICRO™ Catheter, QMODE+™ .....	7
Figure 6 - Ablation Screen.....	12
Figure 7 - Irrigation Control Panel .....	13
Figure 8 - Post Screen: Ablation Procedure Summary .....	16
Figure 9 - Post Screen: Ablation Session Information – Data.....	17
Figure 10 - Post Screen: Ablation Session Information – Graph .....	17
Figure 11 - Output Specifications .....	20
Figure 12 – Output Specifications - QMODE™ and QMODE+™ Ablation Modes .....	20

## Table of Tables

Table 1 - Devices .....	2
Table 2 - Cables .....	2
Table 3 - Catheter Dialog Box .....	4
Table 4 - Cables .....	5
Table 5 - Features on the Presets Screen .....	5
Table 6 - Preset Settings .....	7
Table 7 - Features on the Ablation Screen .....	12
Table 8 - Irrigation Control Buttons .....	14
Table 9 - Post Screen Fields .....	18
Table 10 - QDOT MICRO™ Catheter Alarm Messages .....	21

# 1 Description of This Operation Manual

This manual contains information specifically applicable to use of the nGEN™ Generator with the Biosense Webster QDOT MICRO™ Catheter. The *nGEN™ Generator User Manual: General Features and Functions* contains information about the nGEN™ Generator that applies regardless of the catheters used with the generator.

## 2 Ablation Modes for the QDOT MICRO™ Catheter

The nGEN™ Generator delivers RF energy using two temperature-controlled ablation modes through the QDOT MICRO™ Catheter. QMODE™ is used for ablation power settings up to 50W; QMODE+™ is used for ablation power settings up to 90W.

QMODE™ uses irrigation flow to maximize power delivery without exceeding target temperature during an ablation. Select a maximum power and a target temperature on the nGEN™ Generator. For maximum power settings greater than 35 W, also select a maximum low flow temperature. Depending on the power setting, the generator changes the pump from high to low flow or from low to high flow to reach and maintain the set maximum power to ensure the tip is within the target temperature range.

- For maximum power settings up to 35 W (the pump starts at low flow): The generator changes the pump from low to high flow to maintain the set maximum power when the temperature exceeds the set target temperature during ablation and the power begins to drop.
- For maximum power settings between 36W and 50W (the pump starts at high flow): The generator changes the pump from high to low flow when the temperature drops below or does not rise above the set maximum low flow temperature during ablation. The change to low flow allows the temperature to rise closer to the target temperature.

QMODE+™ delivers high power in a short ablation session. Select a target temperature and maximum power on the nGEN™ Generator. If the target temperature is reached before the set duration expires, the generator modulates the power to maintain the target temperature.

## 3 Warnings

See the *nGEN™ Generator User Manual: General Features and Functions* and the QDOT MICRO™ Catheter instructions for use for applicable warnings.


## 4 Generator Components

For information about the components provided as part of the nGEN™ Generator, see the *nGEN™ Generator User Manual: General Features and Functions*.

## 5 Catheter Setup and Connections

This section describes how to connect the catheter, and identifies the devices and cables shown in the illustration.

**Figure 1 - QDOT MICRO™ Catheter to CARTO™ System**

**Table 1 - Devices**

<b>1</b>	See Note	QDOT MICRO™ Catheter
<b>2</b>	D140102	TX eco Cable
<b>3</b>	FG540000	CARTO™ System
<b>4</b>	M685301	nGEN™ Console

**Note:** Contact a Biosense Webster representative for the QDOT MICRO™ Catheter part number.

**Table 2 - Cables**

<b>A</b>	D135703	TX eco EXT Cable
<b>B</b>	M581007	Cable TX, nGEN™ Console to CARTO™ System, 3 m (10 ft)

## 6 Typical Workflow

### Steps

1. Physically connect all required devices, and selected optional devices, to the generator.
2. Turn on the power supply unit.
3. Turn on the console and monitor(s).
4. At the **Connections** screen, review the device connections.
5. Physically connect all applicable devices to the patient.
6. At the **Presets** screen, confirm the **QDOT MICRO™** selection from the **Catheter** menu, then select or create a preset.
7. At the **Ablation** screen, review the settings. Position the catheter at the first ablation site.
8. Start ablation by pressing the **START** button on either the Primary Monitor or the Secondary Monitor, or by pressing the pedal connected to the Primary Monitor, the Secondary Monitor, or the Console. Press the **STOP** button or release the pedal to stop ablation.


## 7 Generator Screens

For information about the screens that appear regardless of catheter type, see the *nGEN™ Generator User Manual: General Features and Functions*.

### 7.1 Connections Screen: Catheter Dialog Box

The catheter icon on the **Connections** screen indicates only whether a compatible catheter is connected to the nGEN™ Generator. See Figure 1 for details on how to connect the QDOT MICRO™ Catheter to the nGEN™ Generator via the TX eco Cable and CARTO™ System. For information about the **Connections** screen in general, and all other device dialog boxes, see *Section 13.3 Connections Screen, nGEN™ Generator User Manual: General Features and Functions*.

**Figure 2 - Connections Screen: QDOT MICRO™ Catheter**

Press the catheter icon on the **Connections** screen. The **Catheter Type** dialog box appears.

**Table 3 - Catheter Dialog Box**

Button or Field	Description
	<ul style="list-style-type: none"> <li>The diagram shows the connection of <b>Navigation Catheters</b> to the <b>nGEN™ Console</b> and the <b>CARTO™ PIU</b>. See <i>Section 5, Catheter Setup and Connections</i> for more information.</li> </ul>

**Table 4 - Cables**

<b>A</b>	M581007	Cable TX, nGEN™ Console to CARTO™ System, 3 m (10 ft)
<b>B</b>	D135703	TX eco EXT Cable
Not shown	D140102	TX eco Cable

## 7.2 Presets Screen: QDOT MICRO™ Catheter

**Figure 3 - Presets Screen**

**Table 5 - Features on the Presets Screen**

<b>1</b>	<b>Undo</b> button
<b>2</b>	<b>New Preset</b> button
<b>3</b>	<b>Delete</b> button
<b>4</b>	<b>Save Preset</b> button
<b>5</b>	<b>Ablation Settings</b> area
<b>6</b>	<b>Irrigation Settings</b> area
<b>7</b>	<b>Preset</b> Menu
<b>8</b>	<b>Catheter</b> Menu

## 7.2.1 Preset Overview


- When **QDOT MICRO™** is selected from the **Catheter Menu**, the names of the available presets are displayed in the **Preset Menu**. When a preset name is selected, the settings for that preset are displayed in the **Ablation Settings** area and the **Irrigation Settings** area of the screen.
- If the list is longer than the screen, the grey arrow in the **Next** button becomes white. Press the **Next** button to display the next page of presets. Press the **Previous** button to display the previous page of presets.
- The generator comes with one **QDOT MICRO™** preset named **Default**. Additional presets can be created, or presets can be imported from a USB flash drive.
- Values on the **Presets** screen can be changed. Press the value and the background will become orange. Use the **Control** knob, **Plus** button, or **Minus** button to change the value. The changes are automatically updated.

**Caution:** Do not press the control knob. Doing so may inadvertently change the **Presets** values.


## 7.2.2 Preset Settings

The screens below show the fields associated with the QDOT MICRO™ Catheter and are followed by a table that explains each field. The values in the following figures are for illustration only.

**Figure 4 – QDOT MICRO™ Catheter, QMODE™**

**Figure 5 – QDOT MICRO™ Catheter, QMODE+™**

**Table 6 - Preset Settings**

Field	Setting
<b>Ablation Control</b>	
	<ul style="list-style-type: none"> <li>The ablation control displays the ablation mode.</li> <li>Press <b>QMODE</b> for QMODE™ ablation mode or press <b>QMODE+</b> for QMODE+™ ablation mode. The selected ablation mode becomes orange.</li> </ul>
<b>Time</b>	
	<ul style="list-style-type: none"> <li>Change the <b>Time</b> value to set the duration of the ablation session.</li> <li>The value set in this field appears on the <b>Ablation</b> screen.</li> </ul>
<b>Power</b>	
	<ul style="list-style-type: none"> <li>Change the <b>Max.</b> value to the maximum power the generator is to deliver during the ablation.</li> <li>The value set in this field appears on the <b>Ablation</b> screen.</li> </ul>

**Table 6 - Preset Settings**

Field	Setting
-------	---------

### Temperature

- Change the **Target** value to the temperature for the generator to reach and maintain during ablation. The value set in this field appears on the **Ablation** screen.
- Change the **Cutoff** value to set the temperature for the generator to stop ablation.
- The **Target** value must be at least 5°C lower than the **Cutoff** value. When the **Target** value is exactly 5°C lower than the **Cutoff** value, the generator will not allow the **Target** value to be increased unless the **Cutoff** value is increased first to maintain a difference of at least 5°C.



### Impedance

#### Range

- Press the **On/Off** toggle button to turn on **Range** for the generator to stop the ablation and issue an alarm when the maximum or minimum impedance is reached.
- Change the **Max.** and **Min.** values to set the maximum or minimum impedance. It is recommended that impedance during ablation not exceed 250 ohms. The generator will also stop the ablation and issue an alarm if the impedance value becomes infinite (for example, the catheter becomes disconnected from the generator) or the value becomes zero (for example, an internal short).
- When **Range** is **Off**, it is not possible to update the **Max.** and **Min.** impedance values.
- Press the **On/Off** toggle button to turn on **Delta** for the generator to stop the ablation and issue an alarm when the impedance change exceeds the **Max.** value.

**Table 6 - Preset Settings**

Field	Setting
<b>Irrigation</b>	
	<p><b>Auto Pump Control</b></p> <ul style="list-style-type: none"> <li>• <b>Auto Pump Control</b> indicates that the pump is being controlled by the generator (see <i>Section 13.3.6, nGEN™ Pump, nGEN™ Generator User Manual: General Features and Functions</i>).</li> </ul>
 	<p><b>High Flow</b></p> <ul style="list-style-type: none"> <li>• The <b>High Flow</b> values show the irrigation flow rates during ablation and cannot be changed.</li> <li>• When <b>QMODE™</b> is selected, the generator changes between the <b>Max.</b> and <b>Min.</b> flow rates.</li> <li>• Press the <b>Pre-Abl.</b> value to set the time to start the high flow before the ablation starts.</li> <li>• Press the <b>Post-Abl.</b> value to set the time to maintain the high flow after the ablation ends.</li> </ul>
	<p><b>Low Flow</b></p> <ul style="list-style-type: none"> <li>• The <b>Low Flow</b> value shows the irrigation flow rate before the <b>Pre-Abl.</b> time and after <b>Post-Abl. time</b>. This value cannot be changed.</li> </ul>
	<p><b>Max. Low Flow Temp. (QMODE™ only)</b></p> <ul style="list-style-type: none"> <li>• For <b>Max. Power</b> settings between 36W and 50W (the pump starts at high flow), select the <b>Max. Low Flow Temp.</b> for the generator to change the pump from high to low flow when the temperature drops below or does not rise above this value during ablation.</li> </ul>
 	<p><b>Bag Size</b></p> <ul style="list-style-type: none"> <li>• Change the <b>Bag Size:</b> value to set the size of the irrigation solution bag. The value can be changed in increments of 250 ml. See <i>Section 12.7, Changing the Irrigation Solution Bag, nGEN™ Generator User Manual: General Features and Functions</i>.</li> <li>• Press the <b>New Bag</b> button to reset the fluid volume calculator.</li> </ul>
	<p><b>Fluid Warning</b></p> <ul style="list-style-type: none"> <li>• Press the <b>On/Off</b> toggle button to turn on <b>Fluid Warning</b> for the generator to issue the <b>Low Fluid Warning</b> when the <b>Current Fluid Volume</b> is below the set <b>Level</b> value.</li> </ul>

### 7.2.3 Selecting a Preset

#### Steps

1. Press the **Presets** tab.
2. Verify that **QDOT MICRO™** is selected.
3. Select a preset name. The preset name background turns orange when selected.

## 7.2.4 Changing a Preset: Basic Instructions

### WARNING

Modifying settings can significantly affect the performance of the nGEN™ Generator. Use clinical judgment and consider patient conditions when modifying the settings.

#### Steps

1. Press the **Presets** tab.
2. Verify that **QDOT MICRO™** is selected.
3. Select a preset name.
4. Change the settings as desired. When a setting is changed, an asterisk appears after the preset name to indicate that a change was made. If the **Default** preset is selected, settings can be changed but not saved, and the **Save Preset** button will be disabled.
5. Press the **Save Preset** button. The **Save Preset** dialog box appears and displays this message:  
*“Are you sure you want to save changes to this preset?”*
6. Press **Yes** or **No**. If **Yes** is pressed, the changes are saved and the asterisk disappears.
7. If a change does not need to be saved, press the **Undo** button. The **Undo Changes** dialog box appears and display this message:  
*“Are you sure you want to undo the changes to this preset?”*
8. Press **Yes** or **No**.

## 7.2.5 Creating a New Preset

#### Steps

1. Press the **Presets** tab.
2. Verify that **QDOT MICRO™** is selected.
3. Select a preset with settings similar to those that will be created for the new preset.
4. Change the settings as desired. When a setting changes, an asterisk appears after the preset name to indicate that a change was made to the selected preset.
5. Press the **New Preset** button. The keyboard will appear.
6. Preset names can have a maximum of 15 characters. Enter a name for the new preset. If that name already exists, the **Save Preset** dialog box appears and displays this message:  
*“[Preset name] already exists. Please choose as different name.”*  
Use the Backspace key to delete the name, then enter a new one.
7. Press **Save**.

## 7.2.6 Deleting a Preset

### Steps

1. Press the **Presets** tab.
2. Select the preset to delete. If the **Default** preset is selected, the **Delete** button will be disabled.
3. Press the **Delete** button. The **Delete Preset** dialog box appears and displays this message:  
*“Are you sure you want to delete this preset?”*
4. Press **Yes**. After a preset is deleted, the deletion cannot be undone.

## 7.3 Ablation Screen: QDOT MICRO™ Catheter

### Steps

1. After the generator is started, the **Ablation** screen appears.
2. Connect the **QDOT MICRO™** Catheter and the **Ablation** screen displays the settings from the **QDOT MICRO™ Default** preset.
3. If a different preset was previously selected, the **QDOT MICRO™ Connected** dialog box appears and displays this message:  
*“The preset has changed to the Default preset for QDOT MICRO™.”*
4. Press **OK**.

### Notes

- The time, power, and temperature settings on the **Ablation** screen can be changed before and during ablation.
- Changes made in the **Ablation** screen will appear on the **Presets** screen. Changes will be lost if not saved before a procedure ends, or the generator is shut down.
- Data from the ablation sessions is saved in log files that can be exported to a USB flash drive.

### WARNING

Modifying settings can significantly affect the performance of the nGEN™ Generator. Use clinical judgment and consider patient conditions when modifying the settings.


### 7.3.1 Features on the Ablation Screen

The **Ablation** screen displays temperature, impedance, power, and time values as numbers and graphs.


**Figure 6 - Ablation Screen**

**Table 7 - Features on the Ablation Screen**

Feature	Description
 	<ul style="list-style-type: none"> <li>Press the <b>QMODE™</b> or <b>QMODE+™</b> to select the ablation mode. The selected mode turns orange.</li> </ul>
<p><b>Time</b></p> 	<ul style="list-style-type: none"> <li>Two dashes are displayed before and after ablation.</li> <li>The large white number shows the elapsed ablation time.</li> <li>When using irrigated catheters, a minus sign before the number indicates the <b>Pre-Abl.</b> high flow time and counts to <b>0</b> before the ablation starts. A plus sign before the number indicates the <b>Post-Abl.</b> high flow time after the ablation ends.</li> <li>The small white number inside a rectangle shows the ablation duration. This value can be changed before or during ablation by pressing the number. The background will become white. Use the <b>Control</b> knob, <b>Plus</b> button, or <b>Minus</b> button to change the value.</li> </ul>
<p><b>Power</b></p> 	<ul style="list-style-type: none"> <li>Two dashes are displayed before and after ablation.</li> <li>The large yellow number shows the real-time power.</li> <li>The small yellow number inside a rectangle shows the <b>Max.</b> power. This value can be changed before or during ablation by pressing the number. The background will become yellow. Use the <b>Control</b> knob, <b>Plus</b> button, or <b>Minus</b> button to change the value.</li> <li>On the graph, the solid yellow line represents the real-time power and the dashed yellow line represents the <b>Max. power</b>.</li> </ul>

**Table 7 - Features on the Ablation Screen**


Feature	Description
<p><b>Temperature</b></p>  <p>(QMODE™)</p>	<ul style="list-style-type: none"> <li>The large orange number shows the real-time temperature.</li> <li>The small orange number inside a rectangle shows the <b>Target</b> temperature. This value can be changed before or during ablation by pressing the number. The background will become orange. Use the <b>Control</b> knob, <b>Plus</b> button, or <b>Minus</b> button to change the value.</li> <li>On the graph, the solid orange line represents the real-time temperature and the dashed orange line represents the <b>Target</b> temperature.</li> <li>On the graph, if <b>Flow Change Indicator</b> is selected from the <b>Ablation Display</b> screen, the blue arrows and their associated blue dots show when the nGEN™ Pump changes between the <b>Min.</b> and <b>Max. High Flow</b> settings in QMODE™ ablation mode. See <i>Section 13.2.1, Ablation Display, nGEN™ Generator User Manual: General Features and Functions.</i></li> </ul>
<p><b>Impedance</b></p> 	<ul style="list-style-type: none"> <li>The green number shows the real-time impedance.</li> <li>On the graph, a solid green line represents the real-time impedance and the dashed green line represents the impedance value when the ablation session starts.</li> </ul>

**Note:** When a setting is changed, an asterisk appears after the preset name to indicate that a change was made. Setting changes made on the **Ablation** screen also appear on the **Presets** screen.

## 7.4 Irrigation Control Panel

The **Irrigation Control** panel is at the lower right of any screen. Irrigation settings can be changed before ablation, but all buttons are disabled during ablation.

**Figure 7 - Irrigation Control Panel**

4b

The following buttons are displayed.

Press the **Collapse**  button to see fewer buttons or the **Expand**  button to see more buttons.


If a **Bubble** or **Fluid Level Low** alarm occurs, the display expands automatically to show the alarm.

**Table 8 - Irrigation Control Buttons**

1	<b>Flow Rate</b>	<ul style="list-style-type: none"> <li>The <b>Flow Rate</b> button displays the current flow rate. This value is read-only and cannot be changed.</li> <li>When the generator communicates with the pump and there is no irrigation flow, the <b>Flow Rate</b> button displays zero.</li> <li>When the generator cannot communicate with the pump and there is no irrigation flow, the <b>Flow Rate</b> button displays dashed lines.</li> <li>If <b>Manual Pump Control</b> selected, the generator will issue an alarm: <i>"Pump is in manual mode."</i></li> <li>Return to <b>Auto Pump Control</b>.</li> </ul>
	2 ml/min	
	0 ml/min	
	--- ml/min	
	<b>Manual Pump Control</b>	
	M	
2a	<b>Low Flow</b>	<ul style="list-style-type: none"> <li>The <b>Low Flow</b> button indicates when the pump provides irrigation at the set low flow rate.</li> </ul>
		
2b	<b>High Flow</b>	<ul style="list-style-type: none"> <li>The <b>High Flow</b> button indicates when the pump provides irrigation at the set high flow rate.</li> <li>The generator automatically changes from low flow rate to high flow rate when the ablation session starts and from high flow rate to low flow rate when the ablation session ends.</li> </ul>
		
3	<b>Collapse</b>	<ul style="list-style-type: none"> <li>Press the <b>Collapse</b> button to display only the <b>Flow Rate</b> and the <b>Low Flow</b> or <b>High Flow</b> button.</li> </ul>
		
	<b>Expand</b>	<ul style="list-style-type: none"> <li>Press the <b>Expand</b> button to display all irrigation buttons.</li> </ul>
		
		<p><b>Note:</b> The <b>Low Flow</b> button and the <b>High Flow</b> button can also be pressed to collapse and expand the display of buttons.</p>

**Table 8 - Irrigation Control Buttons**

4a	<b>Start Low Flow</b>	<ul style="list-style-type: none"> <li>If this button is not visible, press the <b>Expand</b> button.</li> <li>Press the <b>Start Low Flow</b> button to start low flow on the pump. After the low flow starts, the button changes to the <b>Start High Flow</b> button.</li> </ul>
		
4b	<b>Start High Flow</b>	<ul style="list-style-type: none"> <li>Press the <b>Start High Flow</b> button to start high flow on the pump. After the high flow starts, the button changes to the <b>Start Low Flow</b> button.</li> </ul>
		
5	<b>Stop Flow</b>	<ul style="list-style-type: none"> <li>If this button is not visible, press the <b>Expand</b> button.</li> <li>Press the <b>Stop Flow</b> button to stop the irrigation flow.</li> </ul>
		
6	<b>Flush</b>	<p><b>WARNING</b> Do not flush the catheter when it is in the patient.</p> <ul style="list-style-type: none"> <li>If this button is not visible, press the <b>Expand</b> button.</li> <li>Press the <b>Stop Flow</b> button to enable the <b>Flush</b> button.</li> <li>Press the <b>Flush</b> button to begin flushing the irrigation tubing. The button becomes orange and this message appears:  <i>“To flush, press and hold the <b>Control</b> knob on the monitor. Release the knob to stop flushing.”</i></li> <li>Press and hold the <b>Control</b> knob on the monitor to flush the irrigation tubing.</li> <li>Release the <b>Control</b> knob to stop flushing.</li> <li>Flushing cannot be performed during ablation.</li> <li>When it is not possible to flush, the button is disabled (grey).</li> </ul>
		
7	<b>Current Fluid Volume</b>	<ul style="list-style-type: none"> <li>The <b>Current Fluid Volume</b> displays the amount of fluid in the irrigation solution bag.</li> </ul>
		
8a	<b>Reset Bag</b>	<ul style="list-style-type: none"> <li>If the <b>Reset Bag</b> button is not visible, press the <b>Expand</b> button.</li> <li>Press the <b>Reset Bag</b> button to change the irrigation solution bag. The button changes to orange and the <b>Bag Size:</b> and <b>New Bag</b> button appear. Change the <b>Bag Size:</b> value to the set the size of the irrigation solution bag in increments of 250 ml. Press the <b>New Bag</b> button to reset the fluid volume calculator. See <i>Section 12.7, Changing the Irrigation Fluid Bag, nGEN™ Generator User Manual: General Features and Functions.</i></li> </ul>
	 	
8b	<b>Low Fluid Warning</b>	<ul style="list-style-type: none"> <li>The <b>Low Fluid Warning</b> button appears when <b>Current Fluid Volume</b> drops below the set <b>Level</b> value. If on the <b>Presets</b> screen, the <b>Fluid Warning</b> field is selected <b>On</b>, an alarm message will appear.</li> </ul>
		
		<p><b>Note:</b> The <b>8a</b> and <b>8b</b> buttons appear together on the figure above for illustration; the <b>Irrigation Control</b> panel will display either the <b>8a</b> or <b>8b</b> icon, but not both at the same time.</p>

**Table 8 - Irrigation Control Buttons**

**9 Bubbles**

- The **Bubbles** alarm displays when the pump detects bubbles in the irrigation tubing set.
- If a bubble error occurs, the irrigation tubing must be flushed. Remove the catheter from the patient or close the stopcock to the catheter, then press the **Flush** button. The **Flush** dialog box appears and displays this message:  
*“To flush, press and hold the **Control** knob on the monitor. Release the knob to stop flushing.”*
- Press and hold the knob to flush and release the **Control** knob after the **Bubbles** alarm disappears. Flushing cannot be performed during ablation.

**7.5 Post Screen: QDOT MICRO™ Catheter**

The **Post** screen displays data about the entire ablation procedure or a specific ablation session. The data remains available on the **Post** screen until a selection is made from the **Shut Down** menu. The data is then deleted from the **Post** screen but is saved in a log file.

**Figure 8 - Post Screen: Ablation Procedure Summary**

**Figure 9 - Post Screen: Ablation Session Information – Data**

**Figure 10 - Post Screen: Ablation Session Information – Graph**

**Steps**



1. Press the **Post** tab.
2. Press the **Procedure Summary** button or an **Ablation Session** button.

**Table 9 - Post Screen Fields**

Feature	Description
<p><b>Procedure Summary</b></p> 	<ul style="list-style-type: none"> <li>The number on the <b>Procedure Summary</b> button indicates the total number of ablation sessions in the procedure. Press the button to display the <b>Procedure Summary</b>.</li> <li>The background of the <b>Procedure Summary</b> button changes from black to orange when it is selected, and the procedure information is displayed.</li> </ul>
<p><b>Total Time</b> (Procedure Summary)</p>	<ul style="list-style-type: none"> <li>The <b>Total Time</b> field displays the total duration of all ablation sessions in the procedure in hours/minutes/seconds (00:00:00).</li> </ul>
<p><b>Start Time</b></p>	<ul style="list-style-type: none"> <li>An ablation procedure starts when <b>End Procedure</b> or <b>Reload Software</b> is selected from the <b>Shut Down</b> menu.</li> <li>The <b>Start Time</b> field in the ablation session area displays the time that the ablation session started.</li> <li>The <b>Start Time</b> field in the procedure summary area displays the time that the ablation procedure started.</li> </ul>
<p><b>Date</b> (Procedure Summary)</p>	<ul style="list-style-type: none"> <li>The <b>Date</b> field displays the date that the procedure was started.</li> </ul>
<p><b>Total Irrigation</b> (Procedure summary)</p>	<ul style="list-style-type: none"> <li>The <b>Total Irrigation</b> field displays the amount of fluid that was delivered during the ablation procedure.</li> <li>The number represents the total amount of fluid used during the ablation procedure, including during low flow and during pre-ablation and post-ablation high flow. It does not include fluid delivered during flushing.</li> </ul>
<p><b>Ablation Session</b></p> 	<ul style="list-style-type: none"> <li>The <b>Ablation Session</b> button represents the individual ablation sessions. Press a number to view the information for the corresponding ablation session.</li> <li>The background of the <b>Ablation Session</b> button changes from black to orange when it is selected, and the corresponding ablation session information is displayed.</li> </ul>
<p><b>Next Previous</b></p> 	<ul style="list-style-type: none"> <li>Press the <b>Next</b> button or the <b>Previous</b> button to display ablation numbers that are not visible.</li> </ul>
<p><b>Total Time</b> (Ablation session)</p>	<ul style="list-style-type: none"> <li>The <b>Total Time</b> field displays the duration of the ablation session in minutes/seconds (MM:SS).</li> </ul>
	<ul style="list-style-type: none"> <li>The yellow <b>Maximum</b>, <b>Average</b>, and <b>Total</b> fields display the maximum, average, and total power reached during the ablation session or ablation procedure.</li> </ul>
	<ul style="list-style-type: none"> <li>The orange <b>Maximum</b> and <b>Average</b> fields display the maximum and average temperature reached during the ablation session or ablation procedure.</li> </ul>

**Table 9 - Post Screen Fields**

Feature	Description
	<ul style="list-style-type: none"> <li>The green <b>Maximum</b> and <b>Minimum</b> fields display the maximum and minimum impedance reached during the ablation session or ablation procedure.</li> <li>The <b>Drop</b> field displays the impedance drop during the ablation session.</li> </ul>
<p><b>Irrigation</b> (Ablation session )</p>	<ul style="list-style-type: none"> <li>The <b>Irrigation</b> field displays the amount of fluid delivered during the ablation session.</li> <li>The number represents the total amount of fluid used during the ablation session, including during low flow and during pre-ablation and post-ablation high flow. It does not include fluid delivered during flushing.</li> </ul>
<p><b>Stop Reason</b> (Ablation session)</p>	<ul style="list-style-type: none"> <li>The <b>Stop Reason</b> field displays the reason why the ablation session stopped.</li> </ul>
<p><b>Catheter Type</b> (Ablation session)</p>	<ul style="list-style-type: none"> <li>The <b>Catheter Type</b> field displays the catheter type that was used during the ablation session.</li> </ul>
 <p>(Ablation session)</p>	<ul style="list-style-type: none"> <li>The <b>Data/Graph</b> button is a toggle button. The selected button has an orange background.</li> <li>The orange <b>Data</b> button indicates that the data for the ablation session is displayed.</li> <li>The orange <b>Graph</b> button indicates that a graph of the ablation session is displayed.</li> </ul>
 <p>(Ablation session)</p>	<ul style="list-style-type: none"> <li>The yellow line and number represent the power reached as the ablation session progressed.</li> <li>The orange line and number represent the temperature reached as the ablation session progressed.</li> <li>The solid green line and number represent the impedance reached as the ablation session progressed.</li> <li>The dashed green line represents the impedance value when the ablation starts.</li> <li>Press on any graph line, then use the <b>Control</b> knob to see the power, temperature, and impedance values at any point during the selected ablation session.</li> </ul>


# 8 Technical Data

## 8.1 Output Specifications: QDOT MICRO™ Catheter

**Figure 11 - Output Specifications**

— Power Level    - - - Accuracy Tolerance


Load Impedance: 50 - 250 Ω


**Figure 12 – Output Specifications - QMODE™ and QMODE+™ Ablation Modes**

— Full power    - - - Half power

QMODE™ Ablation Mode

QMODE+™ Ablation Mode

## 9 Troubleshooting

The Alarm Messages table explains each alarm message that may appear when a QDOT MICRO™ Catheter is used.

**Table 10 - QDOT MICRO™ Catheter Alarm Messages**

Message • Reasons	What To Do	Icon	Alarm Tone
<b>CARTO™ - catheter error</b> • The generator has received a message from the CARTO™ System.	• Check catheter error message(s) on CARTO™ System.		✓
<b>CARTO™ - communication error</b> • There is a problem with the communication between the CARTO™ System and the nGEN™ Generator.	• Verify that the Ethernet cable from the CARTO™ System Workstation to the nGEN™ Console is properly connected. • Verify that the Ethernet cable is not damaged. • Reload the software on the monitor(s). • Reload the software on the CARTO™ System. • If the problem persists, turn the console off and then on. Restart the monitor(s). • If the problem persists, replace the cable.		✓
<b>Electrode - delta impedance too high</b> • The impedance change exceeded the maximum impedance delta value set in the preset.	• Verify that the maximum impedance delta setting in the selected preset is appropriate. • Consider possible clinical reasons such as a steam pop or the catheter being pushed suddenly into the tissue. • Restart the console and monitor(s).		✓
<b>Electrode - impedance too high</b> • The impedance exceeded the maximum impedance value set in the preset.	• Verify that the maximum impedance setting in the selected preset is appropriate. • Verify that the catheter is properly inserted in the body and is out of the sheath. • Disconnect the indifferent electrode from the console, then reconnect it. • Disconnect the catheter cable from the console or from the CARTO™ PIU, then reconnect it. • If the catheter cable is connected to the CARTO™ PIU, verify that the cable from the CARTO™ PIU to the console is properly connected. • Improve the contact of the indifferent electrode with the patient's skin or use a different indifferent electrode. • Consider applying an indifferent electrode to a different place on the patient's body to avoid oil, hair, and dirt that increase the impedance.		✓

**Table 10 - QDOT MICRO™ Catheter Alarm Messages**

Message • Reasons	What To Do	Icon	Alarm Tone
<b>Electrode - impedance too low</b> <ul style="list-style-type: none"> <li>The impedance was below the minimum impedance value set in the preset. Rarely, such as in pediatric cases, the patient has low body impedance and there is no problem.</li> </ul>	<ul style="list-style-type: none"> <li>Verify that the minimum impedance setting in the selected preset is appropriate.</li> <li>Restart the console and monitor(s).</li> <li>Replace the sterile catheter cable or the cable between the console and the CARTO™ PIU.</li> <li>Replace the catheter and/or catheter cable(s).</li> </ul>		✓
<b>Electrode - power too high</b> <ul style="list-style-type: none"> <li>During ablation, the generator delivered higher power than the power value set in the preset.</li> <li>The console has an internal error.</li> </ul>	<ul style="list-style-type: none"> <li>Turn off the console. Let it cool down, then turn it on again.</li> <li>Contact Customer Support even if the problem is resolved.</li> </ul>		✓
<b>Electrode - power too low</b> <ul style="list-style-type: none"> <li>During ablation, the generator delivered lower power than the power value set in the preset.</li> </ul> <p>The console has an internal error.</p>	<ul style="list-style-type: none"> <li>Turn off the console. Let it cool down, then turn it on again.</li> <li>Contact Customer Support even if the problem is resolved.</li> </ul>		✓
<b>Electrode – temperature invalid</b> <ul style="list-style-type: none"> <li>When the temperature is invalid and not displayed on the screen.</li> </ul>	<ul style="list-style-type: none"> <li>Verify that a QDOT MICRO™ Catheter is being used and is properly connected.</li> <li>Verify that TX cable M581007 is being used and is properly connected to the Console and the CARTO™ System.</li> </ul>		✓
<b>Electrode - temperature not decreasing</b> <ul style="list-style-type: none"> <li>The temperature of the catheter's electrode is above the target temperature for longer than the time that is predefined in the system.</li> </ul>	<ul style="list-style-type: none"> <li>Verify that the temperature settings in the selected preset are appropriate.</li> <li>Reposition the catheter.</li> <li>Check irrigation flow (lines, catheter tip, pump).</li> <li>Restart the console and the monitor(s).</li> <li>Replace the catheter and/or the catheter cable(s).</li> </ul>		✓
<b>Electrode - temperature slope too high</b> <ul style="list-style-type: none"> <li>The temperature of the ablation electrode increased suddenly.</li> <li>The temperature sensor(s) in the catheter intermittently disconnected and reconnected.</li> <li>The wires in the catheter cable(s) intermittently disconnected and reconnected.</li> </ul>	<ul style="list-style-type: none"> <li>Reposition the catheter.</li> <li>Check irrigation flow (lines, catheter tip, pump).</li> <li>Restart the console and the monitor(s).</li> <li>Replace the catheter and/or the catheter cable(s).</li> </ul>		✓

**Table 10 - QDOT MICRO™ Catheter Alarm Messages**

<b>Message</b> • Reasons	<b>What To Do</b>	<b>Icon</b>	<b>Alarm Tone</b>
<b>Electrode - temperature too high</b> • The temperature of the catheter ablation electrode exceeded the temperature cutoff value set in the preset. Possible reasons are: • The ablation electrode is too hot. • The cutoff value is too low.	<ul style="list-style-type: none"> <li>• Wait a few seconds for the temperature to drop before restarting ablation.</li> <li>• Consider possible clinical reasons for the ablation electrode to overheat.</li> <li>• Restart the console and the monitor(s).</li> </ul>		
<b>Electrode - temperature too low</b> • The generator does not recognize a valid temperature reading.	<ul style="list-style-type: none"> <li>• Verify that the catheter and catheter cable(s) are compatible with the generator.</li> <li>• Restart the console and the monitor(s).</li> <li>• Replace the catheter and/or the catheter cable(s).</li> </ul>		

Biosense Webster, Inc.  
31 Technology Drive, Suite 200  
Irvine, California 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010  
Fax: +1-909-468-2905  
[www.biosensewebster.com](http://www.biosensewebster.com)

### **Customer Support**

Biosense Webster, Inc.  
Tel: +1-866-4SERVCE (+1-866-473-7823)  
Email: [complaints@bwius.jnj.com](mailto:complaints@bwius.jnj.com)

### **Customer Support – Outside USA**

Biosense Webster  
Tel: +32-2-7463-463  
Email: [techsupport@its.jnj.com](mailto:techsupport@its.jnj.com)

BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, QDOT MICRO, QMODE, QMODE+, and nGEN are trademarks of Biosense Webster, Inc. The third party trademarks listed herein are trademarks of their respective owners.

Patent information is available at [www.biosensewebster.com/virtualpatentmarking](http://www.biosensewebster.com/virtualpatentmarking)

© Biosense Webster, Inc. 2020-2022

Revised: 2022-08

EML1384021FU02A

## **Sterile Interface Cables**

Instructions for Use  
**Sterile Interface Cables**

(en) English

Page

4

---

**Symbol Definitions**

Medical Device

Sterilized using ethylene oxide

Catalog number

US catalog number

Batch code

Unique Device Identifier

Manufacturer

Date of manufacture

Authorized Representative in the European Community

Packaging unit

Single sterile barrier system with protective packaging outside

Use-by date

Caution

Consult instructions for use

[www.e-ifu.com](http://www.e-ifu.com)  
e-IFU indicator

Do not use if package is damaged

Do not use if package is opened

Keep away from sunlight

Keep dry

Temperature limit

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Waste Electrical and Electronic Equipment Directive (WEEE)—Dispose of in a manner consistent with required EU Directives in your local jurisdiction.

Pin Connector

**ASSEMBLED MEXICO**

Assembled in Mexico

## Sterile Interface Cables

**Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.**

### INDICATIONS FOR USE

These cables provide a means to connect a Biosense Webster electrophysiology catheter to the appropriate equipment. The cables may be reused subject to the cleaning and sterilization restrictions in this document.

**Table 1 – Cables**

Cable Name	Purpose	Catalog Numbers	Manufacturing Numbers
TX eco EXT Cable	Connects a therapeutic "eco" catheter to a CARTO™ 3 System.	D135703	D-1357-03-S
Interface Cable, CELSIUS™ Catheter to nGEN™ Generator	Connects a CELSIUS™ Catheter to an nGEN™ Generator.	D133703	D-1337-03-S

### CONTRAINDICATIONS

There are no known contraindications for these cables.

### PATIENT TARGET GROUP

These cables are targeted at patients who have been diagnosed with cardiac arrhythmias and are undergoing an electrophysiology procedure.

### INTENDED USERS

The intended users of these cables are appropriately trained personnel in a fully equipped electrophysiology laboratory.

### CLINICAL BENEFIT

These cables provide a means to connect a Biosense Webster electrophysiology catheter, external reference patch, or sheath to the appropriate equipment.

### DIRECTIONS FOR USE

- Do not use the cable after the use-by date on the product label.
- Before using or reusing the cable, visually inspect it for damage such as crushed or elongated sections, cuts, kinks, or nicks. If the cable is damaged, do not use the cable and dispose of it per the "DISPOSAL" section below.
- Connect each cable to the proper equipment by referring to the connector designations at each end of the cable.
- Follow the instructions in the "CLEANING AND STERILIZATION," "STORAGE," and "DISPOSAL" sections below.

### CLEANING AND STERILIZATION

These cables are provided sterile. They may be resterilized as described below. These cables meet the requirements of the ANSI/AAMI EC53-2013.

**Table 2 – Cleaning and Sterilization**

<b>Warnings and Precautions</b>	<ol style="list-style-type: none"> <li>Ensure proper cable connection by referring to the connector designations at each end of the cable.</li> <li>Do not use the cable after the use-by date on the product label.</li> <li>Do not immerse the cable connectors in liquid. Do not wipe the cable connectors with liquid. The presence of liquid inside the connectors may result in adverse effects, including improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.</li> <li>Do not expose the cable to organic solvents because they may damage the cable.</li> <li>Patient or operator injury may result from improper handling, improper connection, or improper use of the cable.</li> <li>To prevent damage, store the cable in a clean and secure area.</li> <li>If a break occurs in the cable wire or if the cable becomes otherwise electrically discontinuous, arcing may occur in the patient-return or active circuit and may burn the patient or create a fire.</li> <li>Failure to properly clean the cable could lead to inadequate resterilization.</li> <li>Do not use radiation sterilization because it will damage the cable.</li> <li>When sterilizing the cables, do not exceed 140°C because temperatures higher than 140°C may damage the cable.</li> <li>Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.</li> </ol>
<b>Limitations on Cleaning and Sterilization</b>	Clean and sterilize the cable after each use. The maximum number of sterilization cycles is specified in Table 3. Exceeding the specified number of sterilization cycles may impair the proper functioning of the cable.

<b>Initial Treatment at Point of Use</b>	There are no requirements for initial treatment at point of use.						
<b>Containment and Transportation</b>	Place the cable in a clean biohazard bag after use for subsequent cleaning, disinfection, and sterilization. Clean, disinfect, and sterilize the cable as soon as it is reasonably practical after use.						
<b>Preparation before Cleaning</b>	Visually inspect the cable for any damage such as crushed or elongated sections, cuts, kinks, or nicks. If damage is observed, do not use the cable and dispose of it per the "Disposal" section below.						
<b>Cleaning: Automated</b>	Automated cleaning of the cable is not recommended because fluid may enter the connectors and leave moisture or residue in or on the connectors. Resulting adverse effects may include improper functioning (such as noisy signals or signal degradation) or arcing of electricity between the connector pins, which can lead to patient injury.						
<b>Cleaning: Manual</b>	<ol style="list-style-type: none"> <li>1. Wear sterile gloves while handling the cable.</li> <li>2. Using a wipe*, use a twisting motion to wipe the entire surface of the cable, but not the cable connectors, to remove soil. Pay particular attention to complex cable features (such as non-smooth areas, joints, and crevices).</li> <li>3. Using another wipe*, use a twisting motion to thoroughly wipe the entire surface of the cable to remove soil.</li> <li>4. Visually inspect the cable to ensure there is no visible soil. Additional wipes may be used to remove any visible soil.</li> <li>5. Leave the cable in ambient conditions until it is visibly dry.</li> <li>6. When the cable is dry, visually inspect the entire surface of the cable in adequate lighting to ensure complete removal of visible soil.</li> </ol> <p>*Recommended Wipes: CaviWipes™ or any other wipes with the chemical composition indicated in the table below may be used.</p> <table border="1"> <thead> <tr> <th>Chemical Composition</th> <th>% of Chemical</th> </tr> </thead> <tbody> <tr> <td>Isopropanol</td> <td>17.2%</td> </tr> <tr> <td>Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride</td> <td>0.28%</td> </tr> </tbody> </table>	Chemical Composition	% of Chemical	Isopropanol	17.2%	Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride	0.28%
Chemical Composition	% of Chemical						
Isopropanol	17.2%						
Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride	0.28%						
<b>Disinfection</b>	<ol style="list-style-type: none"> <li>1. Wear sterile gloves while handling the cable.</li> <li>2. Ensure manual cleaning (see the "Cleaning: Manual" instructions in this table) is performed before disinfection.</li> <li>3. Using a wipe*, use a twisting motion to thoroughly wipe the entire surface of the cable. Pay particular attention to complex cable features (such as non-smooth areas, joints, and crevices).</li> <li>4. Allow the surface of the cable to remain visibly wet for 3 minutes at room temperature (approximately 68°F/20°C).</li> <li>5. Additional wipes may be used to ensure that the cable remains visibly wet for 3 minutes.</li> <li>6. Leave the cable in ambient conditions until it is visibly dry.</li> <li>7. Visually inspect the entire surface of the cable in adequate lighting to ensure the cable is visibly clean.</li> </ol> <p>*Recommended Wipes: CaviWipes™ or any other wipes with the chemical composition indicated in the table below may be used.</p> <table border="1"> <thead> <tr> <th>Chemical Composition</th> <th>% of Chemical</th> </tr> </thead> <tbody> <tr> <td>Isopropanol</td> <td>17.2%</td> </tr> <tr> <td>Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride</td> <td>0.28%</td> </tr> </tbody> </table>	Chemical Composition	% of Chemical	Isopropanol	17.2%	Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride	0.28%
Chemical Composition	% of Chemical						
Isopropanol	17.2%						
Diisobutylphenoxyethoxyethyl dimethylbenzyl ammonium chloride	0.28%						
<b>Drying</b>	Drying is performed as part of the disinfection process (see the "Disinfection" instructions in this table).						
<b>Maintenance, Inspection, and Testing</b>	Before packaging, visually inspect the cable. If visible soil is present, repeat the "Cleaning: Manual" and "Disinfection" instructions in this table. If there is evidence of damage, do not use the cable and dispose of it per the "Disposal" section below.						
<b>Packaging</b>	If the cable is stored after cleaning and disinfection, place it in a breathable pouch that is approved by the FDA and/or by local authorities. Coil the cable into a loop when placing it into the pouch. Do not fold the cable because doing so may break the cable. Use a pouch large enough to prevent stress to the pouch seams and to prevent excessive bending of the cable. A pouch size of 19 cm x 33 cm (7.5" x 13") or larger is recommended. For storage conditions, refer to the "Storage" section below.						
<b>Resterilization</b>	See Table 3 and Table 4 for resterilization methods, parameters, and the maximum number of sterilization cycles. Resterilize the cable after each use.						

The instructions provided above have been validated by the medical device manufacturer as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing as performed (using equipment, materials, and personnel in the processing facility) achieves the desired result. This requires validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences.

Table 3 – Validated Resterilization Methods for Each Cable

Maximum Number of EtO Resterilization Cycles	Maximum Number of Steam Resterilization Cycles	Cable Name	Cable Catalog Numbers
20X	20X	TX eco EXT Cable	D135703
		Interface Cable, CELSIUS™ Catheter to nGEN™ Generator	D133703

Table 4 – Resterilization Parameters for Cables

	Exposure Temperature	Exposure Time	Minimum Drying Times	Gas Concentration	Humidity (RH)	Accommodation/Aeration Time at 55°C
Gravity Steam Resterilization	121°C (250°F)	30 min	15 - 30 min	---	---	---
	132°C (270°F)	15 min	15 - 30 min	---	---	---
	135°C (275°F)	10 min	30 min	---	---	---
100% EtO Gas Resterilization	55°C (131°F)	1 hr	---	725 mg/L	30% - 80%	12 hr

**STORAGE**

Refer to the product label for recommended storage conditions.

**DISPOSAL**

Recycle components or dispose of the product and its residual elements or waste items in accordance with local laws and regulations.

**REPORTING INCIDENTS**

Per regulation 2017/745/EU on medical devices, if this device is used in the European Union (or in a country with an identical regime) and a serious incident occurs during the use of this device or as a result of use of this device, report the incident to the manufacturer and/or its authorized representative and to your national authority.

**ADVERSE REACTIONS**

A number of adverse reactions have been documented for electrophysiology procedures including:

heart block, pulmonary vein stenosis, esophageal fistula and/or injury, stroke (cerebrovascular accident), other arrhythmias (outside diagnosis), life threatening arrhythmias, myocardial infarction, cardiac perforation, pericardial effusion, cardiac tamponade, thrombosis, embolism, pulmonary embolism, air embolism, valvular damage, phrenic nerve injury, vagal nerve injury, pericarditis, coronary artery stenosis, vessel perforation (peripheral and/or central), soft tissue injury, persistent atrial communication, device related infection, embolization of components, device entrapment, surgical intervention (additional), skin burns, and localized skin reaction.

The following complications associated with cardiac catheterization have also been reported in the literature:

major bleed, hematoma, reaction to medications, allergic reaction, vascular access complication, damage to vasculature, implanted device interactions, renal artery stenosis, pneumothorax, ST segment changes, fluid overload, urinary catheter complications, hypotension, sepsis, wound infection, respiratory failure, heart injury, renal injury, heart failure, cardiac arrest and death.

**DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY**

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.

WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

**ELECTRONIC INSTRUCTIONS FOR USE**

This document is available at [www.e-ifu.com](http://www.e-ifu.com).

---

NOTES

**Biosense Webster, Inc.**  
31 Technology Drive, Suite 200  
Irvine, California 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010  
Fax: +1-909-468-2905  
[www.biosensewebster.com](http://www.biosensewebster.com)

BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, CELSIUS, and nGEN are trademarks of Biosense Webster, Inc.

The third-party trademarks used herein are trademarks of their respective owners.

Patent information is available at [www.biosensewebster.com/virtualpatentmarking](http://www.biosensewebster.com/virtualpatentmarking).

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## TX eco Cable

**Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.**

- Do not use if the package is opened or damaged.

### DEVICE DESCRIPTION

The TX eco Cable is used with an extension cable to connect a Biosense Webster therapeutic catheter to the Patient Interface Unit (PIU) of the CARTO™ 3 System (Version 6 and later). The extension cables are listed below:

- For the QDOT™ Micro Catheter: TX eco EXT Cable D135703
- For other Biosense Webster therapeutic catheters: D128603 (catalog # CR3434CT) or D128604 (catalog # CR3425CT)

The TX eco Cable communicates data from a Biosense Webster therapeutic catheter to the CARTO™ 3 System and the RF generator. The information communicated is listed below.

- Force signals
- Location signals
- IC signals from the 3 microelectrodes in the tip of the catheter (for catheters with microelectrodes only)
- Temperature measurements from the 6 thermocouples in the tip of the catheter

### INDICATIONS FOR USE

The TX eco Cable is used to connect a Biosense Webster therapeutic catheter to the Patient Interface Unit (PIU) of the CARTO™ 3 System (Version 6 and later).

### WARNINGS

- Before using this cable, read the *Instructions for Use* for the Biosense Webster therapeutic catheter. If this cable is being used with the QDOT MICRO™ Module, refer to the *Instructions for Use* for the CARTO™ 3 System and the *Instructions for Use and Release Notes* for the QDOT MICRO™ Module before using this cable with the QDOT MICRO™ Module.
- Only personnel who have read and understood the contents of this document may set up and use this cable.
- This cable is only validated for use with a Biosense Webster therapeutic catheter. Do not connect this cable to any other catheter.
- This cable is intended for use only with a Biosense Webster therapeutic catheter, a compatible RF generator, a CARTO™ 3 System, and Biosense Webster cables. Consult your Biosense Webster representative for details.
- This cable must be at least 1 meter (39.4 in) away from the location pad of the CARTO™ 3 System during the procedure.
- This cable must be connected to the PIU of the CARTO™ 3 System for a 5-minute warm-up period before use to ensure accurate temperature readings. After the warm-up period, the temperature readings stabilize. If the cable is used before the end of the warm-up period, a temperature drift might occur.

## INSTALLING AND OPERATING THE CABLE

### Attaching the Cable Holder

- Place the clamp portion of the cable holder (part number M652802) on the side of the CARTO™ System Cart handle or over the bed rail. Then tighten the knob on the cable holder.
- Insert the cable in the cable holder with the cable facing down.

### Connecting the Cable

- Connect the cable's red connector to the MAP socket on the PIU of the CARTO™ 3 System.
- Connect the cable's yellow connector to the QUAD B or DECA socket on the PIU.  
Note 1: Both connectors must be connected to the PIU.  
Note 2: Connecting the connector to QUAD A may result in undesirable ECG noise.
- If using a QDOT™ Micro Catheter, connect the catheter to one end of the TX eco EXT Cable. Connect the other end of the TX eco EXT Cable to the large socket at the end of the TX eco Cable.
- If using another Biosense Webster therapeutic catheter, connect the catheter to one end of the extension cable. Connect the other end of the extension cable to the small round socket at the end of the TX eco Cable.

### Disconnecting the Cable

- Disconnect the red connector from the PIU of the CARTO™ 3 System by pulling the connector from the slide release grip.
- Then disconnect the yellow connector from the PIU by pulling the connector from the slide release grip.  
Note: Both connectors must be disconnected from the PIU. Do not leave one cable connected and the other disconnected.
- After disconnecting the cable, wait 10 seconds before reconnecting the cable to the PIU.

### LED Status

The status of the cable is indicated by an LED. After the cable is properly connected to the PIU of the CARTO™ 3 System and to the Biosense Webster therapeutic catheter, the cable performs a Built-In Test (BIT).

LED	Cable Status	Comments
Green: slow blink	BIT in progress	---
Green: solid	Ready for use	---
Red: fast blink	Error, BIT failed	The temperatures are not sent to the CARTO™ 3 System or RF generator. See the <i>Troubleshooting</i> section in this document.

## CARING FOR THE CABLE

### Cleaning

The cable does not require disinfection or sterilization. Do not steam, autoclave, or otherwise sterilize the cable.

If dust or debris appears on the cable, clean the cable as follows:

- Disconnect the cable.
- Clean the outside of the cable by wiping it with a cloth dampened with alcohol-free hand soap and water.
- Ensure that none of the soap solution enters the cable connectors or sockets.
- Ensure that the cable is dry before connecting it to the CARTO™ 3 System or to a catheter.

### Maintenance

There are no parts in the cable that may be serviced by the user. If the cable fails to perform, contact Customer Support or your Biosense Webster representative for a replacement. The expected useful life of the cable is three years.

### Disposal

Recycle components, or dispose of the product and its residual elements or waste items, in accordance with local laws and regulations.

## TECHNICAL DATA

### Specifications

DC Input	The cable receives power input from the PIU of the CARTO™ 3 System.
Weight	1 kg (2.2 lb)
Dimensions	10 in long x 1.9 in wide x 3 in high (254 mm x 48 mm x 76 mm)
Temperature Accuracy	≤ 2 °C

### Operating, Storage, and Shipping Specifications

	Minimum	Maximum
<b>Operating Specifications</b>		
Ambient Temperature	10°C	30°C
Relative humidity*	25%	75%
<b>Storage and Shipping Specifications</b>		
Ambient Temperature	-30°C	65°C
Relative humidity*	10%	95%

\* In accordance with MIL-STD-1695, relative humidity levels shall be within the range of 30% to 70% in areas where electronic parts and hybrid microcircuits are handled or processed. MIL-STD-1695 requires the same level of relative humidity controls for handling and storage areas, except when items are covered or protected.

### EMC Information

The TX eco Cable is intended for use in the electromagnetic environment as specified in the *Instruction for Use* for the CARTO™ 3 System.

## TROUBLESHOOTING

If the LED on the TX eco Cable is red and is blinking fast, the BIT has failed and has caused an error (see the *LED Status* section in this document). Follow the steps below to correct the problem.

- Disconnect the catheter from the cable.
- Disconnect the cable from the PIU of the CARTO™ 3 System.
- Reconnect the cable to the PIU.
- Reconnect the catheter to the cable.
- If a QDOT MICRO™ Module is being used, follow the directions in the *Instructions for Use and Release Notes* for the QDOT MICRO™ Module.
- If the problem persists, contact Customer Support or your Biosense Webster representative.

In some situations, where 15 kV of air discharge or 8 kV of contact discharge is applied on the system, the CARTO™ 3 System displays 3 error messages simultaneously: Map points cannot be acquired (Error 401), Patient Body Reference has moved (Error 256), and a popup message indicating the Patient Body Reference change. If this occurs, restart the PIU.

## DISCLAIMER OF WARRANTY AND LIMITATION OF LIABILITY

**THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE PRODUCT(S) DESCRIBED HEREIN. UNDER NO CIRCUMSTANCES SHALL BIOSENSE WEBSTER, INC., OR ITS AFFILIATED COMPANIES, BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY APPLICABLE LAW.**

**WITHOUT LIMITING THE FOREGOING, BIOSENSE WEBSTER, INC. OR ITS AFFILIATED COMPANIES, SHALL NOT BE LIABLE FOR ANY SPECIAL, DIRECT, INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES, ARISING OUT OF THE REUSE OF ANY PRODUCT(S) LABELED FOR SINGLE USE OR WHERE REUSE IS PROHIBITED BY APPLICABLE LAW.**

Descriptions and specifications appearing in Biosense Webster, Inc. printed matter, including this publication, are informational only and meant solely to generally describe the product(s) at the time of manufacture and are not made or given as a warranty of the prescribed product(s) in any way.

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BIOSENSE WEBSTER, the Biosense Webster logo, CARTO, and QDOT MICRO are trademarks of Biosense Webster, Inc.