

August 4, 2023

Abbott Laboratories Neha Vatsyayan Regulatory Affairs Project Manager 4551 Great America Pkwy Santa Clara, California 95054

Re: K220031

Trade/Device Name: Alinity h-series System Regulation Number: 21 CFR 864.5220 Regulation Name: Automated Differential Cell Counter Regulatory Class: Class II Product Code: GKZ Dated: March 29, 2023 Received: March 31, 2023

Dear Neha Vatsyayan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Min Wu -S

Min Wu, Ph.D. Branch Chief Division of Immunology and Hematology Devices OHT7: Office of In Vitro Diagnostics Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K220031

Device Name Alinity h-series System

Indications for Use (Describe)

The Alinity h-series System is an integrated hematology analyzer (Alinity hq) and slide maker stainer (Alinity hs) intended for screening patient populations found in clinical laboratories by qualified health care professionals. The Alinity h-series System can be configured as:

• One standalone automated hematology analyzer system.

• A multimodule system that includes at least one Alinity hq analyzer module and may include one Alinity hs slide maker stainer module.

The Alinity hq analyzer module provides complete blood count and a 6-part white blood cell differential for normal and abnormal cells in capillary and venous whole blood collected in K2EDTA or K3EDTA. The Alinity hq analyzer provides quantitative results for the following measurands: WBC, NEU, %N, LYM, %L, MON, %M, EOS, %E, BASO, %B, IG, %IG, RBC, HCT, HGB, MCV, MCH, MCHC, MCHr, RDW, NRBC, NR/W, RETIC, %R, IRF, PLT, MPV, %rP. The Alinity hq analyzer module is indicated to identify patients with hematologic parameters within and outside of established reference ranges. The Alinity hs slide maker stainer module automates whole blood film preparation and staining and stains externally prepared whole blood smears.

For in-vitro diagnostic use.

Type of Use (Select one or both, as applicable)							
Prescription Use (Part 21 CFR 801 Subpart D)							
CONTINUE ON A SEPARATE PAGE IF NEEDED.							
This section applies only to requirements of the Paperwork Reduction Act of 1995.							

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Section 5: 510(k) Summary

This summary of the 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

I. Applicant Name

Abbott Laboratories 4551 Great America Pkwy, Santa Clara, CA 95054 Date Prepared: April 28, 2023

Contact: Neha Vatsyayan Regulatory Affairs Project Manager Abbott Diagnostics Division Phone: (408) 313 4401 E-Mail: <u>neha.vatsyayan@abbott.com</u>

II. Device Information

Trade name (proprietary name): Alinity h-series system Common name (usual name): Automated Hematology Analyzer and Slide Maker Stainer Classification Name: Automated Differential Cell Counter

III. Regulatory Information

Alinity hq (Analyzer) Device Classification: Class II Regulation Description: Automated Differential Cell Counter Governing Regulation: 21 CFR 864.5220 Code: GKZ Alinity hs (Slide maker stainer) Device Classification: Class I Regulation Description: Automated Slide Stainer Governing Regulation: 21 CFR 864.3800 Code: KPA

IV. Predicate Device

Sysmex® XN-Series (XN-10, XN-20) Automated Hematology Analyzers (K112605)

V. Device Description

The Alinity h-series system is a multimodule system that consists of different combinations of one or more of the following modules: a quantitative multi-parameter automated hematology analyzer (Alinity hq) and an automated slide maker stainer (Alinity hs).

The modules are designed to fit together. Each module has an internal conveyor that enables racks of specimen tubes to be transported between modules. The system can move racks between modules to perform different tests on a given specimen (e.g., make slide smears on the Alinity hs).

Principles of Operation

The Alinity hq uses flow cytometry and absorption spectrophotometry technologies to measure, count, and calculate hematological parameters in samples.

Two methods are used to introduce a specimen to the Alinity hq module:

- Closed-tube processing mode
- Open-tube processing mode

The Alinity hs module creates and stains smears from whole blood samples in addition to staining externally prepared smears for morphologic review. The operator selects and may configure staining protocols as needed by the laboratory. The Alinity hs module is configured with the May-Grünwald-Giemsa stain or the Wright-Giemsa stain.

Intended Use

The Alinity h-series system is an integrated hematology analyzer (Alinity hq) and slide maker stainer (Alinity hs) intended for screening patient populations found in clinical laboratories by qualified health care professionals. The Alinity h-series can be configured as:

- One standalone automated hematology analyzer system.
- A multimodule system that includes at least one Alinity hq analyzer module and may include one Alinity hs slide maker stainer module.

The Alinity hq analyzer module provides complete blood count and a 6-part white blood cell differential for normal and abnormal cells in capillary and venous whole blood collected in K₂EDTA or K₃EDTA. The Alinity hq analyzer provides quantitative results for the following measurands: WBC, NEU, %N, LYM, %L, MONO, %M, EOS, %E, BASO, %B, IG, %IG, RBC, HCT, HGB, MCV, MCH, MCHC, MCHr, RDW, NRBC, NR/W, RETIC, %R, IRF, PLT, MPV, %rP. The Alinity hq analyzer module is indicated to identify patients with hematologic parameters within and outside of established reference ranges. The Alinity hs slide maker stainer module automates whole blood film preparation and staining and stains externally prepared whole blood smears.

For *in-vitro* diagnostic use.

Definitions of Reportable Parameters

Definitions of the reportable parameters are presented in Table 5-1.

Abbreviation	Definition			
White Blood Cell Parameters				
BASO	Basophil absolute			
21150	concentration			
%B	Basophil percentage of WBCs			
702	(%BASO)			
EOS	Eosinophil absolute			
205	concentration			
%F	Eosinophil percentage of			
70L	WBCs (%EOS)			
IG	Immature Granulocyte			
10	concentration			
%IG	Immature Granulocyte			
/010	percentage			
IVM	Lymphocyte absolute			
	concentration			
04 I	Lymphocyte percentage of			
70 L	WBCs (%LYM)			
MONO	Monocyte absolute			
MONO	concentration			
04 M	Monocyte percentage of			
%0 IVI	WBCs (%MON)			
NEU	Neutrophil absolute			
NEU	concentration			
% N	Neutrophil percentage of			
701N	WBCs (%NEU)			
WBC	White Blood Cell			
WDC	concentration			

Table 5-1
Peripheral Whole Blood Reportable Parameters

Abbreviation	Definition			
Red Blood Cell Parameters				
НСТ	Hematocrit			
HGB	Hemoglobin concentration			
IRF	Immature Reticulocyte Fraction			
МСН	Mean Cell Hemoglobin			
МСНС	Mean Cell Hemoglobin Concentration			
MCHr	Mean cell hemoglobin of the reticulocyte			
MCV	Mean Cell Volume			
NRBC	Nucleated red blood cell absolute concentration			
NR/W	NRBCs per 100 WBCs			
RETIC	reticulocyte concentration (RETC)			
RBC	Red Blood Cell concentration			
RDW	red blood cell distribution width			
%R	Reticulocyte percentage of RBCs (%RETC)			

Platelet Parameters			
PLT	Platelet concentration		
MPV	Mean Platelet Volume		
%rP	Reticulated Platelet		
/011	percentage		

VI. Comparison of Technological Characteristics

The similarities and differences between the subject device and the predicate device are presented in the Device Similarities and Differences tables below.

	Device Similari	ties
-	Subject Device:	Predicate Device:
Item Intended Use	 Almity h-series System The Alinity h-series system is an integrated hematology analyzer (Alinity hq) and slide maker stainer (Alinity hs) intended for screening patient populations found in clinical laboratories by qualified health care professionals. The Alinity h-series can be configured as: One standalone automated hematology analyzer system. A multimodule system that includes at least one Alinity hq analyzer module and may include one Alinity hs slide maker stainer module. The Alinity hq analyzer module provides complete blood count and a 6-part white blood cell differential for normal and abnormal cells in capillary and venous whole blood collected in K2EDTA or K3EDTA. The Alinity hq analyzer provides quantitative results for the following measurands: WBC, NEU, %N, LYM, %L, MONO, %M, EOS, %E, BASO, %B, IG, %IG, RBC, HCT, HGB, MCV, MCH, MCHC, MCHr, RDW, NRBC, NR/W, RETIC, %R, IRF, PLT, MPV, %rP. The Alinity hq analyzer module is indicated to identify patients with hematologic parameters within and outside of established reference ranges. The Alinity hs slide maker stainer module automates whole blood film preparation and staining and stains externally prepared whole blood smears. 	Sysmex® XN-Series (XN-10, XN-20) The Sysmex XN-10 and XN-20 modules are quantitative multi-parameter automated hematology analyzers intended for <i>in vitro</i> diagnostic use in screening patient populations found in clinical laboratories. The XN-Series modules classify and enumerate the following parameters for whole blood: WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, NEUT%/#, LYMPH%/#, MONO%/#, EO%/#, BASO%/#, NRBC%/#, IG%/#, RDW-CV, RDW-SD, MPV, RET%/#, IRF, IPF, RET-He and has a Body Fluid mode for body fluids. The Body Fluid mode for body fluids. The Body Fluid mode enumerates the WBC-BF, RBC-BF, MN%/#, PMN%/# and TC-BF# parameters in cerebrospinal fluids (CSF), serous fluids (peritoneal, pleural) and synovial fluids. Whole blood should be collected in K ₂ or K ₃ EDTA anticoagulant and serous and synovial fluids in K ₂ EDTA anticoagulant to prevent clotting of fluid. The use of anticoagulants with CSF specimens is neither required nor recommended.

Table 5-2Device Similarities and Differences

Device Similarities (Continued)				
Subject Device:ItemAlinity h-series System		Predicate Device: Sysmex® XN-Series (XN-10, XN-20)		
Test Principle	Test PrinciplePerforms hematology analyses according to flow cytometry method (using Hydro Dynamic Focusing) and absorptionPerforms hematology according to flow cytometry method (using Hydro Dynamic Focusing), and absor 			
Parameters ¹	Whole Blood Mode: WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, MPV, IRF, NEU, %N, LYM, %L, MONO, %M, EOS, %E, BASO, %B, NRBC, NR/W, IG, %IG, RETIC, %R, RDW, MCHr, %rP	Whole Blood Mode: WBC, RBC, HGB, HCT, MCV, MCH, MCHC, PLT, MPV, IRF, NEUT%/#, LYMPH%/#, MONO%/#, EO%/#, BASO%/#, NRBC%/#, IG%/#, RET%/#, RDW-CV, RDW-SD, RET-He#, IPF		
Specimen Type	Whole blood	Whole blood		
Use of Controls/ Calibrators	Yes	Yes		
IPU	Multi-Module connect	Multi-Module connect		
Sample Aspiration/ Fluidic Pathway	Single aspiration pathway	Single aspiration pathway		
Software/Hardware	Rules based rerun / reflex	Rules-based rerun / reflex		

Table 5-2, ContinuedDevice Similarities and Differences

¹ Different names/formats of equivalent parameters are used between the Alinity h-series System and Sysmex[®] XN-series; therefore, equivalent parameters are listed in the same row.

Device Differences				
Item	Subject Device: Alinity h-series System	Predicate Device: Sysmex® XN-Series (XN-10, XN-20)		
Test Principle	The Alinity h-series System uses flow cytometry method with Hydro Dynamic Focusing to analyze whole blood samples including RBC and PLT.	Sysmex XN-Series uses Hydro Dynamic Focusing (Direct Current Detection) for RBC and PLT.		
Parameters	Not Applicable - Body Fluid test selection is not included in this submission.	Body Fluid Mode: WBC-BF, RBC-BF, MN%/#, PMN%/#, TC-BF#		
Specimen TypeNot Applicable - Body Fluid is not included in this submission.		Body Fluids [i.e., cerebrospinal fluids (CSF), serous fluids (peritoneal, pleural) and synovial fluids]		
Reagents	 Diluent HGB Reagent WBC Reagent Retic Reagent 	CELLPACKTM DCL (Diluent) CELLPACKTM DFL (Diluent) LYSERCELL WNR (Lyse) LYSERCELL WDF (Lyse) FLUOROCELL WDF (Lyse) FLUOROCELL WDF (Stain) FLUOROCELL WDF (Stain) FLUOROCELL RET (Stain) FLUOROCELL PLT (Stain) FLUOROCELL WPC (Stain) SULFOLYSER [®] (Lyse)		
Controls/ Calibrators	 Whole Blood: Calibrator - Alinity h-series HemCal Control - Alinity h-series Control 29P No Body Fluids Mode on Alinity hq. 	Whole Blood: • XN-Check - 3 Levels • XN CAL (XN-10/X-20 Calibrator) • XN CAL PF (Platelet F Calibrator) Body Fluid: • XN Check BF – 2 Levels		
Measuring Channels/ Methods Selection	 CBC+Diff (for RBC, WBC, and PLT) CBC+Diff+Retic (for RBC, WBC, PLT and Retic) 	• RET/PLT • WNR, WDF, WNR, WPC (Not available in XN-10) • PLT-F		

Table 5-2, ContinuedDevice Similarities and Differences

Device Differences (Continued)					
Item	Predicate Device: Sysmex® XN-Series (XN-10, XN-20)				
Modules Connected to the Analyzer	Required Water Purification System System Control Center Computer (SCC) <u>Optional</u> Laboratory Automation System (LAS) for automatic sample loading	IPU (Information processing unit) Pneumatic Unit			
Data Transfer Mode	USB, Internet, and Intranet	USB, CD-R, Internet, and Intranet			

Table 5-2, ContinuedDevice Similarities and Differences

VII. Performance Characteristics:

A. Analytical Performance:

1. Method Comparison

The method comparison study was conducted based on guidance from the Clinical and Laboratory Standards Institute CLSI EP09c, 3rd edition² to assess the performance of the Alinity hq when compared to the predicate device, Sysmex (XN-10, XN-20) (K112605). A total of 2,194 unique venous and/or capillary specimens collected in K₂EDTA tubes from pediatric (\leq 21 years) and adult (> 21 years) subjects including a wide variety of disease states (clinical conditions) were tested across 7 clinical sites. For each measurand, each specimen was tested within 8 hours from the time of collection in 1 replicate using either the Closed or Open tube processing mode in the CBC+Diff+Retic test selection on the Alinity h-series System and in 1 replicate on the Sysmex (XN-10, XN-20) System.

Passing-Bablok and Deming regression analyses were performed with the investigational method as the dependent variable (y) and the predicate method as the independent variable (x).

Bias at the critical points (the upper and lower limits of the reference ranges and relevant medical decision levels) were also evaluated for each site individually and for all sites combined. All results were within the predefined acceptance criteria and found to be acceptable. The method comparison results are shown in the table below.

² Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Measurand	N	Sample Range	r (95% CI)	Slope (95% CI)	Intercept (95% CI)
WBC (X 10 ³ /µL)	2002	0.07 - 436.00	1.00 (1.00, 1.00)	1.00 (1.00, 1.00)	0.01 (0.00, 0.03)
%N (%)	1551	7.56 - 98.30	0.99 (0.99, 1.00)	1.00 (1.00, 1.01)	0.22 (-0.10, 0.50)
%L (%)	1640	0.34 - 84.60	1.00 (1.00, 1.00)	1.00 (0.99, 1.00)	0.08 (0.00, 0.16)
%M (%)	1646	0.03 - 49.20	0.98 (0.97, 0.98)	0.99 (0.98, 1.00)	-0.03 (-0.13, 0.06)
%E (%)	1712	0.00 - 37.50	0.99 (0.99, 0.99)	1.02 (1.01, 1.03)	0.02 (0.01, 0.03)
%B(%)	1854	0.00 - 8.37	0.45 (0.41, 0.49)	1.53 (1.47, 1.59)	-0.14 (-0.17, -0.11)
%IG (%) ^a	1545	0.00 - 12.50	0.59 (0.56, 0.62)	0.59 (0.45, 0.73)	-0.30 (-0.44, -0.15)
NR/W ^a	1949	0.00 - 228.00	0.99 (0.99, 0.99)	0.97 (0.93, 1.02)	-0.07 (-0.13, -0.01)
RBC (X 10 ⁶ /µL)	1993	0.60 - 8.03	1.00 (0.99, 1.00)	0.99 (0.99, 0.99)	0.04 (0.03, 0.06)
HGB (g/dL)	2006	1.64 - 23.30	1.00 (1.00, 1.00)	0.99 (0.99, 0.99)	0.24 (0.21, 0.28)
HCT (%)	1999	4.92 - 86.00	0.99 (0.99, 0.99)	1.02 (1.01, 1.02)	-0.49 (-0.68, -0.29)
MCV (fL)	2001	51.40 - 131.00	0.95 (0.94, 0.95)	1.05 (1.03, 1.07)	-4.56 (-6.06, -3.12)
MCH (pg)	1993	15.30 - 47.00	0.98 (0.97, 0.98)	0.97 (0.96, 0.98)	1.25 (0.94, 1.52)
MCHC (g/dL)	1993	25.00 - 39.30	0.66 (0.63, 0.68)	0.97 (0.92, 1.00)	1.51 (0.40, 2.90)
%R (%)	1942	0.12 - 20.80	0.97 (0.96, 0.97)	1.06 (1.04, 1.07)	0.01 (-0.01, 0.03)
IRF	1935	0.00 - 0.70	0.89 (0.88, 0.90)	0.94 (0.92, 0.96)	-0.01 (-0.01, -0.01)
PLT (X 10 ³ /µL)	1933	1.21 - 5144.00	0.99 (0.99, 0.99)	0.97 (0.97, 0.98)	0.27 (-0.41, 1.08)

All Sites Combined - Regression Analysis

^aA summary of the Deming regression model is presented for NRBC, NR/W, IG, and %IG. A summary of the Passing-Bablok regression model is presented for all other measurands.

Measurand	N	Sample Range	r (95% CI)	Slope (95% CI)	Intercept (95% CI)
MPV (fL)	1723	8.04 - 13.30	0.73 (0.71, 0.75)	0.94 (0.91, 0.99)	0.29 (-0.14, 0.65)
%rP (%) ^b	1910	0.55 - 42.10	0.82 (0.81, 0.84)	0.78 (0.76, 0.80)	0.62 (0.55, 0.69)
MCHr (pg)	1933	6.89 - 45.80	0.84 (0.82, 0.85)	1.09 (1.06, 1.12)	-1.28 (-2.27, -0.29)
NEU (X 10 ³ /µL)	1551	0.10 - 55.00	1.00 (1.00, 1.00)	1.01 (1.01, 1.01)	0.00 (-0.02, 0.02)
LYM (X 10 ³ /µL)	1640	0.05 - 27.20	0.99 (0.99, 1.00)	0.99 (0.99, 1.00)	0.02 (0.01, 0.03)
MONO (X 10 ³ /µL)	1646	0.00 - 8.84	0.99 (0.99, 0.99)	1.02 (1.01, 1.03)	-0.02 (-0.03, -0.01)
EOS (X 10 ³ /µL)	1712	0.00 - 4.19	0.99 (0.99, 0.99)	1.02 (1.01, 1.03)	0.00 (0.00, 0.00)
BASO (X 10 ³ /µL)	1854	0.00 - 8.11	0.22 (0.18, 0.26)	1.31 (1.27, 1.37)	0.00 (-0.01, 0.00)
IG (X 10 ³ /µL) ^a	1545	0.00 - 3.15	0.81 (0.80, 0.83)	1.01 (0.85, 1.18)	-0.07 (-0.09, -0.05)
NRBC (X 10 ³ /µL) ^a	1945	0.00 - 17.70	0.91 (0.90, 0.92)	0.88 (0.70, 1.07)	0.01 (0.00, 0.02)
RDW (%)	2003	10.10 - 32.30	0.94 (0.93, 0.94)	0.86 (0.84, 0.87)	2.23 (2.06, 2.45)
RETIC (X 10 ³ /µL)	1935	1.96 - 614.00	0.96 (0.96, 0.97)	1.05 (1.04, 1.06)	0.79 (0.02, 1.64)

All Sites Combined - Regression Analysis (Continued)

^a A summary of the Deming regression model is presented for NRBC, NR/W, IG, and %IG. A summary of the Passing-Bablok regression model is presented for all other measurands.
 ^b %rP (%) on Alinity h-series System is equivalent to the Sysmex XN-10 IPF measurand.

Maagunand	Critical Dainta	Bias		%Bias		
Measurand	Critical Points	Estimate	95% CI	Estimate	95% CI	
	8.00	0.17	0.16, 0.20	2.18	2.03, 2.48	
HGB (g/dL)	12.0	0.14	0.13, 0.17	1.17	1.09, 1.43	
	16.2	0.11	0.09, 0.15	0.65	0.53, 0.90	
	14.0	-0.25	-0.37, -0.11	-1.75	-2.61, -0.81	
	35.4	0.13	0.06, 0.20	0.38	0.18, 0.57	
HCI (%)	46.4	0.33	0.22, 0.45	0.71	0.46, 0.97	
	70.0	0.75	0.51, 1.00	1.07	0.72, 1.42	
	80.0	-0.47	-0.70, -0.29	-0.58	-0.88, -0.36	
MCV (IL)	100	0.55	0.31, 0.82	0.55	0.31, 0.82	
	11.6	0.57	0.54, 0.63	4.93	4.64, 5.42	
RDW (%)	14.0	0.23	0.20, 0.27	1.63	1.46, 1.93	
	1.00	0.01	0.00, 0.02	1.09	0.03, 2.40	
$\mathbf{WDC} (\mathbf{V} 10^{3} \mathbf{I})$	3.54	0.01	0.00, 0.02	0.29	0.02, 0.54	
WBC (X 10 ³ /μL)	9.06	0.01	-0.01, 0.02	0.10	-0.10, 0.21	
	30.0	0.00	-0.08, 0.04	0.01	-0.26, 0.13	
0/ NI (0/)	40.0	0.39	0.25, 0.60	0.99	0.62, 1.50	
%N(%)	70.0	0.52	0.40, 0.60	0.74	0.57, 0.86	
0/1 (0/)	20.0	0.02	-0.03, 0.10	0.09	-0.14, 0.50	
%L(%)	50.0	-0.08	-0.23, 0.10	-0.15	-0.46, 0.20	
	4.00	-0.05	-0.11, -0.01	-1.34	-2.79, -0.25	
%M (%)	8.00	-0.08	-0.12, -0.05	-1.00	-1.53, -0.68	
	0.00	0.02	0.01, 0.03	NA	NA	
%E(%)	6.00	0.16	0.11, 0.21	2.67	1.81, 3.42	
	0.00	-0.14	-0.18, -0.11	NA	NA	
%B(%)	2.00	0.91	0.80, 1.01	45.57	40.05, 50.73	
	0.50	0.01	-0.01, 0.02	1.17	-1.99, 4.10	
	1.42	0.01	0.00, 0.03	0.97	0.03, 1.85	
NEU (Χ 10 ⁹ /μL)	6.34	0.06	0.04, 0.07	0.89	0.66, 1.12	
	25.0	0.22	0.13, 0.30	0.88	0.53, 1.19	
	1.00	0.01	0.00, 0.01	0.77	0.27, 1.28	
LYM (X 10 ³ /µL)	4.00	-0.02	-0.04, 0.00	-0.47	-0.89, 0.05	

All Sites Combined - Estimated Bias at Critical Points

NA = Not applicable since the critical point is zero.

			Bias	%Bias			
Measurand	Critical Points	Estimate	95% CI	Estimate	95% CI		
	0.20	-0.01	-0.02, -0.01	-7.49	-10.18, -5.56		
MONO (X $10^{3}/\mu$ L)	1.00	0.00	0.00, 0.01	0.39	-0.24, 0.95		
	0.00	0.0	0.00, 0.00	NA	NA		
EOS (X 10 ³ /µL)	0.40	0.01	0.01, 0.01	2.72	1.96, 3.50		
	1.50	0.04	0.02, 0.05	2.38	1.56, 3.21		
	0.00	0.00	-0.01, 0.00	NA	NA		
BASO (X 10 ³ /μL)	0.20	0.06	0.05, 0.07	29.22	24.67, 34.63		
	1.00	0.31	0.26, 0.37	30.73	25.66, 36.94		
	0.00	-0.30	-0.44, -0.15	NA	NA		
%IG (%)	1.00	-0.71	-0.74, -0.68	-70.96	-73.67, -68.25		
	2.00	-1.12	-1.26, -0.99	-56.13	-62.86, -49.39		
	0.00	-0.07	-0.13, -0.01	NA	NA		
NR/W	1.00	-0.10	-0.15, -0.04	-9.54	-15.36, -3.72		
DDC (X 106/ L)	4.00	0.00	0.00, 0.01	0.10	-0.03, 0.26		
$RBC(X 10\%\mu L)$	5.60	-0.01	-0.02, 0.00	-0.22	-0.42, -0.02		
	26.7	0.47	0.43, 0.51	1.76	1.60, 1.89		
MCH (pg)	31.9	0.32	0.27, 0.37	1.00	0.85, 1.15		
	32.3	0.43	0.39, 0.50	1.33	1.20, 1.55		
MCHC (g/dL)	35.9	0.31	0.18, 0.40	0.86	0.49, 1.11		
$\mathbf{DETIC} (\mathbf{Y} 10^{3} 10)$	32.0	2.38	1.86, 2.97	7.43	5.80, 9.29		
$KETIC (X 10-7 \mu L)$	129	7.18	6.09, 8.21	5.57	4.72, 6.36		
$0/\mathbf{P}(0/)$	0.80	0.05	0.04, 0.07	6.29	4.79, 8.68		
%K(%)	2.30	0.13	0.12, 0.15	5.82	5.24, 6.68		
IDE	0.04	-0.01	-0.01, -0.01	-25.91	-30.19, -20.67		
IKF	0.37	-0.03	-0.04, -0.02	-8.01	-9.72, -6.41		
	10.0	-0.02	-0.75, 0.85	-0.24	-7.51, 8.53		
$\mathbf{D}\mathbf{I} = (\mathbf{V} + 0^3 0 + 1)$	165	-4.58	-5.06, -3.93	-2.78	-3.06, -2.38		
PLI (Λ 10 ^{-/} μL)	415	-11.94	-13.55, -10.31	-2.88	-3.27, -2.48		
	1000	-29.14	-33.97, -24.59	-2.91	-3.40, -2.46		

All Sites Combined - Estimated Bias at Critical Points (Continued)

NA = Not applicable since the critical point is zero.

Maagunand	Critical Dainta]	Bias	%Bias			
Measurand	Critical Points	Estimate	95% CI	Estimate	95% CI		
	6.40	-0.07	-0.23, 0.06	-1.16	-3.61, 0.93		
MPV (IL)	11.0	-0.34	-0.38, -0.30	-3.06	-3.44, -2.68		
	1.00	0.40	0.35, 0.45	39.80	34.65, 45.18		
%fP(%)	7.00	-0.92	-1.03, -0.80	-13.17	-14.70, -11.50		
MCHr (pg)	29.0	1.30	1.11, 1.46	4.47	3.83, 5.04		
	34.5	1.78	1.62, 1.95	5.17	4.70, 5.64		

All Sites Combined - Estimated Bias at Critical Points (Continued)

NA = Not applicable since the critical point is zero.

2. Sensitivity and Specificity

The Sensitivity and Specificity study was performed based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document CLSI H20-A2³.

Sensitivity and specificity performance of the Alinity h-series System was assessed by identifying distributional abnormalities and morphological flags using blood films from venous and capillary specimens collected in K₂EDTA tubes.

For this analysis, separate 2x2 tables were constructed in order to determine sensitivity for both morphological and distributional abnormalities. The sample size (N) and numbers of true positives (TP), false positives (FP), true negatives (TN), false negatives (FN), sensitivity, specificity, and efficiency are presented in the table below. All results were within the predefined acceptance criteria and found to be acceptable.

	Ν	ТР	FP	FN	TN	Sensitivity (95% CI)	Specificity (95% CI)	Efficiency (95% CI)
Any Morphological Flags and/or Distributional Abnormalities	674	255	84	63	272	80.19% (75.38%, 84.43%)	76.40% (71.64%, 80.72%)	78.19% (74.88%, 81.25%)

All Sites Combined – Sensitivity and Specificity

³ Clinical and Laboratory Standards Institute (CLSI). Reference Leukocyte (WBC) Differential Count (Proportional) and Evaluation of Instrumental Methods; Approved Standard—Second Edition. CLSI Document H20-A2. Wayne, PA: CLSI; 2007.

3. Precision (Repeatability)

Precision was performed based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document H26-A2⁴. Samples from 32 unique donors (16 normal and 16 pathological around medical decision points) were collected in K₂EDTA tubes and tested in a minimum of 32 consecutive replicates for normal samples and a minimum of 10 consecutive replicates for pathological samples. The mean, standard deviation (SD), coefficient of variation (CV), and 95% CI were calculated for each parameter. All results met the predefined acceptance criteria and were determined to be acceptable.

4. System Reproducibility

Reproducibility was performed based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document EP05-A3⁵.

The reproducibility study was performed at three clinical sites using a single lot of Alinity h-series 29P Control (low, normal, high). Each control level was tested for 5 days with 3 runs per day and in a minimum of 2 replicates per run. The within-laboratory %CV or SD for each control level were calculated and presented in the table below. All results met the predefined acceptance criteria and were determined to be acceptable.

⁴ Clinical and Laboratory Standards Institute (CLSI). Validation, Verification, And Quality Assurance Of Automated Hematology Analyzers; Approved Standard - Second Edition. CLSI Document H26-A2. Wayne, PA: CLSI; 2010.

⁵ Clinical and Laboratory Standards Institute (CLSI). Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. CLSI Document EP05-A3. Wayne, PA: CLSI; 2014.

	T. com	N	Maar	Repea	tability	Betwee	en-Run	Betwe	en-Day	Within-L	aboratoryª	Between	n-Device	Reprod	ucibility
Measurand	Level	IN	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
	Low	84	3.04	0.068	2.22	0.000	0.00	0.028	0.93	0.073	2.41	0.027	0.87	0.078	2.56
WBC (X 10 ³ /µL)	Normal	84	7.18	0.138	1.92	0.000	0.00	0.000	0.00	0.138	1.92	0.045	0.63	0.145	2.02
	High	84	16.12	0.185	1.15	0.000	0.00	0.000	0.00	0.185	1.15	0.000	0.00	0.185	1.15
	Low	84	1.40	0.052	3.73	0.014	0.98	0.000	0.00	0.054	3.86	0.034	2.46	0.064	4.58
NEU (X 10 ³ /µL)	Normal	84	3.45	0.097	2.80	0.015	0.42	0.000	0.00	0.098	2.83	0.060	1.72	0.115	3.32
	High	84	8.15	0.187	2.29	0.000	0.00	0.000	0.00	0.187	2.29	0.102	1.25	0.213	2.61
	Low	84	0.88	0.035	4.02	0.000	0.00	0.008	0.93	0.036	4.13	0.012	1.35	0.038	4.34
LYM (X 10 ³ /µL)	Normal	84	1.84	0.054	2.92	0.000	0.00	0.000	0.00	0.054	2.92	0.000	0.00	0.054	2.92
	High	84	3.57	0.082	2.29	0.000	0.00	0.000	0.00	0.082	2.29	0.000	0.00	0.082	2.29
	Low	84	0.33	0.024	7.18	0.000	0.00	0.005	1.65	0.024	7.36	0.007	2.15	0.025	7.67
MONO (X 10 ³ /µL)	Normal	84	0.80	0.051	6.35	0.000	0.00	0.000	0.00	0.051	6.35	0.000	0.00	0.051	6.35
	High	84	1.84	0.073	4.00	0.047	2.57	0.000	0.00	0.087	4.75	0.000	0.00	0.087	4.75
	Low	84	0.07	0.010	13.91	0.000	0.00	0.002	2.76	0.010	14.18	0.000	0.00	0.010	14.18
EOS (X 10 ³ /µL)	Normal	84	0.20	0.014	7.20	0.006	2.77	0.000	0.00	0.015	7.72	0.000	0.00	0.015	7.72
	High	84	0.49	0.025	5.16	0.010	2.13	0.000	0.00	0.027	5.59	0.003	0.55	0.027	5.61

All Sites Combined – Reproducibility Study

NA=Not Applicable; %CVs are not meaningful when measurand result approaches zero. ^a Reproducibility contains repeatability, between-run, between-day, and between-device variance components.

Manager	Land	N	Mean	Repea	tability	Between-Run		Between-Day		Within-Laboratory ^a		Between-Device		Reproducibility	
Measurand	Level	IN	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
	Low	84	0.03	0.008	24.23	0.005	14.55	0.001	2.04	0.009	28.34	0.000	0.00	0.009	28.34
BASO (X 103/µL)	Normal	84	0.06	0.012	18.37	0.000	0.00	0.000	0.00	0.012	18.37	0.000	0.00	0.012	18.37
	High	84	0.15	0.017	10.71	0.007	4.57	0.000	0.00	0.018	11.64	0.003	2.14	0.018	11.84
	Low	84	0.33	0.024	7.40	0.000	0.00	0.000	0.00	0.024	7.40	0.010	2.96	0.026	7.97
IG (X 10 ³ /µL)	Normal	84	0.82	0.047	5.71	0.000	0.00	0.013	1.62	0.049	5.94	0.017	2.06	0.052	6.29
	High	84	1.92	0.079	4.09	0.061	3.16	0.000	0.00	0.099	5.17	0.059	3.08	0.116	6.02
	Low	84	0.00	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA
NRBC (X 103/µL)	Normal	84	0.00	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA
	High	84	2.34	0.062	2.64	0.000	0.00	0.000	0.00	0.062	2.64	0.016	0.69	0.064	2.73
	Low	84	2.67	0.018	0.69	0.000	0.00	0.006	0.23	0.019	0.73	0.000	0.00	0.019	0.73
RBC (X 106/µL)	Normal	84	4.13	0.028	0.68	0.000	0.00	0.000	0.00	0.028	0.68	0.019	0.46	0.034	0.82
	High	84	5.16	0.047	0.92	0.000	0.00	0.008	0.15	0.048	0.93	0.085	1.65	0.098	1.90
	Low	84	7.11	0.041	0.57	0.017	0.25	0.000	0.00	0.044	0.62	0.048	0.68	0.065	0.92
HGB (g/dL)	Normal	84	11.34	0.066	0.58	0.000	0.00	0.000	0.00	0.066	0.58	0.048	0.42	0.082	0.72
	High	84	16.37	0.094	0.57	0.017	0.11	0.020	0.12	0.097	0.60	0.000	0.00	0.097	0.60

All Sites Combined – Reproducibility Study (Continued)

NA=Not Applicable; %CVs are not meaningful when measurand result approaches zero.

Measurand Level		N	Maan	Repeat	ability	Betwee	en-Run	Betwee	en-Day	Within-La	boratoryª	Betwee	n-Device	Reprodu	cibility
Measurand	Level	1	Ivicali	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
	Low	84	22.64	0.179	0.79	0.000	0.00	0.107	0.47	0.209	0.92	0.062	0.27	0.218	0.96
HCT (%)	Normal	84	35.88	0.269	0.75	0.000	0.00	0.038	0.11	0.272	0.76	0.244	0.68	0.366	1.02
	High	84	49.80	0.468	0.94	0.000	0.00	0.155	0.31	0.493	0.99	0.897	1.80	1.023	2.05
	Low	84	84.68	0.120	0.14	0.185	0.22	0.248	0.29	0.332	0.39	0.320	0.38	0.461	0.54
MCV (fL)	Normal	84	86.93	0.212	0.24	0.153	0.18	0.173	0.20	0.313	0.36	0.247	0.28	0.399	0.46
	High	84	96.48	0.092	0.10	0.116	0.12	0.270	0.28	0.308	0.32	0.201	0.21	0.368	0.38
	Low	84	26.58	0.231	0.87	0.076	0.29	0.069	0.26	0.253	0.95	0.154	0.58	0.296	1.11
MCH (pg)	Normal	84	27.48	0.240	0.88	0.000	0.00	0.000	0.00	0.240	0.88	0.041	0.15	0.244	0.89
	High	84	31.72	0.344	1.08	0.000	0.00	0.075	0.24	0.352	1.11	0.495	1.56	0.607	1.92
	Low	84	31.39	0.292	0.93	0.000	0.00	0.154	0.49	0.330	1.05	0.251	0.80	0.415	1.32
MCHC (g/dL)	Normal	84	31.61	0.310	0.98	0.000	0.00	0.000	0.00	0.310	0.98	0.150	0.47	0.344	1.09
	High	84	32.88	0.363	1.10	0.000	0.00	0.112	0.34	0.380	1.16	0.563	1.71	0.679	2.07
	Low	84	13.26	0.084	0.63	0.000	0.00	0.033	0.25	0.090	0.68	0.670	5.05	0.676	5.10
RDW (%)	Normal	84	13.32	0.096	0.72	0.000	0.00	0.021	0.16	0.098	0.74	0.705	5.29	0.711	5.34
	High	84	12.37	0.049	0.39	0.031	0.25	0.035	0.29	0.068	0.55	0.621	5.02	0.625	5.05

All Sites Combined – Reproducibility Study (Continued)

NA=Not Applicable; %CVs are not meaningful when measurand result approaches zero.

Measurand Level		N	Maan	Repea	tability	Betwe	en-Run	Betwo	en-Day	Within-La	aboratoryª	Between	1-Device	Reprod	lucibility
Measurand	Level	IN	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
	Low	84	216.24	3.142	1.45	0.000	0.00	1.132	0.52	3.340	1.54	4.535	2.10	5.632	2.60
RETIC (X 10 ³ /µL)	Normal	84	153.55	3.478	2.26	0.805	0.52	1.374	0.89	3.825	2.49	7.441	4.85	8.367	5.45
	High	84	132.30	4.172	3.15	0.000	0.00	0.000	0.00	4.172	3.15	10.734	8.11	11.517	8.71
	Low	84	0.33	0.011	3.30	0.000	0.00	0.008	2.52	0.014	4.15	0.050	15.31	0.052	15.86
IRF	Normal	84	0.27	0.011	4.23	0.000	0.00	0.005	2.02	0.013	4.69	0.028	10.30	0.031	11.32
	High	84	0.19	0.009	4.56	0.003	1.48	0.000	0.00	0.009	4.80	0.004	2.28	0.010	5.31
	Low	84	73.78	1.977	2.68	0.000	0.00	0.508	0.69	2.041	2.77	1.317	1.79	2.429	3.29
PLT (X 10 ³ /µL)	Normal	84	225.19	3.454	1.53	1.380	0.61	0.216	0.10	3.726	1.65	1.257	0.56	3.932	1.75
	High	84	478.12	7.653	1.60	0.000	0.00	2.889	0.60	8.180	1.71	3.997	0.84	9.105	1.90
	Low	84	9.55	0.070	0.73	0.000	0.00	0.000	0.00	0.070	0.73	0.022	0.23	0.073	0.77
MPV (fL)	Normal	84	9.59	0.041	0.43	0.000	0.00	0.008	0.08	0.042	0.44	0.022	0.23	0.047	0.49
	High	84	9.60	0.031	0.32	0.021	0.22	0.000	0.00	0.037	0.39	0.005	0.05	0.038	0.39
	Low	84	9.77	0.439	4.49	0.000	0.00	0.030	0.31	0.440	4.50	0.000	0.00	0.440	4.50
%rP(%)	Normal	84	8.82	0.140	1.59	0.095	1.08	0.017	0.20	0.170	1.93	0.000	0.00	0.170	1.93
	High	84	9.00	0.105	1.16	0.000	0.00	0.055	0.62	0.119	1.32	0.000	0.00	0.119	1.32

All Sites Combined – Reproducibility Study (Continued)

NA=Not Applicable; %CVs are not meaningful when measurand result approaches zero.

Magurand Level		N	Mean	Repea	tability	Betwee	en-Run	Betwe	en-Day	Within-La	aboratoryª	Between	n-Device	Reprod	ucibility
Measurand	Level	IN	Mean	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
	Low	84	1.04	0.251	24.03	0.154	14.79	0.032	3.03	0.296	28.38	0.000	0.00	0.296	28.38
%B(%)	Normal	84	0.90	0.164	18.32	0.000	0.00	0.000	0.00	0.164	18.32	0.000	0.00	0.164	18.32
	High	84	0.96	0.100	10.39	0.049	5.09	0.000	0.00	0.111	11.57	0.020	2.06	0.113	11.75
	Low	84	10.78	0.780	7.24	0.000	0.00	0.000	0.00	0.780	7.24	0.450	4.17	0.900	8.35
%IG(%)	Normal	84	11.46	0.572	4.99	0.175	1.53	0.078	0.68	0.603	5.26	0.322	2.81	0.684	5.97
	High	84	11.92	0.465	3.90	0.378	3.17	0.000	0.00	0.599	5.03	0.379	3.18	0.709	5.95
	Low	84	0.00	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA
NR/W	Normal	84	0.00	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA	0.000	NA
	High	84	14.53	0.384	2.64	0.000	0.00	0.035	0.24	0.385	2.65	0.093	0.64	0.396	2.73
	Low	84	8.09	0.120	1.49	0.000	0.00	0.000	0.00	0.120	1.49	0.173	2.14	0.211	2.60
%R(%)	Normal	84	3.72	0.078	2.10	0.027	0.73	0.026	0.70	0.087	2.33	0.198	5.32	0.216	5.81
	High	84	2.57	0.084	3.28	0.000	0.00	0.000	0.00	0.084	3.28	0.254	9.90	0.268	10.43

All Sites Combined – Reproducibility Study (Continued)

NA=Not Applicable; %CVs are not meaningful when measurand result approaches zero.

5. Linearity

Linearity was evaluated based on guidance from the Clinical and Laboratory Standards Institute (CLSI) document EP06 2nd edition⁶.

Linearity for RBC, HGB, and NRBC was determined using whole blood to span the analytical measuring interval of each measurand. Linearity for WBC, PLT, and RETIC was determined using commercially available linearity kits. For each measurand, the testing minimally included:

- 9 levels
- 4 replicates of each level
- 1 instrument
- 1 set of reagent lots

A weighted linear regression analysis was used to assess linearity for each measurand. Results are presented in the table below. All results met the predefined acceptance criteria and were determined to be acceptable.

Measurand	Linear Range
RBC	$0.00 - 8.08 \text{ x } 10^{6} / \mu L$
HGB	0.04 - 24.14 g/dL
NRBC	$0.00 - 26.10 \ge 10^3/\mu L$
WBC	$0.00 - 449. \ x \ 10^{3}/\mu L$
PLT	$0.06 - 5325 \ x \ 10^3/\mu L$
RETIC	$0.05 - 644. \ge 10^3/\mu L$

⁶ Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures*. 2nd edition. CLSI Document EP06. Wayne, PA: CLSI; 2020.

6. Carryover

Alinity hq susceptibility to potential carryover was evaluated based on guidance from the Clinical Laboratory and Standards Institute (CLSI) document H26-A2.⁷ Venous whole blood specimens were collected in K₂EDTA tubes. For each measurand, a minimum of 4 carryover runs was completed at each of 4 testing sites for a minimum of 16 total carryover runs per measurand, where each run consisted of testing a high target specimen in 3 replicates followed by testing a low target specimen in 3 replicates. All results met the predefined acceptance criteria and were determined to be acceptable.

7. Potentially Interfering Substances Study

The susceptibility of the Alinity h-series System to interference in the presence of hemoglobin, triglycerides, bilirubin, cholesterol, elevated WBCs, elevated RBCs, elevated PLTs, and microcytic RBCs was tested in samples collected in K₂EDTA tubes. Results are presented in the table below. All results met the predefined acceptance criteria and were determined to be acceptable.

Interferent	Result Level	Measurand Impacted
Hemoglobin	1.0 g/dL	WBC, RBC, MCV, PLT, RETIC
Trightanida	1.15 g/dL	WBC, RBC, HGB, MCV, RETIC
Ingrycende	0.63 g/dL	PLT
Dilimbin unconjugated	0.080 g/dL	WBC, RBC, HGB, PLT, RETIC
Binrubin - unconjugated	0.040 g/dL	MCV
Bilirubin - conjugated	0.080 g/dL	WBC, RBC, HGB, MCV, PLT, RETIC
Chalasteral	0.40 g/dL	WBC, RBC, HGB, MCV, RETIC
Cholesteroi	0.50 g/dL	PLT
Elevated WBCs	99.0 x 10 ³ cells/μL	RBC, PLT, HGB, MPV
Elevated RBCs	No interference was observed acro	oss the measuring range
Elevated PLTs	2840 x 10 ³ cells/μL	WBC, RBC, HGB, MPV
Microcytic RBCs	Microcytosis (MCV < 57 fL)	PLT

⁷ Clinical and Laboratory Standards Institute (CLSI). Validation, Verification, And Quality Assurance Of Automated Hematology Analyzers; Approved Standard - Second Edition. CLSI Document H26-A2. Wayne, PA: CLSI; 2010.

VIII. Other Supportive Instrument Performance Data

1. Limits of Blank, Detection, and Quantitation (LoB, LoD, and LoQ)

Limits of Blank, Detection, and Quantitation were established for the measurands WBC, RBC, HGB, and PLT based on guidance from the Clinical and Laboratory Standards Institute (CLSI) documents EP17-A2⁸ and H26-A2⁹.

Testing was conducted over a minimum of 3 days using a minimum of 2 unique samples per day on each of 2 test selections (CBC+Diff and CBC+Diff+Retic) in a minimum of 5 replicates using each of the 2 sets of reagent lots. The maximum observed limit of blank (LoB), limit of detection (LoD), and limit of quantitation (LoQ) values are summarized in the table below. All results met the predefined acceptance criteria and were determined to be acceptable.

Maagunand	Results								
wieasuranu	LoB	LoD	LoQ						
WBC (x10 ³ /µL)	0.01	0.02	0.03						
RBC (x10 ⁶ /µL)	0.00	0.01	0.01						
HGB (g/dL)	0.08	0.11	0.05						
PLT (x10 ³ /μL)	0.15	0.38	0.29						

Limits of Blank, Detection, and Quantitation (LoB, LoD, and LoQ)

⁸ Clinical and Laboratory Standards Institute (CLSI). Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.

⁹ Clinical and Laboratory Standards Institute (CLSI). Validation, Verification, And Quality Assurance Of Automated Hematology Analyzers; Approved Standard - Second Edition. CLSI Document H26-A2. Wayne, PA: CLSI; 2010.

2. Specimen Stability

For venous specimen stability, a minimum of 10 abnormal and 10 normal venous whole blood specimens K₂EDTA and K₃EDTA tubes. Each specimen was tested in a minimum of 2 replicates. For capillary specimen stability, a minimum of 20 normal whole blood specimens were collected in K₂EDTA and K₃EDTA tubes. Each specimen was tested in a minimum of 1 replicate. All samples were tested within 4 hours (baseline) of specimen collection. Samples stored at Refrigerated Temperature (2-8°C) were tested at up to 24 hours after specimen collection. Samples stored at Room Temperature (18-26°C) were tested at up to 48 hours after specimen collection. The results were used to support the information provided in the system labeling for venous and capillary specimen stability.

3. Anticoagulant Comparability (K₃EDTA versus K₂EDTA)

Anticoagulant Comparability (K₃EDTA versus K₂EDTA) was evaluated based on guidance from CLSI EP35 1^{st} ed¹⁰. A total of 199 unique adult and pediatric donor sets covering relevant medical decision levels and reference ranges and spanning the analytical measurement ranges to the extent possible were tested in 2 replicates for each measurand. The performance between the anticoagulant tube type (K₃EDTA) and anticoagulant tube type (K₂EDTA) was compared.

Comparability between the anticoagulants was assessed based on the mean difference or % difference and a regression analysis using either a Passing-Bablok or Deming regression model. All reportable parameters that were evaluated met their predefined bias acceptance criteria.

4. Microtainer Capillary versus Microtube for Automated Process (MAP)

Comparability between the K₂EDTA Microtainer Capillary tube versus K₂EDTA Microtainer Microtube for Automated Process (MAP) was evaluated

¹⁰ Clinical and Laboratory Standards Institute (CLSI). Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures. 1st ed. CLSI Guideline EP35-A. Wayne, PA: CLSI; 2019.

based on guidance from CLSI document EP35 1st ed¹¹. A total of 44 unique donor sets (normal whole blood specimens) were collected in K₂EDTA Microtainer Capillary and Microtainer Microtube for Automated Process (MAP) blood collection tubes. Each specimen was tested in 1 replicate in the open-tube processing mode for each measurand.

Comparability between the capillary tube types was assessed based on the mean difference or % difference and a regression analysis using either a Passing-Bablok or Deming regression model. All reportable parameters that were evaluated met their predefined bias acceptance criteria.

5. Matrix Comparability (Capillary versus Venous)

Matrix Comparability (Capillary versus Venous) was evaluated based on guidance from CLSI EP35 1st ed.¹² A total of 76 unique venous and capillary donor sets (normal and abnormal whole blood specimens) were collected in Microtainer Microtube for Automated Process (MAP) Microtubes (capillary specimens) and standard K₂EDTA tubes (venous specimens). Each specimen was tested in 2 replicates for each measurand.

Comparability between capillary and venous matrices was assessed based on the mean difference or % difference and a regression analysis using either a Passing-Bablok or Deming regression model. All reportable parameters that were evaluated met their predefined bias acceptance criteria.

6. Sample/Tube Processing Mode Comparability (Open Mode versus Closed Mode)

Sample processing mode comparability was evaluated based on guidance from CLSI EP35 1st ed.¹¹ A total of 226 unique venous specimens covering relevant medical decision levels and reference ranges and spanning the analytical

¹¹ Clinical and Laboratory Standards Institute (CLSI). Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures. 1st ed. CLSI Guideline EP35-A. Wayne, PA: CLSI; 2019.

¹² Clinical and Laboratory Standards Institute (CLSI). Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures. 1st ed. CLSI Guideline EP35. Wayne, PA: CLSI; 2019.

measurement ranges to the extent possible were collected in K₂EDTA tubes. Each specimen was tested in 2 replicates in the closed-tube and open-tube processing modes for each measurand.

Comparability between the sample/tube processing modes was assessed based on the mean difference or % difference and a regression analysis using either a Passing-Bablok or Deming regression model. All reportable parameters that were evaluated met their predefined bias acceptance criteria.

7. Reference Intervals (Expected Values)

The study was performed based on guidance from the Clinical Laboratory and Standards Institute (CLSI) document EP28-A3c¹³ to establish adult (> 21 years old) reference intervals for male and female populations and pediatric (\leq 21 years old) reference intervals for all subgroups (neonate, infant, child, and adolescent). Reference intervals were established by evaluating venous or capillary whole blood specimens collected in K₂EDTA tubes from apparently healthy subjects.

A total of 261 unique venous and 1 capillary whole blood specimens collected from 126 male and 136 female adult subjects were tested in a minimum of 1 replicate to establish adult reference intervals. A total of 360 venous or capillary specimens from pediatric sub-populations: 61 neonates (birth to 1 month); 68 infant (> 1 month to 2 years old), 109 child (> 2 to 12 years old), and 122 adolescents (> 12 to 21 years old) were tested in 1 replicate to establish pediatric reference intervals for each measurand.

¹³ Clinical and Laboratory Standards Institute (CLSI). Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory, 3rd Ed CLSI Guideline EP28-A3c. Wayne, PA: CLSI; 2019.

IX. Conclusion

The results presented in this 510(k) Pre-market Notification demonstrate that the performance of the subject device, Alinity h-series System, is substantially equivalent to the predicate device, Sysmex® XN-Series (XN-10, XN-20).

The similarities and differences between the subject device and predicate device are presented in Section 5-VI.

There is no known potential adverse effect to the operator when using the subject device, Alinity h-series System according to its Operations Manual.