

May 17, 2022

Tianjin Huahong Technology Co., Ltd. % Stuart Situ, Director Landlink Healthcare Technology (Shanghai) Co., Ltd Room 1308, Baohua International Plaza, West Guangzhong Road 555, Jingan District, Shanghai, 200072 China

Re: K220372

Trade/Device Name: Heel Stick Safety Lancet (HHZ-II, HHZ-III) Regulation Number: 21 CFR 878.4800 Regulation Name: Manual Surgical Instrument For General Use Regulatory Class: Class II Product Code: FMK Dated: January 28, 2022 Received: February 9, 2022

Dear Stuart Situ:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Long Chen, Ph.D. Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K220372

Device Name

Heel Stick Safety	Lancet	(HHZ-II,	HHZ-III)
-------------------	--------	----------	----------

Indications for Use (Describe)

Heel Stick Safety Lancet is intended for the collection of capillary blood from the heel of newborn, preemie, and toddler. The lancet has equipped with safety protection features.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) summary-K220372

I Submitter

Tianjin Huahong Technology Co., Ltd. A01, Plant B, No.278, Hangkong Road, Tianjin Pilot Free Trade Zone (Air Port Industrial Park), 300308 Tianjin, China

Establishment Registration Number: 3009498536

Contact person: Ms. Ying Yuan Quality Manager Tel.: +86-18622179097 E-mail: ying.yuan@hh-technology.com

Submission Correspondent:

Ms. Stuart Situ Landlink Healthcare Technology (Shanghai) Co., Ltd. E-mail: <u>stuart.situ@landlink-healthcare.com</u>

Preparation date: Jan 28, 2022

II Proposed Device

Trade Name of Device:	Heel Stick Safety Lancet (HHZ-II,HHZ-III)	
Common name:	Single Use Only Blood Lancet With An Integral Sharps Injury	
	Prevention Feature	
Regulation Number:	21 CFR 878.4850	
Regulatory Class:	Class II	
Product code:	FMK	
Review Panel	General & Plastic Surgery	

III Predicate Devices

510(k) Number:	K210745
Trade name:	Heel Incision Safety Lancet
Classification:	Class I
Product Code:	FMK
Manufacturer	SteriLance Medical(Suzhou) Inc.

IV Device description

The heel stick safety lancet is offered in various depth size (0.65mm, 0.85mm, 1.00mm, 1.50mm) and width size (1.50 mm, 1.75mm, 2.50 mm, 2.80 mm).

It consists of blade base, blade, which is welded with the blade base, spring, top head, bottom, and pushing button. The top head and bottom are intended to provide physical protection to the blade, and PET blister along with the tyvek paper cover the primary sterile barrier system for the device. The product is individually primarily packaged and sterilized to a SAL of 10⁻⁶ by radiation sterilization. It is intended for single use only. The shelf-life of the product is 5 years.

Heel stick safety lancet is intended to be used by professionals. The product is intended for prescription (Rx) only.

V Indication for use

Heel Stick Safety Lancet is intended for the collection of capillary blood from the heel of newborn, preemie, and toddler. The lancet has equipped with safety protection features.

VI Comparison of technological characteristics with the predicate devices

The comparison and discussion between the subject device and the predicate devices are listed in below table 1:

ltem	Proposed device	Predicate device	Discu
	(K220372)	(K210745)	ssion
Product name	Heel Stick Safety Lancet (HHZ-II,HHZ-III)	Heel Incision Safety Lancet	Same
Product Code	FMK	FMK	Same
Regulation No.	21 CFR § 878.4850	21 CFR § 878.4800	Simila r ¹
Class	II	I	Simila r ¹
Prescription/ over-the-cou nter use	Prescription Use	Prescription Use	Same

Table 1 General Comparison of Heel Stick Safety Lancet

Indication for use	Heel Stick Safety Lancet is intended for the collection of capillary blood from the heel of newborn, preemie, and toddler. The lancet has equipped with safety protection features	Heel Incision Safety Lancet is intended for the collection of capillary blood from the heel of newborn, preemie, and toddler. The lancet has equipped with safety protection features	Same
Applicable user	Newborn, preemie, and toddler	Newborn, preemie, and toddler	Same
Safety protection features	Yes	Yes	Same
Reuse durability	Single use	Single use	Same
Sterilization method and SAL	Sterilized by Radiation SAL=10 ⁻⁶	Sterilized by Radiation SAL=10 ⁻⁶	Same
Self-life	5 years	5 years	Same
Component	 Blade Blade base Top head Bottom Pushing button Spring 	 Triggering button Safety button/Screw button Spring Cam Lancet core Shell Blade Swing arm 	Simila r ²
Incision	0.65*1.50 mm	0.65*1.40 mm	Simila
Depth*Width	0.85*1.75 mm	0.85*1.75 mm	r ³
	1.00*2.50 mm	1.00*2.50 mm	
	1.50*2.80 mm	1.14*2.80 mm	
		2.00*3.00 mm	

Materials of	Blade: stainless steel	Blade:304 stainless steel	Same
parts in contact with	Top head: ABS	Shell: ABS	
human body	Bottom: ABS	Triggering button: ABS	
	Pushing bottom: ABS	Safety button: ABS	
Biocompatibil	Conforms to the	Conforms to the	Same
ity	requirements of ISO 10993	requirements of ISO 10993	
	series standards.	series standards.	
Label/Labelin	Complied with 21 CFR part	Complied with 21 CFR part	Same
g	801	801	

¹ The classification and Regulation number are different because FDA issued the final order about Reclassification of Blood Lancets on 11/22/2021.

² The component of proposed devices is different from the predicate devices. However, the performance test for proposed device has been conducted as same as predicate device. The Safety feature performance specification for predicate device is a little difference. However, the safety feature performance test for proposed device has been evaluated and the test result conforms to requirements of ISO 23908:2011 standards. Therefore, the difference on configuration and materials does not affect substantially equivalence.

³ The Incision depth and width of proposed devices is different from the predicate devices. However, the performance test for proposed device has been conducted as same as predicate device. And the depth and width of the proposed devices is covered by the predicated devices. Therefore, the difference on configuration and materials does not affect substantially equivalence.

VII Non-Clinical Testing

The bench testing performed verifies that the proposed device is as safe, as effective, and performs as well as the legally marketed predicate device in terms of critical performance characteristics. These tests are as follow.

Items	Acceptance criteria Results		
	The product color shall be correct,	Meet	the
Appearance	with the surface without burr and	requirement	
	scratches		
Cleanness	The product shall be cleaned	Meet	the
Cleanness	The product shall be cleaned	requirement	
Dimension	HHZ-II :	Meet	the

	length: 32.0 mm ± 1.0 mm; width:	requirement	
	11.5 mm ± 0.5 mm		
	HHZ-III :		
	length: 36.65mm \pm 0.5mm; width:		
	11.6±0.5mm		
Firmness	Blade and blade base shall be firmly	Meet the	е
1 111111055	connected.	requirement	
	Launch performance should be	Meet th	е
Launch performance	good, launch button press	requirement	
	smoothly, no jam		
Launch depth and Launch	Launch depth and launch width	Meet the	е
width	shall meet the requirements.	requirement	
Force to activate the	4 - 10N	Meet the	е
safety feature	4 - TON	requirement	
Testing access to the	The blade shall not touch the	Meet the	е
sharp in safe mode	sphere.	requirement	

Biocompatibility Testing:

The biocompatibility evaluations were conducted in accordance with the 2020 FDA Guidance document Use of International Standard ISO 10993-1 "Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within a Risk Management Process". The tests include the following tests:

Item	Test method	Test results
In Vitro Cytotoxicity	ISO 10993-5: 2009	No Cytotoxicity
Skin Sensitization	ISO 10993-10: 2010	No Skin sensitization
Intracutaneous reactivity	ISO 10993-10: 2010	No irritation
Acute Systemic Toxicity	ISO 10993-11: 2017	No Acute Systemic
		Toxicity
Pyrogenicity	ISO 10993-11: 2017	No thermogenic reaction

Simulated Clinical Use

A simulated clinical use study was performed on 600 device samples each for the Heel Stick Safety Lancet according to FDA Guidance, Guidance for Industry and FDA Staff: Medical Device with Sharps Injury Prevention Feature, issued on August 9, 2005 and ISO 23908 to evaluate the safety mechanism of the proposed device. The results demonstrated that the proposed device met the pre-established criteria.

VIII Clinical Testing

No clinical study is included in this submission.

IX Conclusion

The proposed device has the same indication for use and has similar design features and technological characteristic as the predicate device. Performance testing data demonstrates that the proposed device is safety and effectiveness as the predicated device. Accordingly, the proposed device is substantially equivalent to the predicate device.