

April 27, 2022

Venus Concept USA Inc. William Mcgrail VP, Global Regulatory Affairs & Quality Assurance 1880 N Commerce Pkwy, Suite 2 Weston, Florida 33326

Re: K220592

Trade/Device Name: Venus BlissMAX Regulation Number: 21 CFR 878.4400 Regulation Name: Electrosurgical Cutting And Coagulation Device And Accessories Regulatory Class: Class II Product Code: PBX, PKT, NGX Dated: February 24, 2022 Received: March 1, 2022

Dear William McGrail:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Purva Pandya, D.Eng. Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K220592

Device Name Venus BlissMAX

Indications for Use (Describe)

The Venus BlissMAX device is a diode laser system intended for non-invasive lipolysis of the abdomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less.

In addition, the Venus BlissMAX device is intended for the treatment of the following medical conditions; using the MP2 applicator for delivery of RF energy combined with massage and magnetic field pulses:

- Relief of minor muscle aches and pain, relief of muscle spasm
- Temporary improvement of local blood circulation
- Temporary reduction in the appearance of cellulite

In addition, The Venus BlissMAX device using the FlexMAX applicators is intended for muscle conditioning to stimulate healthy muscles.

The Venus BlissMAX device using the FlexMAX applicators is not intended to be used in conjunction with therapy or treatment of medical diseases or medical conditions of any kind.

The Venus BlissMAX device using the FlexMAX applicators is intended to be operated by a trained professional.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K220592

510(k) Summary (As Required by 21.CFR.807.92)

1. SUBMITTER	
Manufacturer:	Venus Concept USA, Inc.
	1880 N Commerce Pkwy, Suite 2
	Weston, FL 33326, USA
Contact Person:	William H. McGrail
	Vice President, Global RA & QA
	Venus Concept
	Phone: (978) 808-0420
	Email: bmcgrail@venusconcept.com
Date Prepared:	April 27, 2022

2. DEVICE INFORMATION

Trade/Device Name(s):	Venus BlissMAX System
Common or Unusual Name:	Venus BlissMAX System
Regulation Number:	21 CFR § 878.4400 21 CFR § 878.5400 21 CFR § 890.5850
Classification Name:	Electrosurgical Cutting And Coagulation Device And Accessories
<u>Regulation Class:</u> <u>Product Code:</u>	Class II PKT, PBX, NGX,
Review Panel:	General & Plastic Surgery

3. PREDICATE DEVICES

Primary Predicate:	Venus Concept Venus BlissMAX (K213308)
Co-Secondary Predicate:	Cynosure SculpSure (K182741)
Co-Secondary Predicate:	Venus Concept Venus Legacy Pro (K191528)

4. INDICATIONS FOR USE

The Venus BlissMAX device is a diode laser system intended for non-invasive lipolysis of the abdomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less.

In addition, the Venus BlissMAX device is intended for the treatment of the following medical conditions; using the MP2 applicator for delivery of RF energy combined with massage and magnetic field pulses:

- Relief of minor muscle aches and pain, relief of muscle spasm
- Temporary improvement of local blood circulation
- Temporary reduction in the appearance of cellulite

In addition, The Venus BlissMAX device using the FlexMAX applicators is intended for muscle conditioning to stimulate healthy muscles.

The Venus BlissMAX device using the FlexMAX applicators is not intended to be used in conjunction with therapy or treatment of medical diseases or medical conditions of any kind.

The Venus BlissMAX device using the FlexMAX applicators is intended to be operated by a trained professional

5. DEVICE DESCRIPTION

The Venus BlissMAX device is a computerized system comprised of a system console (main unit), four (4) Diode Laser applicators, one (1) MP² (RF+ PEMF+ Vacuum) applicator and four (4) FlexMAX (EMS) applicators. The system delivers laser, bipolar RF and biphasic electrical energies, vacuum pressure, and pulsed electromagnetic fields (PEMF) to the skin and the underlying tissues of the treatment area. The device provides individual adjustment of laser power, EMS intensity level, and RF power, in addition to vacuum levels, for each patient.

The console of the Venus BlissMAX device contains a power supply unit, Laser, RF, and EMS controllers, (power modules, on main board), a suction module (vacuum), a controller unit (on main board), Laser water cooling system (power module, on main board), a touch- screen user interface and display panel.

The applicators are connected to the console via a cable. The RF applicator is comprised of various combinations of RF electrodes, magnetic coils, and vacuum conduits.

The Laser applicators are comprised of a light guide, touch sensors and light-emitting diodes.

The EMS applicator is comprised of two electrodes and a light indicator.

6. TECHNOLOGICAL CHARACTERISTICS

The Venus BlissMAX device delivers laser energy to the subcutaneous tissue layers via the four diode laser applicators connected to the console. The console utilizes diode laser modules as sources of optical energy and the optical output is fiber-coupled through the applicator to the treatment area. The laser applicators are coupled to the patient's body while using a dedicated belt for the entire treatment. Individual adjustment of the laser output power is provided for each applicator to achieve maximum safety and efficiency for the patient. The laser applicators have an integrated contact skin cooling system to enhance safety and comfort of the treatment.

In addition, the Venus BlissMAX device using the MP2 applicator provides RF treatments combined with emitted magnetic fields and vacuum massaging. The Vacuum is mainly used for the massaging of deep tissues by creating mild to deep suction. The vacuum massage improves the contact surface between electrodes and tissue. The RF currents heat the adipose and muscular tissues to trigger tissue level changes leading to temporary reduction in the appearance, of cellulite and temporary relief of muscle pain and muscle spasm. The RF heating effect also improves local blood circulation in the sub dermal layers. The PEMF assists in achieving treatment effect.

Furthermore, The Venus BlissMAX device using the FlexMAX applicators (EMS applicators) contracts muscles by passing electrical currents through electrodes contacting the affected body area. The transcutaneous electrical current is designed to effect underlying, healthy muscles, causing them to contract. The FlexMAX applicators are coupled to the patient's body while using a dedicated belt during the entire treatment. The belt size and number of EMS applicators is determined by the treatment area and its size. Its use on muscles is in accordance with a class II device Powered Electrical Muscle Stimulator (Product Code NGX), based on the FDA guidance document for powered muscle stimulators for muscle conditioning as a special control¹.

7. TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

The following tables compare the Venus BlissMAX (Subject) device to the primary predicate: Venus BlissMAX (K213308) and co-secondary predicate: Cynosure SculpSure Laser System (K182741) and co-secondary predicate: Venus Legacy Pro (K191528) devices with respect to intended use, technological characteristics, and principles of operation, providing detailed information regarding the basis for the determination of substantial equivalence.

¹ The Guidance for Industry, FDA Reviewers/Staff and Compliance Guidance Document for Powered Muscle Stimulator 510(k)s issued on June 9, 1999

Indications for Use Venus BlissMAX (Subject) Venus RlissMAX (K213308) (Primary Predicate) SculpSure Laser System (K182741) (Co-Secondary Predicate) Substantial Equivalence Discussion Manufacturer Venus Concept Venus Concept Cynosure Inc. Same as Primary Predicate Class II, PKT, 21 CFR 878.5400 Class II, PKT, 21 CFR 878.5400 Class II, PKT, 21 CFR 878.5400 Same as both marking the individual system individuals with a BissMAX with a Body Whats Index (BMI) of 30 or Less. The Venus BlissMAX the Venus BlissMAX and thighs in individuals with a Body with a Body Whats Index (BMI) of 30 or Less. The Venus BlissMAX individuals with a Body or less. Same as both Indications for Use Heat-assisted Heat-assisted Body or less. In addition, the device is is intended for non- invasive lipolysis of the abdoment, flanks, back, and thighs in individuals with a Body with a Body Whats Index (BMI) of 30 or Less. In addition, the device is is intended to affect the appearance of visible fat buices in the abdoment, flanks, back, thighs and submental area. When using the petite musk for non-invasive lipolysis of the automatial area. When using the petite musk for non-invasive lipolysis of fluing 4 Applicators - Sapphire light guides - LED contact sensors Same as both Components Console with GUI Console with GUI 4 Applicators - Sapphire light guides - LED contact sensors Console with GUI 4 Applicators - Sapphire light guides - LED contact sensors						
Manuacturer Venus Concept Venus Concept Cynsure inc. Predicate Device Classification Product Code, Class II, PKT, 21 CFR 878,5400 Class II, PKT, 21 CFR 878,5400 Same as both Regulation The Venus BlissMAX device is a idode laser system intended for non- invasive lipolysis of the abdomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less. The Venus BlissMAX device is a idode laser system individuals with a Body mass Index (BMI) of 30 or less. Same as Co- source invasive lipolysis of the abdomen, flanks, back, diffights in individuals with a Body mass Index (BMI) of 30 or less. Same as Co- source invasive lipolysis of the abdomen, flanks, back, diffights in individuals with a Body eless. Same as Co- source is intended for non- invasive lipolysis of the abdomen, flanks, back, diffights in individuals with a BMI of 49 or less. The device is intended a free the appearance of visible flat buges in the abdomen, flanks, back, thighs and submental area. When using the petite mask for non-invasive lipolysis of the abdomen, flanks, back, thighs and submental area. Same as both Lipolysis for Use Heat-assisted Heat-assisted Heat-assisted Same as both Lipolysis for Use Console with GUI Console with GUI 4 Applicators - Sapphire light guides - LED contact sensors Console with GUI 4 Applicators - Sapphire light guides - LED contact sensors Same as both Synt Size 6 As 6 cm ² on each of the Applicator heads 6 x 6 cm ² on each of t			(K213308)	System (K182741) (Co-Secondary	Equivalence	
Product Code, RegulationClass II, PKT, 21 CFR 878.5400Class II, PKT, 21 CFR 878.5400Same as bothRegulationThe Venus BlissMAX device is a diode laser system intended for non- invasive lipolysis of the adbomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less.The Venus BlissMAX the Venus BlissMAX (BMI) of 30 or less.Same as Co- sceendary Predicate and thighs in individuals with a Body Mass Index (BMI) of 30 or less.The Venus BlissMAX thanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less.Same as DothIndications for UseHeat-assistedHeat-assisted Diode LaserIn addition, the device is intended for non- invasive lipolysis of the abdomen, flanks, back, and thighs in individuals with a Body Mass Index (BMI) of 30 or less.In addition, the device is intended for non- invasive lipolysis of the submental area in individuals with a Body Mass Index, (BMI) of 30 or less.In addition, the device is intended for non- invasive lipolysis of the submental area. When using the petite mask for appearance of lax issue in the submental area.Same as bothLipolysis Method DeliveredHeat-assistedHeat-assistedHeat-assistedSame as bothComponentsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicator feadsSame as bothComponents1064 ± 10nm (infrared)1064 ± 10nm (infrared)1064 ± 10nm (infrared)Same as bothSpot Size 0 Peration6 x 6 cm ² on each of the Applicator heads6 x 6 cm	Manufacturer	Venus Concept	Venus Concept	Cynosure Inc.		
device is a diode laser system intended for non- invasive lipolysis of the abdomen, flanks, back, and thigbs in individuals with a Body with a Body Mass Index (BMI) of 30 or less.device is a diode laser system intended for non- invasive lipolysis of the abdomen, flanks, back, nidividuals with a Body Mass Index (BMI) of 30 or less.intended for non- invasive lipolysis of the abdomen, flanks, back, flanks in individuals with a Body Mass Index (BMI) of 30 or less.intended for non- invasive lipolysis of the abdomen, flanks, back, flanks in individuals with a Body Mass Index (BMI) of 30 or less.intended for non- invasive lipolysis of the submental area in individuals with a BMI of 49 or less. The device is intended to affect the appearance of visible fat bulges in the abdomen, flanks, back, thigbs and submental area, the device can also affect the submental area, the device can also affect the submental area.Same as bothLipolysis MethodHeat-assistedHeat-assistedHeat-assistedSame as bothEnergy Used / DeiveredDiode LaserDiode LaserDiode LaserSame as bothComponentsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsSame as bothSpot Size6x 6 cm ² on each of the Applicator heads1060 ± 20m (infrared)Same as bothSpot Size6x 6 cm ² on each of the Applicator heads6x 6 cm ² on each of the Applicator headsSame as poth	Product Code,		21 CFR 878.5400		Same as both	
MethodHeat-assistedHeat-assistedSame as bothEnergy Used / DeliveredDiode LaserDiode LaserDiode LaserSame as bothComponentsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsSame as bothWavelength1064 ± 10nm (infrared)1064 ± 10nm (infrared)1064 ± 10nm (infrared)Same as bothSpot Size $6 \ge 6 \le m^2$ on each of the Applicator heads $6 \ge 6 \le m^2$ on each of the Applicator heads $4 \ge 6 \le m^2$ on each of the Applicator headsSame as bothOperationCWCWCWSame as both		device is a diode laser system intended for non- invasive lipolysis of the abdomen, flanks, back , and thighs in individuals with a Body Mass Index	device is a diode laser system intended for non-invasive lipolysis of the abdomen and flanks in individuals with a Body Mass Index (BMI) of 30 or	intended for non- invasive lipolysis of the abdomen, flanks, back , and thighs in individuals with a Body Mass Index (BMI) of 30 or less. In addition, the device is intended for non- invasive lipolysis of the submental area in individuals with a BMI of 49 or less. The device is intended to affect the appearance of visible fat bulges in the abdomen, flanks, back, thighs and submental area. When using the petite mask for non-invasive lipolysis of the submental area, the device can also affect the appearance of lax tissue in the submental	secondary Predicate and Primary	
DeliveredDiode LaserDiode LaserSame as bothDeliveredConsole with GUIConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensorsSame as bothWavelength1064 ± 10nm (infrared)1064 ± 10nm (infrared)1064 ± 10nm (infrared)Same as bothSpot Size6 x 6 cm² on each of the Applicator heads6 x 6 cm² on each of the Applicator heads4 x 6 cm² on each of the Applicator headsSame as bothOperationCWCWSame as both	Method	Heat-assisted	Heat-assisted	d Heat-assisted		
ComponentsGUI 4 Applicators - Sapphire light guides - LED contact sensorsConsole with GUI 4 Applicators - Sapphire light guides - LED contact sensors4 Applicators - Sapphire light guides - LED contact sensors3 ame as bothWavelength 1064 ± 10 nm (infrared) 1064 ± 10 nm (infrared) 1060 ± 20 nm (infrared)Same as bothSpot Size $6 \times 6 \operatorname{cm}^2$ on each of the Applicator heads $6 \times 6 \operatorname{cm}^2$ on each of the Applicator heads $4 \times 6 \operatorname{cm}^2$ on each of the Applicator headsSame as Primary predicate DeviceOperationCWCWCWSame as both		Diode Laser	Diode Laser		Same as both	
Spot Size 6 x 6 cm ² on each of the Applicator heads 6 x 6 cm ² on each of the Applicator heads 4 x 6 cm ² on each of the Applicator heads Same as Primary predicate Device Operation CW CW Same as both	Components	GUI 4 Applicators - Sapphire light guides	4 Applicators - Sapphire light guides	4 Applicators - Sapphire light guides	Same as both	
Spot Size Applicator heads Applicator heads Applicator heads predicate Device Operation CW CW Same as both	Wavelength	1064 ± 10 nm (infrared)	1064 ± 10 nm (infrared)	1060 ± 20 nm (infrared)	Same as both	
Operation CW CW Same as both	Spot Size					
					Same as both	

Table 1: Substantial Equivalence Table for the Diode Laser Applicators

	Venus BlissMAX (Subject)	Venus BlissMAX (K213308) (Primary Predicate)	SculpSure Laser System (K182741) (Co-Secondary Predicate)	Substantial Equivalence Discussion
Power Density	Up to 1.4 W/cm ²	Up to 1.4 W/cm ²	Up to 1.4 W/cm ²	Same as both
Attachment to Patient	Belt	Belt	Belt	Same as both
Power Requirements	100-240V~50/60 Hz, Single Phase	100-240V~50/60 Hz, Single Phase	100-240V~50/60 Hz, Single Phase	Same as both
Materials	Biocompatible	Biocompatible	Biocompatible	Same as both
Sterility	Non-sterile	Non-sterile	Non-sterile	Same as both

Table 2: Substantial Equivalence Table for the MP2 (RF+ PEMF +Vacuum) Applicators

	Venus BlissMAX (Subject)	Venus BlissMAX (K213308) (Primary Predicate)	Venus Legacy Pro (K191528) (Co-Secondary Predicate)	Substantial Equivalence Discussion
Manufacturer	Venus Concept	Venus Concept	Venus Concept	Same as both
Classification Product Code, Regulation	Class II, PBX, 21 CFR 878.4400 21 CFR 878.5400	Class II, PBX, 21 CFR 878.4400 21 CFR 878.5400	Class II, GEI, PBX 21 CFR 878.4400	Same as both
Indications for Use	The Venus BlissMAX device is intended for the treatment of the following medical conditions; using the MP2 applicator for delivery of RF energy combined with massage and magnetic field pulses: • Relief of minor muscle aches and pain, relief of muscle spasm • Temporary improvement of local blood circulation • Temporary reduction in the appearance of cellulite.	The Venus BlissMAX device is intended for the treatment of the following medical conditions; using the MP2 applicator for delivery of RF energy combined with massage and magnetic field pulses: • Relief of minor muscle aches and pain, relief of muscle spasm • Temporary improvement of local blood circulation • Temporary reduction in the appearance of cellulite.	When used with the Octipolar (LB1) or Diamondpolar (LF1) applicators, the Venus Legacy Pro device is intended for use in dermatologic and general surgical procedures for females for the noninvasive treatment of moderate to severe facial wrinkles and rhytides in Fitzpatrick skin types I- IV. When used with the 4D Body (LB2) and 4D Face (LF2) applicators, the Venus Legacy Pro device is intended for the delivery of non-thermal RF combined with Massage and magnetic field pulses for the treatment of the following medical conditions: • Relief of minor muscles aches and	Same as Primary predicate Device and similar to the Co-Secondary Predicate.

	Venus BlissMAX (Subject)	Venus BlissMAX (K213308) (Primary Predicate)	Venus Legacy Pro (K191528) (Co-Secondary Predicate)	Substantial Equivalence Discussion
			 pain, relief of muscle spasm Temporary improvement of local blood circulation Temporary reduction in the appearance of cellulite 	
Energy Used / Delivered	 RF Energy Pulsed Magnetic Field (PMF) Vacuum 	 RF Energy Pulsed Magnetic Field (PMF) Vacuum 	 RF Energy Pulsed Magnetic Field (PMF) Vacuum 	Same as Primary predicate
Applicator Footprint Dimensions	MP ² : 38.5 cm ²	MP ² : 38.5 cm ²	LB1: 23.7 cm2 LF1: 2.9 cm2 LB2: 38.5 cm2 LF2: 4.9 cm2	Same as Primary predicate
Performance	Frequency: 1MHz Vacuum pressure: - 400mbar Maximal RF output power: up to 150W PMF Power: 15 Gauss (15Hz)	Frequency: 1MHz Vacuum pressure: - 400mbar Maximal RF output power: up to 100W PMF Power: 15 Gauss (15Hz)	Frequency: 1MHz Vacuum pressure: - 400mbar Maximal RF output power: up to 150W PMF Power: 15 Gauss (15Hz)	Same as Co- Secondary Predicate device
Materials	Biocompatible	Biocompatible	Biocompatible	Same as both
Power Requirements	100-120 VAC / 60Hz 220-240 VAC / 50Hz	100-240V~50/60 Hz, Single Phase	100-120 VAC / 60Hz 220-240 VAC / 50Hz	Same as both

	Venus BlissMAX (Subject)	Venus BlissMAX (K213308) (Primary Predicate).	Substantial Equivalence Discussion		
Basic Unit Characteristics					
Classification Produc	tClass II, NGX,	Class II, NGX,	Same		
Code, Regulation	21 CFR 890.5850,	21 CFR 890.5850			
Indications for Use	the FlexMAX applicators is not intended to be used in conjunction with therapy or treatment of medical diseases or medical conditions of any kind. The Venus BlissMAX device using	The Venus BlissMAX deviceusing the FlexMAX applicators is intended for muscle conditioning to stimulate healthy muscles. The Venus BlissMAX deviceusing the FlexMAX applicators is not intended to be used in conjunction with therapy or treatment of medical diseases or medical conditions of any kind. The Venus BlissMAX deviceusing the FlexMAX applicators is	Same		
	intended to be operated by a trained professional.	intended to be operated by a trained professional.			
Components	 The Venus BlissMAX device consists of the following components: Console, including a power supply unit, controller and user interface including an LCD touch 	The Venus BlissMAX device consists of the following components: • Console, including a power supply unit, controller and user interface including an LCD touch	Same		
Mechanism of	 screen. Several types of applicatorsfor different indications (Laser, MP², FlexMAX applicators) connected to theconsole via a cable. Muscle contraction by 	 screen. Several types of applicatorsfor different indications (Laser, MP², FlexMAX applicators) connected to theconsole via a cable. Muscle contraction by 	Same		
Action	electrical pulsing.	electrical pulsing.	Suite		
Power Source	Main Line Frequency (nominal): 50-60Hz Input Voltage (nominal): 100- 240VAC	Main Line Frequency (nominal): 50-60Hz Input Voltage (nominal): 100- 240VAC	Same		
Method of Line Current Isolation	Power Supply isolation	Power Supply isolation	Same		
Patient Leakage	Normal condition:1 µA	Normal condition:1 µA	Same		
Current	Single Fault condition:11 µA	Single Fault condition:11 µA			
Number of Output Modes	1	1	Same		

Number of Output Channels	4	4	Same
Channel Output:	Synchronous	Synchronous	Same
Synchronous or		Synomeneus	Sume
Alternating			
Method of channel	Isolation transformer per channel	Isolation transformer per channel	Same
isolation		1	
Regulated	Regulated voltage with current	Regulated voltage with current	Same
Current or	limit	limit	
Regulated Voltage Software/	Yes	Yes	Same
Software/	i es	res	Same
Microprocessor			
Control			
Automatic EMS Shut	Yes- Automatic EMStreatment	Yes- Automatic EMStreatment	Same
off	shut off	shut off	
		V D 1	
Patient Override Control	Yes- Emergency button	Yes- Emergency button	Same
Control Indicator Display:On /	Vas	Yes	Same
Off Status	1 05	1 es	Same
Indicator Display:	Yes	Yes	Same
Voltage Current	105	105	Same
Level			
Timer Range	Selectable 15, 30, 45, 60 min	Selectable 15, 30, 45, 60 min	Same
Weight	62 Kg / 136.7 lbs	62 Kg / 136.7 lbs	Same
Dimensions	21.7 x 25.6 x 53.2 in	21.7 x 25.6 x 53.2 in	Same
	55 x 65 x 135 cm	55 x 65 x 135 cm	
Output Specifications			
Waveform	Symmetrical Biphasicwaveform	Symmetrical Biphasicwaveform	Same
			-
Shape	Rectangular	Rectangular	Same
Maximum Output	40V@500Ω	40V@500Ω	Same
Voltage	105V@2kΩ	105V@2kΩ	
	160V@10kΩ	160V@10kΩ	
Maximum Output	80 mA @ 500 Ω	80 mA @ 500 Ω	Same
Current	52.5 mA (a) 2 k Ω	52.5 mA @ 2 k Ω 16 mA @ 10 k Ω	
	$16 \text{ mA} @ 10 \text{ k}\Omega$		
Frequency range	1 Hz to 1000 Hz	1 Hz to 1000 Hz	Same
Pulse width range	500 to 2500 [µs]	500 to 2500 [μs]	Same

Net Charge	0 [μC] @ 500Ω	0 [μC] @ 500Ω	Same
	Zero net charge is achieved by	Zero net charge is achieved by	
	using symmetrical biphasic	using symmetrical biphasic	
	waveforms	waveforms	
Maximum Phase	40 [μC] @ 500Ω @ 1000 Hz	40 [μC] @ 500Ω @ 1000 Hz	Same
Charge			
Maximum Current	2.5 mA/cm^2	2.5 mA/cm^2	Same
Density			
Maximum Power	55 [mW/cm ²] @500Ω	55 [mW/cm ²] @500Ω	Same
Density			

8. SUBSTANTIAL EQUIVALENCY AND COMPARISON OF TECHNOLOGICAL SIMILARITIES & DIFFERENCES

As described in the comparison tables above, the Venus BlissMAX subject device has the same intended use and indications for use, similar technological characteristics, and principles of operation as its primary predicate and co-secondary predicate devices. The technological differences between the Venus BlissMAX device and its primary predicate and co-secondary devices do not raise any new issues of safety or effectiveness. The Venus BlissMAX device and its primary predicate device Venus Concept's Venus BlissMAX (K213308) and its co-secondary predicate devices the Venus Concept Venus Legacy Pro (K191528) and Cynosure SculpSure (K182741) are based on the same core technologies of Diode Laser (as in the Venus BlissMAX and Cynosure SculpSure), RF along with PEMF and vacuum massaging (as in the Venus BlissMAX) for the same indications for use. The design and components in the Venus BlissMAX device, including the console and the applicators are similar to the design and components found in the primary predicate and the co-secondary predicate devices.

The technological differences do not alter the device's core technology or performance and have been addressed by the manufacturer through the applicable safety standards (General controls and mitigation measures) and through non-clinical performance testing (Special controls).

9. PERFORMANCE DATA

Verification and Validation tests have been performed to show that the performance specifications of Venus BlissMAX system is the same as the primary predicate and co-secondary predicate devices. As MP² applicator is identical to the primary predicate and co-secondary predicate devices, the Bench test that was conducted in Legacy demonstrates that the use of 150W is considered to be safe. Testing to IEC 60601-1 and IEC 60601-1-2 has been conducted on the Venus BlissMAX System. All performance testing demonstrated that the Venus BlissMAX System performs according to specifications and functions as intended.

The following non-clinical testing was performed on the Venus BlissMAX System. These tests verified that the Venus BlissMAX System was identical to the primary predicate and co-secondary predicate devices.

Electrical Safety and Electromagnetic Compatibility: Electrical safety and electromagnetic

compatibility (EMC) testing for the Venus BlissMAX System was conducted by an independent test laboratory in accordance with IEC 60601-1, Medical electrical equipment, Part 1: General requirements for basic safety and essential performance and with IEC 60601-1-2, Collateral Standard: Electromagnetic Compatibility - Requirements and Tests.

The Venus BlissMAX System was determined to be in conformance with applicable IEC standards (IEC 62366, IEC 60601-1 and IEC 60601-1-2).

The Venus BlissMax was also determined to be in conformance with the following electrical and safety testing:

IEC 60601-2-2 :2017 ed 6 Particular requirements for the basic safety and essential performance of high frequency surgical equipment and high frequency surgical accessories.

EN 60601-2-22:2007+A1:2012 Particular requirements for basic safety and essential performance of laser equipment.

IEC 60601-2-10:2012+A1:2016 Particular requirements for the basic safety and essential performance of nerve and muscle stimulator.

IEC 60825-1:2014 Safety of laser products - Part 1: Equipment classification and requirements.

IEC 60601-1-6: 2013 ed 3 General requirements for basic safety and essential performance - Collateral standard: Usability.

<u>Software</u>: Software verification testing was conducted, and results demonstrated that testing results were found acceptable for software release. Software testing was performed per FDA's guidance document "Guidance for the Content of Premarket Submission for Software Contained in Medical Devices".

Biocompatibility

The patient-contacting components of the device were demonstrated to be biocompatible including evaluation of cytotoxicity, irritation, and sensitization, acute systemic toxicity and materialmediated pyrogenicity per the FDA guidance document "Use of International Standard ISO 10993-1, 'Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process.

Clinical Testing

No Clinical testing was performed.

10. SUBSTANTIAL EQUIVALENCE SUMMARY

The Venus BlissMAX System is the same device and has the same technological characteristics as the primary device: Venus BlissMAX (K213308) and co-secondary predicate: Venus Legacy Pro Legacy Pro (K191528) and has the same intended use and indications as the co-secondary

predicate device: Cynosure SculpSure (K182741). Thus, the Venus BlissMAX System is Substantially equivalent to its predicate devices.

11. CONCLUSION

The performance testing data demonstrates that Venus BlissMAX System is as safe and effective as the legally marketed primary predicate and co-secondary predicate devices. The Venus BlissMAX System did not raise new questions of safety or effectiveness. Therefore, based on the information provided in this Premarket Notification, we conclude that Venus BlissMAX System has demonstrated substantial equivalence to the predicate devices and the performance testing data support the indications for use.