



June 29, 2022

Daxwell, LLC  
% Ray Wang  
General Manager  
Beijing Believe-Med Technology Service Co., Ltd.  
Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd.,  
FangShan District  
Beijing, Beijing 102401  
China

Re: K220834  
Trade/Device Name: Powder Free Nitrile Examination Glove  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LZA  
Dated: May 19, 2022  
Received: May 26, 2022

Dear Ray Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian, M.D., Ph.D.  
Acting Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K220834

Device Name  
Powder Free Nitrile Examination Glove

Indications for Use (Describe)

The Powder Free Nitrile Examination Glove is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary K220834

This 510(k) Summary of 510(k) safety and effectiveness information is being submitted in accordance with requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) Number: K220834

1. Date of Preparation: 06/29/2022

2. Sponsor

**DAXWELL, LLC**

2825 Wilcrest Dr #500

Houston, TX 77042 US

Contact Person: Frank Zhang

Position: Head of Firm

Tel: 01-281-699-0622

Fax: 01-281-699-0617

Email: scm.asst@Daxwell.com

3. Submission Correspondent

**Beijing Believe-Med Technology Service Co., Ltd.**

Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd., FangShan District, Beijing, China,102401

Contact Person: Ray Wang

Position: General Manager

Tel: +86-18910677558

Fax: +86-10-56335780

Email:information@believe-med.com

4. Proposed Device Identification

Trade Name: Powder Free Nitrile Examination Glove

Common Name: Powder Free Nitrile Examination Glove

**Regulatory Information:**

Classification: I

Product Code: LZA

Regulation Number: 21 CFR 880.6250

Review Panel: General Hospital

Indication for Use Statement:

The Powder Free Nitrile Examination Glove is a disposable device intended for medical purposes that is worn on the examiner’s hands to prevent contamination between patient and examiner.

5. Predicate Device Identification

510(k) Number: K212497

Product Name: Nitrile Examination Glove (Powder free, Blue)

Manufacturer: Jiangsu Jinlian Medical Technology Co., Ltd

6. Device Description

The proposed device, Powder Free Nitrile Examination Glove is a disposable device intended for medical purposes that is worn on the examiner’s hands to prevent contamination between patient and examiner.

The proposed devices are Powder Free Nitrile Examination Gloves and include variations of different sizes from X-Small to X-large. The color of the proposed device is Blue and Indigo

The proposed device is not provided as sterilized

The proposed device is made of Nitrile.

Table 1 Device Size Specifications

Size Model	Cuff Thickness (mm)	Palm Thickness (mm)	Finger Thickness (mm)	Width (mm)	Length (mm)	Color
XS	≥ 0.05	≥ 0.05	≥ 0.05	70±10	≥ 220	Blue
S	≥ 0.05	≥ 0.05	≥ 0.05	80±10	≥ 220	
M	≥ 0.05	≥ 0.05	≥ 0.05	95±10	≥ 230	
L	≥ 0.05	≥ 0.05	≥ 0.05	110±10		
XL	≥ 0.05	≥ 0.05	≥ 0.05	120±10		
XS	≥ 0.05	≥ 0.05	≥ 0.05	70±10	≥ 220	Indigo
S	≥ 0.05	≥ 0.05	≥ 0.05	80±10	≥ 220	
M	≥ 0.05	≥ 0.05	≥ 0.05	95±10	≥ 230	
L	≥ 0.05	≥ 0.05	≥ 0.05	110±10		
XL	≥ 0.05	≥ 0.05	≥ 0.05	120±10		

Table 2 Performance and Physical Specifications

Before Aging		After Aging		Pinhole AQL
Tensile Strength	Ultimate Elongation	Tensile Strength	Ultimate Elongation	
15 MPa, min	500 % min	14 MPa, min	500 % min	1.5

**7. Technological Characteristic Comparison Table**

Table 3 General Comparison

ITEM	Proposed Device(K220834) Powder Free Nitrile Examination Glove	Predicate Device (K212497) Nitrile Examination Glove (Powder free, Blue)	Remark
Product Code	LZA	LZA	SAME
Regulation No.	21 CFR 880.6250	21 CFR 880.6250	SAME
Class	I	I	SAME
Intended Use / Indications for Use	The Powder Free Nitrile Examination Glove is a disposable device intended for medical purposes that is worn on the examiner’s hands to prevent contamination between patient and examiner.	The Nitrile Examination Glove (Powder free, Blue) is a disposable device intended for medical purposes that is worn on the examiner’s hands to prevent contamination between patient and examiner.	SAME
Sterility	Non-sterile	Non-sterile	SAME
Single use or reuses	Single use	Single use	SAME
Powdered or Powered free	Powdered free	Powdered free	SAME

Table 4 Device Dimensions Comparison

Proposed Device (K220834) Powder Free Nitrile Examination Glove (Blue, Indigo)	Designation	Size					Tolerance
		XS	S	M	L	XL	
	Length, mm	230	230	230	230	230	±10
	Width, mm	70	80	95	110	120	±10
Thickness, mm:							
	Finger	0.12					±0.03
	Palm	0.10					±0.03
	Cuff	0.09					±0.03
Predicate Device (K212497) Nitrile Examination Glove (Powder free, Purple-Blue, Blue)	Designation	Size				Tolerance	
		S	M	L	XL		
	Length, mm	230	230	230	230	min	
	Width, mm	80	95	110	120	±10	
Thickness, mm:							
	Finger	0.05				min	
	Palm	0.05				min	
	Cuff	0.05				min	
Remark	Analysis 1						

Analysis 1:

The proposed device has different size specification to the predicate device, but all proposed devices meet the specifications of ASTM D 6319.

Table 5 Performance Comparison

ITEM		Proposed Device(K220834) Powder Free Nitrile Examination Glove	Predicate Device (K212497) Nitrile Examination Glove (Powder free, Purple-Blue, Blue)	Remark	
Colorant		Blue, Indigo	Blue	Analysis 2	
Physical Properties	Before Aging	Tensile Strength	15 MPa, min	14 MPa, min	Analysis 3
		Ultimate Elongation	500 % min	500 % min	SAME
	After Aging	Tensile Strength	14 MPa, min	14 MPa, min	SAME
		Ultimate	500 % min	400 % min	Analysis 4

	Elongation		
	Comply with ASTM D6319		Comply with ASTM D6319
Freedom from Holes	Be free from holes when tested in accordance with ASTM D5151 (AQL:1.5)		Be free from holes when tested in accordance with ASTM D5151 (AQL:1.5)
Powder Content	Less than 2 mg per glove when tested in accordance with ASTM D6124		Less than 2 mg per glove when tested in accordance with ASTM D6124.

Analysis 2:

The proposed device has different color to the predicate device, this difference may cause potential biocompatibility risk, for this risk we conducted the biocompatibility testing according to ISO 10993-11 and ISO 10993-10, the test results showed that the proposed devices with blue colorant and indigo colorant did not induce skin irritation, and showed no significant evidence of causing skin sensitization and acute system toxicity reactions.

Analysis 3:

The proposed device has different Tensile Strength before aging specification to the predicate device, but all proposed device meets the specification requirements of ASTM D 6319.

Analysis 4:

The proposed device has different Ultimate Elongation after aging specification to the predicate device, but all proposed device meets the specification requirements of ASTM D 6319.

Table 6 Safety Comparison

ITEM	Proposed Device(K220834) Powder Free Nitrile Examination Glove	Predicate Device (K212497) Nitrile Examination Glove (Powder free, Blue)	Remark
Material	Nitrile	Nitrile	SAME
Biocompatibility	Irritation	Under the conditions of the study, not an irritant	SAME
	Sensitization	Under conditions of the study, not a sensitizer.	
	Systemic toxicity	Under conditions of the study, not a systemic toxicity	Under conditions of the study, not a systemic toxicity
Label and Labeling	Meet FDA’s Requirements	Meet FDA’s Requirements	SAME



### 8.0 Summary of Non-Clinical Testing

Bench tests were conducted to demonstrate that the proposed device complies with the following standards:

ASTM D6319-19, Standard Specification for Nitrile Examination Gloves for Medical Application.

ASTM D5151-19, Standard Test Method for Detection of Holes in Medical Gloves.

ASTM D6124-17, Standard Test Method for Residual Powder on Medical Gloves.

ISO 10993-10:2010, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization.

ISO 10993-11:2017, Biological evaluation of medical devices - Part 11: Tests for systemic toxicity.

Table 7 Performance Test Results Summary

Test Method	Purpose	Acceptance Criteria		Color	Results
ASTM D5151	Testing for Freedom from holes	Freedom from holes (AQL:1.5)		Blue	No water leakage is inspected form 200 samples
				Indigo	
ASTM D6124	Determine the powder residue for powder free gloves	<2.0 mg per glove		Blue	Residual Powder: Average 0. 38 Mg;
				Indigo	
ASTM D6319	Testing for Physical property characteristics	Before Aging	Tensile Strength: 15 MPa, min Ultimate Elongation: 500 % min	Blue	Before Aging: Tensile Strength: ≥ 17 MPa; Ultimate Elongation: ≥ 512%.
		After Aging	Tensile Strength: 14 MPa, min Ultimate Elongation: 500 % min	Indigo	After Aging: Tensile Strength: ≥ 17 MPa; Ultimate Elongation: ≥ 508%.
	Testing For physical dimensions specification	Length: 230 mm ± 10 for all size (XS, S, M, L, XL); Width: 70±10 mm for XS; 80±10 mm for S; 95±10 mm for M; 110±10 mm for L; 120±10 mm for XL. Finger Thickness: 0.12±0.03 mm; Palm Thickness: 0.10±0.03 mm; All acceptance criteria above meet the requirements in Table 2 Dimensions and Tolerances of ASTM D6319		Blue/ Indigo	Length of Size XS: ≥ 229 mm; Width of Size XS: ≥70 (70-73) mm; Palm Thickness of Size XS: ≥0.09 mm; Finger Thickness of Size XS: ≥0.11 mm.  Length of Size S: ≥ 230 mm; Width of Size S: ≥85 (85-86) mm; Palm Thickness of Size S: ≥0.09 mm;

				<p>Finger Thickness of Size S:  <math>\geq 0.11</math> mm.</p> <p>Length of Size M: <math>\geq 230</math> mm;                  Width of Size M: <math>\geq 95</math> (95-96) mm;                  Palm Thickness of Size M:  <math>\geq 0.09</math> mm;                  Finger Thickness of Size M:  <math>\geq 0.11</math> mm.</p> <p>Length of Size L: <math>\geq 230</math> mm;                  Width of Size L: <math>\geq 108</math> (108-110) mm;                  Palm Thickness of Size L: <math>\geq 0.09</math> mm;                  Finger Thickness of Size L:  <math>\geq 0.11</math> mm.</p> <p>Length of Size XL: <math>\geq 230</math> mm;                  Width of Size XL: <math>\geq 120</math> (120-121) mm;                  Palm Thickness of Size XL:  <math>\geq 0.09</math> mm;                  Finger Thickness of Size XL:  <math>\geq 0.11</math> mm.</p>
ISO 10993-11	Evaluate the endpoint of systemic toxicity for biocompatibility	The test article showed “negative” systemic toxicity	Blue	Under the conditions of the study, the test article showed “negative” systemic toxicity.
			Indigo	
ISO 10993-10	Evaluate the endpoint of irritant for biocompatibility	The response of the test article has no skin irritation	Blue	Under the experimental conditions, the test article has no skin irritation on rabbits.
			Indigo	
	Evaluate the endpoint of sensitization for biocompatibility	The test article showed no evidence of causing delayed dermal contact sensitization.	Blue	The test article showed no evidence of causing delayed dermal contact sensitization in the guinea pig.
		Indigo		

### 9.0 Summary of Clinical Testing

Clinical testing is not needed for this device.

**10.0 Conclusion**

Based on the nonclinical tests performed, the subject device is as safe, as effective, and performs as well as the legally marketed predicate device, Nitrile Examination Glove (Powder free, Blue) cleared under K212497.