

June 23, 2023

Medline Industires, Inc. % Joy Gutermuth Senior Specialist (Consultant) Rqm+ 2790 Mosside Blvd. Suite 800 Monroeville, Pennsylvania 15146

Re: K220955

Trade/Device Name: Hudson RCI Variable concentration Large Volume Nebulizer (1770) Regulation Number: 21 CFR 868.5630 Regulation Name: Nebulizer Regulatory Class: Class II Product Code: CAF Dated: May 26, 2023 Received: May 26, 2023

Dear Joy Gutermuth:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Ting Song -S

#### For

Assistant Director DHT1C: Division of Sleep Disordered Breathing, Respiratory and Anesthesia Devices OHT1: Office of Ophthalmic, Anesthesia,

Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number *(if known)* K220955

Device Name Hudson RCI Large Volume Nebulizer

#### Indications for Use (Describe)

The non-prefilled nebulizer is indicated for use when humidity needs to be added to the air flow going to a patient in the form of an aerosol. The large volume nebulizer also allows for the provision of humidified oxygen at select oxygen concentrations from 28% to 98%. This product is for single patient use and is designed to be used in hospital, nursing homes, extended care facilities, outpatient clinics as prescribed by a healthcare professional.

Type of Use (Select one or both, as applicable)	
Rescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

## 510(k) SUMMARY

#### DATE PREPARED

June 23, 2023

#### MANUFACTURER AND 510(k) OWNER

Medline Industries, Inc. 1 Three Lake Dr. Northfield, IL 60093 Phone: 724-640-9680

Contact Person: Nicole Schaffer, Manager Regulatory Affairs, Respiratory

#### **REPRESENTATIVE/CONSULTANT**

Joy Gutermuth, Regulatory Consultant RQM+ Telephone: +1 (412) 810-8172 Email: jgutermuth@rqmplus.com

#### **DEVICE INFORMATION**

Proprietary Name/Trade Name:	Hudson RCI Variable Concentration Large Volume
	Nebulizer
Common Name:	Large Volume Nebulizer
Regulation Number:	§868.5630
Class:	II
Product Code:	CAF

#### PREDICATE DEVICE IDENTIFICATION

The Hudson RCI Variable Concentration Large Volume Nebulizer is substantially equivalent to the following predicates:

510(k) Number	Predicate Device Name / Manufacturer	Primary
		Predicate
K041418	AirLife Misty Finity Large Volume Continuous Neublizer	✓
K141214	Hudson RCI® AquaPak® Prefilled Nebulizer	Reference Device

The predicate devices have not been subject to a design related recall.

## **DEVICE DESCRIPTION**

The large volume nebulizer is a non-prefilled reservoir nebulizer for supplying humidity for inhalation therapy. The device features a wing nut style connector that fits standard flow metered medical gas sources and includes a 500mL capacity jar with minimum and maximum fill lines. Large volume nebulizers utilize an internal venturi nozzle to draw the solution up from the jar through a small plastic pickup tube and into the gas stream to be aerosolized. A rotating collar sets the delivered oxygen concentration by controlling the size of the room air opening around the venturi. A bull-nose style output connector is used to connect 22mm aerosol tubing for delivery to the patient.

#### **INDICATIONS FOR USE**

The non-prefilled nebulizer is indicated for use when humidity needs to be added to the air flow going to a patient in the form of an aerosol. The large volume nebulizer also allows for the provision of humidified oxygen at select oxygen concentrations from 28% to 98%. This product is for single patient use and is designed to be used in hospital, nursing homes, extended care facilities, outpatient clinics as prescribed by a healthcare professional.

#### COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

Medline Industries believes that the large volume nebulizer is substantially equivalent to the predicate devices based on the information summarized here:

The subject device has a similar design and dimensions, and uses similar or identical materials as the devices cleared in K041418 and K141214. The subject device has the same intended use and similar technological characteristics to the devices cleared in both K041418 and K141214. The device has similar instrumentation to the device cleared in K041418 and K141214. These technological characteristics have undergone testing to ensure the device is as safe and effective as the predicates.

Product Features	<u>Proposed</u> Hudson RCI Variable Concentration Large Volume Nebulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	<u>Reference Device</u> Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Classification	Class II	Class II	Class II
Product Code	CAF	CAF	CAF
Regulation Number	§868.5630	§868.5630	§868.5630
<b>Regulation Name</b>	Nebulizer	Nebulizer	Nebulizer

Product Features	<u>Proposed</u> Hudson RCI Variable Concentration Large Volume Nebulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	Reference Device Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Intended Use	To add humidity in aerosol form to a patient's breathing gases.	This device is intended to be used to spray liquids in aerosol form into gases that are delivered directly to the patient for	To add humidity in aerosol form to a patient's breathing gases.
Indications for Use	The non-prefilled nebulizer is indicated for use when humidity needs to be added to the air flow going to a patient in the form of an aerosol. The large volume nebulizer also allows for the provision of humidified oxygen at select oxygen concentrations from 28% to 98%. This product is for single patient use and is designed to be used in hospital, nursing homes, extended care facilities, outpatient clinics as prescribed by a healthcare professional.	This device is intended to be used to spray liquids in aerosol form into gases that are delivered directly to the patient for breathing. Its use is indicated whenever a physician or healthcare professional administers or prescribes medical aerosol products to a patient using a Continuous Large Volume Nebulizer. This product is a single patient use, non-sterile prescriptive device and is designed to be used under medical supervision in hospitals, nursing homes, extended care facilities or outpatient clinics.	The Hudson RCI® AquaPak® Prefilled Nebulizer adds sterile water or saline solution in aerosol form to a patient's breathing gases.

Product Features	<u>Proposed</u> Hudson RCI Variable Concentration Large Volume Nebulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	Reference Device Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Representative		1	
Image	-		
	11		Note that the second seco
		ia ne Marine	

Product Features	Proposed Hudson RCI Variable Concentration Large Volume Nebulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	<b><u>Reference Device</u></b> Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Device	The large volume	The nebulizer is a	The Hudson RCI®
Description	nebulizer is a non-	single patient use	AquaPak <sup>®</sup> Prefilled
	prefilled reservoir	device, which is filled	Nebulizers provide
	nebulizer for	with a fluid, typically	sterile water or sterile
	supplying numidity	respiratory	saline for innalation
	therapy The device	connected to an air	generate aerosol a fine
	features a wing nut	source via flexible	mist of liquid water (or
	style connector that	tubing. The nebulizer	sodium chloride
	fits standard flow	works by having the	solution) that is
	metered medical gas	fluid come into	suspended in the gas to
	sources and includes	contact with the	be inhaled by the
	a 500mL capacity jar	steam of gas. The gas	patient. The Hudson
	with minimum and	shatters the liquid	RCI® AquaPak®
	maximum fill lines.	into small particles.	Prefilled Nebulizers are
	Large volume	impost a haffle that	and provide a fine mist
	internal venturi	further reduces the	of sterile water or
	nozzle to draw the	size of the particles	saline solution to
	solution up from the	The majority of the	inspired gas during
	jar through a small	larger particles settle	aerosol therapy.
	plastic pickup tube	inside the nebulizer	Prefilled sterile
	and into the gas	as a result of gravity	reservoirs for
	stream to be	and inertia, returning	AquaPak <sup>®</sup> Nebulizers
	aerosolized. A	the mist to liquid to	come in three sizes:
	rotating collar sets	repeat the	440ml, 760ml, and
	the delivered oxygen	nebulization process.	10/0ml. Each reservoir
	controlling the size	are then administered	must be used with a suitable adaptor
	of the room air	as the patient inhales	component which
	opening around the	The treatment is	connects the system to
	venturi. A bull-nose	completed when the	a flow-metered gas
	style output	majority of fluid is	source and provides
	connector is used to	nebulized.	nebulizer functionality.
	connect 22mm		
	aerosol tubing for		
	delivery to the		
	patient.		

Product Features	<u>Proposed</u> Hudson RCI Variable Concentration Large Volume Nabulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	Reference Device Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Accessories	Masks (not	Masks (not included)	Nebulizer Adaptors
Accessories	included) tubing	tubing (not included)	masks (no included)
	(not included)	(not moradou)	tubing (not included),
Environment	Hospital or home care setting.	Hospitals, nursing homes, extended care facilities or outpatient clinics.	Hospital or home care setting.
Patient Population	Pediatric (ages 2 years and above) and adults.	Infant, pediatric, and adult	Pediatric (ages 2 years and above) and
Principle of	Iet nebulizer with	Iet nebulizer with	Iet nebulizer with
Operation	adjustable air	adjustable air	adjustable air
• <b>F</b> · · · · · · · ·	entrainment	entrainment	entrainment
Sterilization Method	Non-sterile	Non-sterile	Sterile
Patient Contacting Materials	Polypropylene, polyethylene, LDPE, polystyrene, acrylic,	Thermoplastics	Thermoplastics
	dye		
Single Use	Single patient use, disposable	Single patient use, disposable	Single patient use, disposable
Adjustable Oxygen Settings	FIO <sub>2</sub> 28-98%	FIO <sub>2</sub> 28-98%	FIO <sub>2</sub> 28-98%
Biocompatibility	<ul> <li>External communicating, prolonged/perman ent contact device that indirectly contacts tissue/bone/dentin</li> <li>Indirect gas pathway</li> </ul>	• Unknown from publicly available summary what biocompatibility endpoints were assessed.	<ul> <li>External communicating, prolonged/permanent contact device that indirectly contacts tissue/bone/dentin.</li> <li>Indirect gas pathway</li> </ul>

Product Features	Proposed Hudson RCI Variable Concentration Large Volume Nebulizer	<u>Predicate</u> AirLife Misty Finity Large Volume Continuous Nebulizer (K041418)	<u>Reference Device</u> Hudson RCI® AquaPak® Prefilled Nebulizer (K141214)
Standards Utilized	ISO 10993-1 ISO 10993-5 ISO 10993-10 ISO 10993-11 ISO 10993-18 ISO 18562-1 ISO 18562-2 ISO 18562-3	ISO 10993-1	ISO 10993-1 ISO 10993-5 ISO 10993-10 ISO 10993-11 ISO 10993-18
Non-Clinical Testing	<ul> <li>Packaging</li> <li>Environmental Conditioning (high and low humidity)</li> <li>Aging</li> <li>Oxygen entrainment</li> <li>Lift testing</li> <li>Humidity output</li> <li>Useful life testing</li> <li>Cleaning process</li> </ul>	"Performance evaluation of the proposed and predicated devices consisted of cascade impaction and output rate testing."	<ul> <li>Packaging</li> <li>Environmental Conditioning (high and low humidity)</li> <li>Aging</li> </ul>

## SUMMARY OF NON-CLINICAL TESTING

No FDA performance standards have been established for the Hudson RCI Variable Concentration Large Volume Nebulizer. The following tests were performed to demonstrate safety based on current industry standards:

- Biocompatibility
- Packaging
- Environmental Conditioning (high and low humidity)
- Aging
- Oxygen entrainment
- Lift testing
- Humidity output
- Useful life testing
- Cleaning process

The results of these tests indicate that the Hudson RCI Variable Concentration Large Volume Nebulizer is substantially equivalent to the predicate devices.

#### CONCLUSION

Based on the testing performed, including additional gas pathway biocompatibility according to ISO 18562, it can be concluded that the subject device does not raise new issues of safety or effectiveness compared to the predicate devices. The similar indications for use, technological characteristics, and performance characteristics for the proposed Hudson RCI Variable Concentration Large Volume Nebulizer are assessed to be substantially equivalent to the predicate devices.