

June 16, 2023

Fresenius Kabi % Jason Ma Sr. Manager, Regulatory Affairs Fresenius Kabi USA,LLC 3 Corporate Dr Suite 300 Lake Zurich, Illinois 60047

Re: K221121

Trade/Device Name: Volumat Polyethylene I.V. Administration Set

Regulation Number: 21 CFR 880.5440

Regulation Name: Intravascular Administration Set

Regulatory Class: Class II

Product Code: FPA Dated: May 18, 2023 Received: May 18, 2023

Dear Jason Ma:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

K221121 - Jason Ma Page 2

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For David Wolloscheck, Ph.D.

Assistant Director

DHT3C: Division of Drug Delivery and

General Hospital Devices,

and Human Factors

Porsche Bennet

OHT3: Office of GastroRenal, ObGyn,

General Hospital and Urology Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

K221121				
Device Name Volumat Polyethylene I.V. Administration Set				
Indications for Use (Describe) I.V. Administration Set for the infusion of parenteral fluids and medications from a container into the patient's vascular system through a vascular access device with Agilia VP MC/Volumat MC Agilia pump or gravity only.				
Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)				
CONTINUE ON A SEPARATE PAGE IF NEEDED.				

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

510(k) SUMMARY- K221121

1. Date Prepared: June 16, 2023

2. Submitter Information

Name: Fresenius Kabi

Address: Fresenius Kabi AG

61346 Bad Homburg

Germany

Contact Person: Jason Ma

Sr. Manager, Regulatory Affairs

Fresenius Kabi USA, LLC

Three Corporate Drive, 2nd Floor Lake Zurich, IL 60047 USA

Phone: 224-817-4100 Fax: 847 550 2960

E-mail: Jason.ma01@fresenius-kabi.com

Secondary Contact Person: Keith Dunn

Director, Regulatory Affairs Fresenius Kabi USA, LLC

Three Corporate Drive, 2nd Floor Lake Zurich, IL 60047 USA

Phone: 224-817-2430 Fax: 847 550 2960

E-mail: Keith.Dunn@fresenius-kabi.com

3. Device Name and Classification

Device Trade Name: Volumat Polyethylene I.V. Administration Set

Common Name: I.V. Administration Set

Classification Name: Intravascular administration set

Regulation Number 21 CFR 880.5440

Regulatory Class:Class II **Product Code:**FPA



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

4. Predicate Device

Device Trade Name: Intravascular Administration Sets Common Name: Intravascular Administration Set

Classification Name: 21 CFR 880.5440

Regulatory Class: Class II **Product Code:** FPA **510(k) Number:** K203609

5. Device Description

The Volumat Polyethylene I.V. Administration Sets (PE Sets) are available for dedicated use with the Agilia Infusion Pump and the Agilia VP MC Pump, or they can administer parenteral fluids and medications by gravity flow. Both the Agilia Infusion Pump and the Volumat I.V. Administration Sets were cleared under K121613. The Agilia VP MC Pump was cleared under K210073. The PE Sets include a spike, drip chamber, roller clamp, upstream/downstream clamp, backcheck valve, pump segment, needle-free port, rotating male luer lock, and tubing.

Set Number	Description
M46441395	Volumat Polyethylene I.V. Administration Set with spike, air vent, drip chamber, roller clamp, pump segment, rotating male luer
M46441995	Volumat Polyethylene I.V. Administration Set with spike, air vent, drip chamber, roller clamp, pump segment, downstream clamp, downstream needle-free port, rotating male luer
M46444495	Volumat Polyethylene I.V. Administration Set with spike, air vent, drip chamber, upstream clamp, backcheck valve, upstream needle-free port, roller clamp, pump segment, downstream clamp, downstream needle-free port, rotating male luer
M46445495	Volumat Polyethylene I.V. Administration Set with spike, air vent, drip chamber, upstream clamp, backcheck valve, upstream needle-free port, roller clamp, pump segment, downstream clamp, downstream needle-free port, rotating male luer

6. Principle of Operation

The range of sets provides options for intermittent or continuous delivery of parenteral fluids (solutions, colloids, parenteral nutrition) and medications (including but not limited to diluted drugs, chemotherapy) through clinically accepted intravenous (I.V.) routes of administration.

The device may be used for Adults and Pediatrics Patients.

7. Indication for Use/ Intended Use



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

Indication for Use:

I.V. Administration Set for the infusion of parenteral fluids and medications from a container into the patient's vascular system through a vascular access device with Agilia VP MC/Volumat MC Agilia pump or gravity only.

Intended Use:

Volumat Polyethylene I.V. Administration Sets are intended for the administration of drugs and solutions.

8. Comparison of the Technological Characteristics with the Predicate Device

The technological characteristics of the subject device are substantially equivalent to those of the predicate device in regard to the following technological characteristics:

- Principle of operation and conditions of use of the subject device are the same as those of the predicate device.
- Materials of the proposed device do not raise new questions of safety and effectiveness, as demonstrated by performance testing and biocompatibility evaluation.
- Physical specifications of the subject device are equivalent to those of the predicate device.

A comparison between the predicate device and the subject device is provided in Table 1 and Table 2 below.

Table 1: Summary of Substantial Equivalence Comparison

Technological	Volumat Polyethylene	Baxter Administration	Assessment of
Characteristics	I.V. Administration	Set K203609	Differences
	Sets (Subject Device)	(Predicate Device)	
Indication for Use	IV Administration set	For the administration of	Similar. Both I.V.
	for the infusion of	fluids from a container	administration sets are
	parenteral fluids and	into the patients'	used to administer fluids
	medications from a	vascular system through	from a container to a
	container into the	a vascular access device.	patient. However, the
	patient's vascular		Subject Device
	system through a		interoperates with
	vascular access device		infusion systems that
	with Agilia VP MC		have been cleared in
	pumps/Volumat MC		other 510(k) submissions
	Agilia pump or gravity		including K121613 and
	only.		K210073
Operating	The Volumat	The Baxter	Similar, both devices can
Mechanism	administration sets	administration sets can	be used with an infusion
	deliver the intravenous	be used to deliver	pump or by gravity flow.



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

Technological Characteristics	Volumat Polyethylene I.V. Administration Sets (Subject Device)	Baxter Administration Set K203609 (Predicate Device)	Assessment of Differences
	infusions via an infusion pump through a pumping mechanism ("cassette") or by gravity flow.	solutions for gravity or pump by infusion of I.V. fluids.	The predicate device does not have a pumping mechanism ("cassette"). Difference has been verified through performance testing including flow rate accuracy under various environmental conditions (temperature, pressure, humidity) which
Sterile	Yes	Yes	demonstrate equivalence. Same
Non-Pyrogenic	Yes	Yes	Same
Single Use	Yes	Yes	Same
Length	75-105 inches	69-133.5 inches	Similar. The sets included in the subject device range in length from 75 inches to 105 inches which is within the range of the predicate device. Difference tested according to ISO 8536-4 and flow rate accuracy to demonstrate the subject device performance.
Priming Volume	21-27 ml	6.1 to 21.2 mL	Similar. The Volumat set requires more priming volume than the predicate. Bench testing confirmed that the differences in priming volume do not impact safety or effectiveness.
Internal Tube Diameter	0.122 inches	0.102 to 133 inches	Similar internal and external diameter and
External Tube Diameter	0.164 inches	0.140-0.209 inches	approximate range of the predicate device. Difference tested to ISO 8536-8 has demonstrated equivalence.
Components			
Spike	Yes	Yes	Different. The subject
Drip chamber	Yes	Yes	device includes an



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

Technological Characteristics	Volumat Polyethylene I.V. Administration Sets (Subject Device)	Baxter Administration Set K203609 (Predicate Device)	Assessment of Differences
Backcheck Valve	Yes	Yes	additional pumping
Pump Segment	Yes	No	segment. However, no
Needle Free Port	Yes	Yes	significant difference
Male Luer	Yes	Yes	between the subject and
Roller Clamp	Yes	Yes	predicate administration
Connector Luer	Yes	Yes	sets that would raise new
Lock			issues of safety or
			efficacy.
			The differences have
			been verified in various
			aspects to demonstrate
			the subject device's
			safety and performance
			including:
			Biocompatibility testing
			according to ISO 10993
			collateral standards,
			Microbial Ingress
			Testing, and Particulate
			Testing; Performance Testing according to ISO
			8536-4, ISO 80369-20,
			ISO 8536-8 and ISO
			8436-14.
			0.20111

Table 2: Summary of Material Comparison

Component	Modified Volumat I.V. Administration Set (Subject Device)	Baxter Administration Set K203609 (Predicate Device)	Assessment of Differences
Spike	Polystyrene	Acrylonitrile	Similar, testing to
		Butadiene Styrene	biocompatibility and performance standards demonstrated equivalence
Drip chamber	Polystyrene, styrene- butadiene-copolymer, Polyamide	Polyvinyl Chloride	Different, testing to biocompatibility and performance standards demonstrated equivalence
Backcheck Valve	Methylmethacrylate- Acrylonitrile-Butadiene- Styrene)	Polymethyl methacrylate (acrylic) Silicone Rubber	Similar, testing to biocompatibility and performance standards demonstrated equivalence
Pump Segment	Methylmethacrylate- Acrylonitrile-Butadiene- Styrene	N/A – predicate device does not have a pump segment	Testing to biocompatibility and performance standards demonstrated equivalence



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

Component	Modified Volumat I.V. Administration Set (Subject Device)	Baxter Administration Set K203609 (Predicate Device)	Assessment of Differences
Needle Free Port	Methylmethacrylate-	Polyester	Similar, Testing to
	Acrylonitrile-Butadiene-		biocompatibility and
	Styrene, Copolyester		performance standards
			demonstrated equivalence
Rotating Male Luer	Polyvinylchloride,	Acrylonitrile	Similar, Testing to
Lock	Acrylonitrile butadiene	butadiene styrene	biocompatibility and
	styrene		performance standards
			demonstrated equivalence
Tubing	Polyvinyl chloride,	Polyvinyl chloride	Similar, Testing to
	ethylene-vinyl acetate,		biocompatibility and
	polyethylene		performance standards
			demonstrated equivalence

9. Substantial Equivalence

Intended Use/Indication for Use—Discussion of Differences

Both the subject and predicate device have the same intended use and similar indication for use. Both I.V. administration sets are used to administer fluids.

Technological Characteristics—Discussion of Differences

- Both the subject and predicate device can be used with an infusion pump or by gravity flow. The subject device delivers the intravenous infusions via an infusion pump through a pumping mechanism ("cassette") or by gravity flow, and the predicate device does not have a pumping mechanism ("cassette").
- The subject device ranges in length from 75 inches to 105 inches, which is within the length range of the predicate device.
- The subject device carries more priming volume than the predicate. Bench testing confirmed that the differences in priming volume do not impact safety or effectiveness.
- The subject device includes an additional pumping segment. However, no significant difference between the subject and predicate administration sets that would raise new issues of safety or efficacy.

Conclusion on Substantial Equivalence



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

The Volumat Polyethylene I.V. Administration Set has the same intended use and equivalent indication for use as the predicate device. The subject device has similar technological characteristics to the predicate, and the descriptive and performance information provided within this premarket notification demonstrates that:

- any differences do not raise different questions of safety and effectiveness compared to the predicate device; and
- the proposed device is at least as safe and effective as the legally marketed predicate device.

Based on the comparison of the intended use and the technological characteristics, the subject device is substantially equivalent to the currently marketed predicate device.

10.Performance Testing

Performance Testing—Bench

Functional performance bench testing was conducted to demonstrate that the Volumat Polyethylene I.V. Administration Set performs as intended.

The following performance testing was conducted to support the substantial equivalence determination:

- ISO 8536-8
 - o Leakage
 - o Tensile Strength
 - o Storage Tube Volume
- ISO 80369-20
 - Luer Fittings
 - Stress Cracking
 - o Resistance to Separation from Axial Loading
 - o Resistance to Separation from Unscrewing
 - Resistance to Overriding
- ISO 8536-4
 - Closure-piercing Device Testing
 - o Air-inlet Device Testing
 - Tubing Testing
 - o Flow Requirements of the Infusion Pump
 - o Drip Chamber Testing
 - o Clamp Opening and Closing
 - o Protective Cap Testing
 - Chemical Compatibility Testing



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

- ISO 8536-14
 - Clamps and Flow Regulators Testing
- ISO 11607-1 & 2 and ASTM D4169
 - o Sterile Barrier and Packaging Systems and Simulated Shipping
- Particulate USP<788> Particulate Matter in Injections
 - o Particulate Testing
- Pump Segment Compatibility Testing with Infusion Pump
- Operation under Temperature, Pressure and Humidity
- Microbial Ingress Testing
- Usability Testing

11. Biocompatibility Testing

Following the FDA Guidance: "Use of International Standard ISO 10993-1, Biological evaluation of medical devices-Part 1: Evaluation and testing within a risk management process", the tests selected were for prolonged, blood path indirect, externally communicating devices. The following biocompatibility tests were successfully conducted on the Volumat Polyethylene I.V. Administration Set:

- Acute Systemic Toxicity
- Hemolysis
- Irritation
- Sensitization
- Cytotoxicity
- Pyrogens
- Subacute/subchronic toxicity

12. Sterilization Validation

Sterilization was achieved by ethylene oxide and meets the requirements of DIN EN ISO 11135:2020, which is equivalent to ISO 11135 (2014-7) + AMD 1 (2018-10). The ethylene oxide sterilization method achieves a Sterilization Assurance Level (SAL) of 10⁻⁶.

13. Clinical Testing

No clinical testing was performed as this device does not require clinical studies to demonstrate substantial equivalence with the predicate device.



Three Corporate Drive Lake Zurich, Illinois 60047 T 847-550-2300 T 888-391-6300 www.fresenius-kabi.us

14. Conclusion

The differences between the predicate and the subject device do not raise any new or different questions of safety or effectiveness. The subject device, Volumat Polyethylene I.V. Administration Set is substantially equivalent to the predicate device, I.V. Administration Set, cleared under K203609.