

November 20, 2022

Empatica S.r.l. Alberto Poli Regulatory Affairs & Quality Manager Via Stendhal, 36 Milan, 20144 Italy

Re: K221282

Trade/Device Name: Empatica Health Monitoring Platform Regulation Number: 21 CFR 870.2700 Regulation Name: Oximeter Regulatory Class: Class II Product Code: DQA, DRG, FLL, LEL, GZO Dated: November 14, 2022 Received: November 14, 2022

Dear Alberto Poli:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <u>https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</u>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# James J. Lee -S

James J. Lee, Ph.D. Director DHT1C: Division of Sleep Disordered Breathing, Respiratory and Anesthesia Devices OHT1: Office of Ophthalmic, Anesthesia, Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

### Indications for Use

510(k) Number *(if known)* K221282

Device Name Empatica Health Monitoring Platform

Indications for Use (Describe)

The Empatica Health Monitoring Platform is a wearable device and paired mobile and cloud-based software platform intended to be used by trained healthcare professionals or researchers to remotely monitor physiologic parameters in ambulatory individuals 18 years of age and older in home-healthcare environments.

The device supports continuous data collection for monitoring the following physiological parameters:

- Peripheral skin temperature,
- Electrodermal activity,
- Blood Oxygen Saturation under no motion conditions,
- Activity associated with movement during sleep.

The Empatica Health Monitoring Platform can be used to analyze circadian rhythms and assess activity in any instance where quantifiable analysis of physical motion is desirable.

The Empatica Health Monitoring Platform is not intended for SpO2 monitoring in conditions of motion or low perfusion.

The Empatica Health Monitoring Platform is intended for peripheral skin temperature monitoring, where monitoring temperature at the wrist is clinically indicated.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

### **Empatica Srl Traditional 510(k)** Empatica Health Monitoring Platform



Empatica Health Monitoring Platform – 510(k)

## 510(k) Summary

Version 3.0

The information contained in this document is confidential and is the property solely of Empatica, S.r.l.. Any reproduction in part or in whole without the written permission of Empatica, S.r.l. is prohibited. © Empatica S.R.L, 2022. All rights reserved

**Empatica Health Monitoring Platform** 

### 510(k) Summary

#### I. SUBMITTER

Company Name	Empatica Srl
Establishment Registration Number	3012933969
Contact Person	Alberto Poli, Regulatory Affairs & Quality Manager
Contact Person email	apo@empatica.com
Address	Via Stendhal, 36 - 20144, Milan, Italy
Telephone Number	+39 02 36165068
Date prepared	April 29, 2022

#### II. DEVICE

Trade/Proprietary Name:	Empatica Health Monitoring Platform
Common/Usual Name:	Remote Patient Monitoring System

#### **Primary Product Code:**

Classification	Classification Name	Device	Product	Classification
Regulation		Class	Code	Panel
870.2700	Oximeter	Class II	DQA	Cardiovascular

#### **Secondary Product Codes:**

Classification	Classification Name	Device	Product	Classification
Regulation		Class	Code	Panel
870.2910	Transmitters and Receivers, Physiological	Class II	DRG	Cardiovascular
	Signal, Radiofrequency			
882.5050	Device, Sleep Assessment	Class II	LEL	Neurology
882.1540	Galvanic skin response measurement	Class II	GZO	Neurology
	device			
880.2910	Thermometer, Electronic, Clinical	Class II	FLL	General Hospital

#### III. PREDICATE DEVICES

Predicate Device	Name	Submitter	Product Code(s)	510(k) Number
Primary	Loop System	Spry Health, Inc.	DQA BZQ	K181352
Secondary	Current Wearable Health Monitoring System	Current Health Ltd.	MSX DQA DRG BZQ FLL	K191272
Secondary	ActiGraph CentrePoint Insight Watch	ActiGraph	LEL	K181077
Secondary	Empatica E4	Empatica S.r.l.	GZO	N/A

None of these predicates have been subject to a design-related recall.

No reference devices were used in this submission.

**Empatica Health Monitoring Platform** 

### IV. DEVICE DESCRIPTION

The Empatica Health Monitoring Platform is a wearable device and software platform composed by:

- A wearable medical device called EmbracePlus,
- A mobile application running on smartphones called "Care App",
- A cloud-based software platform named "Care Portal".

The EmbracePlus is worn on the user's wrist and continuously collects raw data via specific sensors. These data are wirelessly transmitted via Bluetooth Low Energy to a paired mobile device where the Care App is up and running. The data received are analyzed by one of the Care App software modules, EmpaDSP, which computes the user physiological parameters. Based on the version of the Care App installed, the user can visualize a subset of these physiological parameters. The Care App is also responsible for transmitting, over cellular or Wi-Fi connection sensors' raw data, device information, Care App-specific information, and computed physiological parameters to the Empatica Cloud. On the Empatica Cloud, these data are stored, further analyzed, and accessible by healthcare providers or researchers via a specific cloud-based software called Care Portal.

The platform is intended to continuously monitor adult patient physiological parameters in homehealthcare environment. It is designed for monitoring patients by trained healthcare professionals or researchers. It is intended to continuously monitor blood oxygen saturation (SpO<sub>2</sub>), peripheral skin temperature (TEMP), and electrodermal activity (EDA). Activity sensors are used to detect sleep periods and to monitor the activity associated with movement during sleep.

#### V. INDICATION FOR USE

The Empatica Health Monitoring Platform is a wearable device and paired mobile and cloud-based software platform intended to be used by trained healthcare professionals or researchers to remotely monitor physiologic parameters in ambulatory individuals 18 years of age and older in home-healthcare environments.

The device supports continuous data collection for monitoring of the following physiological parameters:

- Peripheral skin temperature,
- Electrodermal activity
- Blood Oxygen Saturation under no motion conditions,
- Activity associated with movement during sleep.

The Empatica Health Monitoring Platform can be used to analyze circadian rhythms and assess activity in any instance where quantifiable analysis of physical motion is desirable.

The Empatica Health Monitoring Platform is not intended for SpO2 monitoring in conditions of motion or low perfusion.

The Empatica Health Monitoring Platform is intended for peripheral skin temperature monitoring, where monitoring temperature at the wrist is clinically indicated.

**Empatica Health Monitoring Platform** 

### I. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICES

The application Empatica Health Monitoring Platform is substantially equivalent to the identified predicate devices. The devices have similar Indications for Use, features, technology, and accuracy.

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
Common Name	Oximeter	Oximeter	System, Network and Communication, Physiological Monitors	Sleep assessment device (actigraphy)	Device, Galvanic Skin Response Measurement	N/A
Device Manufacturer	Empatica S.r.l.	Spry Health Ltd.	Current Health Ltd.	ActiGraph, Inc.	Empatica S.r.l.	N/A
Device Classification	I	Ш	Ш	II	ll (510(k) exempt)	N/A
510(k) number	N/A	K181352	K191272	K181077	N/A	N/A
Primary Product Code	DQA	DQA	MSX	LEL	GZO	N/A
Secondary Product Code	DRG, GZO, LEL, FLL	BZQ	FLL, DQA, BZQ, DRG, BZG	-	-	N/A
Intended Use/Indications for Use	The Empatica Health Monitoring Platform is a wearable device and paired mobile and cloud-based software platform intended to be used by trained healthcare professionals or	The Loop System is intended for adult patients in the home environment for passive, noninvasive, intermittent data collection of physiological parameters that will later be transmitted to a web server for	The Current Wearable Health Monitoring System is intended for reusable bedside, mobile and central multi-parameter, physiologic patient monitoring of adult patients in professional healthcare facilities, such as hospitals or	The ActiGraph CentrePoint Insight Watch is a small worn activity monitor designed for documenting physical movement associated with applications in physiological monitoring. The	The Empatica E4 is intended for passive, non- invasive continuous collection of electrodermal activity that will be later transmitted to a web server for remote review by	The subject device indication for use includes the monitoring of a subset of the physiological parameters of all the predicates.

The information contained in this document is confidential and is the property solely of Empatica, S.r.l.. Any reproduction in part or in whole without the written permission of Empatica, S.r.l. is prohibited. © Empatica S.R.L, 2022. All rights reserved Page 5.4

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
	researchers to	remote review by a	skilled nursing	device is intended	clinicians or	
	remotely monitor	clinician. The Loop	facilities, or their own	to monitor the	researchers.	
	physiologic	System measures and	home. It is intended for	activity associated		
	parameters in	records: • arterial	monitoring of patients	with movement		
	ambulatory	oxygen saturation	by trained healthcare	during sleep. The		
	individuals 18 years	(SpO2) • heart rate	professionals. The	Insight watch can		
	of age and older in	(HR) • respiration rate	Current Wearable	be used to analyze		
	home-healthcare	(RR) All of these	Health Monitoring	circadian rhythms		
	environments.	measurements are	System is intended to	and assess activity		
	The device	made when no motion	provide visual and	in any instance		
	supports	is detected by the	audible physiologic	where quantifiable		
	continuous data	System. The Loop	multi-parameter	analysis of physical		
	collection for	System device does	alarms. The Current	motion is		
	monitoring the	not provide	Wearable Health	desirable.		
	following	physiological alarms.	Monitoring System is			
	physiological		intended for			
	parameters:		temperature			
	<ul> <li>Peripheral skin</li> </ul>		monitoring where			
	temperature,		monitoring			
	<ul> <li>Electrodermal</li> </ul>		temperature at the			
	activity,		upper arm is clinically			
	<ul> <li>Blood Oxygen</li> </ul>		indicated. The Current			
	Saturation under		Wearable Health			
	no motion		Monitoring System is			
	conditions,		intended for			
	<ul> <li>Activity associated</li> </ul>		continuous monitoring			
	with movement		of the following			
	during sleep		parameters in			

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
	The Empatica Health Monitoring Platform can be used to analyze circadian rhythms and assess activity in any instance where quantifiable analysis of physical motion is desirable. The Empatica Health Monitoring Platform is not intended for SpO2 monitoring in conditions of motion or low perfusion.		adults: • Pulse rate • Oxygen saturation • Temperature • Movement The Current Wearable Health Monitoring System is intended for intermittent or spot- check monitoring, in adults, of: • Respiration rate • Non-invasive blood pressure • Lung function & spirometry • Weight The Current Wearable			
	The Empatica Health Monitoring Platform is intended for peripheral skin temperature monitoring, where monitoring		Health Monitoring System is not intended for use in high-acuity environments, such as ICU or operating rooms. The Current Wearable Health Monitoring System is not intended			

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
	temperature at the wrist is clinically indicated.		for use on acutely ill cardiac patients with the potential to develop life threatening arrhythmias e.g. very fast atrial fibrillation. For these patients, they should be monitored using a device with continuous ECG. The Current Wearable Health Monitoring System is not a substitute for an ECG monitor. The Current Wearable Health Monitoring System is not intended for SpO2 monitoring in conditions of high			
Target Population	Adult	Adult	motion or low perfusion Adult	Adult	Adult	The subject device and the predicates are identical

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
Anatomical Site	Wrist	Wrist	Upper Arm	Wrist	Wrist	Clinical testing demonstrated the equivalence between the subject device and the predicates. The difference in wearing location on the body does not raise new questions of safety or efficacy.
Over the Counter or Rx	Rx	Rx	Rx	Rx	Rx	The subject device and the predicates are identical
Environment	Home	Home	Professional Healthcare Facilities & Home	Professional Healthcare Facilities & Home	Professional Healthcare Facilities & Home	The subject device includes a subgroup of the predicates, hence this does not raise new questions of safety or efficacy.

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
Alarms	No	No	Yes	No	No	This difference does not raise new questions of safety or efficacy since the Empatica Health Monitoring Platform is not intended, by design, to include alarms to be used in a situation where the presence of alarms is a requirement for the correct patient care.
User Interface	Device screen, Mobile device application, and cloud software platform	Central station	Mobile devices and a central station	Device screen and Mobile device application	Mobile devices and cloud software platform	The differences between the subject device and the predicates do not raise new questions of safety or efficacy

Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences
Energy Source	Battery	Battery	Battery	Battery	Battery	The subject device and the predicates are identical
Battery Type	Rechargeable Lithium-Ion	Rechargeable Lithium- Ion	Rechargeable Lithium- Ion	Rechargeable Lithium-Ion	Rechargeable Lithium-Ion	The subject device and the predicates are identical
Wireless Communication Interface	Bluetooth® Low Energy (device to mobile device) IEEE 802.11 WiFi/cellular to Empatica cloud	Wireless (cellular connection) via charging station to Spry Server.	IEEE 802.11 WiFi	Bluetooth® Low Energy	Bluetooth® Low Energy (device to mobile device) IEEE 802.11 WiFi/cellular to Empatica cloud	All the devices are designed to transmit their data to alternate devices or sites. The different technologies used do not raise new questions of safety or efficacy
Patient contacting materials	Compliant to ISO 10993-1	Compliant to ISO 10993-1	Compliant to ISO 10993-1	Compliant to ISO 10993-1	Compliant to ISO 10993-1	The subject device and the predicates are identical

	Technical and Performance Information for Blood Oxygen Saturation								
Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences			
Technology	SpO2 relies on the principle that hemoglobin at different oxygenation states absorbs light differently based upon the wavelength of light.	SpO2 measured by analyzing reflectance of certain LED frequencies in a photoplethysmogram design. The diodes are mounted in the device such that they are in contact with the skin	SpO2 is measured by analyzing the reflectance of certain LED frequencies in a photoplethysmogram design. The diodes are mounted in the device such that they are in contact with the skin.	N/A	N/A	The subject device and the predicates are identical in that they all use the photoplethysmogram technology			
SpO <sub>2</sub> Range	70-100%	70-100%	70-100%	N/A	N/A	The subject device and the predicates are identical			
SpO <sub>2</sub> Resolution	1%	1%	1%	N/A	N/A	The subject device and the predicates are identical			
SpO₂ Accuracy	2.6% A <sub>rms</sub>	3% A <sub>rms</sub>	± 2 Digits	N/A	N/A	The subject device and the predicates comply with ISO 80601-2- 61 as well as with FDA Guidance for Pulse Oximeters (2013)			

	Technical and Performance Information for Temperature								
Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences			
Technology	high-precision temperature sensor	N/A	In-Built Thermistor	N/A	N/A	The subject device and the predicates are identical			
Temperature Range	0°C to 50°C	N/A	0°C to 50°C	N/A	N/A	The subject device and the predicates are identical			
Temperature Resolution	0.1°C	N/A	0.1°C	N/A	N/A	The subject device and the predicates are identical			
Temperature Accuracy	± 0.1ºC within 30.0ºC - 45.0ºC range	N/A	±0.1°C	N/A	N/A	The subject device and the predicates are identical			

	Technical and Performance Information for Electrodermal Activity								
Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences			
Technology	EDA is measured by analyzing detected changes in the conductivity of the superficial layers of the skin.	N/A	N/A	N/A	EDA is measured by analyzing detected changes in the conductivity of the superficial layers of the skin.	The subject device and the predicates are identical			
EDA Range	0.01 μS – 100 μS	N/A	N/A	N/A	0.01 μS – 100 μS	The subject device and the predicates are identical			
EDA Resolution	1 digit ~ 55 pS	N/A	N/A	N/A	1 digit ~ 900 pS	This difference shall not raise new concerns of device safety or effectiveness			

Technical and Performance Information for Activity and Sleep							
Features	Empatica Health Monitoring Platform (Subject Device)	Loop System	Current Wearable Health Monitoring System	ActiGraph CentrePoint Insight Watch	Empatica E4	Analysis of differences	
Technology	Accelerometer	N/A	N/A	Accelerometer	N/A	The subject device and the predicate are identical	
Accelerometer Type	Microelectromechanical system (MEMS)-based integrated circuit	N/A	N/A	Microelectromechanical system (MEMS)-based integrated circuit	N/A	The subject device and the predicate are identical	
Accelerometer Sampling Rate	Digital method, 26 Hz – 208 Hz	N/A	N/A	Digital method, 32 Hz – 256 Hz	N/A	This difference shall not raise new concerns of device safety or effectiveness	
Accelerometer Dynamic Range	± 16 g	N/A	N/A	± 8 g	N/A	This difference shall not raise new concerns of device safety or effectiveness	
Accelerometer Sensitivity	0.488 milli-g per Least Significant Bit	N/A	N/A	2.4 milli-g per Least Significant Bit	N/A	This difference shall not raise new concerns of device safety or effectiveness	

**Empatica Health Monitoring Platform** 

#### II. PERFORMANCE DATA

#### Non-Clinical testing (Bench testing)

The following non-clinical (bench) testing was conducted to support a determination of substantial equivalence to the predicates and to demonstrate performance. The non-clinical bench tests included:

Test Name	Test Description	Results
Biocompatibility	The wearable device component of the Empatica Health Monitoring	Passed
testing	Platform, called EmbracePlus, was tested in accordance with the FDA	
	Guidance for Industry and Food and Drug Administration Staff "Use	
	of International Standard ISO 10993-1, "Biological evaluation of	
	medical devices - Part 1: Evaluation and testing within a risk	
	management process" September 4, 2020, and International	
	Standard ISO 10993-1:2018 "Biological Evaluation of Medical Devices	
	– Part 1: Evaluation and Testing Within a Risk Management Process,"	
	as recognized by FDA. The battery of testing included the following	
	tests:	
	Cytotoxicity	
	Sensitization	
	Irritation	
	The EmbracePlus wearable device is considered surface contacting	
	for a prolonged duration (>24 hours < 30 days)	
Electrical safety	The wearable device component of the Empatica Health Monitoring	Passed
testing	Platform, called EmbracePlus, was tested in accordance with	
	International Standard IEC 60601-1 for electrical safety	
Electromagnetic	The wearable device component of the Empatica Health Monitoring	Passed
compatibility	Platform, called EmbracePlus, was tested in accordance with	
(EMC) testing	International Standard IEC 60601-1-2 for EMC	
Wireless Radio	Empatica Health Monitoring Platform was tested to ensure it can	Passed
Communication	communicate via wireless radio in its intended environment in	
	compliance with FDA Radio Frequency Wireless Technology in	
	Medical Devices Guidance, issued August 2013	
Usability testing	The Empatica Health Monitoring Platform was assessed with regards	Passed
	to usability for compliance with IEC 62366-1. The EmbracePlus was	
	also tested in accordance with International Standard IEC 60601-1-11	
	for Usability of medical devices.	
Home-Use	The wearable device component of the Empatica Health Monitoring	Passed
testing	Platform, called EmbracePlus, was tested in accordance with	
	International Standard IEC 60601-1-6 for medical devices used in	
Cloaning	home healthcare environments.	Passed
Cleaning validation	The wearable device component of the Empatica Health Monitoring Platform, called EmbracePlus, was tested in accordance with	Passeu
valluation	International Standard ISO 17664 and AAMI TIR 30 to assess device	
	cleaning procedure	
Manual	The wearable device component of the Empatica Health Monitoring	Passed
disinfection	Platform, called EmbracePlus, was tested in accordance with ASTM	1 03360
	E1837:2014 and AAMI TIR 12 to assess device low-level disinfection	
	procedure	
Temperature	The Empatica Health Monitoring Platform was tested to confirm the	Passed
measurement	Skin temperature measurement accuracy and transient time complies	
accuracy	with ISO 80601-2-56 Medical electrical equipment - Part 2-56:	

The information contained in this document is confidential and is the property solely of Empatica, S.r.l.. Any reproduction in part or in whole without the written permission of Empatica, S.r.l. is prohibited. © Empatica S.R.L, 2022. All rights reserved

#### **Empatica Health Monitoring Platform**

Test Name	Test Description					
	Particular requirements for basic safety and essential performance of					
	clinical thermometers for body temperature measurement.					
	[Including: Amendment 1 (2018)] to assess its accuracy.					
Electrodermal	The Empatica Health Monitoring Platform computed electrodermal					
activity	activity (EDA) was tested to determine its equivalence to the					
measurement	predicate device Empatica E4.					
Activity	Bench testing has been performed to demonstrate the equivalence of	Passed				
Counts/Sleep	the Empatica Health Monitoring Platform activity counts and sleep					
	detection with the predicate device.					

#### Software Verification and Validation Testing

Software verification and validation testing were conducted. Documentation was provided as recommended by the FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." All the Empatica Health Monitoring Platform software components were considered a "moderate" level of concern since a failure or latent flaw in the software could result in minor injury to the patient or operator.

#### Cybersecurity

Cybersecurity activities was conducted. Documentation was provided as recommended by the FDA's Guidance for Industry and FDA Staff, " Content of Premarket Submissions for Management of Cybersecurity in Medical Devices." All the Empatica Health Monitoring Platform software components underwent appropriate cybersecurity assessment and testing.

#### Animal study

No animal studies were conducted as part of the submission to prove substantial equivalence.

#### **Clinical Study**

A human clinical investigation study was conducted to demonstrate the performance of the Empatica Health Monitoring Platform.

The clinical study investigated the accuracy of the blood oxygen saturation monitoring in 13 healthy adult subjects with heterogeneous skin types. The Empatica Health Monitoring Platform was compared to the gold standard, arterial blood gas analysis.

This testing was conducted in accordance with ISO 80601-2-61 Medical electrical equipment - Part 2-61: Particular requirements for basic safety and essential performance of pulse oximeter equipment and in accordance with the FDA Guidelines for Pulse Oximeters – Premarket Notification Submissions [510(k)s]: Guidance for Industry and Food and Drug Administration Staff (2013).

The Empatica Health Monitoring Platform was found to be in compliance with both documents. This testing demonstrated an accuracy of  $2.6\% A_{rms}$  across the SpO<sub>2</sub> range of 70-100%. This testing was not conducted in the presence of motion or low perfusion.

No adverse events related to the device were encountered during the execution of both studies. The results of the clinical investigations demonstrate an effectiveness profile similar to the predicate devices.

#### III. CONCLUSION

Based on the information presented in this 510(k) premarket notifications, the Empatica Health Monitoring Platform is substantially equivalent to the predicate devices. The Empatica Health Monitoring Platform is as safe and effective as the currently marketed predicate devices.

**Empatica Health Monitoring Platform** 

Based on testing and comparison with the predicate devices, the Empatica Health Monitoring Platform indicated no adverse indications or results. It is our determination that the Empatica Health Monitoring Platform is safe, effective and performs within its design specifications, and is substantially equivalent to the predicate devices.

The information contained in this document is confidential and is the property solely of Empatica, S.r.l.. Any reproduction in part or in whole without the written permission of Empatica, S.r.l. is prohibited. © Empatica S.R.L, 2022. All rights reserved Page 5.16