

June 26, 2023

Z7 LLC % Linda Braddon Chief Executive Officer Secure BioMed Evaluations 7828 Hickory Flat Highway Suite 120 Woodstock, Georgia 30188

Re: K221488

Trade/Device Name: Z7 Zirconia Implant System Regulation Number: 21 CFR 872.3640 Regulation Name: Endosseous Dental Implant Regulatory Class: Class II Product Code: DZE, NHA Dated: May 25, 2023 Received: May 25, 2023

Dear Linda Braddon:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Andrew I. Steen -S

Andrew I. Steen Assistant Director DHT1B: Division of Dental and ENT Devices OHT1: Office of Ophthalmic, Anesthesia, Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K221488

Device Name Z7 Zirconia Implant System

Indications for Use (Describe)

The Z7 Zirconia Implant System is intended for surgical placement in the patient's upper and lower jaw to provide support for prosthetic devices, such as artificial teeth and in order to restore the patient chewing function. The implants are indicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The ø3.7 mm reduced diameter implants are recommended for central and lateral incisors only.

Type of Use (Select one or both, as applicable)	
⊠ Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K221488 **510(k) SUMMARY** Z7 Zirconia Implant System

Date Prepared	June 26, 2023		
^	Z7, LLC		
G	28 Felmley Road		
Sponsor	Whitehouse Station, NJ		
	908-399-1359		
	Secure BioMed Evaluations		
	Linda Braddon, Ph.D.		
	7828 Hickory Flat Highway		
510(k) Contact	Suite 120		
	Woodstock, GA 30188		
	770-837-2681		
	Regulatory@SecureBME.com		
Trade Name	Z7 Zirconia Implant System		
Common Name	Endosseous Dental Implant		
Code –	DZE, NHA		
Classification	21 CFR 872.3640: Class II		
Primary Predicate	K192053 TAV Medical Ltd. W Zirconia Implants		
Reference Device K172668 TAV Medical Ltd. W Zirconia Implants			
Kelefence Device	K132585 Zibone Ceramic Dental Implant System		
Device Description	The Z7 Zirconia Implant System is an integrated system of endosseous dental		
	implants (Z7 Zirconia Implant One Piece) and PEEK prosthetic parts. The Z7		
	Zirconia Implant One Piece are yttria stabilized tetragonal zirconia (Y-TZP)		
	dental implants composed of a One Piece, monotype implant with an		
	integrated abutment. The implant is manufactured via a ceramic injection		
	molding with the macro and micro surface characteristics of the implant		
	directly structured in the mold. The implant body portion is configured to		
	extend into the bone and osseo-integrate with the alveolar bone. The neck		
	should be positioned 1.8mm above the bone. The implants come in		
	corresponding diameters of 3.7 and 4.3 mm.		
Indications for Use	The Z7 Zirconia Implant System is intended for surgical placement in the		
Statement	patient's upper and lower jaw to provide support for prosthetic devices, such		
	as artificial teeth and in order to restore the patient chewing function. The		
	implants are indicated for immediate loading when good primary stability is		
	achieved and with appropriate occlusal loading.		
	The Q2 Term reduced dispersion implete an according to the term		
	The Ø3.7mm reduced diameter implants are recommended for central and		
	lateral incisors only.		

Comparison of Indications for Use

Subject Device Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053
Intended for surgical placement in the patient's upper and lower jaw to provide support for prosthetic devices, such as artificial teeth and in order to restore the patient chewing function. The implants are indicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading.	Intended for surgical placement in the patient's upper and lower jaw to provide support for prosthetic devices, such as artificial teeth and in order to restore the patient chewing function. The implants are indicated for immediate loading when good primary stability is achieved and with appropriate occlusal loading
The ø3.7 mm reduced diameter implants are recommended for central and lateral incisors only.	The ø3.6 mm reduced diameter implants are recommended for central and lateral incisors only.

There are no significant differences between the subject and predicate device indications for use. The minor differences in wording to align with the provided sizes of the subject device do not change the intended use for demonstration of substantial equivalence.

Characteristic	Subject Device Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053	Reference Device TAV Medical's W Zirconia Implants K172668
Purpose of Device	Subject Device	Same Indications	Similar One-Piece Sizes as Subject Device
Device Classification	Class II	Class II	Class II
Product Code	DZE, NHA	DZE, NHA	DZE, NHA
Material	Yttria-stabilized zirconia (Y-TZP) per ISO 13356:2015	Yttria-stabilized zirconia (Y-TZP) per ISO 13356:2015	Yttria-stabilized zirconia (Y-TZP) per ISO 13356:2015
Surface Topography	Macro and Micro Roughness	Macro and Micro Roughness	Macro and Micro Roughness
Description	Endosseous dental implant	Endosseous dental implant	Endosseous dental implant
Design	One-Piece Ceramic Implants	One-Piece Ceramic Implants	One-Piece Ceramic Implants
	N/A	Two-Piece Ceramic Implants	Two-Piece Ceramic Implants
Diameter	Z7 One-Piece: 3.7mm, 4.3mm	TAV W One-Piece: 3.6mm	TAV W One-Piece: 4.1mm, 4.8mm

Comparison of Technological Characteristics - Implants

Characteristic	Subject Device Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053	Reference Device TAV Medical's W Zirconia Implants K172668
	N/A	TAV W Two-Piece:	TAV W Two-Piece:
		4.1mm, 4.8mm	4.1mm, 4.8mm
	Z7 One-Piece: 10mm,	TAV W One-Piece:	TAV W One-Piece:
	13mm	8mm, 10mm, 12mm,	8mm, 10mm, 12mm,
Longth	1511111	14mm	14mm
Length		TAV W Two Piece:	TAV W Two Piece:
	N/A	8mm, 10mm, 12mm,	8mm, 10mm, 12mm,
		14mm	14mm
Manufacturing	CIM: Ceramic	CIM: Ceramic	CIM: Ceramic
Technology	injection molding	injection molding	injection molding
Sterilization Method	Ethylene Oxide	Gamma irradiation	Gamma irradiation
Intended Use Environment	Dental clinic setting	Dental clinic setting	Dental clinic setting
	Straight	TAV W One-Piece:	TAV W One-Piece:
Abutments		Straight	Straight
Adutments	N/A	TAV W Two-Piece:	TAV W Two-Piece:
		Straight and Angled	Straight and Angled
	Implant device with	Implant device with	Implant device with
Diocompatibility	permanent (>30 days)	permanent (>30 days)	permanent (>30 days)
Biocompatibility	contact with tissue /	contact with tissue /	contact with tissue /
	bone / dentin	bone / dentin	bone / dentin
Shelf-Life	5 years	Not stated in 510(k) Summary	1 year

Comparison of Technological Characteristics – PEEK Healing Caps

Characteristic	Subject Device MABB Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053	Reference Device TAV Medical's W Zirconia Implants K172668
Regulation Number	872.3640	872.3640	872.3640
Classification Product Code	DZE, NHA	DZE, NHA	DZE, NHA
Product Name	PEEK Healing Cap	PEEK Healing Cap for W One Piece Zirconia Implant	PEEK Healing Cap for W One Piece Zirconia Implant
Product Description	Protect the implant during the healing phase up to 180 days.	Protect the implant during the healing phase up to 180 days.	Protect the implant during the healing phase up to 180 days.

Characteristic	Subject Device MABB Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053	Reference Device TAV Medical's W Zirconia Implants K172668
Material	PEEK per ASTM F2026-17	PEEK	PEEK
Diameter (mm)	3.7, 4.3	3.6, 4.1, 4.8	4.1, 4.8
Height (mm)	5	5	5
Angle (°)	0° (Straight)	0° (Straight)	0° (Straight)
Sterility	Non-Sterile	Non-Sterile	Non-Sterile

Comparison of Technological Characteristics – PEEK Temporary Abutments

Characteristic	Subject Device MABB Z7 Zirconia Implant System	Primary Predicate Device TAV Medical's W Zirconia Implants K192053	Reference Device TAV Medical's W Zirconia Implants K172668
Regulation Number	872.3640	872.3640	872.3640
Classification Product Code	DZE, NHA	DZE, NHA	DZE, NHA
Product Name	PEEK Healing Cap	PEEK Healing Cap for W One Piece Zirconia Implant	PEEK Healing Cap for W One Piece Zirconia Implant
Product Description	Serves as a basis for temporary restoration for crown. Up to 180 days.	Serves as a basis for temporary restoration for crown or bridge. Up to 180 days.	Serves as a basis for temporary restoration for crown or bridge. Up to 180 days.
Material	PEEK per ASTM F2026-17	PEEK	PEEK
Diameter (mm)	3.7, 4.3	3.6, 4.1, 4.8	4.1, 4.8
Angle (°)	0° (Straight)	0° (Straight)	0° (Straight)
Sterility	Non-Sterile	Non-Sterile	Non-Sterile

Technological Characteristics

There are no significant technological differences between the subject and predicate devices. Minor differences include the method of sterilization, and the exact size offerings. The subject device uses the same material as the additional predicates and all devices are manufactured using ceramic injection molding techniques. Any technological differences between the subject and predicate devices do not raise new concerns for safety or effectiveness.

Non-Clinical Performance Testing Summary

Test	Test Method Summary	Results
Material Performance Testing	Per ISO 13356:2015 Implants for surgery — Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP)	PASS Device met all predetermined acceptance criteria
Fatigue Testing	Per ISO 14801:2016 Dentistry — Implants — Dynamic loading test for endosseous dental implants	PASS Substantially equivalent fatigue strength to reference device (K132585)
Implant Surface Roughness and	Surfaces macro and micro morphological characterization through 3D Scanning Electron	PASS
Chemical	Microscope (SEM) Imaging. Surface chemical	Device met all
Analysis	analysis through Energy Dispersive	predetermined
Validation	Spectrometry (EDS Analysis).	acceptance criteria

Performance testing for the subject device is summarized in the table below:

Biological evaluation of the subject device was performed according to ISO 10993-1 (cytotoxicity testing per ISO 10993-5). Endotoxin testing on the subject device or suitable test specimens was performed following USP<85> and USP<161> according to the sponsor's endotoxin sampling plan.

The implants are provided sterile to the end user via Ethylene Oxide (EO) sterilization. The applicable cycle was validated in compliance with ISO 11135:2014 "Sterilization of health-care products — Ethylene oxide — Requirements for the development, validation and routine control of a sterilization process for medical devices" to a sterility assurance level of 10⁻⁶ for the worst-case configuration of the subject device. All cycle parameters were within the limits established by the protocol and residue testing yielded acceptable results for EO and ECH (ethylene chlorohydrin) per ISO 10993-7:2008 "Biological evaluation of medical devices — Part 7: Ethylene oxide sterilization residuals". The subject device non-sterile components are end-user steam sterilized. The applicable steam sterilization cycle was validated in compliance with ISO 17665 "Sterilization of health care products -- Moist heat -- Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices" and demonstrated a sterility assurance level of 10⁻⁶.

Accelerated aging has been applied on the final sterile implant packaging and is being followed by real time aging to validate a five-year shelf life.

Conclusions

Based on the similarities of the intended use/indications for use, technological and functional characteristic, and the results of the non-clinical performance testing, the subject device is substantially equivalent to the legally marketed predicate device.