

Manteia Technologies Co., Ltd. % Dandan Chen RA 1903, B Tower, Zijin Plaza No. 1811 Huandao East Road Xiamen, 361001 CHINA

Re: K221706

Trade/Device Name: AccuContour Regulation Number: 21 CFR 892.2050 Regulation Name: Medical image management and processing system Regulatory Class: Class II Product Code: QKB Dated: May 31, 2022 Received: June 13, 2022

Dear Dandan Chen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

March 9, 2023

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Lora D. Weidner, Ph.D. Assistant Director Radiation Therapy Team DHT8C: Division of Radiological Imaging and Radiation Therapy Devices OHT8: Office of Radiological Health Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K221706

Device Name AccuContour

Indications for Use (Describe)

It is used by radiation oncology department to register multi-modality images and segment (non-contrast) CT images, to generate needed information for treatment planning, treatment evaluation and treatment adaptation.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

I . SUBMITTER

Manteia Technologies Co., Ltd. 1903, B Tower, Zijin Plaza, No.1811 Huandao East Road, Xiamen

China Establishment Registration Number: 3016686005

Contact Person: Dandan Chen Position: RA Tel: +86 592-6100813 Email: <u>chendandan@manteiatech.com</u>

Date Prepared: March 9, 2023

II. DEVICE

Name of Device: AccuContour Common or Usual Name: Medical Imaging Software Classification Name: System, Imaging processing, Radiological Regulatory Class: II Product Code: QKB Regulation Number: 21CFR 892.2050 Review Panel: Radiology

III. PREDICATE DEVICE

Device	510(k) Number	Product Name
Predicate Device	K191928	AccuContour TM
Reference Device	K182624	MIM-MRT Dosimetry
Reference Device	K173636	Velocity
Reference Device	K181572	Workflow Box

IV. DEVICE DESCRIPTION

The proposed device, AccuContour, is a standalone software which is used by radiation oncology department to register multi-modality images and segment (non-contrast) CT images, to generate needed information for treatment planning, treatment evaluation and treatment adaptation.

The product has two image processing functions:

(1) Deep learning contouring: it can automatically contour organs-at-risk, in head and neck,

thorax, abdomen and pelvis (for both male and female) areas,

- (2) Automatic registration: rigid and deformable registration, and
- (3) Manual contouring.

It also has the following general functions:

- Receive, add/edit/delete, transmit, input/export, medical images and DICOM data;
- Patient management;
- ➢ Review of processed images;
- ➤ Extension tool;
- Plan evaluation and plan comparison;
- \triangleright Dose analysis.

V. INDICATIONS FOR USE

It is used by radiation oncology department to register multi-modality images and segment (non-contrast) CT images, to generate needed information for treatment planning, treatment evaluation and treatment adaptation.

VI. SUBSTANTIALLY EQUIVALENT (SE) COMPARISION

Table 1 Comparison of Technology Characteristics

ITEM	Proposed Device	Predicate Device	Reference Device	Reference Device	Reference Device
		K191928	K182624	K173636	K181572
Regulatory Inform	nation				
Regulation No.	21CFR 892.2050	21CFR 892.2050	21CFR 892.2050	21CFR 892.2050	21CFR 892.2050
Product Code	QKB	QKB	LLZ	LLZ	LLZ
Indication for	It is used by radiation	It is used by radiation oncology	MIM software is used by trained	Velocity is a software package that	Workflow Box is a software
Use	oncology department	department to register	medical professionals as a tool	provides the physicians a means	system designed to allow users
	to register	multi-modality images and	to aid in evaluation and	for comparison of medical data	to route DICOM-compliant data
	multi-modality images	segment (non-contrast) CT	information management of	including imaging data that is	to and from automated
	and segment	images, to generate needed	digital medical images. The	DICOM compliant.	processing components.
	(non-contrast) CT	information for treatment	medical image modalities	It allows the display, annotation,	Supported modalities include
	images, to generate	planning, treatment evaluation	include, but are not limited to,	volume operation, volume	CT, MR, RTSTRUCT.
	needed information for	and treatment adaptation.	CT, MRI, CR, DX, MG, US,	rendering, registration, and	Workflow Box includes
	treatment planning,		SPECT, PET and XA as	fusion of medical images as an	processing components for
	treatment evaluation		supported by ACR/NEMA	aid during use by diagnostic	automatically contouring
	and treatment		DICOM 3.0. MIM assists in the	radiology, oncology, radiation	imaging data using deformable
	adaptation.		following indications:	therapy planning and other	image registration to support
			• Receive, transmit, store,	medical specialties. Velocity is	atlas based contouring,
			retrieve, display, print, and	not intended for mammography.	re-contouring of the same patient
			process medical images		and machine learning based
			and DICOM objects.		contouring.
			• Create, display and print		Workflow Box is a data routing
			reports from medical		and image processing tool which
			images.		automatically applies contours to

	•	Registration, fusion	data which is sent to one or more
		display, and review of	of the included image processing
		medical images for	workflows. Contours generated
		diagnosis, treatment	by Workflow Box may be used
		evaluation, and	as an input to clinical workflows
		treatment planning.	including, but not limited to,
	•	Evaluation of cardiac left	radiation therapy treatment
		ventricular function and	planning.
		perfusion, including left	Workflow Box must be used in
		ventricular enddiastolic	conjunction with appropriate
		volume, end-systolic	software to review and edit
		volume, and	results
		ejection fraction.	generated automatically by
	•	Localization and definition	Workflow Box components, for
		of objects such as tumors	example image visualization
		and normal tissues in	software
		medical images.	must be used to facilitate the
	•	Creation, transformation,	review and edit of contours
		and modification of	generated by Workflow Box
		contours for applications	component
		including, but not limited	applications.
		to, quantitative analysis,	Workflow Box is intended to be
		aiding adaptive therapy,	used by trained medical
		transferring contours to	professionals.
		radiation therapy	Workflow Box is not intended to
		transferring contours to	automatically detect lesions.

		radiation therapy	
		transferring contours to	
		radiation therapy treatment	
		planning systems and	
		plaining systems, and	
		archiving contours for	
		patient follow-up and	
		management.	
	•	Quantitative and statistical	
		analysis of PET/SPECT	
		brain scans by comparing	
		to other registered	
		PET/SPECT brain scans.	
	•	Planning and evaluation of	
		permanent implant	
		brachytherapy procedures	
		(not including radioactive	
		microspheres).	
	•	Calculating absorbed	
		radiation dose as a result of	
		administering a	
		radionuclide When using	
		device elipically the user	
		device chinically, the user	
		snould only use FDA	
		approved	
		radiopharmaceuticals. If	
		using with unapproved	

			.1. 1 . 1 . 1 .		
			ones, this device should		
			only be used for research		
			purposes.		
			Lossy compressed		
			mammographic images and		
			digitized film screen images		
			must not be reviewed for		
			primary image interpretations.		
			Images that are printed to film		
			must be printed using an		
			FDA-approved printer for the		
			diagnosis of digital		
			mammography images.		
			Mammographic images must be		
			viewed on a display system that		
			has been cleared by the FDA for		
			the diagnosis of digital		
			mammography images. The		
			software is not to be used for		
			mammography CAD.		
Label/labeling	Conform with 21CFR	Conform with 21CFR Part 801	Conform with 21CFR Part 801	Conform with 21CFR Part 801	Conform with 21CFR Part 801
	Part 801				
Operating	Windows	Windows	Windows and MAC system	Windows and MAC system	Windows
System					
Segmentation Fea	1	1		1	
	itures				

					registration based
					re-contouring, machine learning
					based contouring
Compatible	Non-Contrast CT	Non-Contrast CT	Non-Contrast CT	Non-Contrast CT	CT, MR
Modality					
Compatible	No Limitation on	No Limitation on	No Limitation on scanner	No Limitation on scanner	No Limitation on scanner
Scanner Models	scanner model,	scanner model,	model,	model,	model,
	DICOM 3.0	DICOM 3.0 compliance	DICOM 3.0 compliance	DICOM 3.0 compliance required.	DICOM 3.0 compliance
	compliance required.	required.	required.		required.
Compatible	No Limitation on TPS	No Limitation on TPS model,	No Limitation on TPS model,	No Limitation on TPS model,	No Limitation on TPS model,
Treatment	model, DICOM	DICOM	DICOM	DICOM	DICOM
Planning	3.0 compliance	3.0 compliance required.	3.0 compliance required.	3.0 compliance required.	3.0 compliance required.
System	required.				
Unattended	Yes	No	Not stated	No	Yes
workstation					
Registration Feature	ires			•	
Algorithm	Intensity Based.	Intensity Based.	Intensity Based.	Intensity Based.	Intensity Based.
Image	Auto rigid registration	Auto rigid registration	Auto rigid registration and	Auto rigid registration and	Auto rigid registration and
registration	and auto deformable		deformable registration.	deformable registration.	deformable registration.
	registration.				
Compatible	Auto rigid registration:	CT, MRI, PET	CT, MRI, CR, DX, MG, US,	PET/SPECT/CT/MRI	CT, MRI
Modality	CT, MRI, PET		SPECT, PET and XA		
	Auto deformable				
	registration: CT, MRI,				
	CBCT				

Compatible	No Limitation on	No Limitation on	No Limitation on scanner	No Limitation on scanner	No Limitation on scanner
Scanner Models	scanner model,	scanner model,	model,	model,	model,
	DICOM 3.0	DICOM 3.0 compliance	DICOM 3.0 compliance	DICOM 3.0 compliance required.	DICOM 3.0 compliance
	compliance required.	required.	required.		required.
Compatible	No Limitation on	No Limitation on	No Limitation on scanner	No Limitation on scanner	No Limitation on scanner
_	scanner model,	scanner model,	model,	model,	model,
Treatment	DICOM 3.0	DICOM 3.0 compliance	DICOM 3.0 compliance	DICOM 3.0 compliance required.	DICOM 3.0 compliance
Planning	compliance required.	required.	required.		required.
System					
Plan Evaluation F	eature				
Display of	Yes	No	Not stated	Yes	No
DICOM RT					
Plans					
Isodose Line	Yes	No	Not stated	Yes	No
Display					
DVH statistics	Yes	No	Not stated	Yes	No
display					
RT Plans	Yes	No	Not stated	Yes	No
comparison					
Dose Analysis Fea	ature				
Display of	Yes	No	Not stated	Yes	No
DICOM RT					
Doses					
Dose	Yes	No	Not stated	Yes	No
accumulation					
L	1	1	4	1	1

VII. PERFORMANCE DATA

The following performance data were provided in support of the substantial equivalence determination.

1. Non-Clinical Test Conclusion

Deformable registration performance test

The registration performance test was performed on proposed device and reference device (K182624) to evaluate the deformable registration accuracy. All fixed images and moving images are generated in healthcare institutions in U.S. The scanner models covered products from five major vendors. The image registration feature is tested on multi-modality image sets from different patients. The Normalized Mutual Information (NMI) was used for evaluation. NMI values were calculated on two sets of images for both the proposed device and reference device (K182624), respectively. The NMI value of proposed device was compared with that of the reference device. According to the results, it could be concluded that the NMI of proposed device was non-inferior to that of the reference device.

2. Clinical Test Conclusion

No clinical study is included in this submission.

3. Software Verification and Validation Testing

Software verification and validation testing were conducted, and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The software for this device was considered as a "major" level of concern.

Software bench testing in the form of Unit, System and Integration tests were performed to evaluate the performance and functionality of the new features and software updates. Software verification and regression testing have been performed successfully to meet their previously determined acceptance criteria as stated in the test plans.

WI. SUBSTANIALLY EQUIVALENT (SE) CONCLUSION

The proposed device is substantially equivalent to the predicate device AccuContourTM (K191928).