

November 3, 2022

Aegis Spine, Inc. Kihyang Kim Regulatory Affairs 9781 S. Meridian Blvd, Ste 300 Englewood, Colorado 80112

Re: K221719

Trade/Device Name: ABTross ALIF Expandable Cage System Regulation Number: 21 CFR 888.3080 Regulation Name: Intervertebral Body Fusion Device Regulatory Class: Class II Product Code: MAX Dated: October 3, 2022 Received: October 5, 2022

Dear Kihyang Kim:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Brent Showalter, Ph.D. Assistant Director DHT6B: Division of Spinal Devices OHT6: Office of Orthopedic Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number *(if known)* K221719

Device Name

ABTross ALIF Expandable Cage System

Indications for Use (Describe)

ABTross ALIF Expandable Cage System is indicated for intervertebral body fusion of the spine in skeletally mature patients. The device system is designed for use with bone graft (autograft bone and/or allogenic bone graft composed of cancellous and/or corticocancellous bone graft) to facilitate fusion. ABTross ALIF Expandable Cage System is intended for use at either one level or two contiguous levels in the lumbar spine, from L2 to S1, for the treatment of degenerative disc disease (DDD) with up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. Patients should have at least six months of non-operative treatment prior to treatment with a lumbar intervertebral fusion device. The device is intended to be used with supplemental fixation systems that have been cleared for use in the lumbosacral spine (e.g., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems). Hyperlordotic interbody devices ($\geq 20^{\circ}$ lordosis) must be used with at least anterior supplemental fixation.

Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)	
CONTINUE ON A SEPAR	RATE PAGE IF NEEDED.	
This section applies only to requirements	of the Paperwork Reduction Act of 1995.	
DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.		
The burden time for this collection of information is esti time to review instructions, search existing data source and review the collection of information. Send commer of this information collection, including suggestions for	nts regarding this burden estimate or any other aspect	
Department of Hea Food and Drug Ad	alth and Human Services	

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K221719 - 510(k) SUMMARY

The following 510(k) summary is being submitted as required by 21 CFR 807.92(a):

1. SUBMITTER

SUDWITTER	
Submitter's Name:	Aegis Spine, INC.
Submitter's Address:	9781 S. Meridian Blvd, Ste 300
	Englewood, CO 80112
Submitter's Telephone:	+1.303.741.4123
Contact Person:	Kihyang Kim
	+1.303.741.4123
	khkim@aegisspine.com/khkim3747@gmail.com

2. DEVICE NAME

Trade or Proprietary Name	ABTross ALIF Expandable Cage System
Common or Usual Name	Intervertebral Body Fusion Device, Intervertebral cage, Spacer
Regulation class / Number	ClassII, CFR 888.3080
Regulation Name	Intervertebral Body Fusion Device
Product Code	MAX
Classification Panel	Spinal Devices (DHT6B)

3. PREDICATE DEVICE

The subject ABTross ALIF Expandable Cage System is substantially equivalent to the following devices:

Primary Predicate Device

• MAGNIFY TM Spacer System(K142498)

Additional Predicate Devices

- XYPAN Expandable Lumbar Cage (XL, XTP) System(K203531)
- Xenco Medical Lumbar Interbody System(K143158)
- Vu a POD Intervertebral Body Fusion Device (K080822/K121211/K173606)
- AnyPlus PEEK Cages (ALIF) (K100516/K111354)

4. **DESCRIPTION OF THE DEVICE**

The ABTross ALIF Expandable Cage System is interbody fusion devices. This cage system is made of Titanium 6AL-4V Alloy (ASTM F136). And cages are offered in a variety of widths, lengths, heights and lordotic angles designed to adapt to a variety of patient anatomies.

The cage can be expanded in height using the system instrument after being inserted in the unexpanded state. The cages have serrations on the superior endplate and inferior endplate surfaces area to contact vertebrae bone endplate. The device system is designed for use with bone graft (autograft bone and/or allogenic bone graft composed of cancellous and/or corticocancellous bone graft) to facilitate fusion. The device is intended to be used with supplemental fixation systems that have been cleared for use in the lumbrosacral spine (e.g., posterior pedicle screw and rod systems, anterior or lateral plate systems, and anterior screw and rod systems). All implants are provided sterile and intended for SINGLE USE ONLY and should not be reused under any circumstances.

5. INDICATION FOR USE

ABTross ALIF Expandable Cage System is indicated for intervertebral body fusion of the spine in skeletally mature patients. The device system is designed for use with bone graft (autograft bone and/or allogenic bone graft composed of cancellous and/or corticocancellous bone graft) to facilitate fusion. ABTross ALIF Expandable Cage System is intended for use at either one level or two contiguous levels in the lumbar spine, from L2 to S1, for the treatment of degenerative disc disease (DDD) with up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. Patients should have at least six months of non-operative treatment prior to treatment with a lumbar intervertebral fusion device. The device is intended to be used with supplemental fixation systems that have been cleared for use in the lumbosacral spine (e.g., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems). Hyperlordotic interbody devices ($\geq 20^{\circ}$ lordosis) must be used with at least anterior supplemental fixation.

6. **PERFORMANCE DATA**

Performance testing was performed to demonstrate that the subject ABTross ALIF Expandable Cage System is substantially equivalent to other predicate devices.

Static compression, dynamic compression, static and dynamic shear testing according to ASTM F2077, was presented to demonstrate the substantial equivalency of the ABTross ALIF Expandable Cage System to the predicate devices.

- Static Axial Compression Test ASTM F 2077 -18
- Static Compression-Shear Test ASTM F 2077 -18
- Static Expulsion Test ASTM draft F-04.25.02.02
- Static Subsidence Test ASTM F 2267 04 (Reapproved 2018)
- Dynamic Axial Compression Test-ASTM F 2077 -18
- Dynamic Compression-Shear Test ASTM F 2077 -18

Bench testing to evaluate the mechanical properties of the ABTross ALIF Expandable Cage System showed a higher or similar mechanical value than predicate marketed devices.

7. MATERIAL

The ABTross ALIF Expandable Cage System is manufactured from Ti-6Al-4V ELI titanium alloy (ASTM F136). This this is the same material used in the predicate devices.

8. COMPARISON OF TECHNICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

The subject and predicate devices have nearly identical technological characteristics and the minor differences do not raise any new issues of safety and effectiveness. Specifically, the following characteristics are identical between the subject and predicates:

- Instruction for use
- Design
- Expanding Mechanism
- Material
- Approach
- Sterilization & Method

The following technological similarities and differences exist between the subject and predicate devices:

INDICATION FOR USE		
This Subject ABTross ALIF Expandable Cage System / MAX	 with Not indicated for standalone use and must be used with supplemental fixation. Hyper lordotic degree (≥20°lordosis) must be used with at least anterior supplemental fixation. 	
Predicate Devices	Similarities	Differences
MAGNIFY TM Spacers MAGNIFY™ S Spacers K142498 / MAX, OVD	 Below both indications can be used MAGNIFY TM Spacers is not indicated for stand-alone use and must be used with supplemental fixation. 	 MAGNIFY[™] S Spacers is stand-alone devices.
XYPAN Expandable Lumbar Cage (XL, XTP) System K203531 / MAX	 Not indicated for standalone use and must be used with supplemental fixation. 	-
Xenco Medical Lumbar Interbody System K143158/ MAX	 Not indicated for standalone use and must be used with supplemental fixation. 	-
Vu a POD Intervertebral Body Fusion Device K080822/K121211/K173606 /MAX, OVD	 Below both indications can be used Not indicated for standalone use and must be used with supplemental fixation (without use SpinePlate and bone screws) Hyper lordotic degree (≥20°lordosis) must be used with at least anterior supplemental fixation. 	 When used with provided SpinePlate and bone screws, is a stand-alone device.
AnyPlus® Peek Cage K100516/K111354 /MAX	 Not indicated for standalone use and must be used with supplemental fixation. 	-

DESIGN		
This Subject ABTross ALIF Expandable Cage System	 Lordotic rectangular cylindrical shape and has an open space inside for bone graft for Anterior approach Lordotic Min. / Max. angle: 8° ~ 30° Depth / width (mm) / Height :24x32 / 29x38 / 6 ~19 	
Predicate Devices	Similarities	Differences
MAGNIFY TM Spacers MAGNIFY™ S Spacers K142498	 Lordotic rectangular shape and has an open space inside for bone graft for Anterior approach Depth / width / Height: similar 	 Provided screws for stand-alone indication Lordotic Min. /Max. angle: 8° ~ 15°
XYPAN Expandable Lumbar Cage (XL, XTP) System K203531	 Lordotic rectangular cylindrical shape and has an open space inside for bone graft for lateral approach. Depth / width / height: similar / bigger 	 Curved (banana shape) shape and has an open space inside for bone graft for Anterior to Psoas approach. Lordotic Min. /Max. angle: 0° ~ 15°
Xenco Medical Lumbar Interbody System K143158	 Rectangular cylindrical shape and has an open space inside for bone graft for Anterior approach Hyper Lordotic angle: 20° ~ 30° 	-
Vu a POD Intervertebral Body Fusion Device K080822/K121211/K173606	 Rectangular cylindrical shape and has an open space inside for bone graft. Depth / width / Height: similar 	 Provided SpinePlate and bone screws for stand-alone indication Lordotic Min. /Max.angle: 8° ~ 15°
AnyPlus® Peek Cage (ALIF) K100516/K111354	 Rectangular cylindrical shape and has an open space inside for bone graft. Depth / width / Height: similar 	• Lordotic Min. /Max.angle: $0^{\circ} \sim 11^{\circ}$

EXPANDING MECHANISM		
EXPANDING MECHANISM This Subject ABTross ALIF Expandable Cage System Predicate Devices MAGNIFY TM Spacers MAGNIFY™ S Spacers/K142498 XYPAN Expandable Lumbar Cage (XL, XTP) System K203531	 Expandable Height Adjusting Struct with Height adjusting Driver, the an guide block are close to each other, Similarities Expandable Devices - Height adjusting in mm 8~11/10~13/ 2~15/14~17/16~19/18~21 Expandable Devices- Height adjusting mm 6~10/7~11/ 8~12/10~14/12~16 Expandable Height Adjusting Structural Mechanism: tightening 	terior guide block and the posterior
	the bolt with Height adjusting Driver, the anterior guide block and the posterior guide block are close to each other, pushing up the plates.	
Xenco Medical Lumbar Interbody System K143158	-	Non-Expandable Devices
Vu a POD Intervertebral Body Fusion Device K080822/K121211/K173606	-	Non-Expandable Devices
AnyPlus® Peek Cage 100516/K111354	-	Non-Expandable Devices

MATERIAL		
This Subject ABTross ALIF Expandable	Ti-6Al-4V ELI titanium alloy (ASTM F136)	
Cage System		
Predicate Devices	Similarities Differences	
MAGNIFY TM Spacers	Ti-6Al-4V ELI titanium alloy	
MAGNIFY [™] S Spacers/K142498	(ASTM F136)	-
XYPAN Expandable Lumbar Cage (XL,	Ti-6Al-4V ELI titanium alloy	
XTP) System K203531	(ASTM F136)	-
Xenco Medical Lumbar Interbody		PEEK (ASTM F 2026) with
System K143158	-	Tantalum markers (ASTM F560)
Vu a POD Intervertebral Body Fusion Device K080822/K121211/K173606	Ti-6Al-4V ELI titanium alloy	PEEK (ASTM F 2026) with
	(ASTM F136)- Spineplate and	Tantalum markers (ASTM F560)-
	screws	cage
		PEEK (ASTM F 2026) with Ti-6Al-
AnyPlus® Peek Cage K100516/K111354	-	4V ELI titanium alloy makers
		(ASTM F136)

APPROACH		
This Subject ABTross ALIF Expandable	Anterior Approach (ALIF)	
Cage System		
Predicate Devices	Similarities	Differences
		(Provided other approach cages)
MAGNIFY TM Spacers	Anterior Approach (ALIF)	
MAGNIFY [™] S Spacers/K142498	Anterior Approach (ALIP)	
XYPAN Expandable Lumbar Cage (XL,		 lateral approach (LLIF)
XTP) System K203531	-	 Anterior to Psoas approach
		(ATP).
Xenco Medical Lumbar Interbody	Anterior Approach (ALIF)	_
System K143158		
Vu a POD Intervertebral Body Fusion	• Anterior Approach (ALIF)	_
Device K080822/K121211/K173606		_
AnyPlus® Peek Cage K100516/K111354	 Anterior Approach (ALIF) 	Posterior Approach (PLIF)
	Anterior Approach (ALIP)	Transforaminal Approach (TILF)

STERILIZATION & METHOD		
This Subject ABTross ALIF Expandable	Sterile Cage: Gamma Radiation	
Cage System	 Intended for SINGLE USE ONLY 	
Predicate Devices	Similarities	Differences
MAGNIFY TM Spacers	Sterile Cage: Gamma Radiation	
MAGNIFY™ S Spacers K142498	 Intended for SINGLE USE ONLY 	-
XYPAN Expandable Lumbar Cage (XL,	Sterile Cage: Gamma Radiation	Non-Sterile
XTP) System K203531	 Intended for SINGLE USE ONLY 	- Non-Sterne
Xenco Medical Lumbar Interbody	Sterile Cage: Gamma Radiation	
System K143158	 Intended for SINGLE USE ONLY 	-
Vu a POD Intervertebral Body Fusion		Non-Sterile
Device K080822/K121211/K173606	-	
AnyPlus® Peek Cage K100516/K111354	-	Non-Sterile

9. SUBSTANTIAL EQUIVALENCE AND CONCLUSION

The subject ABTross ALIF Expandable Cage System have been found to be substantially equivalent to the predicate devices with respect to technical characteristics, performance, and intended use. The information provided within this premarket notification supports substantial equivalence of the subject

device to the predicate devices. The overall Data lead to the conclusion that the ABTross ALIF Expandable Cage System is substantially equivalent to the predicate devices.