

April 14, 2023

BTL Industries Inc. David Chmel VP of Operations 362 Elm Street Marlborough, Massachusetts 01752

Re: K221865

Trade/Device Name: BTL-094 Regulation Number: 21 CFR 890.5660 Regulation Name: Therapeutic Massager Regulatory Class: Class I Product Code: ISA Dated: March 20, 2023 Received: March 20, 2023

Dear David Chmel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Amber T. Ballard -S

Amber Ballard, PhD Assistant Director DHT5B: Division of Neuromodulation and Physical Medicine Devices OHT5: Office of Neurological and Physical Medicine Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

## **Indications for Use**

510(k) Number *(if known)* K221865

Device Name BTL-094

Indications for Use (Describe) BTL-094 is indicated to be used for:

- Relief of minor muscle aches and pains.
- Temporary increase in local blood circulation.
- Activation of connective tissue.

Type of Use (Select one or both, as	applicable)	

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



## 510(k) Summary

## **General Information**

Sponsor:	BTL Industries, Inc. 362 Elm Street Marlborough, MA 01752 Tel: <u>+1-866-285-1656</u> Fax: +1-888-499-2502
Applicant:	BTL Industries, Inc. 362 Elm Street Marlborough, MA 01752 Tel: <u>+1-866-285-1656</u> Fax: +1-888-499-2502
Contact Person:	David Chmel BTL Industries, Inc. <u>chmel@btlnet.com</u>
Summary Preparation Date:	April 14, 2023

## **Device Name**

Trade/Proprietary Name:	BTL-094
Primary Classification Name:	Therapeutic Massager
Common/Usual Name:	Therapeutic Massager
Classification Regulation:	21 CFR 890.5660, Class I
Classification Product Code:	ISA

## Legally Marketed Predicate Device

The BTL-094 is a state-of-the-art device substantially equivalent to the following products, which are already cleared for distribution in the USA under the following 510(k) Premarket Notification number:

- D-Actor 200 (K173692)
- DolorClast® Radial (K220538)

## **Product Description**

The BTL-094 is a non-invasive therapeutic device, which uses acoustic waves in order to stimulate a local biological response in the treated tissue.

The BTL-094 generates extracorporeally induced massage pulses to stimulate a local biological response in the treated tissue. The main source of acoustic waves is a piezo transducer inside the applicator. The energy transfer from the piezo transducer to the treated tissue is realized via coupling pad.

The main unit is equipped with a color touch screen, which considerably simplifies its operation. The onscreen information guides the user step-by-step through the entire therapy process. A secondary screen is placed on the applicator for the presentation of the actual state of main therapy parameters. The therapeutic parameters are adjustable via touch-screen, buttons and knob. The user is able to set following therapy parameters: intensity, frequency and number of shocks.

The mechanism of action and technological similarities and differences between the BTL-094 device and the predicate devices are described below in the comparison table. The differences do not raise any new types of safety or effectiveness questions.

## **Indications For Use**

BTL-094 is indicated to be used for:

- Relief of minor muscle aches and pains.
- Temporary increase in local blood circulation.
- Activation of connective tissue.



## **Non-clinical Testing**

The BTL-094 device has been thoroughly evaluated for electrical safety. The device has been found to comply with the following applicable medical device safety standards:

IEC 60601-1:2005/A1:2012 FR Recognition Number: 19-4	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2:2014	Medical electrical equipment – Part 1-2: General requirements for
FR Recognition Number: 19-8	basic safety and essential performance - Collateral standard: Electromagnetic disturbances – Requirements and tests
IEC 60601-1-6:2010/A1:2013	Medical electrical equipment - Part 1-6: General requirements for
FR Recognition Number: 5-89	basic safety and essential performance - Collateral standard: Usability
IEC 60601-2-36:2014	Medical electrical equipment - Part 2-36: Particular requirements for
FR Recognition Number: 9-119	the basic safety and essential performance of equipment for extracorporeally induced lithotripsy
IEC 62304:2006/A1:2015	Medical device software – Software lifecycle processes
FR Recognition Number: 13-79	
ISO 14971:2019	Medical devices – Application of risk management to medical
FR Recognition Number: 5-125	devices
ISO 10993-1:2018	Biological evaluation of medical devices – Part 1: Evaluation and
FR Recognition Number: 2-258	testing within a risk management process
ISO 10993-5:2009	Biological evaluation of medical devices – Part 5: Tests for in vitro
FR Recognition Number: 2-245	cytotoxicity
ISO 10993-10:2010	Biological evaluation of medical devices – Part 10: Tests for
FR Recognition Number: 2-174	irritation and skin sensitization

## **Clinical Testing**

Not aplicable

## **Technological Characteristics**

The BTL-094 is a non-invasive therapeutic device, which uses acoustic waves in order to stimulate a local biological response in the treated tissue.

The BTL-094 generates extracorporeally induced massage pulses to stimulate a local biological response in the treated tissue. The main source of acoustic waves is a piezo transducer inside the applicator. The energy transfer from the piezo transducer to the treated tissue is realized via coupling pad.

The mechanism of action and technological similarities and differences between the BTL-094 device and the predicate devices are described below in the comparison table. The differences do not raise any new types of safety or effectiveness questions.

Comparison	with	the	Predicate	Device
------------	------	-----	-----------	--------

510(k)	K221865	K173692	K220538	
number Device name	BTL-094	D-Actor 200	DolorClast® Radial	
Company name	BTL Industries Inc.	Storz Medical	EMS Electro Medical SA	
Indications for Use	Relief of minor muscle aches and pains Temporary increase in local blood circulation Activation of connective tissue	*Relief of minor muscle aches and pains *Temporary increase in local blood circulation *Activation of connective tissue	<ul> <li>Relief of minor muscle aches and pains.</li> <li>Temporary increase in local blood circulation.</li> <li>Activation of connective tissue.</li> </ul>	None
Product Code and Regulation	21 CFR 890.5660 ISA	21 CFR 890.5660 ISA	21 CFR 890.5660 ISA	None
Modes of Action	Extracorporeally induced pressure waves.	Radial pressure waves, or extracorporeal pulse activation respectively	Radial pressure waves, or extracorporeal pulse activation respectively	Not Significantly different Please see the discussion and conclusion below.
Mechanisms of Action	Vibrations generated by electroacoustic technology.	Pneumatically generated vibrations	Pneumatically generated vibrations	Not Significantly different Please see the discussion and conclusion below.
Maximum and Minimum intensity settings	5-30% alternative output power setting scale	1-5bar	1-4bar	Not Significantly different Differences are caused by alternative technology of the pressure wave generation.



Number and size of treatment applicator heads	2; L Pad - 12.2 cm <sup>2</sup> U Pad - 7.3 cm <sup>2</sup>	4; 6mmOD, 15mmOD, 20mmOD, 35mmOD	7; 5mm OD, 10 mm OD, 15 mm OD, 25 mm OD, 40 mm OD	Not Significantly different The application area of the pads in contact with the patient is within the set range of the predicate devices. All materials in contact with the patient have been assessed for biocompatibility
Maximum and minimum displacements of applicator heads	N/A	0.6 – 2.0mm	0.18 - 0.28 mm	N/A The BTL-094 applicator does not have any moving parts that come in contact with the patient.
Type of application (e.g., continuous vibration at a fixed frequency);	Continuous vibration at a fixed frequency and single pulse mode	Continuous vibration at a fixed frequency	Continuous vibration at a fixed frequency	None
Maximum and minimum vibration frequency	1-25 Hz	1-21Hz	1-25 Hz	Not Significantly different
Power Supply	100 V to 240 V AC 50 - 60 Hz 150VA	500VA	Mains power	None
Maximum penetration depth	35 mm	32.3mm	40 mm	Not Significantly different Please see the discussion and conclusion below.
Energy flow density	In the range of 0.01 - 0.11 mJ/mm2	Values of ultrasonic pulse: 5bar/0.284mJ/mm2 3bar/0.176mJ/mm2	0.29 mJ/mm2 at 4 bar 0.14 mJ/mm2 at 2.4 bar	Not Significantly different
Operating mode	Continuous	Continuous	Continuous	None



Projectile mass (g)	N/A	3	3.2	<b>Not applicable</b> The BTL-094 applicator does not contain a projectile.
Pulse repeat rate (1/s)	1-25 Hz	1-21Hz	1-25 Hz	Not Significantly different
Number of pulses (min and max)	Variable in range 0 - 9999	Variable	Variable, Max 5000/treatment	Not Significantly different
Maximum operating temperature	10 - 35 °C	10-40°C	10 - 30 °C	Not Significantly different
Type of acoustic wave generation	Electroacoustic technology with electroacoustic lens system	Pneumatic/ballistic	Pneumatic/ballistic	Not Significantly different Differences are caused by alternative technology of the pressure wave generation.
Positive peak pressure amplitude (MPa)	In the range of 3 - 10 MPa	Values of ultrasonic pulse: 5bar/18.5MPa 3bar/13.4MPa	17 MPa at 4 bar 11.24 MPa at 2.4 bar	Not Significantly different Please see the discussion and conclusion below.
Negative peak pressure amplitude (MPa)	In the range of 3 - 10 MPa	Values of ultrasonic pulse: 5bar/6.8MPa 3bar/5.0MPa	10 MPa at 4 bar 7.2 MPa at 2.4 bar	Not Significantly different Please see the discussion and conclusion below.
Derived focal acoustic pulse energy (mJ)	In the range of 0.1 – 1.4 mJ	Values of ultrasonic pulse: 5bar/6.5mJ 3bar/2.4mJ	5.9 mJ at 4 bar* 2.2 mJ at 2.4 bar* *Averaged across all applicator sizes.	Not Significantly different Please see the discussion and conclusion below.
Derived pulse- intensity integral, integrated over total temporal integration limits (mJ/mm2)	In the range of 0.01 - 0.11 mJ/mm2	Values of ultrasonic pulse: 5bar/0.284mJ/mm2 3bar/0.176mJ/mm2	0.29 mJ/mm2 at 4 bar 0.14 mJ/mm2 at 2.4 bar	Not Significantly different Please see the discussion and conclusion below.
Rise time (ns)	10 ns	Ultrasonic pulse: 2.5µs	3.2 µs	Not Significantly different



		Sonic pulse: 25µs – 2.5ms		Differences are caused by alternative technology of the energy generation.
Compressional pulse duration (µs)	2µs	Ultrasonic pulse: 5.0µs Sonic pulse: 50µs – 5.0ms	2.6 µs (1st peak) 62.7 µs (1st phase)	Not Significantly different

#### Modes of Action, Mechanisms of Action

The technology used in BTL-094 device differs from the predicate device. Nevertheless, the technology creates the same type of acoustic pressure delivered to the tissue. Similar devices that generate a pressure wave based on a different technological principle but with the same or very comparable mechanism of action have already been recognised as substantially equivalent in this device category.

#### Maximum and Minimum intensity settings

BTL-094 has an alternative type of technology from the predicate device and therefore the BTL-094 device has an alternative power setting scale - in percentage and not in bars as a comparative device has. In addition, for comparability with the predicate devices the BTL-094 device intended for the US market is limited to a maximum power level of 30% of maximum driving power.

#### Number and size of treatment applicator heads

The number of application pads are different from predicate device but the area of the subjected device pads that is in contact with the patient is comparable within a set range to the contact area of application tips of predicate devices. Additionally all materials in contact with patient have been assessed for biocompatibility and the maximal depth of pressure wave penetration is comparable to the predicate.

#### Maximum and minimum displacements of applicator heads

The predicate device is generating the pressure wave by the head displacement. Due to an alternative type of technology used for pressure wave generation the BTL-094 device does not have the ability of displacements of applicator heads/pads.

#### Maximum and minimum vibration frequency

BTL-094 device differs from one predicate device in the value of the maximum frequency, where the frequency value of the BTL-094 device can be up to 4Hz higher than the predicate device, but has available the same frequency range as its second predicate device.

#### **Driving Power**

BTL-094 device has an alternative type of technology from the predicate device and therefore the driving power parameter is not applicable. In addition, for comparability with the predicate devices the BTL-094 device intended for the US market is limited to a maximum power level of 30% of maximum driving power.

#### Energy flow density

BTL-094 device uses alternative technology and is able to better pulse targeting, a maximum driving power of the device has been limited at a maximum power level of 30%. The EFD of our device differs from the predicates but it is lower or within setting range of predicat devices.



#### Operating mode

The BTL-094 device has two modes of operation: single therapy mode or sequence mode. The two modes differ based on repetition of therapy sections and settings that can be done during the running therapy. Nevertheless, both the modes have a character of continual therapy.

#### Maximum penetration depth

The BTL-094 device is equipped with two different coupling pad sizes, which allow it to reach up to 35 mm depth. This penetration depth is comparable within setting range of predicate devices.

#### Pulse repeat rate (1/s)

BTL-094 device differs from one predicate device in the value of the maximum frequency, where the frequency value of the BTL-094 device can be up to 4Hz higher than the predicate device, but has available the same frequency range as its second predicate device.

#### Number of pulses (min and max)

BTL-094 device allows to set a variable number of pulses according to the user's needs. This parameter is similar as for the predicate device.

#### Type of acoustic wave generation

The subject device provides the therapy by the application of extracorporeally induced acoustic pressure waves that are generated by means of piezoelectric principle. The device is designed with planar piezoelement and acoustic lens guiding the waves. The generated wave is coupled to biological tissue via the coupling gel pads of the applicator.

This technology creates the same type of acoustic pressure wave with comparable parameters therefore it reach the same mechanism of action as the predicate devices.

#### Positive peak pressure amplitude (MPa)

A maximum driving power of the BTL-094 has been limited at a maximum power level of 30%. The device with this limitation in maximal output power fits with output parameters within or bellow setting range of predicat devices and is not significantly different from the predicate devices.

#### Negative peak pressure amplitude (MPa)

A maximum driving power of the BTL-094 has been limited at a maximum power level of 30%. The device with this limitation in maximal output power fits with output parameters within or bellow setting range of predicat devices and is not significantly different from the predicate devices.

#### Derived focal acoustic pulse energy (mJ)

A maximum driving power of the BTL-094 has been limited at a maximum power level of 30%. The device with this limitation in maximal output power fits with output parameters within or bellow setting range of predicat devices and is not significantly different from the predicate devices.

#### Derived pulse- intensity integral, integrated over total temporal integration limits (mJ/mm2)

A maximum driving power of the BTL-094 has been limited at a maximum power level of 30%. The device with this limitation in maximal output power fits with output parameters within or bellow setting range of predicat devices and is not significantly different from the predicate devices.

#### Rise time (ns)

BTL-094 has an alternative type of technology from the predicate device and therefore the BTL-094 device is able to reach a faster rise time and shorter compressional pulse duration.



#### Compressional pulse duration (µs)

BTL-094 has an alternative type of technology from the predicate device and therefore the BTL-094 device is able to reach a faster rise time and shorter compressional pulse duration.

#### **Substantial Equivalence**

The BTL-094 is substantially equivalent to the predicate devices D-Actor (K173692) and DolorClast® Radial (K220538). The subject device is safe and effective for its intended use.

Any differences between the predicate devices and BTL-094 have no significant influence on the safety or effectiveness of the subject device. Therefore, the BTL-094 is substantially equivalent to the predicate devices.

### Conclusion

Based upon the intended use and known technical information provided in this pre-market notification, the BTL-094 device has been shown to be substantially equivalent to the currently marketed predicate device.