

March 9, 2023

Cerenovus, Inc. Ariell Joiner, Ph.D. Manager, Regulatory Affairs 6303 Blue Lagoon Drive, Suite 315 Miami, Florida 33126

Re: K221934

Trade/Device Name: CEREGLIDE 71 Intermediate Catheter; Cerenovus Aspiration Tubing Set

Regulation Number: 21 CFR 870.1250 Regulation Name: Percutaneous Catheter

Regulatory Class: Class II Product Code: NRY Dated: February 2, 2023 Received: February 3, 2023

Dear Ariell Joiner:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Naira Muradyan -S

Naira Muradyan, Ph.D.
Assistant Director
DHT5A: Division of Neurosurgical,
Neurointerventional
and Neurodiagnostic Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known)
K221934
K2217J4
Device Name
CEREGLIDE™ 71 Intermediate Catheter; Cerenovus Aspiration Tubing Set
Indications for Use (Describe)
CEREGLIDETM 71 Intermediate Catheter
The CEREGLIDE™ 71 Intermediate Catheter, with the Cerenovus Aspiration Tubing Set and a compatible aspiration
pump, is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large
vessel occlusive disease (within the internal carotid, middle cerebral -M1 and M2 segments, basilar, and vertebral arteries
within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or
who fail IV t-PA therapy are candidates for treatment.
Cerenovus Aspiration Tubing Set
The Cerenovus Aspiration Tubing Set is intended to connect the Cerenovus Large Bore Catheter to the canister of the
Nouvag Vacuson 60 Aspiration Pump (or equivalent vacuum pump) and to allow the user to control the fluid flow. The
Cerenovus Aspiration Tubing Set is also intended to connect the CEREGLIDE TM 71 Intermediate Catheter to the canister
of a compatible aspiration pump and to allow the user to control the fluid flow.
of a companion aspiration paint and to anow the aser to contact the flata flow.
Type of Use (Select one or both, as applicable)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary K221934

I. Submitter

Cerenovus, Inc.

6303 Blue Lagoon Drive

Suites 215 & 315 Miami, FL 33126 USA

Contact Person: Ariell Joiner

Tel: (908) 249-0182

Email: ajoiner@its.jnj.com

II. Date Prepared

March 8, 2023

III. Device Information

Table 1. Device Information			
Device Proprietary Name	CEREGLIDE™ 71 Intermediate Catheter; Cerenovus Aspiration Tubing Set		
Common or Usual name	Catheter, Thrombus Retriever		
Classification Name	21 CFR 870.1250 – Catheter, Percutaneous		
Regulatory Classification	П		
FDA Product Code	NRY		

IV. Predicate Device Information

The primary predicate device is listed below in **Table 2** and the reference device is listed in **Table 3**.

Table 2. Primary Predicate Device			
510(k) Number	Date Cleared	Name	Manufacturer
K193380	July 20, 2020	CERENOVUS Large Bore Catheter; CERENOVUS Aspiration Tubing Set	Medos International SARL*
*510(k) was previously held by Codman & Shurtleff, Inc. The 510(k) was transferred to			
Medos International SARL, which is the current holder of the 510(k).			

Table 3. Reference Device			
510(k) Number	Date Cleared	Name	Manufacturer
K173200	June 11, 2018	SOFIA Plus Aspiration Catheter	MicroVention, Inc.

V. Device Description

The CEREGLIDETM 71 Intermediate Catheter is a variable stiffness, single lumen catheter designed to be introduced over a steerable guide wire or microcatheter into the neuro vasculature. The catheter shaft is composed of a stainless steel variable pitch braid with a PTFE inner liner to facilitate movement of guide wires and other devices. The exterior of the catheter shaft is covered with polymer materials, which encapsulate the stainless steel braid construction. The catheter has a stiff proximal shaft which transitions into the flexible distal shaft to facilitate the advancement of the catheter in the anatomy. The distal end of the catheter has a radiopaque marker band to facilitate fluoroscopic visualization and has a hydrophilic coating to provide lubricity to reduce friction during navigation in the vasculature. The proximal end of the catheter has a luer fitting located on the end of the catheter hub which can be used to attach accessories for flushing and aspiration. An ID band is placed at the distal end of the hub over a strain relief.

The catheter is packaged with a Tuohy Borst rotating hemostasis valve (RHV) with a side port and two slit introducers as accessories. The RHV with side port is used for flushing, insertion of catheters, and connection to an external aspiration system. The slit introducers are designed to introduce the catheter into the base catheter and protect the distal tip of the catheter during insertion into the RHV of the base catheter.

The CEREGLIDE™ 71 Intermediate Catheter can be connected to a compatible aspiration pump using the Cerenovus Aspiration Tubing Set.

VI. Indications for Use

The CEREGLIDETM 71 Intermediate Catheter, with the Cerenovus Aspiration Tubing Set and a compatible aspiration pump, is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral – M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

The Cerenovus Aspiration Tubing Set is intended to connect the Cerenovus Large Bore Catheter to the canister of the Nouvag Vacuson 60 Aspiration Pump (or equivalent vacuum pump) and to allow the user to control the fluid flow. The Cerenovus Aspiration Tubing Set is also intended to connect the CEREGLIDETM 71 Intermediate Catheter to the canister of a compatible aspiration pump and to allow the user to control the fluid flow.

VII. Predicate Comparison

A comparison of the similarities and differences of product features between the CEREGLIDETM 71 Intermediate Catheter, Cerenovus Aspiration Tubing Set and the predicate device is presented in **Table 4**.

Table 4. Subject and Predicate Device Comparison Summary			
Description	Subject Device: CEREGLIDE™ 71 Intermediate Catheter; Cerenovus Aspiration Tubing Set	Predicate Device: CERENOVUS Large Bore Catheter (K193380)	
Product Code	NRY	Same	
Regulatory Name	Catheter, Percutaneous	Same	
Classification	Class II – 21 CFR 870.1250	Same	
Basic Design	Variable stiffness single lumen catheter	Same	
Indications For Use	The CEREGLIDE TM 71 Intermediate Catheter, with the Cerenovus Aspiration Tubing Set and a compatible aspiration pump, is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral - M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. The Cerenovus Aspiration Tubing Set is intended to connect the Cerenovus Large Bore Catheter to the canister of the Nouvag Vacuson 60 Aspiration Pump (or equivalent vacuum pump) and to allow the user to control the fluid flow. The Cerenovus Aspiration Tubing Set is also intended to connect the CEREGLIDE TM 71 Intermediate Catheter to the canister of a compatible aspiration pump and to allow the user to control the fluid flow.	The CERENOVUS Large Bore Catheter, with the CERENOVUS Aspiration Tubing Set and NOUVAG Vacuson 60 aspiration pump (or equivalent aspiration pump), is indicated for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral - M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t- PA therapy are candidates for	
<u>Dimensions:</u>			
Length	115 - 137cm	125 - 135cm	
Inner Diameter (ID)	0.071"	Same	
Distal Outer Diameter (OD)	0.082" (0.0825" max)	0.081"	
Proximal OD	0.0823 max) 0.0837" max	0.0825"	
Catheter Coating	Hydrophilic	Same	
Coating Length	55 cm	30 cm	
Tip Configuration	Non-shapeable tip	Same	
Materials:			
Marker Band	Metal Platinum (90%) / Iridium (10%)	Same	
Braid	Stainless Steel	Same	
Liner	PTFE Liner	Same	
Hub	Polyamide	Same	
Strain Relief		Same	
Outer Jacket	Pebax and Urethane	Pebax, Urethane, Nylon	

Table 4. Subject and Predicate Device Comparison Summary, continued					
Description	Subject Device: CEREGLIDE™ 71 Intermediate Catheter; Cerenovus Aspiration Tubing Set	Predicate Device: CERENOVUS Large Bore Catheter (K193380)			
Accessories Included:					
Hemostasis Valve	Included: Tuohy Borst Hemostasis Valve with Side Port Extension Tubing	Same			
Introducer Sheath	Introducer (2)	Peel-Away Sheath Introducer (2)			
Sterilization Method	Ethylene Oxide	Same			
Sterility Assurance Level (SAL)	10^{-6}	Same			
Packaging	Polyethylene Hoop and Mounting Card, Tyvek® Pouch, Carton	Same			
Shelf Life	1 year	1 year			
Required Additional Accessories	Cerenovus Aspiration Tubing Set (K193380) Compatible aspiration pump	Same			
Aspiration Pump Requi					
Minimum Aspiration Pressure	-20 inHg (-68 kPa)	Same			
Maximum Aspiration Pressure	-28 inHg (-95 kPa)	-29 inHg (-98kPa)			
Flowrate (Air)	0 to 60LPM	Same			
Aspiration Tubing Requ	Aspiration Tubing Requirements:				
Tubing ID	0.110 in minimum	Same			
Tubing Length	112 in	Same			
Flow Control Mechanism	Flow Control Switch	Same			

VIII. Non-Clinical Performance Data

<u>Performance Testing – Bench</u>

Appropriate testing was identified based on design, risk analyses and the intended use of the CEREGLIDETM 71 Intermediate Catheter to demonstrate that it is substantially equivalent to the legally marketed predicate device. The following performance data are being provided in support of the substantial equivalence determination. All testing was conducted using sampling methods as required by internal procedure. The bench testing included the following tests:

Table 5. Performance Testing Summary			
Test	Test Summary	Result	
	Design Verification		
Visual Inspection	Confirm that the Catheter meets the visual requirement described in ISO 10555-1 Section 4.4.	PASS: Samples met the established acceptance criteria	
Catheter ID	Verify that the Catheter internal diameters meet the requirements.	PASS: Samples met the established acceptance criteria	
Catheter OD	Verify that the Catheter outer diameters meet the requirements.	PASS: Samples met the established acceptance criteria	
Catheter Working Length	Confirm the working length of the catheter as defined in ISO 10555-1 Section 3.6.	PASS: Samples met the established acceptance criteria	
Catheter Tip Length	Verify the catheter tip length of the Catheter.	PASS: Samples met the established acceptance criteria	
System Air Leakage	Verify that there is no air leak into the hub subassembly.	PASS: Samples met the established acceptance criteria	

VIII. Non-Clinical Performance Data, continued

Table 5. Performance Testing Summary, continued				
Test	Test Summary	Result		
Design Verification				
System Liquid Leakage	Verify that the catheter joint strength meets the freedom from leakage (liquid during pressurization) requirements of ISO 10555-1:2013, section 4.7.	PASS: Samples met the established acceptance criteria		
Delamination of PTFE Liner	Verify that the PTFE has appropriately adhered to the inner lumen of the Catheter with braid reinforcement.	PASS: Samples met the established acceptance criteria		
Kink (Distal & Proximal)	Confirm that the Catheter meets the requirement for the catheter to remain stable and not kink during use.	PASS: Samples met the established acceptance criteria		
Tip Movement	Confirm that the Catheter meets the tip column stiffness requirement.	PASS: Samples met the established acceptance criteria		
Tip Linear Stiffness	Test the tip flexibility of the Catheter, relative to other devices of similar design.	PASS: Samples met the established acceptance criteria		
Coating Lubricity and Durability	Verify the lubriciousness and durability of the Catheter's hydrophilic coating.	PASS: Samples met the established acceptance criteria		
Coating Length	Verify that the Catheter hydrophilic coating length meets the design requirements.	PASS: Samples met the established acceptance criteria		
Peak Tensile Strength	Verify that the Catheter joint strength meets the requirements of Section 4.5 of ISO 10555-1.	PASS: Samples met the established acceptance criteria		
Introducer ID	Verify that the introducer internal diameters meet the requirements.	PASS: Samples met the established acceptance criteria		
Particulate Count	Verify that the coating integrity of the Catheter's outer surface meets the requirements for content of Particle Matter in alignment with USP<788> counting methods and compared to the reference predicate device.	PASS: Samples met the established acceptance criteria		
Burst Pressure	Confirm the maximum hydrostatic pressure the Catheter can withstand using a Crescent Hydraulic Burst-leak Tester.	PASS: Samples met the established acceptance criteria		
Introducer Working Length	Confirm the working length of the introducer.	PASS: Samples met the established acceptance criteria		
Introducer Separation Force	Confirm the force required to separate the introducer.	PASS: Samples met the established acceptance criteria		
Torque Test	To determine the number of revolutions to failure of the Catheter in simulated anatomy.	PASS: Samples exceeded comparator devices in revolutions to failure		
Design Validation				
In Vitro Usability Studies	The in-vitro studies were conducted to demonstrate thrombus retrieval patency/durability, (ancillary) device compatibility, and accessory durability. Studies also included simulated use evaluation of user requirements related to: trackability and tip stability during thrombus removal.	PASS: Samples met the established acceptance criteria		

VIII. Non-Clinical Performance Data, continued

Performance Testing – Animal

No animal studies were required as appropriate verification and validation of the device design were achieved based on the similarities of the proposed device and the predicate device, and from results of bench testing.

Performance Testing – Clinical

Clinical studies were not required as appropriate verification and validation of the device design were achieved based on the similarities of the proposed device and the predicate device, and from results of bench testing.

Sterilization

The CEREGLIDETM 71 Intermediate Catheter, as packaged with included accessories, is sterilized using a validated 100% Ethylene Oxide sterilization process to ensure sterility assurance level (SAL) of 10⁻⁶ in accordance with ISO 11135. The CEREGLIDETM 71 Intermediate Catheter and all accessories meet EO residuals per EN ISO 10993-7 for a limited contact delivery system – externally communicating. The CEREGLIDETM 71 Intermediate Catheter and all accessories are for single use only.

Shelf-Life

The CEREGLIDETM 71 Intermediate Catheter will have a shelf life of one year based on the successful completion of stability testing. Shelf life testing was performed using standard test methods and acceptance criteria. Prior to aging, all samples were exposed to standard transportation conditioning. Results of testing on the subject device all met established acceptance criteria.

Biocompatibility Testing

A biological safety evaluation was conducted on the CEREGLIDE™ 71 Intermediate Catheter in accordance with the FDA guidance, *Use of International Standard ISO 10993-1*, "Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process", (Issued June 16, 2016) based on the changes to the predicate device, CERENOVUS Large Bore Catheter (K193380). Biocompatibility testing was conducted in accordance with International Standard ISO 10993-1, "Biological Evaluation of Medical Devices – Part 1: Evaluation of Testing Within a Risk Management Process" (2018).

The following testing was performed:

Table 6. Biocompatibility Test Summary for CEREGLIDE™ 71 Intermediate				
	Catheter and Introducer			
Test	Test Summary	Results		
Chemical	Determination of extractable species from the test article was	PASS		
Characterization	performed in purified water, isopropyl alcohol, and hexane.			
Cytotoxicity	Cytotoxicity Study Using the ISO Elution Method (compliant with ISO 10993-5:2009/(R)2014 and ISO 10993-12:2012).	PASS		
	1X MEM extract (6 cm²/ mL at 37°C for 72 hours) was evaluated for potential cytotoxic effects using monolayers of L-929 mouse fibroblast cells.			

VIII. Non-Clinical Performance Data, continued

Table 6. Biocompatibility Test Summary for CEREGLIDE™ 71 Intermediate			
Catheter and Introducer, continued			
Test	Test Summary	Results	
Sensitization –	Evaluation of the allergenic potential or sensitizing capacity of a	PASS	
Guinea Pig	test article by screening of contact allergens in guinea pigs and		
Maximization	extrapolating the results to humans.		
Irritation – Rabbit	Determination of any chemicals that may leach or be extracted	PASS	
Intracutaneous	from the test article capable of causing local irritation in the		
Reactivity	dermal tissues of rabbits.		
Acute Systemic	Screening of test article extracts for potential toxic effects as a	PASS	
Toxicity in Mice	result of a single-dose systemic injection in mice.		
Pyrogenicity –	Determination if a saline extract of the test article causes a	PASS	
Materials Mediated	febrile response in rabbits.		
Rabbit Pyrogen	reothe response in radous.		
ASTM Hemolysis	Evaluation of the hemolytic potential of the test articles and test	PASS	
Study: Direct and	article extract according to ASTM method (F756-17).		
Extraction Methods			
SC5b-9	Measure of complement activation in Normal Human Serum	PASS	
Complement	(NHS) indicates whether a test article is capable of generating the		
Activation Assay	activation fragment, SC5b-9, which in turn contributes to the		
	inflammatory immune response in humans.		
ASTM Heparinized	Determination if medical materials exposed to human whole	PASS	
Platelet and	blood would adversely affect the platelet and leukocyte ratios in		
Leukocyte Count	whole blood.		
ASTM Partial	Screening for detection of coagulation abnormalities in the	PASS	
Thromboplastin	intrinsic coagulation pathway. The PTT indirectly measures the		
Time (PTT)	formation of thrombin by its action on fibrinogen, forming the		
	fibrin clot through a method compliant with ISO 10993-4.		
In Vivo	Evaluate the thrombogenic potential of a blood contacting	PASS	
Thromboresistance	medical device in comparison to a predicate device.		

IX. Conclusion

Based upon the intended use, design, materials, function, and side-by-side *in-vitro* and *in-vivo* testing, it is concluded that the subject device, CEREGLIDETM 71 Intermediate Catheter, Cerenovus Aspiration Tubing Set is substantially equivalent to the predicate device, CERENOVUS Large Bore Catheter, CERENOVUS Aspiration Tubing Set (K193380). The differences in materials and design do not raise new questions regarding the safety and effectiveness of the device. The device, as designed, manufactured, packaged and sterilized, is substantially equivalent to the primary predicate and reference devices currently marketed under the Federal Food, Drug and Cosmetic Act.