

April 27, 2023

Guangdong Kingfa Sci. & Tech. Co., LTD. Xiaoge Yu Manager No.28, Delong Ave., Shijiao Town, Qingcheng District Qingyuan, Guangdong 511545 China

Re: K221983

Trade/Device Name: Biodegradable Powder Free Nitrile Examination Glove, Blue Color,

Biodegradable Powder Free Nitrile Examination Glove, Blue Violet Color, Biodegradable Powder Free Nitrile Examination Glove, Green Color

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA Dated: November 8, 2022 Received: April 24, 2023

Dear Xiaoge Yu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian -S

Bifeng Qian., M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number	(if known)
K221983	

Device Name

Biodegradable Powder Free Nitrile Examination Glove, Blue Color

Biodegradable Powder Free Nitrile Examination Glove, Blue Violet Color

Biodegradable Powder Free Nitrile Examination Glove, Green Color

Indications for Use (Describe)

The biodegradable nitrile examination glove is intended to be worn on the hands of examiner's to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.

Type of Use (Select one or both, as applicable)	
Type of dee (defect one of sour, de approache)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K221983-510(k) Summary

I. Submitter

GUANGDONG KINGFA SCI. & TECH.CO., LTD.

No. 28 Delong Avenue, Shijiao Town, Qingcheng District, Qingyuan City, Guangdong Province 511545, China

Contact person: Xiaoge Yu

Position: Manager Tel.: +86-13570952157

E-mail: yuxiaoge@kingfa.com.cn

Revision date: September 7, 2022

II. Proposed Device

Device Trade Name Biodegradable Powder Free Nitrile Examination Glove, Blue

Color

Biodegradable Powder Free Nitrile Examination Glove, Blue

Violet Color

Biodegradable Powder Free Nitrile Examination Glove, Green

Color

Common name: Polymer Patient Examination Glove

Regulation Number: 21 CFR 880.6250

Regulatory Class: Class I Product code: LZA

Review Panel General Hospital

III. Predicate Device

510(k) Number: K162312

Trade name: Biodegradable Powder Free Nitrile Exam Gloves With Aloe Vera,

Green Color

Common name: Patient Examination Gloves

Classification: Class I Product Code: LZA

Manufacturer Shen Wei USA Inc.

IV. Device description

The propose devices is powder free nitrile patient examination gloves, provided as

non-sterile and disposable device. The gloves are offered in six sizes, extra-small (6.5"), small (7"), medium (8"), large (8.5"), extra-large (9"), extra extra-large (9.5"). Three colors are available for all size, includes blue, blue violet and green.

V. Indication for use

The biodegradable nitrile examination glove is intended to be worn on the hands of examiner's to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.

VI. Comparison of technological characteristics with the predicate devices

Table 1 Comparison of Nitrile Examination Gloves

Item Proposed device (K221983)		Predicate device	Discussion	
	(N221903)	(K162312)		
Product name	Biodegradable	Biodegradable	-	
	Powder Free Nitrile	Powder Free		
	Examination Glove,	Nitrile Exam		
	Blue Color	Gloves With Aloe		
	Biodegradable	Vera,		
	Powder Free Nitrile	Green Color		
	Examination Glove,			
	Blue Violet Color			
	Biodegradable			
	Powder Free Nitrile			
	Examination Glove,			
	Green Color			
Product Code	LZA	LZA	Same	
Regulation No.	21 CFR 880.6250	21 CFR	Same	
		880.6250		
Classification	Class I	Class I	Same	
Powder free	Yes	Yes	Same	
Indication for	The biodegradable	The biodegradable	Same	
use	nitrile examination	powder free nitrile		
	glove is intended to be	exam gloves with		
worn on the hands of		aloe vera, green color is a		
	examiner's to prevent			
	contamination between	disposable device		
	patient and examiner.	intended		
	This is a single-use,	for medical		

	powder-free, non-sterile device.	purpose that is worn on the examiner's hand to prevent contamination between patient and examiner. This device is single use only.	
Main Material	Nitrile rubber	Nitrile rubber	Same
Color	Blue\Blue violet\Green	Green	Similar
Size	X-Small (6.5), Small(7), Medium (8), Large(8.5),X-large(9), XX-large(9.5),	Ambidextrous in different sizes per ASTM D6319-10 Dimensions requirement	Similar Only the different standard version. This requirement given in the standard is the same.
Palm width	X- Small(70±10mm) Small (80±10mm) Medium(95±10mm) Large(110±10mm) X-large(120±10mm) XX-large(≧120mm)	Ambidextrous in different sizes per ASTM D6319-10 dimensions requirement	Similar Only the different standard version. This requirement given in the standard is the same.
Length	XS(220mm min) S (220mm min) M (230mm min) L (230mm min) XL (230mm min) XXL (230mm min) XXL (230mm min)	Ambidextrous in different sizes per ASTM D6319-10 Dimensions requirement	Similar Only the different standard version. This requirement given in the standard is the same. Similar

Thickness	Finger: 0.08mm min	different sizes per ASTM D6319-10 dimensions requirement	Only the different standard version. This requirement given in the standard is the same.
Freedom from holes	Meets requirements of the ASTM D6319-19	Meets requirements of the ASTM D6319- 10	Similar Only the different standard version. This requirement given in the standard is the same.
Physical Properties (before aging)	Meets requirements of the ASTM D6319-19	Meets requirements of the ASTM D6319- 10	Similar Only the different standard version. The requirements of physical properties given in the standard are the same.
Physical Properties (after aging)	Meets requirements of the ASTM D6319-19	Meets requirements of the ASTM D6319- 10	Similar Only the different standard version. The requirements of physical properties given in the standard are the same.
Powder residual	≤2.0 mg/gloves	≤2.0 mg/gloves	Same
Sterility	Non-sterile	Non-sterile	Same
Biodegradation Properties	Biodegradable	Biodegradable	Same
For single use	Yes	Yes	Same
Type of use	Over the counter use	Over the counter use	Same
Biocompatibility	Skin Sensitization-ISO 10993-10 Irritation-ISO 10993-10	Skin Sensitization-ISO10993-10 Irritation-ISO 10993-10	Similar

Acute Systemic Toxicity-ISO 10993-11	

As above comparison, the difference in the reference standard version of the subject and predicate device does not raise additional questions for safety and effectiveness of the device. The biocompatibility test and performance test of the subject devices have been performed on the final finished device. The test results shows the predicate device meet the specifications of ASTM D 6319. So we consider this as the proposed device is similar to the predicate device.

VII. Non-Clinical Testing

Non clinical tests were conducted in accordance with following standards to verify that the proposed device met all design specifications.

- ASTM D6319-19, Standard Specification for Nitrile Examination Gloves for Medical Application
- ASTM D3767-03(2020), Practice for rubber-Measurement of Dimensions
- ASTM D5151-19, Standard Test Method for Detection of Holes in Medical Gloves
- ASTM D6124-06(2017), Standard Test Method for Residual Powder on Medical Gloves
- ASTM D573-04(2019), Standard Test Method for Rubber—Deterioration in an Air Oven
- ASTM D412-16, Standard Test Methods for Vulcanized Rubber and Thermoplastic Elastomers—Tension
- ISO 10993-10: 2010 Biological Evaluation Of Medical Devices Part 10: Tests For Irritation And Skin Sensitization.
- ISO 10993-11:2017, Biological evaluation of medical devices Part11:Tests for Systemic Toxicity

Table 2 Summary of Non-Clinical Performance Testing

Test Method	Purpose	Acceptance Criteria	Results
ASTM D6319	Physical	Extra-Small:	Pass
	Dimensions Test	Length: ≥ 220mm Width: 70±10 mm;	
		Small:	
		Length: ≥ 220mm Width: 80±10mm;	
		Medium:	
		Length: ≥ 230mm Width: 95±10mm	
		Large:	
		Length: ≥ 230mm Width: 110±10mm	
		Extra- Large:	
		Length: ≥ 230mm Width: 120±10mm	
		Extra- Extra- Large:	

		Length: ≥ 230mm Width: ≥120mm			
		Thickness (mm): Finger:≥0.05			Pass
		Palm: ≥0.08			
	Physical properties		Tensile Strength	≥14MPa	Pass
	properties		Ultimate Elongation	≥500%	
		After Aging	Tensile Strength	≥14MPa	Pass
			Ultimate Elongation	≥500%	
ASTM D5151	Freedom from pinholes	Meet the requirements of ASTM D5151Test for AQL 2.5			Pass
ASTM D6124	Powder Residue	Meet the requirements of ASTM D6124< 2.0mg			Pass
ISO 10993-10	To determine if the finished device material is an irritant	Non-irritating			Under the conditions of the study not an irritant/ Pass
ISO 10993-10	To determine if the finished device material is a sensitizer	Non- sensitizing			Under conditions of the study, not a sensitizer. / Pass
ISO 10993-11	To determine if the finished device material extracts pose a systemic toxicity concern	Non-acute systemic toxicity		Under conditions of the study, did not show acute systemic toxicity in vivo / Pass	

VIII. Clinical Testing

No clinical study is included in this submission.

IX. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(k) submission K221983, Biodegradable Powder Free Nitrile Examination Gloves, is as safe, as effective, and performs as well as or better than the legally marketed predicate device, Shen Wei Biodegradable Powder Free Nitrile Examination Gloves, cleared under 510(k) K162312.