



February 27, 2023

Saneso, Inc.  
% Parul Chansoria  
Founder and CEO  
Elexes Medical Consulting, LLC  
30 N Gould St Ste R  
Sheridan, WY 82801

Re: K222173  
Trade/Device Name: Saneso Lens Wash Bottle Assembly (SAN-LWB-A)  
Regulation Number: 21 CFR 876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: Class II  
Product Code: OCX  
Dated: January 24, 2023  
Received: January 24, 2023

Dear Parul Chansoria:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any

Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Sivakami Venkatachalam -S

*for*

Shanil P. Haugen, Ph.D.

Assistant Director

DHT3A: Division of Renal, Gastrointestinal,  
Obesity and Transplant Devices

OHT3: Office of GastroRenal, ObGyn,

General Hospital and Urology Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K222173

Device Name  
Saneso Lens Wash Bottle Assembly (SAN-LWB-A)

Indications for Use (Describe)

The Saneso Lens Wash Bottle Assembly (SAN-LWB-A) is intended to be used with FDA-cleared Saneso endoscopes.

It facilitates the following:

1. Sterile water for lens wash function
2. Air / CO2 for luminal insufflation function

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



This 510(k) Summary has been created per the requirements of the Safe Medical Device Act (SMDA) of 1990, and the content is provided in conformance with 21 CFR Part 807.92.

**5.1. Submitter**

Saneso, Inc.  
One Oxford Center, 301 Grant Street Suite 4300  
Pittsburgh, PA 15219

**Contact Person:**

Parul Chansoria, MS, RAC, CQA  
CEO & Founder, Elexes Medical Consulting  
Telephone: +408-475-8091  
Email: [parul@elexes.com](mailto:parul@elexes.com)  
Summary prepared: 07/15/2022

**5.2. Device Information**

Common/Usual Name: Endoscope, Water Bottle with cap and tubing  
Trade Name: Saneso Lens Wash Bottle Assembly (SAN-LWB-A)  
Regulation Name: Endoscope and accessories  
Regulatory Class: Class II  
Classification Panel: Gastroenterology/Urology  
Product Code: OCX  
Regulation Number: 21 CFR Part 876.1500

**5.3. Predicate Device Information**

Saneso Lens Wash Bottle Assembly is substantially equivalent to the following cleared device:

| Company     | Predicate Device Priority | Product           | 510(k) Number |
|-------------|---------------------------|-------------------|---------------|
| Saneso Inc. | Primary                   | Saneso - A series | K210052       |

**5.4. Device Description**

The Saneso Lens Wash Bottle Assembly is intended to be used with Saneso endoscopes.

It facilitates the following:

1. Sterile water for lens wash function
2. Air / CO2 for luminal insufflation function

Saneso lens wash & insufflation system is divided into proximal and distal sections:

Proximal Section: It comprises the Lens Wash Bottle Assembly and terminates at its



connection with the air/water port on the endoscope. It consists of the lens wash bottle, lens wash bottle cap, airline, water line, CO2 line, and air/water connector. This section is the non-patient contacting portion of the system and is covered in this submission.

Distal Section: It starts at the air/water port and terminates at the air/water nozzles on the distal end of the endoscope. It consists of the air/water port, air/water valve, air/water channel, and the lens wash nozzles, which are all located in the endoscope. This section is the patient contacting portion of the system and is part of the Saneso endoscopes, which is FDA 510k approved (K210052).

## **5.5. Intended Use**

### **5.5.1. Subject Device**

The Saneso Lens Wash Bottle Assembly (SAN-LWB-A) is intended to be used with FDA-cleared Saneso endoscopes.

It facilitates the following:

- Sterile water for lens wash function
- Air / CO2 for luminal insufflation function

### **5.5.2. Primary Predicate Device**

#### **Saneso Colonoscope 360-A**

The Saneso Colonoscope 360-A is intended for diagnostic visualization of the lower gastrointestinal tract (including the rectum, colon, and cecum) and Saneso Colonoscope 360-A is not indicated for ileoscopy procedures. The system also provides access to therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Colonoscope 360-A, Processor- A, and other ancillary equipment.

#### **Saneso Gastroscope 360-A**

The Saneso Gastroscope 360-A is intended for diagnostic visualization of the upper gastrointestinal tract (including the esophagus, stomach and duodenum). The system also provides access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Gastroscope 360-A, Processor-A and other ancillary equipment.

#### **Saneso Single Camera Colonoscope-A**

The Saneso Single Camera Colonoscope is intended for diagnostic visualization of the lower gastrointestinal tract (including the rectum, colon and ileocecal valve). The system also provides access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Single Camera Colonoscope, Processor-A and other ancillary equipment.



**Saneso Gastroscope 360- A**

The Saneso Single Camera Gastroscope is intended for diagnostic visualization of the upper gastrointestinal tract (including the esophagus, stomach, and duodenum). The system also provides access to therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Single Camera Gastroscope, Processor-A, and other ancillary equipment.

**Note:** The scope of the submission is not the endoscope (Colonoscope and the gastroscope) but the Lens wash bottle alone. Lens wash bottle is intended to be used only as an accessory to the endoscope system (SAN-A series) which was originally cleared by the FDA with K number K210052. There is no modification in the system, but there is a modification in the accessory where it is changing from single-use to multi-use, and it is supported by the reprocessing validation results.

**5.6. Technological Characteristics**

| <b>Table 1: Substantial Equivalence Table for Saneso Lens Wash Bottle Assembly</b> |  |   |   |
|--|--|---|---|
| <b>Parameter</b>   | <b>Subject Device</b>  | <b>Primary Predicate Device</b>   | <b>Equivalence</b>  |
| <b>510(k) Number</b>   | -  | K210052   | NA  |
| <b>Product code</b>  | OCX  | FDS, FDF  | Equivalent  |
| <b>Regulation No.</b>  | 21 CFR Part 876.1500   | 21 CFR Part 876.1500  | Equivalent  |
| <b>Classification:</b>   | Class II   | Class II  | Equivalent  |
| <b>Regulation Name</b>   | Endoscope and Accessories  | Endoscope and Accessories   | Equivalent  |
| <b>Prescription Use/OTC</b>  | Prescription Use   | Prescription Use  | Equivalent  |
| <b>Intended Use</b>  | The Saneso Lens Wash Bottle Assembly is intended to be used with the FDA cleared Saneso endoscopes. It facilitates the following:<br>a. Sterile water for lens wash function<br>b. Air / CO2 for | <b><u>Saneso Colonoscope 360-A</u></b><br><br>The Saneso Colonoscope 360-A is intended for diagnostic visualization of the lower gastrointestinal tract (including the rectum, colon, and cecum) and Saneso | Equivalent<br><br><b>Note:</b> The scope of the submission is not the endoscope (Colonoscope and the gastroscope) but the Lens wash bottle alone. |

|  |                                      |   |   |
|--|--------------------------------------|---|---|
|  | <p>luminal insufflation function</p> | <p>Colonoscope 360-A is not indicated for ileoscopy procedures. The system also provides access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Colonoscope 360-A, Processor- A and other ancillary equipment.</p> <p><b><u>Saneso Gastroscope 360-A</u></b></p> <p>The Saneso Gastroscope 360-A is intended for diagnostic visualization of the upper gastrointestinal tract (including the esophagus, stomach and duodenum). The system also provides access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Gastroscope 360-A, Processor-A and other ancillary equipment.</p> <p><b><u>Saneso Single Camera Colonoscope-A</u></b></p> <p>The Saneso Single Camera Colonoscope is intended for diagnostic visualization of the lower gastrointestinal tract (including the rectum, colon and ileocecal valve). The system also provides</p> | <p>Lens wash bottle is intended to be used only as an accessory to the endoscope system (SAN-A series) which was originally cleared by the FDA with K number K210052. There is no modification in the system, but there is a modification in the accessory where it is changing from single-use to multi-use, and it is supported by the reprocessing validation results.</p> |
|--|--------------------------------------|---|---|



|                         |   |  |            |
|-------------------------|---|--|------------|
|                         |   | <p>access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Single Camera Colonoscope, Processor-A and other ancillary equipment.</p> <p><b><u>Saneso Gastroscope 360- A</u></b></p> <p>The Saneso Single Camera Gastroscope is intended for diagnostic visualization of the upper gastrointestinal tract (including the esophagus, stomach and duodenum). The system also provides access for therapeutic interventions using standard endoscopy tools. The Saneso system consists of Saneso Single Camera Gastroscope, Processor-A and other ancillary equipment.</p> |            |
| Construction/Components | The Saneso Lens Wash Bottle Assembly consists of a reusable lens wash water bottle (cylindrical bottle), cap, tubing set and air/water connector. | The Saneso Lens Wash Bottle consists of a single-use lens wash water bottle (cylindrical bottle), cap, tubing set and air/water connector.   | Equivalent |
| Sterile/Non Sterile     | Sold Non - Sterile  | Sold Non - Sterile   | Equivalent |
| Accessories             | <ul style="list-style-type: none"> <li>● Cap CO2 outlet</li> <li>● Water line</li> <li>● Cap water outlet</li> </ul>                              | <ul style="list-style-type: none"> <li>● Cap CO2 outlet</li> <li>● Water line</li> <li>● Cap water outlet</li> <li>● CO2 line</li> </ul>   | Equivalent |





|                                  |  |  |            |
|----------------------------------|--|--|------------|
|                                  | <ul style="list-style-type: none"> <li>● CO2 line</li> <li>● Air Line</li> <li>● CO2 clamp</li> <li>● Cap Air inlet</li> </ul> | <ul style="list-style-type: none"> <li>● Air Line</li> <li>● CO2 clamp</li> <li>● Cap Air inlet</li> </ul> |            |
| Usage                            | Reusable   | Single Use   | Different  |
| Method of Application            | Manual   | Manual   | Equivalent |
| Compatible Endoscopes/processors | Saneso - A series (K210052), Saneso Processor - A  | Saneso - A series (K210052), Saneso Processor - A  | Equivalent |

**5.6.1. Similarities**

The Saneso Lens Wash Bottle Assembly SAN-LWB-A is similar to the Lens Wash Bottle which is cleared under K210052 in terms of intended use, construction/components, method of application.

**5.6.2. Differences**

The difference between the Saneso Lens Wash Bottle Assembly SAN-LWB-A and the Lens Wash Bottle which is cleared under K210052 is its usage, ie; from single use to multi-use. However, this difference does not raise new questions of safety or efficacy. Differences were justified using Non-Clinical performance testing - Manual Cleaning Validation and Steam Sterilization Validation and demonstrates that the Lens wash bottle Assembly is safe and effective.

**5.7. Summary of Non-Clinical Tests**

Performance testing was performed as per the relevant standards and test procedures listed in the Declaration of Conformity.

1. Manual Cleaning Validation - Cleaning validation was performed using Saneso’s recommended cleaning procedure. Saneso Lens Wash Bottle Assembly met all the acceptance criteria and demonstrated that the cleaning process does not impact the functionality of the device.
2. Steam Sterilization Validation - Sterilization validation was performed and Saneso Lens Wash Bottle Assembly met all the acceptance criteria. The device is found to be safe and effective.

**5.8. Summary of Clinical Tests**

**5.8.1. Animal Study**

No animal studies have been conducted for the Subject Device.

**5.8.2. Clinical Study**

No clinical studies have been conducted for the Subject Device.

**5.9. Conclusion**

Saneso Lens Wash Bottle Assembly SAN-LWB-A is substantially equivalent to the Lens Wash Bottle which is cleared under K210052 in terms of intended use and technological characteristics. Differences were justified using Non-Clinical performance testing. Thus, Saneso Inc. has concluded that the device does not introduce any significant questions of safety and efficacy and is substantially equivalent.