

April 13, 2023

JR Engineering & Medical Technologies (M) SDN. BHD. % Manoj Zacharias Consultant Liberty Management Group Ltd. 75 Executive Dr. STE 114 Aurora, Illinois 60504

Re: K222350

Trade/Device Name: Sterile Nitrile Surgical Gloves Powder Free

Regulation Number: 21 CFR 878.4460

Regulation Name: Non-Powdered Surgeon's Glove

Regulatory Class: Class I, reserved

Product Code: KGO Dated: March 13, 2023 Received: March 13, 2023

Dear Manoj Zacharias:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K222350			
Device Name			
Sterile Nitrile Surgical Gloves Powder Free			
ndications for Use (Describe)			
A Sterile Nitrile Surgical Gloves Powder Free is a device made of synthetic rubber intended to be worn by operating room personnel to protect a surgical wound from contamination.			
Type of Use (Select one or both, as applicable)			
☐ Prescription Use (Part 21 CFR 801 Subpart D) ☐ Over-The-Counter Use (21 CFR 801 Subpart C)			
CONTINUE ON A SEPARATE PAGE IF NEEDED.			

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510K SUMMARY

K222350

Date of Preparation: April 10, 2023 As required by: 21CFR § 807.92

A. APPLICANT INFORMATION

Submitter Name JR Engineering & Medical Technologies (M) SDN. BHD
Address Lot 8 & 10, Jalan Zurah 3 & Lot 1&3, Jalan Zurah 3A/1, Pusat

Perindustrian 2, 44200 Rasa, Hulu

Selangor, Selangor Darul Ehsan, Malaysia.

Phone +603-60572081 Fax +603-60572181

E-mail <u>ganeshjrmt@gmail.com</u>
Contact Person Mr. Ganesan Subramaniam

Designation Managing Director
Contact Number +6012 224 6677
Contact Email ganeshjrmt@gmail.com

US Agent Manoj Zacharias

US Agent Address Liberty Management Group Ltd.

75 Executive Dr. Ste 114 Aurora, IL 60504 USA

US Agent Telephone 1 (630) 270-2921 US Agent Fax 1 (815) 986-2632

US Agent Email manoj@libertymanagement.us

B. Device Identification

Common Name Surgical Gloves

Device Name Sterile Nitrile Surgical Gloves Powder Free

Product Proprietary or JR MEDIC

Trade Name

Classification Name Surgeon's Gloves

Device Classification 1
Product Code KGO

Regulation Number 21 CFR 878.4460

C. Predicate Device

510k Number K170515 Common Name Surgical Gloves

Device Name Sterile Nitrile Surgical Gloves, Powder Free

Classification name Surgeon's Gloves

Device Classification 1
Product Code KGO

Regulation Number 21 CFR 878.4460

D. Description of the Device

The proposed device, Sterile Nitrile Surgical Gloves Powder Free is a sterile and disposable medical glove intended to be worn by operating room personnel to protect a surgical wound from contamination.

The proposed device is made of synthetic rubber latex, as per standard ASTM D3577-2019, Standard Specification for Rubber Surgical Gloves.

The classification is: Type II - gloves compounded from synthetic rubber latex.

The proposed device is Sterile Nitrile Surgical Gloves Powder Free, and is produced in sizes $6, 6 \frac{1}{2}, 7, 7 \frac{1}{2}, 8, 8 \frac{1}{2}$, and 9. All sizes share the same White color.

The proposed device is sterilized using Gamma Radiation method to achieve the Sterility Assurance Level (SAL) of 10⁻⁶ while packaged to maintain sterility. The shelf life is 3 years.

E. Technological Characteristics

Characteristic	Subject device K222350	Predicate device K170515	Remarks
Product Code	KGO	KGO	Same
Regulation No.	21 CFR 878.4460	21 CFR 878.4460	Same
Class	1	1	Same
	A Sterile Nitrile Surgical Gloves	The powder-free nitrile surgical	
	Powder Free is a device made of	Gloves, is a sterile and single use	
Intended Use	synthetic rubber intended to be	devise made of synthetic rubber	
	worn by operating room	intended to be worn by operation	
	personnel to protect a surgical	room personnel to protect a surgical	Similar
	wound from contamination.	wound from contamination. The	
		gloves do not contain lubricating or	
		dusting powder	
Powdered or Powder free	Powder free	Powder free	Same
Material	Synthetic Rubber	Synthetic Rubber	Same
Color	White	White	Same
Classification as per	Type II - gloves compounded	Type II - gloves compounded from	Same
ASTMD3577-19,	from Synthetic rubber latex	Synthetic rubber latex	
Standard Specification			
for Rubber Surgical			
Gloves			
Sterilization	Radiation, SAL- 10 ⁻⁶	ETO/as well as Radiation, SAL- 10-6	Similar
Label and Labeling	Meet FDA's label Requirements	Meet FDA's label Requirements	Same
Type of use	Over the counter use	Over the counter use	Same
Dimensions Length:-	307mm	265mm min	
Min 265mm			Similar
Size: Width			
6.0 (76±6mm)	78mm	76±6mm	
6.5 (83±6mm)	85mm	83±6mm	Similar
7.0 (89±6mm)	88mm	89±6mm	
7.5 (95±6mm)	97mm	95±6mm	
8.0 (102±6mm)	103mm	102±6mm	
8.5 (108±6mm)	110mm	108±6mm	
9.0 (114±6mm)	116mm	114±6mm	
Cuff, Palm, Finger Tip	Cuff- 0.12mm	Cuff, Palm & Finger has a min	Similar
Thickness	Palm- 0.18mm	0.10mm thickness	
Min 0.10mm	Finger Tip- 0.21mm		
Tensile Strength			
Unaged	24.20MPa	17MPa min	Similar
17MPa minimum			
Ultimate Elongation			
Unaged	871%	650% min	Similar
650% minimum			

Stress at 500%			
Unaged	6.1MPa	7.0MPa max	Similar
7.0 MPa Max			
Tensile Strength	18.53MPa	12 MPa min	Similar
Aged			
12MPa minimum			
Ultimate Elongation			
Aged	653%	490% min	Similar
490% minimum			
Freedom from Holes	AQL 1.5	AQL 2.5	
			Similar
Powder residue for	0.40mg/Glove	< 2mg/Glove	Similar
powder free glove			
Powder content			
< 2 mg/Glove			
Skin Irritation	Under conditions of the testing,	Under conditions of the testing, not	Similar
	not an irritant	an irritant	
Skin Sensitization	Under the conditions of the	Under the conditions of the testing,	Similar
	testing, not a sensitizer	not a sensitizer	
In vitro cytotoxicity	Under the conditions of the	Under the conditions of the testing,	Similar
	testing, not cytotoxic	not cytotoxic	
Material Mediated	Under the conditions of the	No data available	Different
pyrogenicity	testing, non-pyrogenic		
Bacterial Endotoxin Test	< 20 EU/glove	No data available	Different
Systemic Toxicity	Under the conditions of the	No data available	Different
	testing, no acute systemic		
	toxicity		
Sterilization Modality	Radiation	Radiation/EtO	Similar
Sterility Assurance Level	10-6	10-6	Same
Shelf Life	3 years	Unknown	Different

F. Summary of Non-Clinical Performance Testing

Test Method	Purpose	Acceptance Criteria	Result
ASTM D5151-2019	Freedom from holes	AQL 1.5	Pass
	Before aging	3 non-consecutive lots	
ASTM D5151-2019	Freedom from holes	AQL 1.5	Pass
	After accelerated aging	3 non-consecutive lots	
	70±2°C for 166±2 h		
ASTM D5151-2019	Freedom from holes	AQL 1.5	Pass
	After accelerated aging	3 non-consecutive lots	
	50±2°C for 90±1 days		
ASTM D3577-19,	Tensile strength	Minimum 17MPa for all	Pass
Standard	Before aging	sizes	
Specification for		AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM D3577-19,	Tensile strength	Minimum 12MPa for all	Pass
Standard	After accelerated aging	sizes	
Specification for	70±2°C for 166±2 h	AQL 4%	
Rubber Surgical		3 non-consecutive lots	

Gloves			
ASTM D3577-19,	Tensile strength	Minimum 12MPa for all	Pass
Standard	_	sizes	rass
	After accelerated aging		
Specification for	50±2°C for 90±1 days	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			_
ASTM D3577-19,	Ultimate Elongation	Minimum 650% for all	Pass
Standard	Before aging	sizes	
Specification for		AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM D3577-19,	Ultimate elongation	Minimum 490% for all	Pass
Standard	After accelerated aging	sizes	
Specification for	70±2°C for 166±2h	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM D3577-19,	Ultimate elongation	Minimum 490% for all	Pass
Standard	After accelerated aging	sizes	
Specification for	50±2°C for 90±1 days	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM D3577-19,	Stress at 500% elongation	Maximum 7.0MPa for all	Pass
Standard	Before accelerated aging	sizes	1 433
Specification for	Before accelerated aging	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
_		3 non-consecutive fors	
Gloves	G 5000/ 1	70.60	D.
ASTM D3577-19,	Stress at 500% elongation	Maximum 7.0MPa for all	Pass
Standard	After accelerated aging	sizes	
Specification for	70±2°C for 166±2h	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM D3577-19,	Stress at 500% elongation	Maximum 7.0MPa for all	Pass
Standard	After accelerated aging	sizes	
Specification for	50±2°C for 90±1 days	AQL 4%	
Rubber Surgical		3 non-consecutive lots	
Gloves			
ASTM F1929,	Package integrity	No penetration	Pass
Test Method for	After accelerated aging	AQL 0.65%	
detecting seal leaks in	70±2°C for 166±2h	3 non-consecutive lots	
porous medical			
packaging by dye			
penetration			
ASTM F1929,	Package integrity	No penetration	Pass
Test Method for	After accelerated aging	AQL 0.65%	
	50±2°C for 90±1 days	3 non-consecutive lots	
porous medical			
packaging by dye			
penetration			
ASTM D3577-19,	To determine the length of	Min 265mm for all sizes	Size 6 Pass
Standard	the gloves		Size 6 ½ Pass
Specification for			Size 7 Pass
Rubber Surgical			1 ass
reaccor buildioni	l .	1	

C1	T	1		T	_
Gloves				Size 7 ½	Pass
				Size 8	Pass
				Size 8 ½	Pass
				Size 9	Pass
ASTM D3577-19,	To determine the width of	Size 6	76±6mm	Size 6	Pass
Standard	the gloves	Size 6 ½	83±6mm	Size 6 ½	Pass
Specification for		Size 7	89±6mm	Size 7	Pass
Rubber Surgical		Size 7 ½	95±6mm	Size 7 ½	Pass
Gloves		Size 8	102±6mm	Size 8	Pass
		Size 8 ½	108±6mm	Size 8 ½	Pass
		Size 9	114±6mm	Size 9	Pass
ASTM D3577-19,	To determine finger	Minimum 0	.10mm for all	Size 6	Pass
Standard	thickness	sizes		Size 6 ½	Pass
Specification for		AQL 4%		Size 7	Pass
Rubber Surgical				Size 7 ½	Pass
Gloves				Size 8	Pass
				Size 8 ½	Pass
				Size 9	Pass
A CTM D2577 10	T. 1.		10 C 11		
ASTM D3577-19, Standard	To determine palm		.10mm for all	Size 6	Pass
	thickness	sizes		Size 6 ½	Pass
Specification for		AQL 4%		Size 7	Pass
Rubber Surgical Gloves				Size 7 ½	Pass
Gloves				Size 8	Pass
				Size 8 ½	Pass
				Size 9	Pass
ASTM D3577-19,	To determine cuff thickness	Minimum 0	.10mm for all	Size 6	Pass
Standard		sizes		Size 6 ½	Pass
Specification for		AQL 4%		Size 7	Pass
Rubber Surgical				Size 7 ½	Pass
Gloves				Size 8	Pass
				Size 8 ½	Pass
					Pass
ASTM D6124-	Less than 2mg/glove	Maximum 2	2mg/glove for	Pass	
06R17,		all sizes	22		
Standard Test					
Method for Residual					
Powder on Medical					
Gloves					
ISO 10993-10: 2010,	Irritation	Under the co	onditions of	Pass	
Biological evaluation		the testing,	not an irritant		
of medical devices -					
Part 10- Tests for					
irritation and skin					
sensitization	g vi vi	TT 1 '	1:.: 2	D.	
	Sensitization		onditions of	Pass	
Biological evaluation		the testing, not a			
of medical devices - Part 10- Tests for		sensitizer			
irritation and skin					
sensitization					
ISO 10993-5:2009(E)	Cytotoxicity	Under the co	onditions of	Pass	
150 10775-3.2009(E)	CytotoAlCity	onder the C	ondinons of	1 433	

Biological Evaluation of Medical Devices -		the testing, non-cytotoxic	
Part 5-Tests for			
in-vitro Cytotoxicity			
ISO 10993-	Acute systemic toxicity	Under the conditions of	Pass
11:2017(E)		the testing, no acute	
Biological Evaluation		systemic toxicity	
of Medical Devices -			
Part 11- Tests for			
Systemic			
Toxicity and			
Biological Tests			
USP 41	Material mediated	Under the conditions of	Pass
<151>Pyrogen Test	pyrogenicity	the testing, non-pyrogenic	
USP 42 <85>	Bacterial endotoxin limits	Less than 20EU/glove	Pass
Bacterial Endotoxin			
Test			

G. Summary of Clinical Performance Testing

No clinical performance testing was performed in support of this submission.

H. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device Sterile Nitrile Surgical Gloves Powder Free is as safe, as effective and performs as well or better than the legally marketed predicate device "PRIMUS NITRILE GLOVES" Sterile Nitrile Surgical Gloves Powder Free (K170515).