

April 2, 2023

Beijing Rongrui-Century Science & Technology Co., Ltd. % Ray Wang
General Manager
Beijing Believe-Med Technology Service Co., Ltd.
Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd.,
FangShan District
Beijing, Beijing 102401
China

Re: K222370

Trade/Device Name: SpO2 Extension Cable Regulation Number: 21 CFR 870.2900

Regulation Name: Patient Transducer And Electrode Cable (Including Connector)

Regulatory Class: Class II

Product Code: DSA Dated: March 3, 2023 Received: March 3, 2023

# Dear Ray Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

K222370 - Ray Wang Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# James J. Lee -S

James J. Lee, Ph.D.
Director
DHT1C: Division of Sleep Disordered
Breathing, Respiratory and
Anesthesia Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

### DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

# **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K222370				
Device Name SpO2 Extension Cable				
dications for Use ( <i>Describe</i> ) pO2 Extension Cables are intended to be used to connect SpO2 sensors, placed at appropriate sites on the patient to a nonitoring device for SpO2 monitoring by health care professional.				
Type of Use (Select one or both, as applicable)    Prescription Use (Part 21 CFR 801 Subpart D)   Over-The-Counter Use (21 CFR 801 Subpart C)				
CONTINUE ON A SEPARATE PAGE IF NEEDED.				

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

> Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

The assigned 510(k) Number: K222370

# **Tab #3 510(k) Summary**

This 510(k) Summary of 510(k) safety and effectiveness information is being submitted in accordance with requirements of SMDA 1990 and Title 21, CFR Section 807.92.

- 1. Date of Preparation:2023/03/29
- 2. Sponsor Identification

#### Beijing Rongrui-Century Science & Technology Co., Ltd.

3rd Floor, West zone, No.1 Building, No.7 Yard, Fengxian middle Road, HaiDian District, Beijing, 100094, P.R.China.

Contact Person: Calen Chen Position: Quality Manager Tel: +86-769-22243880 Fax: +86-769-22242882 Email: info@r-rui.com

3. Designated Submission Correspondent

Mr. Ray Wang

# Beijing Believe-Med Technology Service Co., Ltd.

Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd., FangShan District, Beijing, 102401, China

Tel: +86-18910677558 Fax: +86-10-56335780

Email: information@believe-med.com

# 4. Identification of Proposed Device

Trade Name: SpO2 Extension Cable

Common Name: Cable, Transducer And Electrode, Patient, (Including Connector)

#### Regulatory Information

Classification Name: Cable, Transducer And Electrode, Patient, (Including Connector)

Classification: II Product Code: DSA

Regulation Number: 21 CFR 870.2900

Review Panel: Cardiovascular

#### Indication for use Statement:

SpO2 Extension Cables are intended to be used to connect SpO2 sensors, placed at appropriate sites on the patient to a monitoring device for SpO2 monitoring by health care professional.

#### Device Description:

SpO2 Extension Cables are intended to be used to connect SpO2 sensors, placed at appropriate sites on the patient to a monitoring device for SpO2 monitoring by health care professional.

The SpO2 Extension Cable is the replacements for similar cables manufactured by Original Equipment Manufacturers (OEM) and other third party after market manufacturers for their respective monitors.

The SpO2 Extension Cable is comprised of Plug, Cable/Leadwires and Connector. It's intended to plug into monitoring device and connect with SpO2 sensor, for transmitting signals which generated by SpO2 sensor to monitoring device.

By using the same types of construction and technological characteristics to the compatible patient monitor and SpO2 sensor, the SpO2 Extension Cable can avoid measured data corrupted.

The compatible sensors and monitors are shown in the following table:

NO.	Model	Compatible Rongrui Sensor	Compatible Monitors	
			& Nellcor OEM Boards	
1	RCT006	RSJ002A, RSA002DN		
		RSJ002I, RSA002DI	NELLCOR: N-600X、N-595	
		RST002A, RSA002DP	NELLCOR: N-000A\ N-393	
		RST002I, RSA002DA		
		RSW002N, RSY002N		
2	RCT018	RSJ001A,RSJ001I,		
		RST001A,RST001I,		
		RSW001N,RSY001N,	Philips: CM100、MP2、X2、MP5	
		RSA001DN,RSA001DI,		
		RSA001DP,RSA001DA		

## 5. Identification of Predicate Device(s)

Predicate Device

K192404

Patient Monitoring Cables (SpO2 Extension Cable)

Shenzhen Coreray Technology Co., Ltd.

#### 6. Non-Clinical Test Conclusion

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

- ➤ ISO 10993-5: 2009 Biological Evaluation Of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity
- ➤ ISO 10993-10: 2010 Biological Evaluation Of Medical Devices Part 10: Tests For Irritation And Skin Sensitization.
- ➤ IEC 60601-1:2012, Medical Electrical Equipment Part 1: General requirements for basic safety and essential performance
- ➤ IEC 60601-1-2:2014, Medical Electrical Equipment-Part 1-2: General Requirements For Basic Safety And Essential Performance Collateral Standard: Electromagnetic Compatibility-Requirements And Tests
- ➤ ISO 80601-2-61:2017 Medical electrical equipment-Part 2-61: Particular requirements for basic safety and essential performance of pulse oximeter equipment
- ➤ ANSI AAMI EC53:2013/(R)2020 ECG trunk cables and patient leadwires

SpO2 Extension Cable has also been evaluated the performance accuracy through Compatibility testing, which proves that no measured data corrupt during communication between SpO2 sensors and host monitors.

#### 7. Clinical Test Conclusion

No clinical study is included in this submission.

#### 8. Substantially Equivalent (SE) Comparison

Table 1 General Comparison

Item	Proposed Device	Predicate Device	Remark
		Patient Monitoring Cables (SpO2	
Product name	SpO2 Extension Cable	Extension Cable)	/
		K192404	
Classification	21 CFR 870.2900	21 CFR 870.2900	SAME

Regulation				
Classification	II	II	SAME	
Product Code	DSA	DSA	SAME	
Common Name	Cable, Transducer And Electrode,	Cable, Transducer And Electrode,	SAME	
Common Name	Patient, (Including Connector)	Patient, (Including Connector)		
	SpO2 Extension Cables are intended	Patient Monitoring Cables are intended		
	to be used to connect SpO2 sensors,	to be used to connect sensors, placed at	SAME	
Indications	placed at appropriate sites on the	appropriate sites on the patient to a		
for use	patient to a monitoring device for	monitoring device for general		
	SpO2 monitoring by health care	monitoring by health care professional.		
	professional.			
Design	Plug, Cable/Leadwires and connectors	Plug, Cable/Leadwires and connectors	SAME	
Usage	Reusable	Reusable	SAME	
Cable lengths	$2m \pm 10\%$	8ft	Difference1	
Wire material	Shielded & Unshielded Copper with	Shielded & Unshielded Copper with	Difference2	
Wife material	PVC or TPE jacket	TPU jacket	Difference2	
Sterility	Non-sterile	Non-sterile	SAME	
Biocompatibility	Meets ISO 10993-5 Cytotoxicity, ISO	Meets ISO 10993-5 Cytotoxicity, ISO	SAME	
	10993-10 Sensitization and Irritation	10993-10 Sensitization and Irritation	SAME	
Electrical safety	IEC 60601-1	IEC 60601-1		
and Performance	ISO 80601-2-61	ISO 80601-2-61	SAME	
and Performance	ANSI AAMI EC53	ANSI AAMI EC53	_	
EMC	IEC 60601-1-2	IEC 60601-1-2	SAME	

### Analysis:

# Difference 1:

There is difference on the cable length, the cable length of the proposed device is  $2m \pm 10\%$  (Approximately 6.6 ft), the cable length of the predicate device is 8ft, but these cables are same function , the cable length difference would not affect its safety and effectiveness.

# Difference 2:

Although the wire material between the proposed device and predicate device have few differences, the proposed device comply with electrical safety standard IEC 60601-1 and the performance standard ISO 80601-2-61 and ANSI AAMI EC53, so the differences would not affect its safety and effectiveness.

# 9. Substantially Equivalent (SE) Conclusion

Based on the comparison and analysis above, the proposed devices are determined to be Substantially Equivalent (SE) to the predicate devices.