

August 4, 2023

Trinity Biotech (Primus Corporation, dba Trinity Biotech)
Kaitlyn Eastman
Regulatory Supervisor
4231 E. 75th Terrace
Kansas City, Missouri 64132

Re: K222635

Trade/Device Name: Premier Resolution System

Regulation Number: 21 CFR 864.7415

Regulation Name: Abnormal Hemoglobin Assay

Regulatory Class: Class II Product Code: GKA Dated: August 30, 2022 Received: August 31, 2022

# Dear Kaitlyn Eastman:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Min Wu, Ph.D.
Branch Chief
Division of Immunology and Hematology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

# **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

K222635
Device Name
Premier Resolution System
Indications for Use (Describe) The Premier Resolution System is an automated High Performance Liquid Chromatography (HPLC) system which performs the separation of hemoglobin species in venous whole blood samples for the quantitative analysis of normal hemoglobin (A, A2, and F), and the qualitative detection of major variant hemoglobin S, C, D-Los Angeles, and E in adult, adolescent, children and infant populations. The assays are performed on venous whole blood samples collected in tubes containing K2EDTA as anticoagulant.
The Premier Resolution System is intended for Professional Laboratory Use only.
The Premier Resolution System is intended for use with analytical components and reagents provided by Trinity Biotech.
The Premier Resolution System is intended to be used in conjunction with other laboratory and clinical findings.
For In Vitro Diagnostic Use.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

# \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



# Section 5: 510(k) Summary

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of 21 CFR §807.92:

#### 1. Submitter / Sponsor

Trinity Biotech 4231 E. 75<sup>th</sup> Terrace Kansas City, MO 64132 USA

Phone: 716.483.7423 Fax: 816.361.1974

Contact Person: Kaitlyn Eastman Date Prepared: July 25, 2023

#### 2. Device

Name of Device: Premier Resolution System Classification Name: Abnormal Hemoglobin Assay

Regulation: 21 CFR § 864.7415

Regulatory Class: Class II Product Code: GKA

### 3. Predicate Device

Bio-Rad Variant II β -thalassemia (K991127)

# 4. Device Description

The Premier Resolution System consists of a high performance liquid chromatographic analyzer, reagents, analytical column and software which allows for the fractionation and quantitation of fetal hemoglobin (Hb F), and hemoglobin A2 (Hb A2), and with fractionation and presumptive identification of abnormal hemoglobin variants. This is accomplished using the principles of ion-exchange (IEX) high performance liquid chromatography (HPLC).

#### 5. Indications for Use Statement and Intended Uses

The Premier Resolution System is an automated High Performance Liquid Chromatography (HPLC) system which performs the separation of hemoglobin species in venous whole blood samples for the quantitative analysis of normal hemoglobin (A, A2, and F), and the qualitative detection of major variant hemoglobin S, C, D-Los Angeles, and E in adult, adolescent, children and infant populations. The assays are performed on venous whole blood samples collected in tubes containing K2EDTA as anticoagulant.

The Premier Resolution System is intended for Professional Laboratory Use only.



The Premier Resolution System is intended for use with analytical components and reagents provided by Trinity Biotech.

The Premier Resolution System is intended to be used in conjunction with other laboratory and clinical findings.

For In Vitro Diagnostic Use.

# 6. Comparison of Technological Characteristics with the Predicate Device

	Predicate					
Aspect	Bio-Rad Variant <sup>TM</sup> II Beta	Premier Resolution System				
or Feature	Thalassemia K991127	Tremier Resolution bystem				
Intended Use	The Variant <sup>™</sup> II B-thalassemia Program is intended for the separation and area percent determinations of hemoglobins A2 and F and as an aid in the identification of abnormal hemoglobins in whole blood using ion-exchange high performance liquid chromatography (HPLC).  The Variant <sup>™</sup> II B-thalassemia Program is intended for use only with the Bio-Rad Variant <sup>™</sup> II Hemoglobin Testing System.  For <i>in vitro</i> diagnostic use only.	The Premier Resolution System is an automated High Performance Liquid Chromatography (HPLC) system which performs the separation of hemoglobin species in venous whole blood samples for the quantitative analysis of normal hemoglobin (A, A2, and F), and the qualitative detection of major variant hemoglobin S, C, D-Los Angeles, and E in adult, adolescent, children and infant populations. The assays are performed on venous whole blood samples collected in tubes containing K2EDTA as anticoagulant. The Premier Resolution System is intended for Professional Laboratory Use only.  The Premier Resolution System is intended for use with analytical components and reagents provided by Trinity Biotech.  The Premier Resolution System is intended to be used in conjunction with other laboratory and clinical findings.  For <i>In Vitro</i> Diagnostic Use.				
Chemistry	Ion Exchange HPLC	Ion Exchange HPLC				
Collection Tubes	K₂EDTA, K₃EDTA	K₂EDTA				
Sample Tube Processing	Aspiration of hemolysate from closed tube	Aspiration of whole blood from closed tube				
Sample Hemolysis	Performed automatically by the system	Performed automatically by the system				
Automated Sample Introduction	Continuous loading with sample racks	Continuous loading with sample racks				
Separation System	Ion-exchange high performance liquid chromatography (HPLC)	Ion-exchange high performance liquid chromatography (HPLC) protein separation				



# **Premier Resolution System**

Premarket Notification

Aspect or Feature	<b>Predicate</b> Bio-Rad Variant™ II Beta Thalassemia K991127	Premier Resolution System					
	protein separation on analytical column based on ionic interaction with the column material and elution by buffer gradient with increasing ionic strength.	on analytical column based on ionic interaction with the column material and elution by buffer gradient with increasing ionic strength.					
Separation Unit	Analytical Column	Analytical Column					
Analysis Throughput	About 6 Minutes	About 4 Minutes - Quick Scan About 8 Minutes - High Resolution					
Calibration	A2 + F Calibrator	A2+F Calibrator					
Control	A2+F Control	A2+F Control					
Use of Other Controls	FASC Position Marker	FASC Position Marker					
Sample Identification	Barcode on sample tube	Barcode on sample tube					
Absorbance wavelength	About 413 nm	413 nm					
Hb Variants Library	Display and Operator Manual	Display and Operator Manual					

#### 7. Performance Data

The following performance data were provided in support of the substantial equivalence determination.

#### **Correlation (Method Comparison)**

To establish equivalence, the method comparison study was conducted at three (3) professional external laboratory sites with both the Premier Resolution System and the Bio-Rad Variant<sup>TM</sup> II Beta Thalassemia. A total of 780 unique patient samples were collected and analyzed in both the Quick Scan and High Resolution assay modes. The following results are presented for the Premier Resolution Quick Scan and High Resolution assay modes:

#### **Premier Resolution Quick Scan Assay**

The Premier Resolution is comparable to the Bio-Rad Variant II for HbF on the Quick Scan Assay by a mean bias of -0.3 between the intervals of 1.1 to 48.9% with 160 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbA on the Quick Scan Assay by a mean bias of 0.7 between the intervals of 2.5 to 89.7% with 682 patient results.





The Premier Resolution is comparable to the Bio-Rad Variant II for HbA2 on the Quick Scan Assay by a mean bias of 0.1 between the intervals of 1.6 to 6.1% with 602 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbS on the Quick Scan Assay by a mean bias of 0.3 between the intervals of 6.8 to 67.1% with 106 patient samples.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbC on the Quick Scan Assay by a mean bias of -1.1 between the intervals of 9.5 to 82.8% with 49 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbD-LA on the Quick Scan Assay by a mean bias of 1.6 between the intervals of 11.6 to 82.7% with 17 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbE on the Quick Scan Assay by a mean bias of -3.0 between the intervals of 5.5 to 70.4% with 25 patient results.

## **Premier Resolution High Resolution Assay**

The Premier Resolution is comparable to the Bio-Rad Variant II for HbF on the High Resolution Assay by a mean bias of -0.4 between the intervals of 1.1 to 46.6% with 158 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbA on the High Resolution Assay by a mean bias of 2.4 between the intervals of 3.5 to 90.5% with 586 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbA2 on the High Resolution Assay by a mean bias of 0.1 between the intervals of 1.6 to 6.0% with 598 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbS on the High Resolution Assay by a mean bias of 1.3 between the intervals of 1.9 to 67.9% with 110 patient samples.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbC on the High Resolution Assay by a mean bias of -1.0 between the intervals of 10.2 to 82.5% with 49 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbD-LA on the High Resolution Assay by a mean bias of 2.7 between the intervals of 11.7 to 84.1% with 17 patient results.

The Premier Resolution is comparable to the Bio-Rad Variant II for HbE on the High Resolution Assay by a mean bias of -4.9 between the intervals of 5.3 to 66.7% with 25 patient results.



# **Single Site Precision**

The Single-Site precision study was conducted following *EP05-A3 Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition.* Samples of varying concentrations of HbA, A2, and F were analyzed on the Premier Resolution using the Quick Scan and High-Resolution assay mode for the quantitation hemoglobin fractions. The retention time and relative retention time precision was also determined for each hemoglobin fraction. The 20x2x2 study design was conducted over 20 days, with two runs per day and two replicates per run. The study yielded 80 data points for each analyte to establish a balanced dataset.

# **Single Site Precision - Quick Scan**

Sample		Mean	Repeatability		ty Between-Run Between-Day			Within- Laboratory		
Description	N	Value	SD	%CV	SD	%CV	SD	%CV	SD	%CV
A High	80	83.56	0.25	0.30%	0.37	0.44%	2.98	3.57%	3.01	3.61%
A Mid	80	51.24	0.35	0.68%	0.15	0.30%	0.25	0.50%	0.46	0.89%
A2 Mid	80	2.57	0.04	1.63%	0.03	1.07%	0.03	1.08%	0.06	2.23%
A2 Low	80	1.75	0.04	2.02%	0.02	0.90%	0.10	5.57%	0.10	5.99%
F High	80	11.38	0.05	0.43%	0.03	0.26%	0.10	0.92%	0.12	1.05%
F Mid	80	1.67	0.04	2.42%	0.00	0.00%	0.04	2.18%	0.05	3.26%
S High	80	74.63	0.33	0.44%	0.18	0.24%	0.55	0.74%	0.67	0.89%
S Mid	80	30.61	0.14	0.47%	0.14	0.45%	0.22	0.73%	0.30	0.98%
C High	80	81.22	0.41	0.50%	0.43	0.53%	0.23	0.28%	0.63	0.78%
C Mid	80	31.91	0.46	1.43%	0.26	0.81%	0.20	0.61%	0.56	1.75%
D High	80	69.98	0.70	1.00%	0.71	1.01%	0.86	1.23%	1.32	1.88%
D Mid	80	38.14	0.30	0.79%	0.37	0.97%	0.61	1.61%	0.78	2.04%
E High	80	81.22	1.38	1.70%	0.00	0.00%	1.85	2.28%	2.31	2.84%
E Mid	80	22.30	0.14	0.62%	0.17	0.75%	0.65	2.90%	0.68	3.06%

# Single Site Precision – High Resolution

Sample		Mean	Repeatability		Between-Run		Between-Day		Within- Laboratory	
Description	N	Value	SD	%CV	SD	%CV	SD	%CV	SD	%CV
A High	80	83.64	0.43	0.52%	0.38	0.45%	2.98	3.56%	3.03	3.63%
A Mid	80	50.60	0.41	0.81%	0.48	0.94%	0.30	0.59%	0.69	1.37%
A2 Mid	80	2.48	0.10	4.13%	0.15	5.92%	0.00	0.00%	0.18	7.22%
A2 Low	80	1.59	0.06	3.92%	0.09	5.98%	0.00	0.00%	0.11	7.15%
F High	80	11.11	0.12	1.07%	0.03	0.27%	0.37	3.31%	0.39	3.49%
F Mid	80	1.34	0.08	5.78%	0.04	3.23%	0.08	6.24%	0.12	9.10%

# **Premier Resolution System**

Premarket Notification

Sample		Mean	Repeatability		Between-Run		Between-Day		Within- Laboratory	
Description	N	Value	SD	%CV	SD	%CV	SD	%CV	SD	%CV
S High	80	76.89	0.57	0.74%	0.42	0.54%	0.74	0.96%	1.02	1.33%
S Mid	80	32.01	0.30	0.92%	0.23	0.73%	0.39	1.22%	0.54	1.69%
C High	80	81.07	0.51	0.63%	0.49	0.61%	0.61	0.75%	0.94	1.15%
C Mid	80	32.84	0.44	1.34%	0.00	0.00%	0.50	1.52%	0.67	2.03%
D High	80	73.83	1.07	1.45%	0.89	1.21%	0.84	1.14%	1.62	2.20%
D Mid	80	39.80	0.34	0.84%	0.23	0.58%	0.33	0.84%	0.53	1.32%
E High	80	77.26	0.93	1.20%	0.94	1.21%	3.33	4.31%	3.58	4.63%
E Mid	80	22.46	0.13	0.57%	0.22	0.98%	0.66	2.93%	0.70	3.14%

# **Multisite Precision**

The Multisite Precision Study was conducted in accordance with CLSI EP05-A3 Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline -Third Edition. A 3x5x5 study design was performed across three (3) external sites over five (5) days with five (5) replicates per day. Three (3) Premier Resolution Systems, one per site and Mobile Phase 1, 2, and Diluent reagents (one lot of each), analytical column, and precision samples for the study were provided to the external sites. The study assessed instrument-to-instrument, operator-to-operator and site-to-site precision.

#### Multisite Precision - Quick Scan

Quick Scan													
Sample		Mean	Repea	Repeatability		Repeatability Between-Day Be			Betwe	Between-Site		Reproducibility	
Description	N	Value	SD	%CV	SD	%CV	SD	%CV	SD	%CV			
FASC F	75	17.16	0.09	0.52%	0.27	1.55%	0.82	4.80%	0.87	5.07%			
FASC A	75	42.89	0.22	0.52%	0.33	0.78%	0.34	0.78%	0.52	1.22%			
FASC S	75	12.87	0.07	0.51%	0.06	0.50%	0.09	0.69%	0.13	0.99%			
FASC C	75	14.14	0.11	0.78%	0.08	0.57%	0.26	1.87%	0.30	2.11%			
A2+F1F	75	2.94	0.04	1.33%	0.04	1.39%	0.12	4.17%	0.14	4.59%			
A2+FIA	75	82.75	0.31	0.38%	1.00	1.21%	1.16	1.40%	1.57	1.89%			
A2+F I A2	75	3.11	0.04	1.17%	0.06	2.04%	0.06	2.08%	0.10	3.15%			
A2+F II F	75	7.23	0.10	1.36%	0.12	1.70%	0.27	3.69%	0.31	4.28%			
A2+F II A	75	45.53	0.24	0.52%	0.24	0.53%	0.27	0.60%	0.44	0.96%			
A2+F II S	75	27.63	0.18	0.65%	0.21	0.77%	0.19	0.69%	0.34	1.22%			
A2 A2	75	1.63	0.03	2.12%	0.03	1.80%	0.05	2.92%	0.07	4.04%			
C-Trait C	75	30.51	0.26	0.85%	0.36	1.17%	0.86	2.82%	0.97	3.17%			
D-Trait D	75	35.03	0.20	0.57%	0.34	0.96%	0.58	1.64%	0.70	1.99%			
E-Trait E	75	23.08	0.34	1.47%	0.31	1.36%	0.73	3.17%	0.87	3.75%			



# **Multisite Precision – High Resolution**

High Resolution										
			Repea	tability	Betwe	en-Day	Betwe	en-Site	Reprod	ucibility
Sample Description	N	Mean Value	SD	%CV	SD	%CV	SD	%CV	SD	%CV
FASC F	75	17.31	0.11	0.64%	0.06	0.35%	1.43	8.28%	1.44	8.31%
FASC A	75	43.97	0.26	0.60%	0.04	0.08%	0.30	0.68%	0.40	0.91%
FASC S	75	13.46	0.10	0.72%	0.06	0.45%	0.07	0.50%	0.13	0.98%
FASC C	75	14.30	0.17	1.17%	0.07	0.47%	0.15	1.05%	0.24	1.65%
A2+F1F	75	2.89	0.07	2.47%	0.04	1.30%	0.07	2.30%	0.10	3.62%
A2+FIA	75	85.21	0.32	0.37%	0.97	1.14%	0.62	0.73%	1.20	1.40%
A2+F I A2	75	3.03	0.06	1.90%	0.05	1.52%	0.02	0.50%	0.08	2.49%
A2+F II F	75	7.29	0.08	1.11%	0.09	1.27%	0.29	3.93%	0.31	4.27%
A2+F II A	75	46.13	0.31	0.67%	0.50	1.09%	0.30	0.65%	0.66	1.44%
A2+F II S	75	28.99	0.19	0.66%	0.32	1.09%	0.37	1.28%	0.52	1.81%
A2 A2	75	1.59	0.05	3.33%	0.03	1.73%	0.03	1.85%	0.07	4.18%
C-Trait C	75	30.53	0.22	0.71%	0.26	0.86%	0.31	1.03%	0.46	1.51%
D-Trait D	75	36.65	0.38	1.03%	0.22	0.59%	0.39	1.07%	0.59	1.60%
E-Trait E	75	22.25	0.27	1.23%	0.30	1.35%	0.46	2.06%	0.61	2.76%

### **Limits of Detection**

The Limit of Detection Study was conducted in accordance with *CLSI EP17-A2 Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition.* This study determined the limit of quantitation, blank, and lower limits of detection, for hemoglobins A, A2, F, S, C, D-Los Angeles, and E on the Trinity Biotech Premier Resolution. For the Limit of Detection evaluation, human whole blood samples containing 4 different levels of hemoglobins A, A2, F, S, C, D-Los Angeles and E were tested on the Premier Resolution System in the Quick Scan and High Resolution modes.

# **Premier Resolution Quick Scan Assay**:

The limit of detection (LoD) for HbF is 0.2% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbF is 1.1% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbA is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.



# Premier Resolution System Premarket Notification

The limit of quantitation (LoQ) for HbA is 2.3% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbA2 is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbA2 is 1.5% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbS is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbS is 1.0% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbC is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbC is 1.0% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbD-Los Angeles is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbD-Los Angeles is 1.5% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbE is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbE is 1.5% of the total area based on 60 determinations of low-level samples.

### **Premier Resolution High-Resolution Assay:**

The limit of detection (LoD) for HbF is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbF is 1.1% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbA is 0.7% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbA is 2.2% of the total area based on 60 determinations of low-level samples.



# Premier Resolution System Premarket Notification

The limit of detection (LoD) for HbA2 is 0.2% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbA2 is 1.5% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbS is 0.3% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbS is 0.9% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbC is 0.3% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbC is 1.7% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbD-Los Angeles is 0.1% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbD-Los Angeles is 1.4% of the total area based on 60 determinations of low-level samples.

The limit of detection (LoD) for HbE is 0.6% of the total area based on 60 determinations of low-level samples where the limit of blank (LoB) is set at 0.0% of the total area.

The limit of quantitation (LoQ) for HbE is 2.7% of the total area based on 60 determinations of low-level samples.

# **Premier Resolution Analytical Specificity**

In studies on the analytical specificity of the quantitation of common hemoglobins by the Premier Resolution Analyzer, the Quick Scan assay and High Resolution assay were evaluated.

The analyte hemoglobins F, A, A2, S, C, D – Los Angeles, and E showed no statistical significance of interference with 20mg/dL of Bilirubin, 90mg/dL of acetylsalicylic, 85 UI concentration of Sodium Heparin anticoagulant,

85 UI concentration of Lithium Heparin anticoagulant, 9 mg/mL concentration of K2EDTA anticoagulant, and 9 mg/mL concentration of K3EDTA anticoagulant when analyzed in Quick Scan and High Resolution assay using Premier Resolution System.

Therefore, there is no interference with the quantitation of hemoglobins; HbF, Hb A, Hb A2, Hb S, Hb C, Hb D-Los Angeles and Hb E from the following interferent and concentrations stated in the table below:

# **Premier Resolution System**

# Premarket Notification

Interferent	Concentration
D-Glucose (HbF Only)	5000mg/dl
Acetaldehyde	20mg/dL
Lipemia (Triglycerides)	4500mg/dL
Icterus (Unconjugated Bilirubin)	20mg/dL
Icterus (Conjugated Bilirubin)	10mg/dL
Acetylsalicylic Acid	90mg/dL
Sodium Heparin	85 UI
Lithium Heparin	85 UI
K2EDTA	9 mg/mL
K3EDTA	9 mg/mL

# 8. Conclusions

Performance testing has demonstrated that the similarities of the assays, the use of the same technology, ion exchange HPLC, and the excellent concordance between the two devices, it can be concluded that the Premier Resolution System is substantially equivalent to the predicate device, the Bio-Rad VARIANT II  $\beta$  Thalassemia. Based on the establishment of substantial equivalence, the safety and effectiveness of the Premier Resolution System has been confirmed.