

April 11, 2023

Shenzhen HugeMed Medical Technical Development Co., Ltd % Kevin Wang
Consultant
Shenzhen Chonconn Medical Device Consulting Co., Ltd.
Room 504, Block C, No. 1029 Nanhai Avenue, Nanshan District Shenzhen, Guangdong 518067
China

Re: K222910

Trade/Device Name: Bronchoscope System Regulation Number: 21 CFR 874.4680

Regulation Name: Bronchoscope (Flexible Or Rigid) And Accessories

Regulatory Class: Class II Product Code: EOQ

Dated: March 7, 2023 Received: March 8, 2023

Dear Kevin Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Joyce C. Lin -S

for Shu-Chen Peng, Ph.D.
Assistant Director
DHT1B: Division of Dental and
ENT Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K222910	
Device Name Bronchoscope System	
Indications for Use (<i>Describe</i>) The Single-use Bronchoscope have been designed to be used wancillary equipment for endoscopy within the airways and track	heobronchial tree.
The Bronchoscope System is for use in a hospital environment. designed for use in adults.	. The Single-use Bronchoscope is a single-use device
Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARA	ATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) Summary

Prepared in accordance with the requirements of 21 CFR Part 807.92

Prepared Date: 2023/04/11Submission sponsor

Name: Shenzhen HugeMed Medical Technical Development Co., Ltd.

Address: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan Community, Yuanshan

Street, Longgang District, Shenzhen, Guangdong Province, 518115, China

Contact person: Cathy
Title: Regulatory engineer
E-mail: cathy.shi@hugemed.cn

Tel: +86-755-22275833

2. Submission correspondent

Name: Shenzhen Chonconn Medical Device Consulting Co., Ltd.

Address: Room 504, Block C, No. 1029 Nanhai Avenue, Nanshan District, Shenzhen,

Guangdong, P.R. China Contact person: Kevin Wang E-mail: kevin@chonconn.com

Tel: +86-755 33941160

3. Subject Device Information

Trade/Device Name	Bronchoscope System		
	Single-use Bronchoscope:BR-M58, BR-M52, BR-M50, BR-		
Model	M40, BR-M32, BR-M22		
	Image Processor: VLM-02		
Common Name	Bronchoscope System		
Regulatory Class	Class II		
Classification	21CFR 874.4680 / Bronchoscope (Flexible or Rigid) / EOQ		
Submission type	Traditional 510(K)		

4. Predicate Device

510(k) number: K173727

Product name: Ambu® aScopeTM 3, Ambu® aViewTM Monitor

Submitter: Ambu Inc.

5. Reference Device

510(k) number: K213782

Product name: Video Bronchoscope System Submitter: Micro-Tech (Nanjing) Co., Lts

6. Device Description

The Bronchoscope System consists of Single-use Bronchoscope (six models shown in below) to be introduced within the airways or tracheobronchial tree and Image Processor (model: VLM-02) for clinical image processing. The Flexible Bronchoscope is inserted through the airways and tracheobronchial tree during Bronchoscopy. The Image Processor provides power and processes the images for medical electronic endoscope.

System name	Component name	Model
Bronchoscope	Single-use Bronchoscope	BR-M58, BR-M52, BR-M50, BR-M40,
System		BR-M32, BR-M22
	Image Processor	VLM-02

The Single-use Bronchoscope is a sterile single used flexible bronchoscope. The Image Processor is a reusable monitor.

The light emitted by the LED cold light source at the distal tip of the Single-use Bronchoscope is irradiated into the body cavity, and the light reflected from the cavity enters the optical system and is captured by the CMOS image sensor. The CMOS acquisition image is controlled by the CMOS drive circuit, and the RGB video signal is output to the Image Processor via the VI circuit. The Image Processor receives video signals from the endoscope, processes the video signals, and outputs the processed video signal to the monitor. The Image Processor also controls the brightness of the LEDs on the endoscope.

The optical components and their arrangement at the distal tip for all models of the Singleuse Bronchoscope are identical.

Single-use Bronchoscope has the following physical and performance characteristics:

- Maneuverable tip controlled by the user
- Flexible insertion cord
- Camera and LED light source at the distal tip
- Sterilized by Ethylene Oxide
- For single use

Image Processor has the following physical and performance characteristics:

- Display the image from the Single-use Bronchoscope on the screen
- Can record screenshots or video of image from the Single-use Bronchoscope
- Can connect to an external monitor
- Reusable device

7. Intended use & Indication for use

The Single-use Bronchoscope have been designed to be used with the Image Processor, endotherapy accessories and other ancillary equipment for endoscopy within the airways

and tracheobronchial tree.

The Bronchoscope System is for use in a hospital environment. The Single-use Bronchoscope is a single-use device designed for use in adults.

8. Comparison to the Predicate Device

Features	Subject Device	Predicate Device	Comparison
K number	K222910	K173727	/
Manufacturer	Shenzhen HugeMed Medical Technical Development Co., Ltd.	Ambu Inc.	/
Model	Single-use Bronchoscope:BR-M58, BR-M52, BR-M50, BR- M40, BR-M32, BR-M22 Image Processor: VLM- 02	Ambu ® aScope TM 3 Broncho System: Ambu® aScope TM 3 Broncho Slim 3.8/1.2 Ambu® aScope TM 3 Broncho Regular 5.0/2.2 Ambu® aScope TM 3 Broncho Large 5.8/2.8 Ambu® aView TM Monitor	
Classification Name	Bronchoscope (flexible or rigid) and accessories	Bronchoscope (flexible or rigid) and accessories	/
Device trade name	Bronchoscope System	Ambu® aScope TM 3 Broncho Slim 3.8/1.2; Ambu® aScope TM 3 Broncho Regular 5.0/2.2; Ambu® aScope TM 3 Broncho Large 5.8/2.8; Ambu® aView Monitor	/
Product Code	EOQ	EOQ	Same
Indication for use	The Single-use Bronchoscope have been designed to be used with the Image Processor,	The aScope 3 Broncho endoscopes have been designed to be used with the aView monitor,	Same

Features	Subject Device	Predicate Device	Comparison
	endotherapy accessories	endotherapy accessories	
	and other ancillary	and other ancillary	
	equipment for endoscopy	equipment for	
	within the airways and	endoscopy within the	
	tracheobronchial tree.	airways and	
	The Bronchoscope	tracheobronchial tree.	
	System is for use in a	The aScope 3 system	
	hospital environment. The	Broncho system are for	
	Single-use Bronchoscope	use in a hospital	
	is a single-use device	environment.	
	designed for use in adults.	The aScope 3 Broncho	
		are single-use devices	
		designed for use in	
		adults.	
Population	Adults	Adults	Same
Anatomic sites	airways and	airways and	Same
	tracheobronchial tree	tracheobronchial tree	
Rx only	Yes	Yes	Same
Technology	The Flexible	The Flexible	Same
	bronchoscope is inserted	bronchoscope is	
	through the airways and	inserted through the	
	tracheobronchial tree	airways and	
	during Bronchoscopy.	tracheobronchial tree	
	Anatomical images are	during Bronchoscopy.	
	transmitted to the user by	Anatomical images are	
	the video processor with a	transmitted to the user	
	CMOS chip at the distal	by the video processor	
	end of the endoscope and	with a CMOS chip at	
	the images showing on a	the distal end of the	
	monitor.	endoscope and the	
		images showing on a	
		monitor.	
Connect to	Monitor	Monitor	Same
devices			
Performance	Complies with: ISO 8600	Complies with: ISO	Same
		8600	
Field of view	120°±10%	85°	Different
(degree)			
Direction of	0°±10°	0°	Same
view (degree)			

Features	Subject Device	Predicate Device	Comparison
Depth of view	3-50mm	8-19 mm	Different
Working length	600mm±10%	600mm	Same
(mm)			
Maximum	BR-M22: 2.2±10%mm	aScope TM 3 Slim	Different
insertion portion	BR-M32: 3.2±10%mm	4.3mm	
width	BR-M40: 4.4±10%mm	aScope™ 3 Regular	
	BR-M50: 5.0±10%mm	5.5mm	
	BR-M52: 5.2±10%mm	aScope TM 3 Large	
	BR-M58: 5.8±10%mm	6.3mm	
Minimum	BR-M22: not applicable	aScope TM 3 Slim	Different
instrument	BR-M32: 1.2mm	1.2mm	
channel width	BR-M40: ≥1.1mm	aScope TM 3 Regular	
	BR-M50: 2.8mm	2.2mm	
	BR-M52: ≥2.1mm	aScope TM 3 Large	
	BR-M58: ≥2.5mm	2.8mm	
Luer/Luer lock	No	Yes	Same as
connection to			K213782
working channel			reference
			device
			without the
			luer lock
			connection
Deflection angle	180° upward and 180°	aScope TM 3 Slim:	Same as
	downward	130°up 130°down	K213782
		aScope TM 3 Regular:	reference
		150°up 130°down	device with
		aScope TM 3 Large:	the 180°
		140°up 110°down	deflection
			angle
Digital video	CMOS	CMOS	Same
technology			
Illumination	LED	LED	Same
source			
Separate	Yes	Yes	Same
monitor			
Energy	Yes	Yes	Same
used/power			
source			
Image/Video	Yes	Yes	Same
capture			

Features	Subject Device	Predicate Device	Comparison
Storage	Yes	Yes	Same
	SD Card	SD Card	
Disposable after	Single-use Bronchoscope:	aScope TM 3: Yes	Same
use	Yes		
	Image Processor:	aView TM : Reusable	Same
	Reusable		
Components in	Yes, EO	Yes, EO	Same
contact with the			
patient is			
delivered sterile			
Suction possible	Yes	Yes	Same
Biocompatibility	Comply with ISO 10993-	Comply with ISO	Same
	1	10993-1	
Shelf life	Single-use Bronchoscope:	aScope TM 3: 3 Years	Same
	3 Years	Comply with ASTM	
	Comply with ASTM	F1980-16	
	F1980-16		
Electrical	Comply with	Comply with	Same
Performance	ANSI AAMI ES60601-	ANSI AAMI ES60601-	
	1:2005/(R)2012 and	1:2005/(R)2012 and	
	A1:2012,	A1:2012,	
	C1:2009(R)2012 and A2:	C1:2009(R)2012 and	
	2010/(R)2012	A2: 2010/(R)2012	
	IEC 60601-1-2: 2014	IEC 60601-1-2: 2014	
	IEC 60601-2-18:2009	IEC 60601-2-18:2009	

9. Performance Data

The following performance data were provided in support of the substantial equivalence determination.

Biocompatibility testing

Biocompatibility of the Single-use Bronchoscope was evaluated in accordance with the FDA guidance "Use of International Standard ISO 10993-1, "Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process" and International Standard ISO 10993-1 "Biological Evaluation of Medical Devices — Part 1: Evaluation and Testing Within a Risk Management Process," as recognized by FDA. The following tests were performed, as recommended:

- Cytotoxicity
- Sensitization

• Intradermal reactivity

The Single-use Bronchoscope is considered surface – mucosal membrane contacting for a duration of less than 24 hours.

Sterilization and shelf life testing

The Single-use Bronchoscope is provided sterile and its shelf-life is 3 years. Sterilization Process has been validated in accordance with ISO 11135:2014.

EO/ECH residual test was performed according to ISO 10993-7:2008.

The shelf life is determined based on optical testing and product performance testing after accelerated aging test according to ASTM F1980-16 Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices.

Package validation was conducted according to ISO 11607-1:2019 and ISO 11607-2:2019, and ASTM F1886/F1886M-16, ASTM F88/F88M-15, ASTM F 1929-15.

Electrical safety and electromagnetic compatibility (EMC)

Electrical safety and EMC testing were conducted on the Bronchoscope System. The system complies with the IEC 60601-1 and IEC60601-2-18 for safety and the IEC 60601-1-2 for EMC.

Software Verification and Validation Testing

Software verification and validation testing were conducted, and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The software for this device was considered as a "moderate" level of concern.

Bench performance testing

The following bench tests were performed:

- 1. Optical performance testing according to ISO 8600 series.
- 2. Color performance (color reproduction), geometric distortion, optical performance (resolution, depth of field and image intensity uniformity), SNR and dynamic range, image frame frequency and system delay test compared with the predicate device.
- 3. Mechanical testings including suction and bending testing.

10. Clinical study

No clinical study is included in this submission

11. Conclusion

Performance testing and compliance with voluntary standards demonstrate that the proposed subject device is substantially equivalent to the predicate device.