

March 15, 2023

Monteris Medical David Mueller Senior Principal Regulatory Affairs Specialist 131 Cheshire Lane; Suite 100 Minnetonka, Minnesota 55305

Re: K222983

Trade/Device Name: NeuroBlate System Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In

Dermatology

Regulatory Class: Class II

Product Code: GEX, HAW, ONO

Dated: February 10, 2023 Received: February 10, 2023

#### Dear David Mueller:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

K222983 - David Mueller Page 2

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reportingcombination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reportingmdr-how-report-medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medicaldevices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatoryassistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Adam D. Digitally signed by Adam D. Pierce -S Date: 2023.03.15

Adam D. Pierce, Ph.D. Assistant Director DHT5A: Division of Neurosurgical, Neurointerventional and Neurodiagnostic Devices OHT5: Office of Neurological and Physical Medicine Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

#### **Indications for Use**

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number <i>(if known)</i>
K222983
Device Name
Monteris Medical NeuroBlate System,
NeuroBlate Fusion-S Software V3.17
Indications for Use (Describe)
The Monteris Medical NeuroBlate System is a neurosurgical tool and is indicated for use to ablate, necrotize, or coagulate intracranial soft tissue, including brain structures (e.g., brain tumor, radiation necrosis, and epileptogenic foci as identified by non-invasive and invasive neurodiagnostic testing, including imaging), through interstitial irradiation or thermal therapy in medicine and surgery in the discipline of neurosurgery with 1064 nm lasers.
The Monteris Medical NeuroBlate System is intended for planning and monitoring thermal therapies under MRI visualization. It provides MRI-based trajectory planning assistance for the stereotaxic placement of MRI compatible (conditional) NeuroBlate Laser Delivery Probes. It also provides near real-time thermographic analysis of selected MRI images.
When interpreted by a trained physician, this System provides information that may be useful in the determination or assessment of thermal therapy. Patient management decisions should not be made solely on the basis of the NeuroBlate System analysis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

## **Section 3: 510(k) Summary - K222983**

## 5a. Device Information:

Category	Comments	
Sponsor:	Monteris Medical	
	131 Cheshire Lane, suite 100	
	Minnetonka, MN 55305	
	866-799-7655	
	www.monteris.com	
Correspondent Contact	David H. Mueller	
Information:	Senior Principal Regulatory Affairs Specialist	
	Monteris Medical,	
	TEL: 763-333-1614	
	Email: DMueller@Monteris.com	
Device Common Name:	Magnetic Resonance Image Guided Laser Thermal	
	Therapy System	
Device Regulation and Name:	21 CFR 878.4810	
	Laser surgical instrument for use in general and	
	plastic surgery and in dermatology	
	21 CFR 882.4560	
	Stereotaxic instrument	
Device Classification &	Class II, GEX, ONO	
Product Code:	Class II, HAW	
Device Proprietary Name:	Monteris Medical NeuroBlate® System	

## **Predicate Device Information:**

Predicate Device:	NeuroBlate System
Predicate Device Manufacturer:	Monteris Medical
Predicate Device Common Name:	Monteris NeuroBlate System
Predicate Device Premarket Notification #	K193375, K201056
Predicate Device Regulation:	21 CFR 878.4810
	Laser surgical instrument for use in general
	and plastic surgery and in dermatology
	21 CFR 882.4560
	Stereotaxic instrument
Predicate Device Classification &	Class II, GEX
Product Code:	Class II, HAW

#### 5b. Date Summary Prepared

February 10, 2023

#### 5c. Description of Device

The Monteris NeuroBlate<sup>®</sup> System is a collection of MRI-compatible laser devices and accessories that create an MRI guided intracranial delivery of precision thermal therapy in the practice of neurosurgery.

As previously described in K201056, the NeuroBlate System is typically used for the minimally invasive ablation of target tissue (tumors, radiation necrosis, epileptic foci) in the brain.

The NeuroBlate System components consist of:

- Families of gas-cooled Laser Delivery Probe (Probe) (SideFire & FullFire) to deliver controlled energy to a target zone.
- Probe Drivers (Advanced Probe Driver, Robotic Probe Driver) which allow the surgeon to precisely position, stabilize and manipulate a probe, endoscope or other device within the target zone.
- An Interface Platform, which attaches to the MRI system patient table and provides supporting electronics for the Probe Drivers and interconnections for the Laser Delivery Probes (e.g., Connector Module).
- A System Electronics Rack and Components, which includes the laser and necessary umbilicals, cables, penetration panels, and small hardware for system mechanical, electrical, and electronic operation.
- A Control Workstation including the *M-Vision*<sup>TM</sup>, *M-Vision Pro*<sup>TM</sup>, M-Vision *Fusion*<sup>TM</sup>, or *Fusion-S*<sup>TM</sup> software, which includes a user interface for procedure planning, interactive monitoring of NeuroBlate procedures, and interfaces to the MRI and hardware subsystems.

The NeuroBlate System is utilized with stereotaxic frames and patient stabilization systems, such as:

- The Monteris Cranial Bolt and Mini-Bolt fixation components, and
- The AtamA Stabilization System, MRI head coils, and other optional accessories, including: drill bits, bolts, thumbscrews, instrument adaptors, accessory host adaptors, MRI trajectory wands, cranial screws, fiducial markers, bone screws, stereotactic manual driver with mandrel and T-handle, and other manual accessory instruments and tools.

There is no change to entire system, with the exception of an additional NeuroBlate Fusion-S<sup>TM</sup> Software Package (V3.17). While the modified software package includes several modifications which do not meet the FDA's "significant change" submission criteria, there is at least one intended change which modifies an existing risk control measure for a hazardous situation and thus, per FDA 510(k) guidance, requires a new 510(k) submission.

#### 5d. Indications for Use

The Monteris Medical NeuroBlate System is a neurosurgical tool and is indicated for use to ablate, necrotize, or coagulate intracranial soft tissue, including brain structures (e.g., brain tumor, radiation necrosis, and epileptogenic foci as identified by non-invasive and invasive neurodiagnostic testing, including imaging), through interstitial irradiation or thermal therapy in medicine and surgery in the discipline of neurosurgery with 1064 nm lasers.

The Monteris Medical NeuroBlate System is intended for planning and monitoring thermal therapies under MRI visualization. It provides MRI-based trajectory planning assistance for the stereotaxic placement of MRI compatible (conditional) NeuroBlate Laser Delivery Probes. It also provides near real-time thermographic analysis of selected MRI images.

When interpreted by a trained physician, this System provides information that may be useful in the determination or assessment of thermal therapy. Patient management decisions should not be made solely on the basis of the NeuroBlate System analysis.

#### 5e. Comparison to Predicate Device

There is no change to entire system, with the exception of an additional (alternative) NeuroBlate Fusion-S Software Package (V3.17). While the modified software package includes several modifications which do not meet the FDA's "significant change" submission criteria, there is at least one intended change which modifies an existing risk control measure for a hazardous situation and thus, per FDA 510(k) guidance, requires a new 510(k) submission

The substantial equivalency comparisons for the NeuroBlate System's Technical Characteristics are described in Table 1.

**Table 1: Substantial Equivalence: NeuroBlate System Comparison Table** 

Table 1. Substantiar	lai Equivalence: NeuroBiate System Comparison Table		
Cl4! -4! -	Predicate:	Intended	
Characteristic	Monteris Medical	Alternative Modified Fusion-S	
	NeuroBlate System	Software (V3.17)	
	(K193375, K201056)	Software (13.17)	
Company	Monteris Medical	Identical	
Regulation Number	21 CFR 878.4810		
	Laser surgical instrument for use in general		
	and plastic surgery and in dermatology	Identical	
	21 CFR 882.4560		
	Stereotaxic instrument		
FDA Product Code	GEX & HAW	Identical	
12111100000000	ODE WINTY	The Monteris Medical NeuroBlate	
Indications for Use	The Monteris Medical NeuroBlate®	System is a neurosurgical tool and is	
1110101101101101101010	System is a surgical tool for brain surgery	indicated for use to ablate, necrotize,	
	(e.g., tumor and epilepsy surgery) and is	or coagulate intracranial soft tissue,	
	indicated for use to ablate, necrotize, or	including brain structures (e.g., brain	
	coagulate intracranial soft tissue,	tumor, radiation necrosis, and	
	including brain structures, through	epileptogenic foci as identified by non-	
	interstitial irradiation or thermal therapy	invasive and invasive neurodiagnostic	
	in medicine and surgery in the discipline of	testing, including imaging), through	
	neurosurgery with 1064 nm lasers.	interstitial irradiation or thermal	
	The Monteris Medical NeuroBlate System is intended for planning and monitoring	therapy in medicine and surgery in the discipline of neurosurgery with 1064 nm lasers.	
	thermal therapies under MRI visualization. It provides MRI-based trajectory planning assistance for the stereotaxic placement of MRI compatible (conditional) NeuroBlate Laser Delivery Probes. It also provides near real-time thermographic analysis of selected MRI images.	The Monteris Medical NeuroBlate System is intended for planning and monitoring thermal therapies under MRI visualization. It provides MRI-based trajectory planning assistance for the stereotaxic placement of MRI compatible (conditional) NeuroBlate	
	When interpreted by a trained physician, this System provides information that may be useful in the determination or assessment of thermal therapy. Patient management decisions should not be made solely on the basis of the NeuroBlate System analysis.	Laser Delivery Probes. It also provides near real-time thermographic analysis of selected MRI images. When interpreted by a trained physician, this System provides information that may be useful in the determination or assessment of thermatherapy. Patient management decisions should not be made solely on the basis of the NeuroBlate System analysis.	
Operating Principle (Technology)	MRI-guided Laser Probe with near-real time thermographic monitoring of target tissue temperature	Identical	
Mode of Action	Laser Interstitial Thermal Therapy (LITT)	Identical	

Characteristic	Predicate: Monteris Medical NeuroBlate System (K193375, K201056)	Intended Alternative Modified Fusion-S Software (V3.17)
Design	<ul> <li>Families of gas-cooled Laser Delivery Probe (Probe) (SideFire &amp; FullFire) to deliver controlled energy to a target zone.</li> <li>Probe Drivers (Advanced Probe Driver, Robotic Probe Driver) which allow the surgeon to precisely position, stabilize and manipulate a probe, endoscope or other device within the target zone.</li> <li>An Interface Platform, which attaches to the MRI system patient table and provides supporting electronics for the Probe Drivers and interconnections for the Laser Delivery Probes, including the supporting cables and interconnections.</li> <li>A System Electronics Rack and Components, which includes necessary umbilicals, cables, penetration panels, and small hardware for system mechanical, electrical, and electronic operation, and</li> <li>A Control Workstation including the M-Vision™, M-Vision Pro, Fusion-S software, which includes a user interface for procedure planning, interactive monitoring of NeuroBlate procedures, and interfaces to the MRI and hardware subsystems.</li> <li>The Monteris Cranial Bolt and MiniBolt fixation components, and</li> <li>The AtamA Stabilization System, MRI head coils, as well as, other optional accessories, including: drill bits, bolts, thumbscrews, instrument adaptors, accessory host adaptors, MRI trajectory wands, cranial screws, bone screws, stereotactic manual driver with mandrel and T-handle, and other manual accessory instruments and tools.</li> </ul>	Identical NeuroBlate System hardware Software is changing as described along with the corresponding labelling updates, i.e., utilizing Fusion-S V3.17 Software
NeuroBlate Software	NeuroBlate M-Vision and Fusion-S Software Version V3.15 and V3.16	NeuroBlate Fusion-S Software Version V3.17
Probe Distal diameter	2.2mm and 3.3 mm	Identical

Characteristic	Predicate: Monteris Medical NeuroBlate System (K193375, K201056)	Intended Alternative Modified Fusion-S Software (V3.17)
3.3mm Probe Size as per Probe Body Length Range	Size #000 Size #00 Size #1 Size #2 Size #3 Size #4 Size #5	Identical
2.2mm Probe Size as per Probe Body Length Range	Size #000 Size #00 Size #1 Size #2 Size #3 Size #4 Size #5	Identical
Probe Diameters and Laser Output Direction	<ul> <li>2.2mm NeuroBlate FullFire (DTP)         Probes     </li> <li>3.3mm NeuroBlate FullFire (DTP)          Probes     </li> <li>3.3mm NeuroBlate SideFire (SFP)         Probes     </li> </ul>	Identical
Energy Source Warnings and	<ul> <li>Donier Medilase D Laser, or</li> <li>NeuroLase Laser</li> <li>Warnings and Cautions included in</li> </ul>	Identical Warnings and Caution statements
Cautions	labeling and IFU content	updated

#### 5f. Summary of Supporting Data

The updated Fusion-S (V3.17) Software development process followed Monteris' documented Quality System and incorporated a design verification and design validation process. This process included an overarching Design Verification and Design Validation Master Plan. This plan describes the design verification and the design validation of the user needs for the Fusion-S (V3.17) Software when used within the NeuroBlate System.

The Design Verification process utilized protocols to detail the associated tests. Each verification test protocol incorporated clearly defined acceptance criteria. The corresponding test reports confirmed (and documented) the design output met the design input for the requirements.

The Design Validation process utilized protocols to detail the associated tests. Each validation protocol described the objective, test method and acceptance criteria. The corresponding test reports confirmed (and documented) the modified NeuroBlate System met the user needs and intended use.

### 5g. Discussion of Performance Testing

The performance of the V3.17 software was tested utilizing automated and manual test methods and was found to be comparable to that of predicate V3.16 software. All tested metrics were measured to be within acceptable limits and improvements were noted in several metrics.

#### 5h. Conclusion

Monteris concludes that the application for the Monteris Medical NeuroBlate System with the alternative V3.17 Software is substantially equivalent to the predicate Monteris NeuroBlate System (K193375, K201056) in intended use, technology, design and physician use.