

March 20, 2023

Integrity Spine % Jennifer Palinchik President JALEX Medical 27865 Clemens Rd, Suite 3 Westlake, Ohio 44145

Re: K223043

Trade/Device Name: The Integrity Spine Core System

Regulation Number: 21 CFR 888.3080

Regulation Name: Intervertebral Body Fusion Device

Regulatory Class: Class II

Product Code: ODP

Dated: February 22, 2023 Received: February 23, 2023

Dear Jennifer Palinchik:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Brent Showalter -S

Brent Showalter, Ph.D.
Assistant Director
DHT6B: Division of Spinal Devices
OHT6: Office of Orthopedic Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K223043						
Device Name The Integrity Spine Core System						
Indications for Use (Describe) When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 disc. DDD is defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies. The device system is designed for use with supplemental fixation and with autograft to facilitate fusion. Patients should have at least six (6) weeks of non-operative treatment prior to treatment with an intervertebral cage.						
Type of Use (Select one or both, as applicable)						
CONTINUE ON A SEPARATE PAGE IF NEEDED.						
This section applies only to requirements of the Paperwork Reduction Act of 1995.						

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary

Submitted By: Integrity Spine

414 W Sunset Rd Ste 205

Dallas, TX 75206

Date: 03/16/2023

Contact Person: Jennifer Palinchik, President

JALEX Medical

Contact Telephone: (440) 541-0060 **Contact Fax:** (440) 933-7839

Device Trade Name: Integrity Spine Core System

Device Classification Name: Intervertebral Body Fusion Device with Bone Graft, Cervical

Device Classification: Class II **Reviewing Panel:** Orthopedic **Product Code:** ODP

Trouder Code. ODI

Predicate Device: K132718 Integrity Spine Core System

Additional Predicate: K133967 Aurora Spine Interbody Fusion System

The predicate devices have not been subject to any design related

recalls.

Device Description:

The Integrity Spine Core System is a cervical interbody fusion system comprised of lordotic cages in two footprints with varying heights designed to accommodate patient anatomy and may be implanted as a single device via an anterior approach.

The Integrity Spine Core System implant components are made of polyether ether ketone (PEEK, Zeniva ZA-500) that conforms to ASTM F2026. The devices are plasma spray coated with titanium that conforms to ASTM F1580. Additionally, the devices contain tantalum markers that conform to ASTM F560 to assist the surgeon with proper placement of the device.

The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899.

Intended Use:

When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 disc. DDD is defined as discogenic pain with degeneration of the disc confirmed by history and radiographic studies. The device system is designed for use with supplemental fixation and with autograft to facilitate fusion. Patients should have at least six (6) weeks of non-operative treatment prior to treatment with an intervertebral cage.



Summary of Technological Characteristics:

The Integrity Spine Core System and the primary predicate have the same intended use and fundamental scientific technology. All devices compare similarly in:

- Design features
- Intended use
- Materials
- Dimensions
- Function

Feature	ological Characteristics Comparison Integrity Spine				
reature	Core System	System (K132718)	Interbody Fusion	Comparison	
	(Subject)	System (K132718)	System (K133967)		
Classification	Intervertebral Body	Intervertebral Body	Intervertebral Body	Equivalent	
Name	Fusion Device with	Fusion Device with	Fusion Device with	Equivalent	
Name					
	Integrated Fixation,	Integrated Fixation,	Integrated Fixation,		
D 1.4	Cervical	Cervical	Cervical	D 1 1	
Regulation	888.3080	888.3080	888.3080	Equivalent	
Product Code	ODP	ODP	ODP	Equivalent	
Device	The Integrity Spine	The Integrity Spine	The Aurora Spine	Equivalent	
Description	Core System is a	Core System is a	Interbody Fusion		
	cervical interbody	cervical interbody	System,		
	fusion system	fusion system	manufactured from		
	comprised of	comprised of parallel	PEEK-Optima®,		
	lordotic cages in	and lordotic cages in	consist of implants		
	two footprints with	two footprints with	available in various		
	varying heights	varying heights	foot prints, heights,		
	designed to	designed to	and lordotic		
	accommodate	accommodate patient	configurations with		
	patient anatomy	anatomy and may be	an open architecture		
	and may be	implanted as a single	to accept packing of		
	implanted as a	device via an anterior	autograft materials.		
	single device via an	approach.	The exterior of the		
	anterior approach.		device has "teeth" or		
		The Integrity Spine	other generally		
	The Integrity Spine	Core System implant	sharp engagement		
	Core System	components are made	members on the		
	implant	of polyether ether	superior and inferior		
	components are	ketone (PEEK, Zeniva	surfaces to help		
	made of polyether	ZA-500) that	prevent the device		
	ether ketone	conforms to ASTM	from migrating once		
	(PEEK, Zeniva	F2026. Additionally,	it is surgically		
	ZA-500) that	the devices contain	positioned. The		
	conforms to ASTM	tantalum markers that	device comes in a		
	F2026. The devices	conform to ASTM	PEEK or PEEK with		
	are plasma spray	F560 to assist the	a plasma-sprayed		



Integrity Spine Core System (Subject)	Integrity Spine Core System (K132718)	Aurora Spine Interbody Fusion System (K133967)	Comparison
coated with titanium that conforms to ASTM	surgeon with proper placement of the device.	commercially pure titanium coating on the superior and	
Additionally, the devices contain tantalum markers that conform to ASTM F560 to assist the surgeon with proper placement of the device.	The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899.		
The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899.			
When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 disc. DDD is defined as	When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 disc. DDD is defined as discogenic pain with degeneration of the disc confirmed by history and	The Aurora Spine Interbody Fusion System is indicated for use as an intervertebral body fusion device in skeletally mature patients with degenerative disc disease (defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies) at one level of the cervical spine with accompanying radicular symptoms. Patients should have six weeks of non-	Equivalent
	Core System (Subject) coated with titanium that conforms to ASTM F1580. Additionally, the devices contain tantalum markers that conform to ASTM F560 to assist the surgeon with proper placement of the device. The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1	Core System (Subject) coated with titanium that conforms to ASTM F1580. Additionally, the devices contain tantalum markers that conform to ASTM F560 to assist the surgeon with proper placement of the device. The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 surgeon with proper placement of the device. The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7- T1 disc. DDD is defined as discogenic pain with degeneration	Core System (Subject) coated with titanium that conforms to ASTM F1580. Additionally, the devices contain tantalum markers that conform to ASTM F560 to assist the surgeon with proper placement of the device. The Integrity Spine Core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 Interbody Fusion System (K133967) commercially pure titanium coating on the superior and inferior surfaces. The Integrity Spine combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. When used as a cervical intervertebral body fusion device, the Integrity Spine Core System is indicated for intervertebral body fusion in skeletally mature patients with degenerative disc disease (DDD) of the cervical spine with accompanying radicular symptoms at one disc level from the C2-C3 disc to the C7-T1 disc to the C3-T3 Interbody Fusion System (K133967) commercially pure titanium coating on the superior and inferior surfaces. The Integrity Spine core System is implanted using a combination of device specific and universal class I instruments manufactured from stainless steel that conforms to ASTM F899. The Integrity Spine Core System is indicated for intervertebral body fusion device in intervertebral body fusion devic



Feature	Integrity Spine	Integrity Spine Core	Aurora Spine	Comparison				
	Core System	System (K132718)	Interbody Fusion					
	(Subject)		System (K133967)					
	discogenic pain	radiographic studies.	operative treatment					
	with degeneration	The device system is	prior to surgery.					
	of the disc	designed for use with	Cervical implants					
	confirmed by	supplemental fixation	are used to facilitate					
	history and	and with autograft to	fusion in the cervical					
	radiographic	facilitate fusion.	spine (C2-T1) and					
	studies. The device	Patients should have at	are placed via an					
	system is designed	least six (6) weeks of	anterior approach					
	for use with	non-operative	using autogenous					
	supplemental	treatment prior to	bone. When used as					
	fixation and with	treatment with an	an interbody fusion					
	autograft to	intervertebral cage.	device,					
	facilitate fusion.		supplemental					
	Patients should		fixation must be					
	have at least six (6)		used.					
	weeks of non-							
	operative treatment							
	prior to treatment							
	with an							
	intervertebral cage.							
Interbody	5-11 mm	5-11 mm	5-12 mm	Equivalent				
Heights								
Interbody	14x11 mm, 17x13	14x11 mm, 17x13 mm	14x12 mm, 16x14	Equivalent				
Footprints	mm		mm					
Lordosis	6°	0° and 6°	5°	Equivalent				
Implant	PEEK Zeniva ZA-	PEEK Zeniva ZA-500	PEEK Optima Per	Equivalent				
Materials	500 Per ASTM	Per ASTM F2026,	ASTM F2026,					
	F2026, Titanium	Tantalum markers Per	Titanium coating per					
	coating Per ASTM	ASTM F560	ASTM F67,					
	F1580, Tantalum		Tantalum markers					
	markers Per ASTM		Per ASTM F560					
	F560							



Non-clinical Testing:

The following non-clinical tests were conducted:

- Static and dynamic compression testing, conducted in accordance with ASTM F2077-11
- Static and dynamic torsion testing, conducted in accordance with ASTM F2077-11
- Subsidence testing, conducted in accordance with ASTM F2267-04
- Expulsion testing, conducted in accordance with ASTM Draft Standard F-04.25.02.02

Conclusion:

Based on the indications for use, technological characteristics, and comparison with the predicate device, the subject device has demonstrated substantial equivalence.