

April 5, 2023

Advanced Medical Solutions Limited Kay McGrath Regulatory Affairs Specialist Premier Park, 33 Road One Winsford Industrial Estate Winsford, Cheshire CW7 3RT United Kingdom

Re: K223310

Trade/Device Name: Antimicrobial Silicone PHMB Foam Wound Dressing Regulatory Class: Unclassified Product Code: FRO Dated: March 6, 2023 Received: March 6, 2023

Dear Kay Mcgrath:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Julie A. Morabito -S

Julie Morabito, Ph.D. Assistant Director DHT4B: Division of Infection Control and Plastic Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K223310

Device Name

Antimicrobial Silicone PHMB Foam Wound Dressing

Indications for Use (Describe)

Antimicrobial Silicone PHMB Foam Wound Dressings are indicated for use in the management of

- post-surgical incisions,
- pressure sores,
- venous stasis ulcers,
- diabetic ulcers,
- donor sites,
- abrasions,
- lacerations,
- superficial and partial thickness burns,
- dermatologic disorders,
- other wounds inflicted by trauma and,
- as a secondary dressing or cover dressing for packed wounds.

Type of Use (Select one or both, as applicable)	
---	--

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Advanced Medical Solutions Limited Premier Park, 33 Road One, Winsford Industrial Estate, Winsford, Cheshire, CW7 3RT, UK **Tel:** +44 (0) 1606 863500 **Fax:** +44 (0) 1606 863600 **Web:** www.admedsol.com Registered in England 2666957 VAT No. GB 636 5551 27

K223310

510(k) Summary

Submitted by:	Advanced Medical Solutions Ltd Premier Park 33 Road One Winsford Industrial Estate Winsford Cheshire CW7 3RT Tel: +44 1606 863500	
Contact Person:	Kay M ^c Grath	
Date of Summary:	April 04, 2023	
Trade Name:	Antimicrobial Silicone PHMB Foam Wound Dressing	
Common Name:	Wound Dressing	
Classification Name:	Dressing, Wound, Drug	
Classification:	Unclassified (Pre-amendment)	
Classification Code:	Product code: FRO	
Predicate Device(s):	Silicone PHMB Foam Wound Dressing (K190819)	



Page 1 of 5



K223310

DeviceAntimicrobial Silicone PHMB Foam Wound Dressing, is a polyurethane foam trilaminate dressingDescription:impregnated with Polyhexamethylene Biguanide (PHMB), an agent that protects the dressing from
bacterial penetration and colonization. The foam in the dressing has a microporous hydrophilic
foam structure that absorbs wound exudate and maintains a moist wound healing environment.

Based on *in vitro* performance data, the Antimicrobial Silicone PHMB Foam Wound Dressing provides a barrier to bacterial penetration through the dressing and the PHMB prevents colonization and proliferation of bacteria, yeast and mold within the dressing for up to 7 days. Antimicrobial Silicone PHMB Foam Wound Dressing, when tested *in-vitro* has demonstrated to be effective against gram positive bacteria, gram negative bacteria, yeast and mold challenge organisms within the dressing.

The perforated wound contact layer contains a gentle silicone adhesive that provides secure, nonirritating adhesion and supports non-traumatic removal during dressing changes.

The device is presented in a border (adhesive) version. The dressing is supplied sterile in a range of sizes between 10.24 in^2 (64cm) to 64 in² (400cm).

Indications forAntimicrobial Silicone PHMB Foam Wound Dressing is indicated for use in the management ofUse:post-surgical incisions, pressure sores, venous stasis ulcers, diabetic ulcers, donor sites, abrasions,
lacerations, superficial and partial thickness burns, dermatologic disorders, other wounds inflicted
by trauma and as a secondary dressing or cover dressing for packed wounds.



Page 2 of 5



K223310

Comparison of Technological Characteristics:

on of Antimicrobial Silicone PHMB Foam Wound Dressing is a modification of the predicate device to include updated product labelling in support of an antimicrobial product offering.

Antimicrobial Silicone PHMB Foam Wound Dressing has substantially equivalent intended use and performance characteristics, identical design, materials and manufacture process to the predicate device Silicone PHMB Foam Wound Dressing (K190819). Antimicrobial Silicone PHMB Foam Wound Dressing is a multi-layer one piece dressing design incorporating an absorbent polyurethane foam pad containing 0.8-1.1 %w/w PHMB. Based on *in vitro* testing, the PHMB within the dressing is efficacious against gram positive bacteria, gram negative bacteria, yeast and mold. The dressing is semi-occlusive allowing the exchange of gases within the dressing such as oxygen and moisture, and has a film that provides a barrier to bacterial penetration through the dressing. The silicone wound contact layer contains a gentle silicone adhesive that provides secure, non-irritating adhesion and supports non-traumatic removal during dressing changes. The following table shows the comparison of technological characteristics between the subject and predicate devices.

	Subject (modified device)	Predicate	Similarities and
Product Name	Antimicrobial Silicone PHMB Foam Wound Dressing	Silicone PHMB Foam Wound Dressing	Differences
Manufacturer	Advanced Medical Solutions Ltd	Advanced Medical Solutions Ltd	Identical
510(k)	K223310	K190819	-
Classification	Unclassified (pre-amendment)	Unclassified (pre-amendment)	Identical
Product code	FRO (Dressing, Wound, Drug)	FRO (Dressing, Wound, Drug)	Identical
Indications for use	Antimicrobial Silicone PHMB Foam Wound Dressings are indicated for use in the management of post-surgical incisions, pressure sores, venous stasis ulcers, diabetic ulcers, donor sites, abrasions, lacerations, superficial and partial thickness burns, dermatologic disorders, other wounds inflicted by trauma and, as a secondary dressing or cover	Silicone PHMB Foam Wound Dressings are indicated for use in the management of post-surgical incisions, pressure sores, venous stasis ulcers, diabetic ulcers, donor sites, abrasions, lacerations, 1st and 2nd degree burns, dermatologic disorders, other wounds inflicted by trauma and, as a secondary dressing or cover	Different - Subject device includes an antimicrobial claim





Advanced Medical Solutions Limited Premier Park, 33 Road One, Winsford Industrial Estate, Winsford, Cheshire, CW7 3RT, UK Tel: +44 (0) 1606 863500 Fax: +44 (0) 1606 863600 Web: www.admedsol.com Registered in England 2666957 VAT No. GB 636 5551 27

K223310

	dressing for packed wounds.	dressing for packed wounds.	
Primary Material	Polyurethane (PU) foam	Polyurethane (PU) foam	Identical
Antimicrobial agent	PHMB (0.8-1.1%w/w)	PHMB (0.8-1.1%w/w)	Identical
Sterilization method (terminal)	Ethylene oxide SAL 10 ⁻⁶	Ethylene oxide SAL 10 ⁻⁶	Identical
Biocompatibility	Biocompatible	Biocompatible	Identical
Storage	Store below 25°C (77°F)	Store below 25°C (77°F)	Identical
Range of available sizes	10.24 in ² (64cm) to 64 in ² (400cm)	10.24 in ² (64cm) to 64 in ² (400cm)	Identical







K223310

Performance Testing Summary:

Page 5 of 5

ance Microbial efficacy performance data submitted in support of this 510(k) includes *in-vitro* testing against a mold challenge organism. Testing was performed on real time aged predicate device in accordance with the well-established modified AATCC TM 100 method previously used for the predicate device.

No other performance tests were conducted for this submission. All performance data leveraged from the predicate device was submitted as part of the original 510(k) submission of the predicate, Silicone PHMB Foam Wound Dressing (K190819) and includes:

Biocompatibility

ISO 10993-1; Biological evaluation of medical devices USP 41-NF36; <151> Pyrogenic Test

Performance testing

BS EN 13726-1; Test methods for primary wound dressings – aspects of absorbency. BS EN 13726-2; Test methods for primary wound dressings – moisture vapour transmission rate of permeable film dressings. BS EN 13726-3; Test methods for primary wound dressings – waterproofness. ASTM D6282-11; Standard Test Method for 90 Degree Peel Resistance of

Adhesives.

Bacterial barrier.

Distribution

ASTM D4169 – Standard Practice for Performance Testing of Shipping Containers and Systems

The subject device, Antimicrobial Silicone PHMB Foam Wound Dressing, is manufactured with the exact same materials and processes as the predicate.

- Rationale for
SubstantialThe modified device, Antimicrobial Silicone PHMB Foam Wound Dressing, is identical to the
predicate, Silicone PHMB Foam Wound Dressing (K190819), with regard to technology, materials,
manufacture process, intended use, and target population. The only difference between the
predicate and subject device is that the subject device has an antimicrobial claim, this minor
modification does not raise any new questions of safety or effectiveness. Therefore, the
Antimicrobial Silicone PHMB Foam Wound Dressing, is substantially equivalent to the predicate,
Silicone PHMB Foam Wound Dressing (K190819).
- **Conclusion:** Antimicrobial Silicone PHMB Foam Wound Dressing is substantially equivalent to the predicate device listed when compared to the technological characteristics such as design, materials, chemical composition, and manufacture and are supplied sterile for single use. Based on the information provided within this 510(k), Advanced Medical Solutions Ltd. concludes that the proposed Antimicrobial Silicone PHMB Foam Wound Dressing is substantially equivalent to the predicate device listed, Silicone PHMB Foam Wound Dressing (K190819).



Certificate No. MD 78010