

August 30, 2023`

Ochsner Clinic Foundation Hakm Murad Bioengineer 1514 Jefferson Hwy New Orleans, Louisiana 70121

Re: K223367

Trade/Device Name: Ochsner Connected Inhaler Sensor

Regulation Number: 21 CFR 868.5630

Regulation Name: Nebulizer Regulatory Class: Class II

Product Code: CAF Dated: July 10, 2023 Received: July 25, 2023

Dear Hakm Murad:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

K223367 - Hakm Murad Page 2

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Ethan L. Nyberg -S

Ethan Nyberg, Ph.D.
Assistant Director
DHT1C: Division of Sleep Disordered
Breathing, Respiratory and
Anesthesia Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (If known)
K223367
Device Name
Ochsner Connected Inhaler Sensor
Indications for Use (Describe)
The Ochsner Connected Inhaler Sensor System includes the Ochsner Connected Inhaler Sensor. The Sensor is an accessory device intended for single-patient use to assist physicians and patients in recording and monitoring actuations of prescribed MDI usage.
The Ochsner Connected Inhaler Sensor Mobile Application records, stores, and transmits usage events from the Sensors to a remote storage system. With the Mobile Application, the user can review information collected from the Sensors and report and review symptoms and other information about their disease management and its impact. The user may also share their information with their caregivers, physicians, and health care providers.
When used with a prescribed MDI, the System can report on information captured during the normal course of use, such as the time between actuations, that can be helpful in assessing MDI technique.
When used under the care of a physician with a prescribed MDI, the System can assist in the management of respiratory health symptoms and exacerbations by providing feedback through reminders, notifications, and self-management education.
The System is intended to be used in populations from Child (>2 years) to adult.
The System can be used both indoors and outdoors; home, work, and clinical settings, as well as on aircraft.
The System may also be used in clinical trials where researchers need to know information about the use of MDI medications by a participant.
The output of the System is not intended to diagnose or replace a diagnosis provided by a licensed physician. The System is not intended for use as an MDI or inhaled medication dose counter, nor is it intended to indicate the quantity of medication remaining in an MDI or inhaled medication.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Subpart D) Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Prepared Date: 8/28/23

Submitter: Ochsner Clinic Foundation

1514 Jefferson Highway New Orleans, LA 70121

Official Contact: Hakm Murad, Ph.D.

Biomedical Engineering Mgr. Ochsner Clinic Foundation 1514 Jefferson Highway New Orleans, LA 70121

5048427785

Hakm.murad@ochsner.org

Trade Name: Ochsner Connected Inhaler Sensor

Common Name: Nebulizer Accessory

Classification Name: Nebulizer

Classification Regulation: 21 CFR §868.5630

Product Code: CAF

Device Description: Electronic MDI Accessory

Predicate Device: Propeller System Model 2 OTC K142516

Intended Use: The Ochsner Connected Inhaler Sensor System includes the Ochsner

Connected Inhaler Sensor. The Sensor is an accessory device intended for

single-patient use to assist physicians and patients in recording and

monitoring actuations of prescribed MDI usage.

The Ochsner Connected Inhaler Sensor Mobile Application records, stores, and transmits usage events from the Sensors to a remote storage system. With the Mobile Application, the user can review information collected from the Sensors and report and review symptoms and other information about their disease management and its impact. The user may also share their information with their caregivers, physicians, and health care providers.

When used with a prescribed MDI, the System can report on information captured during the normal course of use, such as the time between actuations, that can be helpful in assessing MDI technique.

When used under the care of a physician with a prescribed MDI, the System can assist in the management of respiratory health symptoms and exacerbations by providing feedback through reminders, notifications, and self-management education.

The System is intended to be used in populations from Child (>2 years) to adult.

The System can be used both indoors and outdoors; home, work, and clinical settings, as well as on aircraft.

The System may also be used in clinical trials where researchers need to know information about the use of MDI medications by a participant.

The output of the System is not intended to diagnose or replace a diagnosis provided by a licensed physician. The System is not intended for use as an MDI or inhaled medication dose counter, nor is it intended to indicate the quantity of medication remaining in an MDI or inhaled medication.

Technology Comparison and Device Description:

The subject device uses technology similar to the predicate device, Propeller System Model 2 OTC, including Bluetooth wireless connectivity which connects to a similar mobile application. Like the predicate device, the Ochsner Connected Inhaler Sensor has an enclosure to fit to inhalers and a button for detecting use.

Technological characteristics of the Ochsner System and the Comparison: Propeller System are largely equivalent. Similarities include the indications for use, basic principle of operation, data collection information, time of data recording via internal clock, utilization of software for varying types of data review and modification, nonrechargable batteries, and the use of Bluetooth, low energy. The Ochsner System employs these technological characteristics in a similar way to the predicate device.

By reviewing the recorded data displayed by the Ochsner System, the physician or care provider can identify that a patient's state is worsening, and as a result, may choose to take action, such as contacting their patient. These aspects of the device have been verified and validated in order to establish equivalent performance to the equivalent device. This information indicates that the Ochsner System is equivalent to the predicate device in terms of device safety and effectiveness.

Comparison Table

Technology Characteristics	Predicate Device:	Candidate Device:	Comparison
	Propeller System, Propeller Sensor Model 2 OTC	Ochsner Connected Inhaler Sensor	N/A
	510(k) Number: K142516	510(k) Number: K223367	
Prescription/OTC	OTC	OTC	Equivalent
Indications for Use	The Propeller	The Ochsner	Similar: The
	System includes the	Connected Inhaler	Ochsner

Propeller MDI
Model 2 Sensor. The
sensor is an
accessory device
intended for singlepatient use to assist
physicians and
patients in
recording and
monitoring the
actuations of
prescribed MDI
usage.

The Propeller Mobile Application records, stores, and transmits usage events from Propeller Sensors, or via manual user entry, to a remote storage system. With the Propeller Mobile Application the user can review information collected from the MDI sensor, and report and review symptoms and other information about their disease management and its impact. The user may also share their information with their caregivers, physician, and healthcare providers.

The Propeller Web Application is software that, like Sensor System includes the Ochsner Connected Inhaler Sensor. The Sensor is an accessory device intended for single-patient use to assist physicians and patients in recording and monitoring actuations of prescribed MDI usage.

The Ochsner Connected Inhaler Sensor Mobile Application records, stores, and transmits usage events from the Sensors to a remote storage system. With the Mobile Application, the user can review information collected from the Sensors and report and review symptoms and other information about their disease management and its impact. The user may also share their information with their caregivers, physicians, and health care providers.

When used with a prescribed MDI, the System can report on information captured during the normal

Connected Inhaler System does not include a Web Application. The functionalities of the web application are contained entirely within the Ochsner Connected Inhaler Mobile Application since mobile application is necessary for device use. It was determined that there would be no increase in risk for users from this change.

Otherwise the IFUs for both products contain the same purpose for the devices, patient populations, and use environments.

the Propeller Mobile Application, is intended to allow users to review the collected information and characteristics of their MDI and its use, to capture other patientreported information and outcomes, and to allow that information to be shared with their caregivers, physicians, and health care providers.

When together with a prescribed MDI, the system can report on information captured during the normal course of use, such as the time between actuations that can be helpful in assessing MDI technique.

When together with a prescribed MDI, the system can be used to reduce the frequency of respiratory health symptoms and exacerbations by increasing adherence to MDI course of use, such as the time between actuations, that can be helpful in assessing MDI technique.

When used under the care of a physician with a prescribed MDI, the System can assist in the management of respiratory health symptoms and exacerbations by providing feedback through reminders, notifications, and self-management education.

The System is intended to be used in populations from Child (>2 years) to adult.

The System can be used both indoors and outdoors; home, work, and clinical settings, as well as on aircraft.

The System may also be used in clinical trials where researchers need to know information about the use of MDI medications by a participant.

The output of the System is not

medications through the use of feedback such as reminders and notifications, and self-management education.

The Propeller System is intended to be used in populations from Child (>2 years) to Adult.

The Propeller System can be used both indoors and outdoors; home, work, and clinical settings, as well as on aircraft.

The Propeller
System may also be used in clinical trials where researchers need to know information about the use of MDI medication(s) by a participant.

The output of the Propeller System is not intended to diagnose or replace a diagnosis provided by a licensed physician. The Propeller System is not intended for use as an MDI dose counter, nor is it intended to indicate

intended to diagnose or replace a diagnosis provided by a licensed physician. The System is not intended for use as an MDI or inhaled medication dose counter, nor is it intended to indicate the quantity of medication remaining in an MDI or inhaled medication.

	the quantity of medication remaining in an MDI or inhaled medication.		
Design – Attachment to Medication Dispenser	Physically attached to dispenser without inhibiting patient use	Physically attached to dispenser without inhibiting patient use	Equivalent
Principle of Operation	Attaches to the top of the medication canister and performs wireless uploading of usage history of the inhaler	Attaches to the top of the medication canister and performs wireless uploading of usage history of the inhaler	Equivalent
Output Port and Computer Interface	Wireless uploading to database; viewed by PC or other internet-capable devices	N/A	Different: The candidate device is not web-based
Data Collection Technology	Records date and time of MDI usage with button press switch	Records date and time of MDI usage with button press switch	Equivalent
Mobile Platforms	iOS version 7 or higher Android operating	iOS version 14 or higher Android operating	Different: The candidate device supports more
	system version 4.3 or higher	system version 11 or higher	recent versions of mobile platforms
Required Off-the-Shelf Hardware	Apple smartphones or devices with Bluetooth, iOS 7 or higher	Apple smartphones or devices with Bluetooth, iOS 14 or higher	Similar: The candidate device only works on mobile and
	Android smartphones or devices with Bluetooth	Android smartphones or devices with Bluetooth, version	thus only requires apple and android hardware.
	Internet capable device; no processor or memory	11 or higher	

	requirements (see		
Required Browser	Required Browser) Firefox, Chrome,	N/A	Different: The
Required Browser	Safari, Internet	17/1	candidate
	Explorer		device is not
			web-based
Mobile Application	Propeller Health	Ochsner Connected Inhaler Sensor	Equivalent
	Mobile Application records, stores, and	Mobile Application	
	transmits usage	records, stores, and	
	events from the	transmits usage	
	Propeller Health	events from the	
	Sensor via a feature	Ochsner Connected	
	or smartphone and	Inhaler Sensor via a	
	can be used to review	feature or	
	the information	smartphone and can be used to review the	
	captured when using a smartphone	information captured	
	a smartphone	when using a	
		smartphone	
Software	The Propeller Health	N/A, not web based.	Different: The
	Web Application is	But the mobile	candidate
	software intended to	application allows	device is not
	allow users to review the collected	the user to share information with	web based. But its mobile
	information and	their provider and	application
	characteristics of	add other details	allows the
	MDI use, to add	related to their	user to share
	detail associated with	condition for which	information
	a recorded usage	their MDI	with their
	event, and to share	medication(s) are	provider and
	that information with	prescribed.	add other
	their physician in order to provide		details related to their
	additional		condition for
	information		which their
	associated with the		MDI
	condition for which		medication(s)
	their MDI		are prescribed.
	medication(s) are		
Danie Carrettan	prescribed.		T 1
Dose Counter Records Usage	Ma	M _a	
INCLUIUS USARE	No Vac	No Vac	Equivalent
	Yes	Yes	Equivalent
Records Location of Usage (GPS Coordinates)			•

	wireless device if	wireless device if	
	paired with a Sensor	paired with a Sensor	
Keyboard/Input Interface	Dual button	Single button	Similar: The
	interface: primary	interface	candidate
	button and secondary		device also
	button		uses buttons
			as input, but
			only 1.
Digital Display	No	No	Equivalent
Power Source	Single 3V DC Li-ion	Single 3V DC Li-ion	Equivalent
	battery	battery	
Battery Life	1.5 years	1 year	Similar
Low Battery Indicator	Yes, light	Software displays	Similar:
	combination;	battery life, no on	Candidate
	software display of	device indicator	device only
	battery life		uses software
			to display
			battery, while
			the predicate
			also has a
			light-based
			indicator
Patient Reminder	Yes	Yes	Equivalent
Support	Yes	Yes	Equivalent
Patient Data Storage with	Yes	Yes	Equivalent
Software		**	
Patient Data Report	Yes	Yes	Equivalent
Generation with Software	**	**	D 1 1
Patient Data Graphs	Yes	Yes	Equivalent
Generation	37	37	F 1 1
Data Retrieval from Device with Software	Yes	Yes	Equivalent
Case Material – Patient	Lexan polycarbonate	Silicon Rubber	Different: The
Contact by Intact Skin	Lexan poryearoonate	Silicon Rubbel	candidate
(Hands)			device uses a
(Limitus)			different but
			also
			biocompatible
			case material
Electrical Safety	IEC 60601	IEC 60601	Equivalent
Biocompatibility	ISO 10993	ISO 10993	Equivalent
Sterility	Non-sterile	Non-sterile	Equivalent
Stermty	1 voii-sterife	1 NOII-SICITIC	Lquivalent

Test Summary:

Test results indicate that the Ochsner Sensor and its predicate Propeller Sensor Model comply with predetermined specifications. Software verification and validation testing confirms this result.

Non-clinical testing has been carried out to cover functional verification and device performance. This included completion of software verification and validation procedures, with performance testing of the MDI actuation sensor system to ensure data is logged accurately for MDI usage. This established correct functionality of the Ochsner System according to the requirements. Third party testing of the Ochsner System for compliance to IEC 60601 series standards for general safety and electromagnetic compatibility and ISO 10993 series standards for biocompatibility was completed by accredited laboratories prior to this submission. Complete, detailed reports are included in the application for clearance; summary information is included below.

The above testing confirms that the device is substantially equivalent to the predicate device.

Software Testing:

Software and Firmware for the Ochsner System was designed and developed according to a robust software development process aligned with "Design Control Guidance for Medical Device Manufacturers" "The Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices", "Guidance for Off the Shelf Software Use in Medical Devices", and verified and validated using guidance from the "General Principles of Software Validation" as recommended by FDA. Test results indicate that the Ochsner System complies with its predetermined specifications.

Electrical Safety Testing:

The Ochsner Sensor has successfully completed patient safety testing according to IEC 60601-1.

Electrical Compatibility

Testing:

The Ochsner Sensor has successfully completed EMC testing

Compatibility according to IEC 60601-1-2.

Performance Testing – Bench:

The Ochsner System has successfully completed performance testing according to applicable standards and internal testing. Important to highlight in this summary, is the successful performance testing that was completed for wireless/Bluetooth technology in accordance with specifications and also with, "FDA's Guidance on Radio-Frequency Wireless Technology in Medical Devices". In addition, tests required for FCC licensing were successful.

Clinical Testing: 1

No clinical testing was required.

Hazard Analysis for OTC:

Hazard Analysis for OTC included a review of existing hazards as well as how the patient obtains and learns about the system, registers for the system, installs the sensor, uses the Ochsner System to track MDI medication use, shares data with their physician/care team and obtains help & support with OTC labeling. No concerns of safety with the proposed OTC indication were found.

Conclusion:

Hardware testing carried out for the Ochsner System indicates it meets design and performance functional requirements. Software verification demonstrates that device features are effective, and that the system configuration functions equivalently to the predicate device. The Ochsner System also meets standard requirements for electrical safety, electromagnetic compatibility, biocompatibility, and wireless technology in medical devices.

Based upon this comparison of the predicate, and the accompanying testing results for the Ochsner Connected Inhaler Sensor, the Ochsner System is substantially equivalent to the predicate device.