

June 22, 2023

Qingdao Kingon Medical Science and Technology Co., Ltd. Zhang Benrong
Technical Supervisor
RM.1711, Building K, NO.101 Science Ave International
Creative Valley
Qingdao Free Trade Zone, Shandong 266555
China

Re: K223379

Trade/Device Name: Portable Oxygen Concentrator (Model: P2-E7, P2-E)

Regulation Number: 21 CFR 868.5440

Regulation Name: Portable Oxygen Generator

Regulatory Class: Class II Product Code: CAW Dated: May 24, 2023 Received: May 24, 2023

Dear Zhang Benrong:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



For

James J. Lee, Ph.D.

Director

DHT1C: Division of Sleep Disordered Breathing.

Respiratory and Anesthesia Devices

OHT1: Office of Ophthalmic, Anesthesia, Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K223379		
Device Name Portable Oxygen Concentrator (Model: P2-E7, P2-E)		
Indications for Use (Describe) The Portable Oxygen Concentrator (Model: P2-E7, P2-E) is intended to provide supplemental oxygen in a home, institutional, or travel environment.		
Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)		
CONTINUE ON A SEPARATE PAGE IF NEEDED.		

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary as required by section 807.92(c)

I. Date of the summary prepared: 10/30/2022

II. Administrative Information

Manufacturer	Establishment registration number	3014777423
information	Owner/Operator Number	10061814
	Name	Qingdao Kingon Medical Science and Technology Co., Ltd.
	Address	Room 301-302, No.15 HanchengRoad, Qingdao Free Trade Zone, Shandong, China, 266555
	Contact Person	Name: Benrong Zhang Address: Room 301-302, No.15 HanchengRoad, Qingdao Free Trade Zone, Shandong, China, 266555 TEL: +86-18565833539 FAX: +86 532 58792324 Email: augus@kingonmed.com
Submission	Contact Person	Name: Tracy Nanan Address: Room 301-302, No.15 HanchengRoad, Qingdao Free
Correspondent		Trade Zone, Shandong, China, 266555 TEL: +86-18565833539 FAX: +86 532 58792324 Email: ptg2022ptg@163.com

III. Device Information

Type of 510(k)	Special 510(k): Device Modification	
Prior submission	No prior submission	
Common Name Generator, Oxygen, Portable		
Classification	Portable oxygen generator	
name		
Trade Name	Portable Oxygen Concentrator (Model: P2-E7, P2-E)	
Review panel	Anesthesiology	



Product code	CAW
Regulation	868.5540
Number	
Regulation Class	2

IV. Predicate Device Information

Common Name	Generator, Oxygen, Portable
Classification Portable oxygen generator	
name	
Trade Name	Portable Oxygen Concentrator, model: P2-E6
Review panel	Anesthesiology
Product code	CAW
Regulation	868.5540
Number	
Regulation Class	2

V. Device description

Portable Oxygen Concentrator (Model: P2-E7, P2-E) is a portable oxygen generator that is intended to release oxygen for respiratory therapy by means of physical means (a molecular sieve). It supplies a pulsed high concentration of oxygen and is used with a nasal cannula to channel oxygen from the concentrator to the patient. The Portable Oxygen Concentrator is small, portable and may be used in home, institutional, or travel environment.

The portable oxygen concentrator consists of two parts: an oxygen concentrator and accessories. The oxygen concentrator is composed of compressor, battery, solenoid valve, molecular sieve, circuit control system, heat dissipation system, and a flow control device. Accessories include power adapters.

Model difference: The only different between model P2-E7 and P2-E is that P2-E7 has biggest output oxygen flow of 1.4L/min (at 7 gears) and P2-E has biggest output oxygen flow of 1.0L/min (at 5 gears).



VI. Principle of operation

Both model P2-E7 and P2-E has the same principle of operation. The portable oxygen concentrator works by getting use of the molecular sieves character that the internal pressure of a sealed container containing of molecular sieve will increase when injecting air into it. At this time, the molecular sieve will absorb a lot of nitrogen in the air with the increasing of ambient pressure, while the oxygen in the air is still existed in gaseous form, then the oxygen are collected through some pipelines. When the nitrogen absorption process in the container reaches a certain level, then exhaust of the vacuum container and nitrogen will be released from molecular sieve with the ambient pressure decreases. It will detect when the user begins to take a breath and then delivers a pulsed volume of oxygen during the inhalation period. The volume of the oxygen pulse is dependent on the setting value.

VII. Indications for Use

The Portable Oxygen Concentrator (Model: P2-E7, P2-E) is intended to provide supplemental oxygen in a home, institutional, or travel environment.

VIII. Comparison with predicate device

ID	Comparison Items	Subject device	Predicate device	Comparison
1.	510K Number	/	K210371	/
2.	Manufacturer	Qingdao Kingon Medical Science and Technology Co., Ltd.	Qingdao Kingon Medical Science and Technology Co., Ltd.	Same
3.	Device name	Portable Oxygen Concentrator	Portable Oxygen Concentrator	/
4.	Model	P2-E7, P2-E	P2-E6	/
5.	Classification	21CFR 868.5440	21CFR 868.5440	Same
6.	Product Code	CAW	CAW	Same
7.	FDA Class	II	II	Same
8.	Indications for Use	The Portable Oxygen Concentrator (Model: P2-E7, P2-E) is intended to provide supplemental oxygen in a home, institutional, or travel environment.	The Portable Oxygen Concentrator, model: P2-E6 is intended to provide supplemental oxygen in a home, institutional, or travel environment.	Same
9.	Environment of Use	Home, institutional, or travel environment. Home, institutional, or travel environment.		Same
10.	Design	table type table type		Same
11.	Prescriptive	Yes	Yes	Same
12.	Patient Population	Adult	Adult	Same



13.		Wiring cover: PC+ABS	Wiring cover: PC+ABS	
	Material of	Intake hood: PC+ABS	Intake hood: PC+ABS	
	Patient contact	Nozzle fitting: Aluminum alloy	Nozzle fitting: Aluminum alloy	Same
	components	Button panel: PET	Button panel: PET	
		Main housing: PC+ABS	Main housing: PC+ABS	
14.	Duration and	Type of contact: surface device;	Type of contact: surface device;	
	type of contact	Duration: permanent (> 30 d);	Duration: permanent (> 30 d);	Same
	type of contact			
15.		ISO 10993- 5 tested for Cytotoxicity;	ISO 10993-5 tested for Cytotoxicity;	
		ISO 10993-10 tested for Sensitization	ISO 10993-10 tested for Sensitization	
	Complete list of	and Irritation;	and Irritation;	
	all the			Same
	biocompatibility	ISO 18562-2 tested for Particulate	ISO 18562-2 tested for Particulate	
	tests performed	matter;	matter;	
		100 40502 2 + 14 14 14 14	100 10502 0 to 110 1/1 1/1	
		ISO 18562-3 tested for Volatile	ISO 18562-3 tested for Volatile	
1.0	Cin ale Detient	organic Compounds;	organic Compounds;	
16.	Single Patient, multi-use	Yes	Yes	Same
17.	Patient			
17.	Interface	Cannula Port	Cannula Port	Same
18.	Interrace	Pressure Swing Adsorption with	Pressure Swing Adsorption with	
10.	Technology	molecular sieve	molecular sieve	Same
19.	Dimensions	6.30"H*3.35"W*8.70"L	6.30"H*3.35"W*8.70"L	Same
20.	Dimensions	4.34lbs±0.07lbs (with standard	4.34lbs±0.07lbs (with standard	Same
20.	Weight	battery)	battery)	danie
21.	Oxygen	90%-3%/+6% at all settings	Successy	
	Concentration		90%-3%/+6% at all settings	Same
22.		P2-E7: adjustable in 1 increments		
		from 1 to 7.		Different
	Setting	P2-E: adjustable in 1 increments from	adjustable in 1 increments from 1 to 6	(See below note
		1 to 5.		ID_22)
23.		P2-E7: 50mL per breath at setting 5		
	Pulse mode	with 20BPM	50mL per breath at setting 5 with	Como
	bolus size	P2-E: 50mL per breath at setting 5	20BPM	Same
		with 20BPM		
24.	Principle of	by means of molecular sieve	by means of molecular sieve	Same
	operation	by means of molecular sieve	by means of molecular sieve	Jairie
25.	Filters	Input Filter, Patient Filter	Input Filter, Patient Filter	Same
26.	Breath rate	10 - 40 Breath per minute	10 - 40 Breath per minute	Same
27.	User Interface Buttons, LCD Display		Buttons, LCD Display	Same
28.	Power	AC adaptor: 100-240VAC ;50-60 Hz	AC adaptor: 100-240VAC ;50-60 Hz	Same
	requirements	in,19VDC 5.26A out	in,19VDC 5.26A out	



		DC adaptor: 12 - 16V DC in, 19V 6A out	DC adaptor: 12 - 16V DC in, 19V 6A out	
29.	Maximum oxygen discharge pressure	P2-E7: 20.6 PSI (142KPa) P2-E: 20.9 PSI (144KPa)	18.3 PSI (126KPa)	Different (See below note ID_29)
30.	Inspiratory		Same	
31.	Software	Embedded	Embedded	Same
32.	Acoustic Noise	P2-E7: 58.8 dBA at 1.4 LPM	58.2 dBA at 1.2 LPM	Different (See below note
	7.10000010110100	P2-E: 52.0 dBA at 1.0 LPM	30.2 db//dt 1.2 E/W	ID_32)
33.		Battery empty	Battery empty	Same
34.		Low pressure	Low pressure	Same
35.		No pulse	No pulse	Same
36.		High Temp	High Temp	Same
37.	- Alarms	Compressor Failure	Compressor Failure	Same
38.	Alaittis	Fan Failure	Fan Failure	Same
39.		Low Flow	Low Flow	Same
40.		Low Battery	Low Battery	Same
41.		No Breath Detected	No Breath Detected	Same
42.		EEPROM Failure	EEPROM Failure	Same
43.		Flow rates	Flow rates	Same
44.	Status Indicator	Battery Condition	Battery Condition	Same
45.		Alarms	Alarms	Same
46.	Battery Duration	P2-E7: Up to 4.5 hours at 0.21 LPM P2-E: Up to 4.5 hours at 0.21 LPM	Up to 4.5 hours at 0.21 LPM	Same
47.	Operating Environment	Temperature: 41 to 104°F (5 to 40°C) Humidity: 10% to 90%, non- condensing Altitude: 0 to 10,000 ft. (0 to 3048 meters)	Temperature: 41 to 104°F (5 to 40°C) Humidity: 10% to 90%, non- condensing Altitude: 0 to 10,000 ft. (0 to 3048 meters)	Same
48.	Shipping Storage	Temperature: -4 to 158°F (-20 to 70°C) Humidity: 5% to 90%, non-condensing Store in a dry environment	Temperature: -4 to 158°F (-20 to 70°C) Humidity: 5% to 90%, non-condensing Store in a dry environment	Same
49.	Electrical Safety	AAMI ANSI ES 60601-1	AAMI ANSI ES 60601-1	Same



50.	electromagnetic compatibility	IEC 60601-1-2	IEC 60601-1-2	Same
51.	Biocompatibility	4 VOC's less than ambient	4 VOC's less than ambient	Same
52.		ANSI AAMI ES 60601- 1: 2005 / (R)	ANSI AAMI ES 60601- 1: 2005 / (R)	
		2012 and A1: 2012	2012 and A1: 2012	
		IEC 60601-1-2: 2014	IEC 60601-1-2: 2014	
		IEC 60601-1-11: 2015	IEC 60601-1-11: 2015	
		IEC 60601-1-8: 2006+ A1:2012	IEC 60601-1-8: 2006+ A1:2012	
	Standards Met	ISO 80601-2-69: 2014	ISO 80601-2-69: 2014	Same
	Staridards iviet	ISO 80601-2-67: 2014	ISO 80601-2-67: 2014	Same
		ISO 18562-2: 2017	ISO 18562-2: 2017	
		ISO 18562-3: 2017	ISO 18562-3: 2017	
		IEC 62133: 2012	IEC 62133: 2012	
		ISO 10993-5:2009	ISO 10993-5:2009	
		ISO 10993-10:2010	ISO 10993-10:2010	

Note:

- ID_22: The max setting of subject device P2-E7 is 7 and of subject device P2-E is 5 while the predicate device is 6, the extra setting 7 of P2-E7 introduce risks of increased output flow and increased Emission of particulate matter of gas pass way, those risks are mitigated by tests tested according to ISO 80601-2-69: 2014, ISO 80601-2-67: 2014 and ISO 18562-2: 2017, therefore the difference does not raise new questions of safety and effectiveness.
- ID_29: The Maximum oxygen discharge pressure of subject device P2-E7 and P2-E is different with predicate device P2-E6 (K210371), since the subject device has been tested against ISO 80601-2-69: 2014 with positive result, the difference of subject device with predicate device RS-00500 (K162433) do not raise new questions of safety and effectiveness.
- ID_32: The Acoustic Noise of subject device P2-E7 is litter bigger than predicate device P2-E6 (K210371) and subject device P2-E is litter smaller than predicate device P2-E6 (K210371), since the subject device P2-E7 has been tested against ISO 80601-2-69: 2014 and ANSI AAMI ES60601-1 with positive results, the difference of subject device with predicate device P2-E6 (K210371) do not raise new questions of safety and effectiveness.
- IX. Discussion of Non-Clinical Tests Performed for Safety and effectiveness are as

follows

The recognized consensus standards for safety of medical electrical equipment: ANSI AAMI ES60601-1, IEC 60601-1-11 for safety, IEC 60601-1-2 for electromagnetic compatibility, ISO 80601-2-69: 2014 and ISO 80601-2-67: 2014 for performance and IEC 62304 for software verification are complied. See below table for details:



Standards	Standards Name
ANSI AAMI ES60601- 1:2005/(R)2012 and A1:2012	Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance
IEC 60601-1-2: 2014	Medical Electrical Equipment Part 1-2: General Requirements For Basic Safety And Essential Performance Collateral Standard: Electromagnetic Disturbances Requirements And Tests
IEC 60601-1-11: 2015	Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirements And Tests
IEC 60601-1-8: 2006+ A1:2012	Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
ISO 80601-2-69: 2014	Medical electrical equipment. Particular requirements for the basic safety and essential performance of oxygen concentrator equipment
ISO 80601-2-67: 2014	Medical electrical equipment. Particular requirements for basic safety and essential performance of oxygen-conserving equipment
ISO 18562-2: 2017	Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds
ISO 18562-3: 2017	Biocompatibility evaluation of breathing gas pathways in healthcare applications. Tests for emissions of volatile organic compounds (VOCs)
IEC 62133: 2012	Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
IEC 62304:2006+A1:2015	Medical device software - Software life cycle processes

Electrical safety and electromagnetic compatibility (EMC)

Electrical safety and EMC testing were conducted on the subject device P2-E6. The system complies with the AAMI ANSI ES60601-1, IEC 60601-1, IEC 60601-1-8, IEC 60601-1-11, ISO 80601-2-67, and ISO 80601-2-69 standards for electrical safety and the IEC 60601-1-2 standard for EMC.

Software Verification and Validation Testing

Software verification and validation was performed for the subject device in accordance with Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices - Guidance for Industry and FDA Staff, May 2005.



Software Description:

The software for this device was considered as a "moderate" level of concern, since a failure or latent flaw in the software could result in Minor Injury, either to a patient or to a user of the device. The software of the system, on the whole, is accountable for the system scheduler of the device, including fulfilling alarm task, beep task, detect task, display task, monitor task, startup task, breath task, as well as dealing with USART1-interrupt, USART2-interrupt, RTC-interrupt, WDT-interrupt, and controlling the sensor module.

Programming language

STMicroelectronics and Texas Instruments.

- > Hardware platform
- a) Applicable Device Name: Portable Oxygen Concentrator
- b) Applicable model: P2-E7, P2-E
- c) Micro-controller used: STM32F103RCT6, STM32F103C8T6, MSP430G2755
- d) ROM Size: STM32F103RCT6 256K, STM32F103C8T6 64K, MSP430G2755 32K
- e) RAM Size: STM32F103RCT6 48K, STM32F103C8T6 20K, MSP430G2755 -4K
- f) Software Release Version Number: Keil uVision5 V5.25.20, IRA for 430 V7.12.1

X. Discussion of Clinical Accuracy Testing Performed

There was no clinical testing performed.

XI. Conclusions

The Portable Oxygen Concentrator (Model: P2-E7, P2-E) have the same intended use and similar characteristics as the cleared predicate device Portable Oxygen Concentrator, model: P2-E6. Moreover, bench testing contained in this submission supplied demonstrate that the differences existed between Portable Oxygen Concentrator (Model: P2-E7, P2-E) and Portable Oxygen Concentrator, model: P2-E6 (K210371) do not raise any new questions of safety or effectiveness.

The non-clinical tests support the safety of the device and the hardware and software verification and validation demonstrate that the Portable Oxygen Concentrator (Model: P2-E7, P2-E) performs as intended in the specified use conditions are same with predicate device. The performance tests demonstrate that the Portable Oxygen Concentrator (Model: P2-E7, P2-E) performs comparably to the predicate device that is currently marketed for the same intended use. Thus, Portable Oxygen Concentrator (Model: P2-E7, P2-E) is Substantially Equivalent (SE) to the predicate device.