

April 27, 2023

Alimetry Ltd. % Janice Hogan Partner Hogan Lovells US LLP 1735 Market Street, Floor 23 Philadelphia, PA 19115

Re: K223398

Trade/Device Name: Gastric Alimetry Regulation Number: 21 CFR§ 876.1735

Regulation Name: Electrogastrography System

Regulatory Class: II Product Code: MYE Dated: March 27, 2023 Received: March 27, 2023

Dear Janice Hogan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Shanil P. Haugen -S

Shanil P. Haugen, Ph.D.
Assistant Director
DHT3A: Division of Renal, Gastrointestinal,
Obesity and Transplant Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

K223398				
Device Name Gastric Alimetry				
Indications for Use (Describe) The Gastric Alimetry System is intended to record, store, view and process gastric myoelectrical activity as an aid in the diagnosis of various gastric disorders.				
Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)				
CONTINUE ON A SEPARATE PAGE IF NEEDED.				

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY Alimetry's Gastric Alimetry

Submitter's Name, Address, Telephone Number, Contact Person and Date Prepared

Alimetry Ltd.

Phone: +64 27 609 1886

Facsimile: -

Contact Person: Yaara Yarmut, Chief Regulatory Officer.

Date Prepared: November 8, 2022

Name of Device and Name/Address of Sponsor

Trade Name: Gastric Alimetry

Manufacturer: Alimetry

Address: 70 Symonds St.

Grafton

Auckland 1010 New Zealand

Common Name: Gastric Alimetry System

Classification Name:

Electrogastrography system, 21 CFR 876.1735, Product code: MYE

Class II

Predicate and Reference Devices

1. Trade Name: Gastric Alimetry (K213924) (primary predicate) Manufacturer: Alimetry Ltd.

2. Trade Name: Polygraf ID with POLYGRAM NET ElectroGastroGraphy Application

Software (K014269) (reference device)

Manufacturer: Medtronic A/S

Intended Use / Indications for Use

The Gastric Alimetry is intended to record, store, view and process gastric myoelectrical activity as an aid in the diagnosis of various gastric disorders.

Device Description

The Gastric Alimetry System is an electrogastrography (EGG) device, used for non-invasively measuring the myoelectrical activity of the stomach at the surface of the abdomen. The Gastric Alimetry System is intended to record, store, view and process gastric myoelectrical activity as an aid in the diagnosis of various gastric disorders.

The device is used to acquire and digitize gastric myoelectrical data and movement artifacts using an array with recording electrodes on an adhesive patch, which is used for recording the myoelectrical data from the skin surface. An App is used to set up the device and capture patient-reported symptom data. A report is provided to the clinicians at the end of the test which includes myoelectrical signal data for manual analysis, together with computed data summaries and plots. A Supplementary Report is also routinely available to clinicians that includes signal data from all 64 channels on the array.

In the modified Gastric Alimetry System, the following minor updates are introduced in order to provide additional data summaries within the Report:

- Addition of four post-processing data summary metrics:
 - Principal Gastric Frequency showing the frequency of myoelectrical activity occurring within the gastric range
 - 'BMI-Adjusted Amplitude' calculated amplitude for BMI up to the recommended device limit (BMI <35)
 - Gastric Alimetry Rhythm Index a calculated measure of the stability of the gastric rhythm
 - fed:fasted Amplitude Ratio a calculated ratio showing the change in the gastric myoelectrical amplitude after a meal stimulus.
- Addition of a 'Symptom Burden' tracking bar to the front page of the Report. This is a
 calculated average of the patient-reported symptom data already shown in the Report,
 and is provided as a summary next to the spectral plot as a convenience for clinicians.
- Addition of data tables to the Supplementary Report. These provide hour by hour read outs of the symptom logs and summary metrics provided in the main Report, and are made available in Table form as a convenient reference for clinicians.
- Some minor Report rearrangements, with the 'signal strength' and 'best 8 channel' plots moved from the main Report to the Supplementary Report.

The four additional metrics are equivalent to other metrics widely applied in the EGG literature, and are included in the Reference Device, with only minor updates that address recognized inaccuracies that may affect performance.

Technological Characteristics / Substantial Equivalence

The Gastric Alimetry device has the same intended use and indications for use and technological characteristics as the primary predicate. Both devices are intended to record, store, view and process gastric myoelectrical activity as an aid in the diagnosis of various gastric disorders.

Item	Subject Device	Primary Predicate Device	Discussion
Device	Gastric Alimetry System	Gastric Alimetry System	
Manufacturer / K#	Alimetry Ltd.	Alimetry Ltd. / K213924	

Item	Subject Device	Primary Predicate Device	Discussion	
Intended Use and indication for use	To record, store, view and process gastric myoelectrical activity as an aid to the diagnosis of various gastric disorders.	To record, store, view and process gastric myoelectrical activity as an aid to the diagnosis of various gastric disorders.	Same as the predicate.	
User Population	Medical professionals to set up and use the system	Medical professionals to set up and use the system	Same as the predicate.	
	Patients 18 years and older	Patients 18 years and older		
	Specialist GI physicians	Specialist GI physicians		
Accessories				
Setup and Charging Dock	Gastric Alimetry Dock for charging and user convenience	Gastric Alimetry Dock for charging and user convenience	Same as the predicate.	
Electrodes	Disposable; same components as ECG electrodes.	Disposable; same components as ECG electrodes.	Same as the predicate.	
	Peel-and-stick patch.	Peel-and-stick patch.		
Technological Characteristics				
Sampling Frequency	4 Hz	4 Hz	Same as the predicate.	
Low Frequency Range	DC	DC	Same as the predicate.	
High frequency range	2 Hz	2 Hz	Same as the predicate.	
Number of channels	64 (+2 reference)	64 (+2 reference)	Same as the predicate.	
Electrode to recorder interface	Reader located on the Array and directly connected to it.	Reader located on the Array and directly connected to it.	Same as the predicate.	
Screen	Dedicated tablet for system operation	Dedicated tablet for system operation	Same as the predicate.	
Motion Sensor	Accelerometer	Accelerometer	Same as the predicate.	
Data Storage	On data acquisition device and cloud server	On data acquisition device and cloud server	Same as the predicate.	
Weight of device on patient	225 grams – Reader and Array	225 grams – Reader and Array	Same as the predicate.	

Item	Subject Device	Primary Predicate Device	Discussion
Power Source	Battery powered	Battery powered	Same as the predicate.
Software	Setup (including clinical data), device control, data acquisition software and system checks via App running on a tablet.	Setup (including clinical data), device control, data acquisition software and system checks via App running on a tablet.	Same as the predicate.
Patient Symptom Logging	Available on App interface. Presented as plot with times.	Available on App interface. Presented as plot with times.	Same as the predicate.
Skin Preparation	Yes. Shave, measure, skin prep with abrasive conductive gel Additional array template marking step for convenience.	Yes. Shave, measure, skin prep with abrasive conductive gel Additional array template marking step for convenience.	Same as the predicate.
Sterilization	Electrodes are disposable, non- sterile. Reader and Dock are reprocessed, not supplied sterile. Cleaning and disinfection instructions using wipes provided.	Electrodes are disposable, non- sterile. Reader and Dock are reprocessed, not supplied sterile. Cleaning and disinfection instructions using wipes provided.	Same as the predicate.
Safety Features			
System Checks	Impedance monitor and connectivity check displayed prior to recordings.	Impedance monitor and connectivity check displayed prior to recordings.	Same as the predicate.
Instructions to Patient	Yes. User instructs patient to limit movement, talking and sleeping. Additional step of displaying these instructions to patient via App as an added safety measure.	Yes. User instructs patient to limit movement, talking and sleeping. Additional step of displaying these instructions to patient via App as an added safety measure.	Same as the predicate.
Reporting Features			
Visualization of myoelectrical waveforms	Yes, all channels.	Yes, all channels.	Same as the predicate.
Automated Spectral Analysis	Yes, by Fourier-transform	Yes, by Fourier-transform	Same as the predicate.

Item	Subject Device	Primary Predicate Device	Discussion	
Artifact Evaluation	By manual identification of noise in waveforms and with reference to motion sensor.	By manual identification of noise in waveforms and with reference to motion sensor.	Same as the predicate.	
	Additionally provides automated noise detection for convenience.	Additionally provides automated noise detection for convenience.		
Dominant Frequency and Amplitude of Myoelectrical Activity	Yes, data tables, and visualization by spectral graphs and myoelectrical waveforms.	Yes, data tables, and visualization by spectral graphs and myoelectrical waveforms.	Same as the predicate.	
Other Myoelectrical Parameters	Measures of % time in frequency ranges, changes in signal power over time, and cross-channel coupling (wave propagation across the stomach) are available by visual inspection of spectral maps and waveforms.	Measures of % time in frequency ranges, changes in signal power over time, and cross-channel coupling (wave propagation across the stomach) are available by visual inspection of spectral maps and waveforms.	Equivalent. Same plots available for visual inspection as the predicate. Additional spectral metric outputs equivalent to technological reference.	
	Additionally processes these outputs as tables for convenience.	No additional processing of these metrics.	See substantial equivalence discussion.	
Symptom Outputs	Yes. Arranged in time plots for convenience.	Yes. Arranged in time plots for convenience	Same as the predicate.	
Technical / recording quality outputs	Yes. Impedance, motion. Additionally shows Array position and spatial distribution of amplitude for user convenience.	Yes. Impedance, motion. Additionally shows Array position and spatial distribution of amplitude for user convenience.	Same as the predicate.	
Standards with which the Device Complies				
Electrodes	ANSI/AAMI EC12:2000 – compliance with relevant requirements.	ANSI/AAMI EC12:2000 – compliance with relevant requirements.	Same as the predicate.	
Medical Electrical Equipment	IEC 60601-1	IEC 60601-1	Same as the predicate.	
Electromagnetic Compatibility	IEC 60601-1-2	IEC 60601-1-2	Same as the predicate.	

Performance Data

The four data summary metrics included in the modified Gastric Alimetry System are technically equivalent to those used in the Reference Device (Medtronic Polygram NET EGG

System). Minor technical updates were made to the four metrics to address issues that have been recognized to affect the performance of these metrics in the years since the Reference Device was approved. Therefore, additional performance testing was conducted to confirm that substantial equivalence was maintained as compared to the Reference Device.

These analyses were performed using post-market data from a study of 86 subjects, being 43 patients with chronic nausea and vomiting syndromes, gastroparesis, and functional dyspepsia, and 43 healthy matched controls. Each of the four metrics included in the modified Gastric Alimetry System were subjected to direct comparison with the equivalent Reference Device metrics, using this comparison cohort. In all four comparisons, high correlations were demonstrated (r>0.91; p<0.0001), confirming the substantial equivalence of all metrics.

Substantial Equivalence

The modified Gastric Alimetry System has the same intended uses and the same indications, technological characteristics, and principles of operation as its predicate device, the Gastric Alimetry System. The minor technological differences between the Gastric Alimetry and its predicate devices do not raise different questions of safety or effectiveness. In addition, performance testing based on post-market data analysis in patients with various gastric disorders and matched controls showed that the four additional metrics in the modified Gastric Alimetry System remained highly correlated to their equivalent metrics in the Reference Device after minor updates to correct recognized performance issues. Therefore the modified Gastric Alimetry System is shown to be substantially equivalent.