

Xuzhou Kaixin Electronic Instrument Co., Ltd. % Long Yang, COO Shenzhen Hlongmed Biotech Co., Ltd. 1201, Haosheng Business Center, 4096 Dongbin Road, Nanshan Shenzhen, P.R.C Shenzhen, Guangdong 518054 CHINA July 7, 2023

Re: K223448

Trade/Device Name: Bladder Scanner Model: BVT02

Regulation Number: 21 CFR 892.1560

Regulation Name: Ultrasonic pulsed echo imaging system

Regulatory Class: Class II Product Code: IYO, ITX Dated: June 7, 2023 Received: June 7, 2023

Dear Long Yang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

K223448 - Long Yang Page 2

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Julie Sullivan -S

Julie Sullivan, Ph.D.
Director
DHT8C: Division of Radiological Imaging and Radiation Therapy Devices
OHT8: Office of Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known)
K223448
Device Name Bladder Scanner, Model BVT02
Indications for Use (Describe) The BVT02 Bladder Scanner is B-mode pulsed-echo ultrasound device. It is intended as a portable device. The BVT02 Bladder Scanner projects ultrasound energy through the lower abdomen of the patient to obtain images of the bladder which is used to calculate bladder volume noninvasively. The BVT02 Bladder Scanner is intended to be used only by qualified medical professionals. Intended use environment: Professional health care facilities.
Type of Use (Select one or both, as applicable)
☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Diagnostic Ultrasound Indications for Use Form

System: BVT02 Bladder Scanner

Transducer: 2.5S120M2

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track	Specific (Track 1 & 3)	В	M	PWD	CWD	Color	Combined	Other*
1 Only)						Doppler	(Specify)	(Specify)
Ophthalmic	Ophthalmic							
	Fetal							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify)							
Fatal Lauraina	Neonatal Cephalic							
Fetal Imaging & Other	Adult Cephalic							
& other	Trans-rectal							
	Trans-urethral							
	Trans-esoph.(non-Card.)							
	Musculo-skeletal							
	(Conventional)							
	Musculo-skeletal(Superficial)							
	Intravascular							
	Other(Bladder)	N						
	Cardiac Adult							
	Cardiac Pediatric							
Cardiac	Intravascular (Cardiac)							
Caruiac	Trans-esoph.(Cardiac)							
	Intra-cardiac							
	Other (Specify)							
Peripheral	Peripheral vessel							
Vessel	Other (Specify)							

N=new indication; P= previously cleared by FDA; E=added under this appendix

510(k) Summary

(as required by 807.92(c))

The assigned 510(K) number is: <u>K223448</u>

Date of Summary: <u>2023-07-06</u>

1. Submitter information

Manufacturer Name: Xuzhou Kaixin Electronic Instrument Co., Ltd.

Address: Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.

Zip Code: 221004

Contact Person and Title: Chunmei Xu/ Technology Department Manager

Tel: 0086-516-87792599 Fax: 0086-516-87792599 Email: xingkexin@126.com

Zinan. Amgicames 120.ex

2. Contact person

2.1 Primary Contact Person

Long Yang (COO)

Shenzhen Hlongmed Biotech Company

1201, Haosheng Business Center, 4096 Dongbin Road, Nanshan, Shenzhen, P.R.C

Tel: 0086-755-86664986 Fax: 0086-755-86664933

E-mail: yanglong@hlongmed.com

2.2 Secondary Contact Person

Chunmei Xu (Technology Department Manager)

Xuzhou Kaixin Electronic Instrument Co., Ltd.

Kaixin Mansion, C-01, Economic Development Zone, Xuzhou, Jiangsu, China.

Zip Code: 221004

Tel:0086-516-87792599

E-mail: xingkexin@126.com

3. Device Information

Trade/Device Name	Bladder Scanner	
Model	BVT02	
Common Name	Diagnostic Ultrasound System with Accessories	
Classification Name	Ultrasonic Pulsed Echo Imaging System(IYO)/	
	Diagnostic Ultrasound Transducer(ITX)	
Regulatory Class	Class II	
Classification regulation	21CFR 892.1560 / 21CFR 892.1570	
Review Panel	Radiology	
Regulation Medical Specialty	Radiology	
Product Code	IYO/ITX	

4. Predicate Device

510(k) number	K131227
Device name	PadScan HD series Bladder Scanner
Sponsor	Caresono Technology Co., Ltd
Product Code	IYO/ITX

5. Intended Use

The BVT02 Bladder Scanner projects ultrasound energy through the lower abdomen of the patient to obtain images of the bladder which is used to calculate bladder volume noninvasively. The BVT02 Bladder Scanner is intended to be used only by qualified medical professionals.

6. Indications for Use

The BVT02 Bladder Scanner is B-mode pulsed-echo ultrasound device. It is intended as a portable device. The BVT02 Bladder Scanner projects ultrasound energy through the lower abdomen of the patient to obtain images of the bladder which is used to calculate bladder volume noninvasively. The BVT02 Bladder Scanner is intended to be used only by qualified medical professionals. Intended use environment: Professional health care

facilities.

7. Device Description

The BVT02 Bladder Scanner manufactured by Xuzhou Kaixin Electronic Instrument Co., Ltd. provides real - time ultrasound imaging and measuring, and also provides non - invasive volume measurement of the bladder. During image scanning, multiple 2D plane ultrasonic images are acquired in several seconds.

It features:

- Expert operating mode and Lite operating mode.
- Portable.
- Combined power supply with AC adapter and a battery.

8. Contraindications

Do not use the BVT02 Bladder Scanner on following case:

- a) Fetal use or pregnant patients.
- b) Patients with ascites.
- c) Patients with open or damaged skin.
- d) Wounds in the suprapubic region.

9. Comparison to Predicate Device

Xuzhou Kaixin Electronic Instrument Co., Ltd. believes the BVT02 Bladder Scanner described in this submission is substantially equivalent to the predicate devices as follows:

PadScan HD series Bladder Scanner (K131227)

The ultrasound power transmitted from the device is not user adjustable, and BVT02 Bladder Scanner is Track 1 System and meets the FDA's pre-amendment acoustic output limits, so as the predicate devices(PadScan HD 5, PadScan HD 3) are. Although there are some differences such as System Characteristics, Display, Patient Contacting Material, Range and Power, there is no significant differences in technological characteristics that affecting the safety and efficiently. These are evaluated by safety test and acoustic output test.

Table 6 Comparison to the predicate

(PadScan HD 5, PadScan HD 3)

(Fauscan fid 5, Fauscan fid 5)						
Item	Element Of	Proposed Device	Predicate Device-K131227			
	Comparison					
1	Trade Name	Bladder Scanner	PadScan HD series Bladder			
			Scanner			
2	Model	BVT02	PadScan HD 5, PadScan HD 3			
	510k submitter	Xuzhou Kaixin Electronic	Caresono Technology Co., Ltd.			
3		Instrument Co., Ltd.				
4	510(K) Number	K223448	K131227			
	Classifications	21 CFR 892.1560 Ultrasonic	21 CFR 892.1560 Ultrasonic			
	Name &	Pulsed Echo Imaging System	Pulsed Echo Imaging System			
_	Citations	(Product code: IYO)	(Product code: IYO)			
5		21 CFR 892.1570 Diagnostic	21 CFR 892.1570 Diagnostic			
		Ultrasonic Transducer	Ultrasonic Transducer			
		(Product code: ITX)	(Product code: ITX)			
	Intended Use	The BVT02 Bladder Scanner	The PadScan HD series Bladder			
		projects ultrasound energy	Scanner projects ultrasound			
		through the lower abdomen of the	energy through the lower			
		patient to obtain images of the	abdomen of the patient to obtain			
6		bladder which is used to calculate	images of the bladder which is			
0		bladder volume non-invasively.	used to calculate bladder volume			
		The BVT02 Bladder Scanner is	noninvasively. The PadScan HD			
		intended to be used only by	series Bladder Scanner is			
		qualified medical professionals.	intended to be used only by			
			qualified medical professionals.			
	Indications for	The BVT02 Bladder Scanner is	The PadScan HD series Bladder			
	Use	B-mode pulsed-echo ultrasound	Scanner is B-mode pulsed-echo			
		device. It is intended as a portable	ultrasound device. It intended as			
7		device. The BVT02 Bladder	a portable battery-operated			
		Scanner projects ultrasound	device. The PadScan HD series			
		energy through the lower	Bladder Scanner projects			
		abdomen of the patient to obtain	ultrasound energy through the			

Item	Element Of Comparison	Proposed Device	Predicate Device-K131227
	•	images of the bladder which is used to calculate	lower abdomen of the patient to obtain images of the bladder
		bladder volume non-invasively.	which is used to calculate
		The BVT02 Bladder Scanner is	bladder volume noninvasively.
		intended to be used only by	The PadScan HD series Bladder
		qualified medical professionals.	Scanner is intended to be used
		Intended use environment:	only by qualified medical
		Professional health care facilities.	professionals.
		Do not use the BVT02 Bladder	Do not use the PadScan HD
		Scanner on following cases:	series Bladder Scanner on
	Contraindicatio	a) Fetal use or pregnant patients.	following cases:
	ns	b) Patients with ascites.	a) Fetal use or pregnant patients
8		c) Patients with open or damaged	b) Patients with ascites
		skin.	c)Patients with open or damaged
		d) Wounds in the suprapubic	skin
		region.	d) Wounds in the suprapubic
			region
9	Modes of	B mode	B mode
	operation		
		a)Portable	a) Portable
		b)LCD Display	b) LCD Display
10	System	c)Power source: Battery or	c) Thermal Printer
	Characteristics	AD-DC adapter	d) Power source: Battery or
			AD-DC adapter
11	Display	3.5" TFT-LCD	PadScan HD5: 8" TFT-LCD
11			PadScan HD3: 7" TFT-LCD
	Controls for	No	No
12	Change of		
	acoustic output		

Item	Element Of Comparison	Proposed Device	Predicate Device-K131227
	during scan		
13	Transducer Type	Mechanical Sector Probe	Mechanical Sector Probe
14	Measurement localization	Abdomen	Abdomen
15	Transducer Resonant Frequency	2.5 MHz	2.5 MHz
16	Number of elements	1	1
17	Sector Angle	120 degrees	120 degrees
18	No. of Scan Planes	12	12
19	FDA Limits	Track 1	Track 1
20	Product Safety Certification	AAMI / ANSI ES60601-1:2005/(R)2012 And A1:2012, C1:2009/(R)2012 And A2:2010/(R)2012 IEC 60601-2-37:2015	IEC 60601-1:2005 +CORR.1(2006) +CORR.2(2007) IEC 60601-2-37:2007
21	EMC Compliance	IEC 60601-1-2:2014	IEC 60601-1-2:2007
22	Patient Contacting Material	Plastic, Medical PP (Skin Contact) Complies with ISO 10993	Plastic, PE (Skin Contact) Complies with ISO 10993
23	Range Bladder volume range: 20-999ml Accuracy: ±15%, ±15ml		Bladder volume range: 0-999ml Accuracy: ±15%, ±15ml
24	Classification of protection against	Class II equipment Type B equipment	Class II equipment Type B equipment

Item	Element Of Comparison	Proposed Device	Predicate Device-K131227
	electric shock		
26	Real-time scanning	Yes (Pre-scan)	Yes (Pre-scan)
27	Scan time	< 5 seconds	< 5 seconds
27	PC Data Upload	Using USB flash disk	Using USB flash disk
28	Power	AC/DC Adapter: Input:AC100-240V, 50-60Hz Output: DC12.8V 3.0A Battery: Li-ion rechargeable	AC/DC Adapter: Input: AC100-240V, 50/60Hz, Output: DC14V±0.5V Battery: Li-ion rechargeable

Xuzhou Kaixin Electronic Instrument Co., Ltd. believes that the BVT02 Bladder Scanner is substantially equivalent to the PadScan HD series Bladder Scanner of Caresono Technology Co., Ltd.

10. Non-clinical Testing Summary

10.1 Safety

Electrical, mechanical, environmental safety and performance data demonstrates that the device is in compliance with ANSI/AAMI ES60601-1:2005/(R)2012 And A1:2012, C1:2009/(R)2012 And A2:2010/(R)2012 and IEC 60601-2-37:2007+AMD1:2015 CSV.

10.2 EMC

Electromagnetic Compatibility data demonstrates that the device is in compliance with IEC 60601-1-2:2014.

10.3 Performance-Bench Testing

1) The BVT02 Bladder Scanner had been tested as Track 1 device per the FDA Guidance document "Marketing Clearance of Diagnostic Ultrasound Systems and Transducers" issued in June 27, 2019. The acoustic output is measured and calculated per IEC 62359:2010+AMD1:2017 CSV.

2) The BVT02 Bladder Scanner had been tested volume accuracy per the FDA Guidance document "Marketing Clearance of Diagnostic Ultrasound Systems and Transducers" issued in June 27, 2019. All the test results comply with the pre-set acceptability criterion, which is the same as predicate device.

10.4 Biocompatibility

The biocompatibility testing conducted in according with standard Biocompatibility ISO 10993-5:2009 and ISO 10993-10:2010.

11. Substantial Equivalence Conclusion

The BVT02 Bladder Scanner was evaluated with safety (AAMI / ANSI ES60601-1:2005/(R)2012 And A1:2012, C1:2009/(R)2012 And A2:2010/(R)2012 and IEC 60601-2-37:2007+AMD1:2015 CSV), EMC (IEC 60601-1-2:2014), Biocompatibility (ISO 10993-5:2009, ISO10993-10:2010), Acoustic Output (IEC 62359:2010+AMD1:2017 CSV) and volume accuracy. The conclusions drawn from testing of the BVT02 Bladder Scanner demonstrate that the device is as safe and effective as the legally marketed predicate device.