

April 17, 2023

Zhejiang Yinchili Medical Technology Co., Ltd. % Ivy Wang Technical Manager Shanghai Sungo Management Consulting Company Limited 14th Floor, 1500# Central Avenue Shanghai, Shanghai 200122 China

Re: K223454

Trade/Device Name: Custom-made Invisible Aligners Regulation Number: 21 CFR 872.5470 Regulation Name: Orthodontic Plastic Bracket Regulatory Class: Class II Product Code: NXC Dated: January 16, 2023 Received: January 17, 2023

Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bobak Shirmohammadi -S

For Michael E. Adjodha, M. ChE., CQIA Assistant Director DHT1B: Division of Dental and ENT Devices OHT1: Office of Ophthalmic, Anesthesia, Respiratory, ENT and Dental Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K223454

Device Name Custom-made Invisible Aligner

Indications for Use (Describe)

The Custom-made Invisible Aligners is indicated for the treatment of tooth malocclusion in patients with permanent dentition (i.e., all second molars). The Custom-made Invisible Aligners positions teeth by way of continuous gentle force.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) Summary

(As requirement by 21 CFR 807.92)

A. Applicant:

Zhejiang Yinchili Medical Technology Co., Ltd. Address: North 4F, No.239 Yatai Road, Nanhu District, Jiaxing, Zhejiang, CHINA Contact Person: Ms. Xinyan Zhang Tel: +86- 15800780940 Email: <u>zhangxinyan@smartee.cn</u> Date of summary prepared: 2023-03-30

Submission Correspondent:

Primary contact: Ms. Ivy Wang Shanghai SUNGO Management Consulting Co., Ltd. Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China Tel: +86-21-58817802 Email: haiyu.wang@sungoglobal.com Secondary contact: Mr. Raymond Luo Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China Tel: +86-21-68828050 Email: fda.sungo@gmail.com

B. Device:

Trade Name: Custom-made Invisible Aligners Common Name: Sequential Aligners Model: YCL-C <u>Regulatory Information</u> Classification Name: Aligner, Sequential Classification: Class II Product code: NXC Regulation Number: 21 CFR 872.5470 Review Panel: Dental

C. Predicate device:

K203624 (Primary Predicate) Custom-made Invisible Aligners Zhejiang Yinchili Medical Technology Co., Ltd.

Reference device K181739 Invisalign System with Mandibular Advancement Feature Align Technology, Inc.

K223454

D. Indications for use of the device:

The Custom-made Invisible Aligners is indicated for the treatment of tooth malocclusion in patients with permanent dentition (i.e., all second molars). The Custom-made Invisible Aligner positions teeth by way of continuous gentle force.

E. Device Description:

The Custom-made Invisible Aligners System is a series of dental aligners that fabricated of clear, thin thermoformed plastic to progressively reposition the teeth. The product is made from a composite multilayer material with thermoplastic polyurethane and polyethylene terephthalate-1, 4-cyclohexanedimethanol ester (TPU+PETG). Corrective force to reposition the teeth is delivered via minor changes into a position in each subsequent aligner. They are designed to move the teeth to the target position and deliver desired clinical effect.

F. Comparison to predicate device

The Custom-made Invisible Aligners is substantially equivalent in intended use, indications for use, mode of action, mode of use, design, and manufacturing to the predicate device. Only minor differences exist between the subject product and the predicate, which do not affect the safety or effectiveness of the subject device.

Table 1: Comparison to Predicate Device

Device	Subject Device	Predicate Device	Reference Device	Result
510K number	К223454	K203624	K181739	-
Model Name	Custom-made Invisible	Custom-made	Invisalign System	-
	Aligners	Invisible Aligners	with Mandibular	
			Advancement	
			Feature	
Classification	Class II Device, NXC (21	Class II Device, NXC	Class II Device, NXC	Same
	CFR 872.5470)	(21 CFR 872.5470)	(21 CFR 872.5470)	
Classification	Aligner, Sequential	Aligner, Sequential	Aligner, Sequential	Same
Name				
Indications for	The Custom-made	The Custom-made	The Invisalign	Same
use	Invisible Aligners is	Invisible Aligners is	System is intended	
	indicated for the	indicated for the	for the orthodontic	
	treatment of tooth	treatment of tooth	treatment of	
	malocclusion in	malocclusion in	malocclusion.	
	patients with	patients with		
	permanent dentition	permanent dentition		
	(i.e., all second molars).	(i.e., all second		
	The Custom-made	molars). The		
	Invisible Aligner	Custom-made		
	positions teeth by way	Invisible Aligner		
	of continuous gentle	positions teeth by way		

Table 1 provides a comparison of the subject and predicate device.

	force.	of continuous gentle		
		force.		
Mode of	Orthodontic tooth	Orthodontic tooth	Orthodontic tooth	Same
Action	movement occurs	movement occurs	movement occurs	
	through forces applied	through forces	through forces	
	by the device to the	applied by the device	applied by the	
	dentition as each tooth	to the dentition as	device to the	
	follows the	each tooth follows the	dentition as each	
	programmed	programmed	tooth follows the	
	displacement based on	displacement based	programmed	
	a doctor's prescription.	on a doctor's	displacement	
		prescription.	based on a doctor's	
			prescription.	
Anatomical	Oral cavity	Oral cavity	Oral cavity	Same
Site of Use				
Mode of Use	Each aligner is worn by	Each aligner is worn	Aligners are worn	Same
	the patient as	by the patient as	for approximately	
	determined by the	determined by the	1-2 weeks of 20-22	
	treating dental	treating dental	hours of wear per	
	practitioner, generally	practitioner, generally	day, after which it	
	for 22 nrs/day (or full	for 22 nrs/day (or full	is replaced by the	
	time except for eating	time except for eating	This is repeated for	
	and nygiene) for 2	and nygiene) for 2	this is repeated for	
	replaced by the port	replaced by the port	proscribed by the	
	alignor in soquence	alignor in soquence	Dontal Practitionor	
	This is repeated for a	This is repeated for a	Dental Flactitioner	
	duration as prescribed	duration as prescribed		
	hv a Dental	hy a Dental		
	Professional.	Professional.		
Application	Removable	Removable	Removable	Same
Raw Material	Multilayer TPU	Thermoplastic	The Invisalign	Same
Used	(Thermoplastic	copolyester	System uses either:	with the
	polyurethane) + PETG	(polyethylene	1. Multilayer	referenc
	(Polyethylene	terephthalate-ethylen	aromatic	e device.
	terephthalate-1,	e glycol copolyester)	thermoplastic	
	4-cyclohexanedimethan		polyurethane	
	ol ester)		/copolyester.	
			or	
			2. thermoformed	
			polyurethane	
Method of	Thermoforming	Thermoforming	Thermoforming	Same
Manufacturin				
g				
OTC or Rx	KX	KX	KX	Same

K223454

Sterile	No	No	No	Same
Biocompatibili	In compliance with ISO	In compliance with	In compliance with	Similar.
ty	10993, tests including	ISO 10993, tests	ISO 10993	Both are
	Cytotoxicity	including		in
	Oral Mucosa Irritation	Cytotoxicity		complian
	Sensitization	Oral Mucosa Irritation		ce with
	Sub chronic systemic	Sensitization		ISO
	toxicity			10993
	Genotoxicity			
Design			Not available	Similar.
		28		Both are
		A. and		transpar
		Merry Mar		ent
	8			plastic
				films.

G. Non-clinical Test

1) Performance Testing

Bench testing has demonstrated that the device is in compliance with pertinent standards and specifications, the expectations of the dental community and the product labeling. Performance testing was performed to the subject device including thickness, appearance, odor, density, water absorption, dissolution, color stability, tear resistance, wear resistance, Tensile modulus of elasticity to demonstrate its effectiveness.

Manufacturing validation accuracy testing

Manufacturing accuracy validation were conducted to the Custom-made Invisible Aligners. Aligners from 12 different patient case were evaluated at the beginning, middle and end throughout the sequence. The accuracy of 3D molding and aligner molding are checked and meet the pre-established specification. The suitability, function and form of the aligner were checked and comparing it to the treatment design in the software, and the results were complied with the pre-established specifications and acceptance criteria.

Shelf life – 3 years

A 3-year shelf life was determined by accelerated aging testing. Performance testing were conducted after 81 days of accelerated aging under 60° C. The test results showed conformity with the pre-established specifications and acceptance criteria.

2) Biocompatibility Testing

The biocompatibility evaluation for the device was conducted in accordance with "Use of International Standard ISO 10993-1, Biological evaluation and testing within a risk management process –Guidance for Industry and Food and Drug Administration Staff" as recognized by FDA. The aligner is considered mucosal membrane contacting for a duration of greater than 30 days. The testing included the following tests:

- Cytotoxicity
- Irritation

K223454

- Sensitization
- Sub chronic systemic toxicity
- Genotoxicity

The results of the testing met the requirements of the study protocols and the material is considered non-toxic, non-sensitizing and is not an intracutaneous irritant. The results of the studies further support a determination of substantial equivalence.

3) Software Verification and Validation Testing

Software verification and validation testing were conducted on the software that facilitates ordering and processing of the Custom-made Invisible Aligner to support that the device is as safe and effective as the predicates. Documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The software for this device was considered as a "moderate" level of concern.

H. Clinical Test Conclusion

No clinical study is included in this submission.

I. Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device K203624.