

June 13, 2023

Shanghai EA Medical Instruments Co., Ltd. % Breanne Butler Regulatory Affairs Consultant Prime Path Medtech 1321 Upland Dr. Suite 6792 Houston, Texas 77043

Re: K223518

Trade/Device Name: iOrtho

Regulation Number: 21 CFR 872.5470

Regulation Name: Orthodontic Plastic Bracket

Regulatory Class: Class II

Product Code: PNN Dated: May 19, 2023 Received: May 19, 2023

Dear Breanne Butler:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

K223518 - Breanne Butler Page 2

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Michael E. Adjodha -S

Michael E. Adjodha, M.ChE.
Assistant Director
DHT1B: Division of Dental and
ENT Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

K223518

Device Name :0-41Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

IOTINO
Indications for Use (Describe) iOrtho is intended for use as a medical front-end device providing tools for management of orthodontic cases, systematic inspection, detailed analysis, treatment simulation and virtual appliance design options (Export of Models, Indirect Bonding Transfer Media, Sequential aligners) based on 3D models of the patient's dentition before the start of an orthodontic treatment. It can also be applied during the treatment to inspect and analyze the progress of the treatment. It can be used at the end of the treatment to evaluate if the outcome is consistent with the planned/desired treatment objectives. The use of iOrtho requires the user to have the necessary training and domain knowledge in the practice of orthodontics,
as well as to have received a dedicated training in the use of the software.
Type of Use (Select one or both, as applicable) Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

> Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY - K223518

A summary of 510(k) substantial equivalence information for this traditional 510(k) in accordance with the requirements of 21 CFR 807.92.

Submitter: Shanghai EA Medical Instruments Co., Ltd.

No.1619 Huishan Avenue, Huishan Economic Development Zone

Wuxi, Jiangsu Province

China

Company Contact Person: Jessica Luo, Regulatory Affairs, Angel Align

Phone: +86.0510.83591717(181)
Email: luoyuqing@angelalign.com

Submission Correspondent: Breanne Butler, Regulatory Affairs Consultant

Address: 1321 Upland Dr. Suite 6792 Houston, TX 77043

Phone: 860-810-5594

Email: bbutler@primepathmedtech.com

Date Prepared: June 12, 2023

Trade Name Name: iOrtho

Common Name: Orthodontic plastic bracket (Software)

Classification Name: Orthodontic plastic bracket (Software)

Product Code: PNN

Device Classification: Class II, 21 CFR 872.5470

Predicate Device: 3Shape A/S Ortho System (K171634)

Reference Device: Wuxi EA Medical Instruments Technologies Limited Clear Aligner

(K203688)

Reference Device: Dentsply Sirona CEREC Ortho Software (K171122)

Device Description:

iOrtho (hereafter referred to as "Proposed Device") includes modifications to the currently marketed software included in K203688, cleared October 8, 2021 (hereafter referred to as "Reference Device"). The Proposed Device is an orthodontic appliance design and treatment simulation software. This software is for use by dental professionals to aid in diagnosis and design solutions for patients. Digital scans (3D) of a patient's dentition can be loaded into the software and the dental professional can then create treatment plans for each individual patient and their needs. The system can be used to fabricate 3D dental models

using standard stereolithographic (STL) files for use in 3D printers. These models can then be used as a template for thermoforming aligners or retainers by Angel Align technicians.

Indications for Use:

iOrtho is intended for use as a medical front-end device providing tools for management of orthodontic cases, systematic inspection, detailed analysis, treatment simulation and virtual appliance design options (Export of Models, Indirect Bonding Transfer Media, Sequential aligners) based on 3D models of the patient's dentition before the start of an orthodontic treatment. It can also be applied during the treatment to inspect and analyze the progress of the treatment. It can be used at the end of the treatment to evaluate if the outcome is consistent with the planned/desired treatment objectives.

The use of iOrtho requires the user to have the necessary training and domain knowledge in the practice of orthodontics, as well as to have received a dedicated training in the use of the software.

Comparison to Predicate Devices:

The Proposed Device is functionally equivalent to the following device: 3Shape A/S Ortho System (K171634, cleared Jan 17, 2018) (hereafter referred to as "Predicate Device"), and possesses minor differences to the previously cleared Reference Device to allow treating physicians include mandible repositioning and traction accessories during treatment planning and simulation. The following table demonstrates the functional specifications of the Proposed Device are substantially equivalent to the Predicate Device, minorly different to the Reference Device, and raises no new questions regarding safety and effectiveness of the device.

Device Comparison Table

Device Comparison Table						
	Proposed Device:	Predicate Device:	Reference	Reference		
	iOrtho	3Shape A/S Ortho	Device: Wuxi	Device:		
		System (K171634)	EA Medical	Dentsply		
			Instruments	Sirona	Comparison	
Specification			Technologies	CEREC	Result	
			Limited. Clear	Ortho	11000.110	
				Software		
			Aligner (K203688)	(K171122)		
Regulation	21 CFR 872.5470	21 CFR 872.5470	21 CFR 872.5470	21 CFR	Same	
Number				872.5470		
Classification	Orthodontic Plastic	Orthodontic Plastic	Orthodontic	Orthodontic	Same as	
Name	Bracket (Software)	Bracket (Software)	Plastic Bracket	Plastic Bracket	Predicate	
	` '	, ,		(Software)		
Product Code	PNN	PNN	NXC	PNN	Same	
Classification	Class II	Class II	Class II	Class II	Same	
	iOrtho	The 3Shape Ortho	Clear Aligners are	CEREC Ortho		
		System™	indicated for the	Software is		
	is intended for use	is intended for use	alignment of	intended for		
Indication for	as a medical front-	as a medical front-	teeth during	use with image	Similar to	
Use	end device	end device	orthodontic	data acquired	predicate	
	providing tools for	providing tools for	treatment of	from handheld		
	management of	management of	malocclusion	intra oral 3D		
	orthodontic cases,	orthodontic	in patients with	cameras and		

	Branged Davisa:	Predicate Device:	Reference	Reference	
	Proposed Device: iOrtho			Device:	
	1011110	3Shape A/S Ortho	Device: Wuxi	Dentsply	
		System (K171634)	EA Medical	Sirona	
Specification			Instruments	CEREC	Comparison
op comount in			Technologies	Ortho	Result
			Limited. Clear	Software	
			Aligner	(K171122)	
			(K203688)	(K1/1122)	
	systematic	models, systematic	permanent	desktop	
	inspection, detailed	inspection, detailed	dentition (i.e. all	laboratory	
	analysis, treatment	analysis, treatment	second molars).	scanners to	
	simulation and	simulation and		create 3D	
	virtual appliance	virtual appliance		virtual models	
	design options	design options		to be used for	
	(Export of Models,	(Custom Metal		data acquisition	
	Indirect Bonding	Bands, Export of		and modeling	
	Transfer Media,	Models, Indirect		analysis for	
	Sequential aligners)	Bonding Transfer		orthodontic	
	based on 3D models	Media) based on		patients and	
	of the patient's	3D models of the		conditions. The	
	dentition before the	patient's dentition		CEREC Ortho	
	start of an	before the start of		Software 3D	
	orthodontic	an orthodontic		model data can	
	treatment. It can	treatment. It can		be exported to	
	also be applied	also be applied		orthodontic	
	during the	during the		design software	
	treatment to	treatment to		to aid in the	
	inspect and analyze	inspect and analyze		design of	
	the progress of the	the progress of the		orthodontic	
	treatment. It can be	treatment. It can		appliances.	
	used at the end of	be used at the end			
	the treatment to	of the treatment to			
	evaluate if the	evaluate if the			
	outcome is	outcome is			
	consistent with the	consistent with the			
	planned/desired	planned/desired			
	treatment	treatment			
	objectives.	objectives.			
	The use of iOrtho	The use of the			
		Ortho System™			
	requires the user to	requires the user			
	have the necessary	to have the			
	training and domain	necessary training			
	knowledge in the	and domain			
	practice of	knowledge in the			

	Dranged Davises	Duadiasta Davissa	Deference	Doforosco	
Specification	Proposed Device: iOrtho	Predicate Device: 3Shape A/S Ortho System (K171634)	Reference Device: Wuxi EA Medical Instruments Technologies Limited. Clear Aligner	Reference Device: Dentsply Sirona CEREC Ortho Software (K171122)	Comparison Result
	orthodontics, as well as to have received a dedicated training in the use of the software. • Stand Alone Software • Imports Digital	practice of orthodontics, as well as to have received a dedicated training in the use of the software. • Stand Alone Software • Imports Digital	Produce 3D-model file of the PVS impression or an arms.	• Stand Alone Software • Imports	
Technological Features	 Dental Data Useful for as an aid in Diagnosis, and treatment planning Virtual Planning of tooth movement Patient follow-up monitoring and management Cephalometry 	Patient Scans Can be used to design Dental Casts Useful for Diagnosis, treatment planning, and CAD design Virtual Planning of tooth movement	impression or digital scan. Identifies the individual teeth that will require treatment (i.e. repositioning) Creates a treatment plan (i.e. 3-D models that represent the treatment plan). The treating dental practitioner reviews these images using software and has the option to reject or request modifications	Digital Patient Scans Can be used to design Dental Models	Similar to Predicate

Specification	Proposed Device: iOrtho	Predicate Device: 3Shape A/S Ortho System (K171634)	Reference Device: Wuxi EA Medical Instruments Technologies Limited. Clear Aligner (K203688) to the set-up	Reference Device: Dentsply Sirona CEREC Ortho Software (K171122)	Comparison Result
Minimum Hardware/Soft ware Requirements	 OS: Windows 7)/Mac OS 10.12 RAM: 2 GB Video Card: Any that supports webGL Hard Drive Space: 50MB Web browser (iOrtho): Chrome 63+, Firefox 58+, Safari 12+, Edge 79+ Internet Access 	 OS: Windows 7, 8, 10 64-bit RAM: 8 GB Video Card Memory: 1 GB Hard Drive Space: 250 GB CPU: Intel Core i5 or equivalent Mouse: with wheel button 	• OS: Windows 7, 10, Mac OS X • RAM: 2 GB (minimum) • Hard Drive Space: 50 MB • Video Card: Graphics Display Card supporting webGL • Web Browser: Chrome 9+, Firefox 4+, Safari 5.1+, Edge 5.1+, IE11+ • Internet Access	OS: Windows 7, 64-bit RAM: 8 GB Video Card Memory: 1 GB Hard Drive Space: 250 GB CPU: Intel QuadCore 1.6 GHz processor	Similar
Login Method	Username and password	Username and password	Username and password	Unknown	Same
Supported Anatomic Areas	Maxilla/Mandible	Maxilla/Mandible	Maxilla/Mandibl e	Maxilla/Mandib le	Same
Intended Use					

Specification Managing	Proposed Device: iOrtho	Predicate Device: 3Shape A/S Ortho System (K171634)	Reference Device: Wuxi EA Medical Instruments Technologies Limited. Clear Aligner (K203688)	Reference Device: Dentsply Sirona CEREC Ortho Software (K171122)	Comparison Result
Patient and case base data	Yes	Yes	Yes	Yes	Same
Collection of study material	Yes	Yes	Yes	Yes	Same
Segmenting of study material (Segmenting gums/bone)	Yes – auto- segmented with manual adjustment	Yes – manual selection	Yes – manual selection	Yes – auto- segmented with manual adjustment	Similar
Alignment of study material	Yes	Yes	Yes	Yes	Same
Measuring study material	Yes	Yes	Yes	Yes	Same
Measurement point selection	Automatic with manual adjustment	Manual selection	Manual selection	Automatic with manual adjustment	Similar
Analyzing study material	Yes	Yes	Yes	Yes	Same
Treatment simulation	Yes	Yes	Yes	N/A	Same
Treatment Options	Up to unlimited number of stages	No limit on number of stages	Up to unlimited number of stages	N/A	Similar
Virtual appliance design	Yes	Yes	Yes	Yes	Same
Surface scan for intraoral scanner	Yes	Yes	Yes	Yes	Same
Surface scan from STL file	Yes	Yes	Yes	Yes	Same
Ana	alysis and Treatment				
Arch shape	Yes	Yes	Yes	Yes	Same

Specification	Proposed Device: iOrtho	Predicate Device: 3Shape A/S Ortho System (K171634)	Reference Device: Wuxi EA Medical Instruments Technologies Limited. Clear Aligner (K203688)	Reference Device: Dentsply Sirona CEREC Ortho Software (K171122)	Comparison Result		
Overbite/overj et	Yes	Yes	No	Yes	Same as Predicate		
Occlusal map	Yes	Yes	Yes	Yes	Same		
3D treatment simulation	Yes	Yes	Yes	N/A	Same		
Orthodontic Appliance Search	Yes	Yes	Yes	N/A	Same		
Appliance virtual preparation	Yes	Yes	Yes	N/A	Same		
Orthodontic appliance design	Yes	Yes	Yes	N/A	Same		
Orthodontic appliance export	Yes	Yes	Yes	N/A	Same		
<u>Ma</u>	Managing Patient and Case Base Data						
Creating, editing, and copying patient data	Yes	Yes	Yes	Yes	Same		
Creating, editing, and copying case data	Yes	Yes	Yes	Yes	Same		

Comparison of Indications for Use to the Predicate and Reference Devices:

Based on the above comparison, the indications of use of the Proposed Device is nearly identical to that of the Predicate Device and similar to the Reference Device.

Based on the similarity of indication for use, the Proposed Device can be considered substantially equivalent to the Predicate Device.

Comparison of Technological Characteristics to Predicate and Reference Devices:

Based on the above comparisons, the design, construction, and performance characteristics of the Proposed Device are similar to that of the Predicate Device, and similar to the Reference Devices. The differences identified are not substantially different in operation of the device. Thus, the Proposed Device can be considered substantially equivalent to the Predicate Device.

Summary of Performance Data and Substantial Equivalence:

Utilizing FDA Guidance document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" (issued May 11, 2015) the Proposed Device, iOrtho underwent appropriate integration, verification, and validation testing. The software passed the testing and performed per its intended use.

The software has been designed, integrated, verified, and validated in accordance with IEC 62304-Medical device software – software life cycle processes.

Conclusion:

Based on comparison of indications for use, technological features, performance testing, and software validation testing, the Proposed Device, iOrtho, is substantially equivalent to the legally marketed Predicate Device, 3Shape Ortho System (K171634).