

August 31, 2023

Ki Mobility LLC Mark Murphy Vice President of Operations 5201 Woodward Drive Stevens Point, Wisconsin 54481

Re: K223527

Trade/Device Name: Little Wave Clik, Rogue XP, Little Wave XP, Spark Regulation Number: 21 CFR 890.3850 Regulation Name: Mechanical Wheelchair Regulatory Class: Class I, reserved Product Code: IOR Dated: August 29, 2023 Received: August 30, 2023

Dear Mark Murphy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Tushar Bansal -S

for Heather Dean, PhD Assistant Director, Acute Injury Devices Team DHT5B: Division of Neuromodulation and Physical Medicine Devices OHT5: Office of Neurological and Physical Medicine Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Indications for Use

510(k) Number *(if known)* K223527

Device Name

Little Wave Clik; Rogue XP; Little Wave XP; Spark;

Indications for Use (Describe)

The Ki Mobility Little Wave Clik, Rogue XP, Little Wave XP, and Spark manual wheelchairs are intended to provide mobility to pediatrics limited to a seating position.

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) #: K223527	510(k) Summary	Prepared on: 2023-01-11			
Contact Details		21 CFR 807.92(a)(1)			
Applicant Name	Ki Mobility LLC				
Applicant Address	5201 Woodward Drive Stevens Point WI 54481 United States				
Applicant Contact Telephone	(715) 303-6447				
Applicant Contact	Mr. Mark Murphy				
Applicant Contact Email	mmurphy@kimobility.com				
Device Name		21 CFR 807.92(a)(2)			
Device Trade Name	Little Wave Clik; Rogue XP; Little Wave XP; Spark				
Common Name	Mechanical wheelchair				
Classification Name	Wheelchair, Mechanical				
Regulation Number	890.3850				
Product Code	IOR				
Legally Marketed Predi	21 CFR 807.92(a)(3)				
Predicate # Predica	ate Trade Name (Primary Predicate is listed first)	Product Code			
K983639 Zippie	GS	IOR			

Device Description Summary

The Ki Mobility Little Wave Clik manual wheelchair is intended to provide mobility to pediatrics restricted to a sitting position. The Little Wave Clik manual wheelchair can be self-propelled by the occupant with access to hand-rims on the rear wheels or moved by an attendant with access to push handles. It is a rigid frame wheelchair with a standard weight capacity for all models of 165 pounds/75 kg. The Little Wave Clik wheelchair has a welded tubular high strength 1.125" aluminum frame with bend/trim configurations that vary the seat angle, knee angle, front and rear seat heights, seat depth, and caster wing offset. Frame offsets provide more room for the user's legs without increasing the frame length or reducing seat support.

21 CFR 807.92(a)(4)

The seat frame is dimpled to create a consistent indexing detent array for frame mounted components in 0.25" increments. An adjustable clamp is used to set the rear wheel center of gravity position relative to the backrest. A set of clamps and an aluminum cross tube comprise an adjustable cross strut to set the seat width. The strut position may be adjusted for optimal positioning for rigidity for any given seat depth. The front caster mount also contains a button to mate with dimples along the front tubes for vertical caster adjustment.

The Little Wave Clik has several components intended to grow with the pediatric user either with built in adjustment, or growth kits of separate parts. The Clik cross tubes and camber tubes come in growth kits with multiple accommodated chair widths. Growth of the chair is also accommodated from footplates that can handle multiple widths and are adjustable to fit footrest length.

The Little Wave Clik manual wheelchair has a folding backrest frame for ease of storage and transport. It utilizes discrete holes to adjust the backrest angle relative to the seat tube. It also utilizes an adjustable position backrest to adjust the seat depth.

The tubular aluminum frame accommodates both seat slings and removable depth adjustable aluminum seat pans that have overlapping portions that can grow in depth. The pans and slings are available for use with wheelchair seat cushions. The backrest frame accommodates either back upholstery or planar/contoured wheelchair backrests. For adjustment, the tension adjustable back upholstery has straps that accommodate multiple seat widths. Seat pans and slings have overlapping portions that can grow in depth. The frame is also designed to accommodate aftermarket seating.

The seat cushion is a separate medical device adapted for use to the Little Wave Clik and is the primary contact surface to the occupant along with the backrest upholstery. Little Wave Clik components such as armrests and footrests will also have contact to the occupant. The chair may be used as a seat in a motor vehicle and conforms to ISO 7176-19: 2022 as well as ANSI/RESNA WC-4 Section 19: 2017 requirements. Two sets of transit tie down loops are bolted to the frame.

The Clik XP and Clik XPe order forms are order form variations of the same device to meet different market's requirements for pricing and reimbursement. All can be configured identically.

The Ki Mobility Rogue XP manual wheelchair is intended to provide mobility to pediatrics restricted to a sitting position. The Rogue XP manual wheelchair can be self-propelled by the occupant with access to hand-rims on the rear wheels or moved by an attendant with access to push handles. It is a rigid frame wheelchair with a standard weight capacity for all models of 200 pounds/91 kg. The Rogue XP wheelchair has a welded tubular high strength 1.375" aluminum frame with bend/trim configurations that vary the seat angle, knee angle, front and rear seat heights, seat depth. Taper frames provide additional positioning of user's legs without the use of other components or accessories.

The Rogue XP frame has an index of holes along the sides for mounting the backrest. Additional hole can be provided behind the backrest to accommodate adjustable depth. An adjustable clamp is used to set the rear wheel center of gravity position relative to the backrest. The chair utilizes a welded cross strut to maximize the rigidity of the chair. The Rogue XP has a welded caster wing for stiffness with a knuckle housing. An internal pinion allows adjustment of caster angle forwards and backwards. The caster mount takes multiple fork and caster types, including a +3/4" stem option.

The Rogue XP has several components intended to grow with the pediatric user either with built in adjustment, or growth kits of separate parts. The cross tube weldment is split in the center with an insert the can grow the chair up to 3 inches in width the rigidizer bar on the backrest has a similar split with insert to grow the width of the chair. Growth of the chair is also accommodated from footplates that can handle multiple widths and are adjustable to fit footrest length.

The Rogue XP manual wheelchair has a folding backrest frame for ease of storage and transport. It utilizes discrete holes to adjust the backrest angle relative to the seat tube.

The tubular aluminum frame accommodates both seat slings and removable aluminum seat pans. The pans and slings are available for use with wheelchair seat cushions by adherence view AB tape. The backrest frame accommodates either back upholstery or planar/ contoured wheelchair backrests. For adjustment, the tension adjustable back upholstery has straps that accommodate multiple seat widths. The frame is also designed to accommodate aftermarket seating.

The seat cushion is a separate medical device adapted for use to the Rogue XP and is the primary contact surface to the occupant along with the backrest upholstery. Rogue XP components such as armrests and footrests will also have contact to the occupant.

The chair may be used as a seat in a motor vehicle and conforms to ISO 7176-19: 2022 as well as ANSI/RESNA WC-4 Section 19: 2017 requirements. Two sets of transit tie down loops are bolted to the frame.

The Rogue XP, Rogue XPe and Rogue XP TTL order forms are order form variations of the same device to meet different market's requirements for pricing and reimbursement. All can be configured identically.

The Ki Mobility Little Wave XP manual wheelchair is intended to provide mobility to pediatrics restricted to a sitting position. The Little Wave XP manual wheelchair can be self-propelled by the occupant with access to hand-rims on the rear wheels or moved by an attendant with access to push handles. It is a rigid frame wheelchair with a standard weight capacity for all models of 165 pounds/75 kg. The Little Wave XP wheelchair has a welded tubular high strength 1.375" aluminum frame with bend/trim configurations that vary the seat angle, knee angle, front and rear seat heights, seat depth, and front frame taper. Taper frames provide additional positioning of user's legs without the use of other components or accessories.

The Little Wave XP frame has an index of holes along the sides for mounting the backrest. Additional holes can be provided behind the backrest to accommodate adjustable depth. An adjustable clamp is used to set the rear wheel center of gravity position relative to the backrest. The chair utilizes a welded cross-strut to maximize the rigidity of the chair. The front caster mount is a two-part system with a clamp embedded that allows adjustment of caster angle front to back. The caster mount takes multiple fork and caster types, including $a + 3/4^{\prime\prime}$ stem option.

The Little Wave XP has several components intended to grow with the pediatric user either with built in adjustment, or growth kits of separate parts. The cross-tube weldment is split in the center with an insert that can grow the chair up to 3 inches in width. The rigidizer bar on the backrest has a similar split with insert to grow the width of the chair. Growth of the chair is also accommodated from footplates that can handle multiple widths and are adjustable to fit footrest length.

The Little Wave XP manual wheelchair has a folding backrest frame for ease of storage and transport. It utilizes discrete holes to adjust the backrest angle relative to the seat tube.

The tubular aluminum frame accommodates both seat slings and removable aluminum seat pans. The pans and slings are available for use with wheelchair seat cushions by adherence via AB tape. The backrest frame accommodates either back upholstery or planar/ contoured wheelchair backrests. For adjustment, the tension adjustable back upholstery has straps that accommodate multiple seat widths. The frame is also designed to accommodate aftermarket seating.

The seat cushion is a separate medical device adapted for use to the Little Wave XP and is the primary contact surface to the occupant along with the backrest upholstery. Little Wave XP components such as armrests and footrests will also have contact to the occupant.

The chair may be used as a seat in a motor vehicle and conforms to ISO 7176-19: 2022 as well as ANSI/RESNA WC-4 Section 19: 2017 requirements. Two sets of transit tie down loops are bolted to the frame. All available options and accessories are listed on the order form in the labeling section of this submission.

The Ki Mobility Spark manual wheelchair is intended to provide mobility to pediatrics restricted to a sitting position. The Spark manual wheelchair can be self-propelled by the occupant with access to handrims on the rear wheels or moved by an attendant with access to push handles. It is a folding frame wheelchair with a standard weight capacity for all models of 165 pounds/75 kg.

The Spark wheelchair has a welded tubular high strength aluminum frame with configurations that vary the front and rear seat heights, seat depth, rear wheel position and caster housing position. A width adjustable cross brace connects the wheelchair side frames and allows for folding of the frame for storage.

The Spark has several components intended to grow with the pediatric user with built in adjustment within a range of its initial configuration. Arm rests offer multiple options, including a pediatric t-arm with integrated side guard, standard t-arms, angle adjustable locking extendable flip up arm rests and tubular flip up arm rests. Footplates can handle multiple widths and are adjustable to fit footrest length.

The tubular aluminum frame accommodates aluminum seat pans for use with wheelchair seat cushions. The frame is designed to accommodate aftermarket seating systems.

The Spark manual wheelchair has a folding backrest frame for ease of storage and transport. It utilizes discrete holes to adjust the backrest angle relative to the seat tube and to adjust the seat depth.

Seat cushions and backrests are separate medical devices and are primary contact surfaces to the occupant. Spark components such as armrests and footrests will also have contact to the occupant.

The chair may be used as a seat in a motor vehicle and conforms to ISO 7176-19: 2022 as well as ANSI/RESNA WC-4 Section 19: 2017 requirements. Two sets of transit tie down loops are bolted to the frame.

Intended Use/Indications for Use

The Ki Mobility Little Wave Clik, Rogue XP, Little Wave XP, and Spark manual wheelchairs are intended to provide mobility to pediatrics limited to a seating position.

Indications for Use Comparison

The Ki Mobility Little Wave Clik, Rogue XP, Little Wave XP and Spark manual wheelchairs are intended to provide mobility to persons limited to a sitting position. The primary use is by pediatrics in need of manual wheeled mobility. This is consistent with the intended use of the Zippie GS.

Technological Comparison

The Ki Mobility Little Wave Clik manual wheelchair, in comparison to predicate device, shares intended use, similarity of configuration of components, wheels and accessories. The seat depth adjustment paradigms of both devices are similar, moving the backrest relative to the front frame. The Little Wave Clik utilizes an array of dimples to adjust the position of the backrest, whereas the Zippie GS utilizes discrete throughbolts. The center of gravity is adjustable on both products, with the Little Wave Clik utilizing an array of dimples, while the Zippie GS utilizes discrete through holes. The Ki Mobility Little Wave Clik uses a nonfolding welded rear frame to provide rigidity to the backrest. The Zippie GS utilizes a dual seat rail system to provide two attachment points for each back cane to provide rigidity. Footrests are separate from side frames. The Zippie GS offers both folding and nonfolding frames, while the Little Wave Clik only offers a nonfolding frame with integrated footrest tubes. Additionally, the Zippie GS offers a reversible frame, putting the large wheels at the front of the chair and the casters at the rear, while the Little Wave Clik only offers the conventional orientation with the large wheels in the rear. Different seat heights can be achieved on both products by interchanging wheels and casters, as well as by adjusting the height of the rear wheels on its support, or the casters within the caster forks. The backrest angle of both products is adjustable by use of an array of holes to set discrete seat to backrest angles. The backrests of both products fold down for stowage. The backrests of both products are either fixed height single tubes or a set of telescoping tubes to achieve the backrest heights desired. The Zippie GS provides a folding mechanism to fold the chair for stowage, while the Little Wave Clik does not fold laterally. However, the Little Wave Clik aluminum frame does not present a significant change to the safety or effectiveness of the wheelchair function and has been qualified through testing per the applicable standards. The Little Wave Clik components and accessories do not pose new risks, being as safe and effective as the predicate device. The Little Wave Clik is substantially equivalent to the predicate device in intended use, design, materials, and operating principles with no new issues of safety or effectiveness.

The Ki Mobility Rogue XP manual wheelchair, in comparison to predicate device, shares intended use, similarity of configuration of

21 CFR 807.92(a)(6)

21 CFR 807.92(a)(5)

21 CFR 807.92(a)(5)

IOR PREDICATE COMPARISON TABLE							
Device 510(k) Number	Ki Mobility Spark	Ki Mobility Rogue XP K223527	Ki Mobility Little Wave XP K223527	Ki Mobility Little Wave Clik K223527	Zippie GS K983639		
Intended Use	The Ki Mobility Spark manual wheelchair is intended to provide mobility to persons limited to a sitting position. The Ki Mobility Spark manual wheelchair is	The Ki Mobility Rogue XP manual wheelchair is intended to provide mobility to persons limited to a sitting position. The Ki Mobility Rogue XP manual wheelchair is	The Ki Mobility Little Wave XP manual wheelchair is intended to provide mobility to persons limited to a sitting position. The Ki Mobility Little Wave XP manual	The Ki Mobility Little Wave Clik manual wheelchair is intended to provide mobility to persons limited to a sitting position. The Ki Mobility Little Wave Clik manual	The Zippie GS manual wheelchair is intended t provide mobility to persons limited to a sitting position. The Zippie GS manual wheelchair is intended t		
Indications for Use	intended to provide mobility to pediatrics limited to a seating position.	intended to provide mobility to pediatrics limited to a seating position.	wheelchair is intended to provide mobility to pediatrics limited to a seating position.	wheelchair is intended to provide mobility to pediatrics limited to a seating position.	provide mobility to pediatric users limited to a sitting position.		
Design Weight Limit	Ki Mobility Spark 165lb/75kg (with or without transit)	Ki Mobility Rogue XP 200lb/91kg (with or without transit)	Ki Mobility Little Wave XP 165lb/75kg (with or without transit)	Ki Mobility Little Wave Clik 165lb/75kg (with or without transit)	Zippie GS 165lb/75kg (with or without transit)		
rame Material	Aluminum	Aluminum	Aluminum	Aluminum	Aluminum		
Seat Width Seat Depth	10"-16" 12"-18"	10" - 18" 10" - 20"	10" - 16" 10" - 16"	8" – 16" 8" – 16"	10" - 18" 10" - 18"		
Seat Height	(13"-21") Front Seat-to-Floor	13.5" - 20.5" (Front Seat-to-Floor)	11.5" - 20.5" (Front Seat-to-Floor)	11.5" – 20.5" (Front Seat-to-Floor)	14" - 20" (Front Seat-to-Floor)		
-	(11"-18.5") Rear Seat-to-Floor	11.5 – 20" (Rear Seat-to-Floor)	11.5 – 20" (Rear Seat-to-Floor)	11.5" – 20" (Rear Seat-to-Floor)	12.5"-20" (Rear Seat-to-Floor) 11.8" (Lateral Folding Frame Option) 20.5" (Ri		
Stowage Width	12.6" (Folded Laterally)	23.6"	22.6"	21.3"	Frame Option)		
Stowage Depth Stowage Height	29.5" 25.6"	31.5" 24.4"	28.9" 23.6"	28.9" 26.6"	29.5" 24.0"		
Fransit approved?	Yes (optional)	Yes (optional)	Yes (optional)	Yes (optional)	Yes (optional)		
Product weight Turning Radius	28.7lb 28.9″	26lb 26.4"	16lb 23.8"	31lb 20.7"	29lb 32.5"		
Frame Options:	Width Adjustable Folding Depth Adjustable Footrest Swing Away	Straight frame – non-folding Tapered frame – non-folding	Straight frame – non-folding Tapered frame non-folding	Standard – non folding Offset frame - non-folding	Folding Non-folding		
Back Angle Range	Footrest Super Low – Swing Away 85°-115°, 5° increments	80° - 100°, 2° increments	81° - 102°, 3° increments	80° - 100°, 2° increments	85° - 120°, 5° increments		
Sack Aligie Kalige	Stroller Handle Back Post	Height Adjustable	Height Adjustable	Height Adjustable	Height adjustable		
Backrest Styles	Height Adjustable Straight with Push Handle Dynamic Rocker Back	Fixed Height Fixed Height Offset	Fixed Height Fixed Height Offset	Fixed Height Fixed Height Offset	Fixed Height		
Rack Heights	Stroller Handle: 20"-24"	Adjustable height: 9"-18"	Adjustable height: 11"-18"	Height Adjustable	Fixed Height Stroller: 19", 22", 25"		
Back Heights	Height Adjustable: 13"-24"	Fixed Height: 9"-18" Fixed Height Offset: 11"-17"	Fixed Height: 9"-18" Fixed Height Offset: 11"-17"	Fixed Height Fixed Height Offset	Height Adj (Rigid): 16"-22" or 19"-25" Height Adj (Folding): 17"-23" or 20"-26"		
	Std Push Handle fixed on back canes	Std Push Handle fixed on back canes	Std Push Handle fixed on back canes	Std Push Handle fixed on back canes	Fixed push handles on back canes		
	Removeable Stroller Back Handle	Low Profile Push Handle fixed on back canes	Low Profile Push Handle fixed on back canes	Low Profile Push Handle fixed on back canes			
Push Handles		Folding push Handle	Folding push Handle	Folding push handle	Angle adjustable push handle (non-folding only		
Fusil Hallules		Bolt-On Push Handle	Bolt-On Push Handle	Bolt-on Push Handle			
		Ergo Stroller Handle Single Stroller Handle	Ergo Stroller Handle Single Stroller Handle	Ergo Stroller Handle Single Stroller Handle			
		Double Stroller Handle	Double Stroller Handle	Double Stroller Handle			
Rigidizer bars	None	Non-adjustable Height Rigidizer Bar Adjustable Height Rigidizer Bar	Non-adjustable Height Rigidizer Bar Adjustable Height Rigidizer Bar	Adjustable Height Rigidizer Bar	None		
	Single Dect Height Adjustable			Height Adjustable T-Arm (Tall, Std, low and	Single Part (ctd and Jaw)		
	Single Post Height Adjustable	Height Adjustable T-Arm (Tall, Std, low)	Height Adj. T-Arm (Std, low, tall and pediatric)	pediatric)	Single Post (std and low)		
	Height Adjustable T-Arm (Std, Tall, Low)	Angle Adj. Locking Extendable Flip Up Armrest	Tubular Flip Up Armrest	Tubular Flip Back Armrest	T-Post w/foam pad (regular, tall)		
Armrests	Angle Adjustable Locking Extendable Flip Up	Swing away armrest (short and tall)	Locking Extendable Flip Back armrest	Angle Adj. Locking Extendable Flip Back Armrest	Flip Up		
	Pediatric Height Adjustable T-Arm	Desk and Full Arm Pad – Standard and Waterfall	Swing away armrest	Swing away armrest (short and tall)	Adjustable Locking Flip up		
	Tubular Flip up	styles Foam Grip	Desk and Full length arm pads	Desk and Full length arm pads	Length Adjustable locking flip up		
	Arm Pads (Std, Waterfall, Foam Grip)	roam Grip	Foam grip type arm pads	Foam Grip	Length Aujustable locking hip up		
	Push to lock Pull to Lock	Push to Lock Pull to Lock	Push to Lock Pull to Lock	Push to Lock Pull to lock	Push to lock Pull to lock		
	Push to Lock Flush Mount	Push to Lock Flush Mount		Push to lock Flush Mount			
Wheel Lock	Hemi-Wheel Locks	Low Profile Scissor Lock Short Thro Scissor Lock Under Seat Scissor Lock	Low Profile Scissor lock Short Thro Scissor lock	Low Profile Scissor lock Short Thro Scissor lock			
		Drum Brake Dual Drum Brake w/push or pull	Drum brake Dual Drum brake w/push or pull	Drum brake Dual Drum brake w/push or pull	Foot Lock Hub Lock (12 and 16" wheels)		
	Extension handles (push and pull) Grade Aids	Extension handles (push and pull) Grade Aids	Extension handles (push and pull) Grade aids	Extension handles (push and pull) Grade aids	Extension handles (push and pull)		
Anti-Tips	Removable/swing-up rear Anti-tips	Removable/swing-up rear Anti-tips	Removable/swing-up rear Anti-tips	Removable/swing-up rear Anti-tips	Removable swing up rear anti-tips		
		Stainless Steel Spoke (18Ct.): 18", 20", 22", 24",			Front Anti-tips		
	Stainless Steel Spoke (18Ct.): 18", 20", 22", 24"	25" Mag: 20", 22", 24"		Stainless Steel Spoke (18 Ct.): 18", 20", 22", 24"' Mag: 20", 22", 24"	18", 20", 22" and 24" Spoke		
Rear Wheels	Mag: 20", 22", 24"		Mag: 20", 22", 24"		12, 16, 20 , 22 and 24 mag wheels		
	Composite Fiber Spoke (18Ct.): 22" 24"	Composite Fiber Spoke (18Ct.): 20" 22", 24", 25"		Composite Fiber Spoke (18 Ct.): 20", 22", 24" Composite Fiber Spoke (12 Ct.): 22", 24"			
Tires	Pneumatic Pneumatic with Airless Insert	Composite Fiber Spoke (12Ct.): 22", 24", 25" Pneumatic Pneumatic with Airless Insert	Composite Fiber Spoke (12 Ct.): 22", 24" Pneumatic Pneumatic with Airless Insert	Pneumatic Pneumatic with Airless Insert	Pneumatic Pneumatic with Airless Insert		
	Pneumatic with puncture protection Solid Polyurethane	Pneumatic with puncture protection Solid Polyurethane	Pneumatic with Puncture Protection Solid Polyurethane	Pneumatic with Puncture Protection Solid Polyurethane	Solid Polyurethane		
	Aluminum Anodized Aluminum with Plastic Coating	Aluminum Anodized Aluminum with Plastic Coating	Aluminum Anodized Aluminum with Plastic Coating	Aluminum Anodized Aluminum with Plastic Coating	Aluminum Aluminum Plastic Coated		
	Aluminum with Non-Slip Tape	Aluminum with Non-Slip Tape	Aluminum with Non-Slip Tape	Aluminum with Non-Slip Tape			
Handrims	Projection Knob with Plastic Coating Ergonomic Thumb Grip: Flexible Polymer	Projection Knob with Plastic Coating Ergonomic Thumb Grip: Flexible Polymer	Projection Knob with Plastic Coating Ergonomic Thumb Grip: Flexible Polymer	Projection Knob with Plastic Coating Ergonomic Thumb Grip: Flexible Polymer	Projection Knob		
	Ergonomic Thumb Grip: Alum with Grip Coating	Ergonomic Thumb Grip: Alum with Grip Coating	Ergonomic Thumb Grip: Aluminum with Grip	Ergonomic Thumb Grip: Alum with Grip Coating			
Camber	0°, 2°, 4°	0°, 2°, 3° [drive ready], 4°, 6°, 8°	Coating 0°, 2°, 3° [drive ready], 4°, 6°, 8°	0°, 2°, 3° [drive ready], 4°, 6°, 8°	0°		
Caster Wheels	.75" Roller Blade (3", 4")	0.75" x 3", 4"- Rollerblade	0.75" x 3", 4" - Roller Blade	0.75" x 3", 4" - Roller Blade	3", 4" and 5" lighted rollerblade		
(W × Ø)	.75" Lighted Roller Blade (3",4",5") 1" Polyurethane (4", 5",6", 7",8")	0.75" x 3", 4", 5"- Lighted Rollerblade 1" x 4", 5", 6" - Polyurethane	0.75" x 3", 4", 5" - Lighted Roller Blade 1" x 4", 5", 6" - Polyurethane	0.75" x 3", 4", 5" - Lighted Roller Blade 1" x 4", 5", 6" - Polyurethane	5" low poly, 5x1.5 semi-pneumatic 6" poly, 6" pneu, 6x1.5 semi-pneumatic		
	1" Polyurethane Aluminum (4", 5")	1" x 4", 5" - Polyurethane Aluminum 1" x 6" - Pneumatic	1" x 4", 5" - Polyurethane Aluminum 1" x 6" - Pneumatic	1" x 4", 5" - Polyurethane Aluminum 1" x 6" - Pneumatic	8" poly, 8" pneumatic, 8x1.5 Semi-pneu.		
	1" Pneumatic (6",8") 1.5" Polyurethane (4", 5",6")	1.5" x 4", 5", 6" - Polyurethane	1.5" x 4", 5", 6" - Polyurethane	1.5" x 4", 5", 6" - Polyurethane	8x2 pneu, 8x2 pneu/airless		
	1.5" Soft Roll Aluminum (3", 4", 5", 6") 2" Polyurethane (6") 2" Pneumatic w/ Foam Insert (8")	1.5" x 4", 5", 6" - Soft Roll Aluminum	1.5" x 4", 5", 6" - Soft Roll Aluminum	1.5" x 4", 5", 6" - Soft Roll Aluminum			
	Aluminum Fork	Aluminum Fork Single Sided Fork	Aluminum Fork Single Sided Fork	Aluminum Fork	Aluminum Fork		
	Fork Stem (Std, ¾" Longer, 1 ½" Longer)	Single Sided Fork	Single Sided Fork	Single Sided Fork			
Castor Fork-	Check Absorbing (Alumicum, Composite)	Shock Absorbing Fork: (Aluminum, Composite)	Shock Absorbing Fork: (Aluminum, Composite)	Shock Absorbing Fork: (Aluminum, Composite)			
Caster Forks	Shock Absorbing (Aluminum, Composite)				i i i i i i i i i i i i i i i i i i i		
Caster Forks	Shock Absorbing (Aluminum, Composite) Caster Pin Locks Housing (Std, Performance)						
Caster Forks	Caster Pin Locks Housing (Std, Performance) Swing Away, Footrest, Extension Mount	Integral to frame	Integral to frame	Integral to frame	Swing Away Footrest, Extension mount		
Caster Forks Footrest Hangers	Caster Pin Locks Housing (Std, Performance)	Integral to frame	Integral to frame	Integral to frame	Swing Away Footrest, Extension mount Swing Away Footrest, Front Mount Elevating Footrest (ext and front mount)		
	Caster Pin Locks Housing (Std, Performance) Swing Away, Footrest, Etension Mount Swing Away, Footrest, Front Mount Elevating Leg Rest (PRO, PRO Pediatric) Residual Limb Support	Integral to frame	Integral to frame	Integral to frame	Swing Away Footrest, Front Mount		
	Caster Pin Locks Housing (Std, Performance) Swing Away, Footrest, Extension Mount Swing Away, Footrest, Front Mount Elevating Leg Rest (PRO, PRO Pediatric)	Integral to frame Tubular w/o Footplate (Open)	Integral to frame Tubular without Footplate (Open)	Integral to frame	Swing Away Footrest, Front Mount Elevating Footrest (ext and front mount)		
	Caster Pin Locks Housing (Std, Performance) Swing Away, Footrest, Extension Mount Swing Away, Footrest, Front Mount Elevating Leg Rest (PRO, PRO Pedilatric) Residual Limb Support Hanger Release (Classic, 4 Way In)			Integral to frame	Swing Away Footrest, Front Mount Elevating Footrest (ext and front mount)		

1	Swing Away Extensions (Short, Med, Long)	Angle Adj. Flip-Under	Angle Adjustable Flip Under	Angle Adjustable Flip Under	Footplate flip-up
	Front Mount Extension	High Mount Angle Adj. Flip-Under	High Mount Angle Adjustable Flip Under	High Mount Angle Adjustable Flip Under	
Seating and positioning	Removable Aluminum Seat Pan	Solid Seat pan	Solid Seat pan	Depth Adjustable Seat Pan	Solid Seat/Back
		Standard Seat Sling	Standard Seat Sling	Standard Seat Sling	
	Seat Cushions	Seat Cushions	Seat Cushions	Seat Cushions	Jay Cushions
	Backrest	Backrests	Backrests	Backrest	Jay Back
	Lumbar pads	Lumbar pads	Lumbar pads	Lumbar pads	
	Shoulder guides	Shoulder guides	Shoulder guides	Shoulder guides	
	Headrest	Headrests	Headrests	Headrests	
	Chest supports	Chest supports	Chest supports	Chest supports	
	Supports (Chest and Trunk)				
	Laterals	Laterals	Laterals	Laterals	Toe Loop
	Positioning belts	Pelvic Positioning belts	Pelvic Positioning belts	Pelvic Positioning belts	Positioning belts
	Calf straps	Calf straps	Calf straps	Calf straps	Leg Strap
	Heel loops				Heel loops
Accessories	Ki Mobility Spark	Ki Mobility Rogue XP	Ki Mobility Little Wave XP	Ki Mobility Little Wave Clik	Zippie GS
Sideguards	Composite Side guards (adult and pediatric)	Composite Side guards (Std and Low)			
		Aluminum Side guards (adult and pediatric)	Aluminum Side guards (adult and pediatric)	Aluminum Side guards (adult and pediatric)	
Backpack/Pouch	Yes	Yes	Yes	Yes	Yes
Spoke Guards	Yes	Yes	Yes	Yes	Yes
Cane/Crutch Holder	Yes	Yes	Yes	Yes	No
ELR Gel Pads	Yes	No	No	No	No
Impact guards	Yes	Yes	Yes	Yes	Yes
Canopy	Yes	No	No	Yes	Yes
Tool Kit	Yes	Yes	Yes	Yes	Yes
One-Arm Drive	Yes	Yes	No	Yes	Yes

components, wheels and accessories. The seat depth adjustment paradigms of both devices are similar, moving the backrest relative to the front frame. Both the Rogue XP and the Zippie GS utilizes discrete bolts through the backrest into the side frames. The center of gravity is adjustable on both products, with the Rogue XP utilizing clamp on design, while the Zippie GS utilizes discrete bolt positions through side frames. The Ki Mobility Rogue XP uses a nonfolding welded rear frame to provide rigidity to the backrest. The Zippie GS utilizes a dual seat rail system to provide two attachment points for each back cane to provide rigidity. Footrests are separate from side frames. The Zippie GS offers both folding and nonfolding frames, while the Rogue XP only offers a nonfolding frame with integrated footrest tubes. Additionally, the Zippie GS offers a reversible frame, putting the large wheels at the front of the chair and the casters at the rear, while the Rogue XP only offers the conventional orientation with the large wheels in the rear. Different seat heights can be achieved on both products by interchanging wheels and casters, as well as by adjusting the height of the rear wheels on its support, or the casters within the caster forks. The backrest angle of both products is adjustable by use of an array of holes to set discrete seat to backrest angles. The backrests of both products fold down for stowage. The backrests of both products are either fixed height single tubes or a set of telescoping tubes to achieve the backrest heights desired.

The Zippie GS provides a folding mechanism to fold the chair for stowage, while the Rogue XP does not fold laterally. However, the Rogue XP aluminum frame does not present a significant change to the safety or effectiveness of the wheelchair function and has been qualified through testing per the applicable standards. The Rogue XP components and accessories do not pose new risks, being as safe and effective as the predicate device. The Rogue XP is substantially equivalent to the predicate device in intended use, design, materials, and operating principles with no new issues of safety or effectiveness.

The Ki Mobility Little Wave XP manual wheelchair, in comparison to predicate device, shares intended use,

similarity of configuration of components, wheels and accessories. The seat depth adjustment paradigms of both devices are similar, moving the backrest relative to the front frame. Both the Little Wave XP and the Zippie GS utilizes discrete bolts through the backrest into the side frames. The center of gravity is adjustable on both products, with both the Little Wave XP and Zippie GS utilizing discrete bolt positions through the tower clamps into the side frames. The Ki Mobility Little Wave XP uses a non-folding welded rear frame to provide rigidity to the backrest. The Zippie GS utilizes a dual seat rail system to provide two attachment points for each back cane to provide rigidity. Footrests are separate from side frames.

The Zippie GS offers both folding and non-folding frames, while the Little Wave XP only offers a non-folding frame with integrated footrest tubes. Additionally, the Zippie GS offers a reversible frame, putting the large wheels at the front of the chair and the casters at the rear, while the Little Wave XP only offers the conventional orientation with the large wheels in the rear. Different seat heights can be achieved on both products by interchanging wheels and casters, as well as by adjusting the height of the rear wheels on its support, or the casters within the caster forks. The backrest angle of both products is adjustable by use of an array of holes to set discrete seat-to-backrest angles. The backrests of both products fold down for stowage. The backrests of both products are either fixed height single tubes or a set of telescoping tubes to achieve the backrest heights desired. The Zippie GS provides a folding mechanism to fold the chair for stowage, while the Little Wave XP does not fold laterally. However, the Little Wave XP aluminum frame does not present a significant change to the safety or effectiveness of the wheelchair function and has been qualified through testing per the applicable standards. The Little Wave XP components and accessories do not pose new risks, being as safe and effective as the predicate device. The Little Wave XP is substantially equivalent to the predicate device in intended use, design, materials, and operating principles with no new issues of safety or effectiveness.

The Ki Mobility Spark manual wheelchair, in comparison to predicate device, shares intended use, similarity of configuration of components, wheels and accessories. The Zippie GS is available with either a foldable or non-foldable frame. The Ki Mobility Spark is only available as a foldable option. Both the Spark and Zippie GS fold by use of an x-brace that attaches to the frame to become rigid, and also have discrete holes that allow for the width of the chair to grow. Additionally, the Zippie GS offers a reversible frame, putting the large wheels at the front of the chair and the casters at the rear, while the Spark only offers the conventional orientation with the large wheels in the rear. The seat depth adjustment paradigms of both devices are similar, moving the backrest relative to the front frame. Both the Spark and Zippie GS utilize discrete through-bolts to adjust the position of the backrest. Different seat heights can be achieved on both products by interchanging wheels and casters, as well as by adjusting the height of the rear wheels on its support, or the casters within the caster forks. The backrest angle of both products is adjustable by use of an array of holes to set discrete seat-to-backrest angles. The backrests of both products fold down for storage. The backrests of both products are either fixed height single tubes or a set of telescoping tubes to achieve the backrest heights desired. The Spark components and accessories do not pose new risks, being as safe and effective as the predicate device. The Spark is substantially equivalent to the predicate device in intended use, design, materials and operating principles with no new issues of safety or effectiveness.

Non-Clinical and/or Clinical Tests Summary & Conclusions 21 CFR 807.92(b)

The Little Wave Clik, Rogue XP, Little Wave XP, and Spark have been tested to meet recognized consensus standards including: ANSI/RESNA WC-4:2017 Section 19: Wheelchairs used as seats in motor vehicles ISO-7176-1:2014 Determination of static stability ISO-7176-3:2012 Determination of Effectiveness of Brakes ISO 7176-5:2008 Determination of dimensions, mass, and maneuvering space ISO 7176-7:1998 Measurement of seating and wheel dimensions

ISO 7176-8:2014 Requirements and test methods for static, impact and fatigue strengths

ISO 7176-11:2012 Test Dummies

ISO 7176-13:1989 Determination of coefficient of friction of test surfaces

ISO 7176-15:1996 Requirements for information disclosure documentation and labeling

ISO 7176-16:2012 Resistance to ignition of postural support devices

ISO 7176-19:2022 Wheeled mobility devices for use as seats in motor vehicles

ISO 7176-22:2014 Set-up procedures

Design testing of Little Wave Clik, Rogue XP, Little Wave XP and Spark met or passed the recognized test standard requirements for manual wheelchairs, and which would apply to predicate devices.

Not Applicable

The Little Wave Clik, Rogue XP, Little Wave XP and Spark manual wheelchairs have the same intended use and similar technological characteristics as the predicate devices. The non-clinical testing to recognized standards exhibits that the device will perform as intended and risk analysis has documented risk reduction and identified requirements for labeling for safe and effective use of the device. The Little Wave Clik is substantially equivalent to the predicate devices as shown in the product design comparison.

The conclusion from testing and comparison to predicate devices demonstrates that the Little Wave Clik, Rogue XP, Little Wave XP and Spark are as safe, as effective, and performs as well as the legally marketed devices identified as predicate.