

August 12, 2023

Anhui Tianyang Pharmaceutical Co., Ltd. Zhang Shunlin Quality Manager 46 Tiantong Road, Tianchang City, Anhui Province Tianchang, Anhui China

Re: K223584

Trade/Device Name: Pre-Filled Normal Saline Flush Syringe

Regulation Number: 21 CFR 880.5200 Regulation Name: Intravascular Catheter

Regulatory Class: Class II

Product Code: NGT Dated: July 13, 2023 Received: July 14, 2023

#### Dear Zhang Shunlin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/efdocs/efpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/efdocs/efpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Bifeng Qian -S

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

## DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

#### **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K223584
Device Name Pre-Filled Normal Saline Flush Syringe
Indications for Use (Describe) Pre-Filled Normal Saline Flush Syringe is intended for use in flushing compatible intravenous administration sets and indwelling intravenous access devices. Use according to the recommendations of the manufacturer for the appropriate device.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

### 510(k) Summary (K223584)

I 510(k) Submitter

Device Submitter: Anhui Tianyang Pharmaceutical Co., LTD.

No.46, Tiantong Road, Tianchang City, Anhui Province

Contact Person: Zhang Sunlin

**Quality Manager** 

Phone: +86-15805509075 E-mail: 752322508@qq.com

II Device

Trade Name of Device: Pre-Filled Normal Saline Flush Syringe

Regulation Number: 21 CFR 880.5200

Classification Name: Saline, Vascular Access Flush

Product Code: NGT

Regulation Number: 21 CFR 880.5200

Regulatory Class II

Review Panel General Hospital

**III Predicate Devices** 

510k Number K213522

Trade Name of Device: AMSafe® Pre-Filled Normal Saline Flush Syringe

Regulation Number: 21 CFR 880.5200
Regulation Name: Intravascular Catheter

Regulatory Class II
Product Code: NGT

#### **IV Device Description**

Pre-Filled Normal Saline Flush Syringe is a polypropylene syringe filled with 0.9% sodium chloride for injection. It contains 3ml, 5ml and 10ml and consists of tip cap, barrel, piston, and plunger. This is a single use, disposable device(s), provided sterile.

#### V Indications for use

Pre-Filled Normal Saline Flush Syringe is intended for use in flushing compatible intravenous administration sets and indwelling intravenous access devices. Use according to the recommendations of the manufacturer for the appropriate device.

#### **VI Technological Characteristics Comparison**

VI-1: Comparison of Pre-Filled 0.9% Normal Saline Flush Syringe

Device		Predicate Device		
Characteristic	Subject Device (K223584)	(K213522)	Discussion	
	, ,	The AMSafe® PreFilled		
		Normal Saline Flush		
	Pre-Filled Normal Saline Flush	Syringe, is intended for use		
	Syringe is intended for use in	in flushing compatible		
	flushing compatible	intravenous		
Indications for	intravenous administration	administration sets and	Similar	
Use	sets and indwelling	indwelling intravenous	Comment 1	
	intravenous access devices.	access devices. Use		
	Use according to the	according to the		
	recommendations of the	recommendations of the		
	manufacturer for the	manufacturer for the		
	appropriate device.	appropriate device.		
Prescription/				
over-the counter	For Rx only	For Rx only	Identical	
use				
	The Pre-Filled Normal Saline			
	Flush Syringe is flushed after			
	the locking connector is	The AMSafe® PreFilled		
	connected to the medical	Normal Saline Flush Syringe		
	catheter connector for clinical	is a three-piece, sterile,		
Operation	use the liquid medicine (0.9%	single use syringe with a 6%	Identical	
Principle	sodium chloride injection) was	(Luer) connector pre-filled	Identical	
	pushed into the medical	with 0.9% Sodium Chloride		
	catheter and used to close	Injection, USP, and sealed		
	and flush the end of the	with a tip cap.		
	catheter in the gap between			
	different drugs.			
		The device has modified to		
	The subject device with Luer	add an extra thread to the		
Design	lock connection fitting and	plunger rod and inside of	Different	
Design	nonvented, female Luer lock	plunger stopper, the female	Comment 2	
	tip cap.	Luer cap has Changed to		
		screw type.		
Chemical	0.9% Sodium chloride	0.9% Sodium chloride	Identical	
composition	injection, USP	injection, USP	Idontioal	
	Tip Cap, Barrel and Plunger:	Barrel and plunger:		
Syringe material	polypropylene;	polypropylene	Identical	
- Jinigo material	Piston: Bromobutyl rubber	Stopper: Butyl rubber (not		
	. lotoni Bromosatyi idosoi	made with natural rubber		

Device Characteristic	Subject Device	Predicate Device (K213522)	Discussion
		latex)	
		Tip cap: polypropylene with	
		white colorant	
Syringe Size	Fill 3ml, 5ml, 10ml in 10cc	Fill 3ml, 5ml, 10ml in 10cc	
and Fill Volumes	syringe	syringe	Identical
and i iii voidines	Fill 3ml, 5ml in 5cc syringe	Fill 3ml, 5ml in 5cc syringe	
Syringe	BOPP heat sealing film	DD wron	Different
Packaging	BOFF fleat sealing filling	PP wrap	Comment 3
Sterilization method and SAL Level	Terminally sterilized by steam, 10 <sup>-6</sup> SAL	Terminally sterilized by steam, 10 <sup>-6</sup> SAL	Identical
Labeled	Yes	Yes	Identical
nonpyrogenic	162	162	iuerilical
Singe Use Only	Yes	Yes	Identical
Shelf Life	2 Years	3 Years	Similar

Comment 1: The Indications for Use of the predicate and subject device are same and the intended use is same, which does not affect the safety and effectiveness of the product.

Comment 2: The Luer cap of the predicate device has been modified, and the subject device is a standard Luer connector that meets the requirements of ISO 80369-7, so it does not affect the safety and effectiveness of the product.

Comment 3: The packaging materials of the subject and predicate are different, but the packaging of the subject device has been verified, and the sterility of the product can be guaranteed within the claimed shelf life of 2 years. Therefore, the safety and effectiveness of the product will not be affected.

#### VII Summary of Non-clinical Testing (Bench)

The non-clinical testing for Pre-Filled Normal Saline Flush Syringe was performed to demonstrate verification testing in conformance with the acceptance criteria of test methods and recognized consensus standards shown below. Real time aged samples from three non-consequtive lots were tested for all performance criteria.

Table VII-1: Performance testing was conducted on the subject device

ID#	Test	Method	Acceptance Criteria	Conclusion
1	1 Physical Testing of Syringe			

1.1	Lubricant	ISO7886-1	ISO7886-1	Pass	
1.2	Dead Space	ISO7886-1	ISO7886-1	Pass	
1.3	Limits for acidity or alkalinity	ISO7886-1	ISO7886-1	Pass	
1.4	Syringe Luer Performance	ISO 80369-7	ISO 80369-7	Pass	
1.5	Sealing performance	ISO7886-1	ISO7886-1	Pass	
2	Sodium Chloride I	njection, USP Testin	g		
2.1	pH value	USP<791>	PH: 4.5-7.0	Pass	
2.2	Oxidizable substance test	USP6-471	USP6-471	Pass	
2.3	Carbonate	USP<191>	USP<191>	Pass	
2.4	Sulfate	USP<191>	USP<191>	Pass	
2.5	Calcium	USP<191>	USP<191>	Pass	
2.6	Ammonium	USP<191>	USP<191>	Pass	
2.7	Iron test	USP<241>	< 2ppm	Pass	
	Limits of	USP<233>	USP<233>	Pass	
2.8	extractable metals	USP<232>	USP<232>		
3	Particulate	AAMI TIR42:2021	≥10µm, ≤6000	≥10µm, ≤361.5	
3	Contamination	AAIVII TIR42.2021	≥25µm, ≤600	≥25µm, ≤0.0	
4	Biocompatibility Testing				
4.1	Bacterial	LIOD 40 40Es	Bacterial endotoxins≤	Pass	
4.1	Endotoxins Test	USP 43<85>	0.5EU/mL		
4.2	In Vitro Cytotoxicity	ISO 10995-5:2009	Non-cytotoxic	Pass	
4.3	Intracutaneous Reactivity Test	ISO 10995-23:2021	Non-irritant	Pass	
4.4	Skin Sensitization Test	ISO 10993-10:2021	Non-sensitizer	Pass	
4.5	Acute Systemic Toxicity Test	ISO 10993-11:2017	No systemic toxicity	Pass	
4.6	Pyrogen Test	ISO 10993-11:2017	Non-pyrogen	Pass	
4.7	In Vitro Hemolysis Test	ISO 10993-4:2017	Non-hemolytic	Pass	

#### **VIII Clinical Test Conclusion**

No clinical study is included in this submission.

#### IX Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the subject Pre-Filled Normal

Saline Flush Syringe is as safe as effective, and performs as well as or better than the legally marketed device.

Date of Summary: August 8, 2023