

Wuhan Zonsen Medical Products Co., Ltd. Linna Ye Registration Manager No 8 Jinchao Rd, Zhucheng Street Wuhan, Hubei China

Re: K223600

Trade/Device Name: Sterilization Wrap Regulation Number: 21 CFR 880.6850 Regulation Name: Sterilization Wrap

Regulatory Class: Class II

Product Code: FRG Dated: July 7, 2023 Received: July 21, 2023

#### Dear Linna Ye:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

August 25, 2023

K223600 - Linna Ye Page 2

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Digitally signed by Eileen Cadel -S Cadel -S Date: 2023.08.25 13:20:06 -04'00'

for

Colin O'Neill, M.B.E. Assistant Director

DHT6B: Division of Spinal Devices OHT6: Office of Orthopedic Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

## **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

See PRA Statement below.

K223600
Device Name
Sterilization Wrap
Indications for Use (Describe)
Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a healthcare provider via the following:
A) Pre-vacuum Steam: 134°C for 4 minutes, 30 minutes dry time. B) Ethylene Oxide Sterilization: 100% Ethylene Oxide (EO) with a concentration of 360-363 mg/L @ 54°C and 40% relative humidity for 480 minutes, 24 hours aeration time at 42 °C
Model: ZSW1610, ZSW 1611.
Material: 45gms SMS, 54gms SMS.
Color of wrap: Pink/Blue (double layers).
Size (inches): 12*12, 15*15, 18*18, 20*20, 24*24, 30*30, 36*36, 40*40, 45*45, 48*48, 54*54, 72*54.

The wrap is intended to allow sterilization of the enclosed medical device(s) and maintain sterility of the enclosed device(s) up to 90 days post sterilization.

Types of medical devices to be sterilized in the Pre-Vacuum Steam Sterilization and Ethylene Oxide Sterilization: General purpose non-lumened reusable metal and nonmetal devices including devices with stainless steel diffusionr estricted spaces such as the hinged portion of forceps and scissors, as well as other general medical instruments etc. Models validated for both sterilization methods with stainless steel as part of the load, which with a maximum weight of each types, range from 2lbs to 9lbs, which dependent on each model's size.

For size2\*12, 15\*15, 18\*18, maximum weight is 2lb, for 20\*20, 24\*24, 30\*30, 36\*36, 40\*40, 45\*45, 48\*48, maximum weight is 6lb, and for 54\*54, 72\*54, maximum weight is 9lb.

CONTINUE ON A CERABATE DACE IS NEEDED			
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)		
Type of Use (Select one or both, as applicable)			
weight is 6lb, and for 54*54, 72*54, maximum weight is 9lb.			
$xy_0$ of the first $x_0$ and $x_0$			

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510(K) Summary

## 510 (K) Summary

This 510(K) Summary of safety and effectiveness information is being submitted in accordance with requirement of 21 CFR807.92

1. Date of Submission: August, 23, 2023.

## 2. Submitter / 510(K) Holder

## Wuhan Zonsen Medical Products Co., Ltd.

No. 8 Jinchao Road, Zhucheng Avenue, Xinzhou District, Wuhan City, Hubei Province, P.R. China 431400

Contact Person: Ms. Lisa Zhang Position: Registration Manager

Email: registration01@zonsenmed.com

# 3. Proposed Device Name

Trade name: Sterilization Wrap

Common name: Sterilization Wrap

Classification Name: Wrap, Sterilization/Indicator, Physical/Chemical Sterilization

**Process** 

Device Class: Class II

Classification Panel: General Hospital

Product Code: FRG

Regulation Number: 21 CFR 880.6850

#### 4. Predicate Devices

Primary Predicate Device: K160755

Product Name: Reliance® Solo Sterilization Wrap, Reliance® Tandem Sterilization

Wrap

Submitter: Ahlstrom Nonwovens LLC

# 5. Device Description

Sterilization Wrap is square or rectangular three-layer (SMS) non-woven sheet which manufactured with Polypropylene spunbond-meltblown-spunbond (SMS) fabric. Sterilization wrap provides a strong barrier which protects against cuts, tears with particularly device sets. The device is designed to be implemented as an outer sterilization wrap which allow sterilization of the enclosed medical device(s) and maintain sterility of the enclosed device(s) up to 90 days post Pre-Vacuum Steam Sterilization or Ethylene Oxide Sterilization.

Sterilization Wrap is two-layer sheet which ultrasonically bonded blue and pink non-woven sheet together along two edges, thus for convenient simultaneous wrapping.

Sterilization Wrap separated into two distinct types depending on different grams of weight of material, for ZSW 1610, it is two 45gsm SMS non-woven sheet ultrasonically bonded together, and for model ZSW 1611, it is two 54gsm SMS sheet. Each model contains series specification, including 12 Sizes (inches): 12\*12, 15\*15, 18\*18, 20\*20, 24\*24, 30\*30, 36\*36, 40\*40, 45\*45, 48\*48, 54\*54, 72\*54. With non-lumened reusable metal and nonmetal devices as part of the load, which with a maximum weight of each type, range from 2lbs to 9lbs, which dependent on each model' size (table 1)

Table 1 Device Specification

				Size		Maximum
Model	Material	Layers	Color	Long (inches)	Width (inches)	Weight(lbs.)
				12	12	2
			Pink and Blue double layers	15	15	2
	45 gsm			18	18	2
ZSW1610	SW1610 45 gsm SMS 2	2		20	20	6
		double	dodoic layers	24	24	6
				30	30	6
				36	36	6

ZSW1611							
ZSW1611 SMS 2 Pink and Blue double layers 2 Pink and Blue 48 48 6 54 9 9 72 54 9 12 12 2 15 15 15 2 18 18 18 2 20 20 6 24 24 6 54 24 6 54 54 54 9					40	40	6
ZSW1611 SMS 2 Pink and Blue double layers 2 Pink and Blue 440 40 6 45 45 45 6 54 54 9					45	45	6
ZSW1611 SMS 2 Pink and Blue double layers 2 Pink and Blue double layers 2 Pink and Blue double layers 36 36 6 40 40 6 45 45 6 48 48 6 54 54 9					48	48	6
ZSW1611					54	54	9
ZSW1611 SMS 2  Pink and Blue double layers 36 36 6 40 40 6 45 45 6 48 48 6 54 54 9					72	54	9
ZSW1611 SMS 2 Pink and Blue double layers 2 Pink and Blue double layers 36 36 6 40 40 6 45 45 6 48 48 6 54 54 9					12	12	2
ZSW1611 SMS 2 Pink and Blue double layers 20 20 6 24 24 6 30 30 6 40 40 6 45 45 45 6 48 48 6 54 54 9					15	15	2
ZSW1611 SMS 2 Pink and Blue double layers 24 24 6 30 30 6 40 40 6 45 45 45 6 48 48 6 54 54 9					18	18	2
ZSW1611					20	20	6
ZSW1611 SMS 2 double layers 36 36 6 40 40 6 45 45 6 48 48 6 54 54 9					24	24	6
SMS   double layers   36   36   6	7SW1611	54 gsm	2	Pink and Blue	30	30	6
45     45     6       48     48     6       54     54     9	Z5W1011	$ SMS ^2$	2	double layers	36	36	6
48     48     6       54     54     9					40	40	6
54 54 9					45	45	6
					48	48	6
72 54 9					54	54	9
					72	54	9

## 6. Indications for Use

Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a healthcare provider via the following:

- A) Pre-vacuum Steam: 134°C for 4 minutes, 30 minutes dry time.
- B) Ethylene Oxide Sterilization: 100% Ethylene Oxide (EO) with a concentration of 360~363 mg/L @ 54°C and 40% relative humidity for 480 minutes, 24 hours aeration time at 42 °C.

Model: ZSW1610, ZSW 1611.

Material (single layer): 45gms SMS, 54gms SMS.

Color of wrap: Pink/Blue.

Size (inches): 12\*12, 15\*15, 18\*18, 20\*20, 24\*24, 30\*30, 36\*36, 40\*40, 45\*45, 48\*48, 54\*54, 72\*54.

The wrap is intended to allow sterilization of the enclosed medical device(s) and maintain sterility of the enclosed device(s) up to 90 days post sterilization.

Types of medical devices to be sterilized in the Pre-Vacuum Steam Sterilization and Ethylene Oxide Sterilization:

General purpose non-lumened reusable metal and nonmetal devices including devices with stainless steel diffusionr estricted spaces such as the hinged portion of forceps and scissors, as well as other general medical instruments etc.

Models validated for both sterilization methods with stainless steel as part of the load, which with a maximum weight of each type, range from 2lbs to 9lbs, which dependent on each model's size.

For size: 2\*12, 15\*15, 18\*18, maximum weight is 2lb, for 20\*20, 24\*24, 30\*30, 36\*36, 40\*40, 45\*45, 48\*48, maximum weight is 6lb, and for 54\*54, 72\*54, maximum weight is 9lb.

## 7. Technological Characteristics Comparison Table

The Sterilization Wrap is compared with the predicate device Reliance® Solo Sterilization Wrap, Reliance® Tandem Sterilization Wrap (K160755). The results are shown below in the Technological Characteristics Comparison Table:

Table 2 General Comparison between proposed and predicate device

Device	Proposed Device	Predicate Device	Comparis
characteristic	Froposed Device	510(K) No. K160755	on Result
Manufacturer	Wuhan Zonsen Medical	Ahlstrom Nonwovens,	
Manufacturei	Products Co., Ltd.	LLC	
		Reliance® Solo	
		Sterilization	
Product name	Sterilization wrap	Wrap, Reliance®	
		Tandem	
		Sterilization Wrap	
510(k) Reference	K223600	K160755	

Class	Class II	Class II	Same
Product Code	FRG	FRG	Same
Regulation Number	21 CFR PART 880.6850	21 CFR PART 880.6850	Same
Intended Use	Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a healthcare provider.	intended to be used to enclose another medical device that is to be sterilized by a health care Provider.	Same
Indication for use	Sterilization Wrap is intended to be used to enclose another medical device that is to be sterilized by a healthcare provider via the following:  Pre-vacuum Steam: 134°C for 4 minutes, 30 minutes dry time.  B) 100% Ethylene Oxide (EO) with a concentration of 360~363 mg/L @ 54°C and 40% relative humidity for 480 minutes, 24 hours aeration time at 42 °C.	intended to be used to enclose another medical device that is to be sterilized by a healthcare provider via the following:  A) Pre-vacuum Steam 270°F/132°C for 4 minutes, validated for dry time of 20 minutes for Models T100/S100, T200/S200, T300/S300,	Similar

	B) 100% Ethylene Oxide
	(EO) with a concentration
	of 725~735 mg/L @
	131°F/55°C and
	40%~80% relative
	humidity for 60 minutes.
	C) Gravity Steam
	250°F/121°C for 30
	minutes
	D)Advanced Sterilization
	Products (ASP)
	STERRAD® 100S
	E) Advanced Sterilization
	Products (ASP)
	STERRAD® 100NX
	(Standard, Express and
	Flex cycles)
	F) STERIS Amsco® V-
	PRO 1, STERIS Amsco®
	V-PRO 1 Plus, STERIS
	Amsco® V-PRO maX
	Low Temperature
	Sterilization Systems
	(Lumen, Non-Lumen and
	Flexible Cycles)
The wrap is intended	to The wrap is intended to
allow sterilization of	the allow sterilization of the
enclosed medical device	e(s) enclosed medical Similar
and maintain sterility of	the device(s) and maintain
enclosed device(s) up to	90 sterility of the enclosed
1	6/17

K223000		310 (K)	Summary
	days post sterilization.	device(s) until used.	
	Polypropylene	Polypropylene	
Material	spunbondmeltblown-	spunbondmeltblown-	Somo
Composition	spunbond (SMS)	spunbond (SMS)	Same
	fabric	fabric	
		The Reliance® Solo and	
		Reliance® Tandem	
		Sterilization Wraps are	
		square or rectangular.	
		nonwoven sheets	
		produced using a three-	
		layer SMS (spunbond-	
	Zonsen Sterilization Wrap is	meltblown-spunbond)	
	square or rectangular three-	process. The Reliance®	
	layer (SMS) non-woven sheet which manufactured with Polypropylene spunbond-meltblown- spunbond (SMS) fabric. The device has two-layer	SMS Sterilization Wraps	
		are separated into two	
		distinct product offerings:	
		district product of the gen	
Design feature		Reliance® Solo: Consists	Same
Design feature		of two sheets of SMS	Same
	sheet which ultrasonically	wrap, ultrasonically	
	bonded blue and pink non-	bonded together along	
	woven sheet together along		
	two edges, thus for	two edges for convenient	
	convenient simultaneous	simultaneous wrapping of	
	wrapping.	one or a collection of	
		medical devices that will	
		be sterilized following	
		standard healthcare	
		practices.	
		Reliance® Tandem:	
		Consists of single sheets	

1122000			
		of SMS wrap, where two	
		sheets are used together	
		for the sequential	
		wrapping of one or a	
		collection of medical.	
		devices that will be	
		sterilized following	
		standard healthcare	
		practices.	
Size and color	depending on different grams of weight of material, Size (inches): 12*12, 15*15, 18*18, 20*20, 24*24, 30*30, 36*36, 40*40,	Reliance® Tandem and Reliance® Solo Wrap: six type Models depend on various grams of weight,	
Single Use vs.		34, 34, 00 00, 72 34.	
Reusable	Single use	Single use	Same
Disposable vs.			
Non-Disposable	Disposable	Disposable	Same
Maximum	Types of medical devices to	All models of Reliance®	
Recommended	be sterilized in the Pre-	Tandem and Solo	Similar
Wrapped Package	Vacuum Steam Sterilization		
wrapped rackage	vacuum Steam Stermzation	varidated for pre-vacuum	

Content Weights	and Ethylene Oxide	steam Sterilization and	
	Sterilization:	ethylene oxide	
	General purpose non-	sterilization with stainless	
	lumened reusable metal and	steel lumens (3 mm in	
	nonmetal devices including	diameter or larger and	
	devices with stainless steel	400 mm in length or less)	
	diffusionr estricted spaces	as part of the load with a	
	such as the hinged portion	maximum weight of	
	of forceps and scissors, as	25lbs dependent on the	
	well as other general	model.	
	medical instruments etc.		
	Models validated for both		
	sterilization methods with		
	stainless steel as part of the		
	load, which with a		
	maximum weight of each		
	type, range from 2lbs to		
	9lbs, which dependent on		
	each model's size.		
	ASTM D3776/ D3776M-20	ASTM D3776/ D3776M-	
	Size: Marking dimension	20	
	tolerance ±10%; Weight:	Size: Marking dimension	
Size and Gram	Marking gram tolerance	tolerance ±10%; Weight	Same
weight	±2g/m2	: Marking gram	Same
	the inspection result is	tolerance ±2g/m2	
	PASS	the inspection result is	
	1 700	PASS	
Bacterial	ASTM F2101-19 Bacterial	ASTM F2101-19	
Filtration	Filtration Efficiency (BFE)	Bacterial Filtration	Same
Efficiency (BFE)	Using a Biological Aerosol	Efficiency (BFE) Using a	

11223000			<u>Summar y</u>
	of Staphylococcus aureus ≥	Biological Aerosol of	
	90%	Staphylococcus aureus ≥	
	the inspection result is	90%	
	PASS	the inspection result is	
		PASS	
	AATCC 127-03, Water	AATCC 127-03, Water	
II14-4:-	Resistance: Hydrostatic	Resistance: Hydrostatic	
Hydrostatic	Pressure ≥ 50 cm H2O	Pressure ≥ 50 cm H2O	Same
Pressure	the inspection result is	the inspection result is	
	PASS	PASS	
	ACTM D5024 00 Ct 1 1	ASTM D5034-09	
	ASTM D5034-09 Standard,	Standard, Breaking	
	Breaking Strength and	Strength and Elongation	
	Elongation of Textile	of Textile Fabrics (Grab	
Tensile Strength	Fabrics (Grab Test)	Test)	Same
	MD≥ 90N; CD≥ 63N	MD≥ 90N; CD≥ 63N	
	the inspection result is	the inspection result is	
	PASS	PASS	
	ASTM D737-18Standard,	ASTM D737-18Standard,	
	Air Permeability of Textile	Air Permeability of	
Air Permeability	Fabrics ≥ 30 cfm	Textile Fabrics ≥ 30 cfm	Same
	the inspection result is	the inspection result is	
	PASS	PASS	
	ASTM D3786-18	ASTM D3786-18	
	Standard Test Method for	Standard Test Method for	
	Bursting Strength of Textile	Bursting Strength of	
	Fabrics—Diaphragm	Textile Fabrics—	Com-
Bursting Strength	Bursting Strength Tester	Diaphragm Bursting	Same
	Method ≥ 130 kpa	Strength Tester Method ≥	
	(18.86psi)	130 kpa (18.86psi)	
	the inspection result is	the inspection result is	

	PASS	PASS	
Tearing Strength	ASTM D5587-15, Tearing Strength of Fabrics, MD≥ 50N, CD≥ 30N. the inspection result is PASS	ASTM D5587-15, Tearing Strength of Fabrics, MD≥ 50N, CD≥ 30N. the inspection result is PASS	Same
Lint Generation	ISO9073-10, Textiles Lint, Cofficient of linting ≤4.0 the inspection result is PASS	ISO9073-10, Textiles Lint, Cofficient of linting ≤4.0 the inspection result is PASS	Same
Material Compatibility	After Pre-vacuum Steam or 100% Ethylene Oxide (EO) sterilization, the materials were not degraded.  Performance test including Size and Gram weight Bacterial Filtration Efficiency (BFE) Hydrostatic Pressure Tensile Strength Air Permeability Bursting Strength Tearing Strength Lint Generation can meet requirements	(EO) sterilization, the	Same
Biocompatibility	ISO 10993-5 non-cytotoxic post sterilization	ISO 10993-5, non- cytotoxic post	Same

		sterilization	
Maintenance of Sterility	90days	90days	Same
Shelf Life	2 years from date of manufacture	unknown	Similar

The design and technological characteristics of the Sterilization wrap is similar to the predicate chosen. There are minor differences between the devices including Indication for use (sterilization parameter), Size and color, Maximum Recommended Wrapped Package Content Weights,

## Note 1: Indication for use (sterilization parameter)

The proposed device and the predicate device have same sterilization method application, including Pre-vacuum Steam or 100% Ethylene Oxide (EO) sterilization. For those two devices, parameter recommend are different. However, the proposed device has been conducted verification in accordance with ISO11135 and ISO 17665, performance test has been done before and after sterilization, all results can meet standards requirements. Therefore, there are no different questions of safety and effectiveness questions pertaining to performance of the proposed device.

#### Note 2: Size and color

The proposed device and the predicate device have same dimention, but they are differences in color. However, the proposed device meets the requirements of the standard ISO 10993-5, articles after the steam and EO sterilized were used for the biocompatibility validation testing. Test result shows that the device has non-cytotoxicity. Therefore, the different technological specifications of the proposed device do not raise different questions of safety.

#### Note 3 Maximum Recommended Wrapped Package Content Weights

Although the maximum recommended wrap package loading are different from the predicate device, the proposed device has been conducted verification in accordance with ISO11135 and ISO 17665, Types of medical devices to be sterilized in the Pre-Vacuum

Steam Sterilization and Ethylene Oxide Sterilization including general purpose non-lumened reusable metal and nonmetal devices including devices with stainless steel diffusionr estricted spaces such as the hinged portion of forceps and scissors, as well as other general medical instruments etc. All models validated for both sterilization methods with stainless steel as part of the load, which with a maximum weight of each types, range from 2lbs to 9lbs, which dependent on each model's size.

Performance test has been done before and after sterilization, all results can meet standards requirements. Therefore, there are no different questions of safety and effectiveness questions pertaining to performance of the proposed device.

## Note 4 Performance test

Both the predicate device and the proposed device has been tested including, Size and Gram weight, Bacterial Filtration Efficiency (BFE), Hydrostatic Pressure, Tensile Strength, Air Permeability, Bursting Strength, Tearing Strength, Lint Generation. Result shows that there are no different questions of safety and effectiveness questions related to performance.

Based on above analysis mentioned above, we can conclude that subject device substantial equivalence with the predicate device.

#### 8. Performance Data

## Non-Clinical Performance Test Conclusion

Non-clinical tests were conducted to verify that the proposed device met all design specifications was Substantially Equivalent (SE) to the predicate Sterilization Wraps complies with the following standards:

Test Items	Standard and Acceptance Criteria	Test Result
Size and Gram	ASTM D3776M-20: Standard Test Methods for Mass Per Unit Area (Weight) of Fabric	
weight	Size: Marking dimension tolerance ±10%; Weight:	Pass
	Marking gram tolerance ±2g/m <sup>2</sup>	
Bacterial	ASTM F 2101- 19 Standard Test Method for Evaluating	Pass

Filtration	the Bacterial Filtration Efficiency (BFE) of Medical		
Efficiency	Face Mask Materials, Using a Biological Aerosol of		
(BFE)	Staphylococcus aureus		
	BFE ≥ 90%		
Hydrostatic	AATCC 127-03 Test Method for Water Resistance:		
	Hydrostatic Pressure	Pass	
Pressure	Hydrostatic Pressure ≥ 50 cm H <sub>2</sub> O		
Tensile	ASTM D5034-09 Standard Test Method for Breaking		
	Strength and Elongation of Textile Fabrics (Grab Test)	Pass	
Strength	MD≥ 90N; CD≥ 63N		
Air	ASTM D737-18Standard Test Method for Air Permeability		
	of Textile Fabrics	Pass	
Permeability	Air Permeability≥ 30 cfm		
	ASTM D3786-18 Standard Test Method for Bursting		
Bursting	Strength of Textile Fabrics—Diaphragm Bursting Strength	Daga	
Strength	Tester Method	Pass	
	Bursting Strength≥ 130 kpa (18.86psi)		
Tearing	ASTM D5587-15 Standard Test Method for Tearing Strength		
Strength	of Fabrics by Trapezoid Procedure	Pass	
Strength	MD≥ 50N; CD≥ 30N		
	ISO9073-10: 2019 Textiles — Test methods for nonwovens		
Lint Generation	— Part 10: Lint and other particles generation in the dry state	Pass	
	Coefficient of linting ≤4.0		
	ISO 10993- 1 Biological evaluation of medical devices -		
	Part 1: Evaluation and testing within a risk management		
Skin Irritation	process.		
AND Skin Sensitization	ISO 10993-5 Biological evaluation of medical devices - Part	Pass	
	5: Tests for in vitro cytotoxicity of medical devices.	1 433	
	ISO 10993-10:2010 Biological Evaluation of Medical		
	Devices – Part 10: Tests for irritation and skin sensitization.		
	Under the conditions of the study, the test article (pouches		

	and chemical indicator) extract did not show no significant	
	evidence of causing skin sensitization in the guinea pig	
	before and after sterilized.	
	Shelf Life Validation Test;	
Shelf-Life Validation	Conducted the shelf-life testing as real time aging method.  After the shelf life indicated as following,  Shelf Life: 2 Years; Shelf Life after Sterilized: 90 days;  The device performance shall be meet the requirements of the device. SAL=10 <sup>-6</sup>	Pass
Sterilization Process Validation for EO	EO Sterilization Process Validation Test.  Conducted the EO sterilization process validation as the method/principle of ISO11135: 2014 Sterilization of health care products - Ethylene oxide - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices (Healthcare Product Sterilization - Ethylene Oxide-Medical Device Sterilization Process Development. Validation and Routine Control Requirements).  Use the half cycle method to validate the EO sterilization cycle claimed in indication for use is effective.  Ethylene Oxide Sterilization: 100% Ethylene Oxide (EO) with a concentration of 360-363 mg/L @ 54°C and 40% relative humidity for 480 minutes, 24 hours aeration time at 42 °C  The device performance shall be meet the requirements of the device with the worst loading, SAL=10-6	Pass

Sterilization	Steam Sterilization Process Validation Test.	Pass
Process	Conducted the Steam sterilization process validation as the	
Validation for	method/principle of AAMI / ANSI ST79: 2010 & A1:2010 &	
steam	A2:2011 & A3:2012 & A4:2013, Comprehensive Guide To	
	Steam Sterilization And Sterility Assurance In Health Care	
	Facilities and ISO 17665- 1:2006 Sterilization of health care	
	products Moist heat Part 1: Requirements for the	
	development, validation and routine control of a	
	sterilization process for medical devices 's requirements.	
	Pre-vacuum Steam: 134°C for 4 minutes, 30 minutes dry	
	time.	
	The device performance shall be meet the requirements of the	
	device with the worst loading, SAL=10 <sup>-6</sup>	
Material	After Pre-vacuum Steam or 100% Ethylene Oxide (EO)	Pass
compatibility	sterilization, the materials were not degraded.	
	Performance test including:	
	Size and Gram weight	
	Bacterial Filtration Efficiency (BFE)	
	Hydrostatic Pressure	
	Tensile Strength	
	Air Permeability	
	Bursting Strength	
	Tearing Strength	
	Lint Generation	
	should meet requirements.	
EO/ECH residue	ISO 10993-7: 2008 Biological evaluation of medical devices	Pass
	— Part 7: Ethylene oxide sterilization residuals	
	EO≤4mg/d, ECH ≤9mg/d	
	1	

Non-clinical tests including physical characteristics shown in table above, shelf life validation, biocompatibility evaluation, EO Sterilization Process Validation and Steam

Sterilization Process Validation has been taken to shown all results can meet the standards. Aged device were used for material compatibility research, after Pre-vacuum Steam or 100% Ethylene Oxide (EO) sterilization, the materials were not degraded and all physical characteristics test result can still meet requirements.

## **Clinical Test Conclusion**

No clinical study is included in this submission.

## 9. Conclusion

The conclusions drawn from the non-clinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as the predicate device K160755.

17/17