

March 17, 2023

Inari Medical, Inc Ellen Nguyen Regulatory Affairs Specialist 6001 Oak Canyon, Suite 100 Irvine, California 92618

Re: K223613

Trade/Device Name: InThrill Thrombectomy System

Regulation Number: 21 CFR 870.5150 Regulation Name: Embolectomy Catheter

Regulatory Class: Class II Product Code: QEW, KRA Dated: February 14, 2023 Received: February 14, 2023

Dear Ellen Nguyen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Gregory W. Digitally signed by Gregory W. O'connell -S
O'connell -S
Date: 2023.03.17
13:13:36-04'00'

Gregory O'Connell
Assistant Director
DHT2C: Division of Coronary
and Peripheral Intervention Devices
OHT2: Office of Cardiovascular Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K223013
Device Name InThrill Thrombectomy System
Indications for Use (Describe) The InThrill Thrombectomy System is indicated for:
- The non-surgical removal of thrombi and emboli from blood vessels, including arteriovenous fistulae and arteriovenous grafts for dialysis access, and synthetic grafts. - Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel/graft.
The InThrill Thrombectomy System is intended for use in the peripheral vasculature.
The InThrill Thrombectomy System is not intended for use in deep vein thrombosis (DVT) treatment.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY

Date prepared March 16, 2023

Name Inari Medical, Inc.

6001 Oak Canyon, Suite 100

Irvine, CA 92618 877.923.4747

Contact person Ellen Nguyen

Regulatory Affairs Specialist

Trade name InThrill Thrombectomy System

Common name Embolectomy catheter

Regulation name Embolectomy catheter

Classification number 21 CFR 870.5150

Product code QEW
Secondary product code KRA
Regulatory class II

Predicate device(s) Capture Vascular, MegaVac Mechanical Thrombectomy System (K171493)

Hotspur Technologies, Inc., PTA-Plus PTA Balloon Catheter (K100842) Inari Medical, Mini-ClotTriever Thrombectomy System (K220887)

Reference device(s) Intramed Laboratories, Inc., Graft Thrombectomy Instruments (K942457)

Rex Medical, Cleaner Rotational Thrombectomy System (K141617)

Description

The InThrillTM Thrombectomy System is a single-use, sterile medical device designed to remove thrombi and emboli from the peripheral vasculature. The InThrillTM Thrombectomy System consists of the InThrill Sheath ("Sheath") and the InThrill Thrombectomy Catheter ("Catheter"), each packaged separately.

The Sheath is an introducer sheath with a distal self-expanding funnel, proximal aspiration port, and proximal hub. A Dilator is provided to aid insertion and positioning of the Sheath.

Radiopaque markers aid Sheath positioning under fluoroscopic visualization. The Sheath and Dilator tips are radiopaque, and there is a radiopaque marker band at the proximal end of the Sheath funnel.

Indications for Use

The InThrill Thrombectomy System is indicated for:

- The non-surgical removal of thrombi and emboli from blood vessels, including arteriovenous fistulae and arteriovenous grafts for dialysis access, and synthetic grafts.
- Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel/graft.

The InThrill Thrombectomy System is intended for use in the peripheral vasculature.

The InThrill Thrombectomy System is not intended for use in deep vein thrombosis (DVT) treatment.

Summary of substantial equivalence

A tabular comparison of the predicate and subject devices is provided below:

	Subject Device InThrill	Primary Predicate MegaVac Mechanical	Predicate PTA-Plus PTA	Predicate Mini-ClotTriever
	Thrombectomy System	Thrombectomy System	Balloon Cather 5mm X 4cm, PTA-Plus PTA Balloon Catheter 6mm X 4cm	Thrombectomy System
K Number	K223613	K171493	K100842	K220887
Manufacturer	Inari Medical, Inc.	Capture Vascular, Inc.	Hotspur Technologies, Inc.	Inari Medical, Inc.
Regulations	21 CFR 870.5150 Embolectomy catheter	21 CFR 870.5150 Embolectomy catheter	21 CFR 870.5150 Embolectomy catheter 21 CFR 870.1250 Percutaneous catheter	21 CFR 870.5150 Embolectomy catheter
Product Code	• QEW • KRA	 QEW QEX DXE	DXE LIT	• QEW
Indications for Use	The InThrill Thrombectomy System is indicated for: The non-surgical removal of thrombi and emboli from blood vessels, including arteriovenous fistulae and arteriovenous grafts for dialysis access, and synthetic grafts. Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel/graft. The InThrill Thrombectomy System is intended for use in the peripheral vasculature.	The MegaVac Mechanical Thrombectomy System is indicated for: The non-surgical removal of emboli and thrombi from blood vessels. The non-surgical removal of thrombi from synthetic grafts. Use in temporary blood vessel/graft occlusion. Injection, infusion, and/or aspiration of contrast media and other fluids into or from a vessel/graft Catheter placement over a guidewire	The PTA-Plus PTA Balloon Catheter is indicated for use within synthetic arteriovenous dialysis fistulae to remove embolic material (thrombus/debris) and dilate stenosis for treatment of obstructive lesions.	The Mini-ClotTriever Thrombectomy System is indicated for: The non-surgical removal of emboli and thrombi from blood vessels. Injection, infusion, and/or aspiration of contrast media and other fluids into or from a blood vessel. The Mini-ClotTriever Thrombectomy System is intended for use in the peripheral vasculature. The Mini-ClotTriever Thrombectomy System is not intended for use in deep vein thrombosis (DVT) treatment.

PAGE 3 OF 6

	The InThrill			
	Thrombectomy			
	System is not intended for use in			
	deep vein thrombosis			
	(DVT) treatment.			
Device	The InThrill	The MegaVac	The PTA-Plus PTA	The Mini-ClotTriever
Description	Thrombectomy	Mechanical	Balloon Catheter is	Thrombectomy
2 cocraption	System is a single-	Thrombectomy	designed for de-	System is a single-
	use, over-the-wire,	System is a single-	clotting and treating	use, over-the-wire,
	catheter-based system	use, over-the-wire,	stenosis in synthetic	catheter-based system
	for the minimally	catheter-based system	dialysis fistulae. It is	for the minimally
	invasive treatment of	for intravascular	an .035" guide-wire	invasive treatment of
	thromboemboli	mechanical	compatible, PTA	thromboemboli in the
	(including organized	thrombectomy,	balloon catheter with	peripheral
	thrombus and	occlusion, aspiration	a proprietary valve	vasculature. The
	adherent thrombotic	and embolectomy in	system which allows	system comprises
	material) in the	the peripheral	injection of contrast	two main
	peripheral	vasculature,	and an embolectomy	components
	vasculature,	including grafts. The	coil for mechanical	packaged separately:
	including	system comprises	removal of thrombus.	Mini-ClotTriever
	arteriovenous fistulae	two main		Sheath (8 Fr)
	and arteriovenous	components:		• Mini-ClotTriever
	grafts for dialysis	MegaVac Catheter		Catheter (8 Fr)
	access, and synthetic	• ThromboWire clot		The Mini-ClotTriever
	grafts. The system comprises two main	retractor		Sheath is comprised of reinforced
	components	The MegaVac		polymeric coaxial
	packaged separately:	Catheter with		sheath shafts
	• InThrill Sheath (8	SafeSeal technology utilizes a silicone		equipped with a self-
	Fr)	coated nitinol braided		expanding distal
	• InThrill	funnel that expands		mesh funnel, a
	Thrombectomy	to occlude antegrade		flush/aspiration port,
	Catheter (8 Fr)	blood flow proximal		and a proximal
	The InThrill Sheath is	to the target work		hemostatic valve. The
	comprised of	zone creating a static		Mini-ClotTriever
	reinforced polymeric	environment in which		Catheter is comprised
	coaxial sheath shafts	to perform the		of reinforced
	equipped with a self-	intervention, while		polymeric coaxial
	expanding distal	also centering and		shafts terminating in
	mesh funnel, a	securing the catheter		an expandable nitinol
	flush/aspiration port,	tip position within the		coring element
	and a proximal hemostatic valve. The	vessel. The MegaVac		(basket). Other
	InThrill	catheter's large-		accessories provided include a pre-dilator
	Thrombectomy	mouth funnel and		and dilator.
	Catheter is comprised	inner diameter allows		and unawr.
	of reinforced	for strong aspiration		
	polymeric coaxial	while easily being able to pass the		
	shafts terminating in	ThromboWire and		
	an expandable nitinol	other lesion		
	coring element	disruptive or		
	(basket). Other	therapeutic devices		
	accessories provided	through it. The		
	1		1	<u> </u>

PAGE 4 OF 6

Target Vessel	include a pre-dilator and dilator. Peripheral vessels (4-	ThromboWire consists of a nitinol embolectomy element that when expanded by the proximal actuation handle can serve to gather and pull matter towards and through the MegaVac catheter during aspiration. Peripheral vessels	Peripheral vessels	Peripheral vessels (4-
	 10 mm) that include: native vessels arteriovenous fistulae arteriovenous grafts 	that include: native vesselssynthetic grafts	that include:	10 mm) that include: • native vessels
	• synthetic grafts		synthetic grants	
Sterility	SAL 10 ⁻⁶ , EO	SAL 10 ⁻⁶ , EO	SAL 10 ⁻⁶ , EO	SAL 10 ⁻⁶ , EO
Shelf-life	2 years	Unknown	Unknown	2 years *
Guidewire	0.035"	Up to 0.044" OD	0.035"	0.035"
compatibility	0.033	Ор 10 0.044 ОД	0.033	0.033
Sheath	Outer shaft: 0.154"	Up to 9 Fr OD	Unknown	Outer shaft: 0.154"
Dimensions	OD/0.137" ID	Орюэнов	Clikilowii	OD/0.137" ID
	Inner shaft: 0.137" OD/0.110" ID Length: 6 cm	Length: up to 300 cm		Inner shaft: 0.137" OD/0.110" ID Length: 6 cm
Shaft Materials	Inner: Pebax 55D, ProPell, PTFE Liner, Stainless steel coil, Radiopaque marker band Outer: Pebax 63D, ProPell, PTFE Liner	Unknown	Unknown	Inner: Pebax 55D, ProPell, PTFE Liner, Stainless steel coil, Radiopaque marker band Outer: Pebax 40D and 25D, ProPell, PTFE Liner
Hemostasis	8 Fr Garrote valve	Unknown	Unknown	8 Fr Garrote valve
Valve	Rotating swivel hub with side port			Rotating swivel hub with side port
Handle	Yes; Slide actuator enclosed within handle housings	Yes; Operator controlled – linear actuation	Unknown	Yes; Slide actuator enclosed within handle housings
Braided Funnel	OD: 10 mm Length: 0.70" Braided nitinol funnel	OD: 2-12 mm Silicone coated braided nitinol funnel	Unknown	OD: 10 mm Length: 0.70" Braided nitinol funnel
Side Port	Tygon tubing with Loctite adhesive 1-way stopcock with female Luer connector	Unknown	Unknown	Tygon tubing with Loctite adhesive 1-way stopcock with female Luer connector

PAGE 5 OF 6

Dilator	OD: 0.110"	Unknown	Unknown	OD: 0.110"
	Total length: 8.8"			Total length: 8.3"
	Tipped LDPE/HDPE			Tipped LDPE/HDPE
	extrusion			extrusion
	Dilator cap			Dilator cap
	Proximal flush port			Proximal flush port
Pre-Dilator	OD: 0.13" (10 Fr)	Unknown	Unknown	OD: 0.13" (10 Fr)
	Polypropylene,			Polypropylene,
	HDPE			HDPE
	Length: 10.2 cm			Length: 10.2 cm
Catheter	OD: 0.111"	Nitinol	Unknown	OD: 0.111"
	Proximal hub with			Proximal hub with
	Tuohy Borst			Tuohy Borst
	hemostasis Y-valve			hemostasis Y-valve
	and 1-way stopcock			and 1-way stopcock
	Materials:			Materials:
	Outer: PTFE			Outer: PTFE
	Liner, SS304V			Liner, SS304V
	Braid, Radiopaque			Braid,
	marker band, 63D			Radiopaque
	Pebax Jacket, 63D			marker band, 63D
	Pebax Fluoro-safe			Pebax Jacket, 63D
	marker band			Pebax Fluoro-safe
	Middle: Braided			marker band
	polyimide,			Middle: Braided
	Radiopaque			polyimide,
	marker band,			Radiopaque
	Pebax 72D			marker band,
	Inner: Braided			Loctite 3942
	polyimide,			<u>Inner</u> : Braided
	Radiopaque 55D			polyimide,
	Pebax tip with			Radiopaque 55D
	ProPell			Pebax tip with
				ProPell
Length	65 cm	Up to 300 cm	55 cm	65 cm
Coring Element	Laser-cut nitinol	Woven nitinol	Unknown	Laser-cut nitinol
	OD: 18 mm	OD: 2-9 mm		OD: 18 mm
	Length: 88 mm			Length: 88 mm

Summary of substantial equivalence

Biocompatibility

The material changes proposed in this submission have no impact on the established biocompatibility of the device. Therefore, the previous passing results demonstrating that the InThrill Thrombectomy System (K220887) and accessories meet biological safety requirements per ISO 10993-1 are still applicable.

Sterilization

The subject device, including its accessories, is sterilized using EtO to achieve a sterility assurance level (SAL) of 10^{-6} using a validated sterilization process in accordance with the principles of ISO 11135:2014/Amd 1:2018 and AAMI TIR 28:2016. There have been no changes proposed that would affect device

sterilization; therefore, the previous sterilization process per K220887 remains applicable.

Non-Clinical Testing

In accordance with the Design Failure Modes and Effects Analysis, verification and validation testing were identified to support the substantial equivalence of the InThrill Thrombectomy System to the predicate devices. These tests included:

Performance Tests

- Graft Leak Testing (Pre-Simulated Use)
- In-Graft Simulated Use, InThrill Thrombectomy System
- Graft Leak Testing (Post-Simulated Use)
- Graft Visual Inspection (Post-Simulated Use)
- Radial Force Testing
- Prescale Contact Paper Pressure Testing
- Graft Abrasion Testing
- Comparative Adherent Clot Removal Testing
- Visual Inspection of Vein-to-Graft Anastomosis (Pre-Simulated Use)
- Simulated Use, InThrill Catheter through Vein-to-Graft Anastomosis
- Visual Inspection of Vein-to-Graft Anastomosis (Post-Simulated Use)

Test results demonstrated that all acceptance criteria were met; therefore, the device conforms to established product specifications.

Pre-Clinical Study

A GLP Animal study was performed, which did not identify any new questions of safety or effectiveness for the InThrill Thrombectomy System in peripheral vessels within the indicated size range.

Clinical testing was not required for the determination of substantial equivalence.

Conclusion

The non-clinical tests demonstrate that the subject device does not raise new questions of safety or effectiveness and, therefore, is substantially equivalent to the predicate device.