

August 2, 2023

Advance Medical Designs, Inc. % David Mackie QA/RA Manager 1241 Atlanta Industrial Drive MARIETTA GA 30066

Re: K223689

Trade/Device Name: Disposable Needle Guides and Grids

Regulation Number: 21 CFR 892.1570

Regulation Name: Diagnostic Ultrasonic Transducer

Regulatory Class: Class II

Product Code: ITX Dated: June 21, 2023 Received: July 3, 2023

Dear David Mackie:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

K223689 - David Mackie Page 2

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Yanna S. Kang -S

Yanna Kang, Ph.D.
Assistant Director
Mammography and Ultrasound Team
DHT8C: Division of Radiological Imaging
and Radiation Therapy Devices
OHT8: Office of Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023 See PRA Statement below.

K223689
Device Name
Disposable Needle Guides and Grids
Indications for Use (<i>Describe</i>) "Disposable Needle Guides and Grids are used to assist and aid physicians in performing an endocavity diagnosis
ultrasound needle guided procedure using guided intervention by providing fixed guiding for the precise insertion of linear instruments, such as needles. The Needle Guides and Grids are designed to aid adult patient population, in need of a biopsy of an internal organ, or internal delivery or removal of fluid within the body cavity, via the use of a needle, during an ultrasound procedure by retaining the needle tip and barrel within the ultrasound beam."
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Advance Medical Designs, Inc.

Common Name: Needle Guides & Grids

This summary of Traditional 510(k) Submission of safety and effectiveness of information is being prepared in accordance with the requirements of 21 CFR 807.92.

April 19, 2023

Assigned 510(k) NUMBER: <u>K223689</u>.

1. Submitter's Identifications:

510(k) owner's name: Advance Medical Designs, Inc.

Owner/Operator #: 1037885

Establishment: Advance Medical Designs, Inc. **Address:** 1241 Atlanta Industrial Drive

Marietta, GA 30066 USA

Phone Number: (770) 422-3125

Facility Registration #: 1037885

Contact person: David Mackie: QA/RA Manager at Advance Medical Designs, Inc.

Phone Number: (770) 422-3125 ext. 244
e-mail: mackied@advmeddes.com

2. <u>Date 510(k) Summary Prepared</u>: April 19, 2023

3. Name of Subject Device and Classification Information:

Trade name: Disposable Needle Guides and Grids

Regulation Number: 21 CFR 892.1570

510(k) Number: K223689

Common Name: Needle Guides & Grids

Classified Name: Diagnostic ultrasonic transducer accessories

Regulatory Class: Class II Product Code: ITX

4. Information for the Predicate Device:

A) PRIMARY PREDICATE DEVICE:

Trade Name/Device Name: Disposable Endocavity Needle / Biopsy Guide

Manufacturer: CIVCO Medical Instruments Co., Inc.

510(k) Number: K972514

Regulation Number: 21. CFR 892.1570

Classification Name: Diagnostic ultrasonic transducer accessories

Regulatory Class: Class II Product Code: ITX

5. Information for Reference Devices:

A) REFERENCE DEVICE:

Trade Name/Device Name: Disposable Guides KDNG00

Manufacturer: KOELIS 510(k) Number: K180970

Regulation Number: 21. CFR 892.1570

Classification Name: Diagnostic ultrasonic transducer accessories

Regulatory Class: Class II Product Code: ITX

B) REFERENCE DEVICE:

Trade Name/Device Name: VitroPRO / Disposable Endocavity Needle Guide

Manufacturer: CIVCO Medical Instruments Co., Inc.

510(k) Number: K222052

Regulation Number: 21. CFR 892.1570

Classification Name: Diagnostic ultrasonic transducer accessories

Regulatory Class: Class II Product Code: ITX

C) REFERENCE DEVICE:

Trade Name/Device Name: Reusable Guide

Manufacturer: KOELIS 510(k) Number: K141334

Regulation Number: 21. CFR 892.1570

Classification Name: Diagnostic ultrasonic transducer accessories

Regulatory Class: Class II Product Code: ITX

Description of Subject Device:

Advance Medical Designs' disposable needle guides & grids devices used to direct needles or instruments along a fixed path to a target location with an ultrasound traducer. They are provided in a variety of sizes to fit different equipment and situations. The Needle Guides and Grids are packaged separately, or within kits, provided sterile, and are labelled as single use only. Each disposable needle guide & grid contains a bracket and needle adapter. Each kit includes a disposable needle guide, a 20ml packet of ASonic® sterile gel, two latex free elastic bands, and a transducer cover (Latex or Latex-Free). The needle guides are non-invasive and have contact with only intact skin.

Needle Guide/Grid:

ABS

Stainless Steel

Polypropylene

Rolled Latex-Free Probe Cover K011265:

Polyisoprene

Rolled Latex Probe Cover K011265:

Latex Rubber

20ml. Sterile Ultrasound Gel 510(k)# K163050:

Device Characteristics of Advance Medical Designs, Inc. Needle Guides and Grids:

- Hypoallergenic, non-irritating
- No toxic effects
- Produced with completely harmless material.
- Does not damage the probe.
- Non-Invasive
- Contact with only intact skin, or with breached surfaces with duration <60 minutes
- Single Use
- Sterile (EtO Sterilization)

Device Identification-Model Numbers and Components:

The differences among the model numbers are limited to the convenience kits they are packaged into, the probe/ transducer/ ultrasound system they are compatible with, and the size specific dimensions indicating which Needle Gauge is appropriate for use. All Product Labels include the name of the specific compatible ultrasound systems that the Needle Guides and Grids are designed to fit. All Product Models are manufactured with the same ABS Polypropylene and injection molded. All accessories and kits are composed of the same accessories:

- 1) Advance Medical Designs (AMD) 20ml. Sterile Ultrasound Gel 510(k)# K163050:
- 2) Rolled Latex or Latex Free Probe/Transducer Cover 510(k)# K011265:

Indications for Use:

"AMD's (Advance Medical Designs) Disposable Needle Guides and Grids are used to assist and aid physicians in performing an endocavity diagnosis ultrasound needle guided procedure using guided intervention by providing fixed guiding for the precise insertion of linear instruments, such as needles. The Needle Guides and Grids are designed to aid adult patient population, in need of a biopsy of an internal organ, or internal delivery or removal of fluid within the body cavity, via the use of a needle, during an ultrasound procedure by retaining the needle tip and barrel within the ultrasound beam."

Intended Use:

The Subject Device provides fixed guiding of the precise insertion of linear instruments, such as needles through mechanical means with the use of diagnostic ultrasound equipment. The Needle Guide and Grid is attached over the transducer/ probe/ scanning instruments. The device provides a fixed path for the imaging guidelines for visualizing guided instrument placement procedures. Advance Medical

Designs (AMD) Needle Guides and Grids are furnished sterile; single use patient/ procedure, and disposable. The single use, disposable feature helps prevent transfer of microorganisms, body fluid, and particulate material to the patient and healthcare worker during reuse of the transducer. The Needle Guides and Grids are intended to clinicians and assistant clinicians, in clinical and hospital settings, to guide linear instruments during any ultrasound procedure that require the precision use of a needle.

<u>Comparison to Legally Marketed Device (Predicate and Reference Devices):</u>

Subject, Predicate Device and Reference Devices' indications for use place AMD's Needle Guides and Grids and CIVCO (510k# K970514) /KOELIS (510k# K180970) Ultrasound Needle Guides in device body contact categories as follows:

- a) Surface devices, intact skin / mucosal membranes/ breached surfaces, limited contact duration (<24 hours)
- b) External communicating devices, blood path indirect, tissue communicating limited contact duration (<24 hours)

Subject and Predicate Device have a similar Intended Use and Subject, Predicate device, and Reference devices (K180970, K222052) provide mechanical means for performing needle/ instrument guided procedures with the use of diagnostic ultrasound transducers. These devices provide the same fixed path for the needle or instrument that when coupled by the ultrasound system software corresponds to the on-screen imaging guidelines for visualizing guided instrument placement procedures.

Predicate device, Reference Devices, and Subject device are furnished sterile; and the entire guide is single use patient/ procedure, disposable. The single use, disposable feature helps prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during the reuse of the transducer. CIVCO Needle guide kits also provide ultrasound packets and covers which are similar to the subject device kits.

Reference Device (K141334) utilizes a reusable guide bracket (after cleaning/ sterilization by user) and a sterile, single use, disposable cannula.

Substantial Equivalence Table:

Company	Advance Medical Designs (SUBJECT DEVICE)	CIVCO (PREDICATE DEVICE)	KOELIS (REFERENCE DEVICE)	CIVCO (REFERENCE DEVICE)	KOELIS (REFERENCE DEVICE)
Device Name	Disposable Needle Guides and Grids	Disposable Endocavity Ultrasound Needle/ Biopsy Guide	Disposable Guides KDNG00	VitroPro/ Disposable Endocavity Needle Guide	Reusable Guide KRNG.EL1 KRNG.EL4
510(k) Number	K223689	K970514	K180970	K222052	K141334
Regulation Number	21 CFR 892.1570	21 CFR 892.1570	21 CFR 892.1570	21 CFR 892.1570/ 21 CFR 884.6100	21 CFR 892.1570
Device class	Class II	Class II	Class II	Class II	Class II

Device	Ultrasound	Ultrasound	Ultrasound	Ultrasound	Ultrasound
Common/	Transducer	Transducer	Transducer	Transducer	Transducer
-	Needle/Instrument	Needle/Instrument	Needle/Instrument	Needle/Instrument	Needle/Instrument
Usual Name	Guide	Guide	Guide	Guide	Guide
Device	Diagnostic	Diagnostic	Diagnostic	Diagnostic	Diagnostic
Classification	Ultrasonic	Ultrasonic	Ultrasonic	Ultrasonic	Ultrasonic
Name	Transducer	Transducer	Transducer	Transducer	Transducer
	Accessories	Accessories	Accessories	Accessories	Accessories
Classification	ITX	ITX	ITX	ITX, MQE	ITX
Product					
Code					
Indications	Disposable Need	le Guides and Grid	s are used to assis	t and aid physiciar	ns in performing
for Use	an endocavity dia	agnosis ultrasound	needle guided pro	ocedure using guid	led intervention
	by providing fixed	d guiding for the p	recise insertion of	linear instruments	s, such as
			rids are designed t		
			n, or internal deliv	•	
		_	, during an ultraso	•	
		arrel within the ult	_	and procedure by	. ctaning the
Intended Use	-		nd References Dev	icas provida fivad	guiding of the
intended Use				•	-
	•		nts, such as needle	•	
	_		uipment. The Nee		
	over the transducer/ probe/ scanning instruments. The device provides a fixed path for				
	the imaging guidelines for visualizing guided instrument placement procedures.				
	The Needle Guides and Grids are intended to aid clinicians and assistant clinicians, in				
	clinical and hospital settings, to guide linear instruments during any ultrasound				
	procedure that requires the precision use of a needle.				
	Subject, Predicate Device and Reference Devices' indications for use and intended use				
	place AMD's Disposable Needle Guides and Grids and CIVCO (510k# K970514) /KOELIS				
	(510k# K180970) Ultrasound Needle Guides in device body contact categories as follows:				
	a) Surface devices, intact skin / mucosal membranes/ breached surfaces, limited				
	contact duration (<24 hours)				
	b) External communicating devices, blood path indirect, tissue communicating				
	limited contact duration (<24 hours)				
	Subject and Predicate Device have a similar Intended Use and Subject, Predicate device,				
	and Reference devices (K180970, K222052) provide similar mechanical means for				
	performing needle/ instrument guided procedures with the use of diagnostic ultrasound				
Di	transducers.	Indiana, Co. 11	Diagramia and di	Indiana Control	In an C. 1d.
Design	Advance	Integrates the	Plastic guide	Integrates the	Inox Guide
	Medical	mounting	designed to be	mounting	designed to be
	Designs Needle	bracket and	clipped on an	bracket and	clipped on an
	Guide and Grid	cannula into a	endocavity	cannula into a	endocavity
	is a plastic	single	ultrasound	single	ultrasound
	guide designed	disposable	probe. An	disposable	probe. An
	to be clipped	component	entry cone to	component	entry cone to
	on to an	that attaches	easily	that attaches	easily
	ultrasound	externally, over	introduce the	externally, over	introduce the
	probe, with an	the transducer	needle into the	the transducer	needle into the
	entry cone to	with a clip-on	tube.	with a clip-on	tube.
		-	tube.	-	tube.
	easily	action.		action.	

	introduce the					
	needle into the					
	channel					
	Fixation mechanism of the Guide on the probe:					
	A clip to allow	A ring locks the	A clip to allow	A ring locks the	A clip to allow	
	the needle	needle guide	the needle	needle guide	the needle	
	guide stability	around the	guide stability	around the	guide stability	
	on the	probe thanks	on the	probe thanks	on the	
	transducer and	to a lateral	transducer and	to a lateral	transducer and	
	2 pins for	screw	2 pins for	screw	2 pins for	
	attachment in		attachment in		attachment in	
	the notches for		the notches for		the notches for	
	the probe		the probe		the probe	
Materials of	Medical Grade	Thermoplastic	Medical Grade	Thermoplastic	304 Stainless	
Construction	Thermoplastic	ABS	Polycarbonate	ABS	Steel	
	ABS	Polycarbonate	All materials	Polycarbonate	3161 Stainless	
	Polypropylene	304 Stainless	have met the	304 Stainless	Steel	
	304 Stainless	Steel	requirements	Steel	17/4 PH	
	Steel	All materials	of ISO 10993-1	All materials	All materials	
	All materials	have met the	for	have met the	have met the	
	have met the	requirements	biocompatibilit	requirements	requirements	
	requirements	of ISO 10993-1	у.	of ISO 10993-1	of ISO 10993-1	
	of ISO 10993-1	for	,	for	for	
	for	biocompatibilit		biocompatibilit	biocompatibilit	
	biocompatibilit	y.		y.	y.	
	y.	,		,	,	
Safety/	Meets ISO	Meets ISO	Meets ISO	Meets ISO	Meets ISO	
Biocompatibi	10993-1	10993-1	10993-1	10993-1	10993-1	
lity	biocompatibilit	biocompatibilit	biocompatibilit	biocompatibilit	biocompatibilit	
	у	у	у	у	у	
	requirements	requirements	requirements	requirements	requirements	
	for limited	for limited	for limited	for limited	for limited	
	contact	contact				
	00	contact	contact	contact	contact	
	duration:	duration:	duration:	contact duration:	contact duration:	
	duration:	duration:	duration:	duration:	duration:	
	duration: • surface	duration: • surface	duration: • surface	duration: • surface	duration: • surface	
	duration: • surface devices of	duration: • surface devices of	duration: • surface devices of	duration: • surface devices of	duration: • surface devices of	
	duration: • surface devices of breached or	duration: • surface devices of breached or	duration: • surface devices of breached or	duration: • surface devices of breached or	duration: • surface devices of breached or	
	duration: • surface devices of breached or compromised	duration: • surface devices of breached or compromised	duration: • surface devices of breached or compromised	duration: • surface devices of breached or compromised	duration: • surface devices of breached or compromised	
	duration: • surface devices of breached or compromised surface	duration: • surface devices of breached or compromised surface	duration: • surface devices of breached or compromised surface	duration: • surface devices of breached or compromised surface	duration: • surface devices of breached or compromised surface	
	duration: • surface devices of breached or compromised surface • External	duration: • surface devices of breached or compromised surface • External	duration: • surface devices of breached or compromised surface • Externally	duration: • surface devices of breached or compromised surface • Externally	duration: • surface devices of breached or compromised surface • Externally	
	duration: • surface devices of breached or compromised surface • External communicating	duration: • surface devices of breached or compromised surface • External communicating	duration: • surface devices of breached or compromised surface • Externally communicating	duration: • surface devices of breached or compromised surface • Externally communicating	duration: • surface devices of breached or compromised surface • Externally communicating	
	duration: • surface devices of breached or compromised surface • External communicating indirect blood	duration: • surface devices of breached or compromised surface • External communicating indirect blood	duration: • surface devices of breached or compromised surface • Externally communicating tissue	duration: • surface devices of breached or compromised surface • Externally communicating tissue	duration: • surface devices of breached or compromised surface • Externally communicating tissue	
	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin	
	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue contact	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue contact	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated	
	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue contact Demonstrated	duration: • surface devices of breached or compromised surface • External communicating indirect blood path/tissue contact Demonstrated	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated to be non-toxic,	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated to be non-toxic,	duration: • surface devices of breached or compromised surface • Externally communicating tissue bone/dentin Demonstrated to be non-toxic,	

	non-hemolytic,	non-hemolytic,	(DOES NOT	(DOES NOT	(DOES NOT
	and	and	INCLUDE	INCLUDE	INCLUDE
	nonpyrogenic	nonpyrogenic	PYROGEN/	PYROGEN/	PYROGEN/
			HEMO/ OR	HEMO/ OR	HEMO/ OR
			ACUTE	ACUTE	ACUTE
			SYSTEMIC)	SYSTEMIC)	SYSTEMIC)
Effectiveness	Subject, Predicate, and all Reference Devices are designed for secure and aligned fit to				
	the transducer or probe, while not altering the transducer or probe design integrity or				
	function. Positive Registration features of the design ensure accurate needle path and				
	placement in relation with the transducer. The exterior shapes of the guides are				
	contoured for the patient comfort with no sharp edges. Advance Medical Designs				
	needle guides and grids devices have the same intended use and their technological				
	characteristics do not raise any different questions of safety or effectiveness, as				
	compared to the legally marketed device. Therefore, the				
	Advance Medical Designs needle guides and grids are substantially equivalent to the				
	legally marketed disposable endocavity needle guide marketed by CIVCO (K972514).				
Sterilization	Ethylene Oxide	Ethylene Oxide	Ethylene Oxide	Ethylene Oxide	Ethylene Oxide
Shelf-Life	3 Years	3 Years	3 Years	3 Years	3 Years
Accessories	Ultrasound gel	Ultrasound gel	None provided.	None provided.	None provided.
	packet and	packet and		Intended	
	covers.	covers.		user to provide	
				cover and	
				guides which	
				are IVF use	
				cleared.	
Energy Type	None	None	None	None	None
Software	None	None	None	None	None

Design principles of operation are similar between subject device and predicate device(s):

- Placing the Probe Attachment Bracket on the transducer using locating features
- Securing the needle guide onto the transducer by fastening the Attachment Bracket Lock
- Identifying the angle / position to encourage optimal needle trajectory.
- Releasing the Attachment Bracket Lock to remove the needle guide from the transducer.

Subject Device and Predicate Device(s) integrate the mounting bracket and cannula into a single, disposable component that attaches externally, over the transducer; however predicate device uses a clip-on action and subject device uses an Attachment Bracket Lock. The method of maintaining a secure attachment does not change the intended purpose and the devices maintain fixed positions as intended.

Materials of Construction and Manufacturing:

The Advance Medical Designs' Needle Guides, Grids, and Accessories (Kit) have a non-pharmacological, immunological, or metabolic mode of action and both subject and predicate devices have same Shelf-Life's of 3 years.

Predicate Device is fabricated from:

- 1) Injection molded thermoplastic components, bonded to stainless steel cannula with medical grade adhesive.
- 2) Injection insert-molded thermoplastic with integral stainless-steel cannula
- 3) Thermoplastics (ABS and Polycarbonate) and 304 Stainless Steel and packaged in Polyethylene and Tyvek

Subject Device is fabricated from:

- 1) Injection molded thermoplastic components, but with a stainless-steel attachment.
- 2) Injection insert- molded thermoplastic with integral stainless-steel cannula
- 3) Thermoplastics (ABS and Polypropylene) and 304 Stainless Steel and packaged in Polyethylene and Tyvek

Predicate Device and Reference Device (K180970) uses ABS Polycarbonate, and the Subject Device uses ABS Polypropylene.

Chemical Characterization of Polypropylene have been conducted and those studies and results are provided with this application and prove that the materials used to manufacture AMD's Needle Guides and Grids are safe, fully biocompatible, and as effective as Predicate Device and Reference Devices. Materials and manufacturing processing for Predicate Device and Subject Device (including 100 % EtO sterilization w/ SAL 10⁻⁶) affects to the healthcare worker and patient via intended use/ indications for use contact of this device have been biologically evaluated using biocompatibility tests for cytotoxicity, irritation, sensitization, pyrogenicity, acute system toxicity, and hemocompatibility.

Sufficient evidence is provided to validate AMD's Claim that subject materials/ device to be non-toxic, non-sensitizing, non-irritating, non-hemolytic, and non-pyrogenic.

Modified Technological Characteristics Deviating from Predicate Device(s) Design, Effectiveness, and Safety:

- 1. Cannula Guide Channels
 - a. Formulation: The guide channel design is based upon the ultrasound system and encourages the needle to align with variable guidelines generated on the system monitor.
 - b. Composition: A rigid frame containing evenly spaced through-holes and a connection to the Probe Attachment Bracket.
 - c. Functionality: These channels guide a cannula to a set position that aligns with gridlines displayed by the corresponding ultrasound monitor.
- 2. Probe Attachment Bracket
 - a. Formulation: Every bracket is uniquely designed to meet each transducer's specific geometry to maximize the contact surface area and limit the ability for the cannula to shift or rotate about the transducer.
 - b. Composition: A mounting structure attached to the Cannula Guide Channels which contains: a surface that fits snugly onto the transducer, a hinge to anchor the Attachment Bracket Lock, a hook to catch the mobile end of the Attachment Bracket Lock.

c. Functionality: This bracket connects the guide channels to the upper side of the transducer in a reliable, inflexible manner such that the channels maintain their position and orientation relative to the transducer. This bracket mates with the bracket lock to physically secure the connection to mitigate the effects of impulse motion or jerk movements.

3. Attachment Bracket Lock

- a. Formulation: This lock is designed to be pulled into tension around the transducer and meets each transducer's specific geometry to achieve adequate holding pressure.
- b. Composition: A semi-rigid arm containing two hooks, one hook for anchoring to the hinge of the Attachment Bracket Lock and another hook on the mobile end to latch onto the hook protruding from the Probe Attachment Bracket.
- c. Functionality: This lock secures the connection between the Probe Attachment Bracket and the transducer. The user wraps this lock around the underside of the transducer and hooks it to one side of the Probe Attachment Bracket.

Summary of Non-Clinical Tests Performed on Subject:

Biocompatibility: The disposable Needle Guides and Grids meet ISO 10993-1 biocompatibility and ASCA- Pilot Biocompatibility Guidance requirements for limited contact duration for surface devices of breached or compromised surface external communicating tissue and indirect blood path.

- a. Cytotoxicity ISO 10993-5
- b. Sensitization ISO 10993-10
- c. Irritation ISO 10993-10
- d. Acute Systemic Toxicity—ISO 10993-11
- e. Material-Mediated Pyrogenicity- 10993-11 and USP 151
- f. Direct and Indirect Hemolysis/Hemocompatibility ISO 10993-4 and ASTM F756

Summary of Bench Testing Performed on Subject Device:

- 1) Cover breach and probe damage testing Water leak testing was performed during Design Input, Design Output, and Design Validation to demonstrate material attachment of needle guide over a cover did not cause damage to cover or probe.
- 2) Retention and movement testing Force testing was performed on needle guide attachment to ensure a minimum force of 8N would not cause the guide to dislodge.
- 3) Needle drag testing Force testing was performed by passing a cannula through the needle guide to ensure binding would not occur and force was less than a 1.5N threshold.
- 4) Needle path verification testing Needle guides were tested on test fixtures to ensure needle path falls within the design tolerances specified for the design.
- 5) Simulated Usability Testing Simulated use evaluations were performed by customers to ensure the design of the needle guide conforms to the user needs and intended use as well as imaging testing conducted through laboratory evaluations.

Clinical Test Performed:

Clinical tests were not required to demonstrate substantial equivalence.

Conclusions:

The comparisons of the technological and non-clinical performance characteristics indicate Advance Medical Designs, Inc. Disposable Needle Guides & Grids have the same intended use and its technological characteristics do not raise any different questions of safety or effectiveness, as compared to the legally marketed device(s). The improvements made in designs regarding secure attachment prove to be superior and establish the same level, if not greater, of safety and usage. Therefore, AMD's Needle Guides and Grids are substantially equivalent to the legally marketed predicate device disposable endocavity guide marketed by CIVCO. Since the comparison of bench testing to clinical outcomes, documented through images provided, the subject device and predicate device are substantially equivalent. Thus, AMD's (Advance Medical Designs) Needle Guides and Grid's claim demonstrates the device performs comparably to the predicate device and reference devices that are currently marketed for the same intended use.