

June 15, 2023

STERIS Corporation Gregory Land Lead Regulatory Affairs Specialist 5960 Heisley Rd Mentor, Ohio 44060

Re: K223715

Trade/Device Name: Celerity Incubator Regulation Number: 21 CFR 880.2800 Regulation Name: Sterilization Process Indicator Regulatory Class: Class II Product Code: FRC Dated: May 19, 2023 Received: May 19, 2023

Dear Gregory Land:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

## Christopher K. Dugard -S

for Clarence W. Murray, III, PhD Deputy Director DHT4B: Division of Infection Control and Plastic Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

#### Indications for Use

Submission Number (if known)

K223715

Device Name

Celerity Incubator

Indications for Use (Describe)

Use the Celerity Incubator to incubate and automatically read STERIS Celerity Biological Indicators for Steam and Vaporized Hydrogen Peroxide sterilization at 55°C - 60°C for a fluorescent result.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# STERIS<sup>®</sup>

### 510(k) Summary For Celerity Incubator

STERIS Corporation 5960 Heisley Road Mentor, OH 44060 Phone: (440) 354-2600 Fax No: (440) 357-9198

Contact:	Mr. Gregory Land Lead Regulatory Affairs Specialist Tel: 440-392-7424 greg_land@steris.com
Submission Date:	June 13, 2023

K Number: K223715

STERIS Corporation = 5960 Heisley Road = Mentor, OH 44060-1834 USA = 440-354-2600

#### 1. <u>Device Name</u>

Trade Name:	Celerity Incubator
Device Class:	Class II
Common/usual Name:	Incubator
Classification Name:	Indicator, Biological Sterilization Process
Classification Number:	21 CFR 880.2800
Product Code:	FRC

#### 2. <u>Predicate Device</u>

K213881 - Celerity Incubator

#### 3. <u>Device Description</u>

The Celerity Incubator incubates and reads fluorescent Vaporized Hydrogen Peroxide (VHP) and Steam Biological Indicators (BIs). The BI is incubated at 59°C with an acceptable tolerance of -4°C/+1°C. During incubation, the BI is monitored for a potential fluorescence signal generated as a result of the production of  $\alpha$ -glucosidase. When a growth response has been detected or when the required incubation time has elapsed, the incubator indicates the results to the user.

#### 4. <u>Indications for Use</u>

Use the Celerity Incubator to incubate and automatically read STERIS Celerity Biological Indicators for Steam and Vaporized Hydrogen Peroxide sterilization at 55°C - 60°C for a fluorescent result.

#### 5. <u>Technological Characteristics Comparison Table</u>

A comparison of technical characteristics between the proposed and predicate devices is summarized in **Table 1**.

Feature	Proposed Device Universal Biological Indicator Incubator	Predicate Device Celerity Indicator (K213881)	Comparison
Indications for Use	Use the Celerity Incubator to incubate and automatically read STERIS Celerity Biological Indicators for	Use the Celerity Incubator to incubate and automatically read STERIS Celerity Biological Indicators for	Same

#### Table 1. Technological Characteristics Comparison Table to predicate device

Feature	Proposed Device Universal Biological Indicator Incubator	Predicate Device Celerity Indicator (K213881)	Comparison
	Steam and Vaporized Hydrogen Peroxide sterilization at 55°C - 60°C for a fluorescent result.	Steam and Vaporized Hydrogen Peroxide sterilization at 55°C - 60°C for a fluorescent result.	
Basis of Readout	Photodiode detects fluorescence produced by enzymatic activity that results from viable biological indicator organisms.	Photodiode detects fluorescence produced by enzymatic activity that results from viable biological indicator organisms.	Same
Incubation Temperature Range	55 - 60 °C	55 - 60 °C	Same
Voltage Range	100-240 VAC with 12 VDC conversion	100-240 VAC with 12 VDC conversion	Same
Test Capacity	8 wells	8 wells	Same
Calibration	Factory calibration – no calibration by customer	Factory calibration – no calibration by customer	Same
Incubation Time	5 minutes or 20 minutes depending on Biological Indicator	20 minutes	Similar
Fluorescence Detection	UV LEDs are used to excite the fluorescent molecule produced by enzyme cleavage of the fluorogenic substrate contained in the BI's media. The emitted light is detected by a photodiode.	UV LEDs are used to excite the fluorescent molecule produced by enzyme cleavage of the fluorogenic substrate contained in the SCBI's media. The emitted light is detected by a photodiode.	Same
Indication of Results	Positive – audible alarm, visual LED lights and screen Negative – optional alarm, visual indication with LED lights and LCD screen User must acknowledge results	Positive – audible alarm, visual LED lights and screen Negative – optional alarm, visual indication with LED lights and LCD screen User must acknowledge results	Same
System Operation	The reader/incubator wells are arranged in 2 banks of 4 wells and preset to 59°C. The user scans the barcode on the label of an activated BI using the system's barcode scanner	The reader/incubator wells are arranged in 2 banks of 4 wells and preset to 59°C. The user scans the barcode on the label of an activated BI using the system's barcode scanner	Same

	<b>Proposed Device</b>	Predicate Device	
Feature	Universal Biological Indicator Incubator	Celerity Indicator (K213881)	Comparison
Feature	Proposed Device Universal Biological Indicator Incubator and places it into an open well. The system detects the well the BI was placed into and begins measurement of fluorescence; a blinking yellow light indicates the incubation is in process and the read initiated. The System uses information from the barcode to apply the appropriate algorithm to the well. A "positive" reading is interpreted as an indication of a potential sterilization cycle failure. A "positive" finding is indicated to the user by red light on the front panel adjacent to the well, by an audible alarm, and by text displayed on the LCD screen. The alarm must be muted by the operator when a positive result is obtained. The LCD screen provides instructions for the user to turn off the alarm for that specific BI. Should another BI become "positive", the alarm will sound again and the above actions are repeated. If the result is not positive at the end of the incubation time, the result is negative.	Predicate Device Celerity Indicator (K213881) and places it into an open well. The system detects the well the BI was placed into and begins measurement of fluorescence; a blinking yellow light indicates the incubation is in process and the read initiated. The System uses information from the barcode to apply the appropriate algorithm to the well. A "positive" reading is interpreted as an indication of a potential sterilization cycle failure. A "positive" finding is indicated to the user by red light on the front panel adjacent to the well, by an audible alarm, and by text displayed on the LCD screen. The alarm must be muted by the operator when a positive result is obtained. The LCD screen provides instructions for the user to turn off the alarm for that specific BI. Should another BI become "positive", the alarm will sound again and the above actions are repeated. If the result is not positive at the end of the incubation time, the result is negative.	Comparison
	Negative result is negative. Negative results are identified by a green light on the front panel adjacent to the well	Negative result is negative. Negative results are identified by a green light on the front nanel adjacent to the well	
	with the "negative" BI and by text on the LCD. In addition, an optional alarm is available for negative results.	with the "negative" BI and by text on the LCD. In addition, an optional alarm is available for negative results.	

#### 6.

<u>Summary of Non-Clinical Performance Testing</u> Testing was performed to evaluate performance and demonstrate substantial equivalence to the predicate as summarized in Table 2.

Test	Acceptance Criteria	Result
Qualification Testing with intended Biological Indicators	Fluorescent Read meets > 97% alignment with 7-day grow out per FDA guidance on reduced incubation time for intended Biological Indicators for Vaporized Hydrogen Peroxide Celerity 5 Biological Indicators	PASS
Human Factors	Typical users are capable of following the written instructions for use to correctly use the Celerity Incubators.	PASS
Electromagnetic Compatibility	IEC 60601-1-2:2014 General Requirements for Basic Safety and Essential Performance – Collateral Standard: Electromagnetic Compatibility – Requirements and Tests	PASS
Electrical Safety Conformance	IEC 60101-1 Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use Part 1: General Requirements	PASS
Software Validation per IEC 62304	The software that controls the system was validated and determined to operate effectively and as designed.	PASS

#### **Table 2. Performance Testing**

#### 7. <u>Conclusion</u>

The Celerity Incubator has met the established performance criteria. Based on the intended use, technological characteristics and non-clinical performance data, the subject device is as safe, as effective, and performs as well as the legally marketed predicate device (K213881), Class II (21 CFR 880.2800), product code FRC.