

June 28, 2023

Mio Labs Inc.
Mark Qian
Senior RA Engineer
#1023, ZGC Innovation Center, 4500 Great America Pkwy
Santa Clara, CA 95054

Re: K223722

Trade/Device Name: MIO Blood Glucose Monitoring System

Regulation Number: 21 CFR 862.1345 Regulation Name: Glucose Test System

Regulatory Class: Class II Product Code: NBW Dated: May 30, 2023 Received: June 2, 2023

#### Dear Mark Oian:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

K223722 - Mark Qian Page 2

801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Paula V. Caposino -S

Paula Caposino, Ph.D.
Acting Deputy Director
Division of Chemistry and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

# **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

k223722	
Device Name	
MIO Blood Glucose Monitoring System	
Indications for Use (Describe)	
MIO Blood Glucose Monitoring System is comprised of the MIC	Blood Glucose Meter and the MIO Blood Glucose Test
Strips. MIO Blood Glucose Monitoring System is intended to que capillary whole blood samples drawn from the fingertips. It is into monitor the effectiveness of diabetes control. It is not intended to the control of t	ended for use by persons with diabetes at home as an aid for neonatal use or for the diagnosis of or screening for
diabetes. This system is intended for self-testing outside the body single person and should not be shared.	(in vitro diagnostic use), and should only be used by a
engle process and one are the state of	
Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

CONTINUE ON A SEPARATE PAGE IF NEEDED.

# \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.qian@transtekcorp.com Tel: 301-910-0529

# Section E: 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

The Assigned 510(k) number is k223722.

#### Submitter's Identification:

Mio Labs Inc.

#1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054

Tel: 301-910-0529

Date Updated: May 25, 2023

#### **Contact Person:**

Mark Qian Senior RA Engineer Mio Labs Inc. #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054

## **Proprietary Name of the Device:**

MIO Blood Glucose Monitoring System

Common Name: Glucose Test System

#### **Classification Name:**

Class II §862.1345 Glucose Test System

Product Code: NBW

#### **Predicate Device:**

VivaChek Ino Smart Blood Glucose Monitoring System VivaChek Laboratories, Inc. 510(k) Number: k173140

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.qian@transtekcorp.com Tel: 301-910-0529

#### **Device Name:**

Proprietary Name	Meter Model	Classification	Product Code	Description	Common Name
MIO Blood Glucose	TeleBGM	862.1345	NBW	System, Test, Blood	Glucose
Monitoring System	2282-G	Class II		Glucose, Over The Counter	Test System

#### **Description:**

MIO Blood Glucose Monitoring System is designed to quantitatively measure the glucose concentration in fresh capillary whole blood. The glucose measurement is achieved by using the amperometric detection method. The test is based on measurement of electrical current caused by the reaction of the glucose with the reagents on the electrode of the test strip. The blood sample is pulled into the tip of the test strip through capillary action. Glucose in the sample reacts with glucose oxidase and the mediator. Electrons are generated, producing a current that is positive correlation to the glucose concentration in the sample. After the reaction time, the glucose concentration in the sample is displayed.

MIO Blood Glucose Monitoring System contains 4G module, the device complies with US federal guidelines, FCC Part 15 Subpart B, FCC Part 2, FCC Part 24 Subpart E, FCC Part 27 Subpart C, and FCC 47 CFR§ 2.1093 based on the test reports.

#### Intended Use:

MIO Blood Glucose Monitoring System is comprised of the MIO Blood Glucose Meter and the MIO Blood Glucose Test Strips. MIO Blood Glucose Monitoring System is intended to quantitatively measure the glucose concentration in fresh capillary whole blood samples drawn from the fingertips. It is intended for use by persons with diabetes at home as an aid to monitor the effectiveness of diabetes control. It is not intended for neonatal use or for the diagnosis of or screening for diabetes. This system is intended for self-testing outside the body (in vitro diagnostic use), and should only be used by a single person and should not be shared.

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.qian@transtekcorp.com Tel: 301-910-0529

# **Comparison to Predicate Device:**

The MIO Blood Glucose Monitoring System is substantially equivalent to VivaChek Ino Smart Blood Glucose Monitoring System (k173140)

Features	Predicate: VivaChek Ino Smart Blood Glucose Monitoring System (k173140)	Candidate: Blood Glucose Monitoring System
	Similarities	
Intended Use	It is intended to quantitatively measure the glucose concentration in fresh capillary whole blood samples drawn from the fingertips. It is intended for use by persons with diabetes at home as an aid to monitor the effectiveness of diabetes control. It is not intended for neonatal use or for the diagnosis of or screening for diabetes. It is intended for self-testing outside the body (in vitro diagnostic use), and should only be used by a single person and should not be shared.	Same
Operation Principle	Electrochemical biosensor	Same
Detection Method	Amperometric	Same
Strip Chemical Composition	Glucose oxidase	Same
Measurement Result	Plasma Glucose	Same
Sample	Fresh capillary whole blood	Same
Memory	500 records	Same
Unit of Measure	mg/dL	Same
Measurement Range	20-600 mg/dL	Same
Sample Volume	0.8µL	Same
Test Time	5 seconds	Same
Operating Relative Humidity	10-90% (non-condensing)	Same
Operating Temperature	41–113°F	Same
Hematocrit Range	20–70%	Same
Automatic Shutoff	2 minutes after last action	Same
Power Source	Rechargeable 3.7 Volt Lithium Ion battery	Same
	Differences	
Battery Type	Rechargeable, non-serviceable, 250mAh, 3.7 Volt DC nominal, lithium polymer battery (5V input charge voltage)	Rechargeable, 800 mAh, 3.7 Volt DC nominal, lithium polymer battery (5V input charge voltage)

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.qian@transtekcorp.com Tel: 301-910-0529

Dimensions	83 mm x 52 mm x 18.7 mm	95.5 mm x 59.1 mm x 20.5 mm
Display Size	32mm x 32 mm	47 mm x 37.5 mm
Weight	Approximately 53g	Approximately 70g
Data Transmission	Bluetooth	4G

# **Laboratory Testing:**

The performance characteristics of the MIO Blood Glucose Monitoring System were evaluated by performing the following studies:

No.	Test/Validation Item
1	Within-Run Precision Evaluation
2	Intermediate Precision Evaluation
3	Linearity Evaluation Study
4	User Evaluation
5	User Evaluation - Accuracy at Extreme Glucose Values
6	Interference Agents Study
7	Hematocrit Effect Study
8	Oxygen Interference Study
9	Flex Studies- Mechanical Vibration Testing
10	Flex Studies- Shock Testing
11	Flex Studies- Electromagnetic compatibility (EMC) Testing
12	Flex Studies- Electrostatic Discharge/Electromagnetic Interference Testing
13	Accelerated Closed Vial Stability Study _Strip
14	Accelerated Open Vial Stability Study _Strip
15	Accelerated Closed Vial Stability Study _Control
16	Accelerated Open Vial Stability Study _Control
17	Sample Perturbation Study
18	Real Time Closed Vial Stability Study _Strip

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.qian@transtekcorp.com Tel: 301-910-0529

19	Real Time Open Vial Stability Study _Strip
20	Real Time Closed Vial Stability Study _Control
21	Real Time Open Vial Stability Study _Control
22	Operating Conditions Evaluation
23	Altitude Effect Evaluation
24	Error Codes Evaluation
25	Short Sample Detection Study
26	Intermittent Sampling Study
27	Control Range Assignment
28	Test Strip Lot Release Criteria
29	Low Battery Study
30	Shipping Study (Strip)
31	Shipping Study (Control)
32	Firmware (Software) Validations
33	Electrical Safety Test
34	Environmental Temperature Test
35	Testing with Used Test Strips
36	FCC Tests
37	Cybersecurity Control DFMEA
38	Meter Cybersecurity Management Plan
39	Web App Software Validation

To confirm the 4G module has not brought any unexpected functional failure or adverse effect, DFMEA, cybersecurity control and software validation were conducted and summarized in above table.

Address: #1023, ZGC Innovation Center, 4500 Great America Pkwy, Santa Clara, CA 95054 Contact: Mark Qian Email: Mark.gian@transtekcorp.com Tel: 301-910-0529

### **Discussion of Laboratory Studies:**

Above laboratory studies were performed on MIO Blood Glucose Monitoring System in accordance with the applicable guidance or standards, and the test results indicated that the acceptance criteria were met. Therefore, the performances from these laboratory studies were acceptable.

# **Discussion of Clinical Study:**

Clinical study (user evaluation) was conducted with intended users using the MIO Blood Glucose Monitoring System. Study results indicated that non-professional, inexperienced lay persons were able to obtain blood glucose readings when using the MIO Blood Glucose Monitoring System. In addition, the participated lay persons were questioned and responded as satisfied with the ease of operation by following the User Manual and the overall performance of the MIO Blood Glucose Monitoring System.

#### **Conclusion:**

The laboratory studies and user evaluation study results demonstrate that the MIO Blood Glucose Monitoring System is safe, effective and easy-to-use. It also demonstrates that the MIO Blood Glucose Monitoring System meets FDA Guidance SMBG for OTC Use.

Based on the same intended use and work principle, and the same technological characteristics listed in the section Comparison to Predicate Device - Similarities, meanwhile the different technological characteristics listed in the section Comparison to Predicate Device do not raise safety and effectiveness questions according to the completed performance validation reports, therefore the candidate device MIO Blood Glucose Monitoring System is substantially equivalent to the predicate device VivaChek Ino Smart Blood Glucose Monitoring System (k173140).