

August 15, 2023

Paragonix Technologies Nathan Yetton Senior Director of Quality 639 Granite St., Suite 408 Braintree, Massachusetts 02184

Re: K223874

Trade/Device Name: BAROguard Regulation Number: 21 CFR 876.5880

Regulation Name: Isolated Kidney Perfusion and Transport System

and Accessories

Regulatory Class: II Product Code: KDN Dated: July 17, 2023 Received: July 18, 2023

Dear Nathan Yetton:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

K223874 - Nathan Yetton Page 2

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Glenn B. Bell -S

for Gema Gonzalez, MS
Acting Assistant Director
DHT3A: Division of Renal, Gastrointestinal,
Obesity and Transplant Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Indications for Use		See PRA Statement below.
510(k) Number <i>(if known)</i> K223874		
Device Name BAROguard		
ndications for Use (Describe) BAROguard™ is intended to be used for the static hypothermic presentation into a recipient using cold storage solutions indicators intended organ storage time for BAROguard™ is up to 8 hour Donor lungs exceeding clinically accepted static hypothermic preseurgeon to determine transplantability in accordance with accepted the intended recipient.	ted for use with the rs.	lungs.
Note: Partial lungs can be transported via BAROguard™ by packaguidelines.	aging lungs per inst	itutional protocol and UNOS
Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter	er Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATI	E PAGE IF NEEDE	D.
This section applies only to requirements of th	e Paperwork Reduc	tion Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

> Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

005_Updated 510k Summary

Paragonix Technologies' BAROguard Device (K223874)

Submitter: Paragonix Technologies Inc.

c/o Vaughn & Associates

639 Granite Street Braintree, MA02184

Contact Person: Nathan Yetton

Paragonix Technologies, Inc.

222 Third St

Cambridge, MA 02142 413.345.1814 (phone)

nathan@paragonixtechnologies.com

Date Prepared: July 14, 2023

Trade Name: BAROguard

Classification Name: Isolated kidney perfusion and transport system and

accessories

Classification: Class II

Regulation Number: 21 CFR § 876.5880

Product Code: KDN

Predicate Device: Paragonix LUNGguard (previously named, SherpaPak

Lung Preservation System) - K192869

Device Description:

The subject BAROguard device results from modifications made to the cleared LUNGguard (previously named, SherpaPak Lung Preservation System) cleared under K192869. The subject BAROguard device consists of the following components:

1) BAROguard SherpaCool Pouches – Phase Change Material (PCM) pouches (identical to the predicate) to maintain temperature of the cold preservation solution and lung throughout transportation. The BAROguard device maintains the temperature between 4°C to 8°C identical to the predicate with the use SherpaCool pouches throughout preservation and transportation.

005 Updated 510k Summary

Paragonix Technologies' BAROguard Device (K223874)

- 2) BAROguard Lung Containment Assembly Nested lung containment bags for the packaging of donor lungs and preservation solution. BAROguard Lung Containment Assembly includes pneumatic connections to the donor lung in the inner-most bag, a pneumatic connection to the BAROguard Shipper Airway Pressure Management System, endotracheal connectors to connect the trachea of the donor lungs to the BAROguard Lung Containment Assembly, and tools for the secure attachment of the endotracheal connectors and closure of the BAROguard nested lung containment bags.
- 3) BAROguard Shipper– Outer transport shipper which comprises a protective and insulative package. The BAROguard Shipper is a rigid, molded expanded polystyrene (EPS) insulative container and into which the SherpaCool pouches and donor lung within the BAROguard Lung Containment Assembly are placed. The maintenance of temperature of the donor lung between 4°C to 8°C is assisted by the EPS insulation of the BAROguard Shipper, providing insulation from the exterior environment to the interior components and BAROguard SherpaCool.

The BAROguard Shipper incorporates an Airway Pressure Management System. The Airway Pressure Management System maintains the donor lung airway pressure when connected to the BAROguard Lung Containment Assembly.

The BAROguard Shipper includes an off-the-shelf datalogger (connected to a temperature probe and airway pressure sensor) which monitors and displays the temperature of the solution surrounding the donor lungs and the airway pressure of the donor lungs.

Although the BAROguard is based on the predicate LUNGguard design, it also includes two new design elements:

- Incorporation of a sterile Lung Containment Assembly instead of use of offthe-shelf bags in K192869.
- Incorporation of an Airway Pressure Management System to maintain donor lung airway pressure during preservation and transportation.



005 Updated 510k Summary

Paragonix Technologies' BAROguard Device (K223874)

Intended Use:

Donor lung preservation and transportation

Indications for Use:

The Paragonix BAROguard is intended to be used for the static hypothermic preservation of lungs during transportation and eventual transplantation into a recipient using cold storage solutions indicated for use with the lungs.

The intended organ storage time for BAROguard is up to 8 hours.

Donor lungs exceeding clinically accepted static hypothermic preservation times should be evaluated by the transplant surgeon to determine transplantability in accordance with accepted clinical guidelines and in the best medical interest of the intended recipient.

Note: Partial lungs can be transported via BAROguard by packaging lungs per institutional protocol and UNOS guidelines.

Functional Testing:

- Biocompatibility testing of any new materials that contact the body
 - The biocompatibility evaluation for the modified devices was conducted in accordance with International Standard ISO 10993-1: 2018 "Biological Evaluation of Medical Devices Part 1: Evaluation and Testing Within a Risk Management Process," as recognized by FDA. The battery of testing included the following tests: Cytotoxicity, Sensitization, Irritation, Acute Systemic toxicity, Material Mediated Pyrogenicity, and Hemocompatibility
- Electrical Safety and EMC testing in accordance with the following standards:
 - IEC 60601-1:2005, AMD1:2012, AMD2:2020, IEC 60601-1-2:2014
 (ed. 4.1), IEC 60601-1-6:2010, AMD1:2013, AMD2:2020, IEC 62366-1:2015, AMD1:2020, FCC 47CFR Part 15.247:09
- Thermal and Airway Pressure Validation
 - Validation demonstrates the ability of the BAROguard device to maintain hypothermic preservation and airway pressure of the donor lung beyond the intended organ storage time of 8 hours.

005 Updated 510k Summary

Paragonix Technologies' BAROguard Device (K223874)

- BAROguard Shipper Verification
 - Verification demonstrates BAROguard Shipper meets specifications beyond the intended organ storage time of 8 hours following exposure to worst-case shipping and handling.
- Lung Containment Assembly Verification
 - Verification demonstrates BAROguard sterile Lung Containment Assembly meets specifications following exposure to sterilization and accelerated aging simulating real-time aging.
- BAROguard Sterile Packaging Validation
 - Verification demonstrates the ability of the BAROguard sterile barrier system to maintain sterility following exposure to sterilization and accelerated aging simulating real-time aging.
- 0.2 Micron Filter Validation
 - Validation demonstrates the bacterial retention of the filter used within the BAROguard device.
- ISO 18562-2
 - \circ The average total particulate matter for the BAROgaurd device was found to be 10 μ g/m3 for the entire 24 hours of continuous airflow through the system.

Technological Comparison with Predicate and Reference Device

The following table compares the Paragonix BAROguard device with the predicate and reference devices.

Table 1. Sub Characteristic	Table 1. Substantial Equivalence Comparison Table Characteristic Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard (K223874)	SherpaPak Lung Preservation System (K192869)	Auto CPAP System (K211155)	Comparison
Intended Use	Donor lung preservation and transportation.	Organ preservation and transportation.	To support treatment of adult Obstructive Sleep Apnea (OSA)	Identical to predicate (substantially equivalent)
for Use	used for the static hypothermic preservation of lungs during transportation and eventual transplantation into a recipient using cold storage solutions indicated for use with the lungs. The intended organ storage time for BAROguard is up to 8 hours. Donor lungs exceeding clinically accepted static hypothermic preservation times should be evaluated by the transplant surgeon to determine transplantability in accordance with accepted clinical guidelines and in the best medical interest of the intended recipient. Note: Partial lungs can be transported via BAROguard by packaging lungs per institutional protocol and LINOS guidelines.	Preservation System (renamed as LUNGguard) is intended to be used for the static hypothermic preservation of lungs during transportation and eventual transplantation into a recipient using cold storage solutions indicated for use with the lungs. The intended organ storage time for the Paragonix SherpaPak Lung Preservation system is up to 8 hours. Donor lungs exceeding clinically accepted static hypothermic preservation times should be evaluated by the transplant surgeon to determine transplantability in accordance with accepted clinical guidelines and in the best medical interest of the intended	is a CPAP (Continuous Positive Airway Pressure) device designed for the treatment of adult Obstructive Sleep Apnea (OSA) only, either in the hospital or at home. Auto CPAP System is for prescription use only. It is a travel CPAP device intended for single-patient use.	predicate, with the only difference being the respective product names (substantially equivalent)
Regulation Number	876.5880	SherpaPak LPS by packaging lungs per institutional protocol and UNOS guidelines. 876.5880	868.5905	Identical to predicate (substantially equivalent)
Product Code	KDN	KDN	BZD	Identical to predicate (substantially equivalent)
Device Classificatio n Name	Isolated kidney perfusion and transport system and accessories.	Isolated kidney perfusion and transport system and accessories.	Noncontinuous ventilator	Identical to predicate (substantially equivalent)

Table 1. Sub	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard (K223874)	SherpaPak Lung Preservation System (K192869)	Auto CPAP System (K211155)	Comparison
Operating Principle	Static hypothermic storage (i.e., Static Cold Storage) between 4° C to 8° C using FDA-cleared preservation solutions for lung organs, with management of donor lung airway	Static hypothermic storage (i.e., Static Cold Storage) between 4° C to 8° C using FDA cleared preservation solutions for lung organs. Lungs clinically inflated to static	The system provides continuous positive airway pressure from 4 to 20 cmH ₂ O above	Substantially equivalent to predicate The subject device
	pressure throughout preservation including during air transportation with continuous positive airway pressure of 12-15 cmH ₂ O (±1 cmH ₂ O) above ambient atmospheric pressure.	pressure of 12-15 cmH ₂ O prior to storage ¹	the ambient atmospheric pressure.	provides the identical static hypothermic preservation as the predicate and also incorporates an Airway Pressure Management System with positive airway pressures similar to the reference device that is used in the same anatomy for the same physiological purpose. (See discussion following table on Substantial Equivalence in Operating
				Principle)
Intended Storage Time	Thermal qualification demonstrates the device can maintain 4° C to 8° C through the intended organ maximum cold ischemic time (CIT) with high and low temporature excursions (i.e., in to high and low temporature excursions (i.e., in to	Thermal qualification demonstrates the device can maintain 4° C to 8° C through the intended organ maximum cold ischemic time (CIT) with high and low temporature	N/A	Identical to predicate (substantially equivalent)
Single Use	Entire system is single-use/single-patient only.	Entire system is single-use/single-patient	Reusable system for	Identical to predicate
Meets UNOS	Yes	Yes	Single-patient use	(substantially equivalent)

¹ Copeland et al., Donor heart and lung procurement: A consensus statement, J Heart Lung Transplant 2020;39:501–517

² http://www.optn.transplant.hrsa.gov

I dold Ti Can	Table Ti Sassialitial Edutationes combanison rasic			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard	SherpaPak Lung Preservation System (K192869)	Auto CPAP System	Comparison
	(K223874)		(K211155)	
Shipper	Outer molded EPS Shipping container with	Outer molded EPS Shipping container with	N/A	Substantially equivalent to
	wheels, extending handle, integrated Airway Pressure Management System, and off-the-	wheels, extending handle, and off-the-shelf datalogger		predicate
	shelf datalogger	uataivyyei		Identical EPS insulative
				shipper material and
				identical phase change
				(cooling) material are
				used for the subject and
				predicate devices.
				The subject device differs
				from predicate in that it
				incorporates an Airway
				Pressure Management
				System. (See discussion
				following table on
				Substantial Equivalence in
				Operating Principle)
				See "Monitor" row for
				details on the datalogger.

Table 1. Sub	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard (K223874)	SherpaPak Lung Preservation System (K192869)	Auto CPAP System (K211155)	Comparison
Organ Container	Use of proprietary Lung Containment Assembly (including triple nested bags and endotracheal	Use of 3M Steri-Drape isolation bags cleared under K832318 (Product Code KGY,	A/N	Substantially equivalent
	from the lungs to the airway pressure	The instructions for use for the subject		predicate devices each
	management system. The packaging is consistent with standard practice and OPTN	devices recommend triple bagging the lungs which is consistent with standard practice and		use three bags for transporting the lung
	policy.	OPTN policy.		organ. The BAROguard
	The nested bag assembly and all connections	The cleared, off-the-shelf organ isolation		has undergone functional
	are provided sterile.	bags are provided sterile, and neither the		testing and
	Sterile components are gamma radiation sterilized to SAL 10-6	OEM's.		demonstrate that they achieve their intended
	All stocile components that some into discot			effect and meet the same
	contact or indirect via fluid contact with the donor lungs have been tested in accordance			bagging with off-the-shelf bags in the cleared
	with FDA recognized standards.			device.
				(See 018_Updated Laboratory Bench Testing, subsection named Design
				Containment Assembly for testing conducted on the nested bags)
Phase Change Material	Phase change material cold packs to maintain 4°-8° C temperature range for the intended transport time.	Phase change material cold packs to maintain 4°-8° C temperature range for the intended transport time.	N/A	Identical to predicate (substantially equivalent)
	SherpaCool made of phase change material:	SherpaCool made of phase change material:		
	PCM Manufacturer: savERNG	PCM Manufacturer: savERNG		
	Cold Packs by Akuratemp (formerly RGEES), LLC	Cold Packs by Akuratemp (formerly RGEES), LLC		

Table 1. Sub	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard	SherpaPak Lung Preservation System (K192869)	Auto CPAP System	Comparison
	(K223874)		(K211155)	
Pressure	12-15 cmH ₂ O (±1 cmH ₂ O) preset within device,	Standard of care static pressure of 12-15	4-20 cm H2O (in 0.5	Substantially equivalent to
Range	not user adjustable	cmH ₂ Othe lungs should not be over- or	cm H2O increments),	predicate
		underinflated."3	≤30 cm H2O under	
			single fault conditions	Paragonix observed
				pressures between 4-35
				cmH ₂ O during simulated
				air transportation with the
				predicate device and has
				thus designed the subject
				device to address the
				excursions outside of the
				ISHLT-recommended
				range. The subject device
				has a tighter airway
				pressure range than
				reference device. (See
				discussion following table
				on Substantial
				Equivalence in Operating
				Principle)

³ Copeland et al., Donor heart and lung procurement: A consensus statement, J Heart Lung Transplant 2020;39:501–517

Table 1. Subs	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard (K223874)	SherpaPak Lung Preservation System (K192869)	Auto CPAP System (K211155)	Comparison
Air source	Air in the operating room where the retrieved lungs were placed in the BAROguard	Donor lungs inflated and closed within operating room prior to loading into device	Room air delivered through non-sterile device with a blower	Substantially equivalent to predicate
	This air in the internal chamber of the non- sterile BAROguard Shipper is drawn by the			Similar to reference device. Both systems that
	pneumatic pump and delivered to donor lungs through sterile organ containment assembly			provide air to airway are non-sterile. In addition, the
	only when needed.			subject device has a filter with higher filtration
				efficiency than the
				reference device. Further,
				delivered to the lungs by
				the subject device is
				to the reference device.
				Therefore, the air source
				does not raise a different
				question of safety or effectiveness.
Air Filter Efficiency	0.2 Micron Filter has a log-reduction value (LRV) of B. diminuta greater than 8	N/A	Greater than 20% retentive for 10	Substantially equivalent to predicate
			micron particulates	See the prior discussion
Air	4°C-8°C	4°C-8°C	5°C-35°C	Identical to predicate
icilperature				Comparable air temperature range as
				(Substantially equivalent)

Table 1. Sub	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard	SherpaPak Lung Preservation System (K192869)	Auto CPAP System	Comparison
	(K223874)		(K211155)	
Operating	Sea Level to 8000ft (750-1015 hPa)	Sea Level to 8000ft (750-1015 hPa)	760 - 1060 hPa	Identical to predicate
ricssuic				Comparable operating pressure range as
				(Substantially equivalent)

Table 1. Sub	Table 1. Substantial Equivalence Comparison Table			
Characteristic	Subject Device	Predicate Device	Reference Device	Substantial Equivalence
	BAROguard (K223874)	SherpaPak Lung Preservation System (K192869)	Auto CPAP System (K211155)	Comparison
	BlueTooth Low Energy Enabled Data Logger:	BlueTooth Low Energy Enabled Data Logger:		Substantially equivalent
	 Off-the-shelf Data Logger from Onset Computer Corporation 	 Off-the-shelf Data Logger from Onset Computer Corporation 		Monitoring of temperature is identical to predicate.
	 Monitors temperature and pressure BlueTooth capability to transfer 	 Monitors temperature BlueTooth capability to transfer 		Subject device incorporates a different
	temperature reading to iOS or Android device with HoboConnect Ann	temperature reading to iOS or Android device with InTemp App		model of data logger from the same vendor that
	Range: -20° to 70° C (-4° to 158° F)	7 77		allows for monitoring both temperature and pressure.
		Range: -30° to 70° C (-22° to 158° F)		Similar pressure accuracy
	50° C (± 0.4° F from 32° to 122° F)	Temperature Accuracy: ± 0.5° C from 0° to 50° C (± 0.9° F from 32° to 122° F)		(See 018_Updated
	riessure Accuracy: ± 0.3% of reading	Pressure Accuracy: N/A		subsection named Design Verification of BAROquard
Monitoring	Time Accuracy: ±1 minute/month	•		Shipper, for testing conducted on the data
	Optional downloadable mobile app:	Time Accuracy: ±1 minute/month	Pressure Accuracy:	logger.)
	Allows for viewing temperature, pressure, and time elapsed over the mobile device.	Optional downloadable mobile app:	± (0.5 hPa +4%)	
	:	Allows for viewing temperature and time		
	The app does not support, supplement, and/or augment the performance of the parent device.	elapsed over the mobile device.		
	The BAROguard device is always accompanied	The app does not support, supplement,		
	by medical personnel during transport and the	and/or augment the performance of the		
	medical personnel can view information directly	parent device. The predicate is always		
		transport and the medical personnel can view		
		information directly on the datalogger without use of the app.		

005_Updated 510k Summary

Paragonix Technologies' BAROguard Device (K223874)

Conclusion

By design, the subject device operates identically to the predicate when there is not a change in the ambient pressure and addresses a potential limitation associated with the use of the predicate device (and traditional ice coolers used within the Standard of Care) when there is a change in the ambient pressure. Testing confirms that the BAROguard meets the same performance specifications as the predicate device and also meets the requirements of maintaining lung airway pressure during air transportation.

The Airway Pressure Management System in the BAROguard device is similar to the pressure control systems reference the Auto CPAP System (K211155) for the maintenance of continuous positive airway pressure within the donor lung. The reference device is used in the same anatomical location (the airway) and for the same physiological purpose (maintaining a constant, set pressure within that airway above ambient atmospheric pressures) as the BAROguard device.

Based on the design and testing conducted, the subject BAROguard device is substantially equivalent to the predicate device.