

June 16, 2023

Chess Medical Inc. % Christopher Sloan President Sloan Regulatory Consulting 322 Hart Road Gaithersburg, MD 20878

Re: K223890

Trade/Device Name: NAJA Gastrointestinal Catheter

Regulation Number: 21 CFR 876.5980

Regulation Name: Gastrointestinal tube and accessories

Regulatory Class: Class II Product Code: FGD Dated: May 16, 2023

Received: May 16, 2023

#### Dear Christopher Sloan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling

(21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-safety/medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Sivakami Venkatachalam -S

for
Shanil P. Haugen, Ph.D.
Assistant Director
DHT3A: Division of Renal, Gastrointestinal,
Obesity and Transplant Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

### **Indications for Use**

510(k) Number (if known)

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023

Expiration Date: 06/30/2023
See PRA Statement below.

K223890
Device Name NAJA Gastrointestinal Balloon Catheter
11. 11. Gusti Omicistima i Bantoni Gunivici
Indications for Use (Describe) The NAJA Gastrointestinal Balloon Catheter is intended to administer contrast medium to the small intestine.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)
CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

# \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510(k) Summary NAJA Gastrointestinal Balloon Catheter

#### **Submitter**

Name	Chess Medical
	1650 Sherbrooke Street West, Suite 9E
Address	Montreal, Quebec H3H 1C9
	Canada
<b>Contact Person</b>	Dr. Yen-I Chen
Phone	(514) 291-3154
Email Address	cyen33@gmail.com
Date Prepared	December 23, 2022

#### Device

Trade Name	NAJA Gastrointestinal Balloon Catheter
Common Name	Gastrointestinal Balloon Catheter
Classification Name	Gastrointestinal Tube and Accessories
Classification Number	21 CFR 876.5980
<b>Product Code</b>	FGD
Regulatory Class	II

#### **Predicate Device**

Name [510(k) Number]	Maglinte Enteroclysis Balloon Catheter [K884379]
----------------------	--

# **Device Description**

The NAJA Gastrointestinal Balloon Catheter ("NAJA device") is a single-use, single-operator, over-the-wire, double balloon catheter that is compatible with 0.035" guidewires. The device is introduced orally and designed for temporary occlusion of the gastrointestinal (GI) lumen (small intestine), allowing infusion of fluids for enteroclysis studies (i.e., imaging tests of the small intestine). The device is designed to fit inside the working channel of a gastroscope and reaches its intended location by tracking over a 0.035" guidewire.

The NAJA device features a separate hub connector that is attachable and detachable from the catheter. This feature allows removal of the gastroscope from the catheter once the balloons are in the correct location and before the hub connector is attached for balloon inflation and infusion of fluids into the GI lumen. The hub connector is supplied with 1-way stopcocks on the inflation

and infusion ports. The NAJA device is also provided with a 30cc syringe for balloon inflation and infusion of fluids. The balloons are inflated independently using air.

#### **Indications for Use**

The NAJA Gastrointestinal Balloon Catheter is intended to administer contrast medium to the small intestine.

# Comparison of Intended Use and Technological Characteristics with the Predicate Device

The proposed NAJA device has the same intended use and similar technological characteristics as the Maglinte Enteroclysis Balloon Catheter ("Maglinte Catheter"). The differences in the design between the devices do not raise different questions of safety and effectiveness. The table below summarizes the comparison between the NAJA device and Maglinte Catheter.

# **Comparison of NAJA Device to Maglinte Catheter**

	Attribute	NAJA Device (Subject Device)	Maglinte Catheter (Predicate Device)
510(k) Num	ber	To be assigned	K884379
Applicant		Chess Medical	Lafayette Pharmacal
Classificatio	n Regulation;	876.5980 (gastrointestinal tube and	876.5980 (gastrointestinal tube and
Device Class	S	accessories); class II	accessories); class II
Product Cod	e	FGD	FGD
(Device Typ	e)	(Catheter, Retention, Barium Enema	(Catheter, Retention, Barium Enema
		With Bag)	With Bag)
Intended Use	e / Indications for Use	intended to administer contrast	intended to administer contrast
		medium to the small intestine.	medium to the small intestine.
Sterile / Sing	gle Use	Yes	Yes
Sterilization	Method	Sterile (Ethylene Oxide)	Sterile (Ethylene Oxide)
Principle of	Operation	Inflation of one or two compliant	Inflation of one compliant balloon for
		balloons for occlusion of the small	occlusion of the small intestine,
		intestine, followed by infusion of	followed by infusion of contrast
		contrast medium in the small	medium in the small intestine, distal to
		intestine, distal to the proximal	the balloon.
		balloon or in the inter-balloon	
		space.	
Method of D	Delivery	Oral insertion, through a compatible	Nasal or oral insertion.
		working channel (≥ 3.7 mm) of a	
		gastroscope.	
su	Number of Balloons	Two	One
Balloons	Balloon Inflation Media	Air	Air
Ba	Balloon Diameter	20mm – 50mm	30 mm

	Attribute	NAJA Device (Subject Device)	Maglinte Catheter (Predicate Device)
	Balloon Length	36 mm – 48 mm	28 mm
	Inflation Volume	8 ml – 76 ml	15 ml
	Balloon Material	Polyurethane	Latex
	Ballooli Waterial	1 ory dictilatio	Latex
	Number of Marker	Five	One
sp	Bands	(1) Distal tip, (4) proximal and distal ends of the two balloons	(1) Distal tip
Marker Bands	Marker Band Length	Distal tip – 4 mm Proximal/Distal balloon ends – 2 mm	Distal tip – 5 mm
X	Marker Band Material	Polyether Block Amide (Pebax®) 5533 SA 01 MED w/ Tungsten 80%	Metallic (material unknown)
	Working Length	260 cm	155 cm
	Number of Layers	Two Interior – Multi-Lumen Extrusion Exterior – Braided Shaft Extrusion	One Bi-Lumen Extrusion
	Number of Lumens	Four (1) Guidewire lumen (1) Infusion lumen (2) Individual balloon inflation lumens	Two (1) Guidewire lumen, doubles as infusion lumen (1) Balloon inflation lumen
	Number of Inflation Exit Ports	Two One inflation exit port per balloon	Two Two inflation exit ports to inflate a single balloon
	Number of Infusion Exit	Three	Four
	Ports	In between the two balloons	Distal to the single balloon
Shaf	Waterials	Interior – Multi-Lumen Extrusion Polyether Block Amide (Pebax®)  Exterior – Braided Shaft Extrusion Polyether Block Amide (Pebax®)	Radiopaque Polyvinyl Chloride
	Outer Diameter	Distal 30cm (Multi-Lumen  Extrusion)  Outer Diameter: 7.5 Fr  Proximal 206cm (Braided Shaft  Extrusion)  Outer Diameter: 10 Fr	Outer Diameter: 13 Fr
	Distal Tip	Open lumen at distal tip, guidewire can advance	Closed distal tip, guidewire cannot advance
	Proximal End	Enclosed within a removable hub	Enclosed within a permanent hub

Attribute		NAJA Device (Subject Device)	Maglinte Catheter (Predicate Device)
	Detachable	Yes	No
		Able to attach/detach proximal hub	Permanently fixed to proximal end of
dt.		to/from proximal end of shaft	shaft
Proximal Hub	Number of Ports	Four	Two
ma		(1) Guidewire port	(1) Guidewire port, doubles as infusion
rox		(1) Infusion port	port
Ā		(2) Individual balloon inflation ports	(1) Balloon inflation port
	Mechanism to maintain	One-way Stopcock	Check Valve
	inflation volume		
Gu	idewire Compatibility	0.035in	0.065in

# **Biocompatibility Testing**

The following biocompatibility testing was performed on the NAJA device.

Biocompatibility Test	Acceptance Criteria Origin	Result	Conclusion
Cytotoxicity Study Using the ISO Elution Method	ISO 10993-5: 2009, Biological Evaluation of Medical Devices – Part 5: Tests for in vitro cytotoxicity	The test article extract showed no evidence of causing cell lysis or toxicity. The test article extract met the requirements of the test since the grade was less than a grade 2 (mild reactivity).	Pass
ISO Intracutaneous Study	ISO 10993-23: 2021, Biological Evaluation of Medical Devices – Part 23: Tests for irritation	The test article met the requirements of the test since the difference between each test article extract overall mean score and corresponding control extract overall mean score was 0.0 and 0.0 for the SC and SO test article extracts, respectively.	Pass
ISO Maximization Sensitization Study	ISO 10993-10: 2021, Biological Evaluation of Medical Devices – Part 10: Tests for skin sensitization	The test article extracts showed no evidence of causing delayed dermal contact sensitization in the guinea pig. The test article was not considered a sensitizer in the guinea pig maximization test.	Pass

# **Bench and Performance Testing**

The following bench and performance testing were conducted on the NAJA device.

Design Verification Test	esign Verification Test Specification	
Package Performance Testing		
Legibility of Markings	IFU, pouch label, and carton label are legible.	Pass

Design Verification Test	Specification	Test Result
Seal Width	Pouch seal width is ≥ 0.19".	Pass
Pouch Seal Visual Inspection	Pouch seals are free of channels or other damages that could impair sterility.	Pass
Bubble Emission Test	No evidence of leakage through any surface of pouch.	Pass
Pouch Seal Strength	Pouch seal can withstand a force $\geq 3.35$ N over a seal length of 25.4mm.	Pass
Pouch Ease of Opening	Pouch is easily opened with gloved hands, while keeping internal contents stable.  Components inside pouch are easily removable from backboard.	Pass
Component Visual Inspection	Device components are free of damages that could affect performance of product, such as: kinks, cracks, separated components, etc.	Pass
<b>Device Performance Testing</b>		
Visual Inspection		
Surface Defects	All surfaces and open channels of device are free of loose foreign matter, pores, cracks, and remainders of tooling agents.	Pass
Atraumatic Tip	Distal tip is smooth, rounded, and free from sharp edges.	Pass
Atraumatic Infusion Holes	Infusion holes are smooth and free from sharp edges.	Pass
Radiopaque Markers	Five (5) markers are present on device.	Pass
Dimensional Analysis		
0.035" GW Loading – Front/Back	Device accommodates front- and back-loading of a 0.035" GW.	Pass
Proximal Extrusion Dimensions	Dimensions suitable for insertion into hub connector.	Pass
10.5Fr Channel Compatibility	Distal 30cm of device fits through a 10.5Fr channel.	Pass
Working Length	Device working length is $2360 \pm 20$ mm.	Pass
Tip Length	Tip length is $20 \pm 2$ mm.	Pass
Scaffold Length	Scaffold length is $120 \pm 2$ mm.	Pass
Catheter Outer Diameter	Outer diameter is $\leq 3.5$ mm for the catheter shaft, balloon welds, distal tip.	Pass
Simulated Use		
Number of Balloons and Size	Device has two balloons that can be inflated to a diameter of 50mm.	
Hub Connector Detachability	Hub connector is detachable from device and can be attached for fluid delivery.	Pass
Shaft Insertion Feature	Shaft marker is no longer visible once inserted into hub connector.	

<b>Design Verification Test</b>	Specification	Test I	Result
Gastroscope Compatibility	Device fits through working channel of gastroscope.		
Simulated Use Testing	The device remains functional throughout the lifecycle of one clinical procedure.		
Unique Hub Connector Identification	Proximal terminations of hub assembly are appropriately identified. Proximal terminations of the hub assembly are appropriately identified.		
Legibility of Device Markings	Wording on hub assembly labels is legible.		
<b>Functional Testing</b>			
Hub Connector and Hub Assembly Pressure	Pressure decay specification is met.	Pa	ss
Balloon Pressure and Crosstalk	Pressure decay specification is met.	Pa	iss
Balloon OD/L Ratio – unconfined	Unconfined balloon shape specification is met.	Pa	iss
Balloon Diameter and Inflation Volume Relationship	Relationship between balloon diameter and inflation volume shall be determined.	Balloon C Volume [mL] 9 20 44 76	ompliance OD [mm] 20 30 40 50 (MAX)
Balloon OD/L Ratio – confined	Confined balloon shape specification is met.	Pa	iss
Inflation Time	Injection of 30mL of air in ≤ 4s.	Pa	ISS
Deflation Time	Removal of 30mL of air in $\leq$ 17s.	Pa	ISS
Infusion Time	Infusion of 30mL of water in ≤ 19s.	Pa	iss
Channel Occlusion	Device is able to occlude channels 20-50mm.	Pass	
<b>Destructive Testing</b>			
Rated Burst Volume	Balloons have a rated burst volume (RBV) ≥ 81.5mL.	Pa	ss
Bend Radius	No shaft kink and/or fracture when wrapped 180° around a 20mm bend radius.	Pass	
Shaft Buckling Force	Shaft buckle force specification is met.	Pa	iss
Hub to Shaft Slip Force	Hub connector slip force specification is met.	Pa	iss
Tensile – Proximal Shaft to Extrusion	Tensile force of proximal shaft to extrusion is $\geq 15$ N.	Pass	
Tensile – Distal Shaft to Extrusion	Tensile force of distal shaft to extrusion is $\geq 15$ N.	Pa	iss

<b>Design Verification Test</b>	Specification	Test Result
Tensile –Extrusion at Scaffold	Tensile force of scaffold shall be $\geq 15$ N.	Pass
Shear – Balloon Distal Neck	Shear force of distal balloon neck to extrusion bond is $\geq 15$ N.	Pass
Peel – Balloon Proximal Neck	Peel force of proximal balloon neck to extrusion bond is $\geq 15$ N.	Pass
Tensile – Extension Tube to Hub Connector	Tensile force of extension tube to hub connector bond is $\geq 15$ N.	Pass
Tensile – Luer to Extension Tube	Tensile force of luer to extension tube bond is $\geq 15$ N.	Pass

#### **Conclusions**

The NAJA Gastrointestinal Balloon Catheter (subject device) is substantially equivalent to the legally marketed Maglinte Enteroclysis Balloon Catheter (predicate device) for administering contrast medium to the small intestine. The two devices have the same intended use and similar technological characteristics. The differences in the design between the devices do not raise different questions of safety and effectiveness. The results of biocompatibility testing confirm that the materials of construction of the NAJA device are safe, and the results of bench and performance (design verification) testing demonstrate that the NAJA device can function adequately in accordance with its intended use.

Chess Medical Inc. concludes that the NAJA Gastrointestinal Balloon Catheter (subject device) is substantially equivalent to the Maglinte Enteroclysis Balloon Catheter (predicate device).