

April 4, 2023

Neuronetics Amanda Pentecost, Ph.D. Regulatory Affairs Specialist 3222 Phoenixville Pike Malvern, Pennsylvania 19355

Re: K230029

Trade/Device Name: NeuroStar Advanced Therapy System Regulation Number: 21 CFR 882.5805 Regulation Name: Repetitive Transcranial Magnetic Stimulation System Regulatory Class: Class II Product Code: OBP, QCI Dated: January 4, 2023 Received: January 4, 2023

Dear Dr. Pentecost:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <u>https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</u>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Pamela D. Scott -S

Pamela D. Scott
Assistant Director
DHT5B: Division of Neuromodulation and Physical Medicine Devices
OHT5: Office of Neurological and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# Indications for Use

510(k) Number *(if known)* K230029

Device Name NeuroStar Advanced Therapy System (Version 3.7)

Indications for Use (Describe)

The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.

The NeuroStar Advanced Therapy system is intended to be used as an adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary			
Date Prepared:	April 4, 2023		
Applicant:	Neuronetics, Inc. 3222 Phoenixville Pike Malvern, PA 19355		
<u>Contact Person:</u>	Amanda Pentecost, PhD Regulatory Affairs Specialist Phone: 610-640-4202, ext. 1132 Email: <u>amanda.pentecost@neurostar.com</u>		
Secondary Contact:	Robin Fatzinger Sr. Director, Regulatory Affairs Phone: 610-640-4202, ext. 1064 Email: <u>robin.fatzinger@neurostar.com</u>		
<u>Device Trade Name:</u>	NeuroStar NeuroStar TMS Therapy System NeuroStar Advanced Therapy System NeuroStar Advanced Therapy System for Mental Health		
Device Common Name:	Transcranial Magnetic Stimulator		
<b>Classifications:</b>	21 CFR 882.5805, 21 CFR 882.5802		
Product Codes:	OBP, QCI		
Primary Predicate Device:	NeuroStar Advanced Therapy System (K083538, K130233, K133408, K160703, K161519, K201158, and K220127) (Product Code: OBP)		
Secondary Predicate Device:	NeuroStar Advanced Therapy System (K212289) (Product Code: QCI)		





## **Device Description / Technological Characteristics:**

The NeuroStar Advanced Therapy System is a transcranial magnetic stimulation device. Specifically, it is a computerized, electromechanical medical device that produces and delivers non-invasive magnetic fields to induce electrical currents targeting specific regions of the cerebral cortex. Transcranial magnetic stimulation (TMS) is a non-invasive technique used to apply brief magnetic pulses to the brain. The pulses are administered by passing high currents through an electromagnetic coil placed adjacent to a patient's scalp. The pulses induce an electric field in the underlying brain tissue. When the induced field is above a certain threshold and is directed in an appropriate orientation relative the brain's neuronal pathway, localized axonal depolarizations are produced, thus activating neurons in the targeted brain region.

The NeuroStar Advanced Therapy System is an integrated system consisting of a combination of the following components:

- Mobile Console
- System Software
- Treatment Chair
- Ferromagnetic Treatment Coil
- Head Support System
- SenStar® Connect Treatment Link & SenStar® Treatment Link
- Treatment Pack (for use with the SenStar® Connect Treatment Link)
- MT Cap
- TrakStar<sup>™</sup> Patient Data Management System
- D-Tect<sup>™</sup> MT Accessory
- Beam F3 Treatment Pack

There are two proposed changes to the NeuroStar Advanced Therapy System that are the subject of this 510(k). The first proposed change introduces the capability to use wireless communication to transfer data between the NeuroStar System software and the TrakStar Patient Data Management System, as an alternative to using an ethernet cable to facilitate this transfer. The second proposed change allows for the alternative use of the Beam F3 treatment location method (herein: "Beam F3 method"), in addition to the current method that utilizes a location 5 cm away from the motor threshold within the dorsolateral prefrontal cortex (DLPFC) region of the brain. The Beam F3 method requires measuring different skull dimensions and uses these values to calculate the F3 treatment location, which is also located within the DLPFC. This 510(k) introduces software embedded into the NeuroStar System that performs these calculations as well as single-use measuring accessories. Once the F3 treatment location is determined, the treatment parameters and protocols remain the same as those used in the currently marketed NeuroStar Advanced Therapy System.

## **Indications for Use:**

The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in Page 2 of 12



adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.

The NeuroStar Advanced Therapy system is intended to be used as an adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD).

### **Performance Standards:**

The NeuroStar Advanced Therapy System has been tested and conforms to the following recognized consensus standards:

- ISO 10993-1
- ANSI AAMI ES60601-1 / IEC 60601-1
- IEC 60601-1-2
- IEC 60601-1-6

### **Non-clinical Testing:**

The contents of this 510(k) complies with the FDA Guidance Document: "Class II Special Controls Guidance Document: Repetitive Transcranial Magnetic Stimulation (rTMS) Systems - Guidance for Industry and Food and Drug Administration Staff." Non-clinical performance testing was performed according to the standards listed above.

Additionally, usability testing was completed in accordance with IEC 60601-1-6 and also following FDA Guidance Document: "Applying Human Factors and Usability Engineering to Medical Devices." Software was designed, developed, and tested in accordance with relevant FDA guidance documents, IEC 62304, and ISO 14971.

#### **Clinical Testing:**

There is no clinical testing required to support this submission.





<b>Technological Com</b>	parison with Primar	y Predicate (	Product Code: OBP):

	NeuroStar Advanced Therapy System (Subject Device)	NeuroStar Advanced Therapy System (Primary Predicate Device - K083538, K130233, K133408, K160703, K161519, K201158, and K220127)	Explanation of Differences
Indications for Use	The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.	The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.	No Difference
Intended Use	Major Depressive Disorder (MDD) and comorbid anxiety symptoms	Major Depressive Disorder (MDD) and comorbid anxiety symptoms	No Difference
Anatomical Sites	Left dorsolateral prefrontal cortex	Left dorsolateral prefrontal cortex	No Difference
Target Population	Adult patients	Adult patients	No Difference
Clinical Setting	Inpatient and outpatient settings, including physician's offices and clinics, hospitals, and general medical/surgical hospitals	Inpatient and outpatient settings, including physician's offices and clinics, hospitals, and general medical/surgical hospitals	No Difference
Materials	Standard materials commonly used in the manufacture of electrical medical devices	Standard materials commonly used in the manufacture of electrical medical devices	No Difference
Biocompatibility	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	No Difference
Energy Source	Power console with magnetic coil for delivery for magnetic energy	Power console with magnetic coil for delivery for magnetic energy	No Difference





t 610 640 4202 / f 610 640 4206

Electrical Safety & EMC	IEC 60601-1 compliant IEC 60601-1-2 compliant	IEC 60601-1 compliant IEC 60601-1-2 compliant	No Difference
Communication with TrakStar	Wireless (Wi-fi) and Ethernet cable	Ethernet cable	Different The addition of the wireless capability is part of the subject of this 510(k).
Sterility	No parts of the device, accessories or components are required to be sterilized	No parts of the device, accessories or components are required to be sterilized	No Difference
Coil Type	Ferromagnetic Iron core Internal cooling fan	Ferromagnetic Iron core Internal cooling fan	No Difference
Coil Positioning System	Integrated into Head Support System Laser-aided coil placement	Integrated into Head Support System Laser-aided coil placement	No Difference
Treatment Schedule	5 days per week for 6 weeks with taper over 3 weeks (3 sessions first week, 2 sessions second week and 1 session third week) for total of 36 treatment sessions.	5 days per week for 6 weeks with taper over 3 weeks (3 sessions first week, 2 sessions second week and 1 session third week) for total of 36 treatment sessions.	No Difference





t 610 640 4202 / f 610 640 4206

Device Components	<ul> <li>Mobile Console</li> <li>Ferromagnetic Coil for delivering treatment</li> <li>Head Support System for coil positioning</li> <li>MT Cap for coil positioning</li> <li>D-Tect<sup>™</sup> MT Accessory for MT location and level determination</li> <li>Multi-use disposable for contact sensing and magnetic field quality control</li> <li>Single-use treatment pack including disposable hygienic barriers and coil positioning head strap for use with the standard 5 cm method</li> <li>Single-use treatment pack including disposable hygienic barriers and head strap for use with the standard 5 cm method</li> <li>Single-use treatment pack including disposable hygienic barriers and head strap for use with the Beam F3 method for determining treatment location and coil positioning</li> <li>TrakStar System for recording patient data</li> </ul>	<ul> <li>Mobile Console</li> <li>Ferromagnetic Coil for delivering treatment</li> <li>Head Support System for coil positioning</li> <li>MT Cap for coil positioning</li> <li>D-Tect<sup>™</sup> MT Accessory for MT location and level determination</li> <li>Multi-use disposable for contact sensing and magnetic field quality control</li> <li>Single-use treatment pack including disposable hygienic barriers and coil positioning head strap for use with the standard 5 cm method</li> <li>TrakStar System for recording patient data</li> </ul>	Different The addition of the Beam F3 Treatment Pack is part of the subject of this 510(k).
%MT Range	25% to 140% MT	25% to 140% MT	No Difference
Pulses per Second (PPS) Range	For treatment: 1 to 30 PPS For MT determination: 0.1-0.3 PPS	For treatment: 1 to 30 PPS For MT determination: 0.1-0.3 PPS	No Difference
Induced Electrical field at 2 cm at 1.0 SMT	135 V/m (Nominal)	135 V/m (Nominal)	No Difference
Pulse Type	Biphasic Sinusoid	Biphasic Sinusoid	No Difference
Pulse Width	185 μS (Nominal)	185 μS (Nominal)	No Difference



**t** 610 640 4202 / **f** 610 640 4206

NeuroStar.com

Treatment Protocols	Standard Treatment:	Standard Treatment:	
incutinent i rotocolo	Level: 120% MT with allowable adjustments	Level: 120% MT with allowable adjustments	
	Repetition Rate: 10 PPS	Repetition Rate: 10 PPS	
	Stimulation Time: 4 s	Stimulation Time: 4 s	
	Inter-train Interval: As low as 11 s	Inter-train Interval: As low as 11 s	
	Session Duration: As low as 18.75 min	Session Duration: As low as 18.75 min	
	Pulses per Session: 3000	Pulses per Session: 3000	
	Sessions per Week: 5 (Acute phase)	Sessions per Week: 5 (Acute phase)	
	NeuroBurst Treatment:	NeuroBurst Treatment:	
	Level: 80-120% MT with allowable	Level: 80-120% MT with allowable	No Difference
	adjustments	adjustments	
	Stimulation Time: 2 s	Stimulation Time: 2 s	
	Inter-train Interval: 8 s	Inter-train Interval: 8 s	
	Pulses per Burst: 3	Pulses per Burst: 3	
	Interpulse Interval: 20 ms	Interpulse Interval: 20 ms	
	Session Duration: 3.3 min	Session Duration: 3.3 min	
	Pulses per Session: 600	Pulses per Session: 600	
	Bursts per Second: 5	Bursts per Second: 5	
	Amplitude: 0.22-2.08 SMT ( <u>&lt;</u> 5% drop)	Amplitude: 0.22-2.08 SMT ( <u>&lt;</u> 5% drop)	
Treatment Level Range	Standard Treatment:	Standard Treatment:	
freutinent Level Kunge	0.22 to 2.08 SMT	0.22 to 2.08 SMT	
	Calibrated linear output	Calibrated linear output	
	NeuroBurst Treatment:	NeuroBurst Treatment:	No Difference
	0.22 to 1.9 SMT	0.22 to 1.9 SMT	
	80-120% MT	80-120% MT	
	<u>&lt;</u> 5% drop	<u>&lt;</u> 5% drop	
Stimulation Time Pulse Train	Standard Treatment:	Standard Treatment:	
Duration Range	1 PPS: 1 to 600 s	1 PPS: 1 to 600 s	
Duration Range	> 1 PPS: 1 to 20 s	> 1 PPS: 1 to 20 s	No Difference
	NeuroBurst Treatment:	NeuroBurst Treatment:	
	1 to 10 s	1 to 10 s	



**NeuroS** 

herapy fo

Advanced



**t** 610 640 4202 / **f** 610 640 4206

Inter-train Interval Range	Standard Treatment:	Standard Treatment:	
inter-train interval kange	1 PPS: 0 to 600 s	1 PPS: 0 to 600 s	
	>1 PPS: 10 to 60 s	>1 PPS: 10 to 60 s	No Difference
	NeuroBurst Treatment:	NeuroBurst Treatment:	
	1 to 60 s	1 to 60 s	
Pulse per Treatment Session	Standard Treatment:	Standard Treatment:	
Puise per freatment Session	Nominal: 3000	Nominal: 3000	
	Maximum: 5000	Maximum: 5000	No Difference
	NeuroBurst Treatment:	NeuroBurst Treatment:	NO Difference
	Nominal: 600	Nominal: 600	
	Maximum: 2000	Maximum: 2000	
Pulses per Burst (PPB)	NeuroBurst Treatment:	NeuroBurst Treatment:	No Difference
ruises per buist (FFB)	1 to 5	1 to 5	No Difference
	NeuroBurst Treatment:	NeuroBurst Treatment:	
Interpulse Interval	20 to 2000 ms	20 to 2000 ms	No Difference
Bursts per Second (BPS)	NeuroBurst Treatment:	NeuroBurst Treatment:	No Difference
	0.1 to 20.0 Hz	0.1 to 20.0 Hz	No Difference





NeuroStar.com

	NeuroStar Advanced Therapy System (Subject Device)	NeuroStar Advanced Therapy System (Secondary Predicate Device - K212289)	Explanation of Differences
Indications for Use	The NeuroStar Advanced Therapy System is intended to be used as an adjunct for the treatment of adult patients suffering from Obsessive-Compulsive Disorder (OCD)	The NeuroStar Advanced Therapy System is intended to be used as an adjunct for the treatment of adult patients suffering from Obsessive- Compulsive Disorder (OCD)	No Difference
Intended Use	Obsessive Compulsive Disorder	Obsessive Compulsive Disorder	No Difference
Anatomical Sites	Bilateral dorsomedial prefrontal cortex	Bilateral dorsomedial prefrontal cortex	No Difference
Target Population	Adult patients (ages 22-70) with Obsessive- Compulsive Disorder	Adult patients (ages 22-70) with Obsessive- Compulsive Disorder	No Difference
Clinical Setting	Inpatient and outpatient settings including physician's offices and clinics, hospitals, and general medical/surgical hospitals.	Inpatient and outpatient settings including physician's offices and clinics, hospitals, and general medical/surgical hospitals.	No Difference
Materials	Standard materials commonly used in the manufacture of electrical medical devices	Standard materials commonly used in the manufacture of electrical medical devices	No Difference
Biocompatibility	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	No Difference
Energy Source	Power console with magnetic coil for delivery for magnetic energy	Power console with magnetic coil for delivery for magnetic energy	No Difference
Electrical Safety & EMC	IEC 60601-1 compliant IEC 60601-1-2 compliant	IEC 60601-1 compliant IEC 60601-1-2 compliant	No Difference
Communication with TrakStar	Wireless (Wi-fi) and Ethernet cable	Ethernet cable	Different The addition of the wireless capability is part of the subject o this 510(k).
Sterility	No parts of the device, accessories or components are required to be sterilized	No parts of the device, accessories or components are required to be sterilized	No Difference

# Technological Comparison with Primary Predicate (Product Code: QCI):

Page 9 of 12



**t** 610 640 4202 / **f** 610 640 4206

Coil Type	Ferromagnetic	Ferromagnetic	
con type	Iron core	Iron core	No Difference
	Internal cooling fan	Internal cooling fan	
Coil Positioning System	Integrated into Head Support System	Integrated into Head Support System	No Difference
g_,	Laser-aided coil placement	Laser-aided coil placement	
Treatment Schedule	Weeks 1-5: 1 treatment session per day for 5	Weeks 1-5: 1 treatment session per day for 5 days	
	days	Week 6: 1 treatment session per day for 4 days	No Difference
	Week 6: 1 treatment session per day for 4 days	Total of 29 treatment sessions	
	Total of 29 treatment sessions		- 100
	Mobile Console	Mobile Console	Different
	<ul> <li>Ferromagnetic Coil for delivering treatment</li> </ul>	<ul> <li>Ferromagnetic Coil for delivering treatment</li> <li>Head Support System for coil positioning</li> </ul>	<b>T</b> I I.I (.)
	Head Support System for coil positioning	<ul> <li>Multi-use disposable for contact sensing and</li> </ul>	The addition of the
	<ul> <li>Multi-use disposable for contact sensing</li> </ul>	magnetic field quality control	Beam F3 Treatment
	and magnetic field quality control	Single-use treatment pack including disposable	Pack is part of the
	Single-use treatment pack including	hygienic barriers and coil positioning head strap	subject of this 510(k)
Device Commente	disposable hygienic barriers and coil	for use with the standard 5 cm method	and specific to the
Device Components	positioning head strap for use with the		MDD indication.
	standard 5 cm method		
	Single-use treatment pack including		*Note: The Beam F3
	disposable hygienic barriers and head		Treatment Pack is
	strap for use with the Beam F3 method for		not intended to be
	determining treatment location and coil		used as part of the
	positioning*	TrakStar System for recording patient data	OCD treatment
	TrakStar System for recording patient data		protocol.
%MT Range	25% to 140% MT	25% to 140% MT	No Difference
Pulses per Second (PPS)	For treatment: 1 to 30 PPS	For treatment: 1 to 30 PPS	
Range	For MT determination: 0.1-0.3 PPS	For MT determination: 0.1-0.3 PPS	No Difference
Induced Electrical field at 2			
cm at 1.0 SMT	135 V/m (Nominal)	135 V/m (Nominal)	No Difference
Pulse Type	Biphasic Sinusoid	Biphasic Sinusoid	No Difference





t 610 640 4202 / f 610 640 4206

Pulse Width	185 μS (Nominal)	185 μS (Nominal)	No Difference
Treatment Protocols	Level: 100% foot motor threshold level	Level: 100% foot motor threshold level	
Treatment Protocols	Repetition Rate: 20 PPS	Repetition Rate: 20 PPS	
	Stimulation Time: 2 s	Stimulation Time: 2 s	
	Inter-train Interval: 20 s	Inter-train Interval: 20 s	No Difference
	Session Duration: As low as 18.3 min	Session Duration: As low as 18.3 min	No Difference
	Pulses per Session: 2000	Pulses per Session: 2000	
	Sessions Per Week: 5 for Weeks 1-5 and 4 for	Sessions Per Week: 5 for Weeks 1-5 and 4 for	
	Week 6	Week 6	
Treatment Level Range	0.22 to 2.08 SMT	0.22 to 2.08 SMT	No Difference
freatment Lever Kange	Calibrated linear output	Calibrated linear output	No Difference
Stimulation Time Pulse Train	1 PPS: 1 to 600 s	1 PPS: 1 to 600 s	
Duration Range	> 1 PPS: 1 to 20 s	> 1 PPS: 1 to 20 s	No Difference
Duration Range	> 1 FF3. 1 to 20 5	> 1 FF3. 1 to 20 S	
Inter-train Interval Range	1 PPS: 0 to 600 s	1 PPS: 0 to 600 s	No Difference
inter-train interval Kange	>1 PPS: 10 to 60 s	>1 PPS: 10 to 60 s	No Difference
Pulse per Treatment Session	Nominal: 2000	Nominal: 2000	N D'((
	Maximum: 5000	Maximum: 5000	No Difference





## **Conclusion:**

The NeuroStar Advanced Therapy System and the primary/secondary predicate devices have the same intended uses and technological characteristics. The use of the optional wireless communication or Beam F3 treatment location features do not raise any different questions regarding safety or effectiveness.

