

April 28, 2023

T.A.G. Medical Products Corporation, Ltd Shlomi Dines RA/QA Director T.A.G. Medical Products Corporation, Ltd Gaaton, 2513000 Israel

Re: K230058

Trade/Device Name: Bladeless Trocar - Artemis Lap Cannula Regulation Number: 21 CFR 876.1500 Regulation Name: Endoscope And Accessories Regulatory Class: Class II Product Code: GCJ Dated: March 29, 2023 Received: March 29, 2023

Dear Shlomi Dines:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Mark	Digitally signed by
Trumbore	- Mark Trumbore -S Date: 2023.04.28
S	14:38:20 -04'00'

Mark Trumbore, Ph.D. Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230058

Device Name Bladeless Trocar – Artemis Lap Cannula

Indications for Use (Describe)

The Artemis Lap Cannula has applications in abdominal, thoracic, and gynecologic minimally invasive procedures to establish a path of entry for minimally invasive instruments. Artemis Lap Cannula may be used for primary and secondary insertions.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Prepared on: 2023-04-27

Pursuant to CFR 807.92, the following 510(k) Summary is provided:

Contact Details		
	Applicant Name: Applicant Address:	T.A.G. Medical Products Corporation, Ltd. T.A.G. Medical Products Corporation, Ltd Gaaton 2513000 Israel 972-4-9858400 Gaaton 2513000, ISRAEL <u>www.tag-med.com</u>
	Applicant Contact Telephone:	972-4-9858400
	Applicant Contact:	Shlomi Dines
Device Name	Applicant Contact Email:	sdines@tag-med.com
Device Name	Device Trade Name:	Bladeless Trocar – Artemis Lap Cannula
	Common Name:	Endoscope and accessories
	Classification Name:	Endoscope and accessories
	Regulation Number:	876.1500
Legally Marketed	Product Code:	GCJ
Predicate Devices	Predicate #:	K032676
	Predicate Trade Name:	ENDOPATH III Trocar System
	Product Code:	GCJ
Device Description		Artemis Lap Cannula system is a radiolucent, reusable, bladeless laparoscopic trocar, consisting of a cannula, an obturator, a depth limiter, and a disposable standalone seal pack. The trocar is available in two diameters: Ø5mm and Ø12mm, each consists of 4 different length variants. Depth limiter component is available in two diameters and fits either the Ø5mm or Ø12mm cannula regardless of the length. Depth limiter can be used to prevent over penetration during surgical procedures. Artemis Lap Cannula may be used in abdominal, thoracic, or gynecological procedures.

Indications for Use

Comparison of Technological Characteristics The Artemis Lap Cannula has applications in abdominal, thoracic, and gynecologic minimally invasive surgical procedures to establish a path of entry for minimally invasive instruments. Artemis Lap Cannula may be used for primary and secondary insertions.

The proposed device and the predicate device (K032676) have the same basic design, intended use, and sterilization. In comparison to the predicate device, the proposed modifications include differences in the sleeve and obturator material, packaging configuration, sleeve design, and the inclusion of a depth limiter. Differences between the proposed and predicate device do not raise new questions of safety or effectiveness.

Comparison to Predicate			
Device	Subject Device	Predicate Device	Comparison
Device Description Summary	Artemis Lap Cannula system is a radiolucent, reusable, bladeless laparoscopic trocar, consisting of a cannula, an obturator, a depth limiter, and a disposable standalone seal pack. The trocar is available in two diameters: Ø5mm and Ø12mm, each consists of 4 different length variants. Depth limiter component is available in two diameters and fits either the Ø5mm or Ø12mm cannula regardless of the length. Depth limiter can be used to prevent over penetration during surgical procedures. Artemis Lap Cannula may be used in abdominal, thoracic, or gynecological procedures.	The ENDOPATH III Trocars are sterile single patient use instruments consisting of a radiolucent sleeve and obturator in sizes ranging from 5-12 mm in diameter. There are three different obturators Bladeless, Blunt Tip and Dilating Tip. The Bladeless obturator contains a clear, tapered optical element, which when used with an endoscope, provides visibility of individual tissue layers during insertion. The Bladeless obturator accommodates an appropriately sized zero endoscope. The Blunt Tip	The proposed device and the predicate device have the same basic design, intended use, and sterilization. In comparison to the predicate device, the proposed modifications include differences in the sleeve and obturator material, packaging configuration, sleeve design, and the inclusion of a depth limiter. Differences between the proposed and predicate device do not raise new questions of safety or effectiveness.

	1		
		laceration once the abdominal or thoracic cavity has been entered.	
		The trocar sleeve contains two seals, an outer	
		integrated removable self- adjusting seal to	
		accommodate instruments ranging from 5mm to 12mm	
		in diameter where indicated	
		and an internal seal. Together, these seals	
		minimize gas leakage when	
		instruments are inserted or	
		withdrawn through the trocar. The 5mm trocar	
		sleeve does not contain an	
		integrated removable seal and accommodates only	
		5mm instruments. A	
		stopcock valve is compatible with standard luer lock	
		fittings and provides	
		attachment for gas insufflation.	
		The stopcock is in closed	
		position when it is parallel to	
Indications for Use:	The Artemis Lap Cannula has applications in abdominal,	the sleeve. The ENDOPATH III Bladeless Trocar has applications in	The indications for use are identical to the relevant
	thoracic, and gynecologic minimally invasive surgical	abdominal, thoracic, and	part if the predicate device cleared under
	procedures to establish a path	procedures to establish a path	K032676.
	of entry for minimally invasive instruments. Artemis Lap	of entry for endoscopic instruments. The trocar may be	
	Cannula may be used for primary and secondary	used with or without visualization for primary and	
	insertions.	secondary insertions.	
		The ENDOPATH III Dilating Tip	
		Trocar has applications in thoracic, gynecologic	
		laparoscopy and other abdominal procedures to	
		establish a path of entry for endoscopic instruments.	
		The ENDOPATH III Blunt Tip Trocar has applications in	

		thoracic, gynecologic, laparoscopic and other				
		abdominal procedures to				
		establish a path of entry for				
	Characteri	minimally invasive instruments.				
Outer Seal Design	Characteristics/Features:					
Outer Sear Design	Multi-piece (Pacman), overlapping,	Multi-piece (Pacman), overlapping,	Multi-piece (Pacman), overlapping,			
Inner Seal Design	Duckbill design	Design Duckbill	Equivalent			
Obturator Tip Design	Bladeless	Bladeless	Equivalent			
Sleeve Design	Low-profile design	Low-profile design	Equivalent			
Obturator Design	Low-profile design	Low-profile design	Equivalent			
Cannula Design	Low-profile design	Low-profile design	Equivalent			
Outer Seal Material	Polycarbonate and Polyisoprene	Polycarbonate and Polyisoprene				
Inner Seal Material	Polyisoprene	Polyisoprene	Equivalent			
Dimensions (Diameter)	5 mm & 12 mm	5 mm & 12 mm	Equivalent			
Sterilization	Cobalt, irradiation	Cobalt, irradiation	Equivalent			
Dimensions (Length)	75 mm, 100 mm, 150 mm	75 mm, 100 mm, 150 mm	Equivelent			
Sleeve Material	Radel	Polycarbonate	The proposed device offers the sleeve material in Radel rather than polycarbonate.			
Obturator Material	Radel	Polycarbonate	The proposed device offers the obturator material in Radel rather than polycarbonate.			
Packaging	Flexible Film Composite, with lidding film top stock (FMP- 521®)	Copolyester rigid blister, with heat-sealed Tyvek lid	The proposed device incorporates a different packaging configuration compared to the predicate.			
Sleeve Design	locking of seal pack (includes inner and outer seals)	locking of outer seal	The proposed device offers a release button for the seal pack of the inner and outer seals whereas the predicate offers a release button for the outer seal only.			
Depth Limiter	5 or 12 mm in diameter	Not included with system	The proposed device offers a depth limiter in either 5 or 12 mm in diameter whereas the predicate does not have a depth limiter included in the system.			

Nonclinical Testing Discussion

Biocompatibility evaluation was conducted on the proposed device. It was found biocompatible for intended use.

Reprocessing evaluation was conducted on the proposed device. Validated reprocessing instructions are sufficient to clean and sterilize it in healthcare settings.

Nonclinical testing in accordance with ISO 80639-7 was completed. The test data demonstrates success and met the criteria of ISO 80369-7.

Leak testing was conducted on the proposed devices and submitted in this Traditional 510(k). The leak test data demonstrates the proposed devices preform statically equivalent to the predicate device.

Conclusion

Subject and predicate devices have same intended use. Differences in design between the subject and predicate device do not raise new questions of safety or effectiveness. Based on comparison of the technological characteristics, and performance test data, the subject devices is substantially equivalent to the predicate device for requested intended use.