

April 7, 2023

Mazor Robotics Ltd. Marina Shkedy Sr. Regulatory Affairs Specialist 1 HaEshel Street (Building C) Caesarea Business Park, 3079830 Israel

Re: K230064

Trade/Device Name: Mazor X System (Mazor X Stealth Edition) Regulation Number: 21 CFR 882.4560 Regulation Name: Stereotaxic Instrument Regulatory Class: Class II Product Code: OLO, LLZ Dated: December 29, 2022 Received: January 9, 2023

Dear Marina Shkedy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Shumaya Ali -S

Shumaya Ali, M.P.H.
Assistant Director
DHT6C: Division of Restorative, Repair and Trauma Devices
OHT6: Office of Orthopedic Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# Indications for Use

510(k) Number (*if known*) K230064

Device Name

Mazor X System (Mazor X Stealth Edition)

Indications for Use (Describe)

The Mazor X is indicated for precise positioning of surgical instruments or spinal implants during general spinal surgery. It may be used in open or minimally invasive or percutaneous procedures.

Mazor X 3D imaging capabilities provide a processing and conversion of 2D fluoroscopic projections from standard Carms into a volumetric 3D image. It is intended to be used whenever the clinician and/or patient benefits from generated 3D imaging of high contrast objects.

The Mazor X navigation tracks the position of instruments, during spinal surgery, in relation to the surgical anatomy and identifies this position on diagnostic or intraoperative images of a patient.

Type of Use (Select one or both, as a	applicable)
---------------------------------------	-------------

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510(k) Summary

Applicant name:	Mazor Robotics Ltd.
	1 HaEshel Street (Building C)
	Caesarea business Park 3079830, Israel
Contact person(s):	Primary: Marina Shkedy
	Role: Sr. Regulatory Affairs Specialist
	Email: marina.shkedy@medtronic.com
	Phone: +972507848810
	Secondary: Sophie Gorali
	Role : Regulatory Affairs Specialist
	Email: sophie.gorali@medtronic.com

Date Prepared: April 6th, 2023

Name of Device: Mazor X System (Mazor X Stealth Edition)

Classification Name: Stereotaxic instrument

Classification Code: OLO and LLZ

**Device class:** II

**Regulation number:** 882.4560

Panel: Orthopedic

Predicate Devices: Mazor X System (Mazor X Stealth Edition) (K203005)

# Intended Use / Indications for Use

The Mazor X is indicated for precise positioning of surgical instruments or spinal implants during general spinal surgery. It may be used in open or minimally invasive or percutaneous procedures.

Mazor X 3D imaging capabilities provide a processing and conversion of 2D fluoroscopic projections from standard C Arms into volumetric 3D image. It is intended to be used whenever the clinician and/or patient benefits from generated 3D imaging of high contrast objects.

The Mazor X navigation tracks the position of instruments, during spinal surgery, in relation to the surgical anatomy and identifies this position on diagnostic or intraoperative images of a patient.

## **Device Description**

The Mazor X system combines robotic trajectory guidance with navigated surgical instruments (either guided or free hand navigation) to enable the surgeon to precisely position surgical instruments and/or implants according to predefined planning. With the imaging capabilities of the system, the user can also visualize the implants on the patients CT. Same as the predicate device the modified Mazor X consists of a workstation with dedicated software, the surgical system, navigation camera, accessories, instruments and disposable kits. The modified Mazor X, the subject of this 510(k) application, introduces software and hardware modifications to the Mazor X System cleared in 510(k) K203005.

### **Technological Characteristics**

The *modified* Mazor X device is similar in its technological features to its predicate device, the cleared Mazor X. Both the subject and predicate systems combine *robotic* and *navigation* technologies to enable the precise positioning and tracking of surgical instruments or spinal implants during general spinal surgery. Both systems allow registration of image data with patient anatomy by registering the pre-operative CT scan to the intra-operative fluoroscopic scans, and by localizing the position of the 3D patient volume to the robotic system. In both systems, the positioning of surgical instruments and their trajectories are guided by the system in accordance with the planning conducted by the surgeon on the pre-operative CT image. In both systems, the navigation feature provides the option of tracking the surgical instruments during the workstation, the surgical system, navigation camera, accessories, instruments, and disposable kits. The modifications to the cleared Mazor X System, that are the subject of this premarket application.

include software enhancements to enable extended functionality, compatibility with additional, previously cleared, surgical tools and some minor hardware changes.

However, as explained in more detail below, these differences do not raise new or different questions of safety or effectiveness since the principal technology remains very similar and, in both instances, the key question is whether the robotic and navigation functionalities are accurate. The bench performance testing and human factors validation demonstrated that the modified Mazor X is safe and effective as the predicate device.

Comparison of the Proposed Mazor X to the Cleared Mazor X System (K203005) is provided below:

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Product Code,	OLO and LLZ,	OLO and LLZ,	Identical
Class	Class II	Class II	
Indications for	The Mazor X is indicated for precise	The Mazor X is indicated for precise	Identical
Use	positioning of surgical instruments or	positioning of surgical instruments or	
	spinal implants during general spinal	spinal implants during general spinal	
	surgery. It may be used in open or	surgery. It may be used in open or	
	minimally invasive or percutaneous	minimally invasive or percutaneous	
	procedures.	procedures.	
	Mazor X 3D imaging capabilities	Mazor X 3D imaging capabilities	
	provide a processing and conversion of	provide a processing and conversion of	
	2D fluoroscopic projections from	2D fluoroscopic projections from	
	standard C-arms into a volumetric 3D	standard C-arms into a volumetric 3D	
	image. It is intended to be used	image. It is intended to be used	
	whenever the clinician and/or patient	whenever the clinician and/or patient	

 Table 0-1: Comparison of the Proposed Mazor X to the Cleared Mazor X System (K203005)

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
	benefits from generated 3D imaging of	benefits from generated 3D imaging of	
	high contrast objects.	high contrast objects.	
	The Mazor X navigation tracks the	The Mazor X navigation tracks the	
	position of instruments, during spinal	position of instruments, during spinal	
	surgery, in relation to the surgical	surgery, in relation to the surgical	
	anatomy and identifies this position on	anatomy and identifies this position on	
	diagnostic or intraoperative images of a	diagnostic or intraoperative images of a	
	patient.	patient.	
Target	Orthopedic patients	Orthopedic patients	Identical
Population			
Anatomical	Spine	Spine	Identical
Sites			
Environment	Hospital setting (operating room)	Hospital setting (operating room)	Identical
Used			

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Main system	The Mazor X System consists of the	The Mazor X System consists of the	Identical
components	following components:	following components:	
	• Workstation	• Workstation	
	• Surgical System 3Define	Surgical System 3Define	
	Camera: SR300 model	Camera: SR300 model	
	• Bed Connecting Unit (e.g., Bed Frame)	• Bed Connecting Unit (e.g., Bed Frame)	
	• Device accessories for spine application	• Device accessories for spine application	
	Mazor X Navigation Camera	Mazor X Navigation Camera	
	and accessories	and accessories	
Mechanism of	Computer assisted Stereotaxy:	Computer assisted Stereotaxy:	Identical
Action	Instrument position and trajectory	Instrument position and trajectory	
	calculation based on image data &	calculation based on image data &	
	instrument tracking based on optical	instrument tracking based on optical	

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
	navigation. Motorized positioning of the Surgical Arm (spine) with tool	navigation. Motorized positioning of the Surgical Arm (spine) with tool	
	guide through 6 axes.	guide through 6 axes.	
Compatibility with Medtronic instruments and tools	<ul> <li>Mazor X is compatible with the following Medtronic tools and instruments:</li> <li>Medtronic Spine instruments(K182121)</li> <li>Medtronic Navigation NavLock Trackers (K182104)</li> <li>Medtronic Stealth-Midas and Stealth Midas MR8 navigated drill systems (K160713, K183644)</li> </ul>	<ul> <li>Mazor X is compatible with the following Medtronic tools and instruments:</li> <li>Medtronic Spine instruments(K182121)</li> <li>Medtronic Navigation NavLock Trackers (K182104)</li> <li>Medtronic Stealth-Midas and Stealth Midas MR8 navigated drill systems (K160713, K183644)</li> </ul>	Similar. Compatibility with additional, cleared to market, Medtronic instruments and tools using the same navigation technology, enabling tracking the position of the instruments in relation to the surgical anatomy and identify this position on the patient images

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
	<ul> <li>Medifolic havigated surgical instruments - Trials trackers, Inserters and Disc Prep Instruments (K131425,</li> </ul>	<ul> <li>Meditonic havigated surgical instruments - Trials trackers, Inserters and Disc Prep Instruments (K131425,</li> </ul>	additional, cleared to market, tools, does not change the system
	<ul> <li>K163581, K150231) and Medtronic Pedicle Probe (K050438)</li> <li>Medicrea UNiD Analyzer (K212005)</li> <li>Stealth-Midas Rex Drill System</li> </ul>	K163581, K150231) and Medtronic Pedicle Probe (K050438)	indications for use, design or its principles of operation and do not raise new questions of safety and efficacy.
	<ul> <li>Medtronic Navigated Anterolateral Disc Prep Instruments (K211441)</li> </ul>		

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
	<ul> <li>Medtronic Adaptix Interbody System with Titan nanoLOC Surface Technology (K201267)</li> <li>Medtronic CD HORIZO N<sup>TM</sup> Spinal System (K211596)</li> <li>Medtronic Catalyft PL (K214011)</li> <li>Medtronic Anteralign TL and Voyager FNS (K214011)</li> <li>Anteralign LS (K222383), Capstone PTC implants and Enhanced (size 36, length of 36 mm) Capstone implants</li> </ul>		
	(K172199).		

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K205005 (predicate device)	
Features	The Mazor X System consists of the	The Mazor X System consists of the	Identical
	following features:	following features:	
	• Preoperative planning and	• Preoperative planning and	
	operation	operation	
	Advanced 3D Visualization	Advanced 3D Visualization	
	(volume rendering)	(volume rendering)	
	• Robotic guidance and optical	• Mazor X Align	
	navigation of surgical tools	• Robotic guidance of surgical	
		tools	
Imaging	• CT and Fluoro based imaging	• CT and Fluoro based imaging	Identical
Modalities	• X-ray based imaging (Planning)	• X-ray based imaging (planning)	
Registration	CT-Fluoro Merge Registration	CT-Fluoro Merge Registration	Identical
Features	Automatic 3D Image	Automatic 3D Image	
	Registration (Scan and Plan)	Registration (Scan and Plan)	

Technological	Mazor	X System	Mazor	X System	Comparison
Characteristic	Mazor Robotics Ltd.		Mazor R K203005 (pr	obotics Ltd. redicate device)	
Planning	Plan Entry and Target Selection		• Plan Entry	and Target Selection	Identical
Features	• 3D Model I	Building	• 3D Model	Building	
Medical Device	• O- arm Ima	ging System	• O- arm Ima	aging System	Identical
Interfaces	• 2D C-Arm		• 2D C-Arm		
	• 3D C-Arm		• 3D C-Arm		
Dimensions	Workstation	1.85 x 1.20 x 0.64 m (6.07 x 3.94 x 2 ft)	Workstation	1.85 x 1.20 x 0.64 m (6.07 x 3.94 x 2 ft)	Identical
	Surgical System	0.92 x 0.63 x 0.58 m (3.02 x 2.07 x 1.9 ft)	Surgical System	0.92 x 0.63 x 0.58 m (3.02 x 2.07 x 1.9 ft)	
	Mazor X Navigation	$1.9 \times 0.65 \times 0.61$ m (6.2x2.1 x 2 ft)	Mazor X	$1.9 \times 0.65 \times 0.61$ m (6 2x2 1 x 2 ft)	
	Camera (NDI Vega Polaris	III (0.2x2.1 x 2 ft)	Camera (NDI Vega Polaris	III (0.2x2.1 x 2 II)	
Weight	Workstation	n - 185 kg (408 lbs)	Workstatio	n - 185 kg (408 lbs)	Identical
	Surgical Sy	stem – 25 kg (55.1	• Surgical System – 25 kg (55.1		
	lbs)		lbs)		

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K205005 (preuteate device)	
Performance	<ul> <li>Mazor X Navigation Camera (NDI Vega Polaris) - 75Kg (165.3lbs)</li> <li>Mazor X System mean accuracy</li> </ul>	<ul> <li>Mazor X Navigation Camera (NDI Vega Polaris) - 75Kg (165.3lbs)</li> <li>Mazor X System mean accuracy</li> </ul>	Similar. A requirement
	<ul> <li>Mazor A System mean accuracy &lt;1.5mm</li> <li>Navigation Accuracy mean positional error &lt;2mm and mean trajectory error of 2°.</li> <li>Mazor X will also provide Facet Decortication depth accuracy within 1.5mm.</li> </ul>	<ul> <li>Mazor A System mean accuracy &lt;1.5mm</li> <li>Navigation Accuracy mean positional error &lt;2mm and mean trajectory error of 2°.</li> </ul>	for robotic depth mean accuracy (along the trajectory) was added for the Facet Decortication procedure since in the FD procedure the robotic arm serves as a mechanical stopper. The criterion for accuracy is the same as in the cleared device <1.5 mm.

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
		R205005 (predicate device)	
			A series of performance
			bench testing
			demonstrated that the
			absolute robotic depth
			error is smaller than $\pm 1.5$
			mm and that the overall
			system accuracy is
			equivalent to the
			predicate device system
			accuracy. In addition, the
			navigation accuracy was
			tested and found to be
			equivalent to the
			navigation accuracy
			performance of the
			predicate device (mean
			positional error <2mm

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
			and mean trajectory error
			of $2^{\circ}$ ). The performance
			tests demonstrated that
			the modified Mazor X
			device keeps the same
			robotic and navigation
			accuracy as in the
			predicate device and the
			addition of accuracy
			requirement along the
			trajectory doesn't raise
			new questions of safety
			and efficacy.

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Human Factors	The Mazor X System simplifies the	The Mazor X System simplifies the	Identical
	planning and surgical procedure (pre-	planning and surgical procedure (pre-	
	intra-operative planning) and provides	intra-operative planning) and provides	
	the user with additional imaging	the user with additional imaging	
	modalities (3D Define Scan, Auto	modalities (3D Define Scan, Auto	
	Segmentation Process, etc.). Mazor X	Segmentation Process, etc.). Mazor X	
	Planning facilitates performing	Align facilitates performing	
	measurements and calculating angles	measurements and calculating angles	
	within a specific spinal region of	within a specific spinal region of	
	interest, in accordance with standard	interest, in accordance with standard	
	classifications.	classifications.	

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Cybersecurity	Mazor X system contains external	Mazor X system contains external	Identical
	wired communication interface	wired communication interface	
	(ethernet). Mazor X system is	(ethernet). Mazor X system is	
	developed in accordance with the	developed in accordance with the	
	Cybersecurity design controls to	Cybersecurity design controls to	
	ensure Mazor X Cybersecurity and	ensure Mazor X Cybersecurity and	
	maintain safety and effectiveness,	maintain safety and effectiveness,	
	similarly to the currently cleared	similarly to the currently cleared	
	device.	device.	
Standards Met	IEC 62304	IEC 62304	Identical

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Materials	• Ultem 1000HU	• Ultem 1000HU	Identical
	• Carbon Fiber Zacton 350	• Carbon Fiber Zacton 350	
	• Ketron LSG PEEK Natural	• Ketron LSG PEEK Natural	
	Black/ TECAPEEK	Black/ TECAPEEK	
	• Stainless Steel	• Stainless Steel	
	• Aluminum 6061	• Aluminum 6061	
	• Polyethylene	• Polyethylene	
	• MedikoteTM C12 (Me-DLC)	• MedikoteTM C12 (Me-DLC)	
Biocompatibi-	Materials are biocompatible	Materials are biocompatible	Identical
lity			
Compatibility	The Mazor X is compliant with the IEC	The Mazor X is compliant with the IEC	Identical
With the	60601-1-2 (EMC Compatibility)	60601-1-2 (EMC Compatibility)	
Environment	standard.	standard.	

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
and Other			
Devices			
Sterility	No new sterilization methods or	No new sterilization methods or	Identical
	products have been introduced in the	products have been introduced in the	
	Mazor X.	Mazor X.	
	The reusable accessories are sterilized	The reusable accessories are sterilized	
	by steam sterilization, single-use	by steam sterilization, single-use	
	disposable kits by gamma radiation,	disposable kits by gamma radiation,	
	and sterile covers (e.g., Sterile Sleeve	and sterile covers (e.g., Sterile Sleeve	
	for Surgical Arm) by ETO sterilization.	for Surgical Arm) by ETO sterilization.	
	Mazor X navigation tools are sterilized	Mazor X navigation tools are sterilized	
	with the Medtronic Sofamor Danek	with the Medtronic Sofamor Danek	
	spine instruments(K182121) and	spine instruments(K182121) and	
	Medtronic Navigation NavLock tools	Medtronic Navigation NavLock tools	
	(K182104).	(K182104).	

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Electrical Safety	• Power Requirements:	Power Requirements:	Extended power range to
	• 110-120 VAC, 60 Hz	• 110-120 VAC / 60 Hz	support additional
	• 220-240 VAC, 50 Hz	• 220-240 VAC / 50 Hz	markets
	• 100 VAC, 50/60 Hz	• Maximum Power = 1000 VA	
	• Maximum Power = 1000 VA		
	The Mazor X System is compliant with	The Mazor X System is compliant with	
	the IEC 60601-1 (Electrical Safety) std.	the IEC 60601-1 (Electrical Safety) std.	
Mechanical	The Mazor X System is compliant with	The Mazor X System is compliant with	Identical
Safety	the IEC 60601-1 (Electrical Safety)	the IEC 60601-1 (Electrical Safety)	
	standard.	standard.	
Chemical Safety	Not Applicable	Not Applicable	Identical
Thermal Safety	The Mazor X System is compliant with	The Mazor X System is compliant with	Identical
	the IEC 60601-1 standard.	the IEC 60601-1 standard.	

Technological	Mazor X System	Mazor X System	Comparison
Characteristic	Mazor Robotics Ltd.	Mazor Robotics Ltd.	
		K203005 (predicate device)	
Radiation	The Mazor X System is compliant with	The Mazor X System is compliant with	Identical
Safety	the IEC 60601-1-2 (EMC) standard.	the IEC 60601-1-2 (EMC) standard.	
Laser Safety	The Mazor X System is compliant with	The Mazor X System is compliant with	Identical
	IEC 60825-1 (2014) (Safety of laser)	IEC 60825-1 (2014) (Safety of laser)	
	standard.	standard.	

## Substantial Equivalence

The modified Mazor X is as safe and effective as the cleared Mazor X. The modified Mazor X has the same intended uses and indications for use and similar technological characteristics, and principles of operation as its predicate device. The minor technological differences between the modified Mazor X and its predicate device raise no new issues of safety or effectiveness. Performance data demonstrate that the modified Mazor X is as safe and effective as its predicate device. Thus, the Mazor X is substantially equivalent.

## **Non-Clinical Performance Data**

The following testing was conducted to evaluate the device:

- a) Software major level of concern verification and validation testing was conducted as required by IEC 62304 and FDA Guidance on General Principles of Software Validation, January 11, 2002.
- b) Bench testing was conducted to verify that the design outputs meet the design inputs. The functionality of the new or modified features was tested to ensure that system requirements, software requirements and user needs were met as defined. Specifically, the updated Spine Segmentation feature, the addition of Facet Decortication capability, the integration compatibility of the Medicrea UNiD (K212005) and the Planning mode UX / UI improvements were verified and validated through testing.

## Conclusions

The modified Mazor X is substantially equivalent to its predicate, the cleared Mazor X. The modified Mazor X has the same intended use and indications for use and similar technological characteristics, and principles of operation as its predicate device. The minor technological differences between the modified Mazor X and is predicate device, raise no new issues of safety or effectiveness. Non-clinical performance data demonstrate that the Mazor X performs as expected and in a manner that is substantially equivalent to its predicate device. Thus, the Mazor X is substantially equivalent.