

March 31, 2023

Shenzhen Bioeasy Biotechnology Co., Ltd. % Joe Shia
Director
LSI International
504 E Diamond Ave., Suite I
Gaithersburg, MD 20877

Re: K230238

Trade/Device Name: BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup, BIOEASY<sup>TM</sup> U-Catch MAX

Multi-Drug Test Cup Rx

Regulation Number: 21 CFR 862.3650 Regulation Name: Opiate Test System

Regulatory Class: Class II

Product Code: DJG, NFT, NFW, NFY, NGG, NGL, NFV, PTG, PTH, NGM, QAW, QBF

Dated: January 27, 2023 Received: January 30, 2023

#### Dear Joe Shia:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

K230238 - Joe Shia Page 2

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Paula V. Digitally signed by Paula V. Caposino -S
Date: 2023.03.31
20:52:54 - 04'00'

Paula Caposino, Ph.D.
Acting Deputy Director
Division of Chemistry
and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

### **Indications for Use**

510(k) Number (if known)

d-Propoxyphene (PPX)

Nortriptyline (TCA)

Marijuana (THC)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

k230238	
Device Name	
BIOEASY™ U-Catch MAX Multi-Drug Test Cup	
Indications for Use (Describe)	
BIOEASY™ U-Catch MAX Multi-Drug Test Cup tests are competi	tive binding, lateral flow immunochromatographic
assays for qualitative and simultaneous detection of Amphetamine, B	suprenorphine, Secobarbital, Oxazepam, Cocaine, 2-
ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine	e, Methylenedioxymethamphetamine, Morphine,
Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyl	ine and Marijuana in human urine at the cutoff
concentrations of:	
Drug (Identifier)	Cut-off level
Amphetamine (AMP)	500 ng/mL
Buprenorphine (BUP)	10  ng/mL
Secobarbital (BAR)	300  ng/mL
Oxazepam (BZO)	300  ng/mL
Cocaine (COC)	150 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Methamphetamine (MET)	500 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP 300)	300 ng/mL
Methadone (MTD)	300 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for OTC use.

300 ng/mL

1000 ng/mL 50 ng/mL

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Oxazepam, Secobarbital, d-Propoxyphene, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D)	x Over-The-Counter Use (21 CFR 801 Subpart C)	
CONTINUE ON A SEPARATE PAGE IF NEEDED.		

This section applies only to requirements of the Paperwork Reduction Act of 1995.

### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2023 See PRA Statement below.

### indications for US

510(k) Number (if known)		
K230238		
Device Name BIOEASY™ U-Catch MAX Multi-Drug Test Cup Rx		
ladications for the (Describe)	_	

Indications for Use (Describe)

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Oxazepam, Cocaine, 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline Cannabinoids and 6-Acetylmorphine in human urine at the cutoff concentrations of:

Drug (Identifier)	<b>Cut-off level</b>
Amphetamine (AMP)	500 ng/mL
Buprenorphine (BUP)	10 ng/mL
Secobarbital (BAR)	300 ng/mL
Oxazepam (BZO)	300  ng/mL
Cocaine (COC)	150 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Methamphetamine (MET)	500  ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP 300)	300 ng/mL
Methadone (MTD)	300  ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
d-Propoxyphene (PPX)	300 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Cannabinoids (THC)	50 ng/mL
6-Acetylmorphine	10 ng/mL

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for prescription use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Oxazepam, Secobarbital, d-Propoxyphene, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

Type of Use (Select one or both, as applicable)		
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)	
CONTINUE ON A SEPARATE PAGE IF NEEDED.		

This section applies only to requirements of the Paperwork Reduction Act of 1995.

### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

### 510(k) SUMMARY K230238

The purpose of this submission is to add analytes Amphetamine 500, Cocaine 150, Methamphetamine 500, 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP), and 6-acetylmorphine (6-AM) to previously cleared devices (k182530). These five new analytes were evaluated in this submission. For other analytes, please refer to k182530 for Buprenorphine, Secobarbital, Oxazepam, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline and Cannabinoids. In this submission the performance results are presented for the five new analytes, but the lay user study was conducted using the entire panel except of 6-AM.

1. Date: March 31, 2023

2. Submitter: Shenzhen Bioeasy Biotechnology Co., Ltd.

No.2-1 Liuxian 1st Road, Xin'an Sub-District

Shenzhen 518101, China

3. Contact person: Joe Shia

LSI International Inc

504E Diamond Ave., Suite I Gaithersburg, MD 20877 Telephone: 240-505-7880 Email: <a href="mailto:shiajl@yahoo.com">shiajl@yahoo.com</a>.

4. Device Name: BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx

Classification: Class 2

Product Code	Classification	Regulation Section	Panel
DJG	II	21 CFR § 862.3650, Morphine Test	Toxicology (91)
Monoacetylmorphine		System	
DKZ	II	21 CFR § 862.3100, Amphetamine	Toxicology (91)
Amphetamine		Test System	
LDJ	II	21 CFR § 862.3870, Cannabinoids	Toxicology (91)
Cannabinoids		Test System	
DIO	II	21 CFR § 862.3250, Cocaine and	Toxicology (91)
Cocaine		Cocaine Metabolites Test System	
LAF	II	21 CFR § 862.3610,	Toxicology (91)
Methamphetamine		Methamphetamine Test System	
DJG	II	21 CFR § 862.3650, Morphine	Toxicology (91)
Morphine		Test System	
JXM	II	21 CFR § 862.3170,	Toxicology (91)
Oxazepam		Benzodiazepine Test System	
DJG	II	21 CFR § 862.3650, Opiate Test	Toxicology (91)
Oxycodone		System	
DIS	II	21 CFR § 862.3150, Barbiturate	Toxicology (91)
Secobarbital		Test System	
DJG	II	21 CFR § 862.3650,	Toxicology (91)
Buprenorphine		Opiate Test System	

LAF Methylenedioxy- methamphetamine	II	21 CFR § 862.3610, Methamphetamine Test System	Toxicology (91)
LCM Phencyclidine	unclassified	Enzyme Immunoassay Phencyclidine	Toxicology (91)
DJR Methadone	II	21 CFR § 862.3620, Methadone Test System	Toxicology (91)
DJR 2-ethylidene-1, 5- dimethyl-3, 3- diphenylpyrrolidine (EDDP)	II	21 CFR § 862.3620, Methadone Test System	Toxicology (91)
LFG Nortriptyline	II	21 CFR, 862.3910 Tricyclic Antidepressant Drugs Test System	Toxicology (91)
JXN Propoxyphene	II	21 CFR, 862.3700 Propoxyphene Test System	Toxicology (91)

### 5. Predicate Devices: K201630

The Assure Tech Panel Dip Tests/AssureTech Quick Cup Tests

#### 6. Intended Use

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Oxazepam, Cocaine, 2- ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline and Marijuana in human urine at the cutoff concentrations of:

Drug (Identifier) Amphetamine (AMP) Buprenorphine (BUP)	Cut-off level 500 ng/mL 10 ng/mL
Secobarbital (BAR)	300 ng/mL
Oxazepam (BZO)	300 ng/mL
Cocaine (COC) 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	150 ng/mL 300 ng/mL
Methamphetamine (MET)	500 ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP 300)	300 ng/mL
Methadone (MTD)	300 ng/mL
Oxycodone (OXY)	100 ng/mL
Phencyclidine (PCP)	25 ng/mL
d-Propoxyphene (PPX)	300 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Marijuana (THC)	50 ng/mL

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for OTC use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Oxazepam, Secobarbital, d-Propoxyphene, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx tests are competitive binding, lateral flow immunochromatographic assays for qualitative and simultaneous detection of Amphetamine, Buprenorphine, Secobarbital, Oxazepam, Cocaine, 2- ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, Methamphetamine, Methylenedioxymethamphetamine, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline Cannabinoids and 6-Acetylmorphine in human urine at the cutoff concentrations of:

Teorice in an one of .	
Drug (Identifier)	Cut-off level
Amphetamine (AMP)	500  ng/mL
Buprenorphine (BUP)	10 ng/mL
Secobarbital (BAR)	300  ng/mL
Oxazepam (BZO)	300 ng/mL
Cocaine (COC)	150 ng/mL
2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300 ng/mL
Methamphetamine (MET)	500  ng/mL
Methylenedioxymethamphetamine (MDMA)	500 ng/mL
Morphine (MOP 300)	300 ng/mL
Methadone (MTD)	300 ng/mL
Oxycodone (OXY)	100  ng/mL
Phencyclidine (PCP)	25 ng/mL
d-Propoxyphene (PPX)	300 ng/mL
Nortriptyline (TCA)	1000 ng/mL
Cannabinoids (THC)	50 ng/mL
6-Acetylmorphine	10 ng/mL
	_

BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx offers any combinations of the above listed analytes. It is for in vitro diagnostic use only. It is intended for prescription use.

The tests may yield positive results for the prescription drugs Buprenorphine, Nortriptyline, Oxazepam, Secobarbital, d-Propoxyphene, and Oxycodone when taken at or above prescribed doses. It is not intended to distinguish between prescription use or abuse of these drugs. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result.

The tests provide only preliminary results. To obtain a confirmed analytical result, a more specific alternate chemical method must be used. GC/MS or LC/MS is the recommended confirmatory method.

### 7. Device Description

The BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup and BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup Rx are immunochromatographic assays that use a lateral flow system for the qualitative detection of target drug or drug metabolites in human urine. The products are single-use in vitro diagnostic devices. The BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup kit contains a Cup device, a package insert and a urine cup for sample collection. Each test device is sealed with a desiccant in an aluminum pouch.

### 8. Substantial Equivalence Information

A summary comparison of features of the BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup and the predicate devices is provided in following tables.

Table 1: Features Comparison of BIOEASY™ U-Catch MAX Multi-Drug Test Cup and the Predicate Devices

Item	Device	Predicate – K201630
Indication(s)	For the qualitative determination of drugs of	Same (but the number of

for Use	abuse in human urine.	drugs detected is different)	
	Amphetamine (AMP): 500 ng/ml		
	Oxazepam (BZO):300 ng/ml		
	Cocaine (COC): 150 ng/ml		
	11-Nor-△ <sup>9</sup> -Tetrahydrocannabinol-9-COOH (THC):50 ng/ml		
	Methamphetamine (MET): 500 ng/ml		
	Morphine (MOR): 300 ng/mL		
	Oxycodone(OXY): 100 ng/ml		
Calibrator and	Secobarbital (BAR): 300 ng/ml	Same except	
<b>Cut-Off Values</b>	Methadone (MTD): 300 ng/ml	THC at 20 ng/mL	
	Buprenorphine (BUP): 10 ng/ml		
	D,L-Methylenedioxymethamphetamine (MDMA): 500 ng/ml		
	Phencyclidine (PCP): 25 ng/ml		
	Nortriptyline (TCA): 1000 ng/ml		
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP): 300 ng/ml		
	d-Propoxyphene (PPX): 300 ng/ml		
	Competitive binding, lateral flow		
Methodology	immunochromatographic assays based on the	Same	
Wiethodology	principle of antigen antibody	Same	
	immunochemistry.		
Type of Test	Qualitative	Same	
Specimen Type	Human Urine	Same	
<b>Intended Use</b>	For over-the-counter	For prescription use	
Configurations	Cup	Same	

Table 2: Features Comparison of BIOEASY $^{\text{TM}}$  U-Catch MAX Multi-Drug Test Cup Rx and the Predicate Devices

Item	Device	Predicate – K201630
Indication(s) for Use	For the qualitative determination of drugs of abuse in human urine.	Same
Calibrator and Cut-Off Values	Amphetamine (AMP): 500 ng/ml Oxazepam (BZO):300 ng/ml Cocaine(COC): 150 ng/ml 11-Nor-Δ <sup>9</sup> -Tetrahydrocannabinol-9-COOH (THC):50 ng/ml Methamphetamine (MET): 500 ng/ml Morphine (MOR): 300ng/mL Oxycodone(OXY): 100 ng/ml Secobarbital (BAR): 300 ng/ml	Same except THC at 20 ng/mL
	Methadone (MTD): 300 ng/ml	

	Buprenorphine (BUP): 10 ng/ml	
	D,L-Methylenedioxymethamphetamine (MDMA): 500 ng/ml	
	Phencyclidine (PCP): 25 ng/ml	
	Nortriptyline (TCA): 1000 ng/ml	
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP): 300 ng/ml	
	d-Propoxyphene (PPX): 300 ng/ml	
	6-Acetylmorphine (6-AM): 10 ng/mL	
Methodology	Competitive binding, lateral flow immunochromatographic assays based on the principle of antigen antibody immunochemistry.	Same
Type of Test	Qualitative	Same
Specimen Type	Human Urine	Same
Intended Use	For prescription use	Same
Configurations	Cup	Cup

### 9. Test Principle

The BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup tests are rapid tests for the qualitative detection of target drug or drug metabolites in urine samples. The tests are lateral flow chromatographic immunoassays. During testing, a urine specimen migrates upward by capillary action. If target drugs present in the urine specimen are below the cut-off concentration, it will not saturate the binding sites of its specific monoclonal mouse antibody coated on the particles. The antibody-coated particles will then be captured by immobilized drug-conjugate and a visible colored line will show up in the test line region. The colored line will not form in the test line region if the target drug level exceeds its cutoff-concentration because it will saturate all the binding sites of the antibody coated on the particles. A band should form in the control region of the devices regardless of the presence of drug or metabolite in the sample to indicate that the tests have been performed properly.

### 10. Performance Characteristics

### 1. Analytical Performance

### a. Precision

Precision studies were carried out for samples with concentrations of -100% cut off, -75% cut off, -50% cut off, -25% cut off, cut off, +25% cut off, +50% cut off, +75% cut off and +100% cut off. These samples were prepared by spiking drug in negative urine samples. Each drug concentration was confirmed by LC/MS. All sample aliquots were blindly labeled by the person who prepared the samples and didn't take part in the sample testing. For each concentration, tests were performed two runs per day for 25 days per device in a randomized order. The results obtained are summarized in the following tables for Amphetamine 500, Cocaine 150, Methamphetamine 500, 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP), and 6-acetylmorphine (6-AM). The data for Buprenorphine, Methylenedioxymethamphetamine, Secobarbital, Oxazepam, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene,

# Nortriptyline and Cannabinoids were reported in k182530.

# **AMP500**

Concentration by LC/MS (ng/mL) Lot	-100% Cut-off	-75% Cut-off	-50% Cut-off	-25% Cut-off	Cut-off	Cut-off +25%	Cut-off +50%	Cut-off +75%	Cut-off +100%
Number	0	127	250	374	480	590	695	825	915
Lot 1	50-/0+	50-/0+	50-/0+	50-/0+	22-/28+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 2	50-/0+	50-/0+	50-/0+	50-/0+	26-/24+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	26-/24+	50+/0-	50+/0-	50+/0-	50+/0-

## **COC150**

Concentration by LC/ MS (ng/mL) Lot	-100% Cut-off	-75% Cut-off	-50% Cut-off	-25% Cut-off	Cut-off	Cut-off +25%	Cut-off +50%	Cut-off +75%	Cut-off +100%
Number	0	38.2	77.8	113	154	191	230	262	302
Lot 1	50-/0+	50-/0+	50-/0+	50-/0+	28-/22+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 2	50-/0+	50-/0+	50-/0+	50-/0+	29-/21+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	22-/28+	50+/0-	50+/0-	50+/0-	50+/0-

# **MET500**

Concentration by LC/MS (ng/mL)	-100% Cut-off	-75% Cut-off	-50% Cut-off	-25% Cut-off	Cut-off	+25% Cut-off	+50% Cut-off	+75% Cut-off	+100% Cut-off
Number	0	129	240	368	476	605	705	860	1010
Lot 1	50-/0+	50-/0+	50-/0+	50-/0+	27-/23+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 2	50-/0+	50-/0+	50-/0+	50-/0+	23-/27+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	25-/25+	50+/0-	50+/0-	50+/0-	50+/0-

# **EDDP**

Concentration by LC/MS (ng/mL)	-100% Cut-off	-75% Cut-off	-50% Cut-off	-25% Cut-off	Cut-off	+25% Cut-off	+50% Cut-off	+75% Cut-off	+100% Cut-off
Number	0	72.0	146	216	293	360	426	507	582
Lot 1	50-/0+	50-/0+	50-/0+	50-/0+	24-/26+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 2	50-/0+	50-/0+	50-/0+	50-/0+	26-/24+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	25-/25+	50+/0-	50+/0-	50+/0-	50+/0-

# 6-AM

Concentration by LC/MS (ng/mL)	-100% Cut-off	-75% Cut-off	-50% Cut-off	-25% Cut-off	Cut-off	+25% Cut-off	+50% Cut-off	+75% Cut-off	+100% Cut-off
Lot Number	0	2.46	4.88	7.42	9.98	12.3	14.6	16.8	20.1
Lot 1	50-/0+	50-/0+	50-/0+	50-/0+	26-/24+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 2	50-/0+	50-/0+	50-/0+	50-/0+	24-/26+	50+/0-	50+/0-	50+/0-	50+/0-
Lot 3	50-/0+	50-/0+	50-/0+	50-/0+	23-/27+	50+/0-	50+/0-	50+/0-	50+/0-

The following cut-off values are verified.

Drug (Identifier)	Cut-off level

Amphetamine (AMP)	500 ng/mL
Cocaine (COC)	150 ng/mL
Methamphetamine (MET)	500 ng/mL
2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP)	300 ng/mL
6-acetylmorphine (6-AM)	10 ng/mL

### b. Linearity

Not applicable.

### c. Stability and Traceability

The devices are stable at 4-30 °C for 24 months based on real time stability studies. All drug calibrators of the device are traceable to available commercial reference materials.

### d. Interference

Potential interfering substances found in human urine of physiological or pathological conditions were added to drug-free urine and target drugs urine with concentrations at 25% below and 25% above Cut-Off levels. These urine samples were tested using three lots of each device. Compounds that showed no interference at a concentration of  $100\mu g/mL$  are summarized in the following tables.

Acetaminophen Creatinine Ketamine Prednisone Acetophenetidin Deoxycorticosterone Ketoprofen (±)-Propranolol N-Acetylprocainamide Dextromethorphan Labetalol Pseudoephedrine Acetylsalicylic acid Diclofenac Loperamide Ouinine Albumin (100mg/dL) Diflunisal Meperidine Ranitidine Aminopyrine Digoxin Meprobamate Salicylic acid Serotonin (5- Hydroxytyramine) Amoxicillin Diphenhydramine Methoxyphenamine Nalidixic acid Ampicillin 1% Ethanol Sulfamethazine Apomorphine Ecgonine methyl ester Naloxone Sulindac Ascorbic acid β-Estradiol Naltrexone Tetrahydrocortisone 3-(β-Dglucuronide) Tetrahydrocortisone 3-acetate Aspartame Erythromycin Naproxen Atropine Fenoprofen Niacinamide Tetrahvdrozoline Benzilic acid Furosemide Thiamine Nifedipine Benzoic acid Gentisic acid Norethindrone Thioridazine Bilirubin Hemoglobin Noscapine Triamterene Chloral hydrate Hydralazine (±)-Octopamine Trifluoperazine Chloramphenicol Hydrochlorothiazide Oxalic acid Trimethoprim Chlorothiazide Hydrocortisone Oxolinic acid DL-Tryptophan Tyramine Chlorpromazine O-Hydroxyhippuric acid Oxymetazoline Cholesterol 3-Hydroxytyramine Papaverine DL-Tyrosine Clonidine Ibuprofen Penicillin G Uric acid Cortisone Isoproterenol Perphenazine Verapamil Phenelzine Zomepirac (-)-Cotinine Isoxsuprine

### e. Specificity

To test specificity, drug metabolites and other structurally related compounds that are likely to cross-react in urine samples were spiked into negative urine and were tested using three lots of each device. The lowest concentration that caused a positive result for each compound are listed below for Amphetamine 500, Cocaine 150,

Methamphetamine 500, 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP), and 6-acetylmorphine (6-AM). The data for Buprenorphine,

Methylenedioxymethamphetamine, Secobarbital, Oxazepam, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline and Cannabinoids were reported in k182530.

AMP500	Result	%Cross-Reactivity
(Cut-off=500 ng/mL)	Positive at (ng/ml)	
D - Amphetamine	500	100%
L - Amphetamine	10000	5%
DL - Amphetamine	1500	33%
Phentermine	15000	3.3%
Hydroxyamphetamine	4000	12.5%
Methylenedioxyamphetamine (MDA)	10000	5%
d-Methamphetamine	> 100000	< 0.5%
1-Methamphetamine	> 100000	<0.5%
Ephedrine	> 100000	<0.5%
Methylenedioxyethylamphetamine	> 100000	< 0.5%
(MDE)		
3,4-methylenedioxy-methamphetamine	> 100000	< 0.5%
(MDMA)		

COC150 (Cut-off=150 ng/mL)	Result Positive at (ng/ml)	%Cross-Reactivity
Benzoylecgonine	150	100%
Cocaine HCl	375	40%
Cocaethylene	6250	2.4%
Ecgonine	16000	0.9%
Norcocaine	50000	0.3%

MET500	Result	%Cross-Reactivity
(Cut-off=500 ng/mL)	Positive at (ng/ml)	
D(+)-Methamphetamine	500	100%
(+/-)3,4-Methylenedioxy-n-ethylamphetamine (MDEA)	5000	10%
D/L-Methamphetamine	500	100%
p-Hydroxymethamphetamine	5000	10%
D-Amphetamine	> 100000	< 0.5%
L-Amphetamine	>100000	< 0.5%
Chloroquine	25000	2%
(+/-)-Ephedrine	2000	25%
L-Methamphetamine	5000	10%
(+/-)3,4-Methylenedioxyamphetamine (MDA)	> 100000	<0.5%
β-Phenylethylamine	3750	13.3%
Trimethobenzamide	10000	5%
(+/-)3,4- methylenedioxymethamphetamine(MDMA)	5000	10%

EDDP	Result	% Cross-Reactivity
(Cut-off=300 ng/mL)	Positive at (ng/ml)	•

EDDP	300	100%
Methadone	300000	0.1%
Doxylamine	>100000	<0.3%
LAAM HCI	>100000	<0.3%
Alpha Methadol	>100000	<0.3%
EMDP	>100000	<0.3%
Disopyramide	>100000	<0.3%

6-AM (Cut-off=10 ng/mL)	Result Positive at (ng/ml)	% Cross-Reactivity
6-acetylmorphine	10	100%
Acetylcodeine	>10000	<0.1%
Buprenorphine	>10000	<0.1%
Codeine	>10000	<0.1%
Diacetylmorphine	1000	0.01%
Dihydrocodeine	>10000	<0.1%
Ethylmorphine	>10000	<0.1%
Hydrocodone	>10000	<0.1%
Hydromorphone	5000	0.002%
Morphine	10000	0.001%
Morphine-3-glucuronide	>10000	<0.1%
Nalorphine	5000	0.002%
Thebaine	>20000	<0.05%
Dextromethorphan	>100,000	<0.01%
Heroin	100000	0.0001%
Imipramine	>100,000	<0.01%
LAAM (Levacetylmethadol)	>100,000	<0.01%
Levorphanol	>100,000	<0.01%
Meperidine	>100,000	<0.01%
Methadone	>100,000	<0.01%
Mitragynine (kratom)	>20,000	<0.05%
Morphine 6-D-glucuronide	>100,000	<0.01%
Naloxone	>100,000	<0.01%
Naltrexone	>100,000	<0.01%
Naproxen	>100,000	<0.01%
Norbuprenorphine	>10,000	<0.1%
Norbuprenorphine glucuronide	>100,000	<0.01%
Norcodeine	>100,000	<0.01%
Norhydrocodone	>100,000	<0.01%
Normorphine	>100,000	<0.01%
Noroxycodone	>100,000	<0.01%
Noroxymorphone	>100,000	<0.01%
Norpropoxyphene	>100,000	<0.01%
Oxycodone	>100,000	<0.01%
Oxymorphone	>100,000	<0.01%
Oxymorphone-3β-D-glucuronide	>100,000	<0.01%
Tapentadol HCl	>100,000	<0.01%

Tramadol	>100.000	<0.01%
Trainador	100,000	-0.0170

### f. Effect of Urine Specific Gravity and Urine pH

To investigate the effect of urine specific gravity and urine pH, urine samples, with 1.000 to 1.035 specific gravity or urine samples with pH 4 to 9 were spiked with target drugs at 25% below and 25% above Cut-Off levels. These samples were tested using three lots of each device. Results were all positive for samples at and above +25% Cut-Off and all negative for samples at and below -25% Cut-Off.

### 2. Comparison Studies

Method comparison studies for the BIOEASY<sup>TM</sup> U-Catch MAX Multi-Drug Test Cup were performed in-house with three laboratory assistants. Operators ran 80 (40 negative and 40 positive) unaltered clinical samples for each drug. The samples were blind labeled and compared to LC/MS results. The results are presented in the tables below for Amphetamine 500, Cocaine 150, Methamphetamine 500, 2-ethylidene-1, 5-dimethyl-3, 3-diphenylpyrrolidine (EDDP), and 6-acetylmorphine (6-AM). The data for Buprenorphine, Methylenedioxymethamphetamine, Secobarbital, Oxazepam, Morphine, Methadone, Oxycodone, Phencyclidine, d-Propoxyphene, Nortriptyline and Cannabinoids were reported in k182530.

#### **AMP500**

AMI SUU						
U-Catch			Low	Near Cutoff	Near Cutoff	
MAX Multi-		Negative	Negative by	Negative by	Positive by	High Positive
Drug Test			LC/MS	LC/MS	LC/MS	by LC/MS
Cup			(less than	(Between	(Between the	(greater than
			-50%)	-50% and	cutoff and	+50%)
				cutoff)	+50%)	
Viewer	Positive	0	0	3	17	20
A	Negative	7	16	14	3	0
Viewer	Positive	0	0	2	18	20
В	Negative	7	16	15	2	0
Viewer	Positive	0	0	3	18	20
C	Negative	7	16	14	2	0

### **Discordant Results**

Viewer	Sample Number	LC/MS Result	Dip Card Viewer Results
Viewer A	AMPLC063	473	Positive
Viewer C	AMPLC063	473	Positive
Viewer B	AMPLC040	477	Positive
Viewer A	AMPLC039	490	Positive
Viewer B	AMPLC039	490	Positive
Viewer C	AMPLC039	490	Positive
Viewer A	AMPLC004	494	Positive
Viewer C	AMPLC004	494	Positive
Viewer A	AMPLC018	520	Negative
Viewer B	AMPLC018	520	Negative
Viewer A	AMPLC064	525	Negative
Viewer B	AMPLC064	525	Negative
Viewer C	AMPLC064	525	Negative

Viewer A	AMPLC009	540	Negative
Viewer C	AMPLC014	570	Negative

## **COC150**

COC150						
U-Catch			Low	Near Cutoff	Near Cutoff	
MAX Multi-		Negative	Negative by	Negative by	Positive by	High Positive
Drug Test			LC/MS	LC/MS	LC/MS	by LC/MS
Cup			(less than	(Between	(Between the	(greater than
			-50%)	-50% and	cutoff and	+50%)
				cutoff)	+50%)	
Viewer	Positive	0	0	2	19	19
A	Negative	7	16	15	2	0
Viewer	Positive	0	0	2	18	19
В	Negative	7	16	15	3	0
Viewer	Positive	0	0	2	19	19
C	Negative	7	16	15	2	0

## **Discordant Results**

Viewer	Sample Number	LC/MS Result	Easy Cup Viewer Results
Viewer B	COCLC046	139	Positive
Viewer A	COCLC019	142	Positive
Viewer C	COCLC019	142	Positive
Viewer B	COCLC026	142	Positive
Viewer A	COCLC060	145	Positive
Viewer C	COCLC060	145	Positive
Viewer A	COCLC049	154	Negative
Viewer B	COCLC049	154	Negative
Viewer C	COCLC049	154	Negative
Viewer A	COCLC024	157	Negative
Viewer B	COCLC062	157	Negative
Viewer C	COCLC062	157	Negative
Viewer B	COCLC029	159	Negative

### **MET500**

U-Catch			Low	Near Cutoff	Near Cutoff	
MAX Multi-		Negative	Negative by	Negative by	Positive by	High Positive
Drug Test			LC/MS	LC/MS	LC/MS	by LC/MS
Cup			(less than	(Between	(Between the	(greater than
			-50%)	-50% and	cutoff and	+50%)
				cutoff)	+50%)	
Viewer	Positive	0	0	2	17	20
A	Negative	7	16	15	3	0
Viewer	Positive	0	0	3	18	20
В	Negative	7	16	14	2	0
Viewer	Positive	0	0	2	17	20
C	Negative	7	16	15	3	0

**Discordant Results** 

Viewer	Sample Number	LC/MS Result	Dip Card Viewer Results
Viewer B	METL055	427	Positive
Viewer B	METL005	447	Positive
Viewer C	METL005	447	Positive
Viewer A	METL020	465	Positive
Viewer B	METL020	465	Positive
Viewer A	METL063	486	Positive
Viewer C	METL063	486	Positive
Viewer B	METL053	525	Negative
Viewer C	METL053	525	Negative
Viewer A	METL060	530	Negative
Viewer C	METL060	530	Negative
Viewer A	METL036	540	Negative
Viewer B	METL036	540	Negative
Viewer A	METL059	555 Negati	
Viewer C	METL059	555	Negative

# **EDDP**

U-Catch			Low	Near Cutoff	Near Cutoff	
MAX Multi-		Negative	Negative by	Negative by	Positive by	High Positive
Drug Test			LC/MS	LC/MS	LC/MS	by LC/MS
Cup			(less than	(Between	(Between the	(greater than
			-50%)	-50% and	cutoff and	+50%)
				cutoff)	+50%)	
Viewer	Positive	0	0	2	18	20
A	Negative	6	18	14	2	0
Viewer	Positive	0	0	2	18	20
В	Negative	6	18	14	2	0
Viewer	Positive	0	0	3	18	20
C	Negative	6	18	13	2	0

# **Discordant Results**

Viewer	Sample Number	LC/MS Result	Easy Cup Viewer Results
Viewer A	EDDPLC042	244	Positive
Viewer C	EDDPLC042	244	Positive
Viewer B	EDDPLC004	291	Positive
Viewer C	EDDPLC004	291	Positive
Viewer A	EDDPLC051	294	Positive
Viewer B	EDDPLC051	294	Positive
Viewer C	EDDPLC051	294	Positive
Viewer A	EDDPLC074	303	Negative
Viewer B	EDDPLC074	303	Negative
Viewer C	EDDPLC074	303	Negative
Viewer B	EDDPLC052	327	Negative
Viewer A	EDDPLC018	330	Negative

Viewer C	EDDPLC018	330	Negative
----------	-----------	-----	----------

### **6-AM**

U-Catch			Low	Near Cutoff	Near Cutoff	
MAX Multi-		Negative	Negative by	Negative by	Positive by	High Positive
Drug Test			LC/MS	LC/MS	LC/MS	by LC/MS
Cup			(less than	(Between	(Between the	(greater than
			-50%)	-50% and	cutoff and	+50%)
				cutoff)	+50%)	
Viewer	Positive	0	0	2	16	22
A	Negative	6	17	15	2	0
Viewer	Positive	0	0	2	16	22
В	Negative	6	17	15	2	0
Viewer	Positive	0	0	2	15	22
C	Negative	6	17	15	3	0

### **Discordant Results**

Viewer	Sample Number	LC/MS Result	Dip Card Viewer Results
Viewer B	6-AMLC016	9.23	Positive
Viewer A	6-AMLC074	9.68	Positive
Viewer C	6-AMLC074	9.68	Positive
Viewer B	6-AMLC026	9.93	Positive
Viewer A	6-AMLC037	9.95	Positive
Viewer C	6-AMLC037	9.95	Positive
Viewer A	6-AMLC009	10.3	Negative
Viewer C	6-AMLC009	10.3	Negative
Viewer B	6-AMLC027	11.1	Negative
Viewer C	6-AMLC027	11.1	Negative
Viewer A	6-AMLC038	11.4	Negative
Viewer B	6-AMLC038	11.4	Negative
Viewer C	6-AMLC035	11.8	Negative

### Lay-user study

A lay user study was performed at three intended user sites with 140 lay persons. The lay users had diverse educational and professional backgrounds and ranged in age from 20 to > 50 years. Urine samples were prepared at the following concentrations; negative, +/-75%, +/-50%, +/-25% of the cutoff by spiking drugs into drug free-pooled urine specimens. The concentrations of the samples were confirmed by LC/MS. Each sample was aliquoted into individual containers and blind-labeled. Each participant was provided with the package insert, 1 blind labeled sample and a device. Each device was tested. Results are shown below.

Drugs	% of Cut-off	Number of samples	Concentration by LC/MS (ng/mL)	Lay perso	The percentage	
				No. of Positive	No. of Negative	of correct results (%)
AMP	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	130	0	20	100

	1				1	
	-50% Cut-off	20	251	0	20	100
	-25% Cut-off	20	383	0	20	100
	+25% Cut-off	20	635	19	1	95
	+50% Cut-off	20	755	20	0	100
	+75% Cut-off	20	885	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	75.9	0	20	100
BAR	-50% Cut-off	20	150	0	20	100
D/ IIC	-25% Cut-off	20	220	0	20	100
	+25% Cut-off	20	360	19	1	95
	+50% Cut-off	20	429	20	0	100
	+75% Cut-off	20	501	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	37.6	0	20	100
COC	-50% Cut-off	20	74.8	0	20	100
	-25% Cut-off	20	110	1	19	95
	+25% Cut-off	20	183	19	1	95
	+50% Cut-off	20	224	20	0	100
	+75% Cut-off	20	248	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	70.8	0	20	100
BZO	-50% Cut-off	20	148	0	20	100
BEC	-25% Cut-off	20	224	1	19	95
	+25% Cut-off	20	390	19	1	95
	+50% Cut-off	20	452	20	0	100
	+75% Cut-off	20	504	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	124	0	20	100
MET	-50% Cut-off	20	242	0	20	100
IVIL	-25% Cut-off	20	367	0	20	100
	+25% Cut-off	20	610	19	1	95
	+50% Cut-off	20	705	20	0	100
	+75% Cut-off	20	825	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	76.8	0	20	100
MTD	-50% Cut-off	20	147	0	20	100
1,111	-25% Cut-off	20	226	0	20	100
	+25% Cut-off	20	375	19	1	95
	+50% Cut-off	20	441	20	0	100
	+75% Cut-off	20	504	20	0	100
	-100% Cut-off	20	0	0	20	100
МОР	-75% Cut-off	20	79	0	20	100
	-50% Cut-off	20	158	0	20	100
	-25% Cut-off	20	246	1	19	95
	+25% Cut-off	20	389	19	1	95
	+50% Cut-off	20	469	20	0	100
	+75% Cut-off	20	530	20	0	100
OXY	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	24.5	0	20	100

	-50% Cut-off	20	49.3	0	20	100
	-30% Cut-off	20	71.1	0	20	100
	+25% Cut-off	20	118	19	1	95
	+50% Cut-off	20	147	20	0	100
	+75% Cut-off	20	169	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	13	0	20	100
	-50% Cut-off	20	25.3	0	20	100
THC	-25% Cut-off	20	41	1	19	95
	+25% Cut-off	20	65	19	1	95
	+50% Cut-off	20	79	20	0	100
	+75% Cut-off	20	93	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	273	0	20	100
	-50% Cut-off	20	509	0	20	100
TCA	-25% Cut-off	20	809	0	20	100
	+25% Cut-off	20	1190	19	1	95
	+50% Cut-off	20	1510	20	0	100
	+75% Cut-off	20	1680	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	2.57	0	20	100
	-50% Cut-off	20	5.14	0	20	100
BUP	-25% Cut-off	20	6.76	1	19	95
	+25% Cut-off	20	12.8	19	1	95
	+50% Cut-off	20	15.1	20	0	100
	+75% Cut-off	20	17.2	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	6.27	0	20	100
	-50% Cut-off	20	12.5	0	20	100
PCP	-25% Cut-off	20	17.9	1	19	95
	+25% Cut-off	20	30.8	19	1	95
	+50% Cut-off	20	36.4	20	0	100
	+75% Cut-off	20	42.8	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	137	0	20	100
	-50% Cut-off	20	250	0	20	100
MDMA	-25% Cut-off	20	351	0	20	100
	+25% Cut-off	20	600	19	1	95
	+50% Cut-off	20	745	20	0	100
	+75% Cut-off	20	925	20	0	100
	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	76.2	0	20	100
	-50% Cut-off	20	141	0	20	100
EDDP	-25% Cut-off	20	224	1	19	95
	+25% Cut-off	20	375	20	0	100
	+50% Cut-off	20	447	20	0	100
	1 30 /0 Cut-011	20	77/	20		100

	+75% Cut-off	20	522	20	0	100
PPX	-100% Cut-off	20	0	0	20	100
	-75% Cut-off	20	77.4	0	20	100
	-50% Cut-off	20	150	0	20	100
	-25% Cut-off	20	227	1	19	95
	+25% Cut-off	20	351	20	0	100
	+50% Cut-off	20	420	20	0	100
	+75% Cut-off	20	492	20	0	100

Lay-users were also given surveys on the ease of understanding the package insert instructions. All lay users indicated that the device instructions can be easily followed. A Flesch-Kincaid reading analysis was performed on each package insert and the scores revealed a reading Grade Level of 7.

### 3. Clinical Studies

Not applicable.

### 11. Conclusion

Based on the test principle and acceptable performance characteristics including precision, cut-off, interference, specificity, method comparison, and lay-user studies of the devices, it's concluded that the U-Catch MAX Multi-Drug Test Cup is substantially equivalent to the predicate.