

Orthofix US LLC % Shant Aghyarian Regulatory Affairs Program Manager 3451 Plano Parkway LEWISVILLE TX 75056

September 26, 2023

Re: K230252

Trade/Device Name: OFIX MIS App Regulation Number: 21 CFR 892.2050

Regulation Name: Medical image management and processing system

Regulatory Class: Class II Product Code: LLZ

Dated: January 31, 2023 Received: August 30, 2023

Dear Shant Aghyarian:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/efdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/training-and-continuing-education/cdrh-learn) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jessica Lamb, Ph.D.

Assistant Director Imaging Software Team

DHT8B: Division of Radiological Imaging

Devices and Electronic Products
OHT8: Office of Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023
See PRA Statement below.

Submission Number (if known)				
K230252				
Device Name				
OFIX MIS App				
Indications for Use (Describe)				
The OFIX MIS App software assists healthcare professionals in viewing, storing, and measuring images as well as planning orthopedic surgeries. The app allows service providers to plan surgical procedures by making measurements for the placement of surgical implants. Clinical judgement and experience are required to properly use the software.				
Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)				

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) SUMMARY

OFIX MIS App

K230252

510(k) Owner Information

Name: Orthofix US LLC Address: 3451 Plano Parkway

Lewisville, TX, USA

Telephone Number: 214-937-2176 Fax Number: 214-937-3322

Email: shantaghyarian@orthofix.com

Registration Number: 2183449

Contact Person: Shant Aghyarian

Date Prepared: September 21, 2023

Name of Device

Trade Name / Proprietary

Name: OFIX MIS App

Common Name: OFIX MIS App

Product Code(s): LLZ

Classification Name: System, Image Processing, Radiological

Regulatory Class: Class II per 21 CFR 892.2050

Predicate Devices: Nemaris Surgimap 2.0, K141669



Device Description OFIX MIS App

The OFIX MIS App will provide a novel mobile phone solution to supplement the use of calipers when selecting rods used in immobilization and stabilization of pedicle screw type spinal systems. The App will rely on capturing 2D images used in spinal implant procedures through use of mobile phone camera. The App measures the difference between known Orthofix pedicle screw diameters and provides feedback useful in determining appropriate rod length with corresponding Orthofix product number.

Indications for Use

The OFIX MIS App software assists healthcare professionals in viewing, storing, and measuring images as well as planning orthopedic surgeries. The App allows service providers to plan surgical procedures by including tools for measuring anatomical components for placement of surgical implants. Clinical judgment and experience are required to properly use the software.

Technological Comparison

The technological comparison is provided as the Substantial Equivalence table below.



Substantial Equivalence Table

G 15 (1 B)	2 2 4 2 4	211 (7)	D 1 1
Specification/Property	Predicate Device Nemaris Surgimap 2.0 (K141669)	Subject Device	Discussion
	Nemaris Surgimap 2.0 (K141009)	OFIX MIS App	
	The Surgimap software assists healthcare professionals	The OFIX MIS App software assists	The OFIX App doesn't measure anatomical
	in viewing, storing, and measuring images as well as	healthcare professionals in viewing,	components but helps measure distances on
	planning orthopedic surgeries. The device allows service	storing, and measuring images as well	anatomical components between placed markers
	providers to perform generic as well as specialty	as planning orthopedic surgeries. The	(pedicle screws).
	measurements of the images, and to plan surgical	app allows service providers to plan	
	procedures. The device also includes tools for measuring	surgical procedures by making	The predicate is more involved in overall pre-
	anatomical components for placement of surgical	measurements for the placement of	operative planning. This could be determinant of
	implants and offer online synchronization of the database	surgical implants. Clinical judgement	the surgical procedure or approach taken by the
Intended Use / Indications for	with the possibility to share data among Surgimap users.	and experience are required to	surgeon. The OFIX MIS App, is a lower risk tool
Use	Clinical judgment and experience are required to	properly use the software.	used to streamline the process of measuring and
	properly use the software.		choosing the correct rod to be used in the surgical
			procedure. The surgical procedure and approach
			would already be decided at this point. As a
			safeguard, the surgeon has access to traditional means of measurement such as calipers. This is a
			check in case the surgeon wants to confirm the
			OFIX MIS App's measurement.
			Of IA WIIS App sineastrement.
Device Classification Name	Medical image management and processing system	Medical image management and	
Device classification realife	2.20 and and a management and processing system	processing system	Same
Software Functionalities/		processing system	The OFIX MIS App is a simplified software that is
Modalities	Generic, Spine, and Lower Limbs Measurements, Pre-op	Templating (vendor specific	for templating Orthofix specific implants. It is
	Planning, Templating (vendor specific implants and	implants), Database	comparable to the predicate in that manner without
	custom implants), Database, Case Sharing		the additional functionalities that the predicate
			offers. The predicate has previously been used by
			medical professionals to template Orthofix
			implants. After measuring, the OFIX MIS App
			recommends a rod length pulled from its database.
			This doesn't not affect the safety and effectiveness
			of the device and the OFIX MIS App is of lower
			risk than the predicate.
Algorithms	Osteotomy Module	Orthofix Implants Module (Not	The OFIX MIS App is a simplified software that
		machine-learning based)	is for templating Orthofix specific implants. The
			algorithm used is not machine learning based and
			uses information specific to Orthofix implant systems to analyze the captured image (input) and
			recommend the correct implant (output).
User Interface	PC	Mobile Device	This difference is expanded on more in the 510(k)
Oser Interface		Woode Device	Summary in which the OFIX MIS App is
			compared with the Surgimap 2.0 predicate and the
			Nuvaline reference devices. All the devices are
			used outside the sterile field and are accessible
			without impedance (no gloves). Compared the PC
			interface of Monitor/Keyboard/Mouse, The
			Mobile Device touch screen provides the same
			sufficient image handling (moving/zooming in)
			needed to perform the steps to measure and
			template the implant.
			Compared to the predicate, as discussed above,
			the overall risk is lower for the OFIX MIS App
			compared to the predicate in terms of its single
			functionality. The mobile interface makes it more
			accessible in the intraoperative environment. The
			safety and effectiveness of the App is
			demonstrated through user validation and testing, and this difference does not render the OFIX MIS
			App not substantially equivalent.
Obtaining an image	Transferred from other devices	Mobile Device Camera	This difference is expanded on more in the 510(k)
Ostaninis an illiage	Transferred from other devices	172020 Device Califera	Summary in which the OFIX MIS App is
			compared with the Surgimap 2.0 predicate and the
			Nuvaline reference devices. The Nuvaline
			(K162647) uses a Mobile Device Camera to
			obtain the image and predicated against the
			Surgimap 2.0. The use of a Mobile Device
			Camera to obtain the image has been validated to
			demonstrate sufficient input to the OFIX MIS
			App to measure and output an accurate implant
			length. User validation was performed to ensure
			safety and effectiveness. This difference does not
			render the App not substantially equivalent.
DICOM	Yes	N/A	The OFIX MIS App, is a lower risk tool used to
			streamline the process of choosing the correct rod
			to be used in the surgical procedure. It processes
			pictures of a fluoroscopy monitor acquired by a
			cell phone camera.



Non-Clinical Test Conducted for Determination of Substantial Equivalence

Software Verification and Validation Testing

Software verification and validation testing were conducted and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices." The testing demonstrated that the App meets the required specifications.

Usability Testing

The OFIX MIS App demonstrated accuracy and usability with bench top and simulated use testing.

Conclusion

The OFIX MIS App is substantially equivalent to the predicate despite the minor technological differences described. The testing data support the safety and effectiveness of the OFIX MIS App and demonstrate performance as intended in the specified use conditions. This performance is comparable to the predicate's respective software functionality.