



October 27, 2023

Abbott Point of Care Inc.
Brian Ma, Ph.D.
Principal Specialist Regulatory Affairs
400 College Road East
Princeton, New Jersey 08540

Re: K230275

Trade/Device Name: i-STAT CG8+ cartridge with the i-STAT 1 System
Regulation Number: 21 CFR 862.1665
Regulation Name: Sodium Test System
Regulatory Class: Class II
Product Code: JGS, CEM
Dated: September 25, 2023
Received: September 26, 2023

Dear Brian Ma:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,


Paula V. Caposino -S

Paula Caposino, Ph.D.
Acting Deputy Director
Division of Chemistry and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K230275

Device Name

i-STAT CG8+ cartridge with the i-STAT 1 System

Indications for Use (Describe)

The i-STAT CG8+ cartridge with the i-STAT 1 System is intended for use in the in vitro quantification of sodium and potassium in arterial or venous whole blood in point of care or clinical laboratory settings.

The i-STAT CG8+ cartridge with the i-STAT 1 System is intended for use in the in vitro quantification of sodium in capillary whole blood in point of care or clinical laboratory settings.

Sodium measurements are used for monitoring electrolyte imbalances.

Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

The information in this 510(k) summary is being submitted in accordance with the requirements of 21 CFR 807.92.

I. SUBMITTER INFORMATION

Owner Abbott Point of Care Inc.
400 College Road East
Princeton, NJ 08540

Contact Primary: Brian Ma, PhD
Principal Specialist Regulatory Affairs
Phone: + 1 613-688-5949

Secondary: Mojgan Soleimani
Associate Director Regulatory Affairs
Phone: + 1 613-295-0932

Date Prepared October 27, 2023

II. DEVICE INFORMATION

Proprietary Name *i-STAT CG8+* cartridge with the *i-STAT 1 System*

Common Name Chemistry test, analyzer, handheld

510(k) Number: K230275

Product Code	Device Classification Name	Regulation Number	Class	Panel
JGS	Electrode, Ion Specific, Sodium	862.1665	II	Clinical Chemistry
CEM	Electrode, Ion Specific, Potassium	862.1600	II	Clinical Chemistry

III. PREDICATE DEVICE

Proprietary Name *i-STAT CHEM8+* cartridge with the *i-STAT 1 System*

510(k) Number K183688

Product Code	Device Classification Name	Regulation Number	Class	Panel
JGS	Electrode, Ion Specific, Sodium	862.1665	II	Clinical Chemistry
CEM	Electrode, Ion Specific, Potassium	862.1600	II	Clinical Chemistry

IV. DEVICE DESCRIPTION

The *i-STAT CG8+* cartridge is used with the *i-STAT 1* analyzer as part of the *i-STAT 1 System* and contains test reagents to measure sodium (Na) in arterial, venous or capillary whole blood and to measure potassium (K) in arterial and venous whole blood.

The *i-STAT 1 System* is an *in vitro* diagnostic (IVD) medical device intended for the quantitative determination of various clinical chemistry tests contained within *i-STAT* cartridges using whole blood. The *i-STAT 1 System* consists of a portable blood analyzer (*i-STAT 1* analyzer), single-use disposable test cartridges (*i-STAT* cartridges), liquid quality control and calibration verification materials, and accessories (*i-STAT 1 Downloader/Recharger*, *i-STAT Electronic Simulator* and *i-STAT 1 Printer*). The *i-STAT 1 System*, including the *i-STAT CG8+* cartridge, is designed for use by trained medical professionals in point of care or clinical laboratory settings and is for prescription use only.

The *i-STAT CG8+* cartridge contains the required sensors, a fluid pack (calibrant pouch), a sample entry well and closure, fluid channels, waste chamber, and the necessary mechanical features for controlled fluid movement within cartridge. The *i-STAT* cartridge format allows all the tests in the cartridge to be performed simultaneously. All the test steps and fluid movement occur within the *i-STAT CG8+* cartridge. Cartridges require two to three drops of whole blood applied to the cartridge using a transfer device by the trained user before the cartridge is placed within the analyzer.

The *i-STAT 1* analyzer is a handheld, *in vitro* diagnostic analytical device designed to run only *i-STAT* test cartridges. The instrument interacts with the *i-STAT CG8+* cartridge to move fluid across the sensors and generate a quantitative result (within approximately 2 minutes).

V. INTENDED USE STATEMENT

The *i-STAT CG8+* cartridge with the *i-STAT 1 System* is intended for use in the *in vitro* quantification of sodium and potassium in arterial or venous whole blood in point of care or clinical laboratory settings.

The *i-STAT CG8+* cartridge with the *i-STAT 1 System* is intended for use in the *in vitro* quantification of sodium in capillary whole blood in point of care or clinical laboratory settings.

Sodium measurements are used for monitoring electrolyte imbalances.

Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.

VI. SUMMARY COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

Table 1: Similarities and Differences (Test and Instrument): Na and K in Whole Blood						
Feature or Characteristic	Candidate Devices: Na and K Tests in the: <i>i-STAT CG8+</i> cartridge with the <i>i-STAT 1 System</i>	Predicate Device: Na and K Tests in the: <i>i-STAT CHEM8+</i> cartridge with the <i>i-STAT 1 System</i> (K183688)				
Intended Use	<p><i>i-STAT CG8+</i> cartridge</p> <p>The <i>i-STAT CG8+</i> cartridge with the <i>i-STAT 1 System</i> is intended for use in the <i>in vitro</i> quantification of sodium and potassium in arterial or venous whole blood in point of care or clinical laboratory settings.</p> <p>The <i>i-STAT CG8+</i> cartridge with the <i>i-STAT 1 System</i> is intended for use in the <i>in vitro</i> quantification of sodium in capillary whole blood in point of care or clinical laboratory settings.</p> <p>Sodium measurements are used for monitoring electrolyte imbalances.</p> <p>Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.</p>	<p><i>i-STAT CHEM8+</i> cartridge</p> <p>The <i>i-STAT CHEM8+</i> cartridge with the <i>i-STAT 1 System</i> is intended for use in the <i>in vitro</i> quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings.</p> <p>Sodium measurements are used for monitoring electrolyte imbalances.</p> <p>Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.</p>				
Device Classification	Same	Class II				
Product Code	Same	JGS (Na) CEM (K)				
Regulation No.	Same	862.1665 (Na) 862.1600 (K)				
Reportable Range	Same	<table border="1"> <tbody> <tr> <td>Na</td> <td>100 – 180 mmol/L (mEq/L)</td> </tr> <tr> <td>K</td> <td>2.0 – 9.0 mmol/L (mEq/L)</td> </tr> </tbody> </table>	Na	100 – 180 mmol/L (mEq/L)	K	2.0 – 9.0 mmol/L (mEq/L)
Na	100 – 180 mmol/L (mEq/L)					
K	2.0 – 9.0 mmol/L (mEq/L)					

Table 1: Similarities and Differences (Test and Instrument): Na and K in Whole Blood											
Feature or Characteristic	Candidate Devices: Na and K Tests in the: <i>i-STAT CG8+</i> cartridge with the <i>i-STAT 1 System</i>	Predicate Device: Na and K Tests in the: <i>i-STAT CHEM8+</i> cartridge with the <i>i-STAT 1 System</i> (K183688)									
Sample Type	<table border="1"> <tr> <td>Na</td> <td>Arterial, venous or capillary whole blood</td> </tr> <tr> <td>K</td> <td>Arterial or venous whole blood</td> </tr> </table>	Na	Arterial, venous or capillary whole blood	K	Arterial or venous whole blood	Arterial and venous whole blood					
Na	Arterial, venous or capillary whole blood										
K	Arterial or venous whole blood										
Sample Volume	Same	95 µL									
Sample Preparation	Same	Ready to Use									
Sample collection	<table border="1"> <tr> <td></td> <td>Na</td> <td>K</td> </tr> <tr> <td>Without anticoagulant</td> <td colspan="2">Arterial or venous</td> </tr> <tr> <td>With balanced heparin anticoagulant or lithium heparin anticoagulant</td> <td>Arterial, venous, or capillary</td> <td>Arterial or venous</td> </tr> </table>		Na	K	Without anticoagulant	Arterial or venous		With balanced heparin anticoagulant or lithium heparin anticoagulant	Arterial, venous, or capillary	Arterial or venous	With balanced heparin anticoagulant or lithium heparin anticoagulant
	Na	K									
Without anticoagulant	Arterial or venous										
With balanced heparin anticoagulant or lithium heparin anticoagulant	Arterial, venous, or capillary	Arterial or venous									
Traceability	Same	<table border="1"> <tr> <td>Na, K</td> <td>NIST SRM956</td> </tr> </table>	Na, K	NIST SRM956							
Na, K	NIST SRM956										
Calibration	Same	1-point on-board contained within cartridge									
Principle of Measurement	Same	Na, K: Ion-selective electrode potentiometry									
Reagent Format	Same	Cartridge									
Reagent Storage and Stability	Refrigerated at 2-8°C (35-46°F) until expiration date Room Temperature at 18-30°C (64-86°F) for 2 months	Refrigerated at 2 to 8°C (35 to 46°F) until expiration date Room Temperature at 18-30°C (64-86°F) for 14 days									
Analyzer Type	Same	Handheld									

VII. PERFORMANCE CHARACTERISTICS

A. Analytical Performance

a. Precision/Reproducibility:

i. Precision 20 days (Aqueous materials)

The precision of the i-STAT Sodium (Na) and Potassium (K) tests in the *i-STAT CG8+* cartridge on the *i-STAT 1 System* was evaluated using five (5) levels of aqueous material. This 20-day precision testing was based on CLSI document EP05-A3: *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline – Third Edition*. The study was conducted using multiple analyzers and one (1) test cartridge lot over 20 days at one site. Repeatability, between-run, between-day, and within-laboratory precision were estimated for each level. The results of the 20-day precision study for the *i-STAT CG8+* cartridge on the *i-STAT 1 System* are shown in **Table 2**.

Test (units)	Fluid Level	N	Mean	Repeatability		Between-run		Between-day		Within-Laboratory	
				SD	%CV	SD	%CV	SD	%CV	SD	%CV
Na (mmol/L)	CV L1	80	99.3	0.17	0.17	0.09	0.09	0.05	0.05	0.19	0.20
	CV L2	80	121.4	0.20	0.16	0.06	0.05	0.05	0.04	0.22	0.18
	CV L3	80	134.8	0.22	0.16	0.06	0.05	0.06	0.05	0.23	0.17
	CV L4	80	161.3	0.27	0.16	0.12	0.07	0.07	0.04	0.30	0.19
	CV L5	80	181.1	0.32	0.18	0.12	0.06	0.09	0.05	0.35	0.19
K (mmol/L)	CV L1	80	2.09	0.008	0.37	0.003	0.16	0.001	0.04	0.009	0.41
	CV L2	80	2.87	0.007	0.25	0.002	0.09	0.002	0.08	0.008	0.28
	CV L3	80	3.76	0.012	0.31	0.004	0.10	0.006	0.15	0.014	0.36
	CV L4	80	6.41	0.021	0.33	0.008	0.12	0.012	0.18	0.025	0.39
	CV L5	80	7.99	0.027	0.33	0.010	0.13	0.014	0.18	0.032	0.40

ii. Multi-site and operator-to-operator precision (Aqueous materials)

Multi-day precision testing was performed at three (3) sites using a panel of aqueous solutions containing five (5) levels of sodium and potassium. At each site, each level was tested once a day by two (2) operators for five (5) days on six (6) *i-STAT 1* analyzers using *i-STAT CG8+* cartridges. Within-run, between-day, between-operator and within-site (total) variance components were calculated by site. These components were also calculated for all sites combined and provided in the **Table 3** below.

Test (units)	Fluid Level	N	Mean	Within-Run		Between-Day		Between-Operator		Within-Site (Total)		Between-Site		Overall	
				SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV	SD	%CV
Na (mmol/L)	CV L1	91	100.0	0.29	0.29	0.12	0.12	0.08	0.08	0.32	0.32	0.00	0.00	0.32	0.32
	CV L2	90	121.9	0.28	0.23	0.04	0.04	0.00	0.00	0.29	0.23	0.00	0.00	0.29	0.23
	CV L3	97	134.9	0.34	0.25	0.00	0.00	0.06	0.04	0.34	0.25	0.03	0.02	0.34	0.26
	CV L4	90	161.2	0.46	0.29	0.00	0.00	0.00	0.00	0.46	0.29	0.21	0.13	0.51	0.32
	CV L5	90	181.1	0.37	0.21	0.04	0.02	0.00	0.00	0.38	0.21	0.18	0.10	0.42	0.23
K (mmol/L)	CV L1	91	2.10	0.010	0.50	0.000	0.01	0.000	0.00	0.01	0.50	0.000	0.01	0.010	0.50
	CV L2	90	2.81	0.026	0.91	0.000	0.00	0.012	0.44	0.028	1.01	0.008	0.28	0.029	1.05
	CV L3	97	3.70	0.014	0.38	0.000	0.00	0.003	0.09	0.014	0.39	0.000	0.00	0.014	0.39
	CV L4	90	6.32	0.044	0.70	0.000	0.00	0.005	0.08	0.044	0.70	0.008	0.12	0.045	0.71
	CV L5	90	7.89	0.038	0.48	0.000	0.00	0.006	0.07	0.038	0.48	0.011	0.14	0.040	0.50

iii. Precision (Whole Blood)

Whole blood precision of the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge on the *i-STAT 1* System was evaluated using whole blood specimens ¹ collected with lithium heparin. The whole blood precision was assessed using the duplicate test results collected across multiple point of care sites. The results are summarized in **Table 4**.

Table 4: Whole Blood Precision of arterial, venous, and capillary whole blood for i-STAT CG8+ cartridge on the i-STAT 1 Analyzer							
Test (units)	Sample Type	Sample Range	N	Mean	SD	%CV	
Na (mmol/L)	Venous Whole Blood	100-130	17	122.6	0.30	0.24	
		>130-140	99	137.5	0.45	0.33	
		>140-180	67	146.2	0.43	0.30	
	Arterial Whole Blood	100-130	2	128.0	0.00	0.00	
		>130-140	89	137.4	0.42	0.31	
		>140-180	62	142.9	0.37	0.26	
	Capillary Whole Blood	100-130	3	120.8	0.41	0.34	
		>130-140	56	138.1	0.61	0.44	
		>140-180	95	142.4	0.62	0.44	
K (mmol/L)	Venous Whole Blood	2.0-3.5	27	3.22	0.036	1.12	
		>3.5-5.0	135	4.12	0.038	0.92	
		>5.0-9.0	19	6.49	0.032	0.50	
	Arterial Whole Blood	2.0-3.5	23	3.21	0.021	0.65	
		>3.5-5.0	124	4.11	0.032	0.79	
		>5.0-9.0	6	5.67	0.041	0.72	

b. **Linearity/assay reportable range:**

i. Linearity

The study was designed based on CLSI EPO6-Ed2: *Evaluation of the Linearity of Quantitative Measurement Procedures – Second Edition*.

The linearity of the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge with the *i-STAT 1* System was evaluated by preparing whole blood samples of varying analyte levels for each i-STAT test. The i-STAT Sodium and Potassium tests in the *i-STAT CG8+ cartridge* demonstrated linearity over the reportable range for each *i-STAT* test. Regression summary of the response for each i-STAT test versus the concentration of the whole blood samples of varying analyte levels is provided in **Table 5**.

Table 5: Regression Summary for the i-STAT Na and K tests in the i-STAT CG8+ Cartridge on the i-STAT 1 Analyzer						
Test	Units	Reportable Range	Range Tested	Slope	Intercept	R ²
Na	mmol/L	100 – 180	91.3 – 209.8	1.005	-0.525	0.9996
K	mmol/L	2.0 – 9.0	1.79 – 10.04	1.011	0.002	0.9994

¹ The capillary whole blood clinical precision study design involved the performance of two individual fingersticks, collected independently by two operators into two separate capillary tubes and tested on two (2) i-STAT CG8+ cartridges.

c. Detection Limit

i. Limit of Quantitation (LoQ)

The study was based on the CLSI EP17-A2: *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*.

The LoQ of the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge was evaluated on the *i-STAT 1* analyzer using two (2) *i-STAT CG8+* cartridge lots and whole blood that was altered to a low analyte level for each i-STAT test. The LoQ for the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge was determined to be at or below the lower limit of the reportable range for each of the i-STAT tests as shown in **Table 6**.

Test (units)	Lower limit of the reportable range	Determined LoQ
Na (mmol/L)	100	92
K (mmol/L)	2.0	1.6

d. Analytical Specificity

i. Interference

The study was based on CLSI EP07-ED3: *Interference Testing in Clinical Chemistry, Third Edition*.

The interference performance of the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge on the *i-STAT 1* analyzer with the *i-STAT 1 System* was evaluated using whole blood samples based on CLSI EP07-ED3: *Interference Testing in Clinical Chemistry, Third Edition*. The effect of each substance was evaluated by comparing the performance of a control sample, spiked with blank solvent solution, with the test results from a test sample spiked with the potentially interfering substance at the toxic/pathological concentration based on CLSI EP37-ED1: *Supplemental Tables for Interference Testing in Clinical Chemistry, First Edition*, as applicable. A substance was identified as an interferent if the difference between the control and test samples was outside of the allowable error ($\pm Ea$) for the i-STAT test. For an identified interferent, a dose-response was performed to determine the degree of interference as a function of the substance concentration.

Table 7 contains the list of potentially interfering substances tested and the interference results for the *i-STAT CG8+* cartridge.

Substance ²	Test Concentration		i-STAT Test	Interference (Yes/No)	Comments
	mmol/L (unless specified)	mg/dL (unless specified)			
Acetaminophen	1.03	15.6	Na	No	
			K	No	
Acetyl Cysteine	0.92	15	Na	No	

² The compound tested to evaluate the interfering substance is presented in parenthesis.

Table 7: Potentially Interfering Substances and Test Concentrations for the i-STAT Sodium and Potassium tests in the i-STAT CG8+ Cartridge					
Substance ²	Test Concentration		i-STAT Test	Interference (Yes/No)	Comments
	mmol/L (unless specified)	mg/dL (unless specified)			
(N-Acetyl-L-Cysteine)			K	No	
Acetylsalicylic Acid	0.167	3.0	Na	No	
Ammonium (Ammonium Chloride) ³	2.0	10.7	Na	No	
			K	No	
Ascorbic Acid (L-Ascorbic Acid)	0.298	5.25	Na	No	
			K	No	
Benzalkonium (Benzalkonium Chloride) ³	0.03	1.13	K	No	
β -Hydroxybutyric Acid ³	6.0	62.46	Na	No	
			K	No	
Bilirubin	0.684	40	Na	No	
			K	No	
Bromide ³ (Lithium Bromide)	2.5	21.7	Na	No	Use Another Method
			K	No	
	37.5	325.7	Na	No	
			K	No	
Calcium (Calcium Chloride)	5.0	20	Na	No	
			K	No	
Chloride (Lithium Chloride)	3.2	13.6	Na	No	
			K	No	
Cholesterol	11.0	425	Na	Yes	Decreased results > 400 mg/dL
			K	No	
Hemoglobin	10 g/L	1000	Na	No	
			K	No	
Heparin (Sodium Heparin)	3.30 U/mL	330 U/dL	Na	No	
Ibuprofen	1.06	21.9	Na	No	
Intralipid 20%	N/A	2395	Na	No	
		3216	K	No	
Lactate (Lithium Lactate)	10	90	Na	No	
			K	No	
Magnesium (Magnesium Chloride)	4.1	10	Na	No	
			K	No	
Nithiodote (Sodium Thiosulfate) ³	16.7	264.04	Na	Yes	Increased results \geq 2.1 mmol/L
			K	No	
Salicylate (Lithium Salicylate)	0.207	2.86	Na	No	
			K	No	
Triglyceride	16.94	1500	Na	No	
			K	No	
Uric Acid	1.4	23.5	Na	No	

B. Comparison Studies

³ The test concentration for this substance is not included in CLSI guideline EP37 1st edition.

a. Method Comparison with Comparator Device

Method comparison for the *i-STAT CG8+* cartridge with the *i-STAT 1 System* was demonstrated in studies based on CLSI EP09c-ED3: *Measurement Procedure Comparison and Bias Estimation Using Patient Samples – Third Edition*.

Lithium heparin venous and arterial whole blood specimens collected across multiple point of care sites were evaluated using *i-STAT CG8+* cartridges on the *i-STAT 1* analyzer against whole blood specimens tested on a comparative method. For sodium and potassium, the first replicate result from the *i-STAT 1* analyzer was compared to the mean result from the comparative method.

Two (2) capillary whole blood specimens collected from skin puncture with balanced heparin capillary tubes from each study subject across multiple point of care sites were evaluated and analyzed in singlicate on the *i-STAT 1* analyzer against the comparative method. A Passing-Bablok linear regression analysis for sodium was performed using the singlicate result from the *i-STAT 1* analyzer versus the singlicate result of the comparative method.

The venous and arterial data were pooled, and a Passing-Bablok linear regression analysis was performed using the *i-STAT* Potassium results from the *i-STAT CG8+* cartridges on the *i-STAT 1* analyzer versus the comparative method results. Method comparison results comparing the *i-STAT* Potassium performance on the *i-STAT 1* analyzer to comparative methods for arterial and venous are as shown in **Table 8**. In the table, N is the number of specimens in the data set, and r is the correlation coefficient.

Table 8: Method Comparison Results for the i-STAT K test in the i-STAT CG8+ Cartridge with i-STAT 1 System									
Test (units)	Comparative Method	N	Slope	Intercept	r	Xmin	Xmax	Medical Decision Level	Bias at Medical Decision Level
	Arterial/Venous								
K (mmol/L)	i-STAT CHEM8+	340	1.00	0.00	1.00	2.4	8.8	3.0	0.00
								5.8	0.00
								7.5	0.00

The venous, arterial, and capillary whole blood data were pooled, and a Passing-Bablok linear regression analysis was performed using the *i-STAT* Sodium results from the *i-STAT CG8+* cartridges on the *i-STAT 1* analyzer versus the comparative method results.

Method comparison results comparing the *i-STAT* Sodium performance on the *i-STAT 1* analyzer to the comparative method for arterial, venous, and capillary whole blood specimens are shown in **Table 9**. In the table, N is the number of specimens in the data set, and r is the correlation coefficient.

Table 9: Method Comparison Results for for the i-STAT Na test in the i-STAT CG8+ Cartridge with i-STAT 1 System								
Test (units)	Comparative Method		N	Slope	Intercept	r	Medical Decision Level	Bias at Medical Decision Level
	Arterial/Venous	Capillary						
Na (mmol/L)	i-STAT CHEM8+	epoc Blood Analysis System	551	1.00	0.00	0.99	115	0.0
							135	0.0
							150	0.0

The method comparison results for capillary whole blood specimens only for the i-STAT Sodium test are shown in **Table 10**.

Table 10: Results for i-STAT CG8+ Cartridge with i-STAT 1 System – Native and Contrived Capillary Specimens					
Test (units)	N	Slope	Intercept	r	Range
Na (mmol/L)	209	1.00	0.00	0.98	101 - 172

Bias at the medical decision levels for native capillary whole blood specimens only for the i-STAT Sodium test are shown in **Table 11**.

Table 11: Results for i-STAT CG8+ Cartridge with i-STAT 1 System – Native Capillary Specimens Bias at Medical Decision Levels				
Test (units)	N	Medical Decision Level	Bias	
			Estimate	95% CI
Na (mmol/L)	194	115	0.0	(-1.0, 0.0)
		135	0.0	(-1.0, 0.0)
		150	0.0	(-1.0, 0.0)

b. Matrix Equivalence

A matrix equivalence study was conducted to evaluate the performance of the i-STAT Sodium and Potassium tests in the *i-STAT CG8+* cartridge on the *i-STAT 1* System

using non-anticoagulated venous and arterial whole blood specimens. The study design and analysis method were based on recommendations from the Clinical and Laboratory Standards Institute (CLSI) guideline EP35: *Assessment of Equivalence or Suitability of Specimen Types for Medical Laboratory Measurement Procedures, 1st ed.* The matrix equivalence of each test in the *i-STAT CG8+* cartridge was assessed by comparing arterial or venous whole blood specimens collected without anticoagulant (candidate specimen type) to samples collected with balanced heparin or lithium heparin anticoagulant (primary specimen type). Each specimen was tested in duplicate using two (2) *i-STAT CG8+* cartridges with two (2) *i-STAT 1* analyzers. A Passing-Bablok linear regression analysis was performed using the first replicate result from the candidate (y-axis) versus the mean result from the primary specimen (x-axis). The regression analysis results are summarized in **Table 12**. In the table, N is the number of specimens in the data set, and r is the correlation coefficient.

Table 12: Matrix Equivalence Results						
Test (units)	N	Candidate Specimen Range	Primary Specimen Range	r	Slope	Intercept
Na (mmol/L)	295	102-178	102-178	0.99	1.00	0.00
K (mmol/L)	292	2.4-8.7	2.4-8.7	0.99	1.00	0.00

VIII. CONCLUSION

The results of these studies demonstrate that performance of the *i-STAT* Sodium and Potassium tests in the *i-STAT CG8+* cartridge with the *i-STAT 1 System* is substantially equivalent to the predicate device.