

August 16, 2023

Hadleigh Health Technologies % Allison Komiyama Consultant Rqm+ 2251 San Diego Ave Suite B-257 San Diego, California 92110

Re: K230298

Trade/Device Name: Celsi Monitor Regulation Number: 21 CFR 880.2910 Regulation Name: Clinical Electronic Thermometer Regulatory Class: Class II Product Code: FLL Dated: July 17, 2023 Received: July 18, 2023

Dear Allison Komiyama:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Porsche Bennet

Porsche Bennett For David Wolloscheck, Ph.D. Assistant Director DHT3C: Division of Drug Delivery and General Hospital Devices, and Human Factors OHT3: Office of GastroRenal, ObGyn, General Hospital and Urology Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230298

Device Name Celsi Monitor

Indications for Use (Describe)

Celsi Monitor is intended for use in hospitals under a clinician's supervision or at their direction to assist nurses in continuous temperature monitoring of neonates. The device measures abdominal skin temperature through direct contact.

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary K230298

DATE PREPARED

August 16, 2023

MANUFACTURER AND 510(k) OWNER

Hadleigh Health Technologies 200 Gate Five Road, Suite 215 Sausalito, CA 94965 Telephone: +1 (415) 454-3005 Official Contact: Molly McCabe, Managing Director, Projects

REPRESENTATIVE/CONSULTANT

Allison C. Komiyama, PhD, RAC Dulciana Chan, MSE RQM+ Telephone: +1 (412) 816-8253 Email: akomiyama@rqmplus.com dchan@rqmplus.com

DEVICE INFORMATION

Proprietary Name/Trade Name:	Celsi Monitor
Regulation Name:	Clinical Electronic Thermometer
Regulation Number:	21 CFR 880.2910
Regulation Class:	Class II
Product Code:	FLL
Premarket Review:	Drug Delivery and General Hospital Devices, and Human
	Factors (DHT3C)
Review Panel:	General Hospital

PREDICATE DEVICE IDENTIFICATION

The Celsi Monitor is substantially equivalent to the following predicate:

510(k) Number	Predicate Device Name / Manufacturer	Primary Predicate
K083185	Cadi SmartSense Wireless Temperature Monitoring	
	System / Cadi Scientific Pte. Ltd	¥



DEVICE DESCRIPTION

Celsi Monitor is a non-invasive continuous temperature monitor. Celsi Monitor achieves its intended purpose by measuring temperature at the surface of a neonate's abdomen by sensing changes in electrical resistance of a thermistor that is in thermal contact with the skin at the measured point. Celsi Monitor is provided with the Celsi Probe and Celsi Belt. The Celsi Monitor is powered by two rechargeable lithium-ion batteries. Both the Celsi Belt and Celsi Probe are reusable. The device displays measured patient body temperature for up to 12 hours and alarms when the temperature falls outside normothermic levels (<36.5 or >37.5 °C). The device is a skin contacting surface device with duration of contact < 24 hours.

INDICATIONS FOR USE

Celsi Monitor is intended for use in hospitals under a clinician's supervision or at their direction to assist nurses in continuous temperature monitoring of neonates. The device measures abdominal skin temperature through direct contact.

COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

Hadleigh Health believes Celsi Monitor is substantially equivalent to the predicate device based on the information summarized here:

The subject device has a similar design and intended use as the device cleared in K083185. Both devices may be used on neonates, are intended to continuously monitor patient temperature, alert users if temperatures are outside of the normothermic range, and uses thermistors as temperature sensors. Both devices have temperature sensors that are applied directly to the skin on the abdomen. The main differences between the subject and predicate device are the temperature probe application and device display. The subject device uses an adjustable belt to maintain contact between the patient and the temperature probe (Celsi Probe). However, the temperature probe of the predicate device is adhered to the patient with a single use application tape. The temperature probe of the subject device connects to a small monitor (Celsi Monitor) with a screen to display results and alerts to the user. The predicate device's temperature probe wirelessly connects to a receiver which can then send temperature data and alerts to be displayed on a PC or tablet. However, these technological characteristics do not raise different questions of safety and effectiveness and have undergone testing to ensure the subject device is substantially equivalent to the predicate device.

In summary, the subject and predicate devices are based on the following same technological elements:

• Use of thermistor with direct skin contact on the abdomen

The following differences exist between the subject and predicate devices:

- Use of an adjustable belt to maintain contact between the patient and the temperature probe
- Use of a small monitor to display measurements



	Subject Device	Predicate Device	Comparison
	Hadleigh	CADI Scientific Pte.Ltd.	N/A
	Health	Cadi SmartSense Wireless	
	Technologie	Temperature Monitoring	
	s Celsi	System	
	Monitor	K083185	
Indications for	Celsi Monitor is	The Cadi SmartSense	Similar: Both devices
Use	intended for use in	Wireless Temperature	are intended to
	hospitals under a	Monitoring System is	continuously monitor
	clinician's supervision	intended to measure	patient skin
	or at their direction to	abdominal surface	temperature, alert
	assist nurses in	temperature in adults	users to temperatures
	continuous	through neonates. The	outside of the
	temperature	ThermoSENSOR, a reusable	normothermic range,
	monitoring of	temperature sensor, is	and use thermistors as
	neonates. The device	applied to the patient by	temperature sensors.
	measures abdominal	means of single use	Temperature sensors
	skin temperature	application tapes. The	are applied directly to
	through direct contact.	ThermoSENSOR provides	the skin on the
		periodic wireless	abdomen of neonates.
		transmission of temperature	Clinical and non-clinical
		data which is utilized by the	testing have been
		SmartNODE wireless receiver	adequately conducted
		and the SmartSense PC	to demonstrate the
		application program to	differences do not raise
		record, store, and display the	new or different
		temperature information.	questions of safety and
		Warning: This equipment	effectiveness.
		measures and reports	
		abdominal surface	
		temperature - where direct	
		measurements of body core	
		temperature are required, it	
		is recommended to utilize	
		appropriate	
		core temperature	
		monitoring devices for this	
		purpose.	
Product Codes /	FLL / 21 CFR 880.2910	FLL / 21 CFR 880.2910	Same
Regulation			
Number			-
Regulation Name	Clinical Electronic Thermometer	Clinical Electronic	Same
		Thermometer	Como
Monitoring	Thermistor-based	Thermistor-based temperature	Same
Method	temperature	sensor	



	sensor		
Temperature Measurement	Resolution: 0.1 °C Range: 25 to 43.9 °C	Unknown	Similar: Temperature
Measurement	Kalige. 25 to 45.9 C		measurements were evaluated by performance and clinical
			testing to demonstrate differences do not raise
			new or different questions of safety and
			effectiveness.
Contact Duration	12 hours	Unknown	Similar: Safety of contact duration evaluated by biocompatibility and clinical testing to demonstrate differences do not raise new or
			different questions of
Patient	Neonates	Adults	safety and effectiveness. Similar: The subject
Population	Neonates	Pediatrics	device has a patient
ropulation		Neonates	population that is a subset
			of the predicate device.
			Both devices can be used
			on neonates.
Use	Hospital environment	Hospital environment	Same
Environment			
User	Healthcare professionals	Healthcare professionals	Same
Prescription Use	Yes	Yes	Same
Notifications/Al arms	Audio and visual alarms when temperature is <36.0 or >37.9 °C. Visual alarms when temperature is <36.5 or >37.5 °C	Unknown	Similar: Temperature measurement and alarms were evaluated by performance and clinical testing to demonstrate differences do not raise new or different questions of safety and effectiveness.
Power source	Rechargeable lithium-ion battery	Unknown	Similar: Safety and effectiveness of battery evaluated in accordance with IEC 62133.
Measurement site (skin/body)	Skin (abdomen)	Skin (abdomen)	Same
Mode (direct or adjusted)	Direct	Direct	Same
Mode (contact	Contact	Contact	Same



or non contact)			
Device Materials	Silicone, stainless steel, acrylonitrile butadiene styrene, and polyurethane	Unknown	Similar: Safety of materials evaluated by biocompatibility testing. Differences do not raise new or different questions of safety and effectiveness.
Shelf life/ In-use life	2 years	Unknown	Similar: Device use life evaluated by lifetime analysis testing. Differences do not raise new or different questions of safety and effectiveness.

SUMMARY OF NON-CLINICAL TESTING

Biocompatibility:

The subject device was evaluated for cytotoxicity, sensitization, and irritation in compliance to ISO 10993-5:2009, ISO 10993-10:2010, and ISO 10993-23:2010:2021.

Software Verification and Cybersecurity:

The software development and testing were executed with consideration to IEC 62304. Cybersecurity was evaluated per AAMI TIR57.

Electromagnetic Compatibility and Electrical Safety:

The subject device was tested in compliance with IEC 60601-1:2012, IEC 60601-1-2:2014, IEC 60601-1-6:2013, IEC 60601-1-8:2012, ISO 80601-2-56:2017, and IEC 62133:2012.

Performance testing:

Lifetime analysis testing after 100 simulated use cycles ISO 80601-2-56:2017

SUMMARY OF CLINICAL TESTING

A short duration (3 hour) and long duration (12 hour) clinical evaluation was performed using the subject device to evaluate its temperature measurement accuracy and effects on the subject's skin. The subject device was compared to the reference device, the Philips Intellivue MP30 (K061610). Results from both studies demonstrated comparable temperature measurement accuracy between the subject and predicate devices.



CONCLUSION

Based on the performance testing, biocompatibility, software verification and validation, electromagnetic compatibility and electrical safety, and clinical testing, it can be concluded that the differences between the subject and predicate device do not raise new or different questions of safety and effectiveness. Therefore, the subject device, Celsi Monitor is substantially equivalent to the predicate device, Cadi SmartSense Wireless Temperature Monitoring System cleared under K083185.