

August 9, 2023

Puyang Linshi Medical Supplies Co., Ltd. % Boyle Wang General Manager Shanghai Truthful Information Technology Co., Ltd. Room 608, No.738, Shangcheng Rd., Pudong Shanghai, Shanghai 200120 China

Re: K230304

Trade/Device Name: Polyisoprene Surgical Gloves Regulation Number: 21 CFR 878.4460 Regulation Name: Non-Powdered Surgeon's Glove Regulatory Class: Class I, reserved Product Code: KGO Dated: July 8, 2023 Received: July 11, 2023

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Allan Guan -S

For Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control and Plastic Surgery Devices
OHT4: Office of Surgical and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230304

Device Name Polyisoprene Surgical Gloves

Indications for Use (Describe)

The Polyisoprene Surgical Gloves are sterile and single use device intended to be worn on the hands of operating room personnel to protect a surgical wound from contamination.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary (K230304)

This summary of 510(k) safety and effectiveness information is being submitted in accordance with 21 CFR 807.92.

1.0 submitter's information

Name: Puyang Linshi Medical Supplies Co., Ltd.
Address: East of Panjin Road and North of Fumin Road in Puyang County, Puyang City, Henan Province 457001, China.
Phone Number: +86-19839327898
Contact: Catherine Liu
Date of Preparation: 08.01.2023

Designated Submission Correspondent

Mr. Boyle Wang Shanghai Truthful Information Technology Co., Ltd. Tel: +86-21-50313932 Email: Info@truthful.com.cn

2.0 Device information

Trade name:Polyisoprene Surgical GlovesCommon name:Surgeon's GlovesClassification name:Surgeon's GlovesModel(s):6.5, 7, 7.5, 8

3.0 Classification

Production code:KGORegulation number:21CFR878.4460Classification:Class IPanel:General Hospital

4.0 Predicate device information

Manufacturer:Better Care Plastic Technology Co., Ltd.Device:Sterile Polyisoprene Powder Free Surgical Gloves510(k) number:K171047

5.0 Indications for use

The Polyisoprene Surgical Gloves are sterile and single use device intended to be worn on the hands of operating room personnel to protect a surgical wound from contamination.

6.0 Device description

The proposed device is Polyisoprene Surgical Gloves, sterile and disposable devices. The proposed devices are made of polyisoprene. The proposed device is white. The design of proposed device is addressing the standards as ASTM D6124, ASTM D5151, and ASTM D3577.

7.0 <u>Summary comparing technological characteristics with predicate</u> <u>device</u>

ltem	Proposed device	Predicated device	Remark	
510(k) number	K230304	K171047		
Product Code	KGO	KGO	Same	
Regulation No.	21CFR878.4460	21CFR878.4460	Same	
Class	I	I	Same	
Indications	The Polyisoprene Surgical Gloves	This surgeon's glove is a sterile	Same	
for Use	are sterile and single use device	and single use device intended to		
	intended to	be worn on the hands of		
	be worn on the hands of	operating room personnel to		
	operating room personnel to	protect a surgical wound from		
	protect a surgical wound	contamination.		
	from contamination.			
Prescription or Over	Over-The-Counter-Use	Over-The-Counter-Use	Same	
The Counter Use				
Materials	Synthetic polyisoprene rubber	Synthetic polyisoprene rubber	Same	
Design	Single use	Single use	Same	
	Sterile	Sterile	Same	
	Powder-free	Powder-free	Same	
	Hand Specific	Hand Specific	Same	
	Beaded cuff	Beaded cuff	Same	
Color	White	Clear	Difference 1	
Dimensions and	Meets ASTM D3577- 2019	Meets ASTM D3577-09 (2015)	Difference 2	
physical properties				
Sterilization	EO Sterilization	Radiation	Difference 3	
method				
Sterility Assurance	10 ⁻⁶ SAL	10 ⁻⁶ SAL	Same	

Table1-General Comparison

Level (SAL)			
Freedom from	Meets ASTM D3577- 2019 Clause	Meets ASTM D3577- 09(2015)	Difference 2
holes	oles 8.3 (ASTM D5151-19) Ins		
	Inspection level/AQL: GI/AQL 1.5		
Powder-Free	Meets ASTM D 6124-06 The	Meets ASTM D 6124-06 The	Similar
	averaged residual powder content	averaged residual powder	
	for the glove during process	content for the glove	
	validation is 0.1mg per glove	during process validation is	
		0.16mg per glove	
Primary Skin	Under the conditions of the	Under the conditions of the	Difference 4
Irritation	study (per ISO 10993-23), the	study (per ISO 10993-10), the	
ISO 10993- 10:2010	device is not an irritant	device is not an irritant	
Dermal	Under the conditions of the	Under the conditions of the	Same
Sensitization -	study (per ISO 10993-10), not	study (per ISO 10993-10), not	
ISO 10993-10:2010	a sensitizer	a sensitizer	
Acute Systemic	Under the conditions of the study,	Under the conditions of the study,	Same
Toxicity - ISO	there was no mortality or evidence	there was no mortality or	
10993-11: 2006	of Acute systemic	evidence of Acute systemic	
	toxicity	toxicity	
Shelf Life	2 years	3 years	Difference 5

Analysis:

Difference 1: The colors of proposed device are different with those of the predicate, but they all meet the requirements of ASTM D3577-19, so the differences do not raise any new safety or performance questions.

Difference 2: The proposed device and the predicate device meet different requirements ASTM D3577-19 and ASTM D3577-09 as the ASTM D3577-19 is the only Consensus Standard for Rubber Surgical Gloves Recognized by FDA, but the differences do not raise any new safety or performance questions.

Difference 3: The proposed device and the predicate device sterilized by different methods EO Sterilization and Radiation, but the EO validation of proposed device was implemented based on ISO 11135:2014, including sterilizer installation, OQ, PPQ and MPQ of sterilization. And parameters in Sterilization Validation Report can prove these products being SAL of 10-6.

Difference 4: New standard ISO10993-23:2021 replaces the Skin Irritation test in ISO10993-10:2010. Difference 5: The shelf life of proposed device is verified by Product Performance Test Report after 2 Years Accelerated Aging.

8.0 Summary of Non-Clinical Performance Testing

The following performance data has been provided to demonstrate that the subject device meets the acceptance criteria in the standard.

No.	Name of the Test	Purpose	Acceptance Criteria	Results
INO.		Fuipose	Acceptance Chiena	Results
	Methodology / Standard	T 1		
1	ISO 10993-10:2010	This part of ISO	Skin Sensitization Test:	All grades are 0.
	Biological Evaluation Of	10993 assesses	provided	AU
	Medical Devices - Part	possible contact	grades less than 1,	All animals were survived,
	10: Tests For Skin	hazards from	otherwise sensitization.	and no abnormal signs were
	Sensitization.	chemicals released		observed during the study.
		from medical devices,		
2	ISO 10993-23:2021	which may produce	Skin Irritation Test:	The primary irritation index is
	Biological evaluation of	skin and mucosal	If the primary irritation	0.
	medical devices - Part	irritation, eye irritation	index is 0-0,4, the	
	23: Tests for irritation	or skin sensitization.	response category is	The response of the
			Negligible.	proposed device was
			0,5-1,9 means slight	categorized as negligible
			2-4,9 means moderate	under the test condition
			5-8 means severe	
3	ISO 10993-5:2009	This part of ISO	The viab.% of the 100%	Viab.% of 100% test article
	Biological Evaluation Of	10993 describes test	extract of the test article	extract is 6.1%
	Medical Devices - Part	methods to assess	is the final result, and if	
	5: Tests For In Vitro	the in vitro cytotoxicity	viability is reduced to	It means the proposed
	Cytotoxicity	of medical devices.	<70% of the blank, it	device have potential toxicity
			has cytotoxic potential.	to L-929 in the MTT method
4	ISO 10993-11: 2017	To evaluate the	Within the monitoring	There was no evidence of
	Biological evaluation of	potential for medical	period (72 h), if the	systemic toxicity from the
	medical devices — Part	device materials to	toxicosis response of	extract.
	11: Tests for systemic	cause adverse	testing group is not	
	toxicity	systemic reactions.	greater than that of	
	-		control group, the	
			testing sample is	
			regarded as acceptable.	
			J	
		1		

Table 2 Summary of Non-Clinical Performance Testing

5	ISO 10993-7 standards for EO/ECH residual testing	This part of ISO 10993 specifies allowable limits for residual ethylene oxide (EO)	Limit (< 24 h) EO 4 mg ECH 9 mg Prolonged (> 24 h < 30 d) EO 60 mg/30 d ECH 60 mg/30 d	EO residue: ≤ 4 mg in the first 24h; ≤ 60 mg in the first 30d; EO residue shall also \leq 10ug/g; ECH residue: ≤ 9 mg in the first 24h; ≤ 60 mg in the first 30d.
6	ASTM D6124-06 (Reapproved 2017), Standard Test Method for Residual Powder on Medical Gloves	This standard is designed to determine the amount of residual powder (or filter-retained mass) found on medical gloves	powder residue limit of 2.0 mg	0.1mg /glove
7	ASTM D5151-19, Standard Test Method for Detection of Holes in Medical Gloves. 12/23/2019	This test method covers the detection of holes in medical gloves.	Samples number: 200 gloves AQL: 1.5 (ISO 2859) Criterion ≤7 gloves for water leakage	0 glove water leakage found
8	ASTM D5712-15 Standard Test Method for Analysis of Aqueous Extractable Protein in Natural Rubber in Latex, Natural Rubber, and Elastomeric Products Using the Modified Lowry Method	This test method, for the determination of protein levels in latex, is primarily intended to test NR, latex, and elastomeric materials for residual protein content.	Have a recommended aqueous soluble protein content limit of 200 ug/dm ²	52ug/dm ²

9	ASTM	The specification is	Dimensions:	Dimensions:
	D3577-19, Standard	intended as a	-6.5: width 83 ± 6 mm	6.5: width: 83-84 mm
	Specification for Rubber	reference to the	Length \geq 265 mm	Length 280-289 mm
	Surgical Gloves	performance and	-7: width 89 ± 6 mm	Thickness:
	Surgical Gloves	safety of rubber		Finger 0.230-0.260
		surgical gloves. The	Length ≥265 mm	-
			-7.5: width 95 \pm 6mm	mm
		safe and proper use	Length ≥265 mm	Palm 0.188-0.207 mm
		of rubber surgical	-8: width 102 \pm 6mm	Cuff 0.137-0.151 mm
		gloves is beyond the	Length ≥265 mm	7: width 91-93 mm
		scope of this	Thickness:	Length 271-278 mm
		specification.	-Finger ≥0.10 mm	Thickness:
			-Palm ≥0.10 mm	Finger 0.211-0.241
			-Cuff ≥0.10 mm	mm
				Palm 0.181-0.193 mm
			Physical properties:	Cuff 0.133-0.144 mm
			 Before aging 	7.5: width 96-98 mm
			Tensile strength	Length 273-280 mm
			≥17MPa	Thickness:
			Ultimate	Finger 0.218-0.237mm
			● Elongation ≥	
			• 650%	Palm 0.180-0.192 mm Cuff 0.136-0.146 mm
			 Stress at 500% 	8: width 103-105 mm
			Elongation ≤ 7.0 MPa	Length 268-283 mm Thickness:
			•	Finger 0.221-0.272
			After Accelerated	mm
			Aging	Palm 0.189-0.212 mm
			Tensile strength	Cuff 0.143-0.153 mm
			● ≥12MPa	
			● Ultimate Elongation ≥ 490%	Physical properties: Before aging
				-Tensile strength 25.2-30.9 MPa
				-Ultimate Elongation 721% - 777%
				-Stress at 500% Elongation 6.5-7.0 MPa
				After Accelerated Aging
				-Tensile strength 23.0-29.0 MPa
				-Ultimate Elongation 680% - 716%

9. Summary of Clinical Performance Test

No clinical study is included in this submission.

10.0 Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as or better than the legally marketed predicate device.