

August 16, 2023

Shanghai Omni Laser Skinology Co., Ltd.
% Helen Nan
General Manager
New Risen Enterprise Management Consulting Co., Ltd.
Room 302, Building 3, Hangqian Mansion, Hangqian Street
Lucheng District
Wenzhou, Zhejiang 325000
China

Re: K230342

Trade/Device Name: Phototherapy System (OL-PDT950) Regulation Number: 21 CFR 878.4810 Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology Regulatory Class: Class II Product Code: GEX Dated: July 17, 2023 Received: July 17, 2023

Dear Helen Nan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jianting Wang -S

For Tanisha L. Hithe, MS, MHS Acting Assistant Director DHT4A: Division of General Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number *(if known)* K230342

Device Name Phototherapy System (OL-PDT950)

Indications for Use (Describe)

he Phototherapy System use of the red, blue and infrared regions of the spectrum is intended to emit energy to treat dermatological conditions.

The blue light (469nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris.

The red light (633nm wavelength) is generally indicated to treatment of superficial, benign vascular, and pigmented lesions.

The infrared light (835nm wavelength) is generally use for the temporary relief of minor muscle and joint pain, arthritis and muscle spasm; relieving stiffness; promoting the relaxation of muscle tissue; and to temporarily increase local blood circulation where applied.

Type of Use (Select one or both, as applicable	e)	
--	----	--

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov*

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



K230342_510(k) Summary (As required by 21 CFR 807.92)

1.0 Submitter Information

Company:	Shanghai Omni Laser Skinology Co., Ltd.
Address	Floor 3, Building 3, NO. 227, Mingqiang Road, Songjiang
	District, Shanghai, 201612, CHINA
Phone:	+86-021-54847192
E-mail:	avril@omni-laser.com
Contact	Avril Ouyang
Title:	General Manager
Prepared date: 07/17/2023	-

2.0 Device Information

Trade/Device Name:	Phototherapy System
Model:	OL-PDT950
Device:	Powered Laser Surgical Instrument
Review Panel:	General & Plastic Surgery
Product Code:	GEX
Submission Type:	Traditional 510(k)
Regulation Number:	CFR 878.4810
Device Class:	Class II

3.0 Predicate Device Information

K200751		
Photodynamic Therapy System Shangdong Huamei Technology Co., Ltd.		
K222751		
LED Light Therapy Device		
Xuzhou Kernel Medical Equipment Co., Ltd.		
K200104		
Oxylight		
RAJA Trading Company, Inc.		

4.0 Indications for Use

The Phototherapy System use of the red, blue and infrared regions of the spectrum is intended to emit energy to treat dermatological conditions.

The blue light (469nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris.

The red light (633nm wavelength) is generally indicated to treatment of superficial,



benign vascular, and pigmented lesions.

The infrared light (835nm wavelength) is generally use for the temporary relief of minor muscle and joint pain, arthritis and muscle spasm; relieving stiffness; promoting the relaxation of muscle tissue; and to temporarily increase local blood circulation where applied.

5.0 Device Description

The Phototherapy System OL-PDT950 is a vertical device which uses specific wavelengths of light, produced by LEDs (Light emitting diodes), to manage aesthetic conditions.

The device produces light in the red light region of the spectrum (633 ± 15 nm), in the blue light regions of the light spectrum (469 ± 15 nm) and infrared light region of light spectrum (835 ± 15 nm).

Five sets of LEDs panels are available for the device.

6.0 Discussion of Tests Performed

6.1 Clinical Tests

Clinical testing was not performed for the subject device as part of the submission.

6.2 Non-Clinical Tests

The following performance data are provided in support of the substantial equivalence determination:

IEC 60601-1: 2005 + CORR.1:2006 + CORR.2:2007 + A1:2012: Test for Medical Electrical equipment was performed for General Requirements for basic safety and essential performance;

IEC 60601-1-2: 2014: Test for Medical Equipment for General Requirements for basic safety and essential performance: electromagnetic compatibility

IEC60601-2-57:2021: Medical electrical equipment -- Part 2-57: Particular requirements for the basic safety and essential performance of non-laser light source equipment intended for therapeutic, diagnostic, monitoring and cosmetic/aesthetic use;

IEC 62471: 2006: Photobiological safety of lamps and lamp systems.

7.0 Software

Software verification and validation testing were conducted and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices".



8.0 Comparison of Technological Characteristics with the Predicate Device

Device Feature	Subject Device	Predicate Device 1	Predicate Device 2	Predicate Device 3
K Number	K230342	K200751	K222751	K200104
Name	Phototherapy System	Photodynamic Therapy System	LED Light Therapy Device	OxyLight
Product Code	GEX	GEX	GEX	GEX
Indications for use	The Phototherapy System use of the red, blue and infrared regions of the spectrum is intended to emit energy to treat dermatological conditions. The blue light (469nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris. The red light (633nm wavelength) is generally indicated to treatment of superficial, benign vascular, and pigmented lesions. The infrared light (835nm wavelength) is generally use for the temporary relief of minor muscle and joint pain, arthritis and muscle spasm; relieving stiffness; promoting the relaxation of	The Photodynamic Therapy (PDT) Equipment use of the red, blue and infrared regions of the spectrum is intended to emit energy to treat dermatological conditions. The blue light (415nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris. The red light (630nm wavelength) is generally indicated to treatment of superficial, benign vascular, and pigmented lesions. The infrared light (835nm wavelength) is generally use for the temporary relief of minor muscle and joint pain, arthritis and muscle spasm; relieving stiffness;	LED Light Therapy Device use of the red, blue, Yellow and infrared regions of the spectrum is intended to emit energy to treat dermatological conditions. The red light (633±10nm wavelength) is generally indicated to treatment of superficial, benign vascular, and pigmented lesions. The blue light (417±10 nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris. The Yellow light (599±10nm wavelength) is generally indicated to treat dermatological conditions and specifically indicated for treatment of periorbital wrinkles and rhytides.	The Oxylight is intended for dermatological use by physicians and healthcare professionals for the following: LED Technology is intended for: Blue LED – 465nm – to treat dermatological conditions and specifically indicated to treat moderate inflammatory acne vulgaris. Red LED 625nm- for treatment of superficial, benign vascular and pigmented lesions. Yellow LED 590nm - treatment of periorbital wrinkles and rhytides.

Table 1 - General Comparison



Shanghai Omni Laser Skinology Co., Ltd.

muscle tissue; and	o promoting the relaxation of	The infrared light	
temporarily increase loc	al muscle tissue; and to	(835±15nm wavelength) is	
blood circulation whe	e temporarily increase local	generally use for the	
applied.	blood circulation where	temporary relief of minor	
	applied.	muscle and joint pain,	
		arthritis and muscle spasm;	
		relieving stiffness;	
		promoting the relaxation of	
		muscle tissue; and to	
		temporarily increase local	
		blood circulation where	
		applied.	

Table 2 - Performance Comparison

Device Feature	Subject Device (K230342)	Predicate Device 1 (K200751)	Predicate Device 2 (K222751)	Predicate Device 3 (K200104)
W/141-	Red light 633nm±15nm	Red light 630±15	Red light 633nm±15nm	Red light 625nm±5nm
(nm)	Blue light 469nm±15nm	Blue light 415±15	Blue light 417nm±15nm	Blue light 465nm±5nm
(1111)	Infrared light 835nm±15nm	Infrared light 835±15	Infrared light 835nm±15nm	Yellow light 590nm±5nm
	5 panels: 2400 EA LEDs.	5panel: 300 EA LEDs.	RBY Irradiator has 5	Three type, each head type
	The panels may emit the	The panels may emit the	panels:	has only one light.
	three light (red, blue,	three light (red, blue,	Red light:465EA LEDs;	Red, Blue, yellow.
	infrared)	infrared)	blue light: 470EA LEDs;	
	Blue: 800EA,		yellow light: 465EA LEDs;	
	Red: 800EA,		RBI irradiator: has 5	
Panels Type	Infrared: 800EA		panels:	
			Red light:465EA LEDs;	
			blue light: 470EA LEDs;	
			infrared: 465EA LEDs; The	
			panels may emit the three	
			light (red, blue infrared)	
			individual or in	



			combination	
	Each LED lamp bead has 1	Each LED lamp bead has 3	Each panel has three	Unknown
	diodes that emit single	diodes that emit different	different kinds of light-	
	color, the energy power of a	colors, the Energy power of	emitting diodes, and the	
LED power	diode is:	a diode is 3W.	energy power of the diode	
	Red light: 120mw		is 0.5W	
	Blue light: 160mw			
	Infrared light: 100mw			
	Red light: 75mW/cm ²	Red light: 115 mW/cm ²	Red light: $20 \sim 96 \text{ mW/cm}^2$	Red light: 100mW/cm ²
	Blue light: 35 mW/cm ²	Blue light: 120 mW/cm ²	Blue light: 10~120	Blue light: 45 mW/cm ²
Maximum power	IR: 120 mW/cm^2	IR: 70 mW/cm ²	mW/cm ²	Yellow light: 35 mW/cm ²
density			Infrared: $\leq 70 \text{ mw/cm}^2$	_
			Red/IR: 166mW/cm ²	
			Blue/IR: 190mW/cm ²	
Treatment area	1228 cm^2	1410 cm^2	$900 \text{ cm}^2 \pm 10\% \text{ cm}^2$	500cm ² and 860cm ²
Traatmant time	20 minutes (recommended	20 minutes (recommended	20 minutes (recommended	20 minutes (recommended
	Treatment Time)	Treatment Time)	Treatment Time)	Treatment Time)
Working	10~15cm	10~15cm	Unknown	Unknown
distance				
Power supply	AC 100-240V 50/60Hz	AC 100-240V 50/60Hz	AC 100-240V 50/60Hz	Unknown
	440VA	440VA		
Operation	Display Screen	Display Screen	Display Screen	Display Screen
interface				
Electrical	Comply with IEC 60601-1	Comply with IEC 60601-1	Comply with IEC 60601-1	Comply with IEC 60601-1
Safety	and IEC 60601-1-2	and IEC 60601-1-2	and IEC 60601-1-2	and IEC 60601-1-2
Radiation Safety	Comply with IEC 60601-2-	Comply with IEC 60601-2-	Comply with IEC 60601-2-	Comply with IEC 60601-2-
	57	57	57	57
Photobiological	Comply with IEC 62471	Comply with IEC 62471	Comply with IEC 62471	Comply with IEC 62471
safety				



9.0 Conclusion

Based on the comparison and analysis above, the proposed devices are determined to be Substantially Equivalent (SE) to the predicate devices.